Operating manual

VWR® pH Sensors

	EU cat. no	NA cat. no		EU cat. no	NA cat. no
pH SJ 113	662-1382	76460-468	pH SF 224	662-1399	-
pH SJ 114	662-1383	-	pH SF 225	662-1401	-
pH SJ 115	662-1384	-	pH SM 223	662-1402	76460-480
pH DJ 113	662-1385	76460-470	pH SM 224	662-1403	-
pH DJ 114	662-1386	-	pH SM 225	662-1404	-
pH DJ 115	662-1387	-	pH SP 223	662-1405	76460-482
pH SF 113	662-1388	76460-472	pH SP 224	662-1406	-
pH SF 114	662-1389	-	pH SP 225	662-1407	-
pH SF 115	662-1390	-	Redox RD 113	662-1408	76460-484
pH SM 123	662-1391	76460-474	Redox RD 223	662-1409	76460-486
pH SM 124	662-1392	-			
pH SM 125	662-1393	-			
pH SJ 223	662-1395	76460-476			
pH SJ 224	662-1396	-			
pH SJ 225	662-1397	-			
pH SF 223	662-1398	76460-478			

Legal Address of Manufacturer:

Europe

VWR International BV
Researchpark Haasrode 2020
Geldenaaksebaan 464
B-3001 Leuven
+32 16 385011
<http://be.vwr.com>

United States

VWR International LLC
100 Matsonford Rd
Radnor, PA 19087
+1 800-932-5000
www.vwr.com

Importer to UK:

VWR International Ltd
Hunter Boulevard, Magna Park
Lutterworth, Leicestershire, LE17 4XN
<http://uk.vwr.com>

Country of origin: Germany.

Contents

EN English	4
DE Deutsch	5
FR Français	6
ES Español	7
IT Italiano	8
PT Português	9

Supply

All electrodes are supplied tested and ready for use, and meet the following characteristic data at time of delivery.

Electrode slope at 25 deg C	> 95 % of the theoretical value
Offset voltage at the zero point pH 7 25 deg C	< abs. 20 mV

Storage

Unless advised to the contrary, electrodes should be stored in the protective cap or the storage container with their membranes immersed in VWR KCl storage solution 83605.180. Storage temperature should be between 10 deg C and 30 deg C, but not below -5 deg C. Should any KCl solution leak from the protective cap during storage or transit usually seen as a white salt, then it may be easily removed by rinsing the electrode in water. The electrode will not be affected by this leakage.

Calibration

Carefully rinse electrodes with deionized/distilled water and remove air bubbles by shaking as with a mercury clinical thermometer before immersing up to the diaphragm into the buffer solution.

Always rinse the electrode in deionized/distilled water before using other buffer solutions.

Discard the buffer solution after use. Please note that buffer solutions do have a finite shelf life.

If possible calibration buffer solutions should be at the normal working temperature.

Calibrate electrodes in accordance with the instrument's working manual.

Where possible, electrodes should be installed vertically or within 30 degrees to the vertical.

Ageing

All electrodes under go a natural ageing process caused by the exchange reactions of the glass membrane. This usually results in the electrode having a slow response time, decrease of slope, or changes of zero point. High temperature or aggressive samples will reduce the lifetime of the electrodes. Please note that even in storage, electrodes will still be subjected to the ageing process. With refillable electrodes, always check that there is sufficient 3 mol/l KCl electrolyte. The level of the reference electrolyte should be several cm higher than the level of the measuring medium.

Gel-filled combination electrodes are nearly maintenance free and need no refilling of electrolyte.

Cleaning

pH and Redox combination electrodes can be cleaned by rinsing in deionized/distilled water, and if necessary dried with tissue paper. Check that the tissue paper has not left any fibres on the membrane. The Redox electrode's platinum disc or pin may be cleaned by wiping with tissue paper.

Protein contamination may be removed by soaking for 1 hour in cleaning solution Pepsin/HCl.

After use in oily or organic solution, the pH electrode should be cleaned by dipping into acetone or ethanol solution for a short period.

After cleaning the electrode should be rinsed thoroughly in deionized/distilled water before soaking in 3 mol/l KCl. The electrode should always be recalibrated before further use.

Ensure that all electrical terminations and connections are kept clean and dry.

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VWR International AB
VWR International GmbH
VWR International, LLC