

Instruction manual

VWR® NANO SPECTROPHOTOMETER

EU cat. no 76777-152

NA cat. no 76777-150



CE

UK
CA

Legal address of Manufacturer

Europe

VWR International bv,
Researchpark Haasrode 2020
Geldenaaksebaan 464 B-3001
Leuven
+3216385011
be.vwr.com

UK Importer:
VWR International Ltd
Hunter Boulevard, Magna Park
Lutterworth, Leicestershire, LE17 4XN
uk.vwr.com

United States

VWR International LLC
100 Matsonford Rd
Radnor, PA 19087
+1 800-932-5000
www.vwr.com

Country of Origin

Singapore

Table of Contents

Warning	4
Safety information	4
Package Contents	5
Installation	5
Intended Use	5
Product Specifications	6
Description of Buttons and Switches	7
Operation	8
<i>Start-up Interface</i>	8
<i>Main Menu Interface</i>	8
<i>Nucleic Acids Interface</i>	8
<i>Protein A280 Interface</i>	12
<i>Colorimetry Interface</i>	15
<i>UV-Vis Interface</i>	18
<i>OD600 Interface</i>	20
<i>Performance Verification</i>	21
<i>System Settings Interface</i>	23
Troubleshooting	25
Repair and Maintenance	25
User Replaceable Accessories and Spare Parts	25
Technical Service	26
Warranty	26
Equipment Disposal	27

Warning



Any improper operation may cause injury. Please read this manual carefully and operate safely according to the guidelines.



The instrument is a normal indoor instrument which conforms to class I of the GB 4793.1 standard.



This instrument is designed for use in a laboratory environment. The device must be operated by skilled laboratory personnel with appropriate training.



To prevent injury or voiding the warranty, the operator should not attempt to repair the instrument without explicit guidance from VWR. If service is required, please contact VWR service.

Safety Information

Avoiding Electrical Shock

Please follow the guidelines below and read this manual in its entirety to ensure safe operation of the unit.



Before powering on, confirm that the voltage used meets the electrical requirements of the instrument as stated on the rating plate. If the electric cord is damaged, replace it with the same type of cord. Hold the socket firmly before pulling the plug from an outlet. Do not pull the electric cord.



The instrument should be installed in an environment of standard room temperature, low dust, low humidity, and away from direct sunlight, electromagnetic interference, and heat sources. Do not block the vents on the instrument.



Always power off the instrument when you are finished using it. Unplug the power cord and cover the instrument with a cloth or plastic sheet to prevent excessive dust from entering the housing.



Pull the connector plug from the electrical outlet immediately and contact the vendor in the event of:

- Liquid entering the housing.
- Abnormal operation: such as any abnormal sound or smell.
- The instrument is dropped or there is any damage to the housing.
- Any malfunction.

Package Contents

Description	Quantity
VWR NANO SPECTROPHOTOMETER	1
Power Adapter (DC 24V 2A)	1
Power Line	1
Dust Cover	1
Logitech Mouse	1
Printer Paper	2
Glass Cuvette	2
Operation Manual	1
Performance Test Report	1
USB Drive	1

Installation

Required Installation Environment

Environmental Temperature: 5°C ~ 35°C

Relative Humidity: ≤ 70%

Input Voltage: DC 24V, 2A (Adapter CSA, UL, CE marked)

Intended use

For research use only. Not intended for any animal or human therapeutic or diagnostic use.

Product Specifications

Model	VWR Nano Spectrophotometer
Sample Size	0.5µL – 2.0µL
Path Length	0.2mm, 1mm
Light Source / Life	Xenon flash lamp / >10 ⁹ flashes
Detector Type	2048 element linear silicon CCD array
Wavelength Range	200nm — 800nm
Wavelength Accuracy	± 1 nm
Spectral Resolution	≤ 3nm (FWHM@Hg 253.7nm)
Absorbance Precision	0.003Abs (1mm path length)
Absorbance Accuracy	± 1% (7.332Abs, at 260nm wavelength)
Absorbance Range	0.04 — 90 (at 260 wavelength, 10mm equivalent)
Detection Concentration Range	2ng/µL dsDNA ~ 4,500ng/µL dsDNA
Detection Time	< 5 seconds
OD600	Abs range 0~4.000 Abs Abs stability [0,3) ≤0.5% [3,4) ≤2% Abs repeatability [0,3) ≤0.5%, [3,4) ≤2% Abs Precision [0,2) ≤0.005A, [2,3) ≤1%, [3,4) ≤2%
Voltage input	DC 24V, 2A
Power	25W
Dimensions (W×D×H)	20.8cm × 28.0cm ×18.6cm / 8.2in x 11.0in x 7.3in
Weight	3.6 kg / 7.9 lbs

Description of Buttons and Switches

Front View




Fig. 1 Front View

Rear View



Fig. 2 Rear View

Operation

Start-up Interface

Upon powering on the instrument, it will perform a self-check, and the start-up screen will be displayed (**Fig. 3**)



Fig. 3 Start-up Interface

Main Menu Interface




Fig. 4 Main Menu Interface

After start-up, the main menu interface will be displayed. There are 6 options: Nucleic Acids, Protein A280, Colorimetry, UV-Vis, OD600, & System Settings.

Nucleic Acids Interface

Beer-Lambert's Law for DNA/RNA quantitation

The following "Beer - Lambert" equation is used to calculate the concentration of nucleic acids:

$$C = \frac{A * \epsilon}{b}$$

C = Sample DNA concentration, unit : ng/μL

A = Sample absorbance, unit : A

ε = extinction coefficient, unit: ng-cm/μL

b = Path Length, unit: cm

Standard DNA/RNA extinction coefficients :

dsDNA : 50ng-cm/μL

ssDNA : 33ng-cm/μL

RNA : 40ng-cm/μL

When the sample column is used, highly concentrated nucleic acid samples can be measured without dilution using a 1.0mm or 0.2mm path length. The VWR Nano Spectrophotometer will measure and display sample absorbance values of the 10mm pathlength equivalent of up to 90 A. The instrument will accurately measure dsDNA samples up to 4500 ng/μL without dilution. To do this, the instrument automatically detects highly concentrated samples and adjusts the pathlength to 0.2mm to measure the sample absorbance.

Select “Nucleic acid” from the main menu to enter the Nucleic Acids Interface. By default, the system will provide guidance for measuring blanks and samples, as shown below (**Fig. 5**). Touch the screen to dismiss the message displayed. Protocol Guidance can be toggled on or off in the System Settings Interface (see *System Settings Interface*).




Fig. 5 Nucleic Acid Interface – Protocol Guidance




Fig. 6 Nucleic Acid Interface

Fig 5; there are three tabs in the Nucleic Acids Interface: Nucleic Acids, Files, and Help.

- **ID :** **191028_170755** : The sample ID name has a default value of the current date and time. Users can rename the sample ID. One sample ID can contain up to 1000 measurement values.
- **DNA-50** : Select the sample type: DNA-50 for dsDNA, RNA-40 for RNA, ssDNA-33 for ssDNA. For a different nucleic acid type, select “others” and enter the extinction coefficient.
- **Blank** : Perform a blank reading. This step is essential before measurement. Blank absorbance values are typically in the range of 0.004-0.03 Abs and are valid for up to 30 minutes. The instrument will automatically remind the user to perform another blank reading after 30 minutes.
- **Baseline at 340nm** : Spectrum normalization; The baseline is automatically set to the absorbance of the sample at 340nm and can be modified. This feature can remove spectroscopic signals from sample measurements by subtracting the measured absorbance at a specified baseline correction wavelength from the absorbance values at all wavelengths of a measured sample.
 - **Note:** If baseline calibration is not performed, the spectroscopic signals will not be separated from interference/background effects and will lead to inaccurate results.
- **Auto Blank** : When turned on, the blank measurement starts automatically when the upper pedestal is lowered.
 - **Note:** The automatic blank only takes effect when no blank readings have been made in a current run.

- **Auto Sample**  : When turned on, the sample measurement starts automatically when the upper pedestal is lowered.
 - **Note:** Automatic sample detection only takes effect after a blank reading has been performed.
- The  icon appears in the upper right corner to indicate an error in reading blank/sample volumes. Please clean and wipe the pedestal and perform another blank reading. If the problem persists, contact VWR service.

Operation:

1. Set the Sample ID.
2. Clean the upper and lower pedestals with a lint-free wipe, add 2µL buffer solution to perform a blank reading.
3. Clean the buffer solution on the pedestals with a wipe.
4. Measure a 2µL sample volume and click “Measure” to detect the sample.

Note: A blank reading must be performed prior to sample measurements. The sample volume must be equivalent to the volume used to set the blank reading.

5. Clean and dry the pedestal between measurements.

The sample concentration and absorbance ratios will display on the right side of the interface (**Fig. 7**).



Fig. 7 Sample Concentration & Absorbance Ratios

Conc. : Calculated nucleic acid concentration.

A260: The sample absorbance at 260nm (10mm pathlength equivalent).

A280: The sample absorbance at 280nm (10mm pathlength equivalent).

A230: The sample absorbance at 230nm (10mm pathlength equivalent).

A260/A280: The ratio of corrected absorbance values at 260nm to corrected absorbance values at 280nm. An A260/A280 ratio of ~1.8 is generally accepted as pure for DNA (~2.0 for RNA). Acidic solutions will under-represent the A260/A280 ratio by 0.2 – 0.3 units, whereas basic solutions will over-represent the A260/280 ratio by 0.2 – 0.3 units.

A260/A230 : The ratio of corrected absorbance values at 260nm to the corrected absorbance values at 230nm. This ratio can be used as a secondary measure of nucleic acid purity. An A260/A230 ratio within the range of ~1.8 – 2.2 is considered as pure for nucleic acids. Lower ratio values indicate the presence of contaminants that absorb strongly at or near 230nm.

wavelength: **260nm** Absorb: 5.909 : Input a wavelength to view the corresponding sample absorbance value (**Fig. 8**).




Fig. 8 Nucleic Acid Sample Absorbance Curve (220nm – 350nm)

Blank : Perform a blank reading

Sample : Perform a sample reading.

Print : Print the calculated data.

Save Spectral Data : Save the full spectrum absorbance data in 1nm steps.

Enlarge Curve : Zoom into the sample absorbance curve.

Back : Return to the main menu interface.

Nucleic Acid Files Interface

Sample ID	Load #	A260	A280	A230	A260/A280	A260/A230	C(ng/ul)	Time
230725_010914	1	0.005	0.001	0.001	3.09	3.09	0.297	07-25-2023 01:13:25

Fig. 9 Nucleic Acid Files Interface

Select the “Files” tab at the top of the Nucleic Acid Interface (**Fig. 9**). Users can select previously saved results by the file name.

Spectrum Data: View spectral data; To view sample absorbances at a specific wavelength, input a wavelength and the red coordinate line will display the corresponding absorbance value (**Fig. 10**).




Fig. 10 Nucleic Acid Sample Absorbance Curve (220 – 350nm).

Export Load No. : Export the result to a USB flash drive

Delete Load No. : Delete the selected results.

Delete Sample ID: Delete the selected files

Protein A280 Interface

Introduction

The Protein A280 interface can be used to quantify purified proteins that contain amino acids such as tryptophan, tyrosine, or cys-cys disulfide bonds. These amino acids exhibit peak absorbance at 280nm. The following sample types can be selected: "A280", "BSA", "IgG", "Lysozyme", & "Others".

This interface does not require the generation of a standard curve. Sample absorbance values (260nm & 280nm) and A260/A280 ratios can be measured and displayed. A baseline correction can be used for normalization. Like the Nucleic Acids interface, the Protein A280 interface automatically switches the pathlength to 0.2mm when a highly concentrated protein sample is detected and displays 10mm pathlength equivalent data.

Protein A280 Interface

Select "Protein A280" from the main menu interface.




Fig. 11 Protein A280 Interface

Fig. 11; There are three options at the top of the screen, Protein A280, Files, and Help.

- **ID** : The sample ID name has a default value of the current date and time. Users can rename the sample ID. One sample ID can contain up to 1000 measurement values.
- **A280** : Select the sample type: A280, BSA, IgG, Lysozyme. Select "others" and type in the extinction coefficient.
- **Blank**: Perform a blank reading. This step is essential before measurement. Blank absorbance values are typically in the range of 0.004-0.03 Abs and are valid for up to 30 minutes. The instrument will automatically remind the user to perform another blank reading.
- **Baseline at 340** : Spectrum normalization; The baseline is automatically set to the absorbance of the sample at 340nm and can be modified. This feature can remove spectroscopic signals from sample measurements by subtracting the measured absorbance at a specified baseline correction wavelength from the absorbance values at all wavelengths of a measured sample.
 - **Note**: If baseline calibration is not performed, the spectroscopic signals will not be separated from interference/background effects and will lead to inaccurate results.
- **Auto Blank** : When turned on, the instrument will automatically perform a blank detection the first time the pedestal is lifted and closed
 - **Note**: The automatic blank only takes effect when no blank readings have been made in a current run.
- **Auto Blank** : When turned on, the instrument will automatically perform sample measurements once the sample column has formed.
 - **Note**: Automatic sample detection only takes effect after a blank reading has been performed.
- The  icon appears in the upper right corner to indicate an error in reading/sample volumes. Please clean and wipe the pedestal and perform another blank reading. If the problem persists, contact VWR service.

Operation:

6. Set the Sample ID.
7. Clean the upper and lower pedestals with a wipe, add 2μL buffer solution to perform a blank reading.
8. Clean the buffer solution on the pedestals with a wipe.
9. Measure a 2μL sample volume and click "Measure" to detect the sample.

Note: A blank reading must be performed prior to sample measurements. The sample volume must be equivalent to the volume used to set the blank reading.

10. Clean and dry the pedestal between measurements.




Fig. 12 Protein Sample Results

The sample concentration and absorbance ratios will be displayed on the right side of the interface (**Fig. 13**).



Fig. 13 Sample Concentration & Absorbance Ratios


Conc.: Calculated protein concentration.

A260: The sample absorbance under 260nm (10mm pathlength equivalent).

A280: The sample absorbance under 280nm (10mm pathlength equivalent).

A260/A280: The ratio of corrected absorbance values at 260nm to the corrected absorbance values at 280nm. This ratio can be used as a secondary measure of nucleic acid purity. An A260/A280 ratio within the range of ~1.8 – 2.2 is pure for nucleic acids. Lower ratio values indicate the presence of contaminants that absorb strongly at or near 280nm.

wavelength: 280nm Absorb: 2.379 : Input a wavelength to view sample absorbance values (**Fig. 14**).



**Fig. 14 Protein Sample Absorbance Curve
(220nm – 350nm)**

Protein A280 Files Interface




Fig. 15 Protein Files Interface

Note: This interface is similar to the Nucleic Acids detection interface (see **Nucleic Acids Interface**).

Spectrum Data : To view sample absorbances at a specific wavelength, input a wavelength and the red coordinate line will display the corresponding absorbance value (**Fig. 16**).



Fig. 16 Protein Sample Absorbance Curve (220 – 350nm)

Colorimetry Interface

Introduction

BCA, Lowry, & Bradford protein assays are measured using colorimetric methods, which requires the generation of a standard curve.

BCA Protein Assay:

The BCA Assay is an alternative method for determining protein concentration. It utilizes bicinchoninic acid as a detection reagent to measure protein concentrations in crude samples. This assay measures protein solutions at 562nm and uses a standard curve to calculate the concentration. A baseline-correction can be used for normalization. Bicinchoninic acid is used to detect Cu^{+1} , which forms when Cu^{+2} is exposed to an alkaline environment. A purple chelate is formed when two BCA molecules react with 1 Cu^{+1} . The Cu-BCA chelate is measured at 562nm and baseline-corrected using a 750nm absorbance value.

Lowry Assay:

Folin-Ciocalteu is used as the detection reagent to measure protein concentrations in unpurified protein samples. This assay measures the absorbance of protein solutions at 650nm and uses a standard curve to calculate the concentration. A baseline-correction can be used for normalization. When protein solutions are exposed to cupric sulfate in an alkaline environment, tetradentate copper protein complexes form. Folin-Ciocalteu is reduced in proportion to the chelated copper-complexes which results in the formation of a water-soluble blue product that is measured at 650nm and baseline-corrected using a 405nm absorbance value.

Bradford Assay:

The Bradford Assay is an alternative method commonly utilized for determining protein concentration. It is often used for more dilute protein solutions where lower detection sensitivity is needed. Like the BCA and Lowry Assays, the Bradford Assay requires a standard curve before each run. The Bradford uses the protein-induced absorbance shift of Coomassie Blue dye at 595nm to measure protein concentration. The bound protein-dye complex is measured at 595nm and baseline-corrected using a 750nm absorbance value.

Note: Please perform measurements quickly. Coomassie dye-dye & dye-protein aggregates can form the longer it sits, which will result in fluctuations in absorbance readings.

Colorimetry Interface

A standard curve must be generated prior to sample measurement.




Fig. 17 Colorimetry Interface

BCA-562: Select the assay type: BCA-562, Bradford-595, Lowry-650.

Curve : Select the curve.

Operation:


1. Set the colorimetry and curve type.
2. Ensure that the upper and lower pedestals have been wiped clean. Pipette 0.5µL-2µL of buffer solution to perform a blank reading.
3. Wipe the buffer solution off the pedestals with a laboratory wipe.
4. Measure a 0.5µL-2µL sample. Select "Sample" to perform the reading. The sample volume must be equivalent to the volume used to set the blank reading. **Note:** A blank reading must be performed prior to sample measurements.
5. Clean and wipe the pedestals between each sample measurement.

Curve

To generate a standard curve, at least 5 standard sample concentrations are needed. The concentration range of the standards should cover all unknown sample concentrations.

Colorimetry Curve Interface:

Click the “Curve” tab to build a standard curve (**Fig.18**).




The screenshot shows the 'Colorimetry' software interface with the 'Curve' tab selected. The top navigation bar includes 'Colorimetry', 'Curve', 'Files', and 'Help'. Below the navigation bar, there are buttons for 'Blank', 'Sample', 'Save Curve', 'Display Curve', 'Import Curve', and 'Export Curve'. The main area displays a table with columns for 'Sample name', 'g/ml', 'Abs Avg', 'Absorb 1', 'Absorb 2', 'Absorb 3', 'Absorb 4', and 'Absorb 5'. The 'Sample name' column is currently empty, and the 'g/ml' column is set to 'g/ml'. The 'Abs Avg' column is set to '0.000'. The 'Absorb 1' through 'Absorb 5' columns are also set to '0.000'. The bottom status bar shows the VWR logo, 'Save conc.', 'Delete conc.', and the date/time '07/25/2023 01:22:55'.

Fig. 18 Colorimetry Curve Interface

Standard Curve Generation:

Click **New Curve**, input the curve name, and click “Confirm”. The standard curve table will be displayed.



The screenshot shows the 'Colorimetry' software interface with the 'Curve' tab selected. The top navigation bar includes 'Colorimetry', 'Curve', 'Files', and 'Help'. Below the navigation bar, there are buttons for 'Blank', 'Sample', 'Save Curve', 'Display Curve', 'Import Curve', and 'Export Curve'. The main area displays a table with columns for 'Sample name', 'g/ml', 'Abs Avg', 'Absorb 1', 'Absorb 2', 'Absorb 3', 'Absorb 4', and 'Absorb 5'. The 'Sample name' column is currently empty, and the 'g/ml' column is set to 'g/ml'. The 'Abs Avg' column is set to '0.000'. The 'Absorb 1' through 'Absorb 5' columns are also set to '0.000'. The bottom status bar shows the VWR logo, 'Save conc.', 'Delete conc.', 'Import conc.', and the date/time '07/25/2023 01:24:39'.

Fig. 19 Standard Curve Interface

Click to choose the unit for samples and to input the concentration.

Select a standard sample and then click “Blank” to perform a blank reading. Then, click “Measure” to measure the absorbance of the standard sample.

Each standard sample can be measured up to 5 times, and the average value can be used to build the standard curve. Users can delete the standard sample values or delete single values amongst the 5 measurements of an individual sample.

After measuring all standards, click **Save Curve** to save the curve.

Display Curve : View the standard curve as below:




Fig. 20 Standard Curve

- **Import Curve** : Import a standard curve from a USB flash drive.
- **Export Curve** : Export a standard curve to a USB flash drive.
- **Save conc** : Save the input standard sample concentration value.
- **Delete conc** : Delete the input concentration value.

Colorimetry Files Interface




Fig. 21 Colorimetry Files Interface

This interface is similar to the Nucleic Acids detection interface (see **Nucleic Acids Interface**).

BCA-562 : Choose the assay type and the detection data will be displayed.

UV-Vis Interface

Introduction

The UV-Vis interface allows the VWR Spectrophotometer to function as a conventional spectrophotometer. Sample absorbances will display from 200nm to 800nm. Samples with high absorbance values (10mm pathlength equivalent of 90A) can be measured directly.

UV-Vis Measurement




Fig. 22 UV-Vis Interface

Click **Blank** set a blank reading. After setting a blank, click **Enlarge Curve** enter the full wavelength interface of the blank sample.




Fig. 23 Blank Full Wavelength Interface (200 - 800nm)

Users can input up to 5 wavelengths before detection and the corresponding absorbance values will be displayed after detection (**Fig. 24**).

Load #	Wave	Abs.
1	230	0.000
2	260	0.000
3	280	0.000
4	562	0.000
5	595	0.000

Fig. 24 Select Wavelengths for Measurement

After blank detection, click **Sample** to perform a sample reading. After sample measurements, click **Enlarge Curve** to enter the full wavelength interface of the sample.




Fig. 25 Sample Full Wavelength Interface

Operation:

11. Set the Sample ID.
12. Clean the upper and lower pedestals, add 2 μ L buffer solution to perform a blank reading.
13. Clean the buffer solution on the pedestals.
14. Measure a 2 μ L sample volume and click “Measure” to detect the sample.

Note: A blank reading must be performed prior to sample measurements. The sample volume must be equivalent to the volume used to set the blank reading.

15. Clean and dry the pedestal between measurements.

UV-Vis Files Interface

Uv-Vis

Files

Help

Print

Spectrum Data

Export Load No.

Delete Load No.

Delete Sample ID

Export File

Sample ID	☐	Load #	W1	W2	W3	W4	W5	Time
230725_014148	☒	1	230/-0.00	260/-0.00	280/-0.00	562/0.00	595/0.00	07-25-2023 01:44:51

VWR

07/25/2023 01:45:16

Fig. 26 UV-Vis Files Interface

Spectrum Data: Click to enter the UV-Vis Full Spectrum Interface. Input a wavelength to view the corresponding sample absorbance value (**Fig. 27**).




Fig. 27 UV-Vis Full Spectrum Interface

OD600 Interface

Introduction

The OD600 function allows for measuring absorbance values at 600nm. It is commonly used for measuring the growth rate of bacterial cell cultures by measuring absorbance of the culture in growth media at 600nm. The Beer-Lambert equation is used to calculate the concentration (see **Nucleic Acids Interface**). This interface does not require the generation of a standard curve.

OD600 Measurement

Load #	OD600
	0.000

Fig. 28 OD600 Detection Interface

Operation:

1. Set the Sample ID.
2. Clean the upper and lower pedestals; add 2-3mL buffer solution to perform a blank reading.
3. Clean the buffer solution on the pedestals.
4. Measure a 2-3mL sample volume and click “Measure” to detect the sample. **Note:** A blank reading must be performed prior to sample measurements. The sample volume must be equivalent to the volume used to set the blank reading.
5. Clean and dry the pedestal between measurements.

OD600 Files Interface

Sample ID	Load #	OD600	Time
230725_014614	1	0.000	07-25-2023 01:47:21

Fig. 29 OD600 Files Interface

This interface is similar to the Nucleic Acids detection interface (see **Nucleic Acids Interface**).

Performance Verification

Materials Needed for verification:

- Lint Free Laboratory Wipes
- Deionized Water (diH₂O)
- Calibrated Precision Pipettor (2μL)
- PV-1 Solution



Expiration Date: 1 year from date of production

Label	Description
20230411	Production date.
260:2.1221	Absorbance(@260nm with 1mm pathlength):2.1221
302:0.9810	Absorbance(@302nm with 1mm pathlength):0.9810

Step by Step Usage Instructions:

Step #	Step Description	Fig.
Step 01	From the Home Screen, Tap the System icon	Fig 30
Step 02	Tap the PV Button	Fig 30
Step 03	Tap on the Std260(1mm) entry box to display a number keypad. Enter the 260 value of the PV Solution (found on PV-1 tube label)	Fig 31
Step 04	Repeat Step03 for the input value for Std302(1mm) .	Fig 31
Step 05	Pipette 2uL diH ₂ O onto pedestal, lower the arm, and tap Blank .	Fig 31
Step 06	Remove water from upper and lower pedestal surfaces using a clean, dry laboratory wipe.	
Step 07	Ensure PV-1 Solution is thoroughly mixed.	
Step 08	Pipette 2uL of the PV-1 solution onto the lower pedestal. Lower the arm and press Sample to begin the measurement.	Fig 31
Step 09	After the measurement is complete, remove sample from both upper and lower pedestal with a dry laboratory wipe.	
Step 10	Repeat Steps 08 and 09 to measure 9 additional individual replicates of the PV-1 solution. Tap List : the individual results will be displayed on screen. Note: If there is a significant error due pipetting, the individual measurements can be deleted via the List interface	Fig 32
Step 11	Analysis: Select Verification to analyze the measurements. Interpreting results: 0.2mm Avg*5/K1-Std302 ≤ Std302*0.013 , Pass 1.0mm Avg/K2-Std302 ≤ Std302*0.013 , Pass 0.2mm Std. Deviation ≤ 0.013 , Pass 1.0mm Std. Deviation ≤ 0.013 , Pass	Fig 32 Fig 33




Fig. 30




Fig. 31

The PV screen displays the following information:

- Buttons**: Blank, Sample, Export Data, Back
- Table**:

No.	0.2mm	1.0mm
1	0.20173	1.00606
2	0.21144	1.00177
3	0.20735	0.99545
4	0.20335	0.99516
5	0.20956	1.00493
6	0.21162	0.99650
Avg.	0.20873	1.00065
SD	0.00391	0.00504
CV	1.87236%	0.50339%
- Verification Data**:

Std260(1mm)	2.0753
Std302(1mm)	1.0112
- Verification Button**: Verification
- Footer**: VWR logo, 07/25/2023 02:03:04

Fig. 32




Fig. 33

Note: If any of the above (Fig. 33) results in a “FAIL”, please contact VWR Service.

System Settings Interface

Click “System” on the main interface to enter the System Settings Interface (**Fig. 35**).




Fig. 34 System Settings Interface

Date & Time Settings Interface

Click the “Time” icon to enter the date and time settings interface.




Fig. 35 Date and Time Settings Interface

Print

Click the “Print” icon to set the print mode. Select **Automatic Print** to enable automatic printing of results immediately after each measurement.




Fig. 36 Print Settings Interface

Brightness

Click the “Brightness” icon to enter the Brightness Settings Interface. Use the slider to adjust the brightness of the display.




Fig. 37 Brightness Settings Interface

Troubleshooting

Review the information in the table below to troubleshoot operating problems.

Problem	Cause	Solution
Measurement results not precise	Unformed sample column	Recondition the pedestal using the PR-1 kit.
	Calibration required	Perform a calibration verification check using the PV-1 kit. Contact VWR Service.
Instrument does not power on	Power supply is not plugged into instrument	Ensure power supply is plugged into outlet & instrument
	Defective power switch	Replace the power switch
	Defective power adapter	Contact VWR Service.
OD600 module failure	Poor connection between cable and board.	Contact VWR Service.
Insufficient light intensity error	Analysis module defective, optical fiber broken.	Contact VWR Service.
Touch Screen Error	Ineffective grounding of power supply	Provide effective grounding to power supply
Communication Timeout	Communication failure of analysis module	Contact VWR Service

Repair and Maintenance

Encountering problems

For inquiries on the repair and maintenance of the instrument, please contact your local VWR service center for assistance.

User Replaceable Accessories and Spare Parts

Description	Quantity	Cat. No.
VWR PRINTER PAPER SPECTROPHOTOMETER	3	76777-148
VWR CUVETTES GLASS SPECTROPHOTOMETER	2	76777-146

Technical service

Web Resources

Visit the VWR website at www.vwr.com for:

- Complete technical service contact information
- Access to the VWR Online Catalogue, and information about accessories and related products
- Additional product information and special offers

Contact us For information or technical assistance contact your local VWR representative or visit. www.vwr.com.

Warranty

VWR warrants that this product will be free from defects in material and workmanship for a period of two (2) years from date of delivery. If a defect is present, VWR will, at its option and cost, repair, replace, or refund the purchase price of this product to the customer, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear. If the required maintenance and inspection services are not performed according to the manuals and any local regulations, such warranty turns invalid, except to the extent, the defect of the product is not due to such non-performance.

Items being returned must be insured by the customer against possible damage or loss. This warranty shall be limited to the aforementioned remedies. IT IS EXPRESSLY AGREED THAT THIS WARRANTY WILL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND IN LIEU OF THE WARRANTY OF MERCHANTABILITY.

Compliance with local laws and regulations

The customer is responsible for applying for and obtaining the necessary regulatory approvals or other authorisations necessary to run or use the Product in its local environment. VWR will not be held liable for any related omission or for not obtaining the required approval or authorisation, unless any refusal is due to a defect of the product.

Equipment Disposal



This equipment is marked with the crossed out wheeled bin symbol to indicate that this equipment must not be disposed of with unsorted waste.

Instead it's your responsibility to correctly dispose of your equipment at lifecycle -end by handling it over to an authorized facility for separate collection and recycling. It's also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect from health hazards the persons involved in the disposal and recycling of the equipment.

For more information about where you can drop off your waste of equipment, please contact your local dealer from whom you originally purchased this equipment.

By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

Thank you

Our local offices in Europe, AMEA and North America

Austria

VWR International GmbH
Graumannsgasse 7
1150 Vienna
Tel.: +43 1 97 002 0
Email: info.at@vwr.com

Belgium

VWR International BV
Researchpark Haasrode 2020
Geldenaaksebaan 464
3001 Leuven
Tel.: +32 016 385 011
Email: vwr.be@vwr.com

Canada

VWR International
2360 Argentia Road
Mississauga, Ontario L5N 5Z7
Tel.: +1 800 932 5000
Email: Canada_Orders@vwr.com

China

VWR (Shanghai) Co., Ltd
Bld.No.1, No.3728 Jinke Rd,
Pudong New District
Shanghai, 201203- China
Tel.: +400 821 8006
Email: info_china@vwr.com

Czech Republic

VWR International s. r. o.
Veetee Business Park
Pražská 442
CZ - 281 67 Stránská Skalice
Tel.: +420 321 570 321
Email: info.cz@vwr.com

Denmark

VWR International A/S
Tobaksvejen 21
2860 Søborg
Tel.: +45 43 86 87 88
Email: info.dk@vwr.com

Finland

VWR International Oy
Valimotie 17-19
00380 Helsinki
Tel.: +358 9 80 45 51
Email: info.fi@vwr.com

France

VWR International S.A.S.
Immeuble Estréo
1-3 rue d'Aurion
93114 Rosny-sous-Bois cedex
Tel.: 0 825 02 30 30* (national)
Tel.: +33 (0) 1 45 14 85 00
(international)
Email: info.fr@vwr.com
* 0,18 € TTC/ min + prix appel

Germany

VWR International GmbH
Hilpertstrasse 20a
D - 64295 Darmstadt
Free call: 0800 702 00 07
Tel.: +49 (0) 6151 3972 0
(international)
Email: info.de@vwr.com

Hungary

VWR International Kft.
Simon László u. 4.
4034 Debrecen
Tel.: +36 (52) 521-130
Email: info.hu@vwr.com

India

VWR Lab Products Private Limited
No.139. BDA Industrial Suburb,
6th Main, Tumkur Road,
Peenya Post,
Bangalore, India – 560058
Tel.: +91-80-28078400
Email: vwr_india@vwr.com

Ireland

VWR International Ltd
Orion Business Campus
Northwest Business Park
Ballycoolin
Dublin 15
Tel.: +353 01 88 22 222
Email sales.ie@vwr.com

Italy

VWR International S.r.l.
Via San Giusto 85
20153 Milano (MI)
Tel.: +39 02-3320311/02-487791
Email: info.it@vwr.com

Korea

Avantor Performance Materials Korea Ltd
2F ACE Gwanggyo Tower I
Daehak 4ro 17
Yeongtong-gu Suwon
Korea 16226
Tel.: +82 31 645 7250
saleskorea@avantorsciences.com

Mexico

VWR International, S.de R.L. de C.V.
Km. 14.5 Carretera
Tlalnepantla-Cuautitlán
Col. Lechería
Tultitlán Edo. de México
CP 54940
Tel.: +52 (55) 5005 0100
vwrmx@vwr.com

Middle East & Africa

VWR International FZ-LLC
Office 203, DSP Lab Complex,
Dubai Science Park
Dubai, United Arab Emirates
Tel.: +971 4 5573271
info.mea@vwr.com

The Netherlands

VWR International B.V.
Postbus 8198
1005 AD Amsterdam
Tel.: +31 020 4808 400
Email: info.nl@vwr.com

Norway

VWR International AS
Brynsalleen 4
0667 Oslo
Tel.: +47 22 90 00 00
Email: info.no@vwr.com

Poland

VWR International Sp. z o.o.
Limbowa 5
80-175 Gdansk
Tel.: +48 058 32 38 210
Email: info.pl@vwr.com

Portugal

VWR International – Mat. de
Laboratório, Soc. Unipessoal, Lda
Edifício Ramazzotti
Avenida do Forte 6, P-1.09 e P-1.10
2790-072 Carnaxide
Tel.: +351 21 3600 770
Email: info.pt@vwr.com

Singapore

VWR Singapore Pte Ltd
The Metropolis
Tower 1, #05-03
9 North Buona Vista Drive
Singapore 138588
Tel: +65 6505 0760
Email: sales.sg@vwr.com

Spain

VWR International Eurolab S.L.
C/ Tecnología 5-17
A-7 Llinars Park
08450 - Llinars del Vallès
Barcelona
Tel.: +34 902 222 897
Email: info.es@vwr.com

Sweden

VWR International AB
Fagerstagatan 18b
163 94 Stockholm
Tel.: +46 08 621 34 20
Email: info.se@vwr.com

Switzerland

VWR International GmbH
Lerzenstrasse 16/18
8953 Dietikon
Tel.: +41 044 745 13 13
Email: info.ch@vwr.com

UK

VWR International Ltd
Customer Service Centre
Hunter Boulevard
Magna Park
Lutterworth
Leicestershire
LE17 4XN
Tel.: +44 (0) 800 22 33 44
Email: uksales@vwr.com

United States

VWR International, LLC
100 Matsonford Road
Building One Suite 200
Radnor, PA 19087
Tel.: +1 800 932 5000
vwrcustomerService@vwr.com