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1. Safety

1.1 General

- These instructions facilitate the safe and efficient handling of the equipment.
- Before operating the device for the first time the operator must have read and understood the operating manual.
- Adherence to all the specified safety instructions and instructions for actions in this operating manual is a fundamental requirement for working safely.
- In addition, all local accident prevention guidelines and general safety regulations must be observed.
- Store these operating instructions in a well accessible location.

1.2 Safety Symbols

	DANGER - Indicates an imminently hazardous situation, which, if not avoided, will result in death.
	WARNING - Indicates a potentially hazardous situation, which, if not avoided, can result in death or serious injury.
	CAUTION - Indicates a potentially hazardous situation, which, if not avoided, can result in minor or moderate injury.
	NOTICE – Is used to address practices not related to physical injury.
	Safety alert symbol is a symbol that indicates a general hazard. Failure to follow these instructions could lead to personal injury or damage to equipment.
	This symbol indicates a potential danger caused by hot surface.
	This symbol refers to jeopardizing persons by electric power.

1.3 General cause of risks

Product- and media safety

 DANGER	1. Explosion hazard! Danger of death due to explosion! • Do not use the device in an unsafe environment and especially not in explosive atmospheres. It is not EX-protected. • Do not heat explosive samples or objects.
 WARNING	1. Fire Hazard! Burns caused by fire! • The max. working temperature may not be over or just under the flash point of the media. The heating element can act as an ignition source. • Beware of hazards due to: - chemical reactions - flammable materials - combustible media with a low boiling temperature. 2. Risk of accident and health hazard! Risk of injury caused by hazardous substances! • Use a fume hood to avoid coming into contact with toxic or caustic fumes. • All relevant national regulations regarding work safety and accident prevention must be observed when working with hazardous, toxic or pathogenic substances.
 CAUTION	1. Hot surfaces! Contact with the skin causes burns! • Do not touch the upper housing parts or the crucibles during or immediately after operation. • During use the tray openings without crucibles are directly exposed to radiant heat. Avoid any direct contact. • Pay attention to the residual heat after switching off. • Do not directly touch hot media, substances or objects. • Avoid unsafe standing position of crucibles or incorrect container size.

Electrical Safety

Dangerous electrical voltages!

Contact with dangerous voltages can result in life threatening electric shock.

- Only operate the appliance in an environment, where dirt (contamination) is nonconductive.
- Do not use the device if:
 - it is damaged
 - cable is damaged.
- Use the device only with the protective grounding conductor connected.
- Never insert objects through slots, vent holes or the tray openings.
- Only remove the crucibles when the device has been switched off by the time switch.
- Arrange the cables so that the electrical wires do not come in to contact with high-temperature parts.

1.4 Intended use


The device can be used for incineration of **non-explosive**, solid and liquid substances. This also includes the preparatory steps such as pre-heating, drying and fuming-off. The sample material must be placed into a crucible made of quartz glass, porcelain or platinum.

Use the device only as described in these Operating Instructions. Any other usage is considered to be improper and may lead to material damages and even personal injuries. Any selfconversion and alteration undertaken on the device is not permitted. The manufacturer does not bear any liability for damages caused due to improper use.

Use the device in indoor environments like that of a laboratory of research, teaching, industry or trade.

2. Configuration and function

2.1 Overview

Picture: Side view of a rapid incinerator

2.2 Function

The two-storied rapid incinerators are used for time-saving single or series incineration of solid and liquid substances.

The device is suitable for taking up porcelain crucibles of up to approx. 50 mm diameter. Smaller crucibles are always more advantageous, as bright annealing is much quicker. The perforated tray with 8 openings takes up the crucibles for pre-heating, drying and fuming off of also liquid crucible content. Ashing occurs in the lower, 400 mm long triangular heating bed. The uninhibited exposure to oxygen speeds up the incineration process. The inclined position of the crucibles allows the ashing process to be easily observed. The perforated tray and the ashing bed can be simultaneously charged with crucibles.

2.3 Function of the automatic temperature regulation

An integrated NiCr-Ni mantle thermocouple measures the temperature at the heating element and generates the corresponding switching command to the temp. display/control unit. The microprocessor calculates the required heating power and releases the according heating current. As a precaution the heating current is cut-off by opening the flap. Consequently, the temperature decreases and the microprocessor releases more power. If the decrease was only small, the very fast heating may cause superheating of few °C. The controller balances the value immediately by reducing the heating power and the adjustment to the setpoint is 1-2 °C and better. Such superheating of short duration does not damage the samples (inertia of the crucibles).

3. Operating the device

3.1 Set up

The device should be placed into a fume hood or a sufficiently large and dry room with enhanced venting facilities.

Choose a solid ground where radiant heat cannot cause any damage to sensitive surfaces. The distance from walls sensitive to heat or combustible must amount to a minimum of 15 cm. The area needs to be even, solid, nonslip and clean.

Caution: Electric shock hazard! Unwind the power cord completely and make sure that it is safe and do not touch any hot parts of the device.

Compare specified voltage on hotplate label with available mains voltage and insert power plug of device into an earthed socket.

3.2 Switching on

Rapid incinerators are equipped with a time switch up to 120 min running time. The device gets enabled for the corresponding duration by turning the knob clockwise. Further the switch clock is equipped with a controller position "Duration-On" (without time sequence). This is activated by turning the knob counterclockwise.

When the operation flap is opened, current to the heating element is disconnected. Reconnection results automatically on closing the flap. This procedure can be observed by listening to a clicking sound (contactor). Check this safety cut-off function after each switching on by opening and closing the flap several times.

3.3 Setting of the desired temperature

Caution: Risk of burns! Do not touch the upper metal parts or crucibles during use. It can become very hot. During use the tray openings without crucibles are directly exposed to radiant heat. Avoid any direct contact.

If the *-key is not pressed the display shows the actual temperature. Pressing it, the display changes to show the setpoint which increases resp. decreases in 1 °C by pressing the *-key and touch the Up or Down key. Keeping the two keys pressed causes a high-speed adjustment after a few seconds. To store the setting, release the Up or Down key before the *-key. The new value is loaded into the non-volatile memory and remains there even after switching off or power cut.

Reference values for loss of heat within the crucibles			
Set temperature [°C]	Temperature inside lower crucible [°C]	Temperature inside upper crucible [°C]	The temperature was measured in two enamelled crucibles of 45mm dia. and 40mm height. The one in the ashing bed was half filled with fire-clay powder and measured at the bottom by a thermocouple. The
100	60	38	
200	160	76	

300	246	112	bottom of the upper crucible was measured with a special surface sensor. Notice: The loss of heat decreases by covering the 8 upper recesses with crucibles.
400	350	160	
500	440	205	
600	545	262	
700	645	322	
800	748	355	
900	846	436	
950	895	474	
999	944	496	

3.4 Controlling of the adjusted temperature

Press the *-key only. The display shows °C and the adjusted temperature by turns.

3.5 Switching off

Caution: Residual heat! Please note, that the crucibles and the upper housing parts will remain hot after turning it off.

1) Operating as a time switch: The device switches itself off automatically after the entered period of time. A short acoustic signal will follow. Furthermore, the unit can be switched off directly by turning the power controller anticlockwise to the zero position.

2) Operation mode „Duration-On“: Turn the time-switch knob clockwise to zero position.

Notice: The setpoint remains into the non-volatile memory even after switching off or power cut. In order to avoid heating after switching on the device again, the set temperature must be set to 0°C.

4. Failures

Warning: Risk of electric shock and risk of burns! In case of failures or damages, disconnect the device from the mains immediately and allow it to cool.

Because rapid incinerator has been designed in a simple and uncomplicated way, replacing of all spare parts can be done by a qualified electrician. Otherwise please return the device to the supplier.

Error messages

Error	Cause	Correction
The RCD trips.	Humidity in the heating element.	1. Replace the complete heating element (Ref.: EH95) by a qualified electrician. 2. Have the unit repaired by the manufacturer.

1. The device heats up while the operation flap is opened or 2. The device does not heat.	1+2. Bracket of the safety switch has been bent. 2. Heating element or electronics defective.	1+2. When opening or closing the operation flap, it can be heard a click of the contactor. If this is not the case, the metal bracket may be bent/blocked (for e.g. due to improper transport or unpacking). It can be adjusted by carefully bending it. In doing so, it must be ensured that the heating element is still <u>turned off immediately when opening</u> the flap (confirmed by another clicking sound of the contactor). This adjustment may only be done by a qualified electrician and after having the device disconnected from the mains supply. 2. Replace the unit by a qualified electrician or have the device repaired by the manufacturer.
<i>Display flashing:</i> inPt: FAiL	Broken sensor, sensor short circuit	1. Replace the complete heating element incl. thermocouple (Ref.: EH95) by a qualified electrician. 2. Have the device repaired by the manufacturer.
<i>Display flashing:</i> dAtA : FAiL	Non-volatile memory error.	1. De-power briefly. Replace unit by a qualified electrician if problem persists. 2. Have the device repaired by the manufacturer.

If you experience an error that is not listed in this table, please contact your supplier. Please do not try to eliminate it by yourself. Please only send devices in for repair that have been cleaned and are free of materials which might present health hazards.

5. Maintenance and cleaning

Warning: Risk of electric shock and risk of burns! Disconnect device from main power and let it cool down before starting to do maintenance work or cleaning.

Apart from the periodic tests of the safety cut-off function (Chapter 3.2), there is no special maintenance necessary. Please avoid liquid or solid stuff getting into the heating element. This would lead to heavy wastage as well as badly working outlets. Heavy acid exposure in an improperly functioning fume hood can considerably reduce the life span and leaves the devices looking unsightly.


6. Disposal

The meaning of the symbol of the crossed-out bin: For a sustainable and professional disposal of electronical devices the user is obligated to dispose legacy equipment separately from general waste. The device marked with the above symbol must be given to an appropriate facility for recycling.

In case of doubt, obtain information on environmentally sound disposal from the local municipal authority or special disposal companies.

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Harry Gestigkeit GmbH