VWR® Handheld Refractometer

INSTRUCTION MANUAL

European Catalogue Numbers:

635-0629	0-10% Brix
635-0630	0-18% Brix
635-0631	0-32% Brix
635-0632	28-62% Brix
635-0633	45-82% Brix
635-0634	0-55% Brix
635-0635	58-92% Brix
635-0636	0-10% Salinity
635-0637	0-28% Salinity
635-0638	Battery -50-0°C
635-0639	Oechsle 0-140° KMW Babo 0-27°

Version: 1.0

Issued:

15 January 2014

Legal address of manufacturer

Europe

VWR International bvba

Researchpark Haasrode 2020

Geldenaaksebaan 464

B-3001 Leuven

+ 32 16 385011


<http://be.vwr.com>

Country of origin

Made in China

VWR Handheld Refractometers

The handheld refractometers are optical instruments to measure the refractive index and correlated scales. The units have a ATC (automatic temperature compensation) 10 to 30 °C to reference temperature of 20 °C.

Measuring a sample

Step 1

Open daylight plate and place 2 to 3 drops of distilled water or standard solution on the prism. Close the daylight plate. Sample or standard should create a close film between prism and daylight plate without bubbles or dry spots. Wait for at least 30 seconds before going to step 2 to allow the sample to have the same temperature as the prism surface (for a proper temperature compensation and accurate results).

Step 2

Hold daylight plate in the direction of a light source and look into the eyepiece. You will see a circular field with graduations down the centre (you may have to focus the eyepiece to clearly see the graduations). The upper portion of the field should be blue, while the lower portion should be white. You can read the value exactly at the boundary between blue and white field.

Adjustment

We recommend an adjustment before first measurement using distilled water. Look into the eyepiece and turn the calibration screw until the boundary between the upper blue field and the lower white field meets exactly on the zero scale. Adjustment finished. Please pay attention: The adjustment has to be performed at 20 °C room temperature. Once adjusted shifts in ambient temperature within the range 10 to 30 °C are compensated by the ATC function to 20 °C.

Instrument check

Perform step 1 and step 2 using a standard with known value. The unit should show the value of the standard at 20 °C within the main unit and standard accuracy. If the instrument does not show the right value please repeat the adjustment step.

Models

Brix refractometers are used to measure sugar in liquids like fruit juices, honey, soft drinks, wine and other sugar containing samples.

Salinity refractometers are designed for testing the concentration of salt water and brine.

Oechsle (°Oe) refractometer are specially designed for measuring Brix(%), °Oe and KMW of grapes and grape juice.

Battery/anti-freeze refractometers are designed to indicate the status of the battery or motor cooling system in terms of temperature tolerance.

Technical specifications

Measuring range	ATC (to 20 °C)	Resolution	Accuracy	Cat. No.
0 -10% Brix	Yes	0,1%	±0,1%	635-0629
0 -18% Brix	Yes	0,1%	±0,1%	635-0630
0 - 32% Brix	Yes	0.2%	±0,2%	635-0631
28 - 62% Brix	Yes	0,2%	±0,2%	635-0632
45 - 82% Brix	Yes	0,5%	±0,5%	635-0633
0 - 55% Brix	Yes	0.2%	±0.2%	635-0634
58 - 92% Brix	Yes	0,2%	±0,2%	635-0635
0 - 10% Salinity	Yes	0.1%	±0.1%	635-0636
0 - 28% Salinity	Yes	0,2%	±0.2%	635-0637
Battery/anti freeze -50 - 0 °C	Yes	5 °C	±5 °C	635-0638
Oechsle 0 - 140 ° KMW Babo 0 - 27 °	Yes	1Oe 0,2 ° KMW Babo	±1Oe ±0.2° KMW Babo	635-0639

Maintenance

Accurate measurement depends on careful adjustment .Please follow the instructions above. Prism and sample must be at the same temperature for accurate results.

Do not expose the instrument to harsh environment. Do not immerse the instrument in water. If the instrument becomes foggy, water has got inside the unit. We recommend the unit is replaced.

Do not measure abrasive or corrosive chemicals with this instrument. They may damage the prism's coating.

Clean the instrument between each measurement carefully and properly using water and a soft cloth otherwise following measurements maybe incorrect or the prism may cause damage.

This is an optical instrument. It requires careful handling and storage. Any harsh environment or inappropriate handling may cause damage. the optical components.

Technical services

Web resources

Visit the VWR website at **vwr.com** for:

- Complete technical service contact information
- Access to the VWR Online Catalogue and information about accessories and related products
- Additional product information and special offers

Contact us: For information or technical assistance contact your local VWR representative or visit **www.vwr.com**.

Warranty

VWR International warrants that this product will be free from defects in material and workmanship for a period of two (2) years from date of delivery. If a defect is present, VWR will, at its option and cost, repair, replace, or refund the purchase price of this product to the customer, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear. If the required maintenance and inspection services are not performed according to the manuals and any local regulations, such warranty turns invalid, except to the extent, the defect of the product is not due to such non performance.

Items being returned must be insured by the customer against possible damage or loss. This warranty shall be limited to the aforementioned remedies. IT IS EXPRESSLY AGREED THAT THIS WARRANTY WILL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND IN LIEU OF THE WARRANTY OF MERCHANTABILITY.

Compliance with local laws and regulations

The customer is responsible for applying for and obtaining the necessary regulatory approvals or other authorisations necessary to run or use the product in its local environment. VWR will not be held liable for any related omission or for not obtaining the required approval or authorisation, unless any refusal is due to a defect of the product.

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Local VWR offices in Europe and Asia Pacific

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India

VWR Lab Products Private Limited
135/12, Brigade Towers, 2nd Floor
Front wing, Brigade Road,
Bengaluru, India – 560 025
Tel: +91-2522-647911/922 (Mumbai)
Tel: +91-80-41117125/26 (Bangalore)
Fax +91-80-41117120
E-mail: vwr_india@vwr.com

Ireland / Northern Ireland

VWR International Ltd / VWR International
(Northern Ireland) Ltd
Orion Business Campus
Northwest Business Park
Ballycoolin
Dublin 15
Tel.: 01 88 22 222
Fax: 01 88 22 333
E-mail sales@ie.vwr.com

Italy

VWR International PBI S.r.l.
Via San Giusto 85
20153 Milano (MI)
Tel.: 02-3320311/02-487791
Fax: 02-332031307/02-40090010
E-mail: info@it.vwr.com

The Netherlands

VWR International B.V.
Postbus 8198
1005 AD Amsterdam
Tel.: 020 4808 400
Fax: 020 4808 480
E-mail: info@nl.vwr.com

New Zealand

VWR International LP
241 Bush Road
Albany
0632 New Zealand

Norway

VWR International AS
Haavard Martinsens vei 30
0978 Oslo
Tel.: 02290
Fax: 815 00 940
E-mail: info@no.vwr.com

Poland

VWR International Sp. z o.o.
Limbowia 5
80-175 Gdansk
Tel.: 058 32 38 210
Fax: 058 32 38 205
E-mail: labart@pl.vwr.com

Portugal

VWR International - Material de
Laboratório, Lda
Edifício Neopark
Av. Tomás Ribeiro, 43- 3 D
2790-221 Carnaxide
Tel.: 21 3600 770
Fax: 21 3600 798/9
E-mail: info@pt.vwr.com

Singapore

VWR Singapore Pte Ltd
18 Gul Drive
Singapore 629468
Tel: +65 6505 0760
Fax: +65 6264 3780
E-mail: sales@sg.vwr.com

Spain

VWR International Eurolab S.L.
C/ Tecnología 5-17
A-7 Llinars Park
08450 - Llinars del Vallès
Barcelona
Tel.: 902 222 897
Fax: 902 430 657
E-mail: info@es.vwr.com

Sweden

VWR International AB
Fagerstagatan 18a
163 94 Stockholm
Tel.: 08 621 34 00
Fax: 08 621 34 66
E-mail: info@se.vwr.com

Switzerland

VWR International GmbH
Lerzenstrasse 16/18
8953 Dietikon
Tel.: 044 745 13 13
Fax: 044 745 13 10
E-mail: info@ch.vwr.com

Turkey

Pro Lab Laboratuvar Teknolojileri Ltd.Şti.
a VWR International Company
Orta Mah. Cemal Gürsel Caddesi
Ördekcioglu İşmerkezi No.32/1
34896 Pendik - Istanbul
Tel.: +90216 598 2900
Fax: +90216 598 2907
Email: info@pro-lab.com.tr

UK

VWR International Ltd
Customer Service Centre
Hunter Boulevard
Magna Park
Lutterworth
Leicestershire
LE17 4XN
Tel.: 0800 22 33 44
Fax: 01455 55 85 86
E-mail: uksales@uk.vwr.com