**INSTRUCTION MANUAL
MANUEL D'INSTRUCTIONS
BETRIEBSANLEITUNG
MANUALE D'ISTRUZIONI
MANUAL DE INSTRUÇÕES
MANUAL DE FUNCIONAMIENTO**

European Catalogue Numbers:

VisiCam 1.3: 630-1030

VisiCam 3.0: 630-1031

VisiCam 5.0: 630-1032

VisiCam 10.0: 630-1484

Version 6

Issued 21.10.2015

Made in PRC

QC in ITALY

Legal Address of Manufacturer

Europe

VWR International BVBA
Researchpark Haasrode 2020
Geldenaaksebaan 464
B-3001 Leuven
+ 32 16 385011
<http://www.be.vwr.com>

Warning & Safety Information

Do not disassemble the camera. Do not try to cut the USB cable.

Package Contents

Description	Quantity
VisiCam USB camera	1
0.5x eyepiece adapter	1
30mm ring for eyepiece adapter	1
calibration micrometric slide	1
CD-ROM with drivers and software for Windows XP/VISTA/7/8)	1

Product Specifications

- USB 2.0 connection. Plug & Play.
- Dimensions (L) 85mm X (W) 65mm X (H) 42mm.
 - VisiCam 1.3: 1/2.5" CMOS sensor, 1.3 Mpixel. Max resolution: 1280x1024 pixels.
 - VisiCam 3.0: 1/2" CMOS sensor, 3.0 Mpixel. Max resolution: 2048x1536 pixels.
 - VisiCam 5.0: 1/2.5" CMOS sensor, 5.0 Mpixel. Max resolution: 2592x1944 pixels.
 - Visicam 10.0: 1/2.3" CMOS sensor, 10.0 Mpixel. Max resolution: 3856x2764 pixels.
- No need of batteries or power supply. Powered by USB
- White Balance: automatic/manual (switchable by software).
- Real time image and colour quality control.

Recommended System

- IBM PC compatible: Windows 2000 / XP / VISTA / Windows7 / Windows 8
- RAM: 512 MB , HDD: at least 10GB
- USB 2.0 interface
- CD-ROM drive (to install drivers & software)

Intended use

VisiCam cameras can be mounted on all Biological and Stereo Microscopes equipped with trinocular tube.

It can be mounted also on monocular and binocular Microscopes thanks to the eyepiece adapter. It allows an easy connection to any USB 2.0 port of a personal computer to acquire, store and edit images from your Microscope.

Installation

Installation of Camera and Drivers

Remove all the parts from their package. Remove the camera protective cap (leave the black ring in position).

30 mm Ring

0.5x Eyepiece Adapter

Knurled Ring C-CS

If you have a trinocular microscope then screw the C-Mount adapter (not included, please contact a VWR sale responsible) to the camera black ring. Put the C-Mount screwed to the camera onto your microscope trinocular tube.

If you have a binocular or monocular microscope, then screw the 0.5x eyepiece adapter (included) to the camera black ring. Put the eyepieces eyepiece adapter screwed to the camera onto your microscope binocular or monocular tube (you must remove one eyepiece).


NOTE: If the diameter of your microscope eyepiece is 30mm diameter (for example as in stereo Microscopes), then you must use the 30mm ring (included) together with the 0.5x eyepiece adapter.

NOTE: the eyepiece adapter has a knurled ring included. It is not necessary to mount this one between the camera and the adapter if a similar knurled grid is already mounted on the camera.

1. Insert the installation CD-ROM into your CD-ROM driver.
2. Connect the camera to your PC USB port

3. Windows will detect the new hardware and display the installation wizard window. Select "No, not this time" and click "Next".
4. In Windows systems it might happen that the "Compatibility with Windows XP" error during software installation (see picture). Click on continue to continue the installation (this will not affect your system)

5. Select "Install the software automatically" when prompted for the installation CD (if driver files have been copied to another location on your PC, please select "installation from a specific location" and indicate the location where you copied the driver files
6. When the installation is over, click on "Finish". Once this step is terminated, you can remove the CD-ROM from the drive. VisiCam drivers are now installed on your personal computer.

7. Restart the system if requested.

Installation of "IS VisiCam Image Analyser" Software

VisiCam Image Analyser Software allows you working with VisiCam 1.3, VisiCam 3.0, VisiCam 5.0 and VisiCam 10.0 cameras. You can:

- Preview and Acquire images from camera
- Acquire videos
- Control main parameters of camera (colours balance, contrast, white balance...)
- Edit images you acquired (rotate, crop, measure, count....)

Getting started with “IS VisiCam Image Analyser” software

1. Connect the camera and start the software


VisiCam Image Analyser Preview interface

2. The camera is automatically detected and the preview window opens

3. Choose Preview resolution to change preview resolution and Capture resolution to change acquire resolution

4. Click on Wbalance to change white balance options and to make white balance of the image

5. When you are ready to acquire image, Click on Capture to choose acquiring options.

When the image is captured, you can go to the tab “Image” on top of windows and select the image you acquired, so that you can start editing it

Troubleshooting

How to verify if VisiCam camera is properly installed?

Verify That the Device Has Been Loaded in the Devices Manager of your Computer. When the camera is correctly detected by the system, it will loaded in the Devices Manager. To verify that the camera's hardware is being properly detected, please do the following:

- (1) Click "Start", "Settings", and then "Control Panel".
- (2) Double click on the "System" icon.
- (3) Click on the "Device Manager" tab. Windows 2000 users, click on the Hardware Tab, then the Device Manager button.
- (4) Click on the plus (+) sign next to "**Imaging Device**". A properly installed camera will load an entry in this section entitled "Camera Type". If this entry is listed with a yellow exclamation mark (!), proceed to the next step. If there is no entry, please redo camera Drivers installation procedure from the beginning.
- (5) Right Click on the Camera name and select Properties
- (6) Click on "Driver" Tab and Click on "update drivers"
- (7) Driver installation procedure will start again. Follow instructions as in the Driver installation section of this manual.

Camera is connected to the PC but no image is displayed in the capture interface

- (1) Check the VisiCam proper installation (see previous question) and proper PC connection. Check the microscope is turned on.
- (2) If the problem is still present connect the camera to a different USB port if you have more than one. It is important camera is connected to the MAIN USB port and to a secondary one (secondary USB ports are usually in front of computers)
- (3) Install again drivers if system requires it
- (4) If the problem does not resolve restart your computer and try again to launch software

When the capture interface starts with "no video capture device" message

- (1) The camera is not properly connected; check the USB connection and proper driver installation. Check if there is no other camera installed in the computer (if you have a laptop with camera embedded, you must disable it)

Black Preview

- (1) Click in "Exposure" button on software interface and adjust Exposure value or the Brightness and/or contrast settings to see if the image brightens
- (2) Check that your USB port is USB2.0 and not 1.1; if you have a 1.1 try to reduce Preview resolution of camera on Resolution menu of preview window

Colours in Preview window are not properly displayed.

Perform white balance with microscope lamp switched on and without any specimen under microscope

objective

The image in the Preview window appears overexposed even if the camera exposure is set at minimum

- (1) Insert grey filters on the microscope illuminator to reduce the amount of light reaching the camera. If you reduce illumination level of microscope, consider that colour regulation could be altered and you will have to perform again a white balance.

Maintenance

When not in use, please store the camera in a dry place and far from hot sources.

Avoid touching the sensor/filter surface with fingers, and other objects. This may damage the sensor or let dust go over it.

To clean the filter over the sensor it is enough to use a swab stick with 100% alcohol

Do not try to open the camera by yourself, this will void the warranty

User replaceable accessories and spare parts

Description	Qty	Cat. No.
C-mount Adapter for Visiscope TL324	1	630-1688
C-mount Adapter for Visiscope TL524	1	630-1749
C-mount Adapter for Visiscope SZT350	1	630-1793
C-mount Adapter for Visiscope SZTL350	1	630-1793
CS-C Mount adapter	1	630-1122
Calibration Slide	1	630-1123
TV/USB converter kit	1	630-1125
C-MOUNT ADAPTER 0,5x FOR OLYMPUS *MICROSCOPES	1	630-1853
C-MOUNT ADAPTER 0,5x FOR ZEISS* MICROSCOPES	1	630-1854
C-MOUNT ADAPTER 0,5x FOR LEICA* MICROSCOPES	1	630-1855
C-MOUNT ADAPTER 0,5x FOR NIKON* MICROSCOPES	1	630-1856

*not for all models

Technical service

Web Resources

Visit the VWR's website at vwr.com for:

- Complete technical service contact information
- Access to VWR's Online Catalogue, and information about accessories and related products
- Additional product information and special offers

Contact us For information or technical assistance contact your local VWR representative or

visit **vwr.com**.

Warranty

VWR warrants that this product will be free from defects in material and workmanship for a period of two (2) years from date of delivery. If a defect is present, VWR will, at its option and cost, repair, replace, or refund the purchase price of this product to the customer, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear. If the required maintenance and inspection services are not performed according to the manuals and any local regulations, such warranty turns invalid, except to the extent, the defect of the product is not due to such non performance.

Items being returned must be insured by the customer against possible damage or loss. This warranty shall be limited to the aforementioned remedies. IT IS EXPRESSLY AGREED THAT THIS WARRANTY WILL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND IN LIEU OF THE WARRANTY OF MERCHANTABILITY.

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VWR Lab Products Private Limited
VWR Singapore Pte Ltd