VWR

Microscope  **VisiScope®** series 384

INSTRUCTION MANUAL

Model	European Catalogue Number
T384 FL	630-1937

Version: 2
Issued: 19, 03, 2015

Legal Address of Manufacturer

Europe

VWR International bvba
Researchpark Haasrode 2020
Geldenaaksebaan 464
B-3001 Leuven
+ 32 16 385011
<http://be.vwr.com>

Country of origin: **ITALY**

Table of Contents

Warning
Safety Information
Package Contents
Unpacking
Intended use
Symbols and conventions
Product specifications
Overview
Using the microscope
Using the fluorescence
Troubleshooting
User replaceable accessories and spare parts
Technical service
Warranty
Compliance with local laws and regulations
Disposal

Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use.

We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users.

We decline any responsibility deriving from incorrect instrument use uses that does not comply with this manual.

Safety Information

Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position.

Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely.

Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

Package Contents

DESCRIPTION	QUANTITY
Microscope stand with nosepiece, stage, condenser	1
Optical head (trinocular)	1
Objective E-PLAN IOS 4x	1
Objective E-PLAN IOS 10x	1
Objective E-PLAN IOS 20x	1
Objective E-PLAN IOS 40x	1
Objective E-PLAN IOS 100x (oil)	1
Eyepiece WF10x/20mm	2
Fluorescence attachment	1
Fluorescence external power supply	1
Mercury lamp house	1
Dust cover	1
Power supply output 6Vdc	1

Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take care to ensure that the optical items (objectives and eyepieces) do not fall out and get damaged. Using both hands (one around the neck and one around the base), lift the microscope from the container and put it on a stable desk.

Place the observation head onto the top of the neck and tighten the lock-screw. Insert the eyepieces into the eye tubes.

Connect the provided power supply to the power supply input jack on the rear of the microscope.

Intended use

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

Symbols and conventions


The following chart is an illustrated glossary of the symbols that are used in this manual.


	CAUTION This symbol indicates a potential risk and alerts you to proceed with caution
---	---

Product specifications

Head:	Trinocular, 30° inclined, 360° rotating. Interpupillary adjustment 48-75 mm.
Eyepieces:	WF10X/20mm.
Nosepiece:	5-position reversed revolving nosepiece. Ball bearing rotation.
Objectives:	IOS E-PLAN 4x/0.10, 10x/0.25, 20x/0.40, 40x/0.65, 100x/1.25 (oil).
Focusing system:	Coaxial coarse and fine.
Stage:	Double layer mechanical sliding stage, dimensions 216x150mm, moving range 78x54mm.
Condenser:	Belt-drive in X direction..
Illumination:	Abbe condenser, sliding-in, N.A. 1.25 with centering system.
	Transmitted light: P-LED ³ , with manual brightness control.
	Epi-fluorescence: HBO 100W high pressure mercury bulb.
Power supply:	External power supply: Input 100-240Vac 50-60Hz / Output 6Vdc 1A Fluorescence Power Supply 100W. Timer and current display.

Overview

Using the microscope

Adjust the observation head

Loosen the lock-screw, turn the observation head to a comfortable position for observation, and then lock the lock-screw.

Place the specimen on the stage

Lock the specimen slide on the mechanical stage using the slide clamp. Ensure that the specimen is centred over the stage opening by adjusting the coaxial knobs of the stage

Lamp settings

The microscope is fitted with a white LED illuminator. Before turning on the light, read the section about electrical safety precautions. Connect the external power supply into the socket at the base. Turn the brightness adjustment knob to a brightness suitable for observation.

Adjust interpupillary distance

Hold the right and left parts of the observation head by both hands and adjust the interpupillary distance by turning the two parts until one circle of light can be seen.

Focus and diopter adjustment

Turn the dioptic adjustment ring on the left eyepiece to the zero position. Turn the coarse focus knob in order to focus the slide with an objective with low magnification. Adjust the fine focus knob until you obtain a clear and defined picture observing with the right eye, and then act on the left dioptic compensation ring observing with the left eye. When the image appears in focus, choose the necessary objective with the revolving nose-piece. Turn the tension-adjust-knob to get a suitable tension for the focus system.

Condenser adjustment

Turn the condenser adjustment knob to move the condenser up or down for a clear and uniform sample illumination. The two condenser centring screws can be used to centre the condenser.

Set the numerical aperture

Adjust the aperture of the iris diaphragm to set the numerical aperture of the illuminator, thus controlling image contrast and resolution.

Using the fluorescence

Assembling the epi-fluorescence attachment

Extract the attachment from the special packaging and place it on top of the stand of the microscope tightening the locking screw. Then place the optical head above the attachment, setting it with the appropriate locking screw. In the rear part of the attachment slide the lamp housing and connect the cable to the external power supply. Open the lamp housing door and insert the mercury lamp (included in the package) between the two terminals (note that the two poles of the lamp have different sizes so that the insertion in the terminals is unique).

Use and alignment of the lamp

Press the ON/OFF button on the external power supply. Check that the indicator of current after a few seconds shows a value equal to approximately 4A. If the indicator has a value of 0A, press the TRIGGER button. Place on the stage your fluorescence sample and focus. Act on the lamp centering screws until you reach the maximum intensity of emitted fluorescence from the sample. Finally turn the lamp focusing knob in order to achieve uniformity of lighting on the entire field of view.

Warning: for optimal use of the fluorescence wait 10-15 minutes since power-up in order to allow the correct thermal equilibrium of the lamp.

Insertion and replacement of the mercury lamp

If the current drops below 3.5A, replace the lamp. Lifetime is approximately 400 hours.

Take care to cool the lamp for at least 30 minutes before opening its door.

After replacing the lamp, press the reset button to reset the time usage counter.

Unplug any cable from the epi-illuminator, then unscrew the lamp door and open it.

The lamp holder will come out together with the door. Remove the old lamp or the plastic piece (in case this is the first insertion).

It is recommended to use protective glasses because, if not handled with due caution, lamp could crash and injure the operator with glass splinters. Also do not touch the lamp glass during the extraction out of its box, taking care to not exert pressure on any part of it. Plug one lamp end into the flexible metal holder (it fits only in one way), and then insert the other end in the fixed part.

Screw it just enough to keep the lamp in position. If it is too tighten, the lamp could break. Replace the lamp door in its position and connect the cable on the side of the housing to the external power supply.

Troubleshooting

Review the information in the table below to troubleshoot operating problems.

PROBLEM	CAUSE	SOLUTION
LIGHT DOESN'T TURN ON	Power supply not connected	Check that the 6Vdc power supply jack is well inserted on the rear of the microscope
		Rotate the brightness adjustment control and check if there is an increase in the light output
IMAGE CANNOT BE SEEN OR IS DARK.	The iris diaphragm aperture is not completely opened	Completely open the iris diaphragm aperture
	Brightness level is low	Rotate the brightness adjustment potentiometer
	Objective is not aligned with the optical axis	Rotate the nosepiece until an objective is well inserted in the optical path (it "clicks")
IMAGE IS UNCLEAR, BLURRED OR HAS INSUFFICIENT CONTRAST.	Objectives or filters are dirty	Wipe them clean
	The iris diaphragm aperture is not opened correctly	Open the iris diaphragm aperture completely
	Condenser at wrong height	Rotate the condenser knob until you see a uniform illumination

Repair and maintenance

Microscopy environment

This microscope is recommended for use in a clean, dry and shock free environment with a temperature of 0-40°C and a maximum relative humidity of 85 % (non condensing). Use a dehumidifier if needed.

After using the microscope

The microscope should always be kept vertical when moving it so that no moving parts, such as the eyepieces, fall out. Never mishandle or impose unnecessary force on the microscope. Never attempt to service the microscope yourself. Turn down the illumination intensity control and turn the light off. Cover the microscope with the included dust cover, and keep it in a dry and clean place.

Electrical safety precautions

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.

Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

Cleaning the optics

If the optical parts require cleaning first use compressed air.

If that is not sufficient use a soft lint-free cloth with water and a mild detergent.

And as a final option use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.

Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.

Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics. Do not disassemble objectives or eyepieces in an attempt to clean them.

For the best results, use the VWR cleaning kit (see catalogue number below).

If you need to send the microscope to manufacturer for maintenance, please use the packaging if possible.

User replaceable accessories and spare parts

DESCRIPTION	QUANTITY	CAT. NO.
Eyepiece WF10x/20mm.	2	630-2000
Eyepiece WF15x.	2	630-2001
Eyepiece WF20x.	2	630-2002
Micrometer eyepiece WF10x/20mm.	1	630-2003
26x76 mm micrometric slide. Range 1 mm, div. 0,01 mm.	1	630-1650
Objective 4x/0,10 E-PLAN IOS.	1	630-1639
Objective 10x/0,25 E-PLAN IOS.	1	630-1640
Objective 20x/0,40 E-PLAN IOS.	1	630-1641
Objective 40x/0,65 E-PLAN IOS	1	630-1642
Objective PLAN IOS MET 50x/0.75	1	630-2159
Objective 60x/0,80 E-PLAN IOS.	1	630-2049
Objective 100x/1,25 E-PLAN IOS (Oil).	1	630-1643
Complete Phase Contrast Set with IOS PLAN obj. 10x, 20x, 40x, 100x , with Darkfield position	1	630-2059
Polarizing set, filters only	1	630-2285
Rotating table for polarizing set.	1	630-2011
Darkfield condenser for dry objectives.	1	630-2060
Photo adapter for APS-C and Full Frame Reflex cameras	1	630-1645
C-Mount adapter for 1/2" sensor.	1	630-1644
C-Mount adapter for 2/3" sensor.	1	630-2530
Solar battery pack.	1	630-1637
HBO100W high-pressure mercury bulb for fluoresc.	1	630-2158

Technical service

Web Resources

Visit the VWR's website at **www.vwr.com** for:

- Complete technical service contact information
- Access to VWR's Online Catalogue, and information about accessories and related products
- Additional product information and special offers

Contact us For information or technical assistance contact your local VWR representative or visit. **www.vwr.com**.

Warranty


VWR International warrants that this product will be free from defects in material and workmanship for a period of five (5) years from date of delivery. If a defect is present, VWR will, at its option and cost, repair, replace, or refund the purchase price of this product to the customer, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear. If the required maintenance and inspection services are not performed according to the manuals and any local regulations, such warranty turns invalid, except to the extent, the defect of the product is not due to such non-performance.

Items being returned must be insured by the customer against possible damage or loss. This warranty shall be limited to the aforementioned remedies. IT IS EXPRESSLY AGREED THAT THIS WARRANTY WILL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND IN LIEU OF THE WARRANTY OF MERCHANTABILITY.

Compliance with local laws and regulations

The customer is responsible for applying for and obtaining the necessary regulatory approvals or other authorizations necessary to run or use the Product in its local environment. VWR will not be held liable for any related omission or for not obtaining the required approval or authorization, unless any refusal is due to a defect of the product.

Disposal

This equipment is marked with the crossed out wheeled bin symbol to indicate that this equipment must not be disposed of with unsorted waste.

Instead it is your responsibility to correctly dispose of your equipment the end of its life cycle by handling it over to an authorized facility for separate collection and recycling. It is also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect from health hazards the persons involved in the disposal and recycling of the equipment.

For more information about where you can drop off your waste equipment, please contact your local dealer from whom you originally purchased this equipment.

By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

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Ireland / Northern Ireland

VWR International Ltd /
VWR International (Northern Ireland) Ltd
Orion Business Campus
Northwest Business Park
Ballycoolin
Dublin 15
Tel.: 01 88 22 222
Fax: 01 88 22 333
E-mail: sales@ie.vwr.com

Italy

VWR International PBI S.r.l.
Via San Giusto 85
20153 Milano (MI)
Tel.: 02-3320311/02-487791
Fax: 800 152999/02-40090010
E-mail: info@it.vwr.com

The Netherlands

VWR International B.V.
Postbus 8198
1005 AD Amsterdam
Tel.: 020 4808 400
Fax: 020 4808 480
E-mail: info@nl.vwr.com

New Zealand

Global Science - A VWR Company
241 Bush Road
Albany 0632, Auckland
Tel.: 0800 734 100
Fax: 0800 999 002
E-mail: sales@globalscience.co.nz

Norway

VWR International AS
Haavard Martinsens vei 30
0978 Oslo
Tel.: 02290
Fax: 815 00 940
E-mail: info@no.vwr.com

Poland

VWR International Sp. z o.o.
Limbowa 5
80-175 Gdansk
Tel.: 058 32 38 200 do 204
Fax: 058 32 38 205
E-mail: info@pl.vwr.com

Portugal

VWR International -
Material de Laboratório, Lda
Edifício Neopark
Av. Tomás Ribeiro, 43- 3 D
2790-221 Carnaxide
Tel.: 21 3600 770
Fax: 21 3600 798/9
E-mail: info@pt.vwr.com

Singapore

VWR Singapore Pte Ltd
18 Gul Drive
Singapore 629468
Tel.: +65 6505 0760
Fax: +65 6264 3780
E-mail: sales@sg.vwr.com

Spain

VWR International Eurolab S.L.
C/ Tecnología 5-17
A-7 Llinars Park
08450 - Llinars del Vallès
Barcelona
Tel.: 902 222 897
Fax: 902 430 657
E-mail: info@es.vwr.com

Sweden

VWR International AB
Fagerstagatan 18a
163 94 Stockholm
Tel.: 08 621 34 00
Fax: 08 621 34 66
E-mail: kundservice@se.vwr.com

Switzerland

VWR International GmbH
Lerzenstrasse 16/18
8953 Dietikon
Tel.: 044 745 13 13
Fax: 044 745 13 10
E-mail: info@ch.vwr.com

Turkey

Pro Lab Laboratuvar Teknolojileri Ltd.Şti.
a VWR International Company
Orta Mah. Cemal Gürsel Caddesi
Ördekcioglu İşmerkezi No.32/1
34896 Pendik - İstanbul
Tel.: +90216 598 2900
Fax: +90216 598 2907
Email: info@pro-lab.com.tr

UK

VWR International Ltd
Customer Service Centre
Hunter Boulevard - Magna Park
Lutterworth
Leicestershire
LE17 4XN
Tel.: 0800 22 33 44
Fax: 01455 55 85 86
E-mail: uksales@uk.vwr.com

**GO TO VWR.COM FOR THE
LATEST NEWS, SPECIAL OFFERS
AND DETAILS OF YOUR LOCAL
VWR DISTRIBUTOR**

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