Instruction manual

VWR® Advanced Hotplate Stirrer

EU Cat No.

- 444-1502 Ceramic 18 x 18 cm Hotplate Stirrer
- 444-1503 Ceramic 18 x 18 cm Hotplate
- 444-1504 Resin Top 18 x 18 cm Stirrer
- 444-1505 Aluminum 18 x 18 cm Hotplate Stirrer
- 444-1506 Round Top Hotplate Stirrer 135mm
- 444-1507 Ceramic 25 x 25 cm Hotplate Stirrer
- 444-1508 Ceramic 25 x 25 cm Hotplate
- 444-1509 Ceramic 25 x 25 cm Stirrer

Legal address of Manufacturer

Europe

VWR International bv,
Researchpark Haasrode 2020
Geldenaaksebaan 464 B-3001
Leuven
+3216385011
be.vwr.com

UK Importer:
VWR International Ltd
Hunter Boulevard, Magna Park
Lutterworth, Leicestershire, LE17 4XN
uk.vwr.com

United States

VWR International LLC
100 Matsonford Rd
Radnor, PA 19087
+1 800-932-5000
www.vwr.com

Country of Origin

USA

Table of Contents

1. Introduction	
1.1 Safety Information	4
1.2 Intended Use	5
1.3 Package Contents	5
1.4 Installation	5
2. Overview	
2.1 Specifications	6
2.2 Control Panel	8
2.3 Interface	9
2.4 Display	9
3. Operation	
3.1 Getting Ready	10
3.2 Standby Mode	10
3.3 Controlling the Stirrer	10
3.4 Controlling the Top Plate Heater	11
3.5 Using the External Probe	12
4. The Settings Menu	
4.1 Accessing / Exiting	12
4.2 Menu Features	12
4.3 Using the Single Point Calibration (SPC) Feature	13
4.3.1 Start a Single Point Calibration	13
4.3.2 Adjust a Calibrated Temperature	13
4.3.3 Clear a Calibrated Temperature	14
4.4 Enabling / Disabling the Beeper	14
4.5 Changing the Power Recovery Setting	14
4.6 Reset to Factory Default Settings	15
4.7 Software Version	15
5. Maintenance	
5.1 Cleaning	16
5.2 Replacing Power Fuse	16
5.3 Troubleshooting	17
5.4 Service Information	18
6. Technical Data	18
7. Compliance	18
Technical Service	18
Warranty	19
Equipment Disposal	19

1. INTRODUCTION

This manual contains installation, operation and maintenance instructions for the VWR Advanced Hotplate-Stirrer. Please read the manual completely before using.

1.1. Safety Information

Definition of Signal Warnings and Symbols

Safety notes are marked with signal words and warning symbols. These show safety issues and warnings. Ignoring the safety notes may lead to personal injury, damage to the instrument, malfunctions and false results.

WARNING	For a hazardous situation with medium risk, possibly resulting in severe injuries or death if not avoided.
CAUTION	For a hazardous situation with low risk, resulting in damage to the device or the property or in loss of data, or minor or medium injuries if not avoided.
ATTENTION	For important information about the product. May lead to equipment damage if not avoided.
NOTE	For useful information about the product.

Warning Symbols

General hazard

Explosion hazard

Caution, hot surface

Protective conductor terminal

Alternating current

Electrical shock hazard

Safety Precautions

WARNING! The protection provided by the unit may be impaired if used with accessories not provided or recommended by the manufacturer, or used in a manner not specified by the manufacturer.

- Always operate unit on a level surface for best performance and maximum safety.
- DO NOT lift unit by the top plate.
- To avoid electrical shock, completely cut off power to the unit by disconnecting the power cord from the unit or unplugging from the wall outlet.
- Disconnect unit from the power supply prior to maintenance and servicing.
- Spills should be removed promptly, after the unit has cooled down.
- Alkaline spills, hydrofluoric acid or phosphoric acid spills may damage the unit and lead to thermal failure.
- DO NOT immerse the unit for cleaning.
- DO NOT operate the unit at high temperatures without a vessel/sample on the top plate.
- DO NOT operate the unit if it shows signs of electrical or mechanical damage.
- Protective earthing of the equipment is achieved via connection of the provided power cord to a compatible grounded power outlet.

WARNING: unit is not explosion proof. Use caution when unit is on or when heating volatile materials.

WARNING! DO NOT use the unit in explosive atmospheres or with materials that could cause a hazardous environment from processing. Keep in mind the material flash point relative to the target temperature that has been set.

CAUTION! The top plate can reach 500°C, DO NOT touch the heated surface. Use caution at all times. Keep the unit away from explosive vapors and clear of papers, drapery and other flammable materials. Keep the power cord away from the heater plate.

CAUTION! Beware of the following risks when heating.

- Flammable materials
- Low boiling point combustible substances
- Glass breakage as a result of mechanical shaking power
- Incorrect container size
- Too much medium
- Unsafe condition of container

Earth Ground - Protective Conductor Terminal

Alternating Current

1.2. Intended Use

This instrument is intended for use in laboratories, pharmacies, schools, businesses and light industry. It must only be used for processing materials as described in these operating instructions. Any other type of use and operation beyond the limits of technical specifications, without written consent from VWR, is considered as not intended. This instrument complies with current industry standards and the recognized safety regulations; however, it can constitute a hazard in use. If the instrument is not used according to these operating instructions, the intended protection provided by the instrument may be impaired.

1.3. Package Contents

- Hotplate-Stirrer / Hotplate / Stirrer Unit
- Power cord
- Stir Bar (Stirring Models ONLY)

1.4. Installation

Upon receiving the VWR Hotplate-Stirrer, check to ensure that no damage has occurred during shipment. It is important that any damage that occurred in transport is detected at the time of unpacking. If you do find such damage, the carrier must be notified immediately.

After unpacking, place the Hotplate-Stirrer on a level bench or table, away from explosive vapors. Ensure that the surface on which the unit is placed will withstand typical heat produced by the unit and place the unit a minimum of six (6) inches from vertical surfaces. Do not position the equipment such that it is difficult to disconnect the power cord during use. Always place the unit on a sturdy work surface.

The Hotplate-Stirrer is supplied with a 3 conductor, grounded power cord that should be plugged into a matching standard grounded outlet. If the cord supplied does not meet your needs, please use an approved power cord that has ratings equal or exceeding those of the originally provided cord and that complies with the local/national regulations of the country in which the equipment is to be used. Replacement of the plug must be made by a qualified electrician.

2. Product Overview

2.1. Specifications

Model	Advanced 18 x 18 cm Ceramic Hotplate Stirrer	Advanced 18 x 18 cm Aluminum Hotplate Stirrer	Advanced 18 x 18 cm Ceramic Hotplate	Advanced 18 x 18 cm Resin Stirrer
Function	Heating & Stirring		Heating	Stirring
Top Plate Size	18 x 18 cm			
Top Plate Material	Ceramic	Aluminum	Ceramic	Resin
Control	Digital			
Display	75 x 20 mm White Character LCD			
Temperature Range	Ambient +5 °C - 500 °C	Ambient +5 °C - 340 °C	Ambient +5 °C - 500 °C	N/A
Temperature Resolution	1 °C			N/A
Temperature Readability	1 °C			N/A
Temperature Stability of Top Plate*	$\pm 3 \%$ (> 100 °C), $\pm 2 \%$ ($\leq 100 \text{ °C}$)			N/A
Temperature Stability of Temperature Probe**	$\pm 2 \%$ (> 100 °C), $\pm 2 \%$ ($\leq 100 \text{ °C}$)			N/A
SPC Point	1 point for plate temp 1 point for probe temperature			N/A
Probe Connection	PT1000 (Class A)			N/A
Wattage	1050W			35W
Speed Range	80-1600 RPM		N/A	80-1600 RPM
Speed Resolution	10 RPM		N/A	10 RPM
Speed Stability**	+/-2%		N/A	+/-2%
Stirring Capacity	15L H2O		N/A	15L H2O
Weight Capacity	23kg			
IP Rating	IP21			
Housing Material	Painted Aluminum			
Dimensions (L x W x H)	31.8 x 22.3 x 11.2 cm			
Net Weight	3.9 kg	4.0 kg	3.3 kg	3.0 kg
Electrical (50/60Hz)	220-240 V, 7 A			220-240 V, 1 A

Note:

- *The parameter in the specification table is applicable to (5 cm) diameter center of the top plate.
- ** Conditions permitting. Variations in temperature and speed measurement processes, vessels, ambient and sample will impact actual performance. To improve temperature accuracy of the system, use the Single Point Calibration feature.


Model	Advanced 25 x 25 cm Ceramic Hotplate Stirrer	Advanced 25 x 25 cm Ceramic Hotplate	Advanced 25 x 25 cm Ceramic Stirrer	Advanced 135mm Round Top Hotplate Stirrer
Function	Heating & Stirring	Heating	Stirring	Heating & Stirring
Top Plate Size	25 x 25 cm			Ø 13.5 cm
Top Plate Material	Ceramic			Ceramic Coated Stainless Steel
Control	Digital			
Display	115x30mm White Character LCD			75 x 20 mm White Character LCD
Temperature Range	Ambient +5°C – 500°C		N/A	Ambient +5 °C - 380 °C
Temperature Resolution	1°C		N/A	1 °C
Temperature Readability	1°C		N/A	1 °C
Temperature Stability of Top Plate*	±3% (>100°C), ±2°C (≤ 100°C)		N/A	± 3 % (> 100 °C), ± 2 °C (≤ 100 °C)
Temperature Stability of Temperature Probe**	±2% (>100°C), ±2°C (≤100°C)			± 2 % (> 100 °C), ± 2 °C (≤ 100 °C)
SPC Point	1 point for plate temp 1 point for probe temp		N/A	1 point for plate temp 1 point for probe temp
Probe Connection	PT1000 (Class A)		N/A	PT1000 (Class A)
Wattage	2066W		35W	550 W
Speed Range	80-1600 rpm	N/A	80-1600 rpm	
Speed Resolution	10 RPM	N/A	10 RPM	
Speed Stability**	+/-2%	N/A	+/-2%	
Stirring Capacity	22L H2O	N/A	22L H2O	15L H2O
Weight Capacity	23kg			25kg
IP Rating	IP21			
Housing Material	Painted Aluminum			
Dimensions (L x W x H)	40.2 x 28.7 x 11.8 cm			26.4 x 16.3 x 10.1 cm
Net Weight	6.7kg	6.1kg	6.7kg	2.5kg
Electrical (50/60Hz)	220-240V, 10A		220-240V, 1A	220-240 V, 4 A

Note:

*The parameter in the specification table is applicable to (5 cm) diameter center of the top plate.

** Conditions permitting. Variations in temperature and speed measurement processes, vessels, ambient and sample will impact actual performance. To improve temperature accuracy of the system, use the Single Point Calibration feature.

2.2. Control Panel

A	Left Knob	Control temperature and Menu
B	Right Knob	Control speed (and Menu for Stirrer models)
C	Display Screen	Show operation status, parameters, and menu settings
D*	Hot Top Caution Indicator	Hot surface warning (not applicable for Stir only models)
E*	Heater Indicator	Illuminates when heater is working
F*	Speed Indicator	Illuminates when stirrer is working


* The feature D, E, and F are integrated in the LCD display on the Advanced 25 x 25 cm Hotplate Stirrer, Hotplate, and Stirrer.

2.3. Interface

#	Description
A	External Temperature Probe Port
B	Power Entry Module
C	Fuse Drawer
D	Rod Holder

2.4. Display

Heating Display	
A.	Heater Temperature: Switches to external probe temperature when the probe is plugged in and external probe icon is illuminated
B.	Temperature Setting Icon: Switches heater temperature to heat setting when illuminated.
C.	Single Point Calibration Icon
D.	External Probe Icon: Illuminates when the external probe is plugged in.
Stirring Display	
E.	Speed Setting Icon: Illuminates until stirrer reaches the speed setting.
F.	Stir Speed

3. Operation

3.1. Getting Ready

To get ready:

1. Plug the female end of the provided power cord into **Power Entry Module** (PEM) on the rear side of the unit.
2. Plug the male end of the power cord into a matching standard grounded outlet.
3. The unit will beep once, and the screen will illuminate with the following displays:
 - a) The first will display all the icons on the screen. In addition, the **Heater Indicator**, the **Speed Indicator** and the **Hot Top Caution Indicator** on the control panel will light up.
 - b) The second will display the unit type (left) and the software version (right).
 - c) The third will display the unit's electrical voltage (left) and frequency (right).
 - d) The fourth will display the unit's main operating screen.

Note: If the fourth screen displays OFF, the unit is in standby mode.

3.2. Standby Mode

The standby switch on the front side of the unit controls standby mode.

1. When the unit is switched off:
 - a) All heating, stirring, and timing functions will turn off.
 - b) The screen will display OFF.

In addition, If the heater temperature is above 40°C, the on the control panel will remain illuminated as well as the word **HOT** on the screen.

2. When the unit is switched on:
 - a) All heating and stirring functions will remain off.
 - b) The main operating screen will return to display the previous heating and stirring settings.
 - c) The unit is ready for normal use.

3.3. Controlling the Stirrer

1. Rotate the right knob to control the speed setting.
 - a) Clockwise rotation will increase the speed setting.
 - b) Counterclockwise rotation will decrease the speed setting.
2. To turn on the stirrer, press and hold the right knob until the unit beeps and the on the control panel illuminates.
 - a) The unit will beep once to confirm the stirrer has been turned on.
 - b) The indicator light will blink to indicate that the stirrer is on and ramping to the target speed.
 - c) Once the stirrer has reached the target speed, the **Speed Setting Icon** will disappear from the screen, and the **Speed Indicator** will stop blinking and remain illuminated.

Note: the unit will only display the target speed during the process
3. To change the speed setting while the stirrer is on:
 - a) Rotate the right knob clockwise or counterclockwise to the new speed setting.
The speed number on the screen will blink to indicate that the speed setting is not confirmed.
 - b) Briefly press the right knob to confirm the new speed setting.
The speed number on the screen will stop blinking once the new setting is confirmed.

Note: If the speed setting screen remains idle without confirmation for 6 seconds, it will reset to the current setting.
4. To turn off the stirrer, press and hold the right knob until the unit beeps and the lights off.
 - a) The unit will beep once to confirm the stirrer has been turned off.
 - b) The on the control panel will go off.
 - c) The on the screen will illuminate.

STIRRING OPERATING TIPS

The stirrer increases speed at a steady rate until the setpoint is reached. If the stirrer is not reaching its setpoint: 1) the stir bar may be too large, 2) the liquid may be too viscous, 3) the setpoint speed may need to be reduced. Additionally, the magnetic strength of stir bars reduce over time and may need to be replaced. Apart from that, the bottom surface of the vessel can affect the stirring performance of the magnetic stir bar. Adjust the vessel being used or change the stir bar to optimize mixing in the sample. When heating and stirring a reaction vessel within an oil bath

or similar set-up, the stirring function will stir up to approximately one inch (2.5 cm) from the top plate. The stirring speed will vary according to liquid viscosity, spin bar length, and distance from top plate. Adjust one or all of these to achieve the desired stirring speed. For example: the closer the reaction vessel is to the top plate, the stronger the magnetic connection between the unit and the stir bar.

3.4. Controlling the Top Plate Heater

1. Rotate the left knob to control the heat setting.
 - a) Clockwise rotation will increase the heat setting.
 - b) Counterclockwise rotation will decrease the heat setting.
2. To turn on the heater, press and hold the left knob until the unit beeps and the illuminates.
 - a) The unit will beep once to confirm the heater has been turned on.
 - b) The on the control panel will illuminate to indicate that the heater is running.
 - c) During heating, the current heater temperature and the set target temperature along with the **Temperature Setting Icon** will be alternately displayed on the screen, while the **Heater Indicator** will blink.
 - d) When the temperature is reached stably, the screen will display the current heater temperature, and the **Heater Indicator** will stay illuminated.
 - e) When the heater temperature is above 40°C, the **Hot Top Caution Indicator** on the control panel will be illuminated.
3. To change the heat setting while the heater is on:
 - a) Rotate the left knob clockwise or counterclockwise to the new heat setting. The heating temperature on the screen will blink to indicate that the heat setting is not confirmed when heating is active. The **Temperature Setting Icon** will stay illuminated on the screen.
 - b) Briefly press the left knob to confirm the new heat setting. The heating temperature on the screen will stop blinking once the new setting is confirmed. If the setting temperature remains idle without confirmation for 6 seconds, it will reset to the current setting.
4. To turn off the heater, press and hold the left knob until the unit beeps and the **Heater Indicator** lights off.
 - a) The unit will beep once to confirm the heater has been turned off.
 - b) The **Heater Indicator** on the control panel will go off.

CAUTION: THIS DOES NOT MEAN THAT THE TOP PLATE IS SAFE TO TOUCH.

 - c) Once the heater's temperature cools below 40°C, the on the control panel will disappear.

HEATING OPERATING TIPS

Overshoot:

The unit may overshoot the temperature before stabilizing at the setpoint. The two methods to minimize overshoot are:

1. Metal containers minimize overshoot.

CAUTION! When heating metal containers on a ceramic top plate, it is recommended to use the lowest temperature setting possible to limit thermal stress to the ceramic top plate.
2. If a glass vessel is used, anticipate overshoot. Start with a temperature setpoint 5 to 10°C below the desired temperature. When the temperature stabilizes at this lower setting, increase the heater to the final temperature.

The temperature display on the unit represents the estimated top plate temperature, not the sample temperature. When external probe is in use, the temperature display on the unit represents the sample temperature. The vessel contents being heated may be at a lower temperature depending on the size and thermal conductivity of the vessel. It may be beneficial to monitor the temperature of the vessel contents and adjust the setpoint temperature accordingly. If you need precise control, use the VWR External Temperature Probe.

Typical Time to Boil Water

The chart below is an example of an approximate time to boil for the specified amount of water in a specific vessel. These values are only approximate and can vary from unit to unit. Values are based on 23°C water in an ambient environment of 23°C.

Unit Size	Heater Temp. Limit	Volume of Water	Typical Time to Boil
Round Top	380°C	500mL in 1L beaker	≈ 25 min
18 x 18 cm Ceramic	500°C	500mL in 1L beaker	≈ 18 min
18 x 18 cm Aluminum	340°C	500mL in 1L beaker	≈ 16 min
25 x 25 cm	500°C	500mL in 1L beaker	≈ 17 min

3.5. Using the External Probe

1. Connect the VWR External Temperature Probe to the External RTD Probe Port on the rear panel of the unit. Once the VWR External Temperature Probe is connected, the **External Probe Icon** on the screen will illuminate.
2. The display will now show the temperature of the external probe instead of the heater. The **Hot Top Caution Indicator** light will illuminate once the heater temperature reaches 40°C.
Note: When using the VWR External Temperature Probe, the temperature setpoint should be adjusted to the desired sample temperature. If the temperature setpoint is higher than the sample can achieve, an E7 will occur (see the Troubleshooting section). Reduce sample volume or temperature setpoint value. For Example: Water has a theoretical temperature limit of 100°C (boiling). A temperature setpoint greater than 100°C will cause an E7 error.
3. If the external temperature probe is inserted into the external RTD probe port while the heater is running:
 - a) The heater will shut off.
 - b) The unit will display an E7 error.
 - c) The unit will beep 10 times.
4. If the external temperature probe is removed from the external RTD probe port while the heater is running:
 - a) The heater will shut off.
 - b) The unit will display an E4 error.
 - c) The unit will beep 10 times.

Note: To clear an E4 or E7 error, flip the standby switch off and back on. The unit will be ready for normal use.

4. THE SETTINGS MENU

4.1. Accessing / Exiting

1. To access the settings menu, press and hold the left knob until **⌂ U** appears on the screen.
 - a) Continue to hold the left knob after the unit beeps and the **Heater Indicator** on the control panel illuminates.
 - b) The heater will not turn on unless the left knob is released before **⌂ U** appears.
 - c) The settings menu cannot be accessed while the heater or the stirrer is running.
 - d) **⌂ U** will appear briefly then proceed to the top level of the settings menu.
2. Rotate the left knob to navigate the different menu options and briefly press the left knob to select / enter / edit the displayed setting.
3. To exit the menu from the top level, rotate the left knob clockwise until **ESC** is displayed and briefly press the left knob. The unit will return to the main operating screen.
Note: To exit the menu at any time, flip the standby switch off and back on. The unit will be ready for normal use. Turning off the unit will not reset / change the settings.

4.2. Menu Features

The top level of the settings menu has the following features:

1. **SPC** – Single Point Calibration
Single Point Calibration (SPC) improves the accuracy of the heater at user selected temperature points. Up to 1 point (Plate Controlled Heating Temperature) and 1 point (Probe Controlled Heating Temperature) can be stored.
2. **SYS** – System Settings
 System Settings allows the user to set additional features, such as enabling /disabling the **Beeper**, changing the **Power Recovery Setting**, and resetting to **Factory Default Settings**.
3. **ESC** – EXIT
 The unit will save the current settings and return to the main operating screen

4.3. Using the Single Point Calibration (SPC) Feature

4.3.1. Start a Single Point Calibration (SPC)

1. To enter the **Single Point Calibration**, the unit must first be at the top level of the settings menu. Please refer to **Accessing / Exiting** Section for how to enter the top-level menu. Then wait for the screen to show **SPC**.
2. Briefly press the left knob to enter the **SPC** settings menu.
3. Briefly press the left knob again to change the calibration temperature. The temperature will begin to blink to indicate that it can be modified.
4. Rotate the left knob to scroll to the desired temperature.
5. Press and hold the left knob until the unit beeps and the illuminates to begin calibration at the set temperature.
 - a) The unit will begin to heat to the set temperature.
 - b) The **Single Point Calibration Icon** on the display will blink to indicate that the calibration is running.
 - c) The left and right knobs will be disabled until the calibration is completed.
 - d) If the external probe is connected, the stirrer will turn on simultaneously.
 - e) To cancel the calibration while it is running, flip the standby switch off to turn off the unit.
6. Once the unit has reached the calibration temperature, the **Single Point Calibration Icon** and the calibration temperature will blink.
7. With a secondary temperature measurement device, measure the temperature of the top plate, or the heated sample at the location of the external probe if using probe control.
8. Rotate the left knob to scroll to the measured temperature from the secondary temperature measurement device, and then briefly press the left knob to select the new temperature.
9. The unit will begin to regulate temperature with compensated error. When this is complete, **SAVE** will appear and blink on the screen.
10. Measure the temperature at the same location as step 7.
11. Rotate the left knob to select from:
 - a) **SAVE** – to retain calibration (stirring will stop if using probe control).
 - b) **Adj** – to prompt fine tuning of calibration (return to step 8).
 - c) **ABCE** – to cancel the single point calibration and return to the initial calibration menu.
12. If you select **SAVE**, the screen will return to the top level of the settings menu.
13. Rotate the left knob to scroll to **ESC** (EXIT)
14. Briefly press the left knob to return to the main operating screen.
15. If you select **ABCE** in step 11, the screen will return to the top level of the settings menu. Rotate the left knob to **ESC** (EXIT) and briefly press the left knob to return to the main operating screen.

4.3.2. Adjust a Calibrated Temperature

1. To adjust a calibrated temperature, the unit must first be in the top level of the settings menu. Please refer to **Accessing / Exiting** Section for how to enter the top level menu. Then wait for the screen to show **SPC**.
2. Briefly press the left knob to enter the **SPC** settings menu.
3. Rotate the left knob to select which calibration temperature to adjust if two calibration records are stored. The unit can store 1 point (Plate Controlled Heating Temperature) and 1 point (Probe Controlled Heating

Temperature) calibration record. For calibration performed with an external probe, the screen will display **External Probe Icon**.

4. Press and hold the left knob until the unit beeps and the **Heater Indicator** illuminates to begin calibration. Then repeat steps in the previous **Start a Single Point Calibration** Section from step 5 to the end to finish the process.

4.3.3. Clear a Calibrated Temperature

1. To clear a calibrated temperature, the unit must first be at the top level of the settings menu. Please refer to **Accessing / Exiting** Section for how to enter the top-level menu. Then wait for the screen to show **SPC**.
2. Briefly press the left knob to enter the **SPC** settings menu.
3. Rotate the left knob to scroll to **CLR**.
4. Briefly press the left knob to enter the menu. If there is no stored SPC point, select **BACK** to return to the previous screen.
5. Rotate left knob to scroll to the calibration point you want to clear. The calibration points are stored in ascending order by temperature.
6. Briefly press the left knob to confirm your selection.
7. Rotate the left knob to select from:
 - a) **YES** – to clear the selected calibration temperature.
 - b) **NO** – to return to **SPC** settings menu.
8. Briefly press the left knob to confirm your selection and return to the Calibration settings menu.

4.3.4. Use a Calibrated Temperature

1. To heat to a calibrated temperature, the unit should return to the main operating screen.
2. Rotate the left knob to scroll the heat setting to the stored calibration temperature. The **SPC** icon will appear once the temperature is correct.
3. Press and hold the left knob until the unit beeps and the illuminates. The unit will heat to the stored calibration temperature.

4.4. Enabling / Disabling the Beeper

Disabling the beeper will prevent beeps in the following scenarios:

- Starting and stopping the heater
- Starting and stopping the stirrer
- When the heater reaches the set temperature
- Starting Single Point Calibration **SPC**

1. To enter the beeper setting, the unit must first be at the top level of the settings menu. Please refer to **Accessing / Exiting** Section for how to enter the top-level menu.
2. Rotate the left knob to scroll to the **SYS** (System) setting.
3. Briefly press the left knob to enter the System settings menu.
4. The display shows **BEEP**. Briefly press the left knob to enter the beeper setting.
5. Rotate the left knob to scroll to the desired beeper setting. **EN** for enable and **DIS** for disable.
6. Briefly press the left knob to confirm your selection.
7. Rotate the left knob to **BACK**.
8. Briefly press the left knob to return to the top level of the settings menu.
9. Rotate the left knob to scroll to **ESC** (EXIT).
10. Briefly press the left knob to return to the main operating screen.

Note: There is not an icon to indicate that the beeper has been disabled.

4.5. Changing the Power Recovery Setting

Power Recovery is an optional feature that allows the unit to automatically restart heater and stirrer functions when power is returned to the unit after a disconnect. By default, this feature is turned off.

1. To enter the power recovery setting, the unit must first be at the top level of the settings menu. Please refer to **Accessing / Exiting** Section for how to enter the top-level menu.

2. Rotate the left knob to scroll to the **SYS** (System) setting.
 3. Briefly press the left knob to enter the **System** settings menu.
 4. Rotate the left knob to **PLUR** (Power Recovery).
 5. Briefly press the left knob to enter the **Power Recovery** settings.
 6. Rotate the left knob to scroll to the desired Power Recovery setting.
OFF – heating and stirring functions will need to be manually restarted after power restoration.
ON – heating and stirring functions will automatically restart upon power restoration.
 7. Briefly press the left knob to confirm your setting.
 8. Rotate the left knob to **BACK**.
 9. Briefly press the left knob to return to the top level of the settings menu.
 10. Rotate the left knob to scroll to **ESC** (EXIT).
 11. Briefly press the left knob to return to the main operating screen.
- Note:** There is not an icon on the display to indicate that the **Power Recovery** has been activated.

4.6. Reset to Factory Default Settings

Resetting the unit to Factory Default Settings will do the following:

- Clear all Single Point Calibration (SPC) temperatures.
 - Turn off Power Recovery.
 - Re-enable the beeper setting
1. To reset the unit to factory default settings, the unit must first be at the top level of the settings menu. Please refer to **Accessing / Exiting Section** for how to enter the top-level menu.
 2. Rotate the left knob to scroll to the **SYS** (System) setting.
 3. Briefly press the left knob to enter the **System** settings menu.
 4. Rotate the left knob to **RESET** (Reset).
 5. Briefly press the left knob to enter the menu.
 6. Rotate the left knob to scroll to the desired setting. **YES** for resetting and **NO** for not resetting.
 7. Press and hold the left knob until the unit beeps to confirm the **Reset** setting.
 8. Rotate the left knob to **BACK**.
 9. Briefly press the left knob to return to the top level of the settings menu.
 10. Rotate the left knob to scroll to **ESC** (EXIT).
 11. Briefly press the left knob to return to the main operating screen.

Note: There is not an icon on the display to indicate that the unit has been reset to factory default settings.

4.7. Software Version

To check the unit's software version:

1. The unit must first be in the top level of the settings menu. Please refer to **Accessing / Exiting Section** for how to enter the top-level menu.
2. Rotate the left knob to scroll to the **SYS** (System) setting.
3. Briefly press the left knob to enter the settings.
4. Rotate the left knob clockwise until you see numbers, such as **1.04**. That is the software version for your unit
5. To return to the operating screen, rotate the left knob to **BACK**.
6. Briefly press the left knob to return to the top level of the settings menu.
7. Rotate the left knob to scroll to **ESC** (EXIT).
8. Briefly press the left knob to return to the main operating screen.

5. MAINTENANCE

5.1. Cleaning

WARNING: Electric Shock Hazard. Disconnect the equipment from the power supply before cleaning. Make sure that no liquid enters the interior of the instrument.

Attention: Do not use solvents, harsh chemicals, ammonia or abrasive cleaning agents.

The housing may be cleaned with a cloth dampened with a mild detergent if necessary.

The Hotplate-Stirrer / Hotplate / Stirrer is built for long, trouble-free, dependable service. No lubrication or other technical user maintenance is required. It needs no user maintenance beyond keeping the surface clean. The unit should be given the care normally required for any electrical appliance. Avoid wetting or unnecessary exposure to fumes. Spills should be removed promptly after the unit has cooled down. Before using any cleaning or decontamination method except as noted in this section, users should check with the manufacturer that the proposed method will not damage the equipment. Do not use a cleaning agent or solvent on the front panel which is abrasive or harmful to glass, nor one which is flammable. Always ensure the power is disconnected from the unit prior to any cleaning. If the unit ever requires service, contact your VWR representative. The user is responsible for carrying out appropriate decontamination if hazardous material is spilled onto or into the equipment.

CLEANING CERAMIC TOPS:

First remove any burnt-on deposits or spills from the top plate with a scraper (similar to scraping paint off of windowpanes). For your safety, please wear an insulated mitt when using a metal scraper. When the top plate has cooled, apply a few dabs of a non-abrasive cleaner over the surface with a damp paper towel. As a final step, clean with water, and wipe surface with a clean, dry paper towel.


CLEANING ALUMINUM TOPS:

For simple dust and dirt, clean the aluminum top by using a damp cloth with soap and water. For more stubborn deposits, try using a flat edge wooden spatula to scrape off as much as possible. For more stubborn stains, try using a couple of tablespoons of white vinegar to two pints of water and mix well. Dip a clean cloth into the mixture and gently rub the exterior of the aluminum surface. Generally, it is not a good idea to use abrasive pads or cleaners on aluminum, as the metal will scratch easily. If you must use some type of abrasive, try applying baking soda to the surface and then rubbing with a moist cloth. This will work as well as most scouring pads and is less like to create deep scratches in the surface. Be careful not to use steel wool or scouring pads as they can leave the aluminum riddled with little scratches that make it harder to clean in the future. If you feel you must use steel wool, use the finest grade you can find and use as sparingly as possible with as little pressure as possible. Go with the grain rather than using circular motions.

5.2. Replacing Power Fuse


The use of a fuse of a different type or with a different value or bridging or shunting the fuse is not allowed and can possibly cause a hazard to your safety and lead to instrument damage!


WARNING: Electric Shock Hazard. Disconnect the equipment from the power supply before replacing the fuse

To replace the fuse:

1. Put the tip of a screwdriver or a test pen in the circled position below and then pull the fuse holder out.


2. Take out the old fuse and then push the new fuse into the slot. Choose the correct fuse to replace according to the fuse specification listed in the **Dimensions** Section.

Note: the unit comes with a spare fuse in the fuse holder (in the circled position above).

Attention: If the fuse is good and power is available at the outlet, the cord or unit may be defective. Try a new cord. If this does not work, the unit should be sent back for servicing.

3. Push the fuse holder back in to finish the replacement.

5.3. Troubleshooting

The following table lists common problems and possible causes and remedies. If the problem persists, contact VWR or your authorized dealer.

Table 5-3 Troubleshooting

Error*	Cause of Error	How to Fix
Unit fails to power on	Missing or blown fuse	Add or replace fuse as necessary.
E1	Plate RTD open	Not fixable by user, please contact VWR.
E2	Plate RTD short	Not fixable by user, please contact VWR.
E3	No stirring motion / cannot reach speed	Not fixable by user, please contact VWR.
E4	Probe RTD open (Removing the probe while the unit is heating)	Switch unit to standby, then return to normal operating mode.
E5	Probe RTD short (Malfunctioning probe)	Switch unit to standby, remove the probe from the unit, then return to normal operating mode.
E6	A/D lock error	Not fixable by user, please contact VWR.
E7	User Probe Error (Plugging the probe into the unit while it is heating)	Switch unit to standby, then return to normal operating mode.
E8	Plate over temperature	Not fixable by user, please contact VWR.
E9	Plate under temperature	Not fixable by user, please contact VWR.
E10	Triac fault	Not fixable by user, please contact VWR.
E11	Heater damage	Not fixable by user, please contact VWR.
E12	Heater overheating	Switch unit to standby, then return to normal operating mode.
AC Err	Mains frequency out of the range 40~55Hz (rating 50Hz) or 55~70Hz (rating 60Hz)	Regulate mains frequency within the range.

***Note:** Error code instances will stop equipment operation by default.

5.4. Service Information

If the troubleshooting section does not resolve or describe your problem, contact your authorized VWR service agent. For service assistance or technical support in the United States call toll-free 1-888-897-5463 between 8:00 AM and 5:30PM EST. A VWR product service specialist will be available to provide assistance. Outside the USA, please visit our website, <https://vwr.com> to locate the VWR office nearest you.

Serial Number: _____

Date of Purchase: _____

Supplier: _____

6. Technical Data

Equipment Ratings:

Indoor use only

Altitude: 0 to 2000m

Operating temperature: 5°C to 40°C

Storage temperature: -20°C to 65°C

Operating humidity: 20 to 80% relative humidity, non-condensing

Storage humidity: 20 to 80% relative humidity, non-condensing

Voltage fluctuations: Mains supply voltage fluctuations up to $\pm 10\%$ of the nominal voltage.

Overvoltage category

(Installation category): II

Pollution degree: 2

7. COMPLIANCE

Compliance to the following standards is indicated by the corresponding mark on the product.

Mark	Standard
	This product complies with the applicable harmonized standards of EU Directives 2011/65/EU (RoHS), (EU)2015/863 (RoHS), 2014/30/EU (EMC) and 2014/35/EU (LVD). The EU Declaration of Conformity is available online at https://cedclaration.vwr.com/
	This product complies with the EU Directive 2012/19/EU (WEEE). Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. For disposal instructions in Europe, refer to www.VWR.com/
	EN 61326-1, EN 61010-1, EN IEC 61010-2-010, EN IEC 61010-2-051, EN IEC 63000

Global Notice

Warning: This is a Class B product. In a domestic environment this product may cause radio interference in which case the user may be required to take adequate measures.

ISO 9001 Registration

The management system governing the production of this product is ISO 9001 certified.

Technical service

Web Resources

Visit the VWR website at www.vwr.com for:

- Complete technical service contact information

- Access to the VWR Online Catalogue, and information about accessories and related products
- Additional product information and special offers

Contact us For information or technical assistance contact your local VWR representative or visit: www.vwr.com.

Warranty


VWR warrants that this product will be free from defects in material and workmanship for a period of five (5) years from date of delivery. If a defect is present, VWR will, at its option and cost, repair, replace, or refund the purchase price of this product to the customer, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear. If the required maintenance and inspection services are not performed according to the manuals and any local regulations, such warranty turns invalid, except to the extent, the defect of the product is not due to such non-performance.

Items being returned must be insured by the customer against possible damage or loss. This warranty shall be limited to the aforementioned remedies. IT IS EXPRESSLY AGREED THAT THIS WARRANTY WILL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND IN LIEU OF THE WARRANTY OF MERCHANTABILITY.

Compliance with local laws and regulations

The customer is responsible for applying for and obtaining the necessary regulatory approvals or other authorizations necessary to run or use the Product in its local environment. VWR will not be held liable for any related omission or for not obtaining the required approval or authorization, unless any refusal is due to a defect of the product.

Equipment disposal

This equipment is marked with the crossed out wheeled bin symbol to indicate that this equipment must not be disposed of with unsorted waste.

Instead, it's your responsibility to correctly dispose of your equipment at lifecycle -end by handling it over to an authorized facility for separate collection and recycling. It's also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect from health hazards the persons involved in the disposal and recycling of the equipment.

For more information about where you can drop off your waste of equipment, please contact your local dealer from whom you originally purchased this equipment.

By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

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VWR International AB
VWR International GmbH
VWR International, LLC