



393PR-01

G-Biosciences ♦ 1-800-628-7730 ♦ 1-314-991-6034 ♦ technical@GBiosciences.com

A Geno Technology, Inc. (USA) brand name

Protease Screening™

For Testing for Proteases in Protein Samples

(Cat. # 786-137)



think proteins! think G-Biosciences www.GBiosciences.com

INTRODUCTION

Sensitive protease screening is needed to study protease activity present in the sample of interest. Protease screening kit is a simple and quick method for testing the presence of proteases in the protein samples. The screening kit uses dye- labeled protein substrate. The proteases present in the solution will digest the protein substrate and release dye labeled peptides. The absorbance of dye-labeled peptide is measured at 570nm for determination of protease activity. The kit is sufficient for 50 assays in a microtiter format.

ITEM(S) SUPPLIED (Cat. # 786-137)

Description	Size
Incubation Buffer	5ml
Precipitation Agent	5ml
Assay Buffer	6ml
Protease Substrate (lyophilized)	1vial (150µl)

STORAGE CONDITIONS

The kit is shipped at ambient temperature. Store Protease Screening kit at 4°C. When stored and used properly the kit is stable for 1 year.

ADDITIONAL ITEM(S) REQUIRED

96 well titer plates, microfuge tubes, plate reader, spectrophotometer, centrifuge, and micro-plate reader.

PREPARATION BEFORE USE

Substrate

Reconstitute the supplied Protease Substrate (dye-labeled protein) by adding 150µl water into the vial, mix it to dissolve completely. After reconstitution, store the substrate at -20°C and is stable for 3-4 months.

PROTOCOL

NOTE: The reaction set up is designed for 96 well titer plates. The use of 96 well titer plates requires a centrifuge adapted for 96 well titer plates. Alternatively, the assay may be performed in micro-centrifuge tubes and the final reaction product is transferred to 96 well titer plates for the measurement of optical density.

1. Set up the reaction as indicated in the table below. We recommend performing each assay in triplicate:


Blank	Blank	Test Sample
Test Sample	-	1-45µl
Protease Substrate Solution	2.5µl	2.5µl
Incubation Buffer or buffer of your choice	to 50µl	to 50µl

2. Seal the plate or tube and incubate at 37°C for 2-3 hours
NOTE: For slow acting proteases, incubation time may be extended up to 24 hours.
3. After incubation, add 50µl Precipitation Agent, mix the contents and incubate again at 37°C for 10 minutes.
4. Centrifuge the tubes at 12,000xg for 5 minutes or centrifuge 96 well titer plate at 2,000-4,000xg for 15 minutes.
5. Transfer the supernatant to clean tubes or 96-well plate. Add 120µl Assay Buffer in each tube/well and mix. A pink color will develop instantly.
6. If in tubes, transfer to a 96-well plate or micro-cuvette.
7. Read the reaction color at 570nm against the blank. The intensity of the color developed in each tube is proportional to protease activity.

NOTE: The reaction product may be diluted with water to 1ml and read in 1 ml cuvettes, however, the sensitivity will drop by 4-5 fold. Alternatively, the reaction volume may be increased proportionally to give total reaction volume 1ml.

RELATED PRODUCTS

Download our Protease & Phosphatase Inhibitors, Enzymes & Assays Handbook.



<http://info.gbiosciences.com/protease-phosphatase-inhibitors-enzymes-assay-handbook>

For other related products, visit our website at www.GBiosciences.com or contact us.



www.GBiosciences.com