

EN

DO40-5

Dissolved oxygen measuring device

Operating manual

European catalogue number:

EU connection: 664-0234

Version: 1.0

Issued: 07.01.2019

Table of contents

1	Legal address of the manufacturer	4
2	About this documentation	5
2.1	Foreword.....	5
2.2	Purpose of the document.....	5
2.3	Warranty	5
2.4	Compliance with local laws and other legal regulations	5
2.5	Correctness of content.....	6
2.6	Layout of this document.....	6
2.7	Further information	6
3	Safety	8
3.1	Explanation of safety symbols	8
3.2	Foreseeable misuse	8
3.3	Safety instructions	9
3.4	Intended use	10
3.5	Qualified personnel	10
4	Description	11
4.1	Scope of delivery	11
4.2	Job description.....	11
5	The product at a glance	12
5.1	The DO40-K.....	12
5.2	Display elements	12
5.3	Operating elements	12
6	Bases for measurement	14
6.1	Oxygen sensor.....	14
6.1.1	Explanation	14
6.1.2	Design.....	15
6.1.3	Service life	15
6.1.4	Operating position.....	16
6.1.5	Measurement accuracy	16
6.1.6	Residue.....	16
6.2	Instructions for oxygen measurement.....	16
6.2.1	Salinity correction	16
6.2.2	Environmental pressure, water depth and air pressure	17
6.3	Commissioning, filling and maintenance of the sensor.....	17
7	Maintenance	19
7.1	Operating and maintenance notices	19
7.2	Battery	19
7.2.1	Battery indicator	19
7.2.2	Changing battery	19
7.3	Calibration and adjustment	20
7.3.1	Automatic calibration in the air.....	20
8	Operation	22
8.1	Commissioning	22
8.1.1	Explanation	22

8.2	Configuration	22
8.2.1	Explanation	22
8.2.2	Opening the configuration menu.....	22
8.2.3	Configuring parameters of the configuration menu.....	23
8.2.4	Adjustment of the measuring input	24
8.2.5	Configuring parameters of the adjustment menu.....	25
9	Error and system messages.....	27
10	Disposal.....	28
11	Technical data.....	29
12	Service	30
12.1	Online resources.....	30
12.2	Sales partners.....	31

1 Legal address of the manufacturer

VWR International bvba

Researchpark Haasrode 2020

Geldenaaksebaan 464

B-3001 Leuven

Phone: + 32 16 385011

<http://be.vwr.com>



2 About this documentation

2.1 Foreword

Read this document carefully and familiarise yourself with the operation of the product before you use it. Keep this document ready to hand and in the immediate vicinity of the product so that it is available to the personnel/user for reference at all times in case of doubt.

The product was developed according to the state of the art and fulfils the requirements of the applicable European and national Directives. All corresponding documents are available from the manufacturer.

Only technically qualified persons are permitted to carry out commissioning, operation, maintenance and decommissioning. The qualified personnel must have carefully read and understood the operating manual before beginning any work.

2.2 Purpose of the document

- This document describes the operation and maintenance of the product.
- Provides important information for working safely and efficiently with the product.
- In addition to the quick reference guide with all relevant legal and safety content in hard copy, this document is a detailed reference option for the product.

2.3 Warranty

VWR guarantees that this product will remain free from material and production defects for a period of two (2) years from the time of delivery. If there is a defect, VWR shall decide on its own discretion whether to repair or replace the product free of charge or to provide the customer with a refund of the purchase price of the product, insofar as it is sent back during the warranty period. This warranty shall be voided if the product was damaged accidentally or intentionally due to improper use or normal wear. Insofar as the necessary maintenance tasks and inspections are not carried out in accordance with the operating manual and the local requirements, the warranty shall be voided, unless this omission is not causal for the defects of the product.

The customer must insure returned articles against damage and loss. This warranty is limited to the aforementioned rights. IT IS EXPRESSLY AGREED THAT THIS WARRANTY APPLIES INSTEAD OF ANY WARRANTY OF FITNESS AND INSTEAD OF WARRANTY OF GENERAL MERCHANTABILITY.

2.4 Compliance with local laws and other legal regulations

The customer is responsible for obtaining the necessary official approvals and other authorisations that are necessary for operation and use of the purchased product in its place of use. VWR can only be held liable if the customer fails to carry out the necessary activities for this purpose or authorisations are not issued unless a corresponding rejection is based on a defect of the product.

2.5 Correctness of content

The contents of this document were checked for corrected and are subject to a continuous correction and updating process. This does not rule out potential errors. In the event that errors are discovered or in case of suggestions for improvement, please inform us immediately via the indicated contact information in order to help us make this document even more user-friendly.

2.6 Layout of this document

Description

Each chapter is explained at the beginning in the description.

Prerequisite

All mandatory prerequisites are then listed for each step.

Instruction

Tasks to be carried out by the personnel / user are represented as numbered instructions. Adhere to the sequence of the specified instructions.

Representation

Shows an illustrative instruction or a configuration of the product.

Formula

Some instructions include a formula for a general understanding of a configuration, programming or a setting of the product.

Outcome of an action

Result, consequence or effect of an instruction.

Emphases

In order to simplify legibility and provide a clearer overview, various sections / information are emphasised.

- *1234* Display elements
- *Mechanical controls*
- **Product functions**
- *Product labels*
- Cross-reference [► p. 6]
- *Foot notes*

2.7 Further information

Software version of the product:

- V1.2 or later

For the exact product name, refer to the type plate on the rear side of the product.



NOTE

For information about the software version, press and hold the ON button to switch on the product for longer than 5 seconds. The series is shown in the main display and the software version of the product is shown in the secondary display.

3 Safety

3.1 Explanation of safety symbols



DANGER

This symbol warns of imminent danger which can result in death, severe bodily injury, or severe property damage in case of non-observance.



DANGER

This symbol indicates danger for living tissue as well as a variety of materials, which can be damaged or destroyed when coming into contact with this chemical. Caustic effect, protective equipment required!



DANGER

This symbol indicates danger for all life forms, which can result in death or acute or chronic damage to the health after inhaling, swallowing or absorbing this chemical through the skin.



CAUTION

This symbol warns of potential dangers or harmful situations which can cause damage to the device or to the environment in case of non-observance.



NOTE

This symbol indicates processes which can have a direct influence on operation or can trigger an unforeseen reaction in case of non-observance.



NOTE

This symbol instructs the use of eye protection which protects the eyes from harmful influences when working with powerful light, UV radiation, laser, chemicals, dust, splinters or weather influences.



NOTE

This symbol instructs the use of protective gloves which offer protection from mechanical, thermal, chemical, biological or electrical hazards.

3.2 Foreseeable misuse

The fault-free function and operational safety of the product can only be guaranteed if generally applicable safety precautions and the device-specific safety instructions for this document are observed.

If these notices are disregarded, personal injury or death, as well as property damage can occur.



DANGER

Incorrect area of application!

In order to prevent erratic behaviour of the product, personal injury or property damage, the product must be used exclusively as described in the chapter Description [p. 11] in the operating manual.

- Do not use in safety / Emergency Stop devices!
- The product is not suitable for use in explosion-prone areas!
- The product must not be used for diagnostic or other medical purposes on patients!
- The product is not intended to come into direct contact with food. For measurement in foods, samples must be taken and discarded after the measurement!

3.3 Safety instructions

This product has been designed and tested according to the safety requirements for electronic measuring devices.



DANGER

Potassium hydroxide!

The electrode contains potassium hydroxide. This causes burns. All contact with the skin, clothing and eyes should be avoided. Nevertheless, should contact occur, take the following measures.

- Eyes: Flush with flowing water for at least 15 minutes, seek medical attention!
- Skin: Wash with large amounts of water for several minutes!
- Clothing: Remove immediately!
- If swallowed: Drink large amounts of water, do not induce vomiting and seek medical attention!



CAUTION

Erratic behaviour!

On suspicion that the product can no longer be operated without danger, it must be decommissioned and prevented from recommissioning with appropriate labelling. The safety of the user can be impaired by the device if, for example, if it shows visible damage, it no longer works as specified or if it was stored for an extended period of time under unsuitable conditions.

- Visual inspection!
- In case of doubt, send the product to the manufacturer for repair or maintenance!



NOTE

If the product is stored at a temperature above 50 °C, or is not used for an extended period of time, the batteries must be removed. Leaks from the batteries are avoided as a result.



NOTE

This product does not belong in children's hands!

For this purpose, also refer to

 Technical data  29]

3.4 Intended use

The product is designed for analysis of the oxygen concentration and/or oxygen saturation in freshwater and salt water. It is used, for example, for the monitoring of wells, sewer lines and aquaria. A minimum flow to the sensor of approx. 30 cm/sec is necessary for a correct measurement.

See Technical data  p. 29].

3.5 Qualified personnel

For commissioning, operation and maintenance, the relevant personnel must have adequate knowledge of the measuring process and use of the measurements, for which purpose this document makes a valuable contribution. The instructions in this document must be understood, observed and followed.

In order to ensure that no risks arise from the interpretation of the measurements in the concrete application, the user must have additional technical knowledge, because the user is liable in case of damage/danger due to misinterpretation as a result of inadequate technical knowledge.

4 Description

4.1 Scope of delivery

Please check to ensure the completeness of the product after opening the package. You should find the following components:

- Quick reference guide
- Handheld measuring device, ready for operation, including batteries
- Permanently connected oxygen sensor

4.2 Job description

The product offers precision, speed and reliability in a compact, ergonomic housing. Additional impressive features include the dust-proof and waterproof design in accordance with IP 65/67 and the 3-line illuminated display, which offers overhead display at the push of a button. The product can be switched on, switched off and configured and the measurements and parameters can be adjusted and held with the operating elements. The product with maintenance-friendly galvanised O₂ sensor is an entry-level device suitable for everyday use. Concentrations in mg/l or ppm and saturation in percentage can be read directly without using tables. Calibration with environmental air takes place at the push of a button. Use of a GSKA protective cap is recommended for field use in bodies of water in order to protect the membrane.

5 The product at a glance

5.1 The DO40-K



LCD Display



DO40-K



DO40K

5.2 Display elements

Display



Battery indicator

Evaluation of the battery status



Unit display

Display of units, if applicable, with unstable symbol or type of mode, min/max/hold



Main display

Measurement of the current O₂ value or value for min/max/hold

Auxiliary display

Corresponding temperature for the displayed O₂ value with unit.

Bar graph

Progress for calibration and visualisation of the electrode evaluation



NOTE

The unit display shows a rotating circle segment in the first position as long as the measurement is unstable, if the position is unoccupied by the unit display.

5.3 Operating elements



On / Off button

Press briefly

Switch on the product

Activate / deactivate lighting

Long press

Switch off the product

Reject changes in a menu



Up / Down button



Press briefly

Display of the min/max value

Change value of the selected parameter

Long press

Reset the min/max value of the current measurement

Both simultaneously

Rotate display, overhead display



Function key

Press briefly

Freeze measurement

Return to measurement display

Call up next parameter

Long press, 2s

Start menu configuration, CONF appears in the display

Long press, 4s

Start automatic calibration, CAL appears in the display

6 Bases for measurement

6.1 Oxygen sensor

6.1.1 Explanation

The oxygen sensor is an active sensor. It consists of a platinum cathode, a lead anode and potassium hydroxide (KOH) as an electrolyte. If oxygen is present, it is reduced on the platinum cathode and the sensor delivers a signal. If no oxygen is present, no signal is delivered. The anode is consumed by the oxygen measurement. The sensor ages. Furthermore, the sensor loses water through the permeable membrane, in particular, when it is stored in dry air. Therefore, it should be checked and maintained regularly and replaced as necessary.



DANGER

Potassium hydroxide!

The electrode contains potassium hydroxide. This causes burns. All contact with the skin, clothing and eyes should be avoided. Nevertheless, should contact occur, take the following measures.

- Eyes: Flush with flowing water for at least 15 minutes, seek medical attention!
- Skin: Wash with large amounts of water for several minutes!
- Clothing: Remove immediately!
- If swallowed: Drink large amounts of water, do not induce vomiting and seek medical attention!



NOTE

Protective goggles must be worn for all of the following activities!



NOTE


Protective gloves must be worn for all of the following activities!



NOTE

Always store the oxygen sensor damp. It should always be stored in a storage bottle filled with water or in a container filled with water. After storage for an extended period, any potential deposit layers, such as algae, must be cleaned off of the membrane with a soft paper towel prior to measurement.

6.1.2 Design



Platinum electrode	If oxygen is present, it is reduced on the platinum electrode and the sensor delivers a signal. Soiling on the platinum electrode or between the membrane and electrode can influence the measurement.
Storage bottle	The storage bottle is provided to keep the membrane moist. The service life of the sensor is extended as a result. Distilled or deionised water is in the storage bottle; do not add any other liquids!
Membrane head	The membrane head is covered with a thin plastic membrane. Faulty measurements will occur if the membrane is damaged or there are large air bubbles or even an air bubble ring on the membrane. This can also be the cause if a sensor can no longer be calibrated. The GWOK 02 membrane head is a spare part and can be re-ordered separately.
Refill opening	Electrolyte must be filled or added for the initial commissioning of a sensor which is delivered dry, when performing maintenance or after use at high temperatures.

6.1.3 Service life

The sensor signal deteriorates relatively quickly at the end of the service life of the sensors. The electrode evaluation in %, therefore, can only be used as a guide value. A value of 70% does not mean that exactly 70% of the service life is still available, rather that the electrode signal has 70% of a comparison signal.



NOTE

The sensor evaluation is updated by the measuring device after a successfully performed calibration of the oxygen sensor.

The nominal service life can be reduced significantly due to use. Influential factors include:

- Storage / operating temperature
- Contamination of the measured water

- Mechanical stress of the sensor membrane
- Storage in dry air
- Continuous use in elevated carbon dioxide concentrations

6.1.4 Operating position

The oxygen sensor should be arranged vertically upwards with the connecting cable. A slight angle of inclination does not impair the measurement.

6.1.5 Measurement accuracy

The measurement accuracy can be impaired by:

- An inadequate flow below the necessary value of approx. 30cm/sec.
- The water temperature and sensor temperature must be the same. The most accurate measurements are provided when the measuring temperature is calibrated.

6.1.6 Residue

Visible residue collecting as a reaction product in the interior of the membrane cap on the lead-oxide lead anode is red or brown from reaction with oxygen and lead carbonate and is white from reaction with carbon dioxide. These substances can collect on the membrane, but do not normally impair the measuring function and can be mostly removed in the course of maintenance of the sensor. Prior to screwing on the membrane cap, they should be removed as far as possible in order to prevent particles from being trapped between the membrane and the platinum cup. A rapid or excessive formation of lead carbonate after commissioning is an indication of air in the sensor. This is usually due to incomplete filling or a leak due to improper fitting of the cap / fill screw or membrane leak.

6.2 Instructions for oxygen measurement

The following must be observed when measuring dissolved oxygen:

- The storage bottle must be removed before the measurement.
- The sensor must have been calibrated.
- The sensor and the liquid to be measured must have the same temperature. The temperatures of the two must match.
- The sensor must be immersed at least 3 cm into the liquid to be measured.
- A flow speed of at least approx. 30 cm/sec is necessary for exact measurements. Either stir continuously or use an appropriate stirring device
- The measurement is sensitive to jarring! Therefore, make sure that the container is not struck with the sensor while stirring, because this can significantly impair the measured value.
- The oxygen partial pressure, oxygen concentration in mg/l and oxygen saturation in % are calculated from the sensor signal and temperature. In accordance with DIN38408-C22, the measurement is relative to water-vapour-saturated air.

6.2.1 Salinity correction

With increasing salinity *SRL*, which is the value for the salt content in the water, the solubility of oxygen in water decreases, i.e. with the same oxygen partial pressure, fewer mg of oxygen are dissolved per litre of water. To determine this oxygen concentration, therefore, the salinity of the medium must be entered first; see Configuring parameters of the configuration menu [► p. 23]. A salinity correction is not necessary in freshwa-

ter; the value is 0. A salinity of approx. 35 PSU is normal for sea water. The salinity correction is adjusted to aqueous media having a chemical composition corresponding to sea water. The International Oceanographic Tables (IOT) are used as a basis for the correction.

6.2.2 Environmental pressure, water depth and air pressure

The environmental pressure, water depth and air pressure play an important role at the place of measurement for the following points:

- Calculation of oxygen saturation in % ~~5%~~. Clean water can achieve 100 % saturation in air. There must be no oxygen-depleting processes, such as biological decaying processes, chemical effects or oxygen enhancing processes, such as excessive ventilation or photosynthesis. This could result in an oversaturation above 100%.
- Calculation of the oxygen concentration in mg/l
- The valuation of calibration

Adjusting the pressure parameter on the product is recommended prior to calibration. In the scope of measuring accuracy, specifying the current air pressure in the region based on meteorological data or the standard pressure based on sea level is sufficient.

For example:

0 m above NN: 1013 hPa

300 m above NN: 978 hPa

600 m above NN: 943 hPa

1000m above NN: 899 hPa

6.3 Commissioning, filling and maintenance of the sensor

Description	The sensor is delivered dry. Therefore, the sensor is well-suited for storage. The sensor must be filled in good time before the measurement. A wait time of approx. 2 hours after filling should be planned in order to allow the sensor to stabilise. Commissioning with initial filling, filling and maintenance of the sensor are explained in the following chapter.
Prerequisite	– Protective goggles – Protective gloves – A suitable screwdriver – Pipette – Household towel – KOH electrolyte – A spare GWOK 02 membrane head, if applicable
Instruction	1. Unscrew the membrane head. 2. Unscrew the seal screw from the refill opening.



NOTE

If the sensor was already filled and deposits have formed, it should be cleaned with KOH via the fill opening or removed. Loose residue can be rinsed out in the process. The platinum cup on the front of the membrane must be clean. Remove any soiling or electrolyte solution with a paper towel.

3. Fill the pipette with the KOH electrolyte and initially fill the membrane head to $\frac{3}{4}$ full. Rinse off the excess electrolyte.
4. Slowly fill the sensor via the fill opening. In the process, tilt the sensor from side to side and tap the shaft to force out air bubbles. The sensor holds approximately 5 ml. If air bubbles no longer appear and the refill opening is filled up to the rim with KOH, screw the in the seal screw again.
5. Rinse off excess KOH and screw the sensor on with the membrane head. In the process, if there are any recognisable air bubbles under the membrane, KOH must be added.
6. After filling, the sensor should rest for 2 hours before a calibration is started.

Outcome of an action

The sensor has now been refilled. A sensor evaluation during the calibration should produce a result of 100%.



NOTE

Should it no longer be possible to calibrate the sensor or it only delivers unstable measurements, it must be maintained and/or the membrane head must be replaced

7 Maintenance

7.1 Operating and maintenance notices



NOTE

If the product is stored at a temperature above 50 °C, or is not used for an extended period of time, the batteries must be removed. Leaks from the batteries are avoided as a result.



NOTE

If the product is stored at a temperature above 50 °C, or is not used for an extended period of time, the batteries must be removed. Leaks from the batteries are avoided as a result.



NOTE

The electrode should be stored in dry rooms at a temperature between 10 °C and 30 °C. If the storage temperature range is exceeded or undercut, the electrode can be destroyed. It should always be stored wet in distilled or deionised water.

7.2 Battery

7.2.1 Battery indicator

If the empty frame in the battery display blinks, the batteries are depleted and must be replaced. However, the device will still operate for a certain length of time.

If the *bAt* display text appears in the main display, the battery voltage is no longer adequate for operation of the product. Now the battery is fully depleted.

7.2.2 Changing battery



DANGER

Danger of explosion!

Using damaged or unsuitable batteries can generate heat, which can cause the batteries to crack and possibly explode!

- Only use high-quality and suitable alkaline batteries!



CAUTION

Damage!

If the batteries have different charge levels, leaks and thus damage to the product can occur.

- Use new, high-quality batteries!
- Do not use different types of batteries!
- Remove depleted batteries and dispose of them at a suitable collection point!



NOTE

Unnecessary screwing places the water-tightness of the product, among other things, at risk and should be avoided.




NOTE

Read the following handling instructions before replacing batteries and follow them step by step. If disregarded, the product could be damaged or the protection from moisture could be diminished.

Description

Prerequisites

Instruction



Proceed as follows to replace the batteries.

- The product is switched off.
1. Unscrews the Phillips screws and remove the cover.
 2. Carefully replace the two Mignon AA batteries. Ensure that the polarity is correct! It must be possible to insert the batteries in the correct position without using force.
 3. The O-ring must be undamaged, clean and positioned at the intended depth. In order to facilitate assembly and avoid damage, a suitable grease can be applied.
 4. Fit the cover on evenly. The O-ring must remain at the intended depth!
 5. Tighten the Phillips screws.

Outcome of an action

The product is now ready for use again.

7.3 Calibration and adjustment

7.3.1 Automatic calibration in the air

Description

Prerequisite

The following steps describe how to calibrate the product automatically.

- The product is switched on.



NOTE

Calibration takes place in air saturated with water vapour. The GCAL 3610 calibration container or the storage bottle can be used for this purpose. The membrane of the electrode must be dry for the calibration. Prior to performing the calibration of the membrane, dab away water droplets from the membrane with a soft, dry cloth. The following must be observed when using the storage bottle:

Only insert the electrode far enough into the storage bottle that the membrane does not come into contact with the water in the bottle.

	Unscrew the lid of the storage bottle and only position it so that a small exchange of air and pressure equalisation can take place
Instruction	1. Place the electrode in the calibration container. If necessary, wait until the temperature has equalised and a stable value has been achieved. 2. Press the <i>Function key</i> for 4 seconds to open the Calibration menu. The display shows $\overline{RL}$. 3. Release the <i>Function key</i>. 4. The product determines the correct value automatically.
Outcome of an action	After successful completion of the calibration the assessment of the electrode condition is displayed briefly in percent. An aged or contaminated electrode, incorrect adjustment of the pressure, contamination of the platinum electrode or a damaged membrane can be the cause for a lower evaluation. If the calibration is not completed successfully an error message is displayed. $\overline{RL Err.}$ appears in the display See Error and system messages [► p. 27]. Confirm the error message pressing the <i>Function key</i>. The product restarts and the value of the last successful calibration is restored.

8 Operation

8.1 Commissioning

8.1.1 Explanation

Description	The product is switched on with the <i>On/Off button</i> . It may be necessary to configure the product after switching on. See Configuration [► p. 22].		
Prerequisite	– Sufficiently full batteries are inserted in the product.		
Instruction	– Press <i>On/Off button</i> .		
Outcome of an action	Information about the configuration of the product appears in the display.		
	<i>P_{OFF}</i>	Automatic shut-off	Automatic shut-off activated. The product is switched off if no buttons have been pressed after the adjusted time
	<i>t_{0F}</i>	Zero point correction	If a zero point correction of the temperature sensor was made
	<i>t_{5L}</i>	Gradient correction	If a gradient correction of the temperature sensor was made
	<i>S_{RL}</i>	Salinity correction	Links when the salinity correction is active
The product is now ready for measurement.			



NOTE

The product must be calibrated to the electrode prior to starting the measurement. See Calibration and adjustment service.

8.2 Configuration

8.2.1 Explanation

The following steps describe how to adapt the product for your purposes.



NOTE

There are various configuration parameters available depending on the product version and configuration. They can differ depending on the product version and configuration.

8.2.2 Opening the configuration menu

Description	In order to configure the product, you must first open the Configuration menu. The menu is opened as shown in the illustration.
Prerequisite	– The product is switched on.
Instruction	1. Press the <i>Function key</i> for 2 seconds to open the Configuration menu. 2. <i>t_{0F}</i> appears in the display. Release the function key. 3. By briefly pressing the <i>Function key</i>, you can scroll through the parameters. Select the parameter you would like to configure.

4. When you have selected the desired parameter, change the parameter to the desired value with the *Up button* and the *Down button*.
5. The changes are saved after running through the entire **Configuration** menu. **Star** appears in the display. The **Configuration** menu can be exited from any arbitrary parameter by pressing and holding the *Function key* for 2 seconds. The changes made up that point are saved.

Representation

Call up menu	Next parameter	Change value	Save changes	Discard changes
				
2s		Press: Single step Hold: Rapid change	2s	2s

Outcome of an action

The **Configuration** menu is closed after the last parameter.



NOTE

If the product is switched off without saving the configuration, the last save value is reproduced on the next start-up of the product.

8.2.3 Configuring parameters of the configuration menu

Description

The following representation shows the available parameters and various configuration options.

Prerequisite

- The **Configuration** menu is open. See Opening the configuration menu [► p. 22].

Instruction

1. Select the desired parameter you would like to configure.
2. Adjust the desired configuration in the selected parameter with the *Up button* and *Down button*.
3. The available configuration options are listed for each parameter in the following representation.

Representation

Parameter	Values	Meaning
		
Input		
<i>inP</i>		
	<i>SRt</i> %	Oxygen saturation in per cent
	<i>Conc</i> mg/l	Oxygen concentration in mg/l
	<i>Conc</i> ppm	Oxygen concentration in ppm
Pressure		
<i>SEt.P</i>		
	500 .. 4000	Environmental pressure in hPa, corresponding to mbar
Salinity correction		
<i>SRt</i>		
	0 .. 70	Salinity in the measuring medium in PSU, corresponding to g/kg

Shut-off time

*PoFF**oFF*

No automatic shut-off

15 30 60 120 240

Automatic shut-off after a selected time in minutes, during which no buttons have been pressed

Backlighting

*L, tE**oFF*

Backlighting deactivated

15 30 60 120 240

Automatic shut-off of the backlighting after a selected time in seconds, during which no buttons have been pressed

on

No automatic shut-off of the backlighting

Temperature unit

*Unit**°C*

Temperature display in °C

°F

Temperature display in °F

Factory settings

*Init**no*

Use current configuration

*YES*Reset product to factory settings. *Init done* appears in the display

Outcome of an action

The changed value is saved and the **Configuration** menu is closed. *Star* appears in the display. If necessary, the product is restarted automatically in order to adopt the changed values.



NOTE

The configuration is closed if no button is pressed for 2 minutes. Any changes made up to that point are not saved. *End* appears in the display.

8.2.4 Adjustment of the measuring input

Description

The temperature input can be adjusted with the zero point correction and the gradient correction. If an adjustment is made, you change the pre-adjusted factory settings. This is signalled with the *t.oF* or *t.5L* when the product is switched on. The standard settings of the zero point value and the gradient value is *0.00*. It signals that no correction is made.

In order to adjust the product, you must first open the **Adjustment** menu. The menu is opened as shown in the illustration.

Prerequisites

- Sufficiently full batteries are inserted in the product.
- The product is switched off.

Instruction

1. Press and hold the *Down button*.
2. Press the *On/Off button* to switch on the product and open the **Configuration** menu. Release the *Down button*. The display shows the first parameter.
3. By briefly pressing the *Function key*, you can scroll through the parameters. Select the parameter you would like to configure.

4. When you have selected the desired parameter, change the parameter to the desired value with the *Up button* and the *Down button*.
5. In order to save the new parameter value, press and hold the *Function key* for longer than 1 second.

Representation

Call up menu



Hold

Release

Outcome of an action

The **Configuration** menu is closed after the last parameter.

NOTE

If the product is switched off without saving the configuration, the last save value is reproduced on the next start-up of the product.

8.2.5 Configuring parameters of the adjustment menu

Description

The following representation shows the available parameters and various configuration options.

Prerequisites

The **Adjustment** menu is open. See Adjustment of the measuring input [► p. 24].

Instruction

1. Select the desired parameter you would like to configure.
2. Adjust the desired configuration in the selected parameter with the *Up button* and *Down button*.
3. The available configuration options are listed for each parameter in the following representation.

Representation

Parameter	Values	Meaning
Zero point correction		
$\pm 0F$		
	0.00	No zero point correction
	-5.00 .. 5.00	Zero point correction in °C. and/or at °F -9.00 .. 9.00
Gradient correction of the temperature		
$\pm 5L$		
	0.00	No gradient correction
	-5.00 .. 5.00	Gradient correction in %

Formula

Zero point correction:

Displayed value = measured value – $\pm 0F$

Gradient correction °C:

Display = (measured value – $\pm 0F$) * (1 + $\pm 5L$ / 100)

Gradient correction °F:

Display = (measured value – 32 °F – $\pm 0F$) * (1 + $\pm 5L$ / 100) + 32 °F

Example calculation

- Zero point correction $\pm 0F$ to 0.00
- Gradient correction $\pm 5L$ to 0.00
- Display unit *Unit* to °C
- Display in ice water -0.2 °C

- Display in ice water setpoint $t_{\text{dF}} = 0.0 \text{ }^{\circ}\text{C}$
- Display in water bath $36.6 \text{ }^{\circ}\text{C}$
- Display in water bath setpoint $t_{\text{5L}} = 37.0 \text{ }^{\circ}\text{C}$
- t_{dF} = display zero point correction – setpoint zero point
- $t_{\text{dF}} = -0.2 \text{ }^{\circ}\text{C} - 0.0 \text{ }^{\circ}\text{C} = -0.2 \text{ }^{\circ}\text{C}$
- $t_{\text{5L}} = (\text{setpoint gradient correction} / (\text{display gradient correction} - t_{\text{dF}}) - 1) * 100$
- $t_{\text{5L}} = (37.0 \text{ }^{\circ}\text{C} / (36.6 \text{ }^{\circ}\text{C} - (-0.2)) - 1) * 100 = 0.54$

Outcome of an action



The changed value is saved and the **Configuration** menu is closed.

NOTE

If the product is switched off without saving the configuration, the last save value is reproduced on the next start-up of the product.

9 Error and system messages

Display	Meaning	Possible causes	Remedy
No display, unclear characters or no response when buttons are pressed	Battery depleted	Battery depleted	Replace battery
	System error	Error in the product	Send in for repair
	Product is defective	Product is defective	
<i>bAt</i>	Battery depleted	Battery depleted	Replace battery
<i>bAt Lo</i>	Battery depleted	Battery depleted	Replace battery
<i>Cal Err.2</i>	Slope is too low Incorrect oxygen reference	Electrode contaminated or defective	Perform calibration in damp environmental air Maintain the electrode
<i>Cal Err.3</i>	Slope is too high Incorrect oxygen reference	Electrode contaminated or defective	Perform calibration in damp environmental air Maintain the electrode
<i>Cal Err.4</i>	Incorrect calibration temperature	Temperature too low or too high	Range of 5..40 °C
<i>Cal Err.5</i>	Time exceeded during automatic calibration	Unstable electrode signal Contaminated electrode Temperature not equalised	Use calibration container Maintain the electrode Restart calibration
<i>Err.1</i>	Measuring range exceeded	Measurement too high Electrode or product defect Faulty calibration	The measurement is above the permissible range Check electrode Perform calibration Send in for repair
<i>Err.2</i>	Measuring range is undercut	Measurement too low Electrode or product defect	Check electrode Send in for repair
<i>SYS Err</i>	System error	Error in the product	Switch product on/off Replace batteries Send in for repair

10 Disposal



NOTE

This product bears the "recycle batteries" symbol. It means that this product must not be disposed of with household waste. Instead, it is your responsibility to dispose of the product correctly at the end of its service life by dropping it off at an authorised disposal facility which collects the product separately for recycling. It is also your responsibility to decontaminate the equipment in the event of biological, chemical and/or radioactive pollution in order to avoid putting any persons assigned with the disposal and recycling at risk. For further information about the location where you can drop off your product, contact your local dealer or the place where you originally purchased the product.

If you proceed as described above, you are contributing to the protection of natural resources and ensuring that your product is recycled so that no one's health is put at risk.

Thank you!

11 Technical data

Measuring range	O ₂ concentration	O ₂ saturation	Temperature
	0.0 .. 20.0 mg/l	0 .. 200 %	0 .. 50 °C
	0.0 .. 20.0 ppm		32 .. 122 °F
Accuracy (at nominal temperature)	±1.5 % of measured value	±1.5 % of measured value	± 0.3 °C
	± 0.2 mg/l	± 0.2 %	
Temperature compensation		0 .. 50 °C (or 32 .. 122 °F)	
Nominal temperature		25°C	
Measuring cycle		approx. 2 measurements per second	
Connections		Permanently connected oxygen sensor	
Display		3-line segment LCD, additional symbols, illuminated (adjustable white, permanent illumination)	
Additional functions		Min/Max/Hold	
O ₂ calibration		Automatic calibration in the air	
Housing		Break-proof ABS housing	
	Protection rating	IP65 / IP67	
	Dimensions L*W*H [mm] and weight	108 * 54 * 28 mm without electrode 130 g, incl. battery, without electrode 190 g, incl. battery and electrode	
Connections		Permanently connected oxygen sensor	
Operating conditions	Device	-20 to 50 °C; 0 to 95 % r.h. (temporarily 100 % r.h.)	
	Electrode	0 .. 40 °C	
Storage temperature		0 .. 40 °C	
Current supply		2*AA battery (included in the scope of delivery)	
	Current requirement/ battery life	approx. 0.8 mA, approx. 2.7 mA with lighting Service life > 3000 hours with alkaline batteries (without backlighting)	
	Battery indicator	4-stage battery status indicator, Replacement indicator for depleted batteries: "BAT"	
Auto-power-OFF function		The device switches off automatically if this is activated	
Directives and standards		The devices conform to the following Directives of the Council for the harmonisation of legal regulations of the Member States: 2014/30/EU EMC Directive 2011/65/EU RoHS Applied harmonised standards: EN 61326-1:2013 Emission limits: Class B Immunity according to Table 2 Additional errors: < 0.5 % FS EN 50581:2012 The device is intended for mobile use and/or stationary operation in the scope of the specified operating conditions without further limitations.	

12 Service

12.1 Online resources

The following information is provided on the VWR website at www.vwr.com:

- All contact information of the technical customer service
- VWR online catalogue and information about accessories and compatible products
- Additional product information and special offers

Contact If you require information or technical support, please contact your VWR sales centre or visit our website: www.vwr.com.

12.2 Sales partners

Your sales partner		
Australia VWR International GmbH Graumannsgasse 7 1150 Vienna Tel.: +43 01 97 002 0 Email: info.at@vwr.com Belgium VWR International bvba Researchpark Haasrode 2020 Geldenaaksebaan 464 3001 Leuven Tel.: +32 016 385 011 Email: vwr.be@vwr.com China VWR International China Co., Ltd Shanghai Branch Room 256, No. 3058 Pusan Road Pudong New District Shanghai 200123 Tel.: +86-21-5898 6888 Fax: +86-21-5855 8801 Email: info_china@vwr.com Czech Republic VWR International s. r. o. Veetee Business Park Pražská 442 CZ - 281 67 Stříbrná Kalice Tel.: +420 321 570 321 Email: info.cz@vwr.com Denmark VWR International A/S Tobaksvejen 21 2860 Søborg Tel.: +45 43 86 87 88 Email: info.dk@vwr.com Finland VWR International Oy Valimotie 9 00380 Helsinki Tel.: +358 09 80 45 51 Email: info.fi@vwr.com	Germany VWR International GmbH Hilpertstraße 20a D - 64295 Darmstadt Toll-free: 0800 702 00 07 Tel.: +49 (0) 6151 3972 0 (international) Email: info.de@vwr.com Hungary VWR International Kft. Simon László u. 4. 4034 Debrecen Tel.: +36 (52) 521-130 Email: info.hu@vwr.com India VWR Lab Products Private Limited No.139. BDA Industrial Suburb, 6th Main, Tumkur Road, Peenya Post, Bangalore, India – 560058 Tel.: +91-80-28078400 Email: vwr_india@vwr.com Ireland / Northern Ireland VWR International Ltd / VWR International (Northern Ireland) Ltd Orion Business Campus Northwest Business Park Ballycoolin Dublin 15 Tel.: +353 01 88 22 222 Email sales.ie@vwr.com Italy VWR International S.r.l. Via San Giusto 85 20153 Milano (MI) Tel.: +39 02-3320311 Email: info.it@vwr.com	Poland VWR International Sp. z o.o. Limbowa 5 80-175 Gdansk Tel.: +48 058 32 38 200 Email: info.pl@vwr.com Portugal VWR International – Material de Laboratório, Lda Centro Empresarial de Alfragide Rua da Indústria, nº 6 2610-088 Amadora Tel.: +351 21 3600 770 Email: info.pt@vwr.com Singapore VWR Singapore Pte Ltd 18 Gul Drive Singapore 629468 Tel: +65 6505 0760 Email: sales.sg@vwr.com Spain VWR International Eurolab S.L. C/ Tecnología 5-17 A-7 Llinars Park 08450 - Llinars del Vallès Barcelona Tel.: +34 902 222 897 Email: info.es@vwr.com Sweden VWR International AB Fagerstagatan 18a 163 94 Stockholm Tel.: +46 08 621 34 00 Email: info.se@vwr.com Switzerland VWR International AG Lerzenstrasse 16/18 8953 Dietikon Tel.: +41 044 745 13 13 Email: info.ch@vwr.com

France VWR International S.A.S. Le Périgares – Bâtiment B 201, rue Carnot 94126 Fontenay-sous-Bois cedex Tel.: 0 825 02 30 30* (national) Tel.: +33 (0) 1 45 14 85 00 (international) Email: info.fr@vwr.com * €0.18 TTC/min	Netherlands VWR International B.V. Postbus 8198 1005 AD Amsterdam Tel.: +31 020 4808 400 Email: info.nl@vwr.com Norway VWR International AS Haavard Martinsens vei 30 0978 Oslo Tel.: +47 22 90 00 00 Email: info.no@vwr.com	Great Britain VWR International Ltd Customer Service Centre Hunter Boulevard Magna Park Lutterworth Leicestershire LE17 4XN Tel.: +44 (0) 800 22 33 44 Email: uksales@vwr.com
---	--	---