Vacuum Aspiration System (Hereinafter, the “Device”)

INSTRUCTION MANUAL

North American Catalog Number

75870-734

Version: 1
Issued: 17th, Jan, 2017

Disclaimer: The name and address of the Manufacturer of this Device is **confidential.**

Contents

Safety information.....	1
Package contents.....	1
Unpacking.....	2
Installation.....	2
Intended use.....	2
Device specifications.....	2
Instruction for use.....	3
Description of the main parts.....	3
Operation.....	4
Troubleshooting.....	5
Chemical compatibility at 20°C.....	6
Technical service.....	7
Warranty.....	7
Compliance with local laws and regulations.....	9
Device disposal.....	9

Safety Information

Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

The Manufacturer recommends the following guidelines for the use of the Device:

- The inside of the Device is a high-voltage area. Adequate Precaution to be taken and the Device handled by trained personnel only.
- Do not disassemble the Device.
- Disconnect the Device from the mains supply before carrying out maintenance work or changing the fuses.
- Do not use the Device if it is damaged, especially if the main power cable is in any way damaged or defective.
- Repairs may only be carried out by the service technicians from your local VWR office and authorized service partners.
- The Device must be connected to a power outlet that has a protective ground connection.
- If the Device is used in a manner not specified by the manufacturer, the protection provided by the Device may be impaired.
- Do not allow any liquid to enter into the Device control unit.
- Do not operate the Device in a hazardous location or potentially explosive environment.

Package Contents

Main unit	1 pc
Bottle	1 pc
Lid	2 pcs
Operating handle with tube	1 pc
Syringe needle(60mm)	2 pcs
Syringe needle (120mm)	2 pcs
8-channel syringe needle	1 pc
Single-channel tip extrusion adapter	1 pc
8-channel tip extrusion adapter	1 pc
Vacuum tube	1 pc
Power adapter	1 pc
User manual	1 pc
Warranty card	1 pc

Unpacking

Open the package and check the contents per the packing list, if any part is found missing or is damaged please contact your local VWR office / authorized partners.

Installation

Placement-Place the Device on a level surface carefully.

Link the Power Cord - Check to confirm that the Device's power switch is turned off. Plug the power cord into the power supply socket.

Intended Use



The Device is designed for laboratory waste recovery and separation of liquid and solid. It can be widely used in cell culture, DNA extraction, micro-plate waste removal, and any other liquid separation or recovery.

Device Specifications

Vacuum range	0 ~600 mbar (adjustable)
Pump flow rate (rated value)	15 L/min (air)
Aspiration rate of liquids	15 mL/s (aspiration pipette)
Dimensions (W×D×H)	180 ×240×450 mm
Weight	3 kg
Input voltage	100 ~240 VAC, 50/60 Hz
Bottle volume	4 L

Instruction for Use

Overview

Description of Parts

Main unit	Main part of the device, has a vacuum pump inside.
Bottle	Bottle for liquid collection.
Lid	Lid with opening for tubing has the fast-inserting joints, level sensor interface and the air release valve. Note: Please make sure every part is put into position. When you want to disconnect the tubes from the fast-inserting joint,

	hold the metal tab and pull out.
Hand operator	Aspirate liquid by pressing the button of the hand operator.
Liquid level sensor	Used to prevent the overflow of liquids. When an overflow is detected the pump will switch off and the alarm will sound.
Tubing for vacuum	Connect the main unit to the bottle, providing negative pressure environment.
Tubing for aspirating liquid	Connect the hand operator to bottle, aspirates the liquid to the bottle.
Bottle handle	Install the bottle handle between the lid and the bottle, in order to hold the hand operator and as a handle when moving the bottle.
Air release valve	Unscrew the valve to release the vacuum. Note: when you want to open the lid, please release the vacuum first.
Vacuum pressure indicator	Indicates the degree of vacuum pressure. When indicator is flashing the vacuum pressure is building, when indicator stays lit the desired vacuum pressure has been reached.
Vacuum regulation knob	Regulates the setting of the vacuum pressure.
Filter	Connect between the main unit and the bottle to prevent liquid from entering the pump.
Power interface	Connect the power cable.
Power switch	Power on/off.

Operation

- 1) Ensure that all parts are in position.
- 2) Make sure the air release valve is closed.
- 3) Power on.
- 4) Adjust the vacuum regulation knob clockwise to increase the vacuum pressure, vacuum pressure indicator will turn on in sequence. Adjust the vacuum regulation knob counter-clockwise to decrease the vacuum pressure; vacuum pressure indicator will turn off in sequence.
- 5) When the vacuum pressure indicator remains lit, the desired vacuum pressure has been reached. You can now aspirate liquid by pressing the button on the hand operator.

Note: The liquid level sensor is used to prevent the overflow of the liquids. When an overflow is detected the pump will switch off and the alarm will sound. Dispose of the

liquid in the bottle or replace with an empty bottle immediately.

6) Power off.

7) Unscrew the valve to release the vacuum.

Note: when you want to open the lid, please release the vacuum first.

8) Remove the bottle, open the lid, and dispose the liquid.

Troubleshooting

Review the information in the table below to troubleshoot operating problems.

Operating problem	Possible causes	Counter measure
Power switch is on, LED is not lit.	1) Power Failure 2) Faulty Power Adapter 3) Damaged Power Cable	1) Remove the power adapter cable from mains power socket and measure and it should be 100 – 230V AC. if it is not as above, one needs to check the power source. 2) Visually inspect the power adapter and its cable for any physical damage. if there is any damage observed, replace the power adapter. 3) Measure the output voltage from the power adapter and it should be 24V DC. If it is not measured as mentioned above, replace the power adapter.
Cannot adjust the vacuum pressure.	1) Bad Contact from Vacuum Regulating Encoder 2) Faulty Vacuum Regulating Encoder	1) Check if the cable connection from the vacuum regulating encoder is intact. If it is not intact connect properly. 2) Replace the encoder.
Cannot aspirate liquid.	1) Main Board Failure 2) Pump Failure	1) Check if the output voltage from main board to pump is 24V. if the value is abnormal, we need replace the main board. 2) If the output voltage from the PCB is 24V and still there is no aspiration, replace the pump.

Aspiration speed is slow.	1) Liquid barrel cover is not tight or air leak. 2) Air valve is not closed 3) Tube leakage.	1) Tighten the barrel cover or replace the cover seal ring. 2) Close the air valve. 3) Check all the tubes and connectors, replace leaking parts.
---------------------------	--	---

Chemical compatibility at 20°C

1. For biological reagents, the Device has good compatibility to almost all the liquids.
2. For chemical reagents, please refer to the below table.

Notice : Compatible "+", Incompatible "-"			
Reagent	Compatibility	Reagent	Compatibility
Acetic acid	+	Formic acid	-
Acetaldehyde	-	Gallic acid	-
Acetamide	+	Glycol	+
Acetic anhydride	-	Hydrogen peroxide	-
Acetone 70°F (21.1°C)	-	Hydrazine	+
Acetic acid glacial	-	Hydrobromic acid	-
Acetic acid	+	Hydrochloric acid	-
Acrylonitrile	-	Hydrofluoric acid	-
Adipic acid	+	Hydrofluosillicic acid	-
Alcohols	+	Hypochloric acid	-
Aliphatic acid	+	Isopropyl acetate	-
Alkali liquor	+	Isopropyl ether	+
Aluminium nitride	-	Lactic acid	+
Aluminium nitrate	+	Linoleic acid	+
Aluminium potassium sulfate	-	Malic acid	+
Ammonium carbonate	-	Methanol	+
Ammonium hydroxide	-	Methyl acetate	-
Aniline	-	Methyl ethyl ketone	-
Anone	-	N-Amyl acetate	-
Arsenic acid	+	N-Butyl acetate	-
Aqua regia	-	Nitronic acid	-
Benzoic acid	-	Nitroxanthic acid	+
Boric acid	+	Nitric acid	-
Butyl amine	-	Nitramine	+
Butyric acid	-	Oleic acid	-

Calcium hydroxide	+	Palmitic acid	+
Calcium hypochlorite	+	Perchloric acid	-
Calcium nitrite	+	Potassium hydroxide	+
Carbinol	+	Potassium hypochlorite	-
Carbonic acid	+	Propyl acetate	-
Copper sulfate	+	Propanol ethyl carbinol	+
Cupric cyanide	+	Propylene	-
Diacetone alcohol	-	Salicylic acid	+
Diethylamine	+	Sodium acetate	+
Diethylene glycol	+	Sodium hydroxide	+
Diethyl ether	-	Sulfuric acid	-
Dimethylaniline	-	Tetrachloroethylene	-
Dimethylformamide	+	Trichloroacetic acid	-
Ethanolamine	+	Tung oil	+
Ethyl benzenecarboxylate	-	Zinc sulfate	+
Ether	-		

Notes:

Please note that the table is just a directional guide not the manufacturer's commitment. Please read the user manual carefully before use. Good laboratory practice would be to rinse out the liquid handling unit at the end of each day with distilled water to prevent corrosive liquids being left in contact with the parts for too long.

Technical Service

Web Resources

Visit the VWR website at www.vwr.com for:

- Complete technical service contact information.
- Access to the VWR Online Catalogue, and information about accessories and related products.
- Additional product information and special offers.

Contact Us

For information or technical assistance contact your local VWR representative or visit. www.vwr.com.

Warranty

VWR agrees that the device shall be covered by such warranty as is offered by the

manufacturer, which VWR shall pass through to customer, without recourse to VWR. The liability of VWR under this limited warranty will not extend to any Device which is abused, altered or misused by the customer or any other persons or entities or which become defective or non-conforming through the actions or inaction of the Customer or any other persons or entities.

VWR HEREBY DISCLAIMS ALL OTHER WARRANTIES OR GURANTEES WITH RESPECT TO THE SUBJECT MATTER OF THIS AGREEMENT, WETHER STATUTORY, WRITTEN, ORAL, EXPRESS OR IMPLIED INCLUDING, WITHOUT LIMITATION, ANY WARRANTY OF MERCHANT ABILITY, SUITABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

Waiver

Customer acknowledges that VWR is not the manufacturer of any of the Devices. Customer hereby waives any and all actions, claims, suits and demands of any type whatsoever (including, without limitation, claims based on strict liability, product liability, tort or contract) against VWR, its subcontractors or agents and its or their respective employees, for personal injury, wrongful death or property damage arising out of or in any way connected with customer's receipt, storage, possession, purchase or use of the Device. Notwithstanding any other provision of this Instruction Manuel, VWR shall not have any liability for claims based upon the death or bodily injury to any person, or for the loss of, damage to, or destruction of any property so long as VWR or its subcontractors or agents were acting in compliance with customer policies, procedures and specifications of which VWR had been given notice, provided compliance with any such policy, procedure, or specification was the proximate cause of such claim.

Limitation of Liability

IF ANY DECIE WARRANTED HERE UNDER PROVES DEFECTIVE OR NON-CONFORMING,VWR'S SOLE LIABILITY AND CUSTOMER'S SOLE REMEDY HERE UNDER SHALL BE FOR VWR,TO REPLACE,AT NO COST TO CUSTOMER, ANY SUCH DEFECTIVE OR NON-CONFORMING DEVICES WITH A NON-DEFECTIVE OR CONFORMING DEVICE(AS APPLICABLE),OR,AT VWR'S OPTION,CREDIT CUSTOMER'S ACCOUNT FOR ALL AMOUNTS PAID WITH RESPECT TO THE DEFECTIVE OR NON-CONFORMING DEVICE UPON VWR'S RECEIPT OF THE DEFECTIVE OR NON-CONFORMING DEVICE.FOR PURPOSES OF THIS INSTRUCTION MANUAL, A DEFECTIVE OR NON-CONFORMING DEVICE IS DEFINED ONLY AS A DEVICE WHICH IS OUTSIDE OF THE MANUFACTURER'S DEFINED DEVICE SPECIFICATIONS AND RELEASE TOLERANCES,AND SHALL NOT INCLUDE DEVICES THAT FAIL TO MEET ANY FITNESS OF USE BY CUSTOMER OR ANY UNIQUE CUSTOMER OPERATING

CONDITIONS OR APPLICATIONS.

Compliance with Local Laws and Regulations

The customer is responsible for applying for and obtaining the necessary regulatory approvals or other authorizations necessary to run or use the Device in its local environment. VWR will not be held liable for any related omission or for not obtaining the required approval or authorization, unless any refusal is due to a defect of the Device.

Device Disposal

This Device is marked with the crossed out wheeled bin symbol to indicate that this Device must not be disposed of with unsorted waste.

Instead it's your responsibility to correctly dispose of your Device at life-cycle -end by handling it over to an authorized facility for separate collection and recycling. It's also your responsibility to decontaminate the Device in case of biological, chemical and/or radio logical contamination, so as to protect from health hazards the persons involved in the disposal and recycling of the Device.

For more information about where you can drop off your waste of Device, please contact your local dealer from whom you originally purchased this Device.

By doing so, you will help to conserve natural and environmental resources and you will ensure that your Device is recycled in a manner that protects human health.

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