



VWR® Large Capacity Tube Rotator

Instruction Manual



Version: 1 Issued:
2022-11-22



Table of Contents



Section 1 Inspection and Installation	1-1
Packing List	1-1
Connections	1-2
Structure Diagram	1-3
Section 2 Overview	2-1
Specifications	2-1
Environmental Conditions	2-2
Safety Instructions	2-3
Section 3 Operation	3-1
Operations of the Control Panel	3-1
Settings	3-2
Tray Installation and Tilt Angle Adjustment	3-4
Section 4 Safety Tips and Maintenance	4-1
Section 5 TroubleShooting	5-1
Section 6 Optional Accessories/Spare Parts	6-1
Section 7 Warranty Card	7-1

Section 1 Inspection and Installation


- 1. Inspect package and contents upon receipt of the instrument and accessories. If the package is severely damaged or if there are any missing pieces, please contact VWR Customer Service.
- 2. Unpack the instrument, ensure all parts of the instrument and accessories are not missing or damaged. Make sure to take out all the components before discarding the packaging. If any parts or accessories are missing or damaged, please contact VWR Customer Service.
- 3. Place the instrument on a level and firm surface to avoid vibration and noise.

1.1 Packing List


Table -1 Packing List

Model	76595-830	
15x 50ml rotator assembly	1	
General Power Adaptor	1	
US Plug	1	
Screw Driver	1	

1.2 Connections



1.3 Structure Diagram



Section 2 Overview

2.1 Specifications

Speed

Speed Range	5~80rpm
Display Accuracy	1rpm
Speed Accuracy	±1rpm
Display	LED

Note: Maximum speed may vary with heavy or unbalanced loads ■

Angle

Angle Adjustment Range	95°~180°
Step Angle	17°
Adjustment Method	Manual Adjustment

Load

Maximum Capacity	64×Φ18.5mm culture tubes (drum rotor)
Maximum Load (Centered on tray)	5 kg (Contain rotator drum weight)

Time

Range	0 min. ~ 99h59min.
Display	LED

Size

Overall Dimensions	440×296×466mm (17.3×11.7×18.3in)
Standard Match Dimensions	15 x 50ml Rotator Assembly (diameter Φ290mm/11.4 in)
Package Dimensions	486×381×470mm (19.2×15.0×18.5in)

Weight

Net Weight	8.1Kg(17.8lb)
Gross Weight	10.6Kg(23.3lbs)

Power Supply

Power Supply	AC100~240V, 50Hz/60Hz, 0.2A
--------------	-----------------------------

2.2 Environmental Conditions

Application Environmental Conditions: indoor use

Temperature	5° to 40°C (41° to 104 °F)
Humidity	20% to 85%,
Altitude	≤2,000 m
Voltage Fluctuation	±10% of the nominal voltage

Storage Environmental Conditions

Temperature	0° to 60°C (32° to 140°F)
Humidity	20% to 90%, non-condensing

2.3 Safety Instructions

Please read the entire instruction manual before operating the Large Capacity Tube Rotator.



WARNING: DO NOT use the Large Capacity Tube Rotator in a hazardous atmosphere or with hazardous materials for which the unit was not designed. Also, the user should be aware that the protection provided by the equipment may be impaired if accessories used are not provided or recommended by the manufacturer, or if accessories are used in a manner not specified by the manufacturer.

CAUTION ! To avoid electrical shock, completely cut off power to the unit by disconnecting the power cord from the unit or unplug from the wall outlet.

Disconnect unit from the power supply prior to maintenance and service. Any spills should be removed promptly. Biohazard spills should be cleaned using approved liquid promptly. Solvent spills are a fire hazard. Stop the unit immediately, and **DO NOT** operate until clean up is complete and vapors have dissipated.

DO NOT immerse the unit for cleaning.

DO NOT operate the unit if it shows signs of electrical or mechanical damage.


Note: The rotary plate may have a low vibration below 11 rpm, but the vibration will relieve or disappear when a load is added.

Various independent clamps are available for use in convenient combinations with the VWR High Profile Tube Rotator. It is applicable for 1.5 to 50 ml microtubes to mix liquid efficiently and gently. It has an angle range of 95° to 180° and can maintain suspension state of biological samples such as blood mixing. Applications include blood coagulation prevention, emulsion diagnosis, immunoprecipitation, and other related areas.

Section 3 Operation

This Chapter covers the control panel and its operation.

3.1 Operations of the Control Panel



The front panel contains all the controls needed to operate the unit.

1. **Start/Stop button:** Start or stop the instrument.
2. **RUN indicator light:** The light is on when the instrument is running and off when the instrument is in standby.
3. **TIME display window:** The window shows cumulative time (in continuous mode) or remaining time (in timer mode). The range of time displayed is 0 to 99 hours and 59 minutes. The accuracy is 1 minute.
4. **SPEED display window:** The window shows set speed (when the instrument is in standby) or current speed (when the instrument is running).
5. **Set Time Buttons:** UP/DOWN Arrow buttons are used to increase/decrease the set time of the instrument.
6. **Set Speed Buttons:** UP/DOWN Arrow buttons are used to increase/decrease the set speed of the instrument.

Preparation

1. Connect all the components according to the figures shown on page 1-3 of this manual. Make sure to use a grounded power outlet.
2. Press power button on the back right side of the instrument and put it into the "I" position; then the instrument is in standby and ready for use.

3.2 Settings

Time Settings

1. Continuous mode

Press the “” or “” arrow button below the TIME display window. When the number shown on the display window starts flashing, press “” arrow button to decrease the time to 00:00 and then release the button. The time setting is finished after the number shown on the display window has flashed twice.

2. Timer mode

Press the “” or “” arrow button below the TIME display window. When the number shown on the display window starts flashing, press “” or “” arrow button to increase or decrease the time value. Release the button when the time shown on the display window reaches the desired set value. The timer is set after the number shown on the display window has flashed twice.

Speed Settings

Press the “” or “” arrow button below the SPEED display window. When the number shown on the display window starts flashing, press “” or “” arrow button to increase or decrease the speed value. Release the button when the speed shown on the display window reaches the desired set value. The speed setting is finished after the number shown on the display window has flashed twice.

Note: press the “” or “” arrow button for a longer time to accelerate the setting changes.

Run and Stop

1. Continuous Mode

Press “” button and the instrument will start running with the specified settings and the RUN indicator light will be on. The TIME display window will show the cumulative time and the SPEED display window will show the current speed. Press “” button again and the instrument will slow down until it stops. The instrument will then be in standby and the two display windows will show the set values.

2. Timer Mode

Press “” button and the instrument will start running with the specified settings and the RUN indicator light will be on. The TIME display window will show the remaining time and the SPEED display window will show the current speed. Press “” button again and the instrument will slow down until it stops. The instrument will then be in standby and the two display windows will show the set values.

Finish Operation

After the operation is finished, please press the power switch at the back right side the instrument and put it into the “O” position. Unplug the instrument and store the instrument according to the storage guide.

Note:

To ensure shaking operation is smooth and steady, it may take 1 minute for the microprocessor control system to accelerate the tray to the set speed.

Alarm System

In Timer mode, the instrument will stop operation automatically and the alarm will sound when the timer goes off.


Power Recovery

If the power supply is cut off suddenly while the instrument is in operation, the unit will automatically run at the previously set parameters upon power restoration. The display window will flash. Press any button to stop flashing.

3.3 Tray Installation and Tilt Angle Adjustment

Angle Adjustment Method

1. Pull out the black bolt




Pull bolt upward


Note: If bolt does not pull out, please operate according to the following instructions

- 1) Turn motor house back and forth gently by hand and pull bolt upward at the same time.
- 2) If still cannot pull out bolt, please use vise grip plier to pull bolt upward and turn motor house at the same time.
- 3) Please contact us immediately if the problem persists.


2. Rotate the motor house gently by hand to adjust the angle.




3. Pull out bolt and turn 90 degrees and insert the top of the bolt into the long slot of the hexagon block. Push the motor house gently until the bolt fits well.




4. Angle adjustment range is shown in the figures below.



Mixing angle 95°



Mixing angle 135°



Mixing angle 180°

Single Tray Installation

1. Insert the stud on the rotation shaft into the center hole of the tray and make sure the tray is firmly touching the fixing block of the rotary plate, as shown in figure 1 (indicate on next page 3-6).
2. Screw in the fixed handle until the tray fastens, as shown in figure 2 (indicate on next page 3-6).
3. Align the tip of the screw driver with the hole at the side of the fixing block of the rotary plate. Make sure the fixing block of the rotary plate is not moving and continue screwing the fixing handle until the tray is fastened tightly (shown in figure 3, shown on the next page: 3-6).




figure 1.




figure 2.






figure 3.

4. Insert the test tube from the outside of the tube holder until it reaches the appropriate position.



Double Tray Installation

1. Insert the stud on the rotation shaft into the center hole of the tray and make sure the tray is firmly touching the fixing block of the rotary plate, as shown in figure 4.

2. Insert the special-purpose handle into the center hole and screw in the stud until the tray is fastened (shown in figure 5).




figure 4







figure 5

3. Fasten the screw on the tray with the screw driver to keep the tray is fixed well, as shown in figure 6.



Caution
Failure to load or unload the tray as listed above may result in uneven tray rotation or possible stalling of the tray

4. Insert the test tube into the hole. Make sure to place the test tubes close to the center or symmetric around the axis.



Tray, Front Facing

Note
After continued and extended operation, it is possible for the tray to exhibit uneven or erratic rotation. This does not necessarily indicate equipment malfunction and a cause for returning the unit to the manufacturer.

Loading and Unloading the Tray:
Maintaining Tray Balance
Because the rotator is a finely balanced unit, loading and unloading the tray with culture tubes or other containers needs to maintain the overall balance of the tray, especially when it is moving.

Starting from the center and working to the outside edge of the tray, insert or remove tubes that are diametrically opposed—that is, from the top, bottom, left, right. This way the container weight is distributed evenly, allowing optimum tray rotation.

Section 4 Safety Tips and Maintenance

Safety Tips

1. Use an independent power supply.
 2. Check if the local power supply voltage is suitable for use.
 3. Do not drag the power supply cable when unplugging.
 4. Do not use a non-specified power cable nor damaged cable.
 5. Service should only be performed by a qualified professional.
 6. The power supply must be unplugged under the following situations:
 - a. when the unit is moved
 - b. when the electrical cabinet or the moving component is opened
 - c. when the equipment is malfunctioning
 - d. when the equipment is not in use
 7. Wear appropriate protective equipment during operation to prevent spillage of liquid or leakage of flammable or poisonous gas.
- Reagent must be sealed inside the test tube during operation to avoid spillage while rotating.
 - Avoid vibration ensure precise operation of the instrument.
 - Balance loads if possible to improve service life of the instrument. (It is recommended to reduce load or rotation speed if the loads are unbalanced.)
 - Make sure to pull out the test tube horizontally from the holder to avoid side flip of the instrument.
 - Fixed handle must be fastened to prevent rotary plate from falling during operation.

Maintenance

- a. This instrument uses a high-accuracy motor. It is maintenance free and has a long service time, high quality, and low noise level.
- b. Surfaces can be cleaned with a mild detergent and water.

Warning: Avoid dripping detergent or water into the inside of the instrument during cleaning.

Clean Spill

If accidental spillage of liquids occurs on the surface of the instrument due to mishandling or broken containers, please shut down the instrument and clean up the liquid immediately.

If the liquid has already spilled into the unit, cut off the power supply first and immediately clean up the liquid at the surface of the instrument. Place the instrument in a ventilated and dry environment for 24 hours before reuse. If the instrument is not functioning after drying for 24 hours, please contact VWR Customer Service.









Warning: Disassembling/Assembling without a qualified professional's guidance may cause malfunctioning of the instrument.

Section 5 Troubleshooting

Please refer to the following table to troubleshoot if any malfunction occurs.
If the problem still exists, contact your local sales representative.

Error	Cause	Solution
Cannot start machine, LED display window off	Power disconnected	Connect the power
	Power adaptor failure	Replace power adaptor
	Switch off	Switch on
No shaking of the tray	Electrical malfunction	Contact VWR
	Mechanical malfunction	Contact VWR
	Over-weighted or unbalanced load	Adjust the weight and position of load, decrease rotation speed
Loud noise	Component loose	Fasten screws
	Rotary plate loose	Fasten fixed handle
	Uneven or resonant surface	Place instrument at a level and firm surface
Other	Keep record for maintenance	

Section 6 Optional Accessories/Spare Parts

Description	Cat. No.	Parameter	Figure
50 x 1.5/2.0ml rotator assembly	76595-858	50×1.5/2ml	
30 x 15ml rotator assembly	76595-860	30×15ml	
15x 50ml rotator assembly	76595-862	15×50ml	
10 x 10/15ml +20 x 5/7ml assembly	76595-864	10×15ml 20×5ml	
64 place drum for 14mm culture tubes Dia.14mm	76595-866	64×Φ14mm	
64 place drum for 18.5mm culture tubes Dia.18.5mm	76595-868	64×Φ18.5mm	
General Power Adaptor	/	AC 100~240V, 50/60HZ	
US Plug	/	250VAC,10A, 1.8m	

WARRANTY CARD

Save this card carefully please. When you need service, contact the service company with card & invoice.

Warranty Card

User Name		Telephone	
Address			
Zip Code		E-mail	
Product Name		Model	
Serial number		Purchase Date	

- Warranty Scope:
 1. Warranty period is 24 Months from purchase date.
- Warranty is not valid under the following circumstances:
 1. Warranty term is expired.
 2. Damage due to improper use, maintenance or handling.
 3. Unmatched serial numbers, or modification of serial numbers.

Customer service hotline: please contact VWR's customer service department at 1-800-932-5000.
 User: the above information should be completed by the seller.

Maintenance records

Completed by maintenance organizations, and they should help customers to paste maintenance records in the appropriate position!

The 1 time	Maintenance company's name		Maintenance Date	
	Address		Tel	
	Maintenance documents No.		Maintenance personnel to sign	
The 2 time	Maintenance company's name		Maintenance Date	
	Address		Tel	
	Maintenance documents No.		Maintenance personnel to sign	
The 3 time	Maintenance company's name		Maintenance Date	
	Address		Tel	
	Maintenance documents No.		Maintenance personnel to sign	

