TLC Explorer

Documentation System

1.52610.0001

Getting started

MERCK

Content

1	Introduction	4
1.1	Intended Use of TLC Explorer	4
1.2	Safety and Handling Instructions	5
1.2.1	Do NOT Operate the TLC Explorer if	5
1.2.2	Installation Requirements	7
1.2.3	Ultraviolet Radiation Exposure	8
1.2.4	General Notes on Handling	9
1.2.5	User Management and Information Security	10
1.3	Getting started and Operational User Manual	11
2	Get to Know Your TLC Explorer	12
2.1	TLC Explorer Front	12
2.2	TLC Explorer Back	14
2.3	TLC Explorer Baseplate	16
3	First Use	17
3.1	Unpacking and Scope of Delivery	17
3.2	Setting up for Safe Operation	18
3.3	Powering Up and a 'First Look' on Your Plate	19
3.4	Connecting Your Laptop	20
3.5	Creating a First User	23
3.6	Powering Down	25
4	Operational Guide	26
4.1	Device Preparation	26
4.2	Visual Inspection	26
4.3	Image Capture	27
5	Troubleshooting	31
5.1	Identifying Errors	31
5.2	Possible Remedies and Solutions	33
5.3	Resetting the TLC Explorer	37
5.3.1	Resetting Network connection	37
5.3.2	Factory Reset	37
6	Cleaning, Maintenance, and Decommissioning	38
6.1	Cleaning	38
6.2	Contact Technical Support	39
6.3	Ordering Accessories and Replacement Parts	40
6.4	Decommissioning, Storage, and Disposal	42
7	Appendix	44
7.1	Device Measurements and Specifications	44
7.2	Table of Hazardous Substances' Name and Concentration (only for China)	46
7.3	China SRRC-approved ID	47

1 Introduction

1.1 Intended Use of TLC Explorer

The TLC Explorer is designed to document and evaluate thin-layer chromatographic plates. It offers a safe, fast and flexible way to visually inspect and digitally record images of TLC plates after their chromatographic development workflow in the laboratory. For detecting different substance classes, the device contains three different LED light sources which illuminate the plates from above in direct-light mode:

- **White light** is used for the detection of visible, colored substances on thin-layer chromatograms.
- **UV-A** (366 nm) is typically used to stimulate fluorescence. Under this illumination, fluorescent substances and objects appear as bright, different colored spots on a dark background.

It can also be used for plates which are pre-coated with the special fluorescence indicator F366. The substances, in this case, appear as dark or colored spots on a bright background. Other fluorescence indicators, like F254, do not fluoresce under this light and therefore do not interfere with the measurement.

- **UV-C** (254 nm) illumination is used almost exclusively for thin-layer chromatograms. Substances which absorb this light appear as dark spots on a bright background, if the plate is pre-coated with an appropriate F254 fluorescent indicator.

Recorded plate images can be evaluated to generate reports or export the data for further processing with other systems.

The TLC Explorer was developed for use in the laboratory and to be operated by trained and skilled users to perform thin-layer chromatography. Only fully dried TLC plates with a maximum size of 20 x 20 cm should be used.

Please note, that any other use than the intended one described above invalidates claims under warranty.

1.2 Safety and Handling Instructions

Operational safety has been a primary objective when developing the TLC Explorer. However, it can only be guaranteed if the device is operated strictly according to the instructions given in this manual. Furthermore, any operation conflicting with the safety instructions given in this manual may invalidate claims under warranty. All personnel responsible for installation, use and maintenance must therefore read the operating manual (see section 1.3) carefully before working with the TLC Explorer and keep it at hand for later reference.

1.2.1 Do NOT Operate the TLC Explorer if ...

Before connecting the TLC Explorer to a power supply, ensure none of the following applies:

- The device has been stored under adverse conditions over a lengthy period of time.
- The device has been damaged in transport or by incorrect handling (e.g. dropped during transport or contact with liquids).
- Any part of housing, front unit, back cover, drawer or Baseplate is wet, damaged or deformed.
- The UV protection inspection window is missing, damaged or deformed.
- The device does not stand on a dry, flat and stable surface.
- The drawer cannot be extended fully or does not close correctly.
- The back cover is not mounted properly, and all four screws fastened.
- Any connector (USB-C, USB-A, Ethernet and WiFi Antenna) or cable is wet, damaged or deformed.

- The power supply is none of the recommended two adapters (USB-C power adapter DELL HA65NM170 or Power bank DELL PW7018LC).
- The device is significantly colder than its environment, such that humidity could condensate and damage electric circuits (please wait until device reaches room temperature).
- The device environment temperature lies outside the specified 15 - 40°C or the relative humidity outside the specified 20 - 80%.
- The device is or has been exposed to acidic, alkaline, solvent vapors or other gasses that may cause corrosion.
- Liquids or other materials have been spilled on or in the device. If so, remove them immediately. Make sure that the device is not damaged and completely dry before using it again.

If you expect any of the issues above apply to your device, do not connect, power up or operate. For questions on how to proceed, please contact the supplier of your TLC Explorer or review the TLC Explorer service site using the link in Figure 1.1:

Figure 1.1: Find the TLC Explorer service page at:
www.sigmaaldrich.com/TLCservice

1.2.2 Installation Requirements

For a safe and reliable operation, please follow the following guidelines when setting up the device in your laboratory:

- The device must be installed on a flat and dry lab bench suitable to support the weight of about 10 kg of the device.
- Do not expose the instrument to any strong vibration or shock.
- Protect the instrument from excessive dust.
- Locate the TLC Explorer, cables, and, if applicable, the power bank, at a safe distance from heating sources (e.g. gas burners, electric heaters or ovens), chemical substances (e.g. solvents) or other devices which could cause damage to them.
- To power the device, use only the USB-C power adapter DELL HA65NM170 or Power bank DELL PW7018LC, and connect them to the USB-C port at the back of the device.
- Connect the power source to a grounded outlet located within 2 meters of the device, following the regulations for the use of electronic devices in your laboratory.
- The line voltage at the user location must meet the specifications stated on the power adapter.
- Do not place the instrument near equipment that generates intense magnetic fields such as electric welding equipment, high frequency furnaces, pole transformers, etc..
- Connect the instrument to power lines that are free from sudden changes or voltage fluctuations.
- In case power motor driven equipment (such as a stirrer or shaker) must be connected to the same power line as your TLC Explorer, ensure that a noise reduction unit is installed.

1

2

3

4

5

6

7

- 1 • Ensure that cables are arranged in a way that they neither pose a risk to the staff working with, and around, the device, nor that they are accidentally drawn from their respective sockets.
- 2 • There must be enough space in front of the device to open the drawer completely and handle the Baseplate in a safe and convenient way.
- 3 • Place the device on a height that allows users to conveniently look through the inspection window. Avoid reflections of external light sources on the glass front that disturb the view through the inspection window.
- 3 • The device should only be used in an environment and operated in a manner that meets all official lab safety requirements.

4 1.2.3 Ultraviolet Radiation Exposure

5 Ultraviolet radiation, used in TLC plate evaluation to stimulate fluorescence, can cause damage to eyes and skin. Safety regulations in most countries therefore demand that eyes and skin are protected from hazardous UV radiation. The TLC Explorer incorporates several safety measures to ensure that no UV-A or UV-C light is leaking from the device when used as described in this guide. Under no circumstances try to circumvent these safety measures!

- 6 • Only look at TLC plates under UV illumination through the UV-protection inspection window.
- Ensure the inspection window is not damaged or deformed.
- Before activating UV illumination, check that the drawer and back cover are correctly closed and not damaged.
- 7 • Never operate the illumination unit outside the TLC Explorer.
- Follow general laboratory safety regulations for working with UV light.

1.2.4 General Notes on Handling

- Ensure the device has at least 30 seconds available to completely shut down before unplugging the power supply – similar to your desktop computer.
- Only TLC plates without solvent residue can be documented – please ensure the TLC plate is completely dry after the previous chromatographic development or derivatization step.
- Do not place objects higher than 8 mm on the Baseplate, as it might damage the device when closing the drawer and would be out of focus for image capture.
- Never lift or carry the device by holding the drawer handle. Instead, hold it on the bottom and grab the rubber feet to prevent it slipping out of your hands.
- Close the drawer whenever the device is not in use to avoid damages and contamination.
- When the TLC Explorer is being transported, the TLC Explorer must be closed, the Baseplate removed. No TLC plate should be inserted.
- Only attach FAT32 formatted USB sticks to USB-A ports for backups and updates. Do not use external HDDs nor any other media connected via a cable.
- Except for the replacement of the illumination unit, the TLC Explorer shall only be opened, adjusted, or repaired by specialist personnel authorized by the manufacturer. Noncompliance invalidates any claim under warranty.
- Retain the original packaging, including the inner packaging, to protect the instrument against physical impact if it has to be transported, especially for warranty claims. Please note, that damage caused by improper transport voids all warranty claims.

1

2

3

4

5

6

7

1.2.5 User Management and Information Security

It is assumed that the devices are placed in laboratories where physical access is limited to authorized staff only. To give priority to a smooth and efficient laboratory workflow, the TLC Explorer, therefore, features no means to restrict access to specific users or protect data with passwords or encryption. User accounts mainly serve to structure recorded data. There is no administration account with special privileges.

Any person with physical access or network connection to the device can

- create new users
- delete or login as an existing user
- access and delete all stored data
- initiate a restart or shut down of the instrument
- change device settings (incl. connectivity parameters)
- conduct a factory reset.

Depending on how you choose to integrate the TLC Explorer in your IT environment (see Operational User Manual), a large group of users may have access to your device(s). If the TLC Explorer is connected to a LAN or a WiFi network, usually everyone with access to that network can also open the user interface and login as an existing or new user (see section 3.4). If confidentiality is a concern, copy the respective data to your desktop or tablet, and then remove it from the TLC Explorer.

- Plan and test carefully the TLC Explorer connectivity settings to ensure only authorized and trained personnel have access.
- Download or copy your data to a memory stick and then save it at an external storage location. Do not use the TLC Explorer for permanent data storage.

- Only apply verified updates from the manufacturer.
- Only attach USB memory sticks (flash drives) you are certain do not contain malicious software.

1.3 Getting started and Operational User Manual

This 'Getting started' guide is only meant to help in setting up the TLC Explorer for operation and getting started with image capture. Details on processing images of TLC plates and configurations are provided in the Operational User Manual, which can be downloaded by scanning the QR code in Fig. 1.2 or on the official TLC Explorer support site <https://sigmaaldrich.com/TLCservice>. Additionally, the Operational User Manual can be opened from the TLC Explorer's software by clicking on 'Help' in the main menu (three horizontal bars in the upper right corner).

Figure 1.2: Download the Operational User Manual at: [www.sigmaaldrich.com/TLCservice](https://sigmaaldrich.com/TLCservice)

1

2

3

4

5

6

7

2 Get to Know Your TLC Explorer

2.1 TLC Explorer Front

Figure 2.1: Identifying elements on the device front

ID	Item	Description
F1	Power button	Place a finger for at least two seconds on icon to start or switch off the device.
F2	Illumination 'Vis'	Switch on/off white light source.
F3	Illumination '366'	Switch on/off UV-A light source.
F4	Illumination '254'	Switch on/off UV-C light source.
F5	Drawer	Open drawer to place TLC plate on Baseplate.
F6	Drawer open indicator	A red light indicates that the drawer is not closed correctly. While the red light is on, no image can be captured, and UV light sources are deactivated.
F7	Inspection window	Allows to inspect illuminated TLC plate by eye (protecting the eye from UV radiation).
F8	Shutter slider	Gently push down the slider to open inspection window.

1

2

3

4

5

6

7

2.2 TLC Explorer Back


Figure 2.2: Identifying elements on the device back

ID	Item	Description
B1	USB-C port	This USB-C port solely serves to connect the power supply and features no data transfer.
B2	USB-A port	Port to connect memory sticks (flash drive) for data export, backups or software updates.
B3	USB-A port	Port to connect memory sticks (flash drive) for data export, backups or software updates.
B4	Ethernet port	RJ45 port to connect to LAN.
B5	WiFi antenna connector	To mount the WiFi antenna, straighten up the joint at the antenna's base and screw it clockwise into the connector.
B6	Serial number	A sticker with the TLC Explorer's serial number. The QR code below the serial number includes a link to the full device label, which is physically situated on the bottom of the device.
B7	Back cover	Only open to clean device inside or exchange the illumination unit (for details see Operation User Manual).
B8	Illumination unit	Installed behind back cover, the exchangeable illumination unit houses the LED light sources for the different illumination types.

2.3 TLC Explorer Baseplate


Figure 2.3: Identifying elements on and related to the Baseplate

ID	Item	Description
P1	Plate area	Place TLC plates directly on the black coating.
P2	Frame	No object on the plate area should be higher than the frame height when inserting the Baseplate with the drawer.
P3	Alignment pin	Before closing the drawer, ensure that the pin has entered the socket on the drawer, so that the Baseplate lies flat on the drawer.
P4	Alignment socket	

3 First Use

3.1 Unpacking and Scope of Delivery

Your TLC Explorer is built and inspected according to the relevant guidelines and norms for electronic instruments (for details see Operational User Manual). It left the factory in a safe and secure technical condition. In case you note any damage to the packaging that could have impacted the device, document these and contact your supplier. Do not operate the TLC Explorer if you suspect any of the criteria listed in section 1.2.1 to apply.

It is strongly recommended to retain the packaging material for later transport and/or storage of the device and especially for warranty claims.

After opening the cardboard box, first take out all accessories, including the Baseplate, and the upper protective carton. Ensure all loose parts are removed before proceeding. Turn the box, so that the glass front is facing you, and prepare a suitable space to place the device, for when you have taken it out of the box. Only then lift the TLC Explorer (which has a weight of about 10 kg) by holding it on the bottom under both sides. Remove all transportation fixtures – like the adhesive tape fixing the drawer – and protective foils.

Verify that your delivery includes all the following:

- One TLC Explorer with installed illumination unit (to check illumination unit, see section 3.3 below)
- One Baseplate
- One WiFi antenna
- One USB-C power adapter with eight power cables for different socket types
- One dust cover
- Two stickers with two QR codes – used to connect to the device more easily (see Fig. 3.1)
- This printed 'Getting started' guide and safety instructions

1

2

3

4

5

6

7

3.2 Setting Up for Safe Operation

Please carefully read section 1.2.2 about installation requirements before selecting a suitable location for your device. Check that the drawer **F5** can be completely opened, the Baseplate inserted, and TLC plates be seen through the inspection window **F7** in a safe, convenient and ergonomic way for all planned users. Especially consider different body heights and possibly disruptive light reflections on the device's glass front.

Keep an easily accessible copy of this operational user guide near the TLC Explorer, e.g. as a digital document on the tablet or desktop computer you plan to use to control the device.

3.3 Powering Up and a 'First Look' on Your Plate

Choose a power cord **A3** that fits the power outlet, connect it to both the power adapter **A2** and the power grid, and then plug the USB-C connector into the port **B1**. On the front side, you should see the power icon **F1** light up. Gently hold a finger on the power icon for about two seconds, until the icon begins to flash, indicating that the device is booting. After about 30 seconds the icon stops flashing, shines brighter, and the device is ready to operate.

Open the drawer **F5** and insert the Baseplate. The side with the wide edge and logo should be facing outward, i.e. be next to the drawer handle. Verify that the alignment pin **P3** has entered the socket **P4**, and that the Baseplate lies flat on the drawer. To close the drawer, push against the middle of the handle until you feel a slight resistance. The drawer is designed to close the last few centimeters on its own. A special damper system takes care to correctly close the drawer without shifting the plates.

Once the red drawer open light **F6** is off, briefly hold your finger against the circle next to the label 'Vis' **F2** to switch on the white light source inside. Pulling down the shutter slider **F8** and looking through the inspection window **F7**, you should see the illuminated Baseplate.

Switch off the light (touching **F2**), open the drawer and place a TLC plate on the plate area **P1**. Close the drawer again and observe the plate through the inspection window under the featured illumination types (touching **F2**, **F3** or **F4**).

In case your TLC Explorer acts differently than described above, refer to section 5 to identify and resolve the issue before operating the device.

1

2

3

4

5

6

7

3.4 Connecting Your Laptop

There are different ways to connect various types of user devices (laptops, desktop computer, tablets or even smart phones) to the TLC Explorer, which can be configured as described in the Operational User Manual. A recommended basic approach is to connect your laptop to the TLC Explorer's own WiFi Network. This approach is quickly set up and allows to download images and data you generate with just a few clicks directly to the device you are using.

Figure 3.1: Easily connect by scanning these two QR codes on the label of your device

On a laptop (or desktop PC featuring WiFi) open the network settings and search for a network name (SSID) starting with 'tlc-explorer' followed by the serial number **B6**. Choose 'WPA/WPA2' as the security setting and enter the password shown on the label for Step 1 (see Fig. 3.1) to establish the connection. It is even simpler to connect tablets or smart phones by scanning the QR code 1 on the sticker provided in the package. You may see a warning that the chosen WiFi has no internet connection – this is correct and by design.

Once the WiFi connection is established, the user interface can be opened in your web browser. Although all common, up-to-date browsers can be used, it is recommended to use Chrome or Firefox for the best user experience.

To open the TLC Explorer user interface, you may either scan the QR code 2 on the sticker or enter “http://tlc” in the browser address bar.

In case the browser has difficulties to find “http://tlc” (e.g. because the device is connected to the internet via a mobile data network in parallel), you can also enter “http://192.168.200.1” instead. Some browsers enforce encrypted connections and may automatically change the entered address to begin with “https://”, to then produce a warning the site may not be secure. The browser cannot know that you are accessing the TLC Explorer in a protected local network, posing far less risk than opening a site from the internet. You may either ignore the warning and proceed, or disable “https-only” in the browser settings, and then manually enter again “http://tlc” in the browser address bar.

1

2

3

4

5

6

7

3 First Use - 3.4 Connecting Your Laptop

You should see a welcome screen similar to the following (Fig. 3.2):


Figure 3.2: First opening of TLC Explorer user interface

Setting a bookmark and/or selecting the menu point 'Add to Home Screen' from the browser menu (often marked as three dots in a column at the upper right browser corner) will significantly simplify opening the site in future.

3.5 Creating a First User

Begin personalizing your TLC Explorer by selecting a language and entering a first username. Since usernames are internally used to name files and directories, not all special characters or symbols are allowed. In case you entered an invalid character, the line under the field will turn red and a warning message will appear underneath.

The full name, which you can enter in the next field, is less restrictive and optional. It may be used to identify the user e.g. as the author of reports.


Figure 3.3: Adding a first user

Keep 'Device and Network Settings' for now at their default values. They can be modified any time later at the 'Settings' page.

Expand and carefully read the legal terms before clicking in the box to the left of 'Accept Legal Terms'. Your agreement to these terms is needed for you to operate the instrument and software.

After you have clicked on 'Add User', the initial setup is complete and the device ready to be used. You will be automatically forwarded to the standard welcome screen seen at figure 3.4. Select the user you created and click 'Login' to start capturing images of your TLC plates.


Figure 3.4: Welcome screen with user selection

3.6 Powering Down

Like any desktop computer, the TLC Explorer needs to be properly shut down before unplugging the power supply. Failing to do so may damage the device.

In order to shut down the device correctly, click on the main menu icon (three horizontal bars in the upper right corner), and select 'Shutdown' from the menu that appears. From the dialog depicted in figure 3.5 you may choose to shut down the device, restart (i.e. reboot) the device, or just restart the software application. You may also touch the power icon **F1** for about 2 seconds.

If you chose to restart the device or application, just wait for about a minute until the welcome screen shown in figure 3.4 reappears, and the device is ready again.


Figure 3.5: Shutdown or restart device

4 Operational Guide

Whether you just want to take a quick look at a TLC plate, capture and export an image of a plate (maybe with some annotations) for some document, or need a quantitative analysis with a complete report – the TLC Explorer lets you do just that.

However, most of the software's features to process and document TLC plates are out of the scope of this 'Getting started' guide. A detailed description of the software can be found in the Operational User Manual.

4.1 Device Preparation

After connecting the TLC Explorer to the power grid (or power bank), a dim white light under the power icon **F1** shows that the device is ready to boot. Place your finger on the icon for about 2 seconds and allow the instrument some 30 seconds – while the icon is flashing – to get ready for operation.

Place a TLC plate on the Baseplate and close the drawer **F5**. Once the red open drawer indicator **F6** turns off, you will note a brief white flash of the circle next to 'Vis' **F2** – the system has just taken a snapshot to show a preview of your plate in the software.

4.2 Visual Inspection

Prerequisites:

- The TLC reader is on; the power icon **F1** is lit up brightly (not flashing nor dimmed).
- The drawer is closed correctly and the red open drawer indicator **F6** is off.
- No user is capturing images at the same time (see also section 4.3).

Select the illumination type you would like to use by briefly holding your finger on the respective button (**F2**, **F3**, or **F4**) until the circle lights up. Then gently push down the shutter slider **F8** to open the inspection window **F7**. Depending on your observation angle, you can push the shutter slider fully or only partly down to see the relevant part of the TLC plate. You can switch to another illumination type by touching the respective button or turn off the illumination by pressing the same button again.

4.3 Image Capture

Prerequisites:

- The TLC reader is on; the power icon **F1** is lit up brightly (not flashing nor dimmed).
- Drawer and inspection window are closed correctly, and the red open drawer indicator **F6** is off.
- One or more TLC plates with a minimum size of 2.5 x 7.5 cm and a maximum size of 20 x 20 cm lie flat on the plate area **P1**.
- Only one user at a time can capture images. While other users may load and process Image Sets in parallel, competing attempts to capture images will cause error messages.

Assuming you have connected a tablet or mobile phone as recommended in section 3.4, you can just click on the desktop icon, or bookmark you saved in your browser, to open the welcome screen seen in figure 3.5. Select your user-name and click on 'Login'.

The page that opens after each login (see Fig. 4.1) lets you either capture new images or load an existing Image Set. At the left – or, for small screen width, on the top – you will see the low-resolution image of the Baseplate that was captured when you closed the drawer. You can refresh it by simply clicking the image.


Figure 4.1: Capture images or load existing Image Set

Select the option 'Capture new Image Set' from the drop-down menu and enter a name for the new Image Set in the text field. Because this name will be used as part of image filenames, not all characters are allowed. Furthermore, the name needs to differ from other Image Set names stored for the user. In case an Image Set name is invalid, the software informs you about the reasons.

Checkboxes let you choose which illumination types you wish to use for the image capture. Expanding the 'Use manual Exposure' section allows you to set the exposure time in milliseconds for each of the chosen illumination types. Click on the camera icon to let the system take a low-resolution image with the selected exposure time and shows it as a preview. With the 'Auto' button, you request the system to take a series of images and determine automatically the optimal exposure time for your plate. Please be aware that this process might take several seconds, and only proceed when you see the preview updated.

With a click on 'Capture Now', the instrument starts capturing images. You can see which image is currently being taken, both in the text above the progress bar and on the respective illumination button (**F2** - **F4**) on the device lighting up.

The resulting images of all illumination types are stored together as an 'Image Set'. The system then searches for plates in the images. As soon as all plates are detected, rotated (aligning their bottom edges horizontally), and extracted as separate images for further processing, the user is automatically forwarded to the next step for plate editing.

In case you only intend to use the captured images (or the extracted plate images), you can bypass all further processing by using the file manager to download the images you have just taken from the file manager, which can be opened via the main menu (three horizontal bars in the upper right corner). Alternatively, you can simply click the 'Download image' button at the bottom of the 'Plate' page to save the image currently shown in the preview area. The Image Set can still be loaded and processed any time later on.

1

2

3

4

5

6

7


Figure 4.2: Manually set exposure time for selected illumination types.

5 Troubleshooting

5.1 Identifying Errors

An unexpected behavior or error message is sometimes the consequence of a previous unnoticed error. To speed up the search to find an error’s origin and potential solutions, the chart in Fig. 5.1 lists a sequence of checks. Once the issue has been identified, please follow the possible remedies referenced by number to the right of the issue.

Please DO NOT unplug the system from the power supply without properly shutting it down (see section 3.6). Doing so may seriously harm the computer system inside the TLC Explorer and cause additional problems. A safe alternative to shut down and reboot the device is to press and hold the power icon **F1**.

1

2

3

4

5

6

7

5 Troubleshooting - 5.1 Identifying Errors


Figure 5.1: Identifying issues and how to address them

5.2 Possible Remedies and Solutions

- 1) To reboot the device, hold a finger on the power icon **F1**, until the icon starts flashing (indicating the device is shutting down). After the icon has stopped flashing and remains in a dimmed state, wait a few seconds and press the icon again to restart the device.
- 2) If the power icon does not light up, the device may not be connected to a working power supply. A simple way to check the power source is to connect a tablet or mobile phone to the USB-C connector that is plugged in the USB-C port at the device back. However, this test only applies to the USB-C power adapters recommended for the TLC Explorer (see Operational User Manual) – third party power adapter and cables may charge a phone but have specifications that are incompatible with the TLC Explorer. In case the power adapter and cable work as desired with other devices, but not the TLC Explorer, you may need to contact the reseller of your device for technical support (see section 6.2).
- 3) In case the drawer does not close completely, verify that the alignment pin **P3** underneath the Baseplate has entered the socket **P4** on the drawer, so that the Baseplate lies flat on the drawer. There should be no objects higher than 8 mm on the Baseplate. If the drawer seems blocked by an object inside the device, you can open the back cover **B7** and remove the obstacle.
- 4) The red drawer open light **F6** is on whenever sensors indicate that the drawer is not closed correctly. While the drawer is open, light sources cannot be switched on and no images can be taken. Check if the drawer really is correctly closed. In rare cases, a software or sensor issue might cause the systems to incorrectly detect the drawer as open. If you suspect this to be the case, reboot the device.

1**2****3****4****5****6****7**

- 1 5) For safety reasons, the TLC Explorer only allows the switching on of light sources if the drawer is correctly closed and the back cover **B7** correctly installed. Also, while another user is capturing images, manual light switches **F2**, **F3**, **F4** are deactivated.
- 2 6) If one or all light sources are not working, or the illumination appears to be dark on one side, the illumination unit might need to be replaced. The procedure takes only a few minutes and is explained in detail in the Operational User Manual. In case a replaced illumination unit remains dark after the replacement procedure, verify that the unit's power connector is plugged in (see Operational User Manual), and the back cover **B7** has been correctly closed.
- 3 7) Issues with the WiFi configuration may also be resolved by connecting to the device using an Ethernet connection. You can then access the settings page and correct the configuration via the Ethernet connection.
- 4 8) To reset all network settings to their default values, please follow the instructions given below in section 5.3.1. You should then be able to connect to the TLC Explorer's own WiFi network by following the steps described in section 3.4.
- 5 9) The network administrator may have installed security measures on the network configuration, or your device settings, preventing your device(s) from connecting to the TLC Explorer. If in doubt, please ask your network administrator.
- 6 10) Some mobile devices will give preference to network connections that have internet access, which the TLC Explorers own WiFi network does not have. It might, therefore, be necessary to open your mobile device's WiFi settings and disable the 'Auto-connect' for saved networks that have internet access – those which are also available apart from the TLC Explorer.
- 7

- 11) Please note that each TLC Explorer broadcasts a WiFi signal with a different name (usually including its serial number) which is possibly protected by a different password. In case your laboratory works with multiple TLC Explorer devices, it is necessary to use the QR code on the sticker for the specific device to connect (as long as the default network settings have not been changed).
- 12) The TLC Explorer tries to inform all devices in the same network that it can be reached with the URL "http://" followed by the device name (see section 3.4). However, some networks may not support this feature and require users to address the device by its IP, e.g. "http://192.168.200.1". This may also be necessary, if mobile devices in parallel feature an internet connection via mobile data and try to find the URL "http://" followed by the device name in the world wide web.
- 13) Some browsers enforce encrypted connections and may automatically change the entered address to begin with "https://", to then produce a warning that the site may not be secure. The browser cannot know that you are accessing the TLC Explorer in a protected local network, posing far less risk than opening a site from the internet. Alternatively, you may disable "https-only" in the browser settings (see also section 3.4).
- 14) The 'Capture Now' button is only enabled if a unique and valid Image Set name has been entered and the system is not busy with image capture. An extended delay may occur if a user has pressed an 'Auto' button multiple times in the 'Use manual Exposure' section to find the optimal exposure time. This process may take a while, especially for UV-light sources, because a series of images has to be captured and processed.

1

2

3

4

5

6

7

- 1 15) To ensure no stray light is impacting the image quality – e.g. due to an open inspection window **F7** –, the camera checks if the Baseplate is completely dark with no light source on before starting the image capture process. This check can be disabled on the 'Settings' page, though it is recommended to enable it.
- 2 16) Only one user at a time can capture images. While other users may load and process Image Sets in parallel, competing attempts to capture images will cause error messages.
- 3 17) If the available storage space is insufficient to save captured images, an error message will appear. Please transfer all data not currently needed for processing to other devices and delete their copy on the TLC Explorer. At the bottom right corner of the file manager (reached via the main menu in the upper right corner), you can see the 'Remaining free space'. Please have at least 250 MB free to ensure the system runs smoothly.
- 4 18) The software searches for TLC plates as bright rectangular shapes with a minimum size. To get multiple plates detected automatically, place them a few millimeters apart, so that the device detects at least 2 mm of the black coating around each plate.
- 5 19) If the plate appears correctly illuminated through the inspection window **F7**, but in homogeneously illuminated in the plate image, an incorrect illumination correction might have been applied. It is possible to deactivate the correction step, or to recalibrate the calculation. For details, please refer to the Operational User Manual.
- 6 20) As an ultima ratio, all user data, settings and software updates can be deleted by performing a factory reset, as described below in section 5.3.2.
- 7 21) Contact the reseller of your TLC Explorer for technical support.

5.3 Resetting the TLC Explorer

5.3.1 Resetting Network Connection

In case some changes to network settings, or a forgotten WiFi password, hinder you from connecting to the TLC Explorer, the network settings can always be reset to the factory default.

To reset the network connection, touch and hold any **two** (not all three!) illumination buttons (**F2**, **F3**) while the TLC Explorer is starting, i.e. as long as the power icon **F1** is flashing. You should then be able to connect to the device using the QR codes shown in Fig. 3.1.

5.3.2 Factory Reset

Initiating a factory reset should be done with care, as it **erases all users with their data as well as all settings**, including network settings. A factory reset also re-installs the software version with which the device was delivered and erases all other software versions from the device.

To start a factory reset, touch and hold all three illumination buttons (**F2**, **F3**, **F4**) while the TLC Explorer is starting, i.e. as long as the power icon **F1** is flashing.

1

2

3

4

5

6

7

6 Cleaning, Maintenance, and Decommissioning

6.1 Cleaning

It is recommended to clean the TLC Explorer – and especially the Baseplate – with a damp lint-free tissue on a regular basis. To remove more resilient dirt, the cleaning tissue may also be used damp with a mixture of alcohol and water at the ratio of 1:1.

Before starting the cleaning process, shut down the device (see section 3.6) and unplug the USB-C power supply **B1**. Only then remove the Baseplate from the drawer to wipe it clean on all sides. In case silica particles from TLC plates, or dust, are on the drawer, or in the housing, both drawer and housing can be cleaned in the same manner – be sure not to touch the illumination unit. Alternatively, the back cover **B7** can be opened, and the dirt inside the device removed by carefully using a hand vacuum cleaner at low power and a tissue.

Please note: Spilled liquid or other material should be removed immediately (see also section 1.2.1)!

6.2 Contact Technical Support

For questions about warranty, warranty related issues, or questions specific to your device, please contact the reseller of your TLC Explorer. Similarly, if in doubt about how to proceed when facing one of the situations listed under 'Do NOT Operate the TLC Explorer if ...' in section 1.2.1, please contact the reseller of your device or a local supplier via www.sigmaaldrich.com.

For generic technical questions and software updates, please visit the official TLC Explorer support site using the link provided in Fig. 6.1 below.

Figure 6.1: Find the TLC Explorer service page at:
www.sigmaaldrich.com/TLCservice

1

2

3

4

5

6

7

6.3 Ordering Accessories and Replacement Parts

Currently, the TLC Explorer product portfolio includes the following products which can be ordered separately as additional equipment or replacement parts:

Cat. No.	Name	Description
1.52610.0001	TLC Explorer	TLC Documentation system, including the parts listed in section 3.1 as scope of delivery.
1.52613.0001	Illumination Unit	Replacement part which is only needed in case a light source in your TLC Explorer has broken. The unit houses LEDs and optics for all illumination types, i.e. visible light as well as UV light, and can only be exchanged as a whole.
1.52612.0001	Baseplate	An additional Baseplate can simplify and speed up the workflow, since one Baseplate can be prepared while a second one is recorded in the TLC Explorer (see section 2.3).

Cat. No.	Name	Description
1.52611.0001	Power bank	Powers the TLC Explorer independent of a power grid – e.g. for the use away from a mains power source.
1.52618.0001	Dust cover	Replacement cover to protect your TLC Explorer against dust.

Common TLC consumables like TLC plates, solvents and stains, as well as tools used in the TLC process (e.g. developing chambers) are available on the website of Sigma Aldrich at www.sigmaaldrich.com.

Figure 6.2: Find TLC consumables and tools at: www.sigmaaldrich.com

6.4 Decommissioning, Storage, and Disposal

1

If you decide to decommission your TLC Explorer, you can either store it for later use or return it to the reseller or distributor of your TLC Explorer.

2

There are no batteries in the TLC Explorer. You, therefore, do not need to remove any batteries before storing, shipping or returning the device.

3

However, there might still be sensitive data on the internal storage. It is strongly recommended to initiate a factory reset (see section 5.3.2) to remove any user data before decommissioning your TLC Explorer.

4

The device should also be cleaned of all visible residue and encrusted material. In case the device was in contact with hazardous chemicals or infectious/bio-hazardous materials, it must be decontaminated accordingly. Please note that the reseller can refuse to collect a decommissioned device if it has not been correctly decontaminated.

5

To store the device, remove any TLC plates and check that no USB stick is plugged in. Do not leave the Baseplate or any accessories inside the housing, but rather pack them separately next to the housing. Ideally, use the original packaging, including the inner packaging, to protect the device while in storage. Ensure that the storage place is dry, with a relative humidity between 20 - 80%, and the temperature remains within a range of 15 - 40 degrees Celsius.

6

In accordance with European Union directive 2012/19/UE on the management of waste electrical and electronic equipment (WEEE), it is strictly forbidden to dispose the TLC Explorer in unsorted municipal or household waste. Instead, it should either be brought to an official collection and recycling center for electrical and electronic waste or returned to the reseller of your TLC Explorer.

7

For further information on take-back programs offered in your region, please visit www.sigmaaldrich.com/EN/en/life-science/quality-and-regulatory-management/recycling.


Figure 6.3: Only dispose device at official collection and recycling center or return it to your reseller.

1

2

3

4

5

6

7

7 Appendix

7.1 Device Measurements and Specifications

1	(External) Power supply	USB-C Power Supply Unit / 65W
2	Power supply	External Power Supply Unit (HA65NM170) is provided with separate Power supply cord(s) (1.8 m long) fitting type B, G, N, I, D, E/F plugs (the plug must comply with local regulations)
3	Power requirements	100 V – 240 V, 50/60 Hz for the external power supply unit
4	Power consumption	Standard working condition: 18.8 W; standby mode: 0.83 W
5	Measuring technology	Documentation system for TLC plates by using video densitometric measurement.
6	Wavelength range	The TLC Explorer provides the standard illumination for detecting TLC plates (see also section 12.2): • Visible light (VIS) • UV-C - WL: 254nm • UV-A – WL: 366nm
7	Camera	RGB sensor with 3280 × 2464 active pixels; one pixel captures an area of approx. 85 × 85 µm on the TLC plate
	Communication interfaces	• USB: 2 × USB-A (for the direct connection of a USB Stick) • Ethernet: LAN • WLAN IEEE 802.11 b/g/n (2.4 GHz)
	Protection (IP) class	IP2X (for the main enclosure)
	Ambient Condition: Temperature	• The device can be operated in a temperature range from 15°C to 40°C. • The temperature for storage and transportation should be between 0°C and 40°C.

7 Appendix - 7.1 Device Measurement and Specifications

Ambient Condition: Relative humidity	• The device can be operated in a humidity range from 20%RH to 80%RH, non-condensing conditions shall be ensured. • The humidity for storage and transportation should be in a range from 15%RH to 95%RH, non-condensing conditions shall be ensured.
Ambient Condition: Altitude	< 2000 m
Dimensions	374 x 312 x 290 mm length x width x height
TLC plates sizes	Usage of TLC plates with a size of 20 cm x 20 cm and smaller
Weight	approx. 10.4 kg
Instrument safety	Complies with Low voltage Directive 2014/35/EU IEC 61010-1 ANSI/UL EN 61010-1 CSA-C22.2 No. 61010-1-12
EMC	Complies with EMC Directive 2014/30/EU EN/IEC 61326-1 / EN 301 489-1 FCC
RF	Complies with RED 2014/53/EU EN 300 328, FCC Contains: FCC ID 2ABCB-RPICM4 IC:20953-RPICM4

1

2

3

4

5

6

7

7.2 Table of Hazardous Substances' Name and Concentration (only for China)

Component Name	Hazardous substances' name					
	(Pb)	(Hg)	(Cd)	(Cr ⁶⁺)	(PBB)	(PBDE)
Housing	O	O	O	O	O	O
Optics	O	O	O	O	O	O
Electronics	O	O	O	O	O	O

O: Indicates that this toxic or hazardous substance contained in all of the homogeneous materials for this part is below the limit requirement in SJ / T11363-2006.

X: Indicates that this toxic or hazardous substance contained in at least one of the homogeneous materials used for this part is above the limit requirement in SJ / T11363-2006.

- Data listed in the table represents best information available at the time of publication
- Product to be operated under Environmental Conditions as defined in the Product Manual to maintain the declared EFUP
- Consumables or parts may have their own label with an EFUP value less than the system. Periodic replacement of those consumables or parts to maintain the declared EFUP is done in accordance with the declared EFUP is done in accordance with the Product Maintenance Procedures

7.3 China SRRC-approved ID

CMITT ID: 2023DJ15480

1
2
3
4
5
6
7

Supelco®

Analytical Products

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

Merck Life Science KGaA, 64271 Darmstadt, Germany,
Tel. +49(0)6151 72-2440

www.sigmaaldrich.com

© 2023 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.
Merck and Supelco are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

The Merck logo, consisting of the word "MERCK" in a bold, green, sans-serif font.

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