

GENERAL DESCRIPTION

J.T.Baker® Cell Lysis Solution Biotech Reagent is a cGMP product supplied as a ready to use reagent containing a proprietary blend of biodegradable nonionic surfactants to solubilize the cell membrane lipid bilayer to release enclosed viral particles. The innovative formulation provides viral particle protection and higher recovery during the lysis step. J.T.Baker® Cell Lysis Solution is not listed on REACH SVHC (substances of very high concern) list, free of alkylphenol ethoxylates, and is available in a variety of sizes to support processes from development to large scale commercial batches. J.T.Baker® Cell Lysis Solution may be stored at ambient temperature. The product is a ready to use solution in 100X concentration that has a light yellow to dark amber color. J.T. Baker® Workflow Solutions offer consistent performance, high quality, and cost-effective products, with assured availability due to our global presence. Visit [Biopharma solutions | Avantor \(avantorsciences.com\)](https://www.biopharmasolutions.com) for additional information on Avantor™ J.T.Baker® Workflow Solutions for CGT Workflows.

APPLICATION

J.T.Baker® Cell Lysis Solution is designed for the extraction of viral particles, such as Adeno-associated Virus (AAV) particles from HEK293 and SF9 cells and prevention of viral particle damage during the cell lysis step. J.T.Baker® Cell Lysis Solution has been shown to not impact enzymatic activity of endonuclease activity when used in combination with J.T.Baker® Endonuclease, enabling robust and efficient cell lysis and DNA/RNA degradation in a single step. J.T.Baker® Cell Lysis Solution does not contain Mg^{2+} . When used in combination with J.T.Baker® Endonuclease, or other enzymes for AAV vector bioprocessing, addition of Magnesium Chloride Hexahydrate, part number JT2449, to a target concentration ≥ 2 mM is recommended.

FEATURES

- Readily biodegradable, free of alkylphenols
- Comparable performance to Triton™ X-100
- Formulation protects viral particles
- Ease of use
- Reduced foaming and viscosity
- Quality based product specifications
- 0.1 μ m filtered solution
- Global cGMP manufacturing and distribution facilities
- Proven performance with J.T. Baker® Workflow solutions

BENEFITS

- Addresses environmental concerns with current industry standard
- Robust and efficient, providing rapid lysis over a range of cell densities
- Improved process yield
- Able to be cleared in typical downstream processing steps, detectable by routine analytical methods
- Minimize processing challenges
- Consistent performance for research use to commercial scale
- Elimination of aseptic manipulation to reduce chance of contamination
- Reliable supply and availability in most countries worldwide
- End to end solutions to streamline workflows and supply

SPECIFICATIONS

Appearance	Tested in 3 – 5 mL volume in scintillation vial, clear and colorless or yellow solution*
Refractive Index	1.1 – 1.5
Viscosity at 25°C	For Information Only (Target 7.5 to 9.8 cP)
Density	0.9 – 1.2 g/mL
Target pH	For information only
Total Aerobic Microbial Count	Max 100 CFU/mL
Total Yeast and Mold Count	Max 10 CFU/mL
UV-Vis Absorption maximum	For Information Only – record at 380-450nm

*Product may have a light yellow to dark amber color when observed in larger volume.

Avantor Performance Materials, LLC | 100 Matsonford Road, Suite 200 | Radnor, PA 19087, United States

PROPERTIES

J.T.Baker® Cell Lysis Solution exhibits comparable cell lysis efficiency to Triton™ X-100 at typical usage concentration. Cell lysis efficiency is shown by measurement of cell viability in HEK293 cells, initial cell density approximately 5.0×10^6 vc/mL, following lysis using J.T.Baker® Cell Lysis Solution at 1:100 dilution for one-hour, ambient temperature. (Figure 1). J.T.Baker® Cell Lysis Solution is effective over a broad range of operating conditions supporting variable incubation times (Figure 1), a wide range of cell densities (Figure 2) and shows robustness across a wide range of temperatures (Figure 3). J.T.Baker® Cell Lysis Solution offers protection of viral vectors from sheer stress caused by agitation during cell lysis, enabling higher recovery compared to traditional cell lysis solution. In one study, both Triton™ X-100 and J.T.Baker® Cell Lysis Solution were added to viral vector solutions at typical use-dilutions of 0.1% and 1%, respectively, and agitated for one-hour at 250 RPMs. Recovery was determined using qPCR. Results indicate viral vector recoveries increased 15% using J.T.Baker® Cell Lysis Solution compared to use of Triton™ X-100. (Figure 4). J.T.Baker® Cell Lysis Solution is effective on HEK293 and SF9 cells, demonstrated by incubating initial densities of approximately 6.0×10^6 cells/mL with 1:100 Cell lysis solution, resulting in 4.5% and 1.5% cell viability, respectively, after a 5-minute incubation.

Figure 1: J.T.Baker® Cell Lysis Solution efficacy comparing various cell lysis solutions and incubation times at 37°C




Figure 2: J.T.Baker® Cell Lysis Solution effect on various HEK293 cell concentrations after 60 minutes at 16 - 24°C




Figure 3: J.T.Baker® Cell Lysis Solution effectiveness on HEK293 cells at an initial density of 5×10^6 cells/mL after a one-hour incubation across various temperatures.





Figure 4: J.T.Baker® Cell Lysis Solution viral vector recovery versus traditional cell lysis solution demonstrating a 15% increase in recovery over Triton™ X-1



INSTRUCTIONS FOR USE

Following transfection and production of AAV, harvest adherent cells and pool into a suitable container. It may be possible to add directly to culture vessel used for adherent cell line. For suspension cells, J.T.Baker® Cell Lysis Solution may be added directly to culture. J.T.Baker® Cell Lysis Solution is typically used at a concentration of 1:100 (v/v). Incubate under agitation or on a shaker at 4 - 37°C for 60 - 120 minutes. Specific process conditions should be confirmed by end-user with in-house testing. Post-incubation, transfer cell lysate to an appropriate container for centrifugation and further processing. Store J.T.Baker® Cell Lysis Solution at 15 - 25°C (ambient temperature).