



KOVA® Plastics System for Standardized Urinalysis

KOVA SYSTEM FOR STANDARDIZED URINALYSIS

Urinalysis, as currently performed in many laboratories, is carried out using various nonstandard procedures. These procedures vary from laboratory to laboratory and often the actual technique within the laboratory varies depending on the person performing the tests.

Sources of variation in conventional urinalysis:

- variable urine volumes
- different centrifugal conditions creating varying amounts of sediment for microscopic examination
- different amounts of sediment collected and suspended under the cover glass
- technique variations among individuals performing the procedure.

To standardize the urinalysis procedure, a constant specimen volume, centrifugal force and sediment volume must be maintained and a consistent method of microscopic examination and reporting of results should be used. Kova's KOVA System achieves this standardization by reducing variation, including technique differences among technicians.

INTENDED USE

The **KOVA System** offers a procedure and products that can be used to produce standardized results during routine urinalysis. Volume control, consistency and hygiene are provided from collection and transport to microscopic analysis of urine sediment. **KOVA-Trol™** can be used for complete quality control of physical, chemical and microscopic examination test procedures.

ADVANTAGES

If the described procedure is followed consistently, one can use the values obtained in urinalysis with confidence. Clinicians can follow the progress and treatment of patients with certainty; any changes that occur outside the narrower limits that this system allows can be considered significant.

Laboratories may be compared and patients under observation can have their urinalysis done at different laboratories with comparable results.

KOVA SYSTEM AND SYSTEM COMPONENTS

Product Number	Product Description	Determinations Per Package
87153	KOVA Plastics System Super Pac 1000 w/Caps 100 KOVA Glasstic® Slide 10 (10 chambered), 1000 KOVA Petters, 1000 KOVA Super Tubes, 1000 KOVA Caps	1000
87154	KOVA Plastics System Super Pac 1000 100 KOVA Glasstic® Slide 10 (10 chambered), 1000 KOVA Petters, 1000 KOVA Super Tubes,	1000
87162	KOVA Plastics System Super Pac 1000 with Grids 100 KOVA Glasstic® Slide 10 (10 chambered) w/Grids, 1000 KOVA Petters, 1000 KOVA Super Tubes,	1000
87155	KOVA Plastics System Pac II 100 KOVA Slide II (4 chambered), 400 KOVA Petters, 400 KOVA Super Tubes,	400
87159	KOVA Plastics Mini Pack 5 KOVA Glasstic® Slide 10, 20 KOVA Petters, 20 KOVA Tube, 20 KOVA Caps	50
87156	KOVA Plastics System Value Pac 500 50 KOVA Glasstic Slide 10 with grids 500 KOVA Economy Tubes, 100 KOVA Caps	500
87158	KOVA Plastics System Value Pac 500 with Grids 50 KOVA Glasstic Slide 10 (10 chambered), 500 KOVA Petters, 500 Economy Tubes	500
87141	KOVA Plastics KO-LEC-PAC® 500 KOVA Super Tubes, 500 KOVA Caps, 500 KOVA Cups, 500 Labels and 5 Transport Racks	500
87100	KOVA Plastics Slide II with Grid for quantitation; 100 x 4 well slides; with each 1mm x 1mm grid square	400
87118	KOVA Plastics Slide II (without grid) 100 x 4 well slides	400
87146	KOVA Plastics Glasstic® Slide IO 100 x 10 well slides in crystal clear Acrylic	1000
87157	KOVA Plastics Glasstic® Slide IO 50 x 10 well slides in crystal clear Acrylic	500
87144	KOVA Plastics Glasstic® Slide IO with Grid 100 x 10 well slides in crystal clear Plexiglas® with quantitation grids; each chamber contains 6.6 µl and has a 3 mm x 3 mm grid with fine divisions of 0.33 mm x 0.33 mm. The test procedure includes a method for quantitating cells per µl of patient samples.	1000

KOVA SYSTEM AND SYSTEM COMPONENTS - CONTINUED

Product Number	Product Description	Determinations Per Package
87137	KOVA Plastics Super Tube Graduated non-sterile disposable collection and centrifuge tubes made of high impact, unbreakable plastic to eliminate cracking or breaking during centrifugation.	500
87138	KOVA Plastics Economy Tube As above but in economical, break-resistant styrene plastic.	500
87135	KOVA Plastics Petter Disposable plastic transfer pipette designed to retain 1.0 ml of urine after centrifugation. The unique lock tip provides a one-step contamination-free decanting method.	500
87139	KOVA Plastics Cap Recommended to prevent spillage during transportation, as well as aerosol contamination during centrifugation.	500
87136	KOVA Plastics Decanting Rack Rack for decanting up to 10 specimens.	1 Rack
87116	KOVA Stain Modified Sternheimer Malbin stain, facilitates differentiation of formed elements.	3 x 25 mls

KOVA SYSTEM FOR QUALITY CONTROL - KOVA Controls

Human urine controls for physical, chemical and microscopic quality control.

Product Number	Description KOVA-Trol	Packaging
87329	KOVA-Trol I - High Abnormal Without Abnormal Urobilinogen Value Assignment 4 x 15 mls	
87325	KOVA-Trol I - High Abnormal Without Abnormal Urobilinogen Value Assignment 4 x 60 mls	
87326	KOVA-Trol I - High Abnormal Without Abnormal Urobilinogen Value Assignment 10 x 60 mls	
87332	KOVA-Trol I - High Abnormal With Urobilinogen Value Assignment	4 x 60 mls
87333	KOVA-Trol I - High Abnormal With Urobilinogen Value Assignment	10 x 60 mls
87334	KOVA-Trol I - High Abnormal With Urobilinogen Value Assignment	4 x 15 mls
87130	KOVA-Trol II - Low Abnormal	4 x 15 mls
87128	KOVA-Trol II - Low Abnormal	10 x 60 mls
87331	KOVA-Trol III - Normal with hCG	4 x 15 mls
87327	KOVA-Trol III - Normal with hCG	4 x 60 mls
87328	KOVA-Trol III - Normal with hCG	10 x 60 mls

Simulated urine controls for physical, chemical and microscopic quality control.

Product Number	Description KOVA Liqua-Trol with Microscopics	Packaging
87122	KOVA Liqua-Trol Level II (Normal with hCG)	2 x 120ml
87176	KOVA Liqua-Trol Level I (Abnormal)	2 x 120ml
87123	KOVA Liqua-Trol Level II (Normal with hCG)	4 x 120ml
87177	KOVA Liqua-Trol Level I (Abnormal)	4 x 120ml

Product Number	Description KOVA Liqua-Trol without Microscopics	Packaging
87112	KOVA Liqua-Trol Level II (Normal with hCG) Control	3 x 15ml
	KOVA Liqua-Trol Level I (Abnormal) Control	3 x 15ml

SPECIMEN COLLECTION AND TRANSPORT

The KOVA Plastics System KO-LEC-PAC® is recommended for use in the following manner:


1. Label the KOVA Plastics Tube and give the patient a 3 ½ oz. KOVA Plastics Cup.
2. Instruct the patient to collect the voided urine in the KOVA Plastics Cup.
3. Transfer the urine specimen from the KOVA Plastics Cup to the KOVA Plastics Tube, filling it to the 12 ml graduation.
4. Secure the KOVA Plastics Cap on the KOVA Plastics Tube and place it in the KOVA Plastics Transport Rack for transportation and storage.
5. Deliver to the laboratory as soon as possible, preferably within two hours, but no more than four hours following specimen collection.

KOVA SYSTEM TEST PROCEDURE

- Check the specific gravity by placing one or two drops of urine in a temperature-compensated refractometer or use a chemistry test strip containing a specific gravity parameter and record the results.
- Using reagent test strips, perform chemical testing according to the manufacturer's instructions. Record the observed results. KOVA-Trol I, II and III should be included in each batch to assure proper quality control of physical, chemical and microscopic test procedures.
Note: KOVA-Trol I, II or III may be used individually or as a set to evaluate the specific range or application.
- Centrifuge the KOVA Plastics Tubes (each containing 12ml of urine specimen or KOVA-Trol) at a relative centrifugal force (rcf) of 400 for five minutes; approximately 1500 revolutions per minute (rpm) with a 6-inch radius rotor. Formula used:

$$rcf = 28.38 (R) \left(\frac{N}{1000} \right)^2 \quad \begin{matrix} R = \text{radius of rotor in inches} \\ N = \text{revolutions per minute} \end{matrix}$$

The rotating radius is the distance measured from the rotor axis to the tip of the liquid inside the tubes at the greatest horizontal distance from the rotor axis.



- Remove the KOVA Plastics Tubes from the centrifuge being careful not to disturb or dislodge the sediment.
- Insert a KOVA Plastics Petter into the KOVA Plastics Tube. Push the KOVA Plastics Petter to the bottom of the KOVA Plastics Tube until it seats firmly (at the 1ml graduation).
- Decant and discard 11ml from the KOVA Plastics Tube while the KOVA Plastics Petter is locked in position in the KOVA Plastics Tube. This will retain 1ml of urine sediment at the bottom of the KOVA T Plastics ube.
- Withdraw the KOVA Plastics Petter from the KOVA Plastics Tube.
- Add one drop of KOVA Stain to the 1ml of urine sediment.
Note: KOVA Stain is an aid to assist in the cellular differentiation of elements and is optional.
- Using the KOVA Plastics Petter, gently resuspend the sediment and stain until a homogeneous mixture is obtained.
- Withdraw a small sample of the urine sediment stain mixture by squeezing the bulb of the KOVA Plastics Petter.
- Transfer the sediment mixture to the KOVA Plastics Slide by placing one drop in the cut-out notch of each chamber. When chambers 1-5 are on the top row, the notch is at the top left corner of chambers, when chambers 6-10 are on the top row, the notch is at the top right corner of chambers. The chamber will fill by capillary action. Avoid touching the V-shaped barrier between the chambers while dispensing fluid. Incorrect positioning in dispensing may cause overflowing from one chamber to the next.
- Remove any excess specimen remaining on the open recessed area by touching the open edge with absorbent material.
- Place the KOVA Plastics Slide on a microscopic stage under the objective lens.
- Scan the slide chamber under low power magnification (10X eyepiece/10X objective) to enumerate casts. Enumerate all other formed elements under high power magnification (10X eyepiece/40X objective). Do not reuse KOVA Plastics products.

EXPECTED VALUES - MICROSCOPY†

HPF = High Power Field 400X

LPF = Lower Power Field 100x

- 1 + = Occasional form noted
2 + = Noted in every field
3 + = Large amounts in every field
4 + = Full field

Analyte	Normal	Abnormal	Reporting Results
WBC	0-5/HPF	> 5/HPF	Numbers/HPF
RBC	0-3/HPF	> 3/HPF	Numbers/HPF
Epithelial Cells	0	Any (other than squamous)	Numbers/HPF
Crystals	0-3/HPF (non-pathogenic)	> 3 Any abnormal	Numbers/HPF
Yeasts	0	Any	1 + to 4 +/HPF
Trichomonads	0	Any	1 + to 4 +/HPF
Casts	0	Any especially > 1 hyaline cast/LPF	Numbers/LPF
Bacteria	0-5/HPF	> 5/HPF	1 + to 4 +/HPF
Fat	0	Oval fat bodies or free fat	1 + to 4 +/HPF

† Bernard Statland, MLO. p 13-14; Jan. 1985

References for General Information

- Bradley, G.M., Benson, E.S., Todd-Sanford Clinical Diagnosis by Laboratory Methods, 15th Edition, Phila. Saunders, 1974.
- Kurtzman, N.A. and Rogers, P.W. (1974). A Handbook of Urinalysis and Urinary Sediment. Chas. C Thomas, Springfield, IL.
- Little, P.J. (1962). Urinary white-cell excretion. Lancet. pp. 1149-1151.
- Little, P.F. (1964). A comparison of the urinary white cell concentrations with the white cell excretion rate. Brit. J. Urol. 36, 360-363.
- Thomas, M.(1971). A rapid slide method of urine cell counts. Med. Lab Technol. 28, 38-39.
- Moore, T., Hira, N.R., and Stirland, R.M. (1965). Differential urethrovessical urinary cell count. Lancet. pp. 626-627.
- Siegle, M.D., Lab Med., 12:781, 1981.
- Sternheimer, R. and Malbin, B. (1951). The clinical recognition of pyelonephritis with a new stain for urinary sediments. Am. J. of Med., 11:312-323.
- Muschetta, P.A. and Waters, Jr. F.O. (1962). Manual of Medical Laboratory Techniques. Herbert-Spence, Inc. New York, N.Y., Second Edition, pp 44-45.
- Lippman, R. W. (1957). Urine and the Urinary Sediment. Chas. C Thomas, Springfield, IL.
- Dudas, H.C., Lab Med. 12:765. 1981.
- Weller, J.M. and Greene, J.A. (1966). Examination of the Urine. Meredith Publishing Co., New York.

U.S. PATENT NUMBERS

4,022,576
4,563,332
4,299,441
D258,388
4,637,693
4,937,415
4,997,266
RE 33,826



Symbols

	Batch/Lot Code		Catalog Number		Unique Device Identifier		Manufactured in Country
	Expiration/Use By		Contains Quantity		Invitro Diagnostics Use		Storage Limits
	Manufacturer		Do not reuse		Instructions for Use/Electronic Instructions for Use		

www.kovaplastics.com



Alltrista Plastics, LLC
20 Setar Way Reedsville, Pa 17084 United States
Customer Service: +1 864-879-8100

EC REP

Advena Ltd. Tower Business Centre, 2nd Flr.
Tower Street, Swatar, BKR 4013 Malta