



Laboratory Centrifuge VWR Micro Star 30R



INSTRUCTION MANUAL

European Catalogue Numbers:

521-2574 Micro Star 30R inclusive angle rotor RA 24/2
and aerosol tight rotor lid

521-2575 Micro Star 30R without rotor

Version: 2

Issued: 21.10.2020



Legal Address of Manufacturer

Europe

VWR International bvba
Researchpark Haasrode 2020
Geldenaaksebaan 464
B-3001 Leuven
+ 32 16 385011
<http://be.vwr.com>

Country of origin: SLOVENIA

PREFACE

Dear customer,

Congratulations for purchasing our centrifuge Micro Star 30R. You have selected a reliable device which combines many advantages.

A wide selection of programming options and electronic control allow trouble-free operation. The brushless motor provides quiet, maintenance-free operation without any carbon dust pollution.

Your device is equipped with user-friendly functions, which makes the operation and settings easier for you. Built-in safety functions prevent incorrect usage and they control the complete operation of the centrifuge.

The centrifuge has the functionality to save programs. You can save up to 100 different data sets in the memory. The centrifuge always keeps the last run program in its memory for an unlimited amount of time, allowing the program to be restarted at any time - even if the centrifuge was turned off in-between. All important operation parameters can be seen at a glance on the displays.

The settings are executed via a knob and the keys on the control panel.

The centrifuge and its interior are easy to clean.

We thank you for your confidence and we hope you enjoy your new centrifuge.

TABLE OF CONTENTS

Page

PREFACE	2
TABLE OF CONTENTS.....	3
1. SAFETY INFORMATION.....	4
1.1 WARNING SYMBOLS.....	4
1.2 SAFETY INSTRUCTIONS	4
2. INTRODUCTION.....	5
2.1 INTENDED USE.....	5
2.2 DELIVERY CONTENT AND OPTIONAL EQUIPMENT	5
2.3 UNPACKING	5
2.4 INSTALLING THE CENTRIFUGE	6
2.5 OVERALL VIEW.....	6
3. OPERATION	7
3.1 CONTROL PANEL	7
3.2 KNOB AND KEY FUNCTIONS	8
3.3 TURNING ON THE CENTRIFUGE.....	9
3.4 OPENING AND CLOSING OF THE LID.....	10
3.5 ROTOR MOUNTING AND DISMOUNTING	10
3.6 IMBALANCE INDICATOR.....	11
3.7 PRE-SELECTING TEMPERATURE.....	12
3.8 PRE-SELECTING SPEED/RCF.....	12
3.9 PRE-SELECTING RUN TIME, ACCELERATION AND BRAKING LEVEL	13
3.10 PROGRAM SETTING	14
3.11 STARTING AND STOPPING THE CENTRIFUGE	15
3.12 ROTOR RECOGNITION	15
3.13 DISPLAY OF THE SET VALUES DURING THE RUN	15
3.14 FAST SPIN.....	15
3.15 FAST COOL FUNCTION	16
3.16 RESET OF LID LATCH	16
3.17 AUXILIARY LID RELEASE	16
3.18 AUTOMATIC BRAKING AFTER POWER RESTORATION.....	16
3.19 TEMPERATURE CORRECTION.....	17
4. MAINTENANCE AND CARE.....	17
4.1 REGULAR SERVICE	17
4.2 CLEANING THE CENTRIFUGE	17
4.3 CLEANING THE ROTOR.....	18
4.4 ROTOR STERILIZATION	18
4.5 SHIPPING	18
5. TROUBLESHOOTING.....	19
5.1 ERROR MESSAGES	19
5.2 EXIT FROM ERROR DISPLAY	20
5.3 USER'S PARAMETERS	21
5.4 REPLACING THE CENTRIFUGE FUSES.....	22
6. TECHNICAL DATA.....	22
7. CALCULATIONS.....	23
7.1 CALCULATION OF CENTRIFUGAL FORCE.....	23
7.2 CALCULATION OF MAXIMUM PERMITTED ROTOR SPEED	23
8. TECHNICAL SERVICE.....	23
9. WARRANTY.....	24
10. COMPLIANCE WITH LOCAL LAWS AND REGULATIONS	24
11. EQUIPMENT DISPOSAL	24
EU DECLARATION OF CONFORMITY.....	25
YOUR DISTRIBUTOR	27

1. SAFETY INFORMATION

Before using the machine, make sure to read and understand this manual thoroughly. Keep the manual close to the machine, easily accessible to all the users.

Improper operation can cause injury to persons or damage to the equipment.

1.1 WARNING SYMBOLS

The following are the warning symbols that are used in this manual.

	This symbol indicates a potential risk and alerts you to proceed with caution.
---	--

1.2 SAFETY INSTRUCTIONS

In the interest of your own personal safety, always observe the following safety instructions:

- Do not use rotors and buckets, which show clear signs of corrosion or mechanical defects. Please check the accessories at regular intervals.
- Always load the rotor with the same test tubes on all positions or symmetrically with the pairs of the same test tubes. To prevent negative consequences of unbalanced rotor, like damaged bearings and motor axle, or inadequate results of centrifugation, equal or equally loaded test tubes must be symmetrically arranged according to the rotation axis.
- Please use only the original accessories for centrifugation.
- Do not move or knock the centrifuge during operation!
- Repairs must only be performed by an authorized service technician.
- The centrifuge may only be used for specified applications. It may not be used in a hazardous or potentially flammable environment or for centrifugation of explosive or highly reactive substances.
- When handling toxic, aggressive or radioactive materials, observe national regulations or regulations defined by World Health Organization.
- Fluids or materials used for cleaning and disinfecting should be disposed of in accordance with approved laboratory regulations.
- If any liquids are spilled in the rotor chamber, on the rotor or accessories, the surfaces must be cleaned immediately. You can use a damp cloth and mild soap solution. This is particularly important for the cleaning of the bores of the fixed-angle rotors.
- Density of the liquid must not be exceeded 1.2 g/ml at maximum rotational speed.
- During longer spin times, test tubes may heat up. Observe the requirements and regulations specified by test tube manufacturer.
- The use of organic solvents and reagents may have adverse effect on the stability of plastic test tubes.
- Rotors are high-grade components which are subjected to extreme mechanical strain. Aluminium rotors are protected against corrosion.
- Please ensure that the rotors are protected from mechanical damage. Even slight scratches and cracks can cause severe inner damage to the rotor material.
- Please clean your rotors regularly using a neutral cleaning liquid (e.g. Extran). This will protect the rotors and maintain their service life.

2. INTRODUCTION

2.1 INTENDED USE

Centrifuge Micro Star 30R is a universal laboratory refrigerating centrifuge. It is designed for use in medical, scientific and industrial laboratories, for separating the substances with different specific densities, by centrifugal force. In particular, it is intended for preparation of human samples (body fluids), in approved test tubes, before further analysis. Not intended for in vitro diagnostic use. The maximum rotational speed of 18,000 RPM gives the centrifugal force of 31,150 xg.

2.2 DELIVERY CONTENT AND OPTIONAL EQUIPMENT

The following is enclosed with the centrifuge:

- 1 Instruction manual
- 1 Hexagonal rotor key (magnetic key holder on the back of the device, see page 7)
- 1 Power cord
- 1 Angle rotor RA 24/2 with aerosol tight lid for 24 x 1.5/2 ml tubes (only with ECN 521-2574)

Optional rotors and accessories:

- 521-2576 Angle rotor RA 24/2, 24 x 1.5/2 ml micro tubes with aerosol tight lid, max. RCF 31,150 xg, max. speed 18,000 RPM
- 521-2577 Angle rotor RA 16/5, 16 x 5 ml micro tubes with aerosol tight lid, max. RCF 21,000 xg, max. speed 14,000 RPM
- 521-3377 Angle rotor RA 8/15, 8 x 15 ml conical, max. RCF 4,400 xg, max. speed 6,000 RPM
- 521-3379 Adapter 5 ml, 13 x 75 mm (PK 1) for rotor RA 8/15, max. RCF 3,566 xg
- 521-3380 Angle rotor RA 6/50, 6 x 50 ml conical, max. RCF 4,400 xg, max. speed 6,000 RPM
- 521-2585 Adapter 15 ml conical (set of 6), for rotor RA 6/50
- 521-2586 Swing-out rotor for 4 microplates with 2 carriers, max. RCF 2,200 xg, max. speed 5,000 RPM
- 521-2587 Swing-out rotor RS 6/12 with buckets, 6 x 12 ml, max. RCF 2,500 xg, max. speed 4,400 RPM

2.3 UNPACKING

The weight of the centrifuge is 50 kg. To prevent possible injuries, at least two people should lift and carry the centrifuge by holding it at the bottom from opposite sides.

Open the carton box. Take out the accessories and remove the packaging material. Reach with your hands under the centrifuge and lift it from the box together with another person.



When lifting the centrifuge, never hold it by the front or top plastic part of the housing or by the lid, as the appliance may get damaged!

Retain the packaging material for any subsequent transport or storage, which are allowed only in the original packaging.

2.4 INSTALLING THE CENTRIFUGE

The centrifuge should only be operated indoors. Place the centrifuge on a stable, solid, horizontal and clean surface, without vibrations. Make sure that the centrifuge is not exposed to direct sunlight. To ensure sufficient ventilation, there should be enough clearance on all sides of the centrifuge. It must be far enough away from the wall and other devices. According to recommendations of the EN 61010-2-020 standard, a safety clearance of 30 cm should be observed around the centrifuge during operation. Please remove all objects from this area. If the centrifuge isn't levelled, imbalances can occur and the centrifuge can be damaged. Do not place anything under the centrifuge feet to level the centrifuge.



After installation you have to wait for at least four hours, before connecting the centrifuge to the mains power supply. This prevents possible damage to the compressor, which may result from improper transport.


This also prevents damage to electronic components due to condensation, which can occur, when you bring the device from a cold environment to a warm environment.


Before connecting the centrifuge to the mains power supply, check that the mains voltage and frequency correspond to the specifications on the name label of the centrifuge. The mains cable of the centrifuge may only be connected to a properly grounded wall socket.

To disconnect the mains supply from the centrifuge in the event of malfunction, an emergency switch separate from the centrifuge must be available. This switch should be outside the room, where the centrifuge is installed, or next to the entrance to the room.

2.5 OVERALL VIEW

Main parts of the centrifuge are designated on the following pictures:






3. OPERATION

Before using the centrifuge for the first time, familiarize yourself with the control panel and the functions of the knob and keys/icons.


3.1 CONTROL PANEL



1	Temperature display
2	TEMPERATURE key
3	Indicator for selection of rotational speed (RPM)
4	Display for rotational speed (RPM) or centrifugal force (RCF)
5	Indicator for selection of centrifugal force (RCF)
6	SPEED key
7	Display for time (TIME), acceleration level (▲) and braking level (▼)
8	Indicator for selection of time (TIME)
9	Indicator for selection of acceleration level (▲)

10	TIME key
11	Indicator for selection of braking level (▲)
12	Program display
13	PROGRAM key
14	Knob
15	FAST COOL key
16	LID key
17	START key
18	STOP key

3.2 KNOB AND KEY FUNCTIONS

 (Knob)	By rotating the knob, you change the values of the parameters on display. By rotating the knob clockwise, the value increases, and by rotating it counter-clockwise, the value decreases. By pressing on the knob, you confirm the new value for the parameter.
 (TEMPERATURE)	By pressing this key, you can set the temperature. The key lights up. Set new temperature by rotating the knob and confirm it by pressing on the knob. The key light will then turn off.
 (SPEED)	By pressing this key, you can set the rotational speed (RPM) or centrifugal force (RCF). The key lights up. By repeated pressings on this key, you select which parameter will be set. When indicator RPM is lit, you can set the rotational speed, and when indicator RCF is lit, you can set the centrifugal force. Set new values of both parameters by rotating the knob and confirm them by pressing on the knob. The key light will then turn off.
 (TIME)	By pressing this key, you can set the time (TIME), acceleration level (▲) and braking level (▲). The key lights up. By repeated pressings on this key, you select which parameter will be set. When indicator TIME is lit, you can set the run time, when indicator ▲ is lit, you can set the acceleration level, and when indicator ▲ is lit, you can set the braking level. Set new values of all parameters by rotating the knob and confirm them by pressing on the knob. The key light will then turn off.
 (PROGRAM)	By pressing this key, you can set the programs. The key lights up. Select new program by rotating the knob and confirm it by pressing on the knob. The key light will then turn off. If you want to change the values of parameters of an existing program, set the values of individual parameters as described above, when this key is lit.
 (FAST COOL)	By pressing this key, you turn on the fast cool function. The key lights up. Rotor starts to run automatically with 8,000 RPM. Fast cool function cools down the centrifugal chamber to the set temperature. When the temperature is reached, display for rotational speed shows message READY, and the key light will then turn off.

 (LID)	<u>SHORT PRESS:</u> By pressing this key, you open the lid of the centrifuge. When the lid opens, the key light turns off. When closing the lid, press it down with your hands, until it locks up. When the lid is closed, the key lights up. The key light turns off during the run of the centrifuge and thus indicates that you can't open the lid during the run. <u>LONG PRESS:</u> By pressing and holding this key, you set the lid latch to initial position. This is needed in case of power failure during the opening of the lid, when lid latch motor stops in undefined position. When power returns, it could happen that you can't open or close the lid. If this happens, press the key for approximately two seconds, until you hear the sound of lid latch motor, then release it. After that, you will be able to open and close the lid normally.
 (START)	<u>SHORT PRESS:</u> By pressing this key, you start the run of the centrifuge. The key lights up and thus indicates that the centrifuge is running. At the end of the set running time or manual stopping of the centrifuge, the braking procedure of the rotor is activated and the centrifuge stops. The next run is possible, when the rotor stops completely. <u>LONG PRESS:</u> By pressing and holding this key, you perform the fast spin. The centrifuge runs for as long as you hold the key pressed. Time of run in seconds is displayed on time display. The acceleration and braking levels are fixed to 9 (the highest), and they can't be changed. When you release the key, the braking procedure of the rotor is activated and the centrifuge stops.
 (STOP)	By pressing this key, you stop the run of the centrifuge. The key lights up and thus indicates that the centrifuge is stopping. The braking procedure of the rotor is activated and the centrifuge stops. Until the centrifuge stops completely, the START key remains lit, thus indicating that the rotor is still turning.

3.3 TURNING ON THE CENTRIFUGE

Use the enclosed power cord to connect the centrifuge to the mains power supply.

Turn on the main switch, which is located on the front of the right side of the centrifuge.

All segments (eights) are displayed on all displays at first, and all the keys and indicators light up simultaneously. Then centrifuge model (C-30R) and program version (X.XX) are displayed, then lines, and finally the values of operation parameters are displayed.

The values of operation parameters (temperature, rotational speed, centrifugal force, run time, acceleration and braking level) are automatically set to the last used values.

3.4 OPENING AND CLOSING OF THE LID

	By pressing the LID key, you open the lid of the centrifuge. When the lid is open, the key light turns off. When closing the lid, press it down with your hands, until it locks up. When the lid is closed, the key lights up. The key light turns off during the run of the centrifuge and thus indicates that you can't open the lid during the run.
---	---



Always open the centrifuge lid completely, to prevent it from falling!
When closing the lid, make sure to always place your fingers on the top side of the lid and never in the gap between the lid and the housing of the centrifuge, otherwise you could crush your fingers!

3.5 ROTOR MOUNTING AND DISMOUNTING

Before attaching the rotor on the motor axle, make sure that axle and rotor are clean and undamaged. Wipe all fixing surfaces (motor axle and rotor cone) with clean soft cloth. Thus you will avoid potential damages to the axle and motor.

- Mount the rotor on the motor axle and firmly tighten the rotor nut by turning it clockwise, using the supplied hexagonal rotor key.
- To dismount the rotor, turn the rotor nut counter-clockwise, using the hexagonal rotor key and then remove the rotor.



Do not use the rotors, rotor lids and test tubes, which are mechanically or chemically damaged or with visible corrosion defects!



The rotor and the rotor lid must always be securely fastened. Do not begin with centrifugation before the rotor has been securely fastened!



3.6 IMBALANCE INDICATOR

The centrifuge is fitted with an imbalance detector, to ensure safety.

Imbalance at high speed may indicate a tube breakage or leak or rotor crash. Therefore additional care should be taken depending on the samples loaded.

You must always load the rotor with test tubes symmetrically. You may only use approved test tubes. Weight difference of the samples in test tubes should be as low as possible in order to avoid potential damages to the motor and to minimize running noise and vibration.

Following are examples of correctly and wrongly loaded rotors:



If an imbalance is detected at speeds higher than approx. 300 RPM, the centrifuge will automatically stop and an error message IMBAL OUR will be displayed.

3.7 PRE-SELECTING TEMPERATURE

The temperature can be set from -9 °C to 40 °C, in steps of 1 °C.


If the difference between actual and set temperature is more than ± 3 °C, the temperature display will flash.

	By pressing the TEMPERATURE key, you can set the temperature. The key lights up and the current set temperature is displayed. Set new temperature by rotating the knob and confirm it by pressing on the knob. The key light turns off and the actual temperature is displayed.
---	---

NOTE:

To achieve very low temperatures, it is necessary to appropriately reduce the rotational speed of the rotor (e.g.: If you set the temperature to -9 °C, you may not be able to use the max. speed of 18,000 RPM, if you want the centrifuge to reach this very low temperature).

3.8 PRE-SELECTING SPEED/RCF

	By pressing the SPEED key, you can set the rotational speed (RPM) and centrifugal force (RCF). The key lights up. By repeated pressings on this key, you select which parameter will be set. When indicator RPM is lit, you can set the rotational speed, and when indicator RCF is lit, you can set the centrifugal force. Set new values of both parameters by rotating the knob and confirm them by pressing on the knob. The key light will then turn off. • Rotational speed (indicator RPM is lit) Rotational speed can be set from 200 to 18,000 RPM in steps of 10 RPM. The maximum rotational speed is automatically set according to the rotor type in use.
	• Relative Centrifugal Force (indicator RCF is lit) Relative centrifugal force can be set from 5 to 31,150 xg in steps of 10 xg. Relative centrifugal force is calculated according to the radius of used rotor, so the minimum and maximum forces are dependent on the rotor type.

3.9 PRE-SELECTING RUN TIME, ACCELERATION AND BRAKING LEVEL

By pressing the TIME key, you can set the time (TIME), acceleration level (▲) and braking level (▼). The key lights up. By repeated pressings on this key, you select which parameter will be set. When indicator TIME is lit, you can set the run time, when indicator ▲ is lit, you can set the acceleration level, and when indicator ▼ is lit, you can set the braking level. Set new values of all parameters by rotating the knob and confirm them by pressing on the knob. The key light will then turn off.

	• Run time Run time can be set between 0.10 and 99.5 minutes. Time setting from 0.10 to 9.59 minutes is possible in 1 second steps. Between 10.0 and 99.5 minutes the setting is possible in 10 second steps. By rotating the knob either clockwise until after 99.5 is displayed, or counter-clockwise until after 0.10 is displayed, you can set continuous operation (HLd on display). NOTE: 0.10 = 10 s, 9.59 = 9 min 59 s, 10.0 = 10 min, 99.5 = 99 min 50 s
	• Acceleration level (AcX) Acceleration level can be set from 0 to 9. Level 0 means very slow acceleration, level 9 very fast acceleration.
	• Braking level (brX) Braking level can be set from 0 to 9. Level 0 means stopping without braking, level 9 maximum braking.

NOTE:

All the operation parameters can be changed during the operation of the centrifuge, and it starts to work with the new settings after that.

Changing of run time during the operation:

You should be aware, that when you extend run time, the difference between newly set and originally set time is added to current time, and when you shorten run time, the difference is subtracted.

Example: The centrifuge started with time set to 10 minutes. It has been running for 3 minutes. Then you changed the time to 5 minutes. The centrifuge will run for another 2 minutes.

3.10 PROGRAM SETTING



You can save 100 programs, with different operational parameters, in the centrifuge memory.

NOTE:

You can cancel the procedure for program setting by pressing the STOP key and return to stand-by mode without confirming the data.

If you want to use the existing program, follow the next procedure:

By pressing the PROGRAM key, you can set and select programs. The key lights up. Select new program by rotating the knob and confirm it by pressing on the knob. The key light turns off.

Operation parameters set in selected program are displayed on display.

If you want to change the existing program, follow the next procedure:

By pressing the PROGRAM key, you can set the programs. The key lights up. Select the program that you want to change, by rotating the knob. While the PROGRAM key is lit, set the values of individual operation parameters (temperature, rotational speed, centrifugal force, run time, acceleration level and braking level) using the procedures, described in the previous chapters.

If you want to change another program, select it by rotating the knob, and set the values of individual operation parameters for this program. Otherwise skip this step.

By rotating the knob, select a program that you want to use and confirm it by pressing on the knob. The light for the PROGRAM key will then turn off.

Operation parameters set in selected program are displayed on display.

The following table shows factory preset values of operation parameters for all 100 programs:

Program	Temperature (°C)	Rotational speed (RPM)	Run time (min)	Acceleration level (AcX)	Braking level (brX)
0	4	13,000	10	6	6
1	4	10,000	2	9	9
2	4	1,500	10	6	6
3	4	3,000	10	9	9
4	10	4,000	3	6	6
5	20	15,000	3	9	9
6	4	4,000	3	6	6
7	20	8,000	5	9	9
8	20	5,000	10	6	6
9	4	16,000	10	9	9
10	25	1,200	10	5	5
11	20	3,000	5	5	5
⋮	20	3,000	5	5	5
99	20	3,000	5	5	5

3.11 STARTING AND STOPPING THE CENTRIFUGE

Before starting the centrifuge, turn on the main switch, open the centrifuge lid, insert the rotor and load it with test tubes, close the centrifuge lid and set the operation parameters or select a suitable program.

	By pressing the START key, you start the run of the centrifuge. The key lights up and thus indicates that the centrifuge is running. After the expiration of set running time or after manual stopping of the centrifuge, the braking procedure of the rotor is activated and the centrifuge stops. The next run is possible, when the rotor stops completely.
	By pressing the STOP key, you stop the run of the centrifuge. The key lights up and thus indicates that the centrifuge is stopping. The braking procedure of the rotor is activated and the centrifuge stops. Until the centrifuge stops completely, START key is lit, thus indicating that the rotor is still turning.



Do not move or knock the centrifuge during operation!

3.12 ROTOR RECOGNITION

Rotor recognition is executed automatically, every time when the centrifuge starts. After rotor change, the centrifuge stops and display shows ROTOR CHG. Press the STOP key to clear this message. New maximum rotational speed is automatically set, to correspond with the inserted rotor. To start the centrifuge again, press the START key.

3.13 DISPLAY OF THE SET VALUES DURING THE RUN

During the run, the centrifuge shows the current values of the operation parameters on displays. If you wish to check the set values of the parameters, select the desired parameter by pressing the corresponding key. The key lights up. Display shows the set value of the parameter. To exit the display of the set values, press on the knob. The key light will then turn off.

3.14 FAST SPIN

	By pressing and holding the START key, you perform the fast spin. The centrifuge will run as long as you hold the key pressed. Time of run is displayed on time display, and the program display shows FS (fast spin). The acceleration and braking levels are fixed to 9 (the highest), and they can't be changed. When you release the key, the braking procedure of the rotor is activated and the centrifuge stops.
---	--

3.15 FAST COOL FUNCTION

	By pressing the FAST COOL key, you turn on the fast cool function. The key lights up. START key lights up also. The program display shows FC (fast cool). Rotor RA 24/2 and RA16/5 starts to run automatically with 8,000 RPM. Fast cool function cools down the centrifugal chamber to the set temperature. When the temperature is reached, display for rotational speed shows message READY, and the light for the FAST COOL key turns off. Rotor keeps on running and START key is still lit, until you stop the centrifuge by pressing the STOP key. If you don't stop the centrifuge manually, fast cool function turns off automatically after 1 hour.
---	--

NOTE:

If maximum speed of inserted rotor is lower than 8,000 RPM, the rotor will be running with the maximum allowed speed for that rotor.

3.16 RESET OF LID LATCH

	By pressing and holding the LID key, you set the lid latch to initial position. This is needed in the case of power failure during the opening of the lid, when lid latch motor stops in undefined position. When power returns, it could happen that you can't open or close the lid. If this happens, press the key for approximately two seconds, until you hear the sound of lid latch motor, then release it. After that, you will be able to open and close the lid normally.
--	--

3.17 AUXILIARY LID RELEASE

In case of power failure, the lid can be opened manually. If power failure occurs during the operation of the centrifuge, the rotor may continue rotating for several minutes, before it stops.



Before auxiliary lid release, turn off the main switch of the centrifuge and wait until rotor fully stops. Check this by looking through the lid window. Otherwise, injury from rotating rotor may occur!

On the underside of the centrifuge, by the front left foot, there is a plastic plug, which you pull out of the hole. There is a string fastened to the plug. Pull the string to open the lid of the centrifuge. Then insert the string and the plug back in the hole.

3.18 AUTOMATIC BRAKING AFTER POWER RESTORATION

If power failure occurs during the operation of the centrifuge, the rotor may continue rotating for several minutes, before it stops. When power returns, the centrifuge checks if the rotor is still rotating. If it is still rotating then the display shows message Error brX (X is a number from 0 to 9 and represents the set braking level). The centrifuge begins with the automatic rotor braking, which is performed with preset braking level. When rotor stops, display shows error message MAINS INT and thus warns you, that power failure occurred during operation and the centrifugation was interrupted. By pressing the STOP key, you exit from error display and go back to stand-by mode. You can use the centrifuge again.

3.19 TEMPERATURE CORRECTION

If you notice a difference between preset temperature on display and actual measured temperature, the centrifuge has temperature correction functionality. Temperature correction is performed by entering the temperature difference in the service parameters of the centrifuge.

Temperature correction may only be performed by an authorized service technician!

Measured temperature depends on a number of factors, like environment temperature, rotational speed, run time, etc.

4. MAINTENANCE AND CARE

4.1 REGULAR SERVICE

We recommend having the centrifuge and associated rotors checked by authorized service at least once a year. You must thoroughly clean and disinfect the centrifuge prior to service.

4.2 CLEANING THE CENTRIFUGE

After every centrifugation, please remove any condensed water from the rotor chamber using a soft, absorbent cloth. After use, leave the centrifuge lid open, so the inside can dry out.

For regular cleaning of the outside surface of the centrifuge and the rotor chamber, use mild neutral detergent. Make sure that no liquid penetrates the inside of the housing.

Open the lid of the centrifuge and turn off the main switch. Disconnect the power cord from mains socket. Remove the rotor with rotor key. Clean all accessible surfaces of the device and accessories at least once a week and every time, when contaminated.

After cleaning with detergent, the rubber seal around the rotor chamber should be thoroughly cleaned with water and lubricated with glycerine, to prevent it from becoming brittle.

For cleaning and disinfection, use only neutral cleaners and disinfectants. Before cleaning or decontaminating the centrifuge, using means and methods not recommended in this manual, you should consult with the manufacturer, in order to avoid the damage to the centrifuge.

To ensure safe and long operation of the centrifuge, please avoid the use of aggressive chemicals, which can damage the centrifuge, rotor and accessories. Please check them regularly for damage caused by corrosion.

Cleaning intervals:

For the sake of personal, environmental, and material protection, it is your duty to clean and if necessary disinfect the centrifuge on a regular basis.

Maintenance	Recommended interval
Clean rotor chamber	Daily or when contaminated
Clean rotor	Daily or when contaminated
Accessories	Daily or when contaminated
Cabinet	Once per month
Ventilation holes	Every six months



Refrain from using any other cleaning or decontamination procedure than those recommended here, if you are not entirely sure that the intended procedure is safe for the equipment.
Use only approved cleansers.
If in doubt, contact VWR.

4.3 CLEANING THE ROTOR

The rotor and accessories must be regularly cleaned to prevent contamination and corrosion caused by residue. Check the rotor and the centrifuge housing monthly. This applies in particular to the rotor bores. For cleaning the rotor use a neutral cleaning liquid. This will protect the rotor and extend its service life.



Do not use damaged rotors and accessories for centrifugation!

To avoid the damage to the rotor, replace the sealing rings regularly.

4.4 ROTOR STERILIZATION

The rotors are autoclavable at 121 °C, for 20 minutes. After the rotor has been autoclaved for a maximum of twenty times, seals of the rotor must be replaced (this is valid for rotors with seals).

4.5 SHIPPING

Before shipping the centrifuge, please bear the following in mind:
The centrifuge must be clean and decontaminated.

The decontamination must be confirmed with a decontamination certificate. Contact customer service for more details.

5. TROUBLESHOOTING

5.1 ERROR MESSAGES

If an error occurs during the operation of the centrifuge, an error message appears on the display, and the centrifuge stops automatically. Temperature display and program display show message Er, while speed display and time display show the error message. For the list of errors, see below table.

DISPLAY		PROBLEM	SOLUTION	WHO REPAIRS
SPEED	TIME			
ROTOR	CHG	Rotor change	Repeat run	User
	HI	Rotor speed too high	Reduce speed Check rotor sensor	User Service-SP
IMBAL	SEN	Imbalance sensor	Check imbalance sensor	Service-SP
			Check imbalance sensor	Service
			Electronics error	Service
	OUR	Imbalance too high	Check rotor loading arrangement	User
			Check the samples weight in the rotor	User
			Check if the rotor and rotor lid are fastened	User
			Check the rotor and the lid	User
			Check rotor sensor	Service-SP
			Repeat the balancing procedure	Service
LID	OPN	The lid of the centrifuge open	Close the lid of the centrifuge	User
	SEN	Lid latch not engaged	Open and close the lid of the centrifuge again	User
			Check lid sensors	Service-SP
MOTOR	SEN	Speed sensor error	Check speed sensor	Service-SP
			Check speed sensor	Service
			Electronics error	Service
	SPD	Speed deviates for more than ± 500 RPM / 5 s	Check rotor, motor and frequency regulator	Service
			Electronics error	Service
T CHA	SEN	Chamber temperature sensor error	Check chamber temperature sensor	Service-SP
			Check chamber temperature sensor	Service
			Electronics error	Service
	OUR	Chamber temperature over 50 °C	Check refrigerating system	Service

DISPLAY		PROBLEM	SOLUTION	WHO REPAIRS
SPEED	TIME			
T CON	SEN	Condenser temperature sensor error	Check condenser temperature sensor	Service-SP
			Check condenser temperature sensor	Service
			Electronics error	Service
	OUR	Condenser temperature over 80 °C	Check refrigerating system	Service
			Check refrigerating parameters	Service-SP
T REC	TI	Set temperature not reached in 2 hours	Check the suitability of centrifuge installation and environmental conditions	User
			Check refrigerating system	Service
			Check refrigerating parameters	Service-SP
DRIVE	HIV	Voltage overload on the DC link	Reduce braking level	User
			Error on frequency regulator	Service
			Error on braking resistor	Service
	LOV	Voltage too low on the DC link	Check power supply	Service
	OC	Current overload of the motor	Repeat run	User
			Reduce acceleration level	User
			Check start-up parameters	Service-SP
			Check motor	Service
			Error on frequency regulator	Service
	HOT	Driver temperature too high	Reduce speed	User
			Check motor	Service
			Error on frequency regulator	Service
MAINS	INT	Power failure during the run	Repeat run	User
READY		Fast cool function reached the set temperature	Stop the centrifuge by pressing the STOP key	User
-----	---	After 1 hour of standstill centrifuge goes to sleep mode (lines on all displays)	Press any key	User

Note: SP = service parameters

5.2 EXIT FROM ERROR DISPLAY



By pressing the STOP key, you exit from error display and go back to stand-by mode.

If error is still displayed, turn off the main switch of the centrifuge and turn it on again.

If the centrifuge still doesn't return to stand-by mode, call service!

5.3 USER'S PARAMETERS

With user's parameters, you can check and set some of the operational properties of the centrifuge.

For entering to user's parameters, the centrifuge must be in stand-by mode. Simultaneously press the TEMPERATURE and PROGRAM keys and hold them for about 5 seconds. When the first parameter is displayed and the TIME key lights up, release the keys.


You can select the parameters by rotating the knob. If you want to change the setting of the parameter, press the SPEED key, which lights up. The light for the TIME key will then turn off. Then change the setting of the parameter by rotating the knob and confirm it by pressing on the knob. The light for the SPEED key will turn off and the TIME key will light up.

To exit from user's parameters, press the START key.

DISPLAY		PARAMETER DESCRIPTION
SPEED	TIME	
dIS	bSE	Block the changing of operation parameters during the operation of the centrifuge. (dIS = disable, Enb = enable)
dECrE	tIn	Decreasing or increasing of time display. (dECrE = decreasing, InCrE = increasing)
Enb	bEP	Turn on the beeper. (dIS = disable, Enb = enable)
1	SLE	Time, after which the centrifuge goes to sleep mode and fast cool function turns off. (0 - 99 hours)
Enb	COL	Start of cooling at the closing of the lid. (dIS = disable, Enb = enable)
dIS	tSS	Time starts running after the speed is reached, instead at motor start. (dIS = disable, Enb = enable)
dIS	tSC	Time starts running after the temperature is reached, instead at motor start. (dIS = disable, Enb = enable)
Enb	SAP	Temperature compensation (display of sample temperature instead of chamber temperature). (dIS = disable, Enb = enable)

Note: Values in SPEED column of the above table are factory default values of user's parameters.

5.4 REPLACING THE CENTRIFUGE FUSES



The following fuses are required for Micro Star 30R:
2 x 16AT 250V

- Unplug mains plug from the mains socket.
- By pressing the locking device on the left side of the fuse compartment, the fuse holder is released and you can pull it out.
- Replace fuses.
- Insert fuse holder and push it, until it locks.

6. TECHNICAL DATA

Power supply:	230 V $\pm$ 10%, 50 / 60 Hz
Power consumption:	Centrifuge: 700 W Compressor: 350 W
Electric current:	5 A
Fuses:	2 x 16AT 250V
Refrigerant:	R-449A; 0.26 kg; 17 bar; GWP: 1397; CO ₂ e: 0.363 t
Protection class:	I
Rotational speed:	200 to 18,000 RPM
Maximum centrifugal force:	31,150 xg
Maximum load:	6 x 50 ml
Maximum kinetic energy:	7,300 Nm
Max. density of material to be centrifuged:	1.2 g/ml
Noise level at max. speed:	< 56 dB(A)
Run time:	10 s to 99 min 50 s, continuous operation (HOLD)
Number of programs:	100 programs
Acceleration:	levels from 0 to 9
Deceleration:	levels from 0 to 9 (0 - no braking)
Temperature of centrifugation:	-9 °C to 40 °C
Fast cool function (FC):	from 20 °C to 4 °C in approx. 5 minutes
Ambient temperature:	2 °C to 35 °C
Maximum relative humidity:	75%, non-condensing
Dimensions (W x D x H):	360 x 675 x 288 mm
Weight:	50 kg

We reserve the right to alter specification details without prior notice or liability!

7. CALCULATIONS

7.1 CALCULATION OF CENTRIFUGAL FORCE

For the calculation of the centrifugal force (RCF), stated as a multiple of the gravitational force "g", use the following formula:

$$\text{RCF} = 11.18 \times r \times (n / 1000)^2$$

RCF Relative centrifugal force (xg)

r Radius of the rotor (cm)

n Rotational speed (RPM)

7.2 CALCULATION OF MAXIMUM PERMITTED ROTOR SPEED

Users are responsible and must consider the limitations for maximum permitted rotor speed and about correct rotor load.

The maximum permitted speed for each type of rotor is marked on each rotor. It is defined for the use of samples with maximum density of 1.2 g/cm³.

If you need to use higher density samples, maximum permitted rotor speed must be reduced according to the following formula:

$$M = (1.2 \times n^2 / S)^{1/2}$$

M Reduced maximum permitted rotor speed

n Maximum permitted rotor speed for samples with density of 1.2 g/cm³

S Density of used sample

8. TECHNICAL SERVICE

Web Resources

Visit the VWR's website at www.vwr.com for:

- Complete technical service contact information.
- Access to VWR's Online Catalogue, and information about accessories and related products.
- Additional product information and special offers.

Contact us

For information or technical assistance contact your local VWR representative or visit www.vwr.com.

9. WARRANTY

VWR International warrants that this product will be free from defects in material and workmanship for a period of two (2) years from date of delivery. If a defect is present, VWR will, at its option and cost, repair, replace, or refund the purchase price of this product to the customer, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear. If the required maintenance and inspection services are not performed according to the manuals and any local regulations, such warranty turns invalid, except to the extent, the defect of the product is not due to such non-performance.

Items being returned must be insured by the customer against possible damage or loss. This warranty shall be limited to the aforementioned remedies. IT IS EXPRESSLY AGREED THAT THIS WARRANTY WILL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND IN LIEU OF THE WARRANTY OF MERCHANTABILITY.

10. COMPLIANCE WITH LOCAL LAWS AND REGULATIONS

The customer is responsible for applying for and obtaining the necessary regulatory approvals or other authorizations necessary to run or use the Product in its local environment. VWR will not be held liable for any related omission or for not obtaining the required approval or authorization, unless any refusal is due to a defect of the product.

11. EQUIPMENT DISPOSAL



This equipment is marked with the crossed out wheeled bin symbol to indicate that this equipment must not be disposed of with unsorted waste.

Instead it's your responsibility to correctly dispose of your equipment at lifecycle-end by handling it over to an authorized facility for separate collection and recycling. It's also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect from health hazards the persons involved in the disposal and recycling of the equipment.

For more information about where you can drop off your waste of equipment, please contact your local dealer from whom you originally purchased this equipment.

By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

Thank you



VWR International bvba
Haasrode Researchpark 2020, Geldenaaksebaan 464, B-3001 Leuven
<http://be.vwr.com>

EU Declaration of Conformity

Manufacturer: VWR International bvba, Haasrode Researchpark 2020, Geldenaaksebaan 464, B-3001 Leuven

The undersigned states that this declaration is issued under the sole responsibility of the manufacturer, VWR International bvba, and that the object(s) of the declaration described in annex 1 are in conformity with the relevant Union harmonization legislation.

In the event of unauthorized modification of any of the products listed in the annex 1, this declaration becomes invalid.

Object Name: VWR Laboratory centrifuges

Relevant EC Directives:	Low Voltage Directive	2014/35/EU as last amended
	EMC Directive	2014/30/EU as last amended
	Machinery Directive	2006/42/EC as last amended
	RoHS	2011/65/EU as last amended

Harmonized Standards:	EN 61010-1:2010
	EN 61010-2-020:2006
	EN 55011:2009
	EN 55014-1:2006
	EN 55014-2:2015
	EN 61000-3-2:2014
	EN 61000-3-3:2013
	EN 61000-6-1:2007
	EN 61326-1:2013
	EN ISO 12100:2010
EN 50581:2012	

Document Number: 28429.2018.01

Place & date of issue: Leuven, 25/06/2018

Signed for and on behalf of VWR International bvba:


Kris Nijs
Managing Director

Tel: 016 385 011 – Fax: 016 385 385 Email : vwr.be@vwr.com
BTW.TVA; BE0403593.343



VWR International bvba
Haasrode Researchpark 2020, Geldenaaksebaan 464, B-3001 Leuven
<http://be.vwr.com>

Annex 1

Cat. No.	Description
521-2574	Centrifuge Micro Star 30 R with rotor RA 24/2
521-2575	Centrifuge Micro Star 30 R without rotor

Tel: 016 385 011 – Fax: 016 385 385 Email : vwr.be@vwr.com
BTW.TVA; BE0403593.343



Your Distributor

Austria

VWR International GmbH
Graumannsgasse 7
1150 Vienna
Tel.: +43 01 97 002 0
Email: info.at@vwr.com

Belgium

VWR International bvba
Researchpark Haasrode 2020
Geldenaaksebaan 464
3001 Leuven
Tel.: +32 016 385 011
Email: vwr.be@vwr.com

China

VWR International China Co., Ltd
Shanghai Branch
Room 256, No. 3058 Pusan Road
Pudong New District
Shanghai 200123
Tel.: +86-21-5898 6888
Fax: +86-21-5855 8801
Email: info_china@vwr.com

Czech Republic

VWR International s. r. o.
Veetee Business Park
Pražská 442
CZ - 281 67 Stříbrná Skalice
Tel.: +420 321 570 321
Email: info.cz@vwr.com

Denmark

VWR International A/S
Tobaksvejen 21
2860 Søborg
Tel.: +45 43 86 87 88
Email: info.dk@vwr.com

Finland

VWR International Oy
Valimotie 9
00380 Helsinki
Tel.: +358 09 80 45 51
Email: info.fi@vwr.com

France

VWR International S.A.S.
Le Périgares – Bâtiment B
201, rue Carnot
94126 Fontenay-sous-Bois cedex
Tel.: 0 825 02 30 30* (national)
Tel.: +33 (0) 1 45 14 85 00
(international)
Email: info.fr@vwr.com
* 0,18 € TTC/min

Germany

VWR International GmbH
Hilpertstraße 20a
D - 64295 Darmstadt
Freecall: 0800 702 00 07
Tel.: +49 (0) 6151 3972 0
(international)
Email: info.de@vwr.com

Hungary

VWR International Kft.
Simon László u. 4.
4034 Debrecen
Tel.: +36 (52) 521-130
Email: info.hu@vwr.com

India

VWR Lab Products Private Limited
No.139, BDA Industrial Suburb,
6th Main, Tumkur Road, Peenya Post,
Bangalore, India – 560058
Tel.: +91-80-28078400
Email: vwr_india@vwr.com

Ireland / Northern Ireland

VWR International Ltd / VWR
International (Northern Ireland) Ltd
Orion Business Campus
Northwest Business Park
Ballycoolin
Dublin 15
Tel.: +353 01 88 22 222
Email sales.ie@vwr.com

Italy

VWR International S.r.l.
Via San Giusto 85
20153 Milano (MI)
Tel.: +39 02-3320311
Email: info.it@vwr.com

The Netherlands

VWR International B.V.
Postbus 8198
1005 AD Amsterdam
Tel.: +31 020 4808 400
Email: info.nl@vwr.com

Norway

VWR International AS
Haavard Martinsens vei 30
0978 Oslo
Tel.: +47 22 90 00 00
Email: info.no@vwr.com

Poland

VWR International Sp. z o.o.
Limbowa 5
80-175 Gdansk
Tel.: +48 058 32 38 200
Email: info.pl@vwr.com

Portugal

VWR International –
Material de Laboratório, Lda
Centro Empresarial de Alfragide
Rua da Indústria, nº 6
2610-088 Alfragide
Tel.: +351 21 3600 770
Email: info.pt@vwr.com

Singapore

VWR Singapore Pte Ltd
18 Gul Drive
Singapore 629468
Tel: +65 6505 0760
Email: sales.sg@vwr.com

Spain

VWR International Eurolab S.L.
C/ Tecnología 5-17
A-7 Llinars Park
08450 - Llinars del Vallès
Barcelona
Tel.: +34 902 222 897
Email: info.es@vwr.com

Sweden

VWR International AB
Fagerstagatan 18a
163 94 Stockholm
Tel.: +46 08 621 34 00
Email: info.se@vwr.com

Switzerland

VWR International AG
Lerzenstrasse 16/18
8953 Dietikon
Tel.: +41 044 745 13 13
Email: info.ch@vwr.com

UK

VWR International Ltd
Customer Service Centre
Hunter Boulevard
Magna Park
Lutterworth
Leicestershire
LE17 4XN
Tel.: +44 (0) 800 22 33 44
Email: uksales@vwr.com