[EN] Hematology controls

D0000910
Version: 19.0

In vitro diagnostic product, intended for professional use only.

Intended use

Clinical hematology laboratories require material for quality control of automated, semi-automated and manual procedures that measure whole blood parameters. J.T.Baker® Parameter or Retic Controls are hematology controls for these procedures. When handled like a patient sample and assayed on a properly calibrated and functioning instrument the control will provide values within the expected range indicated on the assay sheet. Daily use of these controls provides quality control data to confirm the precision and accuracy of instrument operation.

Reagents

Parameter Controls contain stabilized human RBC's, platelet components and fixed RBC's and or simulated WBC's for partial differential analysis to simulate WBC's, in a plasma like fluid with preservatives.

Product no., pack size and open vial stability

See page 5.

Instructions for use

- Remove the control from the refrigerator and allow vials to warm at room temperature (18 to 30°C) for 20 minutes before mixing.
- Place the control on a mechanical mixer for 20 minutes or follow next steps to mix the control manually. Important: do not place the vial on a Vortex-mixer.
- Hold the vial horizontally between the palms of the hands. Roll the vial back and forth for 30 seconds and gently invert the vial 10 times. Avoid vigorous shaking. Continue to mix in this manner until the cells are completely and uniformly suspended.
- After mixing let the vial rest undisturbed about 15 seconds to allow small air bubbles to disperse. Gently invert the vial 10 times immediately before sampling. Analyze the control using the same technique used for a patient sample.
- After sampling screw cap vials, carefully wipe the vial ring and cap with lint-free gauze and replace the cap immediately after cleaning.
- Place vials back in the refrigerator within 30 minutes after measuring the controls. Store in upright position.

Storage and stability

Parameter Controls are to be stored upright at 2-8°C when not in use. Stored at this temperature, the Controls are guaranteed stable until the expiry date.

Procedures

Instrument procedure: make dilutions and assay according to manufacturer's instructions for patient samples. Refer to assay values and range for the system in use.

Manual procedure: reference methods can be applied to 8 parameter and 3-Diff Controls. Refer to a manual of clinical laboratory procedures.

Expected results

The mean assay values and standard deviation for each Parameter Control are derived from replicate analyses on whole blood calibrated instrumentation as well as by manual reference methods. The values obtained on Parameter Controls prior to its expiry date should be within the expected range. The expected ranges listed represent estimates of instrument or inter-laboratory variation for each parameter. Inter-laboratory variation is usually accounted for by instrument calibration, maintenance and operating technique or reagent brand. For this reason, the assay values given are guide-numbers useful for control but are not absolute assays for calibration.

Values and expected ranges for instruments not listed on the Assay Information sheet must be established by the user. It is recommended that at least 5 consecutive analyses will be performed on a properly calibrated instrument for each level to establish the assay mean and standard deviation.

Warning and precautions

Warning: Potential bio hazardous material.

Parameter Controls are intended solely for use by trained, qualified personnel. Human blood components used in the Parameter Controls were found to be non-reactive for HBsAg and antibody to HIV when tested with licensed reagents. No known test methods can provide complete assurance that products derived from human blood will not transmit infectious diseases. Follow the same precautions as with patient samples when handling or disposing of vials. Do not inject or consume by mouth. Avoid direct mouth pipetting of samples.

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8-Parameter Control tube

3463	low	2.5 ml	4weeks
3464	normal	2.5 ml	4weeks
3465	high	2.5 ml	4weeks

8-Parameter Control VET

(for veterinary use only)

2904	low	2.5 ml	4weeks
2905	normal	2.5 ml	4weeks
2906	high	2.5 ml	4weeks

8 Parameter Control extended vial

3633	low	2.5 ml	4weeks
3634	normal	2.5 ml	4weeks
3635	high	2.5 ml	4weeks

3-Diff Control vial

3433	low	2.5 ml	3weeks
3434	normal	2.5 ml	3weeks
3435	high	2.5 ml	3weeks

3-Diff Control

3502	low	4.5 ml	3weeks
3503	Normal	4.5 ml	3weeks
3504	High	4.5 ml	3weeks

3-Diff Control extended

3421	low	2.5 ml	4weeks
3422	normal	2.5 ml	4weeks
3423	high	2.5 ml	4weeks

BC-Diff 5 Control

3613	low	3 ml	2weeks
3614	normal	3 ml	2weeks
3615	high	3 ml	2weeks

CD-Diff Control

3452	low	3 ml	2weeks
3453	normal	3 ml	2weeks
3454	high	3 ml	2weeks

K-Diff Control

3455	low	2.5 ml	3weeks
3456	Normal	2.5 ml	3weeks
3457	high	2.5 ml	3weeks

XE-Diff Control

3731	low	4.5 ml	2weeks
3732	Normal	4.5 ml	2weeks
3733	high	4.5 ml	2weeks

WBC Reduced RBC Control

3698	Low	3 ml	4weeks
3699	High	3 ml	4weeks

Effective date: 08 October 2020
Update date: 08 October 2020

Avantor Performance Materials Poland S.A.
ul. Sowińskiego 11, 44-101 Gliwice, Poland
Tel: +48 32 2392 000

The devices as mentioned in this sheet comply with
the In Vitro Medical Devices Directive 98/79/EC

In the tables below the product number, level, volume and minimum open-vial stability are listed.