
Laboratory Pump - VP 820

INSTRUCTION MANUAL

European Catalogue Numbers:

Euro Plug: 181-0071

UK Plug: 181-0073

Swiss Plug: 181-0072

Version: 2

Issued: 04.02.2013

Legal Address of Manufacturer

United States

VWR International
1310 Goshen Parkway
West Chester, PA 19380
800-932-5000
<http://www.vwr.com>

Europe

VWR International bvba
Researchpark Haasrode 2020
Geldenaaksebaan 464
B-3001 Leuven
+ 32 16 385011
<http://be.vwr.com>

Country of origin

Germany

Table of Contents

English

1. About this document.....	4
2. Use.....	4
3. Safety.....	5
4. Technical Data.....	6
5. Design and Function.....	7
6. Installation and connection.....	7
7. Operation.....	8
8. Servicing.....	9
9. Troubleshooting.....	11
10. Spare parts and accessories.....	12
11. Technical service.....	13
12. Warranty.....	13
13. Equipment disposal.....	13

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Package contents

- 1 Pump
- 1 Electric plug
- 1 Instruction manual

1. About this document

1.1 Using the Operating Instructions

The operating instructions are part of the pump.

- ➔ Carefully study the operating instructions before using a pump.
- ➔ Always keep the instruction manual handy in the work area.
- ➔ Pass on the operating instructions to the next owner.

1.2 Symbols and Markings

Warning

A danger warning is located here.

Possible consequences of a failure to observe the warning are specified here. The signal word, e.g. Warning, indicates the danger level.

WARNING ➔ Measures for avoiding the danger and its consequences are specified here.

Danger levels

Signal word	Meaning	Consequences if not observed
DANGER	warns of immediate danger	Death or serious injuries and/or serious damage are the consequence.
WARNING	warns of possible danger	Death or serious injuries and/or serious damage are possible.
CAUTION	warns of a possibly dangerous situation	Minor injuries or damage are possible.

Tab. 1

Other information and symbols

- ➔ An activity to be carried out (a step) is specified here.

1. The first step of an activity to be carried out is specified here. Additional, consecutively numbered steps follow.

This symbol refers to important information.

2. Use

2.1 Proper use

The pumps are exclusively intended for transferring gases and vapors.

Owner's responsibility

Operating parameters and conditions

Only install and operate the pumps under the operating parameters and conditions described in chapter 4, Technical data.

Make sure that the installation location is dry and the pump is protected against rain, splash, hose and drip water.

Requirements for transferred medium

Before using a medium, check whether the medium can be transferred danger-free in the specific application case.

Before using a medium, check the compatibility of the materials of the pump head, diaphragm and valves with the medium.

Only transfer gases which remain stable under the pressures and temperatures occurring in the pump.

Accessories

Laboratory equipment or additional components connected to a pump have to be suitable for use with the pneumatic capabilities of the pump.

2.2 Improper use

The pumps may not be operated in an explosive atmosphere.

The pumps are not suitable for transferring dusts.

The pumps are not suitable for transferring liquids.

The pumps must not be used to create vacuum and overpressure simultaneously.

An overpressure must not be applied to the suction side of the pump.

3. Safety

i Note the safety precautions in chapters 6. *Installation and connection*, and 7. *Operation*.

The pumps are built according to the generally recognized rules of technology and in accordance with the occupational safety and accident prevention regulations. Nevertheless, dangers can result during their use which lead to injuries to the user or others, or to damage to the pump or other property.

Only use the pumps when they are in a good technical and proper working order, in accordance with their intended use, observing the safety advice within the operating instructions, at all times.

Personnel

Make sure that only trained and instructed personnel or specially trained personnel work on the pumps. This especially applies to assembly, connection and servicing work.

Make sure that the personnel has read and understood the operating instructions, and in particular the "Safety" chapter.

Working in a safety-conscious manner

Observe the accident prevention and safety regulations when performing any work on the pump and during operation.

Do not expose any part of your body to the vacuum.

Open housing parts with notice sticker (see fig. 1) only after separating mains plug from power source.


Fig. 1: Notice sticker

Handling dangerous media

When transferring dangerous media, observe the safety regulations when handling these media.

Handling flammable media

Be aware that the pumps are not designed to be explosion-proof.

Make sure the temperature of the medium is always sufficiently below the ignition temperature of the medium, to avoid ignition or explosion. This also applies for unusual operational situations.

Note that the temperature of the medium increases when the pump compresses the medium.

Hence, make sure the temperature of the medium is sufficiently below the ignition temperature of the medium, even when it is compressed to the maximum permissible operating pressure of the pump. The maximum permissible operating pressure of the pump is stated in the technical specifications (see chapter 4).

If necessary, consider any external sources of energy, such as radiation, that may add heat to the medium.

In case of doubt, consult the VWR customer service.

Environmental protection

Store all replacement parts in a protected manner and dispose of them properly in accordance with the applicable environmental protection regulations. Observe the respective national and international regulations. This especially applies to parts contaminated with toxic substances.

English

Standards

The pumps conform to the Directive 2011/65/EU (RoHS2).

The pumps conform to the safety regulations of the EC Directive 2004/108/EC concerning Electromagnetic Compatibility and the EC Directive 2006/42/EC concerning Machinery.

The following harmonized standards have been used:

- DIN EN 61010-1
- DIN EN 61000-3-2/-3
- DIN EN 55014-1/-2

The pumps correspond to IEC 664:

- the overvoltage category II
- the pollution degree 2.

Customer service and repairs

Only have repairs to the pump carried out by the VWR Customer Service responsible.

Only authorized personnel should open those parts of the housing that contain live electrical parts.

Use only genuine parts from VWR for servicing work.

4. Technical Data

i All pumps are secured against overheating with thermal switches and are equipped with a mains fuse.

Pump materials

Component	Pump material*
Pump head	PTFE
Diaphragm	PTFE coated
Valves/Sealings	FFPM

Tab. 2

*according to DIN ISO 1629 and 1043.1

<i>Pneumatic performance</i>	
Max. permissible operating pressure [bar g]	1.0
Ultimate vacuum [mbar abs.]	≤ 8
Delivery rate at atm. pressure [l/min]*	max. 20
<i>Pneumatic connections</i>	
Hose connection [mm]	ID 10
<i>Ambient and media temperature</i>	
Permissible ambient temperature	+ 5 °C to + 40 °C
Permissible media temperature	+ 5 °C to + 40 °C
<i>Other parameters</i>	
Weight [kg]	9.3
Dimensions: L x H x W [mm]	312 x 207 x 144
Maximum permissible ambient relative humidity	80 % for temperatures up to 31 °C, decreasing linearly to 50 % at 40 °C
Max. altitude of site [m above sea level]	2000
<i>Electrical Data</i>	
Voltage [V]	230
Frequency [Hz]	50
Max. operating current [A]	0.7
Power consumption pump [W]	120
Maximum permitted mains voltage fluctuations	+/- 10 %

Electrical Data	
Fuse pump** (2x) T [A]	1.6
Protection class motor	IP44


Tab. 3

* Liters in standard state (1,013 mbar)

** For spare part-No. see chapter 10

5. Design and Function


Design

- | | |
|--------------------------|------------------------|
| 1 Connection piece | 5 Inlet (suction side) |
| 2 Pneumatic connection | 6 Power switch |
| 3 Pump head | |
| 4 Outlet (pressure side) | |

Fig. 2

Function Diaphragm Pump

- | | |
|--------------------|------------------|
| 1 Outlet valve | 4 Diaphragm |
| 2 Inlet valve | 5 Eccentric |
| 3 Transfer chamber | 6 Connecting rod |
| | 7 Pump drive |

Fig. 3: Pump head

Diaphragm pumps transfer, compress (depending on pump version) and evacuate gases and vapors.

The elastic diaphragm (4) is moved up and down by the eccentric (5) and the connecting rod (6). In the downward stroke it aspirates the gas to be transferred via the inlet valve (2). In the upward stroke, the diaphragm presses the medium out of the pump head via the outlet valve (1). The transfer chamber (3) is hermetically separated from the pump drive (7) by the diaphragm.

6. Installation and connection

Only install and operate the pumps under the operating parameters and conditions described in chapter 4, Technical data.

Observe the safety precautions (see chapter 3).

6.1 Installation

- ➔ Before installation, store the pump at the installation location to bring it up to room temperature.
- ➔ See chapter 4, Technical data, for the dimensions of pump.
- ➔ Install the pump so that the motor fan can intake sufficient cooling air.
- ➔ Make sure that the installation location is dry and the pump is protected against rain, splash, hose and drip water.
- ➔ Choose a safe location (flat surface) for the pump.
- ➔ Protect the pump from dust.
- ➔ Protect the pump from vibration and jolt.

6.2 Connection

- ➔ Only connect components to the pump which are designed for the pneumatic data of the pump (see chapter 4).

- If the pump is used as a vacuum pump, safely discharge the pump exhaust at the pump's pneumatic outlet.

i A marking on the pump head shows the direction of flow.

1. Remove the protective plugs from the pneumatic connectors of the pump.
2. Connect the suction line and pressure line.
3. Lay the suction and pressure line at a downward angle to prevent condensate from running into the pump.
4. Insert the power cable's plug into a properly installed shockproof socket.

7. Operation

7.1 Preparing for Start-up

Before switching on the pump, observe the following points:

	Operational requirements
Pump	▪ All hoses attached properly ▪ Fan openings not blocked ▪ Specifications of the power supply correspond with the data on the pump's type plate. ▪ The pump outlet is not closed or constricted.

Tab. 4

7.2 Starting

- Only operate the pump under the operating parameters and conditions described in chapter 4, Technical data.
- Make sure the pump is used properly (see chapter 2.1).
- Make sure the pump is not used improperly (see chapter 2.2).
- Observe the safety precautions (see chapter 3).

WARNING

Hazard of the pump head bursting due to excessive pressure increase

- Do not exceed max. permissible operating pressure (see chapter 4).
- Monitor pressure during operation.
- If the pressure exceeds the maximum permissible operating pressure, immediately shut down pump and eliminate fault (see chapter 9. Troubleshooting).
- Only throttle or regulate the air or gas quantity in the suction line to prevent the maximum permissible operating pressure from being exceeded.
- If the air or gas quantity in the pressure line is throttled or regulated, make sure that the maximum permissible operating pressure of the pump is not exceeded.

Pump standstill

- With the pump at a standstill, open pressure and suction lines to normal atmospheric pressure.

WARNING

Automatic starting can cause personal injury and pump damage

When the operation of the pump is interrupted by the thermal switch, the pump will restart automatically after cooling down.

- After triggering of the thermal protection or in the event of power failure, remove the pump's mains plug from the socket so that the pump cannot start uncontrollably.
- Attempt work on the pump only if the pump is separated from mains power.

7.3 Switching pump on and off

Switching pump on

i The pump may not start up against pressure during switch-on. If a pump starts against pressure, it may block. This activates the thermal switch, and the pump switches off.

- Make sure that no pressure is present in the lines during switch-on.

→ Switch on pump with mains switch (see fig. 2)

Switching off the pump/removing from operation

- When transferring aggressive media, flush the pump prior to switch-off to increase the service life of the diaphragm (see chapter 8.2.1).
- Switch off pump with mains switch (see fig. 2).
- Open pressure and suction lines to normal atmospheric pressure.
- Disconnect the power source.

8. Servicing

8.1 Servicing Schedule

Component	Servicing interval
Pump	Regular inspection for external damage or leaks
Diaphragm and valve plates/sealings	Replace at the latest, when pump output decreases

Tab. 5

8.2 Cleaning

i When cleaning, make sure that no liquids enter the inside of the housing.

8.2.1 Flushing Pump

- Before switching off the pump, flush it with air (if necessary for safety reasons: with an inert gas) for about five minutes under atmospheric conditions (ambient pressure).

8.2.2 Cleaning Pump

- Only clean pump with a damp cloth and non-flammable cleaning agents.

8.3 Changing Diaphragm and Valve Plates/Sealings

Conditions

- Pump is switched off and mains plug is removed from the socket
- Pump is clean and free of hazardous materials
- Tubes removed from pump's pneumatic inlet and outlet

Spare part/tool

Spare part/tool
Service Set (according to chapter 10)
Philips-head screwdriver No. 2
Soft pencil

Tab. 6

Information on procedure

- Always replace diaphragm and valve plates/sealings together to maintain the pump performance.

With multi-head pumps, parts of the individual pump heads can be confused.

- Replace the diaphragm and valve plates/sealings of the individual pump heads consecutively.

Health hazard due to dangerous substances in the pump!

Depending on the substance transferred, caustic burns or poisoning are possible.

WARNING

- Wear protective clothing if necessary, e.g. protective gloves.
- Flush pump before replacing the diaphragm and valve plates/sealings (see chapter 8.2.1).

Removing pump head

1. On the pneumatic head connection (3), loosen the union nuts (2) by hand. Then slightly loosen the angle-fitting (1) in the pump head (4) by turning it anticlockwise, so that the connecting tube can be pulled out.
2. Mark the position of top plate (fig. 7/5), head plate (fig. 7/6), intermediate plate (fig. 7/8) and adapter relatively to each other by a drawing line (for two-headed pumps: at both pump heads) with a felt-tip pen. This helps to avoid incorrect assembly later.
3. Loosen the outer screws (fig. 4/5) on the pump head/pump heads.
4. Carefully remove pump head / pump heads.

Change diaphragm

- i** Replace the diaphragms of both pump heads consecutively in order to ensure that the same number of diaphragm spacers is used.
1. Push down one diaphragm until other diaphragm is pushed upwards to its highest position.
 2. Carefully unscrew this diaphragm (1) anti-clockwise using both hands.
 3. Replace spacer thick (2) and spacers thin (3) onto the screw thread of the new diaphragm (same number and order).
 4. Screw in the new diaphragm and tighten it by hand.
 5. Complete steps 1 through 4 for the second pump head.

Change valve plates/sealings

- i** Replacing the valve plates/sealings of both pump heads consecutively.
1. At one pump head: Loosen screw(s) (1) in the center of the top plate (2).
 2. Remove top plate (2) and head plate (3) from intermediate plate (5). Valve plates/sealings (4) are visible
 3. Remove old valve plate/sealings.
 4. Clean intermediate plate (5) carefully (if there should be deposits in the recesses in the intermediate plate).
 5. Insert new valve plates/sealings (4) in the recesses in the intermediate plate (5).
 6. Carry out the steps 1 to 5 for the second pump head.
 7. Dispose of the old diaphragms and valve plates/sealings properly.

Refitting pump head

1. At one pump head: Apply pressure all around the edge of the diaphragm.
2. Place the intermediate plate (8) with valve plates/sealings on the adapter in accordance with the felt-tip pen marking.
3. Place the head plate (6) on the intermediate plate (8) in the position indicated by the guide pin (7).
4. Place the top plate (5) on the head plate (6) in accordance with the felt-tip pen marking.
5. Gently tighten screws (4) in diagonal order.


Fig. 4: Removing pump head


Fig. 5: Change diaphragm


Fig. 6: Change valve plates/sealings

6. Insert screw(s) (1) with disk springs (2, 3) in the center of the top plate (5). In doing so, make sure that the disk springs are arranged properly.
7. Screw in the screw/screws (1) in the centre of the pump top plate (5) until it is flush with the top plate (they are flush with the top plate); then screw one final half turn to tighten.
8. Carry out steps 1 to 7 for the second pump head.
9. Refit the pneumatic head connection: Place tube onto the connecting part of the angle fitting, turn angle fitting to a straight position and tighten the union nut.

Final steps

1. Reconnect system tubing.
2. Reconnect the pump to the electricity supply.


Fig. 7: Refitting pump head

9. Troubleshooting

Extreme danger from electrical shock!

- ➔ Disconnect the pump power supply before working on the pump.
- ➔ Make sure the pump is de-energized and secure.

DANGER

➔ Check the pump (see Tab. 7 to 10).

Pump produces no flow	
Cause	Fault remedy
No voltage in the power source	➔ Check room fuse and switch on if necessary.
Thermal switch has operated following to over-heating.	➔ Disconnect pump from mains. ➔ Allow pump to cool. ➔ Trace cause of over-heating and eliminate it.
Connections or lines blocked.	➔ Check connections and lines. ➔ Remove blockage.
External valve is closed or filter is clogged.	➔ Check external valves and filters.
Condensate has collected in pump head.	➔ Detach the condensate source from the pump. ➔ Flush pump (see chapter 8.2.1).
Diaphragm or valve plates/sealings are worn.	➔ Replace diaphragm and valve plates/sealings (see chapter 8.3).

Tab. 7

Flow rate, pressure or vacuum too low	
<i>The pump does not achieve the output specified in the Technical data or the data sheet.</i>	
Cause	Fault remedy
Condensate has collected in pump head.	➔ Detach the condensate source from the pump. ➔ Flush pump (see chapter 8.2.1).
There is gauge pressure on pressure side and at the same time vacuum or a pressure above atmospheric pressure on suction side.	➔ Change the pressure conditions.
Pneumatic lines or connection parts have an insufficient cross section.	➔ Disconnect pump from system to determine output values. ➔ Eliminate throttling (e.g. valve) if necessary. ➔ Use lines or connection parts with larger cross section if necessary.

Flow rate, pressure or vacuum too low <i>The pump does not achieve the output specified in the Technical data or the data sheet.</i>	
Cause	Fault remedy
Leaks occur on connections, lines or pump head.	→ Check that tubes sit correctly on hose nozzles. → Replace leaky tubes. → Eliminate leaks.
Connections or lines completely or partially jammed.	→ Check connections and lines. → Remove the jamming parts and particles.
Head parts are soiled.	→ Clean head components.
Diaphragm or valve plates are worn.	→ Replace diaphragm and valve plates (see chapter 8.3).
After diaphragm and valve plates have been replaced	→ Check that the spacers have been replaced onto the diaphragm screw thread. → Check head connection and hose connections for leaks. → Possibly carefully tighten the outer screws (fig. 7/4) of the top plate crosswise.

Tab. 8

Pump is switched on, but does not run, the on/off-switch on the pump is not lit	
Cause	Fault remedy
Pump is not connected with the power source.	→ Connect pump to mains power.
No voltage in the power source	→ Check room fuse and switch on if necessary.
Fuse in the pump is defective.	→ Remove pump's mains plug from the socket. → Loosen marked lid on underside of the pump. → Select and replace suitable fuse (see chapter 4).

Tab. 9

Pump is switched on, but does not run, the on/off-switch on the pump is lit	
Cause	Fault remedy
The thermal switch has opened due to overheating	→ Remove pump's mains plug from the socket. → Allow pump to cool. → Trace cause of over-heating and eliminate it.

Tab. 10

Fault cannot be rectified

If you are unable to determine any of the specified causes, send the pump to VWR Customer Service.

1. Flush the pump to free the pump head of dangerous or aggressive gases (see chapter 8.2.1).
2. Clean the pump (see chapter 8.2.2).
3. Send the pump to VWR with a filled out decontamination declaration (see appendix) and specification of the medium transferred.

10. Spare parts

A Service Set contains all spare parts needed for one complete service:
2 diaphragms and 4 valve plates/sealings.

Service Set for pump type:	Order-No.:
VP 820	057358

Tab. 11

11. Technical service

Web Resources

Visit the VWR's website at www.vwr.com for:

- Complete technical service contact information
- Access to VWR's Online Catalogue, and information about accessories and related products
- Additional product information and special offers

Contact us For information or technical assistance contact your local VWR representative or visit www.vwr.com.

12. Warranty

VWR International warrants that this product will be free from defects in material and workmanship for a period of two (2) years from date of purchase. If a defect is present, VWR will, at its option, repair, replace, or refund the purchase price of this product at no charge to you, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear.

For your protection, items being returned must be insured against possible damage or loss. This warranty shall be limited to the replacement of defective products.

IT IS EXPRESSLY AGREED THAT THIS WARRANTY WILL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND IN LIEU OF THE WARRANTY OF MERCHANTABILITY.

13. Equipment disposal

This equipment is marked with the crossed out wheeled bin symbol to indicate that this equipment must not be disposed of with unsorted waste.

Instead it's your responsibility to correctly dispose of your equipment at lifecycle -end by handling it over to an authorized facility for separate collection and recycling. It's also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect from health hazards the persons involved in the disposal and recycling of the equipment.

For more information about where you can drop off your waste of equipment, please contact your local dealer from whom you originally purchased this equipment.

By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

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Molte grazie!

Decontamination Declaration

i The condition for the repair of a pump by VWR is the certification of the customer on the transferred media and on the cleaning of the pump (decontamination declaration).

- ➔ Copy this page.
- ➔ Enter the pump model, the Serial No. and the transferred media in the form below and sent the signed form together with the flushed and cleaned pump to VWR Customer Service.

Europe

VWR International bvba
Researchpark Haasrode 2020
Geldenaaksebaan 464
B-3001 Leuven
+ 32 16 385011

<http://be.vwr.com>

Customer decontamination declaration for repair order

We herewith confirm that the following media have been pumped with the pump listed below, and that the pump has been flushed and cleaned.

Pump model	
Serial No.	
Fed media	

The pump contains neither aggressive, biological, radioactive, poisonous nor other dangerous media.

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Ireland

VWR International Ltd
Orion Business Campus
Northwest Business Park
Ballycoolin
Dublin 15
Tel.: 01 88 22 222
Fax: 01 88 22 333
e-mail sales@ie.vwr.com

Northern Ireland

VWR International Ltd
A10 Harbour Court, 7 Heron Rd
Sydenham Business Park
Belfast BT3 9HB
Tel.: 028 9058 5800
Fax: 028 9080 7812
e-mail sales@ie.vwr.com

Italy

VWR International s.r.l.
Via Stephenson 94
20157 Milano (MI)
Tel.: 02 332 03 11
Fax: 800 152 999
E-mail: info@it.vwr.com

The Netherlands

VWR International B.V.
Postbus 8198
1005 AD Amsterdam
Tel.: 020 4808 400
Fax: 020 4808 480
E-mail: info@nl.vwr.com

Norway

VWR International AS
Kakkelovnskroken 1
P.B. 45, Kalbakken
0901 Oslo
Tel.: 02290
Fax: 815 00 940
E-mail: info@no.vwr.com

Portugal

VWR International - Material de
Laboratório, Lda
Edifício Neopark
Av. Tomás Ribeiro, 43- 3 D
2790-221 Carnaxide
Tel.: 21 3600 770
Fax: 21 3600 798/9
E-mail: info@pt.vwr.com

Spain

VWR International Eurolab S.L.
Ronda Can Fatjó, nº 11
Edifici Tecnopark, 3
Parc Tecnològic del Vallés
08290 Cerdanyola del Vallés
(Barcelona)
Tel.: 902 222 897
Fax: 902 430 657
E-mail: info@es.vwr.com

Sweden

VWR International AB
Fagerstagatan 18a
163 94 Stockholm
Tel.: 08 621 34 00
Fax: 08 621 34 66
E-mail: info@se.vwr.com

Switzerland

VWR International AG
Lerzenstrasse 16/18
8953 Dietikon
Tel.: 044 745 13 13
Fax: 044 745 13 10
E-mail: info@ch.vwr.com

UK

VWR International Ltd
Customer Service Centre
Hunter Boulevard
Magna Park
Lutterworth
Leicestershire
LE17 4XN
Tel.: 0800 22 33 44
Fax: 01455 55 85 86
E-mail: uksales@uk.vwr.com