## **VWR® Thermometer**

### **INSTRUCTION MANUAL**

**North American Catalog Numbers: 61161-283, 46610-024,  
15551-004, 12777-842, 89172-082, 89172-084**

**European Catalogue Numbers: 620-2017, 620-2073, 620-2092,  
620-2099, 620-1826**

Version: 1  
Issued: 3<sup>rd</sup> November 2014

## Legal Address of Manufacturer

### United States

VWR International  
Radnor Corporate Center  
Building One, Suite 200  
100 Matsonford Road, P.O. Box 6660  
Radnor, PA 19087  
800-932-5000  
<http://www.vwr.com>

### Europe

VWR International bvba  
Researchpark Haasrode 2020  
Geldenaaksebaan 464  
B-3001 Leuven  
+ 32 16 385011  
<http://be.vwr.com>

Country of Origin – Hong Kong

### INTENDED USE

This item is intended for general laboratory use to monitor temperature.

### VWR TRACEABLE® STAINLESS-STEEL PROBES THERMOMETER INSTRUCTIONS - ENGLISH

### SPECIFICATIONS

| Cat. No.                       | Range                           | Resolution                                     | Accuracy                           | Probe                                                                                 |
|--------------------------------|---------------------------------|------------------------------------------------|------------------------------------|---------------------------------------------------------------------------------------|
| <b>61161-283/<br/>620-2017</b> | –58 to 572°F<br>(–50 to 300°C)  | 0.1° from –<br>20° to 200°,<br>1° otherwise    | ±1.0°C<br>between<br>–20 and 100°C | stainless steel-probe<br>diameter 0.14 inch, length 5<br>inches                       |
| <b>46610-024/<br/>620-2073</b> | –58 to 572°F<br>(–50 to 300°C)  | 0.1° from –<br>20° to<br>200°, 1°<br>otherwise | ±0.5°C at<br>tested points         | stainless steel-probe<br>diameter 0.14 inch, length 5<br>inches                       |
| <b>15551-004/<br/>620-2092</b> | –58 to 572°F<br>(–50 to 300°C)  | 0.1° from –<br>20° to 200°,<br>1° otherwise    | ±1.0°C<br>between<br>–20 and 100°C | Stainless-steel 0.14-inch;<br>stem length, 8 inches; overall<br>length, 10-1/4 inches |
| <b>12777-842<br/>620-2099</b>  | –58 to 572°F<br>(–50 to 300°C)  | 0.1° from –<br>20° to<br>200°, 1°<br>otherwise | ±0.4°C at<br>tested points         | Stainless-steel 0.14-inch;<br>stem length, 8 inches; overall<br>length, 10-1/4 inches |
| <b>89172-082/<br/>620-1826</b> | –22 to 122°F and<br>–30 to 50°C | 0.1°                                           | ±1.0°C                             | Sealed bottle Probe                                                                   |
| <b>89172-084</b>               | –22 to 122°F and<br>–30 to 50°C | 0.1°                                           | ±0.4°C at<br>tested points         | Sealed bottle Probe                                                                   |

## ULTRA™ THERMOMETER ACCURACY

Ultra™ thermometers are tested at selected test points to be within tighter than normal tolerances to assist in providing improved accuracy. Other points will not necessarily fall within the same accuracy as those found at the selected test points, but will be within an accuracy of  $\pm 1^{\circ}\text{C}$ .

**CAT. NO. 61161-283/620-2017; 46610-024/620-2073 (CABLE WITH STAINLESS-STEEL PROBE)**

**CAT. NO. 15551-004/620-2092; 12777-842/620-2099 (STAINLESS-STEEL PROBE)**

**CAT. NO. 89172-082/620-1826; 89172-084 (SEALED BOTTLE PROBE FOR REFRIGERATOR/FREEZER APPLICATIONS)**

## OPERATION

1. Press the ON/OFF button to turn the unit on. Press the °F/°C button to switch between Fahrenheit and Celsius.
2. Insert the probe into the material to be measured and read the display.
3. Press the ON/OFF button to turn the unit off. To conserve battery life, always turn the unit off when not in use.

## MAX/MIN MEMORY

1. To view the minimum temperature reached since activating the unit, press the MAX/MIN button. "MIN" appears on the display to indicate the minimum temperature recorded.
2. Press the MAX/MIN button a second time, within 3 seconds, to view the maximum temperature reached since activating the unit. "MAX" appears on the display to indicate the maximum temperature recorded.

**Note:** Three seconds after pressing the MAX/MIN button, the unit automatically returns to the current temperature reading. ("MIN" and "MAX" are no longer displayed.)

3. Turn the unit off to reset the minimum and maximum memories.

**Cat. No. 89172-082/620-1826; 89172-084 IMPORTANT: The intended use for this thermometer is to monitor temperatures inside refrigerators and freezers.**

The entire unit is completely sealed and waterproof. Do not attempt to remove the thermometer from the bottle.

Velcro® and double sided foam tape are supplied and may be used to mount the thermometer or the base side to any surface. In a refrigerator application, this prevents the thermometer from taking up valuable shelf space.

## **Cat. No. 89172-082/620-1826; 89172-084 USING THE MEMORY TO MONITOR A REFRIGERATOR /FREEZER**

Following is an example of how to use the memory to monitor the temperature inside a refrigerator or freezer. This example is provided only as a helpful guide and is not intended to replace existing facility requirements or procedures.

## Unit Setup Example

1. Turn the unit on.
2. Select the desired unit of measure.
3. Place the thermometer into the supplied base.
4. Place the thermometer and base inside the refrigerator/freezer.

At this point, allow sufficient time for the bottle thermometer to reach equilibrium with the true current temperature inside the refrigerator/freezer.

5. Once the thermometer has reached equilibrium, clear the minimum and maximum temperature memory. (See the “Minimum/Maximum Memory” section.)

The memory will provide a record of the single lowest and highest temperature achieved. When the temperatures are recorded into the manual log, the memory will allow the user to see if the temperature inside the refrigerator/freezer has gone outside of the acceptable range.

## MONITORING PROCEDURE EXAMPLE

Keep a notebook or spreadsheet as a manual log.

1. At the same time every day, record the following into the manual log:
  - Current Date and Time
  - Current Temperature Reading
  - Minimum Temperature Reading (MIN)
  - Maximum Temperature Reading (MAX)
2. Once the above items have been manually recorded, clear the minimum and maximum temperature memory. (See the “Minimum/Maximum Memory” section.)

By clearing the memory each day, the minimum and maximum temperature memory will provide a record of the minimum and maximum temperature that has been achieved inside the refrigerator/freezer over the past 24 hour monitoring period. The memory will also allow the user to see if the temperature inside the refrigerator went outside of the acceptable range.

## ALL OPERATIONAL DIFFICULTIES

If this thermometer does not function properly for any reason, replace the battery with a new high quality battery (see the “Battery Replacement” section). Low battery power can occasionally cause any number of “apparent” operational difficulties. Replacing the battery with a new fresh battery will solve most difficulties.

## BATTERY REPLACEMENT

Erratic readings, a faint display, no display, or [==] appearing on the display are all indications that the battery must be replaced. To replace the battery, remove the three battery cover screws located on the back of the unit. Remove the battery cover. Replace the exhausted battery with a new 1.5-volt, silver oxide, #357 size battery. Make certain the positive (+) side is visible. Replace the battery cover and battery cover screws. Do not overtighten the screws as this may damage the unit.

**NOTE:** When replacing the battery cover, ensure the screws are tightened securely to maintain the unit's water-proof seal.

## Technical service

### Web Resources

Visit the VWR's website at [www.vwr.com](http://www.vwr.com) for:

- Complete technical service contact information
- Access to VWR's Online Catalogue, and information about accessories and related products
- Additional product information and special offers

Contact us for information or technical assistance contact your local VWR representative or visit. [www.vwr.com](http://www.vwr.com).

### Warranty

VWR International warrants that this product will be free from defects in material and workmanship for a period of two (2) years from date of purchase. If a defect is present, VWR will, at its option, repair, replace, or refund the purchase price of this product at no charge to you, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear. For your protection, items being returned must be insured against possible damage or loss. This warranty shall be limited to the replacement of defective products. IT IS EXPRESSLY AGREED THAT THIS WARRANTY WILL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND IN LIEU OF THE WARRANTY OF MERCHANTABILITY.

**GO TO VWR.COM FOR THE LATEST NEWS, SPECIAL OFFERS AND DETAILS OF YOUR LOCAL VWR DISTRIBUTOR**

### Equipment disposal

This equipment is marked with the crossed out wheeled bin symbol to indicate that this equipment must not be disposed of with unsorted waste. Instead it's your responsibility to correctly dispose of your equipment at lifecycle -end by handling it over to an authorized facility for separate collection and recycling. It's also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect from health hazards the persons involved in the disposal and recycling of the equipment. For more information about where you can drop off your waste of equipment, please contact your local dealer from whom you originally purchased this equipment.

By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

### Compliance with local laws and regulations

The customer is responsible for applying for and obtaining the necessary regulatory approvals or other authorizations necessary to run or use the Product in its local environment. VWR will not be held liable for any related omission or for not obtaining the required approval or authorization, unless any refusal is due to a defect of the product.

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# Local VWR offices in Europe and Asia Pacific

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## India

VWR Lab Products Private Limited  
135/12, Brigade Towers, 2nd Floor  
Front wing, Brigade Road,  
Bengaluru, India - 560 025  
Tel: +91-2522-647911/922 (Mumbai)  
Tel: +91-80-41117125/26 (Bangalore)  
Fax +91-80-41117120  
E-mail: vwr\_india@vwr.com

## Ireland / Northern Ireland

VWR International Ltd / VWR International  
(Northern Ireland) Ltd  
Orion Business Campus  
Northwest Business Park  
Ballycoolin  
Dublin 15  
Tel.: 01 88 22 222  
Fax: 01 88 22 333  
E-mail sales@ie.vwr.com

## Italy

VWR International PBI S.r.l.  
Via San Giusto 85  
20153 Milano (MI)  
Tel.: 02-3320311/02-487791  
Fax: 02-332031307/02-40090010  
E-mail: info@it.vwr.com

## The Netherlands

VWR International B.V.  
Postbus 8198  
1005 AD Amsterdam  
Tel.: 020 4808 400  
Fax: 020 4808 480  
E-mail: info@nl.vwr.com

## New Zealand

VWR International LP  
241 Bush Road  
Albany  
0632 New Zealand

## Norway

VWR International AS  
Haavard Martinsens vei 30  
0978 Oslo  
Tel.: 02290  
Fax: 815 00 940  
E-mail: info@no.vwr.com

## Poland

VWR International Sp. z o.o.  
Limbowa 5  
80-175 Gdansk  
Tel.: 058 32 38 210  
Fax: 058 32 38 205  
E-mail: labart@pl.vwr.com

## Portugal

VWR International - Material de  
Laboratório, Lda  
Edifício Neopark  
Av. Tomás Ribeiro, 43- 3 D  
2790-221 Camaxide  
Tel.: 21 3600 770  
Fax: 21 3600 798/9  
E-mail: info@pt.vwr.com

## Singapore

VWR Singapore Pte Ltd  
18 Gul Drive  
Singapore 629468  
Tel: +65 6505 0760  
Fax: +65 6264 3780  
E-mail: sales@sg.vwr.com

## Spain

VWR International Eurolab S.L.  
C/ Tecnología 5-17  
A-7 Llinars Park  
08450 - Llinars del Vallès  
Barcelona  
Tel.: 902 222 897  
Fax: 902 430 657  
E-mail: info@es.vwr.com

## Sweden

VWR International AB  
Fagerstagatan 18a  
163 94 Stockholm  
Tel.: 08 621 34 00  
Fax: 08 621 34 66  
E-mail: info@se.vwr.com

## Switzerland

VWR International GmbH  
Lerzenstrasse 16/18  
8953 Dietikon  
Tel.: 044 745 13 13  
Fax: 044 745 13 10  
E-mail: info@ch.vwr.com

## Turkey

Pro Lab Laboratuvar Teknolojileri Ltd.Şti.  
a VWR International Company  
Orta Mah. Cemal Gürsel Caddesi  
Ördekcioglu İşmerkezi No.32/1  
34896 Pendik - İstanbul  
Tel.: +90216 598 2900  
Fax: +90216 598 2907  
Email: info@pro-lab.com.tr

## UK

VWR International Ltd  
Customer Service Centre  
Hunter Boulevard  
Magna Park  
Lutterworth  
Leicestershire  
LE17 4XN  
Tel.: 0800 22 33 44  
Fax: 01455 55 85 86  
E-mail: uksales@uk.vwr.com