Instruction Manual

VWR® E-BURETTE Motor Operated Burette with 3 Calibrated Pre-set Speeds

EU cat. no 613-7002

NA cat. no 76319-574
76319-576

CE

UK
CA

Table of Contents

EN

Instruction Manual

01 - 33

FR

Mode d'emploi

67 - 99

DE

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34 - 66

Legal address of Manufacturer

Europe

VWR International bv,
Researchpark Haasrode 2020
Geldenaaksebaan 464 B-3001 Leuven
+3216385011
be.vwr.com

UK Importer

VWR International Ltd
Hunter Boulevard, Magna Park
Lutterworth, Leicestershire, LE17 4XN
uk.vwr.com

United States

VWR International LLC
100 Matsonford Rd
Radnor, PA 19087
+1 800-932-5000
www.vwr.com

Country of Origin

India

Contents

Safety information
Package contents
Unpacking
Installation
Intended use
Symbols and conventions
Product specifications
Instructions for use
Troubleshooting
General maintenance
Accessories and spares
Technical service
Warranty
Disposal
Appendix

Safety Instructions

VWR E-Burette may sometimes be used for operations involving hazardous materials and equipment. It is beyond the scope of this manual to address all of the potential risks associated with its use in such applications. It is the sole responsibility of the user to consult and establish appropriate safety and health practices and determine the applicability of regulatory limitations, prior to the usage of this instrument.

Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

1. Follow general and safety instructions for hazard prevention. E.g., always wear protective clothing, protective gear for the eyes and hands.
2. Carefully observe the specifications provided by reagent manufacturers.
3. When dispensing inflammable fluids, avoid the built up of static charge. Make sure that you do not dispense into plastic vessels and do not wipe the instrument/equipment with a dry cloth.
4. Use the instrument only for dispensing liquids, with strict regard to the defined limitations of use and operating constraints.
5. Observe operating exclusions. When in doubt, contact the manufacturer or supplier.
6. Always use the instrument in such a way that neither the user nor any other person is in danger. While dispensing, the discharge tube must always point away from you or any other person. Only dispense into appropriate vessels and avoid splashes.
7. Do not press the piston when the discharge tube closure is attached.
8. Do not remove the discharge tube while the dispensing cylinder is being filled.
9. Reagents can accumulate in the cap of the discharge tube. Thus, clean it regularly.
10. Do not carry the mounted instrument by the cylinder sleeve or the valve block. Breakage or loosening of the cylinder may lead to personal injury.
11. Use only original accessories and spare parts.
12. Do not attempt to make any technical alterations. Do not dismantle the instrument any further than is described in the Operation Manual.
13. Always check the instrument for visual damage before use.
14. If there is any sign of a potential malfunction (e.g. piston difficult to move, stuck valve or leakage), immediately stop dispensing. Consult the 'Troubleshooting' section of this Operation Manual and contact the manufacturer if needed.

Package Contents

Description	Quantity
E-Burette-50 ml	1
Telescopic tube	1
Control Panel	1
Bottle Adapters	5
Dispenser Thread Reducer	1
Charging AC Adaptor	1
Control Panel Stylus	1
USB Drive for Software & Driver	1
Mini USB Cable	1
Extra Amber Colored Window	1
14.4 V Li-ion Rechargeable Battery	1
Calibration Certificate	1
Operation Manual	1

Unpacking

- Connecting the batteries

1. Open the battery panel at the back of the control panel.
(Fig. 1)

Fig. 1

2. Connect the batteries to the panel wire by inserting the battery connector into the panel connector firmly. It locks with a 'click' sound.
(Fig. 2)

Fig. 2

3. Close the panel.
(Fig. 3)

Fig. 3

4. Press the ON/OFF button.
(Fig. 4)

Fig. 4

5. The screen will display 'Welcome'.
(Fig. 5)

Fig. 5

NOTE:

In order to replace the batteries, press the lock button to disengage the battery lock and remove the old batteries.

Installation

1. Adjust the length of the Telescopic tube (FEP tube) to fit your particular reservoir. If you require a longer tube, it is provided on request. (Fig. 6)


Fig. 6

2. Fix the Telescopic tube as shown. The tube's ends have different diameters – fix it from the wider end.

(Fig. 7)

Fig. 7

3. Choose the correct adapter for the bottle.

The threaded platform base of E-BURETTE has a 30 mm screw thread. Five adapters are supplied to suit containers with a 28, 32, 38, 40, 45 mm and 30 mm (inbuilt adapter) screw neck.

(Fig. 8)

Fig. 8

4. Fix the adapter on the reservoir bottle by screwing it in clockwise direction. (Fig. 9)

Fig. 9

5. Mount the instrument. Screw it to the reservoir by applying gentle hand torque applied to the threaded platform base only. In case of removal, apply the same technique to the base, in the opposite direction.

6. Connect the cable from the instrument housing to the port at the back of the control panel with a firm push, as shown. Check for the cable position. The flat surface of the connector with an arrow mark will face down.

(Fig. 10)

Fig. 10

7. If charging is required, insert the charging pin at the back of the control panel and connect the power cord with the AC adaptor firmly and then plug the two pin power cord to the 220 V power supply.
(Fig. 11)

Fig. 11

8. The battery indicator will show the charging symbol.
(Fig. 12)

Fig. 12

9. The E-Burette is now ready to be switched on. Press the ON/OFF button at the back of the control panel to switch ON the control panel.
(Fig. 13)


Fig. 13

Transporting the instrument

- When mounted to a reagent bottle, always carry the instrument as shown and always store it in the upright position.
(Fig. 14)

Fig. 14

Computer Connectivity

Insert the pen drive provided with the package in your PC/Laptop for driver files and software.

Please note that there are two set up files in the pen drive.

- CDM21228
- E-Burette Setup Files

Procedure to install the driver & software:

- Installing the CDM21228 driver files:

1. Double click the CDM21228 setup icon.
(Fig. 15)


Fig. 15

2. Press the Extract tab from the FTDI CDM driver window.
(Fig. 16)


Fig. 16

3. Click to accept the licence agreement.
(Fig. 17)


Fig. 17

4. Then, click the Finish button on the device driver installation wizard.
(Fig. 18)


Fig. 18

- Installing the instrument software

1. Double click the E-BURETTE setup icon.
(Fig. 19)


Fig. 19

2. Click Next to start the Install Shield Wizard for Communication window.
(Fig. 20)


Fig. 20

3. Click Next to choose the destination folder of the software files.
(Fig. 21)


Fig. 21

4. Click Install for the ready to install program window.
(Fig. 22)


Fig. 22

5. Click Next for the user account control window.
(Fig. 23)


Fig. 23

6. Click Finish for the Install Shield Wizard completion window.
(Fig. 24)


Fig. 24

7. Communication icon will then appear on the desktop.
(Fig. 25)

Fig. 25

Intended Use

VWR E-BURETTE is a state-of-the-art Motor Operated Burette with sophisticated features and functionality. It is intended for use in laboratories and the field for easy, safe and accurate titrations.

Symbols and Conventions

CAUTION This symbol indicates a potential risk and alerts you to proceed with caution

Product Specifications

Performance Data	50 ml
Unit Type	Microprocessor-controlled motorized piston lift system with recirculation valve
Volume Ranges	0.01 ml – 500 ml (Max piston lift 50 ml, automatic addition)
Volume Accuracy	A= 0.05% / CV = 0.05%
Increment	10 µl
Operating Temperature	+10 °C to +40 °C (from 50 °F to 104 °F)
Power Supply	18 V, 3 A (Output supply : 110 V, 230 V)
Battery Operation	One complete battery charge gives 3 hours (approx.) of continuous operation, (Battery capacity indicated on TFT display and charging time 6.5 hrs.)
Quality Control	In compliance with DIN EN ISO 8655-3
Touchscreen	4" TFT Screen
PC Interface	Mini USB Cable
Languages	Available in 4 languages (English, Spanish, German and French)

Error Limits

The error limits (Accuracy and Coefficient of Variation) mentioned in Product Specifications are in accordance with the nominal capacity (or maximum volume) indicated on the instrument. These are obtained by using the instrument with distilled water at equilibrium, ambient temperature of 20 °C, while operating the device smoothly and steadily.

The error limits are well within the limits of DIN EN ISO 8655-3.


Instructions for Use

Overview

VWR E-BURETTE is designed for dispensing liquids directly from the reservoir bottle.

It is calibrated in accordance with the guidelines of the DIN EN ISO 8655 – 3.

When the instrument is correctly used, the dispensed liquid comes into contact with only the following chemically resistant materials: PTFE, FEP and Borosilicate glass.

Getting Started

Switching on the control panel

1. Switch on the instrument using the ON/OFF switch at the back of the control panel.
2. When it boots, the instrument will warn the user to switch the knob to recirculation mode.

(Fig. 26)

Fig. 26

3. Turn the knob to re-circulation mode and click the tick mark (✓) on the control panel screen when ready.

(Fig. 27)

Fig. 27

4. The instrument will reset automatically. It is now ready to use.

Purging

- This step is recommended before using the Burette to ensure bubble free dispensing.

1. Click on the purging mode from the Home Screen on the control panel. (Fig. 28)


Fig. 28

2. The screen on the control panel will warn the user to turn the knob to re-circulation mode.
(Fig. 29)

Fig. 29

3. Turn the knob to re-circulation mode and click the tick mark (✓) on the control panel screen when ready.

(Fig. 30)

Fig. 30

4. The E-BURETTE will purge automatically. If the device is still not properly purged, go to step 1 and repeat the procedure till no large air bubble is visible below the piston.

- A few air bubbles up to 1 mm in size are permissible.

Important Guidelines

Physical Limits

This E-BURETTE is designed for titrating liquids, observing the following physical limits:

- Maintain instrument and reagent between +10 °C to +40 °C (from 50 °F to 104 °F).
- Maintain reagent vapor pressure up to 600 mbar (maximum).
Aspirate slowly above 300 mbar in order to prevent the liquid from boiling.
- Kinematic viscosity upto 500 mm²/s.
(dynamic viscosity [mPas] = kinematic viscosity [mm²/s] x density [g/cm³])
- Relative air humidity between 20% and 90%.
Use fluids with density upto 2.2 g/cm³.

Operating Limitations

- Chlorinated and fluorinated hydrocarbons or chemical combinations which form deposits may make it difficult for the piston to move or jam it. If the piston becomes too difficult to move, immediately clean the instrument.
- While dispensing inflammable fluids, avoid the built up of static charge. Make sure that you do not dispense into plastic vessels and do not wipe instruments with a dry cloth.
- E-BURETTE is designed for general laboratory applications and complies with relevant standards, like DIN EN ISO 8655-3. Please check the compatibility of the instrument for a specific application (e.g. trace material analysis, food sector etc.). Approvals for specific applications, like production and administration of food, pharmaceuticals and cosmetics are not available.

Operating Exclusions

Do not use the instrument with:

- Liquids that attack FEP, PFA and PTFE (e.g. dissolved Sodium Azide*).
- Liquids that attack borosilicate glass (e.g. Hydrofluoric acid).
- Hydrochloric acid > 40% | Tetrahydrofuran or Oxolane | Trifluoroacetic acid.
- Explosive liquids (e.g. Carbon Disulfide).
- Suspensions (e.g. of charcoal) as the solid particles may clog or damage the instrument.
- Liquids that attack PP (cap).
- An aggressive atmosphere, e.g. HCl fumes.

* The permissible concentration of dissolved Sodium Azide is not more than 0.1%.

Storage Conditions

Store the instrument and accessories in a clean, cool and dry place. Storage temperature from -20 °C to +50 °C (from -4 °F to 122 °F) is advised with relative humidity between 5% and 95%.

Sensitive Media

For light-sensitive media (e.g., Iodine, Potassium Permanganate and Silver Nitrate solutions), we recommend the use of 'Amber Coloured Windows' or the inspection windows provided in the box.

1. To replace the inspection window, unclip the default inspection window by pressing it from its top and then removing it.

(Fig. 31)

Fig. 31

2. Place the Amber Coloured Window in the slot provided at the bottom of the housing and press it down until you hear a click sound.

(Fig. 32)

Fig. 32

General Operating Instructions

Control panel brightness setting

- User can change brightness from setting screen on the control panel.

1. Click on the settings icon from the home screen.

(Fig. 33)

Fig. 33

2. In the settings windows, click the BRIGHTNESS icon.

(Fig. 34)

Fig. 34

3. From the BRIGHTNESS window, select the level of brightness and tick mark to save it.

(Fig. 35)


Fig. 35

Language Settings

- User can select the preferred language from the language settings.

1. Click on the settings icon from the home screen

(Fig. 36)

Fig. 36

2. Select the preferred language from the tab by clicking on it and tick mark (✓) to save it.

E-BURETTE supports four languages English, German, Spanish and French.

(Fig. 37)

Fig. 37

Operation

Titration

Wear protective clothing, protective gear for the eyes and hands. Liquid may accumulate in the cap. Follow all safety instructions and observe limitations of use and operating constraints.

1. Remove cap from the discharge tube.

(Fig. 38)

Fig. 38

2. Click on the Titration mode on the home screen.

(Fig. 39)

Fig. 39

3. The screen on the control panel will warn the user to turn the knob to titration mode.

4. Turn the knob clockwise to set at titration mode and click the tick (✓) mark on the screen of the control panel when ready.

(Fig. 40)

Fig. 40

5. Place the discharge tube orifice against the inner wall of a suitable receiving vessel.

6. Now fill the barrel from the titration screen by clicking the FILL button (Fig. 41 - left part). By pressing the STOP button, the filling volume can be stopped as and when you need to.

(Fig. 41 - right part)

Fig. 41

7. Three dispensing speeds are provided, including drop-wise dispensing, which allows the user to achieve the end point very accurately.

- **Fast mode**

In this mode, the dispensing speed is 3.8 ml/sec.

(Fig. 42)

Fig. 42

- **Medium mode**

In this mode, the dispensing speed is 1.3 ml/sec
(Fig. 43)

Fig. 43

- **Dropwise mode**

In this mode, the dispensing speed is 0.010 ml/click (Fig. 44)

Fig. 44

8. Dispensing can be performed by touching the screen and stops immediately when the touch is removed.

9. Press the TARE button to reset to Zero & repeat.
(Fig. 45)

Fig. 45

Data Transfer Setup

- **Connecting the USB Mini Cable**

1. Connect the Mini USB cable from the control panel to your PC/Laptop for transferring and storing the control panel's data on the latter. (Fig. 46)

Fig. 46

2. Insert the wider end of the USB cable in the USB port of the PC/Laptop and the smaller end at the back of the control panel.
(Fig. 47)

Fig. 47

● Software Initialization

1. Double click on the Communication icon to run the software.
(Fig. 48)

Fig. 48

2. From the Ports List button, select the port from the drop down menu.
(Fig.49)


Fig. 49

3. Click on OPEN PORT.
(Fig. 50)


Fig. 50

4. The window will capture the model number, serial number and the volume range of the E-BURETTE.
(Fig.51)


Fig. 51

5. Click ADD from the control panel to save the readings from the control panel to the PC/Laptop.
(Fig. 52)

Note:

The ADD button shall blink and the data shall be transferred to the PC/Laptop.

Fig. 52

6. The readings will also be saved in the internal memory. Please note that all fields are captured in the file.

Memory

1. Click ADD from the control panel to save the readings from the control panel to the internal memory.
2. To recall the previous recorded reading (s), press the RECALL button.
(Fig. 53)

Fig. 53

3. Please note that the recall grid displays the last 20 recorded readings. The LAST saved reading will be displayed on the number 1 position of the table.
(Fig. 54)

Fig. 54

4. For quick recall of the last reading, press the LAST button.
(Fig. 55)

Fig. 55

Please note that the last saved readings will be displayed till the time LAST button is pressed.

Creating a File Name

1. Select the port from the port list dropdown menu.
2. Click Open port.
3. The Software captures the Model No., Sr. no. and volume range.
4. After transferring the data from control panel to the software window we can save the file by clicking the "Save as Excel File".
(Fig. 56)


Fig. 56

Printing the data

1. The data recorded while operating the E-BURETTE can be printed directly from the software by clicking the PRINT tab.
(Fig. 57)
2. For this, select the printer and print the readings.
You can also access the saved .xls file to print.


Fig. 57

Refresh

1. To start afresh for new readings on the software, click on the REFRESH tab.
(Fig. 58)
2. The readings on the grid shall be deleted and new readings can now be transferred in it.


Fig. 58

Exit

1. To end working with the E-BURETTE, press the Exit button
(Fig. 59).
2. The panel will warn the user to turn the knob in the re-circulation mode
3. Click on the tick (✓) button to proceed. The piston will empty the barrel and go to the home position.
Press the ON/OFF button to switch OFF the control panel.

Fig. 59

Troubleshooting

Problem	Possible cause	Solution
Piston difficult to move	Formation of crystals or deposition of dirt	Perform a cleaning cycle.
Filling not possible	Filling valve stuck	Clean the filling valve. If the valve ball is stuck, use a 200 µl plastic tip to loosen it. (Fig. 60)  Fig. 60
Air bubbles in the instrument	Purging not complete Filling tube loose or damaged Filling tube does not dip into the liquid	Purge the instrument again. Fasten the telescopic filling tube firmly. If necessary, cut the tube off by approx. 1 cm from the top or replace it. Fill up the bottle, or correctly adjust the length of the telescopic filling tube.
The volume delivered is smaller than that indicated	The instrument has not been completely primed	Prime the instrument again.
The instrument doesn't indicate any function	Internal error	Perform a restart.

General Maintenance

Whenever cleaning is required, run the dispenser under distilled water. Fill the barrel completely with distilled water and press the fast dispense button to dispense the water completely.

- Procedure to Disassemble the Delivery Pipe:

1. Unscrew the chuck nut by rotating it in anticlockwise direction and pull out the delivery pipe.

(Fig. 61)

Fig. 61

2. Clean the delivery pipe with de-ionized water.

- Procedure to Assemble the Delivery Pipe:

1. First push the delivery pipe into the lower housing till it stops going in further.

(Fig. 62)

Fig. 62

2. Screw the chuck nut to complete the assembly by turning it in clockwise direction.

(Fig. 63)

Fig. 63

Note: Use rubbing alcohol on a cloth or cotton to clean the external body.

User Calibration

VWR E-BURETTE has been laboratory calibrated at its nominal volume. However, due to changes in environmental conditions and the viscosity of the media which is dispensed, we recommend gravimetric testing every 3-12 months. Gravimetric volume testing according to DIN EN ISO 8655-6 is performed as follows:

1. Click the settings icon on the bottom right of the Home Screen.

(Fig.64)

Fig. 64

2. From the settings menu, choose Calibration.

(Fig. 65)

Fig. 65

3. Fill the nominal volume of double distilled de-ionized water by clicking the FILL button on the calibration screen of the operating panel.

(Fig. 66)


Fig. 66

4. Dispense the filled liquid by clicking on the DISPENSE button on the calibration screen of the control panel.

(Fig. 67)

Fig. 67

5. Measure the dispensed liquid on a balance and enter the value in mg using the provided electronic keyboard. Click the tick mark (✓) to proceed.

Note: CAL label will start appearing on top of all screens.

6. Repeat this procedure until the nominal volume is achieved on the electronic balance.

Factory Reset

- The factory reset is for the Manufacturer factory settings which is password protected and intended to be used at the manufacturing facility.

(Fig. 68)

Fig. 68

- If the reset button is pressed, the user calibrated settings shall be reset with the default factory settings & the CAL label will disappear from the screens.

Accessories and Spares

Descripton	Quantity	NA Cat. No.	EU Cat. No.
VWR ADAPTER ASSORTMENT DISENSER PK4	1	76319-594	612-6865
VWR ADAPTER BOTTLE TOP DSPN 38MM PK3	1	76319-602	612-6866
VWR ADAPTER BOTTLE TOP DSPN 45MM PK3	1	76319-598	612-6867
VWR ADAPTER BOTTLE TOP DSPN 32MM PK3	1	76319-596	612-6868
VWR ADAPTER BOTTLE TOP DSPN 40MM PK3	1	76319-600	612-6869
VWR ADAPTER BOTTLE TOP DSPN 28MM PK3	1	76319-604	612-6870
VWR ADAPTER BOTTLE TOP DSPN 25MM PK3	1	76319-606	612-6871

Technical service

Web Resources

Visit the VWR website at www.vwr.com for:

- Complete technical service contact information
- Access to the VWR Online Catalogue, and information about accessories and related products
- Additional product information and special offers

Contact us For information or technical assistance contact your local VWR representative or visit www.vwr.com

Warranty

VWR warrants that this product will be free from defects in material and workmanship for a period of two (2) years from date of delivery. If a defect is present, VWR will, at its option and cost, repair, replace, or refund the purchase price of this product to the customer, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear. If the required maintenance and inspection services are not performed according to the manuals and any local regulations, such warranty turns invalid, except to the extent, the defect of the product is not due to such non-performance.

Items being returned must be insured by the customer against possible damage or loss. This warranty shall be limited to the aforementioned remedies. IT IS EXPRESSLY AGREED THAT THIS WARRANTY WILL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND IN LIEU OF THE WARRANTY OF MERCHANTABILITY

Compliance with Local Laws and Regulations

The customer is responsible for applying for and obtaining the necessary regulatory approvals or other authorisations necessary to run or use the Product in its local environment. VWR will not be held liable for any related omission or for not obtaining the required approval or authorisation, unless any refusal is due to a defect of the product.

Equipment Disposal

This equipment is marked with the crossed out wheeled bin symbol to indicate that this equipment must not be disposed of with unsorted waste. Instead it's your responsibility to correctly dispose of your equipment at lifecycle -end by handling it over to an authorized facility for separate collection and recycling.

It's also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect from health hazards the persons involved in the disposal and recycling of the equipment.

For more information about where you can drop off your waste of equipment, please contact your local dealer from whom you originally purchased this equipment.

By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

Thank you

Appendix

Chemical Compatibility Chart

Chemicals from A to Z

The following tables enlist the most frequently used reagents.

They provide useful information for the safe and efficient use of dispenser.

All the safety precautions and recommendations within this Operation Manual must be followed carefully.

Code explanations:

A = Good resistance

B = Acceptable with limitations

C = Not recommended

1 = Acid vapours (better resistance with lower concentration).

Do not leave the instrument on the bottle.
Rinse with distilled water after use.

2 = Risk of damage, softening or discoloration of external parts through vapours.

Do not leave the instrument on the bottle.
Rinse with distilled water after use.

3 = Chemical degradation of glass parts (plunger/barrel).

List of Reagents

Chemicals A - Z	
A	
Acetaldehyde (Ethanal)	A
Acetic acid 96%	A
Acetic acid 100% (glacial)	B/2
Acetic anhydride	B/2
Acetone (Propanone)	B/2
Acetonitrile (MECN)	A
Acetophenone	B/2
Acetyl Chloride	B/2
Acetylacetone	A
Acrylic acid	A
Acrylonitrile	B/2
Adipic acid	A
Allyl alcohol	A
Aluminum chloride	A
Amino acids	A
Ammonia 20%	B/2
Ammonia 20-30%	B/2
Ammonium chloride	A
Ammonium fluoride	A
Ammonium molybdate	A
Ammonium sulfate	A
Amyl alcohol (Pentanol)	A
Amyl chloride (Chloropentane)	B/2
Aniline	A
Ascorbic acid	A
n-Amyl acetate	B/2
B	
Barium chloride	A
Benzaldehyde	A
Benzene	B/2
Benzine	A
Benzoyl chloride	B/2
Benzyl alcohol	A
Benzyl chloride	B/2
Bis(2-ethylhexyl) phthalate	B/2
Boric acid 10%	A
Bromine	C/2
Bromobenzene	B/2
Bromonaphtalene	A
Butanediol	A
Butanol	A
Butyl acetate	B/2
Butyl methyl ether	B/2
Butylamine	B/2
Butyric acid	B/2
C	
Calcium carbonate	A
Calcium chloride	A
Calcium hydroxide	A
Calcium hypochlorite	A

List of Reagents

Chemicals A - Z	
C	
Carbon disulfide	B/2
Carbon tetrachloride	B/2
Chlorine dioxide	B/2
Chlorine water	B/2
Chloro naphthalene	B/2
Chloroacetaldehyde 45%	A
Chloroacetic acid	A
Chloroacetone	B/2
Chlorobenzene	B/2
Chlorobutane	B/2
Chloroethanol	B/2
Chloroform (Trichloromethane)	B/2
Nitro-hydrochloric acid (Aqua regia)	B/2
Chlorosulfonic acid	B/2
Chlorosulfuric acid 100%	B/1/2
Chromic acid 100%	B/1/2
Chromosulfuric acid 100%	C/1/2
Citric acid	A
Copper fluoride	A
Copper sulfate	A
Covi-Ox-T70/ Mixed Tocopherol	A
Cresol	A
Cumene (Isopropylbenzene)	B/2
Cyanoacrylate	A
Cyclohexane	B/2
Cyclohexanone	B/2
Cyclopentane	B/2
D	
1,2-Diethylbenzene	B/2
1,4-Dioxane (Diethylene dioxide)	B/2
1-Decanol	A
Decane	A
Di-(2-ethylhexyl) peroxydicarbonate	B/2
Dibenzyl ether	B/2
Dichloroacetic acid	A
Dichlorobenzene	A
Dichloroethane	A
Dichloroethylene	B/2
Diesel oil (Heating oil)	A
Diethanolamine	A
Diethylamine	B/2
Diethylene glycol	A
Diethylether	B/2
Dimethylacetamide	A
Dimethyl sulfoxide (DMSO)	B/2
Dimethylaniline	A
Dimethylformamide (DMF)	B/2

List of Reagents

Chemicals A - Z	
E	
Ethanol	A
Ethanolamine	B/2
Ether	B/2
Ethyl acetate	B/2
Ethylbenzene	B/2
Ethylene chloride	B/2
Ethylene diamine	A
Ethylene glycol	A
F	
Fluoroacetic acid	B/2
Formaldehyde (Formalin)	A
Formamide	A
Formic acid	A
G	
Gamma-butyrolactone	A
Gasoline	B/2
Glycerin <40%	A
Glycolic acid 50%	A
H	
Heating oil (Diesel oil)	A
Heptane	A
Hexane	A
Hexanoic acid	A
Hexanol	A
Hydriodic acid	B/2
Hydrobromic acid	A
Hydrochloric acid 20% (HCl)	A
Hydrochloric acid 37% (HCl)	B/1
Hydrofluoric acid (HF)	C/3
Hydrogen peroxide	A
I	
Iodine	A
Iodine bromide	C/2
Iodine chloride	C/2
Isoamyl alcohol	A
Isobutanol	A
Isooctane	A
Isopropanol	A
Isopropyl ether	B/2
Iso-propylamine	B/2
K	
Kerosene	A
L	
Lactic acid	A
M	
2-Methoxyethanol	A
Methanol	A
Methoxybenzene (Anisol)	B/2
Methyl benzoate	B/2
Methyl chloride (Chloromethane)	B/2

List of Reagents

Chemicals A - Z	
M	
Methyl ethyl ketone (MEK/Butanone)	B/2
Methyl formate	A
Methyl iodide (Iodomethane)	B/2
Methyl methacrylate (MMA)	B/2
Methyl propyl ketone (2-Pentanone)	A
Methyl tert-butyl ether	B/2
Methylene chloride (Dichloromethane) (DCM)	B/2
Methylpentanone	A
Mineral oil (engine oil)	A
Monochloroacetic acid	A
N	
Nitric acid 100%	A
Nitric acid 30-70%	A
Nitric acid dil. <30%	A
Nitrobenzene	B/2
Nitromethane	B/2
N-methyl-2-pyrrolidone (NMP)	A
O	
Octane	A
Octanol	A
Oil (vegetable, animal)	B/2
Oil of turpentine	B/2
Oleic acid	A
Oleum (Fuming Sulfuric acid)	A
Oxalic acid	A
P	
Pentane	B/2
Peracetic acid	A
Perchloric acid 100%	B/2
Perchloric acid diluted	A
Perchloroethylene	B/2
Petroleum	B/2
Petroleum ether / spirit	B/2
Phenol	A
Phenylethanol	B/2
Phenylhydrazine	B/2
Phosphoric acid 100%	A
Phosphoric acid 85%	A
Piperidine	B/2
Potassium chloride	A
Potassium dichromate	A
Potassium dihydrogen phosphate	A
Potassium hydroxide	A
Potassium iodide	A
Potassium permanganate (persulfate)	A
Potassium peroxydisulfate	A
Potassium sulfate	A
Propionic acid (Propanoic acid)	A
Propylene glycol (Propane-1,2-diol)	A

List of Reagents

Chemicals A - Z	
Propylene oxide	A
Picric acid (Trinitrophenol)	B/2
Pyridine	B/2
Pyruvic acid	A
R	
Resorcin	A
S	
Salicylaldehyde	A
Scintillation fluid	A
Silver acetate	A
Silver nitrate	A
Sodium acetate	A
Sodium chloride (kitchen salt)	A
Sodium dichromate	A
Sodium fluoride	A
Sodium hydroxide 30%	A
Sodium hypochlorite	A
Sodium thiosulfate	A
Sulfonitric acid 100%	B/2
Sulfur dioxide	B/2
Sulfuric acid 100%	B/2
Sulfuric acid <10%	A
Sulfuric acid (10-75%)	B/1
Sulfuric acid (Cold conc.)	A
Sulfuric acid (Hot conc.)	B/2
T	
1,1,2-Trichlorotrifluoroethane	B/2
Tartaric acid	A
Tetrachlorethylene	B/2
Tetrahydrofuran (THF)	B/2
Tetramethylammonium hydroxide	A
Toluene	B/2
Trichlorethylene	B/2
Trichloroacetic acid	B/2
Trichlorobenzene	B/2
Trichloroethane	B/2
Triethanolamine	A
Triethylamine	A
Triethylene glycol	A
Trifluoroacetic anhydride (TFAA)	B/2
Trifluoromethane (Fluoroform)	B/2
U	
Urea	A
X	
Xylene	B/2
Z	
Zinc chloride 10%	A
Zinc sulfate 10%	A

Our local offices in Europe, AMEA and North America

Austria
VWR International GmbH
Graumannsgasse 7
1150 Vienna
Tel.: +43 1 97 002 0
Email: info.at@vwr.com

Belgium
VWR International BV
Researchpark Haasrode 2020
Geldenaaksebaan 464
3001 Leuven
Tel.: +32 016 385 011
Email: vwr.be@vwr.com

Canada
VWR International
2360 Argentia Road
Mississauga, Ontario L5N 5Z7
Tel.: +1 800 932 5000
Email: Canada_Orders@vwr.com

China
VWR (Shanghai) Co., Ltd
Bld.No.1, No.3728 Jinke Rd,
Pudong New District
Shanghai, 201203- China
Tel.: +400 821 8006
Email: info_china@vwr.com

Czech Republic
VWR International s. r. o.
Veetee Business Park
Pražská 442
CZ - 281 67 Strábrná Skalce
Tel.: +420 321 570 321
Email: info.cz@vwr.com

Denmark
VWR International A/S
Tobaksvejen 21
2860 Søborg
Tel.: +45 43 86 87 88
Email: info.dk@vwr.com

Finland
VWR International Oy
Valimotie 17-19
00380 Helsinki
Tel.: +358 9 80 45 51
Email: info.fi@vwr.com

France
VWR International S.A.S.
Immeuble Estréo
1-3 rue d'Aurion
93114 Rosny-sous-Bois cedex
Tel.: 0 825 02 30 30* (national)
Tel.: +33 (0) 1 45 14 85 00 (international)
Email: info.fr@vwr.com
* 0,18 € TTC/ min + prix appel

Germany
VWR International GmbH
Hilpertstrasse 20a
D - 64295 Darmstadt
Freecall: 0800 702 00 07
Tel.: +49 (0) 6151 3972 0 (international)
Email: info.de@vwr.com

Hungary
VWR International Kft.
Simon László u. 4.
4034 Debrecen
Tel.: +36 (52) 521-130
Email: info.hu@vwr.com

India
VWR Lab Products Private Limited
No.139, BDA Industrial Suburb,
6th Main, Tumkur Road,
Peenya Post,
Bangalore, India – 560058
Tel.: +91-80-28078400
Email: vwr_india@vwr.com

Ireland
VWR International Ltd
Orion Business Campus
Northwest Business Park
Ballycoolin
Dublin 15
Tel.: +353 01 88 22 222
Email sales.ie@vwr.com

Italy
VWR International S.r.l.
Via San Giusto 85
20153 Milano (MI)
Tel.: +39 02-3320311/02-487791
Email: info.it@vwr.com

Korea
Avantor Performance Materials Korea Ltd
2F ACE Gwanggyo Tower I
Daehak 4ro 17
Yeongtong-gu Suwon
Korea 16226
Tel.: +82 31 645 7250
saleskorea@avantorsciences.com

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VWR International AB
VWR International, LLC