



CaptairBio

PCR workstations

Complete protection for
RNA/DNA amplification

SAFETY

Protection of RNA/
DNA samples. UV
safety switched with
SMART technology
for guaranteed safety
from harmful UV rays

PERFORMANCE

H14 HEPA filtration
that meets EN 1822
standards.

ADAPTABILITY

Increase sample pro-
tection with the op-
tion of adding carbon
filtration, combining
molecular and partic-
ulate filtration.


SIMPLICITY

Delivered completely
knocked down (CKD)
for ease of installation
in any setting. Sets up
in minutes.

CONNECTIVITY


SMART technology
for real-time perfor-
mance monitoring
ensuring peak oper-
ation during critical
handlings.





Life in the laboratory becomes simpler and safer

CaptairBio PCR Workstations are designed to protect sensitive biological applications against both environmental pollution and cross-contamination. They feature a high efficiency filtration system that provides a particle-free workstation around the manipulation. High energy UV light is used to decontaminate the worktop from biological cross-contamination between two operations.



Particulate free workstation*

- Protection against external contamination
- Internal air quality achieved by high efficiency particulate filter(s) (HEPA H14 / ULPA U16)
- Carbon filter (optional) to protect handlings from VOCs present in the laboratory atmosphere



UV decontamination

- Protect your samples from cross-contamination
- Powerful UV decontamination (254 nm lamp power)
- Adjustable timer
- Automatic UV lamp of switch in case the sash is opened when the UV light is on



Easy to clean

- Work surface is easy to clean
- Seamless worktop with smooth corners (available in TRESPA®TopLab^{PLUS} or Stainless steel (304L)
- Low porosity material




Ergonomic design

- 3 models available for your handlings with large openings for easy access to your work
- Slanted sash provides an ergonomic position for comfort and productivity
- High luminosity, internal LED lighting (daylight, light intensity > 800 lux)

Simpler to use*

- SMART Technology informs users about their protection using light and sound.
- Light and sound pulses provide real time information indicating that:




- The sash is open when the UV light is on.
- Air face velocity is compromised: check sash, pre-filter or HEPA / ULPA filter.
- Fan failure has occurred.

Safer to operate*

- ULPA U16 filters guarantee 99.99995% filtration efficiency for particles larger than 0.1µm.
- HEPA H14 filters guarantee 99.995% filtration efficiency for particles larger than 0.1µm.
- Add a molecular filter for additional protection from fumes and vapors present in the laboratory air.
- The anemometer monitors a drop in pressure indicating that pre-filter or filter replacement is required.


320*



AVT Cat. #52428-824
Supplier #321PCR



AVT Cat. #52428-826
Supplier #391PCR



Model	320*	52428-824	52428-826		
Safety standards	NF EN 61010 - CE Marking	NF EN 61010 - CE Marking - EN 1822:1998 (HEPA H14 & ULPA U16 Filters)			
External width (in-mm)	(31 ^{7/8}) 810	(31 ^{7/8}) 810	(39 ^{5/8}) 1006		
External depth (in-mm)	(24 ^{1/4}) 615	(24 ^{1/4}) 615	(24 ^{1/4}) 615		
External height min-max (in-mm)	(28 ^{1/8}) 714	(38 ^{1/8} - 42) 968-1067	(38 ^{1/8} - 42) 968-106730		
Internal width (in-mm)	(30 ^{7/8}) 784	(30 ^{7/8}) 784	(38 ^{1/8}) 968		
Internal depth min-max (in-mm)	(19 ^{5/8} - 20 ^{1/8}) 498-511	(19 ^{5/8} - 20 ^{1/8}) 498-511	(19 ^{5/8} - 20 ^{1/8}) 498-511		
Internal height (in-mm)	(22 ^{1/4}) 565	(23 ^{3/8}) 594	(23 ^{3/8}) 594		
Voltage / Frequency (V-Hz)	220-240 / 50-60 or 100-110 / 60	100-240 / 50-60			
Air flow (m³/h-CFM)	-	1P	1C1P	1P	1C1P
		200 / 118	245 / 144	200 / 118	245 / 144
Power consumption (Watts)	25	40	45	40	55
Decibel level (dBA)	< 40	54	57	55	57
Side and front panels	10 mm thick synthetic glass is designed to protect users from harmful UV rays and β (Bêta) emitted from radioactive isotopes such as: T(3H), 14C, 32P				
Structure	Corrosion resistant electro-galvanized steel coated with anti-acid polymer				
Filtration module	Polypropylene				

Filtration

Model	320*	52428-824	52428-826
Particulate filter (1P)	-	HEPA H14: This filtration technology traps particles larger than 0.1 µm with 99.995% efficiency according to the MPPS method set forth in the EN 1822-1 standard. ULPA U16: This filtration technology traps particles larger than 0.1 µm with 99.99995% efficiency according to the MPPS method set forth in the EN 1822-1 standard.	
Carbon filter (optional) (1C)	-	Adding a carbon filter to your enclosure allows protection of your samples from VOCs. AS filter : For organic vapors	
Particulate pre-filter	-	Protects particulate filters from dust in the laboratory environment (only for 1P version)	

Features

Model	320*	52428-824	52428-826
Worktop	Stainless steel 304 L / TRESPA® TopLab PLUS		
Bactericidal UV Lights	15W - Wavelength : 254 nm		
	0.08 mJ/ s/cm²	0.08 mJ/ s/cm²	0.13 mJ/ s/cm²
Internal lighting	Compact fluorescent lighting - 18 Watts - IP 67	LED - IP 44-6000K	
	500 Lux	900 lux	950 lux
Connectivity	RJ45 Cable connection to view and change workstation settings (cable included)		
Anemometer	Monitors a drop in pressure that indicates pre-filter replacement is required.		

Accessories

Model	320*	52428-824	52428-826
Benches	Rolling cart (Mobicap) or fixed bench (Benchcap)		
Molecode S	-	Automatic detection of VOC filter breakthrough	

*320 does not include air flow