VWR

Stereomicroscope  series 300

INSTRUCTION MANUAL

Model	European Catalogue Number
SZT360-6	630-2982

Version: 1
Issued: 08, 03, 2019

Legal Address of Manufacturer

Europe

VWR International bvba

Researchpark Haasrode 2020

Geldenaaksebaan 464

B-3001 Leuven

+ 32 16 385011

<http://be.vwr.com>

Country of origin: **ITALY**

Table of Contents

1. Warning	3
2. Safety Information	3
3. Package Contents	3
4. Unpacking	3
5. Intended use	3
6. Symbols and conventions	3
7. Product Specifications	3
8. Overview	4
9. Assembling	5
9.1 Assembling procedure	5
10. Use of the microscope	6
10.1 Adjust interpupillary distance	6
10.2 Focus the specimen	6
10.3 Adjusting the tension of the focus knob	6
10.4 Dioptic compensation	6
10.5 Magnification	7
10.6 Use of additional lens	7
10.7 Use of contrast disc	8
10.8 Use of illumination	8
10.9 Use of eyeshields	8
10.10 Microphotography	8
11. Troubleshooting	9
12. Repair and maintenance	10
13. Technical service	10
14. Warranty	10
15. Compliance with local laws and regulations	10

1. Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use.

We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users.

We decline any responsibility deriving from incorrect instrument use that does not comply with this manual.

2. Safety Information

Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.

Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

3. Package Contents

DESCRIPTION	QUANTITY
Stereo head 0.67x - 4.5x	1
Stand with double arm LED illuminator	1
Eyepieces WF10x	2
External power supply 12V	1

4. Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (head and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

5. Intended use

For research use only. Not intended for any animal or human therapeutic or diagnostic use.

6. Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

	CAUTION This symbol Indicates a potential risk and alerts you to proceed with caution
--	---

7. Product Specifications

Model	Head	Eyepieces	Objectives	Working Distance	Stand	Illumination
SZT360-6	Trinocular, 360° rotating, 45° inclined	Extra Wide Field 10x/23	0.67.....4.5x	110 mm	Pillar with rack and pinion system	Transmitted: 100 mm LED disc (60 chip LEDs) Incident: double arm P-LED ³ .

8. Overview

9. Assembling

9.1 Assembling procedure

1. Insert the microscope into the head holder and lock the clamping knob ①. (Fig. 1)
2. Remove the eyepieces dust caps and insert the eyepieces in the empty eyepiece sleeve. (Fig. 2)
3. Lock the eyepieces using the locking knob ②. (Fig. 3)
4. Connect the jack of the power supply to the socket in the rear part of the microscope base. (Fig. 4)

10. Use of the microscope

10.1 Adjust interpupillary distance


Hold the right and left eyepiece tube with both hands and adjust the interpupillary distance by moving the two parts until one circle of light can be seen (Fig. 5).

- Should two circles appear, the interpupillary distance is too big.
- Should two overlapped circles appear, the interpupillary distance is too small.

10.2 Focus the specimen


Put the sample to be observed on the stage plate and focus the sample using the focusing knobs of the stand ①. If necessary adjust the height of the microscope head along the vertical stand. Remember to lock the lock screw ② after setting the height of the microscope. (Fig. 6)

10.3 Adjusting the tension of the focus knob

Grab the knobs ③ with both hands and, holding right knob rotate the left one. Tension increases or decreases according to the sense of rotation of the left focus knob. (Fig. 7)


- **This adjustment allows to increase or decrease the tension of the knob by avoiding an involuntary descent of the microscope body under its own weight. Adjust the tension just above the point where the focus is stable.**

10.4 Dioptic compensation

This compensation makes possible for people wearing glasses to adjust the microscope to their eyes and use the microscope without glasses.

1. Put the zoom knob ④ down to the lowest magnification and focus the specimen with the focusing knobs ①. (Fig. 8)
2. Put the zoom knob to the maximum magnification and repeat the focusing.
3. Return to the lowest magnification: the specimen will be out of focus.

- Adjust the diopter compensation ring of the right eyepiece ⑤ until the image of the right eyepiece is clear and sharp. Repeat the procedure for the left eyepiece ⑥. (Fig. 9)
- Then, check the focus of the image for the whole zoom range. It should now be perfectly parfocal (focus is always maintained during the change of magnification).

10.5 Magnification

Select the desired magnification by adjusting the zoom knob ④. (Fig. 8)

The microscope body is equipped with a “click-stop” function that allows to obtain a precise setting of the desired magnification.

Click stop can be activated or deactivated by operating with the provided allen wrench in the hole placed in the right side of the microscope body. (Fig. 10)

To increase or decrease the total magnification, change the eyepieces and/or add an appropriate additional lens if necessary.

10.6 Use of additional lens

Screw the desired additional lens on the microscope body. (Fig. 11)

Each additional lens has a specific Working Distance.

The stroke of the focus adapter cannot compensate the different working distances of the several additional lenses.

Should the focus adapter not be able to focus properly the specimen, move the whole microscope body accordingly.

Total magnification used can be calculated as:

Eyeiece magnification * Zoom magnification * Objective lens magnification

Eyeiece	10x		15x		20x		25x	
Field Number (mm)	23		16		12		9	
Objective	Total mag.	FOV (mm)	Total mag.	FOV (mm)	Total mag.	FOV (mm)	Total mag.	FOV (mm)
0.3X	2.01x-13.5x	114.42-17.04	3.02x-20.25x	52.98-7.90	4.02x-27x	29.85-4.44	5.02x-33.75x	17.93-2.67
0.5X	3.35x-22.5x	68.66-10.02	5.02x-33.75x	31.87-4.74	6.7x-45x	17.91-2.67	8.37x-56.25x	10.75-1.6
0.75X	5.02x-33.75x	45.82-6.81	7.54x-50.62x	21.22-3.16	10.05x-67.5x	11.94-1.78	12.56x-84.37x	7.17-1.07
1X	6.7x-45x	34.32-5.11	10.05x-67.5x	15.24-2.37	13.4x-90x	8.96-1.33	16.75x-112.5x	5.37-0.80
1.5X	10.05x-67.5x	22.89-3.41	15.07x-101.25	10.62-1.58	20.1x-135x	5.97-0.89	25.12x-168.75x	3.58-0.53
2X	13.4x-90x	17.16-2.56	20.1x-135x	7.96-1.19	26.8x-180x	4.48-0.67	33.5x-225x	2.69-0.40

10.7 Use of contrast disc

User can use the black/white disc in order to enhance the contrast of the image when using the incident illumination.


1. Remove the opaline transparent disc used for transmitted light.
2. Insert black/white disc with the desired side facing up. (Fig. 12)

If user wants to look at bright opaque objects, place the disc with the black side up.

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10.9 Use of eyeshields

• Use with eyeglasses

Fold rubber eyeshields with both hands. Folded eyeshields avoid scratching the lenses of eyeglasses. (Fig. 14a)

• Use without eyeglasses

Raise eyeshields and observe at the microscope placing eyes to the shields, avoiding external light to disturb the observation. (Fig. 14b)

10.10 Microphotography

Trinocular head can be connected to cameras via an adapter.

Please refer to the adapter and camera manuals for further details.

11. Troubleshooting

Review the information in the table below to troubleshoot operating problems.

PROBLEM	CAUSE	SOLUTION
I. Optical Section:		
LED operates, but field of view remains dark	Power supply is unplugged.	Connect .
	Brightness is too low	Set brightness to a proper level
The edge of the field of view is vignetted or the brightness is asymmetric.	The incident illuminator is not correctly oriented	Change the angle of the incident illuminator.
Dirt or dust is visible in the field of view	Dirt/dust on the specimen	Clean the specimen
	Dirt/dust on the eyepieces	Clean the eyepieces
Visibility is poor. • Image is not good. • Contrast is poor. • Details are indistinct. • Image glares	The lenses (eyepieces, additional lens) are dirt	Thoroughly clean all the optical system
One side of the image is out of focus	The specimen is out of place (tilted)	Place the specimen flat on the stage.
II. Mechanical Section:		
Coarse focus knob is hard to turn	Tension adjustment ring is too tight	Loosen tension adjustment ring
Focus is unstable	Tension adjustment ring is too loose	Tighten tension adjustment ring
III. Electrical Section:		
LED doesn't turn on.	Power supply not connected	Check for proper connection
Brightness is not enough	Brightness setting is too low	Adjust brightness
Light blinks	Power supply not well connected	Check for proper connection
IV. Observation Tube:		
Field of view of one eye does not match that of the other.	Interpupillary distance is incorrect.	Adjust interpupillary distance.
	Incorrect diopter adjustment.	Adjust diopter.
	Your view is not accustomed to microscope observation.	Upon looking into eyepieces, try looking at overall field before concentrating on specimen range. You may also find it helpful to look up and into distance for a moment before looking back into microscope.
V. Microphotography:		
Bright spots appear on the image	Stray light entering in the microscope through eyepieces or camera viewfinder	Cover eyepieces and viewfinder with a dark cloth

12. Repair and maintenance

Microscopy environment

This microscope is recommended to be used in a clean, dry and shock free environment with a temperature of 0-40°C and a maximum relative humidity of 85 % (non condensing). Use a dehumidifier if needed.

To think about when and after using the microscope

The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out. Never mishandle or impose unnecessary force on the microscope.

Never attempt to service the microscope yourself.

After use, turn off the light immediately, cover the microscope with the included dust-cover, and keep it in a dry and clean place.

Electrical safety precautions

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.

Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

Cleaning the optics

If the optical parts need to be cleaned try first to: use compressed air. If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent. And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether. Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.

Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics. Do not disassemble objectives or eyepieces in attempt to clean them.

For the best results, use the VWR cleaning kit (see catalog number below).

If you need to send the microscope to manufacturer for maintenance, please use the original packaging.

13. Technical service

Web Resources

Visit the VWR's website at **www.vwr.com** for:

- Complete technical service contact information
- Access to VWR's Online Catalogue, and information about accessories and related products
- Additional product information and special offers

Contact us For information or technical assistance contact your local VWR representative or visit. **www.vwr.com**.

14. Warranty

VWR International warrants that this product will be free from defects in material and workmanship for a period of five (5) years from date of delivery. If a defect is present, VWR will, at its option and cost, repair, replace, or refund the purchase price of this product to the customer, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear. If the required maintenance and inspection services are not performed according to the manuals and any local regulations, such warranty turns invalid, except to the extent, the defect of the product is not due to such non-performance.

Items being returned must be insured by the customer against possible damage or loss. This warranty shall be limited to the aforementioned remedies. IT IS EXPRESSLY AGREED THAT THIS WARRANTY WILL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND IN LIEU OF THE WARRANTY OF MERCHANTABILITY.

15. Compliance with local laws and regulations

The customer is responsible for applying for and obtaining the necessary regulatory approvals or other authorizations necessary to run or use the Product in its local environment. VWR will not be held liable for any related omission or for not obtaining the required approval or authorization, unless any refusal is due to a defect of the product.

Disposal

This equipment is marked with the crossed out wheeled bin symbol to indicate that this equipment must not be disposed of with unsorted waste.

Instead it is your responsibility to correctly dispose of your equipment the end of its life cycle by handling it over to an authorized facility for separate collection and recycling. It is also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect from health hazards the persons involved in the disposal and recycling of the equipment.

For more information about where you can drop off your waste equipment, please contact your local dealer from whom you originally purchased this equipment.

By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

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Ireland / Northern Ireland

VWR International Ltd /
VWR International (Northern Ireland) Ltd
Orion Business Campus
Northwest Business Park
Ballycoolin
Dublin 15
Tel.: 01 88 22 222
Fax: 01 88 22 333
E-mail: sales@ie.vwr.com

Italy

VWR International PBI S.r.l.
Via San Giusto 85
20153 Milano (MI)
Tel.: 02-3320311/02-487791
Fax: 800 152999/02-40090010
E-mail: info@it.vwr.com

The Netherlands

VWR International B.V.
Postbus 8198
1005 AD Amsterdam
Tel.: 020 4808 400
Fax: 020 4808 480
E-mail: info@nl.vwr.com

New Zealand

Global Science - A VWR Company
241 Bush Road
Albany 0632, Auckland
Tel.: 0800 734 100
Fax: 0800 999 002
E-mail: sales@globalscience.co.nz

Norway

VWR International AS
Haavard Martinsens vei 30
0978 Oslo
Tel.: 02290
Fax: 815 00 940
E-mail: info@no.vwr.com

Poland

VWR International Sp. z o.o.
Limbowa 5
80-175 Gdansk
Tel.: 058 32 38 200 do 204
Fax: 058 32 38 205
E-mail: info@pl.vwr.com

Portugal

VWR International -
Material de Laboratório, Lda
Edifício Neopark
Av. Tomás Ribeiro, 43- 3 D
2790-221 Carnaxide
Tel.: 21 3600 770
Fax: 21 3600 798/9
E-mail: info@pt.vwr.com

Singapore

VWR Singapore Pte Ltd
18 Gul Drive
Singapore 629468
Tel.: +65 6505 0760
Fax: +65 6264 3780
E-mail: sales@sg.vwr.com

Spain

VWR International Eurolab S.L.
C/ Tecnología 5-17
A-7 Llinars Park
08450 - Llinars del Vallès
Barcelona
Tel.: 902 222 897
Fax: 902 430 657
E-mail: info@es.vwr.com

Sweden

VWR International AB
Fagerstagatan 18a
163 94 Stockholm
Tel.: 08 621 34 00
Fax: 08 621 34 66
E-mail: kundservice@se.vwr.com

Switzerland

VWR International GmbH
Lerzenstrasse 16/18
8953 Dietikon
Tel.: 044 745 13 13
Fax: 044 745 13 10
E-mail: info@ch.vwr.com

Turkey

Pro Lab Laboratuvar Teknolojileri Ltd.Şti.
a VWR International Company
Orta Mah. Cemal Gürsel Caddesi
Ördekcioglu İşmerkezi No.32/1
34896 Pendik - İstanbul
Tel.: +90216 598 2900
Fax: +90216 598 2907
Email: info@pro-lab.com.tr

UK

VWR International Ltd
Customer Service Centre
Hunter Boulevard - Magna Park
Lutterworth
Leicestershire
LE17 4XN
Tel.: 0800 22 33 44
Fax: 01455 55 85 86
E-mail: uksales@uk.vwr.com

**GO TO VWR.COM FOR THE
LATEST NEWS, SPECIAL OFFERS
AND DETAILS OF YOUR LOCAL
VWR DISTRIBUTOR**

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