VWR Thermal Shake *lite*

INSTRUCTION MANUAL

EUROPEAN CATALOGUE NUMBER:

460-0249 Thermal Shake *lite*, 230V

VERSION: 1

ISSUED: 01-2015

Legal Address of Manufacturer

Europe

VWR International bvba
Researchpark Haasrode 2020
Geldenaaksebaan 464
B-3001 Leuven
+ 32 16 385011
<http://be.vwr.com>

Country of Origin

United States

TABLE OF CONTENTS

SECTION 1.0	INTRODUCTION
SECTION 2.0	SYMBOLS AND CONVENTIONS
SECTION 3.0	SAFETY PRECAUTIONS
SECTION 4.0	SPECIFICATIONS
SECTION 5.0	EQUIPMENT OVERVIEW
SECTION 6.0	CONTROLS OVERVIEW
SECTION 7.0	OPERATION GUIDE
SECTION 8.0	SERVICE AND TROUBLESHOOTING
SECTION 9.0	PARTS LIST
SECTION 10.0	WARRANTY
SECTION 11.0	EQUIPMENT DISPOSAL

INTRODUCTION

The VWR Thermal Shake *lite* has microprocessor controlled cooling & heating capabilities. The incubator can be equipped with various block sizes to accommodate several tube sizes and microplates. It can perform a wide range of applications including, sample storage, storage and reaction of various kinds of enzyme, DNA amplification, pre-denaturation of electrophoresis, and serum solidification. The VWR Thermal Shake *lite* includes the following features:

- ◆ VFD display with high-luminosity;
- ◆ Temperature setting value and practical value displayed simultaneously;
- ◆ Display setting time and practical time simultaneously;
- ◆ Aluminum block protecting the sample from contamination;
- ◆ Aluminum block that is easy to replace, clean, sterilize and suitable for different type of tubes;
- ◆ Internal over-temperature protection improving reliability;
- ◆ Warning function after running time ends;
- ◆ Temperature calibration

SYMBOLS AND CONVENTIONS

Your VWR Thermal Shake *lite* uses internationally accepted graphic symbols to help convey information to the user and to call the users attention to important safety precautions and guides for using this equipment.

Indicates that user should consult manual further description or discussion.

Failure to heed this information may present a risk of damage or injury to persons or equipment.

Indicates "Heating"

WARNING: To avoid accidental burning be very careful touching the heater block of the unit. It can be very hot after it is used at high temperatures. Allow the heater block to cool down before handling

SAFETY PRECAUTIONS

During the operation, maintenance, or repair of the VWR Thermal Shake *lite*, the following safety measures should be taken. Otherwise, the safeguards provided by the VWR Thermal Shake *lite* are likely to be damaged, the rated safety level to be reduced, and the rated operation conditions to be affected.

VWR International shall not be in any way responsible for the consequences resulting from operator's not observing the following requirements.

This instrument is intended for indoor use only.

Grounding

A.C. power's grounding should be reliable to safeguard against an electric shock. The 3-pin plug supplied with the VWR Thermal Shake *lite* 's power cable is a safety device that should be matched with a suitable grounded socket. Never allow the third ground pin to be floating. If the 3-pin plug cannot be inserted, it is recommended to ask an electrician to install an appropriate power socket

Keep Away from Electric Circuits

The operator should not open the VWR Thermal Shake *lite* without first consulting VWR's Service Department. Changing components or adjusting certain parameters inside the device must be performed by certificated professional maintenance personnel only. Do not change elements while the power is still on.

A.C. Power Considerations

Before turning on the power, always make sure that the mains voltage is within the range of required power supply ($\pm 10\%$ difference is allowed) and the rating current of the power socket meets the required specification as shown in below contents.

A.C. Power Line Considerations

As an accessory of the VWR Thermal Shake *lite*, the A.C. power supply should be the default device. Once it is damaged, the A.C. power line may not be repaired, but must

be replaced by a new one. The power supply should be free of heavy objects during the VWR Thermal Shake *lite*'s operation. Keep the power supply away from high traffic areas.

Connect the A.C. Power Line

While connecting or disconnecting the power line, user should insert the plug firmly to ensure good contact between the plug and socket. Pull the plug, but not the cable, when the plug needs to be disconnected from the mains.

Design Environments

The VWR Thermal Shake *lite* should be placed in a low-humidity, dust-free, well-ventilated room without caustic gas or powerful magnetic interference. The VWR Thermal Shake *lite* should not be operated in close proximity to water sources, such as pools and water pipes. Never cover or obstruct the openings of the VWR Thermal Shake *lite*, which are designed for ventilation and to prevent the device's interior from becoming too hot. When a single device is running, the shortest distance between its openings and the nearest object should be at least 50cm (20 inches); when two devices or more are running simultaneously, the shortest distance is 100 cm (40 inches) among these machines. Do not place the device on a soft surface. Doing so will result in adverse ventilation near the device's bottom openings.

Operating conditions that cause a high temperature environment will result in degraded performance or failure of the The VWR Thermal Shake *lite*. Additionally, the device should be protected against any kind of heat sources such as sunlight, ovens, or central heating equipment. If the VWR Thermal Shake *lite* will sit idle for extended periods of time, it is recommended to disconnect the power line from the mains and cover the device with a piece of soft cloth or plastic to prevent against dust.

Once one of the following events occurs, the Operator is advised to disconnect the power cable from the mains, and contact the distributor or ask a certified maintenance engineer for assistance.

- Liquid enters into the device;
- The device malfunctions, giving off an abnormal sound or odor;
- The device falls onto the floor or the housing is damaged
- Significant changes in the device's performance

SPECIFICATIONS

NORMAL OPERATING CONDITIONS

Ambient temperature:	10°C ~ 30°C
Relative humidity:	≤70%
Power supply:	AC220-240V 50/60Hz 150W
Placement table-board:	stable and horizontal

PERFORMANCE

Temperature Control Range	0° ~ 105°C
Temperature Range	Ambient ~ 100°C
Timing Range	1 min ~ 99h59min
Cooling Time	From ambient to ambient -10° ≤ 8min From 100° to ambient + 10° ≤ 15min
Heating Time	≤ 12min (from 20°C to 100°C)
Temperature Accuracy	≤ ± 0.5 °C
Temperature Uncertainty	≤ ± 0.5 °C
Heating Rate	about 6°C /min (from 20°C to 100°C)
Block Temperature Uniformity	≤ ± 0.5 °C
Mixing Rate	300~1500rpm
Amplitude	3mm

Before power-on, please check whether the above operating conditions are satisfied. Pay special attention to the power cable's grounding and make sure whether unit is placed on stable and horizontal surface.

TRANSPORTATION AND STORAGE CONDITIONS

Ambient temperature : -20°C ~+55°C

Relative humidity: ≤80%


BASIC PARAMETERS

Model Parameters	
Blocks	460-0250, 40 x 1.5mL 460-0251, 54 x 0.5mL 460-0252, 96 x 0.2mL 460-0253, 24 x 15mm 460-0254, Water bath (ID: 115X73X38mm) (not suggested to use when mixing function is on) 460-0255, 26 x 0.5mL and 24 x 1.5mL 460-0256, 40 x 2.0mL 460-0257, 96 Well ELISA Plate
Dimension Inches/mm	12.9" x 6.5" x 9.8" / 328×166×249mm(L×W×H)
Weight Lbs. (kg)	18.75 lbs (8.5Kg)

EQUIPMENT OVERVIEW

This chapter describes the VWR Thermal Shake *lite*'s mechanical structure, keyboard and key functions, and overview before powering-on. It should be read carefully before the VWR Thermal Shake *lite* is first operated.

1 Structure Description

① Block

② Case(Main Body)


③ Display Panel and
Keyboard

④ Sticker


⑤ Power Connector, Fuse and Power Switch

CONTROLS OVERVIEW

Diagram of display panel

Button explanation

Temperature set up button. Press “▲” or “▼” to adjust the value and set the requested temperature. Keep pressing “▲” or “▼” may help to speed up the setting time.

Vibration time set up button. Press “▲” or “▼” to adjust the value and set the requested vibration time. Keep pressing “▲” or “▼” may help to speed up the setting time.

Vibration speed set up button. Press “▲” or “▼” to adjust the value and set the requested vibration speed.

Instant spot vibration button. Press it to start spot vibration. The running speed is the highest vibration speed. It is also a function

Start button or stop button.. Press this button to start or stop the mixing of the device.(valid for only mixing)

OPERATION GUIDE

This chapter gives a through instruction to the temperature set up, vibration time set up, vibration speed set up of this instrument. It also explains how to correct the discrepancy between the actual temperature and displayed temperature, and how to set up the buzzer.

1 Pre-check before turn on power

Before turning on power, check and confirm the following items first:


1. Unit is on stable and level surface;
2. The voltage should comply with the spec; (Refer to Chapter 2 for the power supply spec);
3. The power plug should be tightly plugged into the outlet;
4. The ground connection should be reliable;
5. The module connection should be reliable.

If the display is abnormal when the power is turned on, please shut down immediately and contact the supplier.

2 Set up of temperature, vibration time and vibration speed.

a) Approximately 5 seconds after the power is turned on, the monitor will display a temperature of 28.5, which is the instant temperature of the module. The display of “TEMP 100.0” is the set up temperature of last run, the display of “TIME 80:00” is the set up vibration time of last run, and the display of “SPEED 1000” is the vibration speed of last run.



At same time, the device will automatically heat or cool, to achieve the setting temp. point.

a) The digits of the TEMP display may be reduced or increased by pressing “Temp” or “Temp” respectively.

The digits of the TIME display may be reduced or increased by pressing “Time” or “Time” respectively.

Press the above mentioned buttons for more than 3 seconds to set up the exact value.

If the temperature should be set at 40° and the time should be set up at 60:00, then press the “Temp” and not release to see the number reducing; When the number shows 40, release the button, and it will be automatically saved in 3 seconds. Next press the “Time” and not release to see the number reducing; when the number shows 60:00, release the button, and it will be automatically saved in 3 seconds.

c) Press “Speed” or “Speed” to set up for the speed. The number will increase or reduce at a rate of 50rpm, release the button when the number reaches exact value, and it will be automatically saved in 3 seconds.

When set up is done, press “Start/Stop” to start mixing.

After the setting value of temperature is saved successfully, the device will automatically heat or cool, to achieve the setting temp. point.

3 PULSE vibration

Press “Short Prog” to start instant vibration or short term vibration. Press the button to start it and release the button to stop it. The running speed of pulse vibration is the highest vibration speed.

4 Correct temperature discrepancy

The temperature of the instrument has been calibrated before shipment. However, due to various conditions, there might be discrepancy between the actual temperature and displayed temperature. On such case, you may use the calibration button to correct the discrepancy.

1. To ensure the accuracy of the temperature, please wait for at least 30 minutes after the instrument is set at a constant temperature, and then perform the calibration.
2. Please use certified standard class two mercury thermometer to calibrate this instrument.
3. Calibration point: Center hole of the module. Please fill the hole with paraffin oil and immerge the thermometer bulb in it.


The method of the temperature calibration of this instrument is two-point linear calibration. The calibration point may be set freely, and by setting the second calibration point the same as the first calibration point, you can adjust the calibration point to one temperature point.

For two-point calibration, if set the two points at 40°C and 100°C, then other temperature points will be automatically adjusted per the linear relationship of the two calibration points.

Please read the following operation instruction:


a) Press both “Temp” and “Temp” at the same time, and the display will show the set up value of last time (for example 20°C) and the cursor of the first calibration point; C-1;

b) Press the temperature set up button “Temp” or “Temp” to set up the first calibration temperature such as 40°C.



c) Perform the same operation again to set up the second calibration temperature such as 100°C;

Note: when setting up the temperature points, an order of low temperature point first and high temperature after is recommended. If on the case high temperature point is set first and low temperature is set after, then when performing the following actual calibration, please follow the order of low temperature first and high temperature after.

d) Press both “Temp” and “Temp” at the same time to start the calibration. The display will show 40.0 and the temperature of the instrument will automatically increase to 40°C. Once the temperature is kept constant, timing will start. 30 minutes later, press “Short Prog” for one time and the display will show the blink of “TEMP.” At this time, the display will show the actual temperature reading of the thermometer. For example, if the thermometer reads 40.5°C, you can adjust the temperature display to 40.5°C by pressing the temperature set up button. Next press “Short Prog” for one time and the display will show 100.0°C. Following the same instruction, when the temperature of the instrument increases to 100°C, wait for 30 minutes,

input the calibrated value and press



“Short/Prog” to save it.

When the calibration of 40°C is done, press

“Start/Stop” to exit; at this time the calibration of 40°C is ineffective.

When the calibration is done for once at both 40°C and 100°C, check for the discrepancy between the two calibration points and the actual temperature points. It should be within 0.5°C. If it's more than 0.5°C, please follow the above instruction and recalibrate until the request is met.

5 Buzzer set up

A Upon failure's occurrence or ending, the instrument will buzz “DI, DI...” You may disable it if you do not need this function. However, the default setting is to have it on.

B This instrument has keyboard beep; once the button is pressed, it will beep “DI.” You may disable it if you do not need this function. However, the default setting is to have it on.

Please read the following operation instructions:



a) Press both “” and “” at the same time, And the display will show On, Beep;

b) Press “” for once and change On to OFF;

c) Press “” and “” at the same time to save the setting, and buzz will be disabled upon failure’s occurrence or ending. The display will then show the keyboard beep status. If the display will show On K_S. Press “” once and change On to OFF;

d) Then press both “” and “” at the same time to save the setting, and keyboard beep will be disabled.

If you want to exit during set up, you may press “”, and the setting will be ineffective.

TROUBLESHOOTING AND REPAIR

This chapter describes possible faults of this instrument along with instructions to the cause analysis and teaches you how to resolve the fault.

Troubleshooting and Repair

No.	Problem	Probable Cause	Remedy
1	No display shown in the monitor when the power is on.	Power is not connected.	Check the power and connect it.
		Fuse has burned out.	Change the fuse Item no. 460-0266
		Switch is broken.	Change the switch.
		Others.	Contact VWR
2	During run, the housing shakes abnormally and makes abnormal noise.	Use improper block	Change the block. Contact VWR
		Unit is placed on uneven surface	Place unit on even, stable surface
		Round feet badly worn	Change round feet.
		Others	Contact VWR
3	Bad temperature uniformity of block.	Improper installation of block.	Re-install the block.
		Bad cooling components.	Contact VWR
4	The system no longer operates after changing block.	Bad grounding of machine.	Check power cable to insure reliable grounding; Do remember to turn off the power when changing block.

No.	Problem	Probable Cause	Remedy
5	Serious discrepancy between the actual temperature and displayed temperature.	Broken sensor or bad contact.	Contact VWR
6	The cooling of the module drastically slows down or the temperature cannot reach below room temperature.	Thermoelectric cooling module is broken.	Contact VWR
		Fan is broken or does not work.	
7	The module is neither able to heat nor cool.	Temperature sensor is broken.	Contact VWR
		Thermoelectric cooling module is broken.	
8	Vibration stops.	Drive is broken.	Contact VWR
		Motor is broken.	Contact VWR
		Fuse has burned out.	Change the fuse. Item no. 460-0266

User shall not open the cover of this instrument when under warranty; If housing does need to be opened to perform certain repair per the above instruction, please contact VWR.

Technical service

Web Resources

Visit the VWR's website at www.vwr.com for:

- Complete technical service contact information
- Access to VWR's Online Catalogue, and information about accessories and related products
- Additional product information and special offers

Contact us For information or technical assistance contact your local VWR representative or visit. www.vwr.com.

PARTS LIST

User replaceable accessories and spare parts

Description	Cat. No.
L socket head wrench	460-0264
Cushioning balls for foot	460-0265
Fuse- used for 230V model (250V1.5A Φ 5*20mm)	460-0266

WARRANTY

VWR International warrants that this product will be free from defects in material and workmanship for a period of two (2) years from date of delivery. If a defect is present, VWR will, at its option and cost, repair, replace, or refund the purchase price of this product to the customer, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear. If the required maintenance and inspection services are not performed according to the manuals and any local regulations, such warranty turns invalid, except to the extent, the defect of the product is not due to such non-performance.

Items being returned must be insured by the customer against possible damage or loss. This warranty shall be limited to the aforementioned remedies. IT IS EXPRESSLY AGREED THAT THIS WARRANTY WILL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND IN LIEU OF THE WARRANTY OF MERCHANTABILITY.

Compliance with local laws and regulations

The customer is responsible for applying for and obtaining the necessary regulatory approvals or other authorisations necessary to run or use the Product in its local environment. VWR will not be held liable for any related omission or for not obtaining the required approval or authorisation, unless any refusal is due to a defect of the product.

EQUIPMENT DISPOSAL

This equipment is marked with the crossed out wheeled bin symbol to indicate that this equipment must not be disposed of with unsorted waste.

Instead it's your responsibility to correctly dispose of your equipment at lifecycle -end by handling it over to an authorized facility for separate collection and recycling. It's also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect from health hazards the persons involved in the disposal and recycling of the equipment.

For more information about where you can drop off your waste of equipment, please contact your local dealer from whom you originally purchased this equipment.

By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

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Ireland / Northern Ireland

VWR International Ltd /
VWR International (Northern Ireland) Ltd
Orion Business Campus
Northwest Business Park
Ballycoolin
Dublin 15
Tel.: 01 88 22 222
Fax: 01 88 22 333
Email: sales@ie.vwr.com

Italy

VWR International PBI S.r.l.
Via San Giusto 85
20153 Milano (MI)
Tel.: 02-3320311/02-487791
Fax: 800 152999/02-40090010
Email: info@it.vwr.com

The Netherlands

VWR International B.V.
Postbus 8198
1005 AD Amsterdam
Tel.: 020 4808 400
Fax: 020 4808 480
Email: info@nl.vwr.com

Norway

VWR International AS
Haavard Martinsens vei 30
0978 Oslo
Tel.: 02290
Fax: 815 00 940
Email: info@no.vwr.com

Poland

VWR International Sp. z o.o.
Limbowa 5
80-175 Gdansk
Tel.: 058 32 38 200 do 204
Fax: 058 32 38 205
Email: info@pl.vwr.com

Portugal

VWR International -
Material de Laboratório, Lda
Edifício Neopark
Av. Tomás Ribeiro, 43- 3 D
2790-221 Carnaxide
Tel.: 21 3600 770
Fax: 21 3600 798/9
Email: info@pt.vwr.com

Spain

VWR International Eurolab S.L.
C/ Tecnología 5-17
A-7 Llinars Park
08450 - Llinars del Vallès
Barcelona
Tel.: 902 222 897
Fax: 902 430 657
Email: info@es.vwr.com

Sweden

VWR International AB
Fagerstagatan 18a
163 94 Stockholm
Tel.: 08 621 34 00
Fax: 08 621 34 66
Email: kundservice@se.vwr.com

Switzerland

VWR International GmbH
Lerzenstrasse 16/18
8953 Dietikon
Tel.: 044 745 13 13
Fax: 044 745 13 10
Email: info@ch.vwr.com

Turkey

VWR International Laboratuvar Teknolojileri Ltd.Şti.
Orta Mah. Cemal Gürsel Caddesi
Ördekcioglu İşmerkezi No.32/1
34896 Pendik - İstanbul
Tel.: +90216 598 2900
Fax: +90216 598 2907
Email: info@pro-lab.com.tr

UK

VWR International Ltd
Customer Service Centre
Hunter Boulevard - Magna Park
Lutterworth
Leicestershire
LE17 4XN
Tel.: 0800 22 33 44
Fax: 01455 55 85 86
Email: uksales@uk.vwr.com

Australia

VWR International, Pty Ltd.
Level 1, Unit 1a/60 Enterprise Place
Tingalpa, Queensland, 4173
Tel.: 1300 727 696
Fax: 1300 135 123

China

VWR International China Co., Ltd
Rm.219, 2100 Dongming Road
Pudong New District
Shanghai 200123
Tel.: +86-21-5898 6888
Fax: +86-21-5855 8801
Email: info_china@vwr.com

India

VWR Lab Products Private Limited
135/12, Brigade Towers, 2nd Floor
Front wing, Brigade Road,
Bengaluru, India – 560 025
Tel.: +91-80-41117125/26 (Bengaluru)
Tel.: +91-2522-647911/922 (Mumbai)
Fax: +91-80-41117120
Email: vwr_india@vwr.com

New Zealand

Global Science - A VWR Company
241 Bush Road
Albany 0632, Auckland
Tel.: 0800 734 100
Fax: 0800 999 002
Email: sales@glbalscience.co.nz

Singapore

VWR Singapore Pte Ltd
18 Gul Drive
Singapore 629468
Tel: +65 6505 0760
Fax: +65 6264 3780
Email: sales@sg.vwr.com

**GO TO VWR.COM FOR THE LATEST
NEWS, SPECIAL OFFERS AND
DETAILS OF YOUR LOCAL VWR
DISTRIBUTOR**

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