VWR®

Portable Turbiditymeter TIR 210

INSTRUCTION MANUAL

**European Catalogue Number:
705-1166**

Version: 1
Issued: 24. December 2020

Legal Address of Manufacturer

Europe

VWR International bvba

Researchpark Haasrode 2020

Geldenaaksebaan 464

B-3001 Leuven

+ 32 16 385011

<http://be.vw3r.com>

Country of origin

Made in Germany

Table of Contents

Warning	5
Safety Information	5
Delivery contents	6
Unpacking	6
Set up	6
Replacement of batteries	7
Intended use	7
Symbols and conventions	8
Product Specifications	8
Overview	9
Description of Buttons	9
Operating principles	9
Operational keys	10
Display backlight	10
Recall of stored data	10
Menu selections	10
diS – Recall of stored data	10
Transmission of stored data to a PC	11
Setting date and time (24-hour-format)	11
Calibration	12
User calibration	12
When to calibrate	13
Factory calibration reset	13
Guidelines for turbidity measurements	14
Cleaning of vials & sampling containers	14
Correct positioning of the vial (Ø 24 mm)	15
Correct filling of the vial	15
Abbreviations	15
TSK standards	15
Storing TSK standards	16
Using < 0.1 NTU standards	16
Preparing and using standards – infrequent use	16
Preparing and using standards – frequent use	16
Formazin Standards	16
Preparation of dilution water	16
Preparation of 4000 NTU Formazin Stock Solution	17
Preparation of dilutions from the 4000 NTU Formazin Stock Solution	17
Indexing and matching sample vials	18
Indexing a single sample vial	18
Indexing a set of sample vials	19

Measurement techniques	20
Degassing – removal of bubbles	20
Measurement of high turbidity values.	21
Measurement of low turbidity values	21
Troubleshooting	22
Operating messages.	22
Error codes.	22
Repair and maintenance	22
User replaceable accessories and spare parts	23
Technical service	23
Warranty.	24
Compliance with local laws and regulations	24
Equipment disposal.	24

Warning

CAUTION

The accuracy of the instrument is only valid if the instrument is used in an environment with controlled electromagnetic disturbances according to DIN 61326. Wireless devices, e.g. wireless phones, must not be used near the instrument.

Safety Information

Read this entire manual carefully before use!

Carefully read and follow the SAFETY INSTRUCTIONS at the beginning of this manual!

The system may only be used by qualified persons!

Keep this manual for later reference!

Also observe important notices in order to avoid malfunctions and faults.

If possible, save the complete transport packaging for later transport.

Follow the safety instructions listed here for your own safety!

The safety instructions draw your attention to potential dangers. They also contain information for how you can prevent danger with appropriate conduct.

Use is only permitted for qualified persons.

 WARNING!	The device may not be operated in explosion-prone rooms. Health hazard! Observe the required protective measures when handling the sample! Health hazard!! Do not handle depleted batteries or leaking material with bare hands! Wear protective gloves! Avoid contact with the eyes and skin!
 WARNING!	Breakage of glass. Danger of cutting injuries. Despite protection from the transport packaging, breakage of glass can occur during transport. The resulting broken shards may have sharp edges and could cause cutting injuries when handled carelessly. Look out for broken glass when unpacking, and wear protective gloves. When handling the system, always ensure a firm grip in order to avoid breakage of glass! Do not strike the test bottles against edged during transport of the system!
 WARNING!	Turbidity Standards and TSK Standards are formulated for chemical analysis and must not be used for any other purpose. Reagents must not get into the hands of children! Some of the reagents contain substances which are not entirely harmless environmentally. Be aware of the ingredients and take proper care when disposing of the test solution.

WARNING!

Ensure proper disposal of reagent solutions.
Material Safety Data Sheets are available on request

Delivery contents

Standard contents for TIR 210:

- ☐ 1 Turbiditymeter in plastic case
- ☐ 1 Cap for Turbiditymeter
- ☐ 1 Instruction manual
- ☐ 1 Test Certificate
- ☐ 1 9V battery
- ☐ 1 battery compartment cover
- ☐ 2 screws
- ☐ 1 screwdriver
- ☐ 1 set of 2 sample cuvettes with black lid, height 55 mm
- ☐ 1 cleaning cloth for cuvettes
- ☐ 1 set of Turbidity Standards TSK
- ☐ 1 USB cable

Different Refill Packs available on request.


Unpacking

Carefully inspect all items to ensure that every part of the list below is present and no visible damage has occurred during shipment. If there is any damage or something is missing, please contact your local distributor immediately.

Set up

Before working with the photometer insert the batteries (delivery contents).

Replacement of batteries:

CAUTION:

To ensure that the instrument is water proof:

- **seal ring (C) must be in position**
- **battery compartment cover (D) must be fixed with the two screws**

If the batteries are removed for more than one minute the date and time menu starts automatically when the photometer is switched on the next time.

Intended use

The TIR 210 is a portable turbidity meter, based on the requirements of ISO 7027 (Water Quality – Determination of Turbidity). The instrument features auto ranging over the range of 0.01 to 1100 NTU/FNU.

The turbidimeter is supplied in a case complete with accessories and spares. Calibration Standards guarantee stability and reproducibility of the results. Always replace the cover on the sample chamber to protect against dust.

Always adhere to the limits specified in chapter "Technical Data". To ensure correct and intended usage, the instructions, especially the "Important Safety Information", have to be read and understood.

The instrument is not used as intended if

- it is used in application areas not stated in these instructions.
- it is used in conditions which deviate from the ones described in these instructions.

The safety instructions in this operating manual must be observed.

The instrument may only be used under the conditions and for the purposes for which it was designed.

The instrument must be handled with care (do not throw, drop, etc.) and used in accordance with the technical data. It must be protected against soiling.

Symbols and conventions

	CAUTION! This symbol indicates a potential risk and alerts you to proceed with caution
---	---

Product Specifications

Light source:	LED, Filter ($\lambda_1 = 860 \text{ nm}$)
Principle:	Nephelometric (Non Ratio)
Range:	0.01 – 1100 NTU ¹⁾
Resolution:	0.01 – 9.99 NTU = 0.01 NTU 10.0 – 99.9 NTU = 0.1 NTU 100 – 1100 NTU = 1 NTU
Accuracy:	$\pm 2.5 \%$ of reading, or $\pm 0.01 \text{ NTU}$, whichever is bigger, in the range of 0.01 - 500 NTU ; $\pm 5 \%$ of reading from 500 to 1100 NTU,
Repeatability:	$\pm 1 \%$ of reading or $\pm 0.01 \text{ NTU}$ which ever is greater
Interface:	Micro USB
Data transmission:	CSV format
Data storage:	125 data records
Battery:	9 V-block battery (Life approx. 550 tests without PC connection)
Auto-OFF:	Automatic switch off 10 minutes after last keypress
Ambient conditions:	5–40°C 30–90% rel. humidity (non-condensing)
CE:	Certificate for Declaration of CE-Conformity at www.vwr.com

To ensure maximum accuracy of test results, always use the reagent systems supplied by the instrument manufacturer.

Subject to technical modification!

¹⁾FNU is applicable to "Non Ratio" instruments.

Overview

Reagent	Quantity	cat. no.
Turbidity Standards TSK: Standard < 0.1 NTU, 20 NTU, 200 NTU, 800 NTU	1 Set	84209.600

Description of Buttons

Switching the Turbiditymeter on or off


Confirm a selection with the [MODE] key.

Perform a measurement.

Press the [!] key to turn the display backlight on or off


Operating principles


The instrument measures turbidity in the range 0.01 to 1100 NTU/FTU featuring auto ranging. The light source is an infrared LED (light emitting diode) with a wavelength of 860 nm.

The emitted light is reflected by turbidity in the sample. The scattered light will be detected at an angle of 90° by a photodiode.

This principle is part of ISO 7027.

The international Reference Standard for turbidity is a Formazin solution. Results related to these standards are indicated as FNU (Formazine Nephelometric Units).

Operational keys

Performing turbidity measurement

Switch the unit on using the [ON/OFF] key.

The display shows the following:

Fill a clean vial with the water sample up to the mark, screw the cap on and place the vial in the sample chamber making sure that the X marks are aligned.

Press the [READ] key.

The "Method" symbol flashes for approx. 8 seconds.

The result is shown in **NTU**.

Repeating the test:

Press the [READ] key again.

Display backlight

Press the [!] key to turn the display backlight on or off. The backlight is switched off automatically during the measurement.

Recall of stored data

If the instrument is switched on, press the [!] key for more than 4 seconds to access the recall menu.

Menu selections

Press the [MODE] key and **hold**.

Switch the unit on using the [ON/OFF] key.

Allow the 3 decimal points to be displayed before releasing the [MODE] key.

The [!] key allows for selection of the following menu points:

- ▲ a) and b);
 - a) recall stored data ([!] press once) and
 - b) transfer to a PC ([!] press twice)

- ▲ ▼ setting the date and time

- ▼ user calibration



The selected menu is indicated by an arrow in the display.

Confirm the selection with the [MODE] key.

▲ a) diS – Recall of stored data

The photometer shows 125 data sets in the following format (automatically proceeds every 3 seconds until result is displayed):

Number n xx (xx: 16...1)
 Year YYYY (i.e. 2014)
 Date mm.dd (month:month:day)
 Time hh:mm (hour:hour:minute:minute)
 Test Method
 Result x,xx

The [READ] key repeats the current data set.

The [MODE] key scrolls through all stored data sets.

Quit the menu by pressing [!].

1 b) Transmission of stored data to a PC

From the TIR 210 data can be transmitted to a PC. PC data recording software (Available on www.vwr.com → 705-1166 → Related Information) enables receipt of the data stored on the turbidity meter by a PC.

The data can be exported to an Excel® sheet and used directly by the user as an evaluation tool in the accustomed manner. If Excel® is not available, the data can be saved as a .txt file and is then available for later evaluation.

The data is transferred from the turbidity meter in .csv format.

Connect the turbidity meter to the PC with the USB cable included in the scope of delivery.

Start the data recording software on the PC and activate the connection to the turbidity meter.

After menu selection (see above) by pressing the [MODE]-button, start the transfer, the device displays "PrtG" (printing) for about 1 seconds. Then the number of the first data record is displayed and the data is transmitted. All stored data records are transmitted in succession. The device switches to measuring mode when the process is finished.

The transmission process can be cancelled by pressing the [On/Off] button. The device switches off.

Mode

PrtG

On
Off

2 3 Setting date and time (24-hour-format)

After confirming the selection with the [MODE] key the value to be edited will be shown for 2 sec.

The setting starts with the year (YYYY) followed by the actual value to be edited. Same applies for month (mm), day (dd), hour (hh) and minutes (mm). Set the minutes first in steps of 10, press the [!] key to continue setting the minutes in steps of 1.


Increase the value by pressing the [MODE] key.

Decrease the value by pressing [READ] key.

Proceed to the next value to be edited by pressing [!] key.

After setting the minutes and pressing the [!] key the display will show "IS SET" and instrument returns to the measurement mode.

ATTENTION: If the battery is removed for more than one minute the date and time menu starts automatically when the photometer is switched on the next time.

Mode

SET

DATE

YYYY

(2. sec)

Mode

Read

!

Calibration

4 User calibration

0.10

StAn

Read

01:00

0.10

20

StAn

Read

01:00

20

200

StAn

Read

01:00

200

800

StAn

Read

01:00

800

After confirming the selection with the [MODE] key **0.10/StAn** appears alternating in the display.

Place the < 0.1 NTU standard in the sample chamber making sure that the marks are aligned. Place the cover on the instrument.

Press the [READ] key.

The countdown starts.

After the countdown is finished the reading starts automatically.

The display flashes for approx. 8 seconds.

After the measurement **20/StAn** appears alternating in the display.

Invert the 20 NTU standard and place it in the sample chamber, making sure that the marks are aligned. Place the cover on the instrument.

Press the [READ] key.

The countdown starts.

After the countdown is finished the reading starts automatically.

The display flashes for approx. 8 seconds.

After the measurement **200/StAn** appears alternating in the display.

Invert the 200 NTU standard and place it in the sample chamber, making sure that the marks are aligned. Place the cover on the instrument.

Press the [READ] key.

The countdown starts.

After the countdown is finished the reading starts automatically.

The display flashes for approx. 8 seconds.

After the measurement **800/StAn** appears alternating in the display.

Invert the 800 NTU standard and place it in the sample chamber, making sure that the marks are aligned. Place the cover on the instrument.

Press the [READ] key.

The countdown starts.

After the countdown is finished the reading starts automatically.

The display flashes for approx. 8 seconds.

USEr	After the measurement "USEr" alternating with Stor is shown in the display.
Stor	
!	Press [!] key to store the calibration.
Stng	First appears Stng (Storing), followed by StEd (Stored).
StEd	


When to calibrate

The turbidimeter was calibrated at the factory with Formazin Primary Standard and does not require user calibration before use.

As the electronic and optical design of this instrument provides long-term stability the need for frequent calibration is minimised. Perform calibration with TSK every 3 months or if required, more often.

Factory calibration reset

The instrument can be reset to factory calibration. This could be necessary if compromised standards have been used to calibrate the instrument. Then the default settings can be used until new standards have been delivered. During this time the instrument may measure with a higher tolerance.

To reset the calibration press both the [MODE] and [READ] key and **hold**.

Switch the unit on using the [ON/OFF] key.
Release the [MODE] and [READ] keys after approx. 1 second.

The following messages will appear in turn on the display:

SEL	The factory setting is active. (SEL stands for Select)
CAL	

or:

SEL	Calibration has been set by the user. (If the user calibration is to be retained, switch the unit off using the [ON/OFF] key).
cAL	


Calibration is reset to the factory setting by pressing the [MODE] key.
The following messages will appear in turn on the display:

Switch the unit off using the [ON/OFF] key.

Guidelines for turbidity measurements

- a) Vials and caps should be cleaned thoroughly after each test to avoid interferences. Minor residuals can cause errors.
- b) The outside of the vial must be clean and dry, before starting the test. Wipe the vials with a smooth cloth to remove fingerprints or waterdrops.
- c) The vials must be positioned in the sample chamber with the mark on the vial aligned with the mark on the instrument.

- d) Always perform the test with securely capped vials.
- e) Bubbles on the inside of the vial lead to errors.
See chapter "Removing bubbles (Degassing)".
- f) Avoid spillage of water in the sample chamber. If water should leak into the instrument, it can destroy electronic components and cause corrosion.
- g) Contamination of the lens in the sample chamber can result in errors. Check at regular intervals – and if necessary – clean the light entry surfaces of the sample chamber using a moist cloth or cotton buds.
- h) Large temperature differences between the instrument and the environment can lead to errors – e.g. due to the formation of condensation in the area of the lens or on the vial. For best results, perform tests with sample temperatures between 20°C (68°F) and 25°C (77°F).
- i) To avoid errors caused by stray light do not use the instrument in bright sunlight.
- j) Use the instrument in a clean, dust-free environment on a table that is free from vibration / agitation.

Cleaning of vials & sampling containers

Vials, caps and sampling containers should be cleaned thoroughly after each test to avoid interferences. Minor residuals can cause errors.


Residuals:

The vial must be cleaned according to the type of sample measured.


- Replace scratched vials immediately.
- Rinse vials thoroughly with deionised water after each measurement.
- Clean all glassware thoroughly with laboratory detergent and rinse with deionised water.

- Clean heavy contamination by filling the vials with 1:1 HCl followed by multiple rinses with distilled or deionised water.
- Allow vials to air dry.
- Touch vials only at the top to minimise dirt and fingerprints.
- Wipe the vials with the delivered cleaning cloth to remove waterdrops and fingerprints.

Correct positioning of the vial (Ø 24 mm):

Correct filling of the vial:

Abbreviations

Abbreviation	Definition
NTU	Nephelometric Turbidity Unit
FTU	Formazine Turbidity Unit
FNU	Formazine Nephelometric Unit
FAU	Formazine Attenuation Unit

TSK standards

The TSK standards purchased with the instrument should be used to calibrate the instrument. It is possible to perform calibration with fresh dilutions of a Formazin Stock Solution. Using dilutions of a Formazin Stock Solution it is recommended to calibrate and measure with an indexed single vial or a set of matched vials.

84209.600	Set of Turbidity Standards TSK
-----------	--------------------------------

For optimum results when using TSK standards, adhere to the following recommendations:

Storing TSK standards

- TSK standards should always be stored in their original vial.
- Preferably store the vials upright.
- Store standards between 5 and 25°C.
- Avoid prolonged exposure to temperatures exceeding 35°C.
- Store away from direct sunlight.
- Always allow the standards to acclimatise to ambient instrument temperature before use (not to exceed 35°C).
- TSK standards have a shelf life of 12 months if stored correctly.

Using < 0.1 NTU standards

Caution: Never shake or invert this < 0.1 NTU standard.

- After delivery allow the vial to stand for at least 24 hours before use.
- If the standard has been shaken it can take several hours for all the bubbles to dissipate and we suggest the vial is left for at least 24 hours as outlined above.
- If the standard has been inverted accidentally wait at least 15 minutes before using.

Preparing and using standards – infrequent use

Note: These instructions apply to all standards except the < 0.1 NTU standard. If you are using the < 0.1 NTU standard follow the instructions given above. These instructions apply to standards that have been sitting undisturbed for more than one week or to new standards.

1. Shake the standard vigorously for 2–3 minutes.
2. Let the vial stand undisturbed for 5 minutes.
3. Invert the vial between 5–10 times.
4. Immediately place the vial in the sample chamber and wait for 1 minute (countdown).

Preparing and using standards – frequent use

Note: These instructions apply to all standards except the < 0.1 NTU standard. If you are using the < 0.1 NTU standard follow the instructions given above. These instructions apply to frequently used standards (daily or weekly usage).

1. Invert the vial 10 times.
2. Immediately place the vial in the sample chamber and wait for 1 minute (countdown).

Preparation of dilution water

Obtain at least 1000 ml of high quality water (e.g. distilled, demineralised or deionised water). Check the turbidity of the dilution water before use. If the turbidity is greater than 0.5 NTU (FNU) the water should be filtered with membrane filter (0.1 µm). Clean the required glassware with 1:1 hydrochloric (37%) acid and rinse several times with the dilution water.

20252.290	Hydrochloric acid (37%)
-----------	-------------------------

Preparation of 4000 NTU Formazin Stock Solution

In lieu of preparation we recommend using a 4000 NTU Stock Solution that is available from specialised suppliers because these standards are high quality and you do not have to handle the raw materials.

84210.180	FORMAZIN STOCK SOLUTION (4000 NTU), 100ML
84211.230	FORMAZIN STOCK SOLUTION (4000 NTU), 250 ML

Caution: Please observe the handling instructions of the MSDS. Gloves, goggles and breathing equipment must be worn!

To prepare a Formazin Stock Solution from raw materials:

1. Dissolve 0.5 g Hydrazine sulfate ($\text{NH}_2)_2\text{H}_2\text{SO}_4$ in 40 ml dilution water.
2. Dissolve 5.0 g Hexamethylentetramine in 40 ml dilution water.
3. Transfer both solutions quantitativ to a 100 mL volumetric flask and fill to the mark with dilution water.
4. Mix the solution thoroughly.
5. Allow this solution to stand for at least 24 hours at $25 \pm 3^\circ\text{C}$ ($77 \pm 5^\circ\text{F}$) in a dark place (brown glass bottle).
6. The turbidity is developed over this period.

The Formazin Stock Solution has a maximum shelf life of 1 year if stored correctly (in a dark place). The Production of a Formazin Stock Solution is described in "EN ISO 7027" and "Standard Methods for Examination of Water and Wastewater".

24696.186	Hydrazine Sulfate, ACS reagent, 100 g
24560.260	Methenamine $\geq 99.0\%$ AnalR NORMAPUR® analytical reagent, 500 g

Preparation of dilutions from the 4000 NTU Formazin Stock Solution

To prepare the dilutions from a 4000 NTU Formazin Stock Solution and dilution water:

Standard	Step 1	Step 2	Step 3
			
20 NTU	Add 100 ml of dilution water to a clean 200 ml volumetric flask.	Pipette 1.00 ml well mixed 4000 NTU Formazin Stock Solution to this 200 ml flask.	Fill to the mark with dilution water. Close and mix the flask.
200 NTU	Add 50 ml of dilution water to a clean 100 ml volumetric flask.	Pipette 5.00 ml well mixed 4000 NTU Formazin Stock Solution to this 100 ml flask.	Fill to the mark with dilution water. Close and mix the flask.

Standard	Step 1	Step 2	Step 3
800 NTU	Add 50 ml of dilution water to a clean 100 ml volumetric flask.	Pipette 20.00 ml well mixed 4000 NTU Formazin Stock Solution to this 100 ml flask.	Fill to the mark with dilution water. Close and mix the flask.


Use class A volume pipette and glass flasks.

Use dilution water for the < 0.1 NTU standard.

Indexing and matching sample vials

Matched sample vials are required to minimise the effects of optical variation from glass vial to glass vial. Alternatively an indexed single sample vial can be used for every measurement. Once vial orientation in the sample chamber is located, always use this new orientation mark for all measurements independent of the white pre-printed triangle mark.

Indexing a single sample vial

1. Fill a clean and dry vial up to the mark with dilution water (see chapter: Preparing Dilution Water).
2. Cap the vial.
3. Hold the vial by the cap and clean it with the supplied cleaning tissue to remove water spots and finger prints.
4. Switch the turbidimeter on.
5. Place the vial in the sample chamber.
6. Make sure that the marks are aligned.
7. Place the cover on the instrument.
8. Press [Read] key.
9. Note the displayed result.
10. Rotate the vial in the sample chamber about 45°.
11. Place the cover on the instrument.
12. Press [Read] key.
13. Note the displayed result.
14. Repeat this procedure until the smallest NTU value is obtained.
15. Mark the vial at this point. Use this orientation mark for all further measurements.

Note:

Assume the indication "Underrange" is equal to 0.00 NTU.

Indexing a set of sample vials

1. Fill several clean and dry vials up to the mark with dilution water.
2. Cap the vials.
3. Hold the vial by the cap and clean it with the supplied cleaning tissue to remove water spots and finger prints.
4. Switch the turbidimeter on.
5. Place the first vial in the sample chamber making sure that the marks are aligned.
6. Place the cover on the instrument.
7. Press [Read] key.
8. Note the displayed result.
9. Rotate the vial in the sample chamber about 45°.
10. Place the cover on the instrument.
11. Press [Read] key.
12. Note the displayed result.
Repeat this procedure until the smallest NTU value is obtained.
13. Mark the vial.
14. To index more vials, repeat point 1. to 13. for each vial.
15. Repeat this procedure until the reading matches the first vial reading within ± 0.01 NTU.
16. Mark the vial.
17. Repeat this procedure to match other vials.

Note:

It may not be possible to match all vials due to variation in the glass.

Measurement techniques

Degassing – removal of bubbles

Note: Do not use with TSK Standards

If the Turbidity is low it is important to remove air bubbles from the sample, using one or a combination of the following methods:

- Addition of a surfactant
- Application of a partial vacuum
- Application of heat
- Use of an ultrasonic bath

Note:

This procedure can influence the nature of the sample and therefore the turbidity reading.

Type of sample	Method	Description of the method	Notes:
Samples oversaturated with air	Addition of a surfactant	Surfactants minimise the surface tension of a sample, allowing entrained gases to disappear.	Particles in the sample settle more rapidly, so the sample must be swirled before being measured. Avoid vigorous shaking as this causes the surfactant to foam.
Liquid samples without readily volatile components	Use of a partial vacuum	A vacuum can be created with the help of a clean, oil-free syringe or pump fitted onto a vial. The vacuum reduces the atmospheric pressure, so that trapped air bubbles can be removed.	Volatile components can escape from the sample. The vacuum may compound the air bubble problem in viscous samples.
Viscous samples	Use of an ultrasonic bath	The ultrasonic waves excite the sample, effectively removing air bubbles from most samples.	Ultrasonic waves can change the particle size in the sample, therefore changing the turbidity.
Very viscous samples	Heating the sample	Heating the sample makes it less viscous, air bubbles can disappear more easily. The sample has to cool to its original temperature.	Volatile components can disappear from the sample. The attributes of suspended particles change, therefore changing the turbidity.

Measurement of high turbidity values

High turbidity samples with more than 1100 NTU ("overrange") may be diluted. The dilution water should be a water with very low turbidity as described in chapter "Preparing dilution water".

For accurate dilution proceed as follows:

Mix the water sample well and pipette x ml of the water sample (see table below) into a 100 ml volumetric flask. Fill with low turbidity water up to the mark and mix gently.

Water sample (x ml)	Multiplication factor
10	10
25	4
50	2

Fill the diluted water sample into the vial, perform reading and multiply the displayed result with the multiplication factor.

Note:

The dilution of the water sample may alter the characteristics of the suspended particles and produce erroneous results.

Measurement of low turbidity values

Accurate and repeatable measurements of low turbidity values depend on exact measurement techniques.

- Use a clean, unscratched and indexed vial.
- Rinse the vial three times with the sample.
- Add the water sample to the mark.
- Allow the vial to stand for 1-5 minutes so that bubbles can disappear.
- Carefully invert the vial (so that settled particles disperse in the sample).
- Place the vial in the sample chamber and press the Read key.
- Perform multiple measurements, until a reproducible value is displayed (leave the vial in the sample chamber).

Note the smallest consistent and reproducible value.

Troubleshooting

Operating messages

Hi	Measuring range exceeded or excessive turbidity.
Lo	Result below the lowest limit of the measuring range.
	Replace 9 V battery, no further tests possible.
btLo	Battery capacity is too low for the display backlight; measurement is still possible.

Error codes

E 30	The measured result is far outside the defined tolerances. Reasons: e.g. dirty optics.
E 140/ E 160	Too much light reaching the detector. Reasons: e.g. extraneous light. Place the cap on the sample chamber and repeat measurement.
E 173	The standards have been measured in the wrong sequence. Please repeat calibration.
E 177	User calibration incorrect / erase Factory calibration is activated.
E 178	Factory calibration incorrect / erase

Repair and maintenance

If it is suspected that the instrument cannot be used without possibly imposing a danger, it should be turned off immediately and the potential danger be identified before the equipment is used again. The safety of the user may be diminished if the instrument

- exhibits visible damages.
- no longer works as specified.
- was stored for an extended period in unsuitable conditions.

In case of doubt, send the instrument to the manufacturer for repair or maintenance.

User replaceable accessories and spare parts

Accessories	cat. no.
BATTERY HIGH ENERGY 9V-BLOCK	VART4922121411
CLEANING BRUSH, 11 CM LENGTH	705-1095
PLASTIC STIRRING ROD, 10 CM LENGTH, LIGHT GREY	705-1096
MEASURING BEAKER, 100 ML	705-1098
SET OF 10 PLASTIC STIRRING RODS, 10 CM LENGTH, LIGHT GREY	705-1099
CLEANING CLOTH FOR CUVETTES	705-1102
SET OF 12 SEALING RINGS FOR ROUND VIALS 24 MM	705-1103
SET OF 12 SAMPLE CUVETTES WITH BLACK LID	705-1104

Technical service

Web Resources

Visit the VWR website at **www.vwr.com** for:

- Complete technical service contact information
- Access to the VWR Online Catalogue, and information about accessories and related products
- Additional product information and special offers

Contact us For information or technical assistance contact your local VWR representative or visit **www.vwr.com**.

Warranty

VWR warrants that this product will be free from defects in material and workmanship for a period of two (2) years from date of delivery. If a defect is present, VWR will, at its option and cost, repair, replace, or refund the purchase price of this product to the customer, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear. If the required maintenance and inspection services are not performed according to the manuals and any local regulations, such warranty turns invalid, except to the extent, the defect of the product is not due to such non-performance.

Items being returned must be insured by the customer against possible damage or loss. This warranty shall be limited to the aforementioned remedies. IT IS EXPRESSLY AGREED THAT THIS WARRANTY WILL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND IN LIEU OF THE WARRANTY OF MERCHANTABILITY.

Compliance with local laws and regulations

The customer is responsible for applying for and obtaining the necessary regulatory approvals or other authorisations necessary to run or use the Product in its local environment. VWR will not be held liable for any related omission or for not obtaining the required approval or authorisation, unless any refusal is due to a defect of the product.

Equipment disposal

This equipment is marked with the crossed out wheeled bin symbol to indicate that this equipment must not be disposed of with unsorted waste.

Instead it's your responsibility to correctly dispose of your equipment at lifecycle -end by handling it over to an authorized facility for separate collection and recycling. It's also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect from health hazards the persons involved in the disposal and recycling of the equipment.

For more information about where you can drop off your waste of equipment, please contact your local dealer from whom you originally purchased this equipment.

By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

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Obrigado

Short manual

Routine Measurement 1. Press ON key. 2. Rinse out a clean vial three times with the sample to be tested. 3. Place the vial in the sample chamber and align correctly X. 4. Put on the sample chamber cover (light shield). 5. Press Read key to start measurement. 6. Record the NTU value.	User-Calibration Press the MODE key and hold. Switch the unit on using the ON key. Allow the 3 decimal points to be displayed before releasing the MODE key. The I key allows for selection of the following menu point: User calibration (CAL). Confirm the selection with the MODE key. 1. Place the <0.1 NTU standard in the sample chamber, making sure that the marks X are aligned. Put the sample chamber cover on and press Read key. Reading starts automatically after count down. 3. Invert the 20 NTU standard and place it into the sample chamber, making sure that the marks X are aligned. Put the sample chamber cover on and press Read key. Reading starts automatically after count down. 4. Invert the 200 NTU standard and place it into the sample chamber, making sure that the marks X are aligned. Put the sample chamber cover on and press Read key. Reading starts automatically after count down. 5. Invert the 800 NTU standard and place it into the sample chamber, making sure that the marks X are aligned. Put the sample chamber cover on and press Read key. Reading starts automatically after count down. 6. Press I key to store.
--	--

GB

TBR-200_TBR-210_RRev_1.12.2020

Correct Storing and Handling Standards • TSK standards should always be stored in their original vial. • Store standards between 5 and 25°C (avoid prolonged exposure to temperatures exceeding 35°C). • Store away from direct sunlight. • Always allow the standard to acclimatise to ambient instrument temperature before use (not to exceed 35°C). • TSK Standards have a shelf life of minimum 12 month if stored correctly.	Using <0.1 NTU Standard After delivery allow the vial to stand for least 24 hours before use. Never shake or invert this <0.1 NTU standard. If the standard has been shaken it can take several hours for all the bubbles dissipate and we suggest the vial is left for at least 24 hours as outlined above. If the standard has been inverted accidentally wait at least 15 minutes before using. Preparing and Using Standards – infrequent use Note: These instructions apply to all standards excepted the <0.1 NTU standard. If you are using the <0.1 NTU standard follow the instructions above. These instructions apply to standards that have been sitting undisturbed for more than one week or to new standards. 1. Shake the standard vigorously for 2-3 minutes. 2. Let the vial stand undisturbed for 5 minutes. 3. Invert the vial between 5-10 times. 4. Immediately place the vial in the sample chamber and wait for 1 minute (count down). Preparing and Using Standards – frequent use Note: These instructions apply to all standards excepted the <0.1 NTU standard. If you are using the <0.1 NTU standard follow the instructions given above. 1. These instructions apply to frequently used standards (daily or weekly usage). 2. Invert the vial 10 times. 3. Immediately place the vial in the sample chamber and wait for 1 minute (count down).
--	--

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VWR International Ltd