




Comprehensive solutions Simplified workflows

Thermo Scientific Insight Pro Software
for Evolution UV-Vis Spectrophotometers

Accelerate your analysis

No matter how you use your UV-Visible spectrophotometer, from simple measurements to the most sophisticated research studies, Thermo Scientific™ Insight™ Pro Software helps you reach answers consistently and with unprecedented speed. Build your own QA/QC procedures to remove costly manual operations and

increase your laboratory productivity. Delivered in the language you prefer with the simple, effective tools you need to achieve current regulatory compliance or accelerate your research studies; Insight Pro Software enhances your Thermo Scientific™ Evolution™ UV-Vis Spectrophotometer.



Traditional UV-Visible applications are redefined with Insight Pro Software, providing simple and powerful tools to tackle even your most challenging assignments.

Applications	Fixed	Scan	Quant	Rate	Live display
Key features	Sample correction	Automated peak/valley picking	Automated pass/fail analysis	Extend run times	Walk-up simplicity
	Data trending	Advanced spectral processing	Simplified report generation	Temperature ramping	One-touch measurements
	Automated result calculations			Multi-stage rate calculations	Simplified scan

Simplify your workflow with automated calculations and result analysis


Minimize errors and increase productivity in your laboratory by letting Insight Pro Software perform method calculations and analyses for you.

- Calculate standard concentration values from measured weights/volumes
- Correct results for standard and sample properties and preparations steps
- Apply user-defined equations to sample results
- Perform simple statistical analyses including mean, standard deviation and %RSD
- Determine Pass/Fail results according to your method specifications

Personalize your home screen the way you work

A customizable home screen conveniently puts your assays and methods right in front of you for immediate access. Organize your laboratory with custom user groups that allow you to hide or


display applications and methods to match each group's needs and proficiency. Users can go directly to the methods they use every day and start collecting data instantly.



Move your procedures from paper to automation

Setting up your experiments is easy with the step-by-step design of Insight Pro Software. Organized to think the way you do, Insight Pro Software intuitively guides you through the set-up process.


- Sequentially displayed settings tabs move you quickly from start to finish
- Experimental options and settings are logically grouped so you can find them easily
- Built-in intelligence pre-fills dependent selections for added simplicity
- Informational tips available throughout the software provide you with added assistance when you need it



Customize reports to fit your needs

Composing detailed reports is effortless with Insight 2 report configuration options. Select the information you want to include and a customized report is generated automatically. Or, export your results for later retrieval in a format that is convenient for you.

- Define global report settings to automatically include your company's name and logo, time and date stamps, footers and page numbering on every report
- Choose to display your sample and method information, tables, graphs and experimental notes as desired for added flexibility
- Quickly adjust your tables and graphs to maximize space with a convenient Fit to Page Width option
- Export your results quickly and conveniently as *.xml, *.csv or *.tsv files for use in other software platforms



Integrate accessories into your workflow

Complete control of automated accessories provides superior performance for advanced sample handling requirements. Calibration and alignment tools ensure your accessories are performing optimally every time.

- Control and monitor the temperature of your reactions in up to 8 positions with supported temperature control accessories
- Gain full confidence in your results with hands-free performance verification and accessory calibration modules for your Thermo Scientific™ Evolution™ Smart Accessories.
- Perform your high-throughput activities with ease using a programmable sipper, cell changer or autosampler accessory




Choose the best fit for your application from QA/QC to research

Powerful analytical tools embedded in Insight Pro Software empower users to interpret analytical information instantly. Specialized display options let you view results the way you want for more effective data analysis.


Streamline QA/QC analysis

- Automated result calculation options eliminate your manual calculations and reduce user error
- Detailed graphical displays of calibration curves and sample run charts provide quick, visual checks of your system suitability and data trends
- Control limits and automated Pass/Fail interpretation provide you immediate compliance results
- Perform multiple calculations and extract parameters in scanning quant to create a complete analytical picture, improving your method development experience




Simplify your research applications

- A multitude of data collection option and graphing features are available to you for maximum flexibility in your experiments
- Pre-programmed peak analysis and convenient spectral processing tools aid in your sample identification and quantitation for scanning measurements
- Multi-staged curve fitting and consecutive reaction modes provide you with comprehensive data fitting options for your kinetics experiments
- Sophisticated math analysis functions including derivatives and smoothing let you display data in the format you need



Quickly measure and preview your sample

The Live Display feature in Insight Pro Software offers walk-up simplicity for real-time single wavelength measurements or quick identification of a sample peak. Display your results in absorbance or transmittance mode and print them for your experimental records.



Ensure consistency of your instrument performance

System Performance Verification Report

Company name: Thermo Scientific
Test Name: Wavelength Accuracy (Helium cad)

Measurement Description	High Limit	Low Limit	Measured	Result
Wavelength of 660.31 nm line	661.21	659.21	660.21	Pass
Wavelength of 544.81 nm line	545.71	543.71	544.71	Pass
Wavelength of 431.41 nm line	432.31	430.31	431.31	Pass
Wavelength of 365.31 nm line	366.21	364.21	365.21	Pass
Wavelength of 279.76 nm line	280.66	278.66	279.66	Pass
Wavelength of 241.31 nm line	242.21	240.21	241.21	Pass

Company name: Thermo Scientific
Test Name: Wavelength Reproducibility (Helium cad)

Measurement Description	High Limit	Low Limit	Measured	Result
Standard deviation of 561 nm peak	0.01	0.00	0.00	Pass

Company name: Thermo Scientific
Test Name: Resolution (Toluene/Benzene)

Measurement Description	High Limit	Low Limit	Measured	Result
Resolution (240 nm peak)	2.11	1.61	2.11	Pass

Company name: Thermo Scientific
Test Name: Stray Light (XL, 220 nm)

Measurement Description	High Limit	Low Limit	Measured	Result
Stray Light at 220 nm (%T)	0.05	0.00	0.00	Pass

Company name: Thermo Scientific
Test Name: Fluorescence Accuracy (Dibenzanthracene)

Multiple results files included in report







Ensure the accuracy and reliability of your QA/QC data while improving the efficiency of your laboratory with hands-free performance verification. Combine one of our Thermo Scientific™ Calibration Validation Carousels (CVCs) with Insight Pro Software to minimize errors and instrument downtime.

- Automated testing can save more than four hours of your analyst's time, improving the productivity of your laboratory
- Automatic CVC identification and file matching eliminate your transcription activities and return results that are ready for sign-off when tests are complete
- Traceable standards provide you with compliance to industry guidelines and regulations for both industrial and pharmaceutical laboratories

Enforce data security of your electronic records

If your laboratory requires 21 CFR Part 11 compliance, Thermo Scientific™ Insight Pro Security Software is here to make your life easier. Using the same easy-to-use interface, Insight Pro Security Software combines security and data integrity assurance with the versatility required for your multi-user laboratory.

- Seamless integration with Windows security features ensures that any changes to files associated with Insight Pro Security Software are monitored and logged, even when the software is not running
- Simple, straightforward user maintenance with Thermo Scientific™ Security Administration (TSA) allows administrators to manage your users' access and privileges from one central server location



Learn more at thermofisher.com/insightpro

thermo scientific