Model No.

N° de modèle

Modellnummern

Modelo No.

Modello No.

**MFLX77200-12**

**MFLX77200-20**

OPERATING MANUAL:

**MASTERflex® L/S®**  
**Compact Pump Drives**

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**MASTERflex® L/S®**  
**Entrainements de**  
**Pompes Compacts**

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**MASTERflex® L/S®**  
**Compact Pumpenantriebe**

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
**MASTERflex® L/S®**  
**Propulsores de Bombas**  
**Compactos**

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# TABLE OF CONTENTS

Title	Page
SAFETY PRECAUTIONS .....	2
INTRODUCTION AND GENERAL DESCRIPTION .....	3-4
SETUP AND DRIVE OPERATION .....	5
MAINTENANCE .....	5
SPECIFICATIONS .....	5
WARRANTY .....	6
PRODUCT RETURN .....	6
TECHNICAL ASSISTANCE.....	6

## PUMP FOR LIQUIDS ORIGINAL INSTRUCTIONS

### SAFETY PRECAUTIONS

**CAUTION:** Use of Pump Heads other than those specified will result in unsatisfactory performance, damage to the drive unit and voiding of applicable warranties.

The drive system can only be used with one Pump Head at a time. Use of this drive unit with multiple Pump Heads will result in unsatisfactory performance, damage to the drive unit and voiding of applicable warranties.

When mounting the Pump Head to the drive unit, the mounting screws should only be tightened fingertight. Excessive tightening of the mounting screws will result in binding and unsatisfactory performance.

**WARNING:** Be sure the Pump Drive is turned off before proceeding. Loose clothing could be caught in the drive.

### WARNING: PRODUCT USE LIMITATION

These products are not designed for, nor intended for use in patient-connected applications, including, but not limited to, medical and dental use and, accordingly, have not been submitted for FDA approval.

## INTRODUCTION AND GENERAL DESCRIPTION

The **MASTERFLEX® L/S® COMPACT Pump Drives** are designed as an inexpensive variable speed drive which accepts **MASTERFLEX L/S STANDARD, QUICK LOAD®** and **EASY-LOAD®** Pump Heads. These 200 rpm drives accommodate one Pump Head and one tube at controlled speed as low as 30 rpm. All units operate from an external Power Supply. The AC models are supplied with a Universal Power Supply which provides a DC output for connection to the Pump Drive.

The single-turn adjustable Speed Control, *Figure 2*, provides variable flow operation. The green **POWER** On indicator lights whenever the pump is operating. The Power On/Direction Switch turns power on when either clockwise or counterclockwise pump rotor direction is selected. The **MAX** Button is used for priming and purging and operates the pump at maximum speed while depressed.


Figure 1. MASTERFLEX L/S Compact Pump Drive


Figure 2. Pump Drive—Front Panel

The rear panel, *Figure 3*, contains a DC power input jack for connection of primary power and a 4-terminal barrier strip for connection of remote start/stop and for a DC backup supply.

**Remote Start/Stop:** Terminals 3 and 4 on the rear panel terminal strip, (*Figure 3*) are used for remote start/stop operation. Pump direction and speed are not remotely controllable. In non-remote operation, these terminals are connected together by a shorting bar. For remote control by switch closure, remove the shorting bar and connect the two terminals of an isolated remote control switch or relay to terminals 3 and 4. A closure of the remote switch contacts will turn the pump system on. Opening the contact will turn the pump system off.

**WARNING:** *Be sure the Pump Drive is turned off before proceeding. Loose clothing could be caught in the drive.*

**Pump Heads Accepted:** Only the MASTERFLEX L/S STANDARD, QUICK LOAD and EASY-LOAD are acceptable for use with this pump drive. Only one Pump Head at a time can be used with this drive unit.

**CAUTIONS:** *Use of Pump Heads other than those specified will result in unsatisfactory performance, damage to the drive unit and voiding of applicable warranties.*

*The drive system can only be used with one Pump Head at a time. Use of this drive unit with multiple Pump Heads will result in unsatisfactory performance, damage to the drive unit and voiding of applicable warranties.*

*When mounting the Pump Head to the drive unit, the mounting screws should only be tightened fingertight. Excessive tightening of the mounting screws will result in binding and unsatisfactory performance.*

**Tubing Accepted:** The Pump Heads accepted for use with drive can be used with MASTERFLEX precision tubing sizes L/S® 13, L/S® 14 and L/S® 16 in all formulations. For L/S® 17 and L/S® 25 size tubing, use silicone or C-FLEX® formulations only.

Use only MASTERFLEX precision tubing with MASTERFLEX pumps to ensure optimum performance. Use of other tubing may void applicable warranties.


Figure 3. Pump — Rear Panel

## SETUP AND DRIVE OPERATION

1. Unpack the drive: save packaging material until proper product operation has been verified.
2. Select and mount Pump Head and load tubing (see Pump Head manual).
3. For 77200-20 plug in Universal Power supply and connect to the back of the drive. (See Figure #3.)
4. For 77200-12 supply 11.0 to 15V DC to terminal strip positions 1 and 2. (See Figure #3.)
5. The drive system can be set to operate in either a clockwise or a counterclockwise direction. The same control used to select direction also turns power on or off. To operate in a clockwise direction place direction control switch in the "up" position. The POWER indicator should light and the fluid will flow from right to left. To operate in a counterclockwise direction, place direction switch in the "down" position. The POWER indicator should light and the fluid will flow from left to right. To turn pump off, place switch in the "center" position.
6. The pump speed is controlled by the variable speed control. Turning the control clockwise increases the speed. Tubing life is decreased with increased operating speed. Speed can be controlled up to a 6:1 range.

## MAINTENANCE

**Cleaning:** Keep the drive enclosure clean with mild detergents. Never immerse nor use excessive fluid.

**Motor/ Brush Replacement:** The motor used in this inexpensive drive is not serviceable or repairable and has brushes which cannot be replaced. The motor is rated for a continuous duty life of 3000 operating hours minimum. Operation at other than continuous duty cycles will result in longer motor life. If motor failure for whatever reason occurs, the motor must be replaced.

## SPECIFICATIONS

### Output:

<b>Operating Speed:</b>	35 to 200 rpm
<b>Number of Pump Heads:</b>	1
<b>Torque Load:</b>	20 oz-in, maximum
<b>Direction of rotation:</b>	Clockwise and Counterclockwise

### Input:

<b>Operating Voltage/Frequency:</b>	
Model 77200-20	115V AC nominal, 50/60 Hz (90–130VAC) @500 mA AC 230V AC nominal, 50/60 Hz (190–260VAC) @300 mA AC
Model 77200-12	13.5V DC nominal (11.0–15.0V DC) @ 2.4 A DC
<b>Installation Category:</b>	
Model 77200-12 & -20	Category I per IEC664 (Signal Level)
Remote Start/Stop:	Contact closure connection at terminal strip contacts 3 and 4

### Construction:

Dimensions (L x W x H):	5.750 in x 5.14 in x 4.14 in (14.6 cm x 13.0 cm x 10.5 cm)
Weight:	1.88 pounds (0.85 kg)
Color:	Black
Material:	Painted Steel
Enclosure Rating:	IP22 per IEC529

### Environment:

Operating Temperature:	32°F to 104°F (0°C to 40°C)
Storage Temperature:	-49°F to 149°F (-45°C to 65°C)
Humidity (Non-cond):	10% to 90%
Altitude:	6600 ft (2000 m)
Pollution Degree:	Pollution Degree 2 per IEC664 (Indoor usage — lab, office)

### Compliance (For CE Mark):

EN61326-1/A2: 2001 (EMC Directive)  
 Converter is UL listed and CSA approved.  
 Regulatory agency specifications not applicable to the balance of the unit due to low voltage.

## **WARRANTY**

***Use only MASTERFLEX precision tubing with MASTERFLEX pumps to ensure optimum performance. Use of other tubing may void applicable warranties.***

The Manufacturer warrants this product to be free from significant deviations from published specifications. If repair or adjustment is necessary within the warranty period, the problem will be corrected at no charge, if it is not due to misuse or abuse on your part, as determined by the Manufacturer. Repair costs outside the warranty period, or those resulting from product misuse or abuse, may be invoiced to you.


***The warranty period for this product is one (1) year from the date of purchase.***

## **PRODUCT RETURN**

To limit charges and delays, contact the seller or Manufacturer for authorization and shipping instructions before returning the product, either within or outside of the warranty period. When returning the product, please state the reason for the return. For your protection, pack the product carefully and insure it against possible damage or loss. Any damages resulting from improper packaging are your responsibility.

## **TECHNICAL ASSISTANCE**

If you have any questions about the use of this product, contact the Manufacturer or authorized seller.

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Masterflex LLC