ADP400 Series Polarimeters & ADS400 Saccharimeters

Thank you for purchasing this instrument. In order to ensure that this product provides many years of service please follow the guidance in this document. To register your instrument with the manufacturer and print a warranty certificate visit: www.bellinghamandstanley.com

Unpacking

Carefully remove all of the packing material. It is recommended that the box and other packing materials are retained so that, should the need arise, the instrument can be safely returned to the supplier for service. Check that all parts listed below are present and that no transit damage has occurred. If any are damaged or missing, contact the supplier immediately.

Contents list

Quantity	Item	Code
1	ADP430 Polarimeter Or ADS435 Saccharimeter Or ADS438 Saccharimeter Or ADP450 XPC Peltier Polarimeter Or ADS455 XPC Peltier Saccharimeter Or ADS458 XPC Peltier Saccharimeter	37-430 37-425 37-485 37-450 37-445 37-495
1	Mains leads	55-109
1	Power supply	55-105
1	Quick Start Guide	37-234
1	Key-fob type ID tags (pack of 3)	22-071
1	Polarimeter Tube	See below
1	XPC Adaptor	See below

Part code numbers:

Order Code	Tube	XPC Adaptor
37-25 / 37-30 / 37-85	35-47 Centre fill tube (200mm)	n/a
37-26 / 37-86	36-57 Funnel flow type tube (100mm)	n/a
37-27 / 37-87	36-58 Funnel flow type tube (200mm)	n/a
37-45 / 37-50 / 37-95	35-91 Stainless steel tube (200mm)	35-96 200mm
37-51	35-90 Stainless steel tube (100mm)	35-95 100mm
37-46 / 37-96	35-92 Stainless tube, funnel filler (100mm)	35-95 100mm
37-47 / 37-97	35-93 Stainless tube, funnel filler (200mm)	35-96 200mm
37-98	35-13 Glass tube, PEEK end (200mm)	35-96 200mm

XPC Models: Place the adaptor on the instrument rails so that the flat surface is in contact with the temperature controlled plate.

Positioning the system

- Dry and indoors.
- Capable of holding the instrument's weight.
- Away from draughty or hot equipment like fans or heaters.
- Out of direct sunlight or strong ambient light.
- Away from potential sources of interference, such as RFI generating equipment.
- Within reach of a power point.
- Not using a power circuit that also has large motors or noise generating equipment connected to it.

Do not block the air vents.

Mains connection

The power supply adapter is supplied with a moulded mains cord and plug to suit one of several socket types. For UK leads, replace the fuse only with the type indicated on the plug.

Voltage	110 to 230 V ~ ±10%
Frequency	50 to 60 Hz
Maximum current	2 A

RISK OF ELECTRIC SHOCK:

- For indoor use only.
- Must be kept dry.
- Disconnect the equipment from the mains supply before unplugging the mains lead from the instrument.

- Do not cover, designed to operate with free air convection.
- No cleaning required.
- Do not jet wash.

Switching on and off

Plug in the power supply and switch on the mains supply, after startup the instrument will show a loading screen indicating the software version. Once loaded, the instrument will be in standby mode. To switch on, press the Power softkey.

Allow the instrument to stabilise for at least 30 minutes after switching on before taking accurate readings.

To switch to standby mode press and hold the Power softkey for several seconds and select Yes on the confirmation screen. It is then safe to turn off the mains supply and disconnect the power lead.

Navigating through menus

The softkey functions change depending upon what inputs are required at that time. If a softkey has no text then that button is disabled.

Menu items can be selected by either:


- Pressing the Up/Down softkeys until it is highlighted and then pressing Select.
- Pressing the item number key.

Alternatively, a USB keyboard, available from Bellingham + Stanley, can be used to select items and navigate menus.

If the menu item contains a list of options to choose from then either:

- Select the menu item and then choose the option from the displayed list.
- Press the Forward/Back arrow keys to step through the options.

Menu items that can be modified by using the arrow keys are identified by arrows being displayed next to the menu items value.

Return to the Measurement Display by pressing Quit to step back through the menus.

Setup Wizard

When the instrument is switched on for the first time, a Setup Wizard will step through a number of questions so that the system can be configured for a particular application. First select the desired language by pressing the associated number key. Press the Next softkey to access each setting screen:

Copy user guides and PC software to USB: Makes copies of the full User Guides and PC software compatible with the instrument to a removable USB storage device. Insert a USB removable storage device before selecting. Full User Guides provide detailed instructions on instrument operation.

Set the standard measurement conditions: Make changes to the Measurement Settings.


Set clock: Set the instrument to local time.

Set the printer options: Specify how reading results will be printed.

When the desired conditions are set in each screen press the Quit softkey to proceed.

ADP400 Series Polarimeters & ADS400 Saccharimeters

Measurement Display
From the Measurement Display readings can be viewed and saved.

1. Reading result: Readings will be displayed here. Warning and alert messages will also appear here. To take a reading lift the chamber covers and place the sample on the adaptor/trails. Close the chamber.

The instrument will take continuous readings. The reading result will be displayed in grey until stable, then turn black.

- 2. Temperature control (XPC models):** the stability status, set temperature and current temperature.
- 3. Time:** The current time.
- 4. Progress bar:** Demonstrates the reading response time. Does not appear if stability active.
- 5. Temperature:** The temperature when the reading was taken.
- 6. Optical Density (OD):** The optical density is a measure of the light and the sample optical density.
- 7. Reading Configuration:** Displays the measurement settings for the currently displayed reading.

Softkeys

- Zero button:** Zeros the instrument.
- Record button:** Pushing the button will record the result to the Reading Log or printer depending on the print settings.
- Mode button:** Gives access to edit all Basic and Advanced reading settings e.g. scale and temperature compensation.
- Menu button:** Opens the instrument's Main Menu from where Calibration can be performed, saved results downloaded and all instrument settings can be changed.

This symbol is an internationally agreed indicator that the product bearing it should not be disposed of as general waste or garbage which might end up in landfill sites, but should instead be sent for special processing and/or recycling in those countries where appropriate legislation and facilities are in place.

Full User Guides can be downloaded from the Instrument's Help Menu or Setup Wizard

Caution or warning.

Electrical hazard.

Menu Structure

1. Calibration	1. Span (Top Calibration)	1. Add new...	1. Settings...	1. Use Methods	3. Selection Type	1. Add new...	2. Settings...	1. Protection	1. basic 5. cfr	2. low	3. medium	4. high	1. Setup Wizard	2. Information	3. Contact	4. administrator	3. calibrator	4. operator	1. fast (2s)	2. medium (20s)	3. slow (40s)	1. From Library	1. Add new...	2. Scales	1. From Data	2. From Data	3. Temp. Comps.	1. Add new...	1. From Data	4. Recording	1. Mode	1. none 5. statistics	2. single 6. multi	3. save	4. single + save	1. none 5. network	2. serial	3. usb	4. pdf	1. Fields	1. serial	2. none	3. csv	4. xml	1. Location	2. Sample ID	1. off	2. number	3. numberinc	4. text	5. dateinc	1. off	2. number	3. text	4. login	1. none	2. batch	3. user	4. batch and user	1. Baud Rate	2. Word Length	2. Network	1. DHCP	2. IP Address	3. Subnet	Mask	4. Default	Gateway	1. Language	2. Time / Date	1. dd/mm/yy	2. mm/dd/yy	1. no	2. Summer Time	2. Yes	3. Time Server	4. Set Clock...	4. UTC Offset	1. Brightness	2. Remote Display	3. Setup	4. Clone From Backup	5. Temperature Control (XPC models only)	1. Mode	1. off	2. on	2. Tick Criteria
-----------------------	---------------------------	---------------	----------------	----------------	-------------------	---------------	----------------	---------------	-----------------	--------	-----------	---------	-----------------	----------------	------------	------------------	---------------	-------------	--------------	-----------------	---------------	-----------------	---------------	-----------	--------------	--------------	-----------------	---------------	--------------	--------------	---------	-----------------------	--------------------	---------	------------------	--------------------	-----------	--------	--------	-----------	-----------	---------	--------	--------	-------------	--------------	--------	-----------	--------------	---------	------------	--------	-----------	---------	----------	---------	----------	---------	-------------------	--------------	----------------	------------	---------	---------------	-----------	------	------------	---------	-------------	----------------	-------------	-------------	-------	----------------	--------	----------------	-----------------	---------------	---------------	-------------------	----------	----------------------	--	---------	--------	-------	------------------

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