# Whatman™ Puradisc™ 25 GF/F Disposable Filter Device

## Product Information sheet

### Warning

**For research use only.**

**Not recommended or intended for diagnosis of disease in humans or animals.**

**Do not use internally or externally in humans or animals.**

The Puradisc 25 GF/F disposable filtration device was designed to provide a disposable filter with a high flow rate and long life, for fine filtration. This filter is ideal for removing fine particulates down to 0.7 µm. Unlike membrane filters with a comparable retention value, this product has a rapid flow rate and a high loading capacity.

Puradisc 25 GF/F has a polypropylene (PP) housing with Glass Microfiber filter medium. This filter medium is very effective with fine particulates such as precipitated proteins. It can also be used as a prefilter in conjunction with smaller pore size filtration devices.

Disposable filtration devices provide great labor saving efficiency while ensuring consistent filtration when compared to hand assembled reusable filter housings.

This document provides general information on the products listed. The specifications in the Technical Data section are intended to provide the basis for establishing functional use, as well as setting the quality assurance test performance levels.

- Whatman Brand GF/F Filter Media
- Extremely High Flow Rate
- Long Life - High Loading Capacity
- Chemical & Physical Stability
- Binder Free Glass Microfiber
- Particulate Retention Down to 0.7 µm
- Autoclavable
- Inlet: Female Luer Lock (FLL)  
Outlet: Male Slip Luer (ML)
- No adhesives used in construction

## Puradisc 25 GF/F - 25 mm Filters

Catalog Number	Product	Pore Size	Qty./Pkg.
6825-2517	Puradisc 25 GF/F	0.7	50
6753-2510	Puradisc 25 GF/F	0.7	200

## Typical Applications for Puradisc 25 GF/F Filter Devices

- TCLP Samples
- Prefiltration: Extends Life of Other Filter Media or Membranes
- Clarification of 'Difficult' Biochemical Solutions and Fluids
- Fine Filtration: Liquids requiring particulate Removal Down to 0.7 µm For Gases Retention Efficiency is Significantly Greater.
- Filtration of Precipitated Proteins
- Liquid Scintillation Counting
- Radioimmunoassay
- Water and Wastewater Analysis
- Sample Preparation

## Operating Instructions:

**Safety:** When considering the special factors of your application, consult the Technical Data to determine correctness of use. Do not exceed the pressure, temperature or chemical compatibility recommendations. High pressures can be obtained when using syringes. The smaller the syringe the higher the pressure that can be generated. As a general guide, the following pressures can be obtained by hand with the syringes indicated: 20 mL, 80 psi; 10 mL, 140 psi; 5 mL, 180 psi; 3 mL, 200 psi; 1 mL, 250 psi. Each user should determine the pressure they can generate by hand with a specific size syringe and take appropriate safety precautions not to exceed the recommended rating for the device used. If these limitations are exceeded, bursting of the device may occur resulting in loss of sample or personal injury.

**GF/F Media Considerations:** The glass microfiber filter medium is produced entirely from the finest grades of borosilicate glass microfiber and contains no binders either as manufacturing aids or as wet strengthening additives. Glass microfiber filters are generally resistant to weakening or disruption of the fibrous matrix by inorganic or organic liquids.

**Autoclaving:** Autoclave at 121°C (131°C max.) for 20 minutes.

### To use with a syringe: (3 mL, 5 mL, 10 mL)

- 1) Fill the syringe with the solution to be filtered.
- 2) Secure the filled syringe to the FLL on the inlet with a twisting motion.
- 3) With outlet pointed upward, gradually apply thumb pressure to the syringe plunger to initiate flow.
- 4) Change filters when flow becomes too slow or resistance becomes excessive.

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

### Technical Data: Puradisc 25 GF/F Disposable Filter Devices

<b>Dimensions:</b>	28.3 mm (1.1 in) x 22.9 mm (0.9 in)
<b>Weight:</b>	2.7 grams
<b>Filtration Area:</b>	4.2 cm <sup>2</sup>
<b>Maximum Pressure:</b>	5.2 bar (75 psi)
<b>Housing:</b>	Polypropylene
<b>Hold Up Volume:</b>	Full Housing 0.16 mL With air purge <0.1 mL
<b>Filter Media:</b>	Whatman GF/F Glass Microfiber
<b>Flow Direction:</b>	Flow should enter from the inlet
<b>Connectors:</b>	Inlet - Female Luer Lock (FLL) Outlet - Male Slip Luer (ML)
<b>Autoclaving:</b>	Autoclave at 121°C (131°C max) at 15 psi for 20 minutes.

### Water Flow Rates

### Liquid Flow Rate (mL/minute)

Catalog No.	Product	Media	Nominal Pore Size Microns	Water Flow Rate <sup>(1)</sup> at 1.0 bar (14.5 psi)	
6825-2517	Puradisc 25 GF/F	Glass Microfiber	0.7	270 mL/min	50
6825-2527	Puradisc 25 GF/F	Glass Microfiber	0.7	270 mL/min	200

(1) Typical values

## Legal.

GE, imagination at work and GE Monogram are trademarks of General Electric Company.

Whatman and Puradisc are trademarks of GE Healthcare companies.

© 2008–2013 General Electric Company – All rights reserved.  
Previously published Jul 2008.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.

For your local office contact information, visit  
[www.gelifesciences.com/contact](http://www.gelifesciences.com/contact)

GE Healthcare UK Limited  
Amersham Place  
Little Chalfont, Buckinghamshire,  
HP7 9NA, UK

<http://www.gelifesciences.com>

### **GE Healthcare offices:**

GE Healthcare Bio-Sciences AB  
Björkgatan 30, 751 84 Uppsala,  
Sweden

GE Healthcare Europe GmbH  
Munzinger Strasse 5, D-79111 Freiburg,  
Germany

GE Healthcare Bio-Sciences Corp.  
800 Centennial Avenue, P.O. Box 1327,  
Piscataway, NJ 08855-1327,  
USA

GE Healthcare Japan Corporation  
Sanken Bldg. 3-25-1, Hyakunincho,  
Shinjuku-ku, Tokyo 169-0073,  
Japan