

Superdex™ 75 Increase, new generation size exclusion chromatography (SEC) columns for recombinant proteins


GE expands the range of new generation SEC columns

Our new generation columns called “Increase” offer higher resolution and shorter run times than their predecessors.

Superdex 75 Increase is designed for small-scale preparative purification and analysis of recombinant tagged or other proteins with molecular weights (M_r) from 3000 up to 70 000.

Comprising smaller, more rigid beads with a narrower particle size distribution and a higher selectivity, Superdex 75 Increase columns offer:

- runtime reduced down to one third with maintained resolution for results faster
- up to 50% higher resolution for improved purity and analysis results
- higher lot-to-lot consistency




Superdex 75 Increase enables reduced runtime down to one third compared with Superdex 75.



With the addition of Superdex 75 Increase, you will have three choices of the new generation SEC media for small-scale preparative purification and analysis (sample volume 4 to 500 µl):



Superdex 75 Increase
Recombinant tagged proteins
Fractionation range
 M_r 3000 to 70000




Superdex 200 Increase
MAb and other antibodies
Fractionation range
 M_r 10 000 to 600 000




Superose™ 6 Increase
Larger proteins and protein complexes*
Fractionation range
 M_r 5000 to 5 000 000

The fractionation ranges of the three different chromatography media complement each other

Superose 6 Increase 10/300 GL



Superdex 200 Increase 10/300 GL




Sample:


- 1 IgM (M_r 970 000)
- 2 Thyroglobulin (M_r 669 000)
- 3 Ferritin (M_r 440 000)
- 4 BSA (M_r 66 000)
- 5 Myoglobin (M_r 17 600)
- 6 Vitamin B₁₂ (M_r 1300)

Comparison between Superose 6 Increase and Superdex 200 Increase. Superose 6 Increase gives higher resolution for proteins with $M_r \geq 440\,000$, while proteins of lower molecular weight are better resolved on Superdex 200 Increase.

Superdex 200 Increase 10/300 GL



Superdex 75 Increase 10/300 GL



Sample:

- 1 Thyroglobulin (M_r 669 000)
- 2 Aldolase (M_r 158 000)
- 3 Conalbumin (M_r 75 000)
- 4 Ovalbumin (M_r 44 000)
- 5 Carbonic anhydrase (M_r 29 000)
- 6 Cytochrome C (M_r 12 300)
- 7 Aprotinin (M_r 6500)
- 8 Vitamin B₁₂ (M_r 1300)

Comparison between Superdex 200 Increase and Superdex 75 Increase. Superdex 75 Increase gives higher resolution of proteins with $M_r \leq 44\,000$. Superdex 200 Increase shows better resolution for proteins above this molecular weight.

Ordering information

Dimensions (mm)	Superdex 75 Increase VWR Cat. No.	Superdex 200 Increase VWR Cat. No.	Superose 6 Increase VWR Cat. No.
10 × 300	75799-300	89497-272	10192-228
5 × 150	75799-302	89497-274	10192-230
3.2 × 300	75799-304	89497-276	10192-226

For additional information, please visit vwr.com

GE, GE monogram, Superdex, and Superose are trademarks of General Electric Company.
© 2016 General Electric Company. First published Jan. 2016.
GE Healthcare UK Limited, Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA, UK.
29186446 AC 06/2016



Distributor
GE Healthcare



1.800.932.5000
vwr.com

Prices and product details are current when published; subject to change without notice. Certain products may be limited by federal, state, provincial, or local regulations. VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in US dollars unless otherwise noted. Offers valid in US, void where prohibited by law or company policy, while supplies last. VWR, the VWR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks of VWR International, LLC and its related companies. All other marks referenced are registered by their respective owner(s). Visit vwr.com to view our privacy policy, trademark owners and additional disclaimers. ©2016 VWR International, LLC. All rights reserved.