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


# 2

## 50 50 T EN

For more detailed information, refer to the meter documentation.

**Caution**  
 ⚠ Personal injury hazard. Glass components can break. Handle with care to prevent cuts.  
 ⚠ Chemical exposure hazard. Refer to the current material safety data sheets (MSDS) for safety protocols.

**Introduction**  
 This probe is a combination gel-filled electrode with a built-in temperature sensor for pH measurements in general aqueous samples.

**Technical information**

Range	0 to 14 pH
Junction	Ceramic
Reference	Ag/AgCl
Temperature	Pt 1000
Minimum depth	18 mm (0.7 in.)
Operating temp.	0 to 80 °C (32 to 176 °F)
Limits	Low-conductivity, dirty or viscous solutions and reducing sample types
Connector type	MP-5

**Connect the probe**  
 ⚠ Do not turn or twist the probe connector.

**1 Calibration**  
 ⚠ Air bubbles under the probe tip when submerged can cause slow stabilization or error in measurement. If bubbles are present, gently shake the probe until bubbles are removed.

For 1 or 2-point pH calibration: push ✓ after 1 or 2 standards.  
**2 Sample measurement**  
 Continuous: push ✓ 2 times.  
 Change the parameter: push and hold ✓ during stabilization.

**Maintenance and storage**  
 Cleaning: soak the probe in a cleaning solution. Rinse with deionized water.  
 Hydration: soak the dry probe in a dilute HCl solution for several hours.  
 Storage: put the probe in the protector with KCl 3 M solution.

**Troubleshooting**  
 Be sure that the correct probe type is used. Do not hold the tubes by hand to prevent temperature changes.

**Warranty**  
 6 months for manufacturing faults. This warranty does not cover improper use or wear.

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