Leica DM IL LED

Leica DM IL LED Cellfactory

Leica DM IL LED Fluo

Leica DM IL LED Fluo Cellfactory

Leica DM IL M LED

Instructions for use • Bedienungsanleitung • Mode d'emploi

Leica Microsystems CMS GmbH, Instructions for uses 11 933 984, Revision 2.1, 2015-08-15

Living up to Life

*Leica*  
MICROSYSTEMS

Revision 2.1, published 2015-08-15 by:

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Leica Microsystems CMS GmbH

Ernst Leitz-Straße 17-37

D-35578 Wetzlar (Germany)

<http://www.leica-microsystems.com>

Responsible for contents:

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*Leica*  
MICROSYSTEMS

# Copyrights

All rights to this documentation are held by Leica Microsystems CMS GmbH. Reproduction of text or illustrations (in whole or in part) by print, photocopy, microfilm or other method (including electronic systems) is not allowed without express written permission from Leica Microsystems CMS GmbH.

The instructions contained in the following documentation reflect state-of-the-art technology. We have composed the texts and illustrations with great care. Still, we are always grateful for comments and suggestions regarding potential mistakes within this documentation.

The information in this manual is subject to modification at any time and without prior notification.

# Contents

<b>1. Important Information about this Manual</b>	<b>6</b>
1.1 Text symbols, pictograms and their meanings	6
<b>2. Intended Purpose of the Microscope</b>	<b>8</b>
<b>3. Safety Notes</b>	<b>10</b>
3.1 General safety notes	10
3.2 Electrical safety	10
3.3 Transport and storage	12
3.4 Notes on handling light sources	12
3.5 Notes on handling immersion oil	13
3.6 Notes on handling acids and bases	13
3.7 Disposal	13
<b>4. Overview of the Instrument</b>	<b>14</b>
<b>5. Unpacking</b>	<b>21</b>
<b>6. Assembling the Microscope</b>	<b>24</b>
6.1 Stand	24
6.2 Attaching the condensers	25
6.3 Inserting the transmitted light illumination carrier	25
6.4 Electrical connection of the transmitted light illumination carrier	27
6.5 Inserting the tubes	27
6.6 Eyepieces and graticules*	29
6.7 Objectives	29
6.8 Inserting the filter	30
6.9 Specimen stage	30
6.10 Light source for the transmitted light axis	31
6.11 Electrical connection of the microscope	31
<b>7. Assembling the Options</b>	<b>32</b>
7.1 Assembling the fluorescence filter blocks*	32
7.2 Inserting the phase contrast slider on the transmitted light illumination carrier*	33
7.3 Inserting the IMC slit diaphragm slider* on the transmitted light illumination carrier	34
7.4 Inserting the IMC modulator*	34
7.5 Assembling the 106z lamp housing* with mercury lamps	35
7.6 Assembling the lamp housing LH115 LED38	35
7.7 Leica EL6000*	39
7.8 Adapting cameras to binocular photo tubes*	40
7.9 Inserting the multiple discussion attachment*	41
7.10 Inserting the ErgoModule*	42
<b>8. Operation</b>	<b>43</b>
8.1 Basic settings for transmitted light	43
8.2 Objectives	46
8.3 Transmitted light	47
8.4 Phase contrast	49
8.5 Integrated Modulation Contrast (IMC)	50
8.6 Incident light fluorescence	54
8.6.1 Switching on and adjusting the halogen and mercury lamps in the 106z lamp housing*	54
<b>9. Troubleshooting</b>	<b>58</b>
<b>10. Care of the Microscope</b>	<b>62</b>
10.1 Dust cover	62
10.2 Cleaning	62
10.3 Handling acids and bases	63
<b>11. Technical Description</b>	<b>64</b>
<b>12. Index</b>	<b>72</b>
<b>13. EU Declaration of Conformity</b>	<b>73</b>

# 1. Important Information about this Manual

## Caution!

This Instructions for use is an essential component of the product; it must be read carefully before the product is assembled, put into operation or used, and must be kept for later reference.

This operating manual contains important instructions and information for the operational safety and maintenance of the microscope and accessories. Therefore, it must be kept in a safe place.

## 1.1 Text symbols, pictograms and their meanings

(1.2)

→ p.20

Numbers in parentheses such as “(1.2)”, correspond to illustrations (in the example, Figure 1, Item 2).

Numbers with an arrow, e.g. → p. 20, refer to a specific page of the manual.

WARNING indicates a hazard with a medium degree of risk that, if not avoided, can result in death or serious injury.

CAUTION indicates a hazard with a low degree of risk that, if not avoided, can result in minor or moderate injury.

## Caution!

Special safety instructions within this Instructions for use are indicated with the triangle symbol shown here, and have a gray background.

Notice! The microscope and accessories can be damaged when operated incorrectly.

## 1. Important Information about this Manual

Warning of hazardous electrical voltage! Risk of electrical shock!

Warning of optical radiation! Never look directly into the light beam! Wear safety goggles!

Warning of electromagnetic field

Danger due to hot surface.

Instructions for disposing of the instrument, its accessories and consumables.

Connection to ground!

Explanatory note.

\*

Item not contained in all configurations.

Device for in vitro diagnostics (IVD).

IVD manufacturing date, example 11 / 2011 for November 2011.

Manufacturing date, example 11 / 2011 for November 2011.

China RoHS 50 years EFUP  
(Environmentally friendly use period)

# 2. Intended Purpose of the Microscope

The microscopes Leica DM IL LED, Leica DM IL LED Cellfactory, Leica DM IL LED Fluo and Leica DM IL LED Fluo Cellfactory are inverted light microscope and intended for use as a general lab microscope for routine examinations of biological specimens.

Directions for use:

Inspection, numbers, classification, identification and monitoring of cell and tissue cultures, examining biological specimens, liquids and sediments that lead to information relating to physiological and/or pathological states, inborn anomalies and the monitoring of therapeutic measures.

The resulting measures are based on the interpretation of the medical professional.

Contraindications:

Not suitable for testing potentially infectious specimens.

### IVD

The Leica DM IL LED, Leica DM IL LED Cellfactory, Leica DM IL LED Fluo and Leica DM IL LED Fluo Cellfactory microscope complies with the EU Directive 98/79/EEC concerning in vitro diagnostics.

The DMIL M LED microscope is intended for materials science, geological or mineralogical examinations and applications in the laboratory area.

The requirements in EU Directives 2006/95/EC and 2014/35/EU regarding the safety of electrical equipment, and 2004/108/EC and 2014/30/EU regarding electromagnetic compatibility have been met for all the instruments mentioned above for use in industrial environments.

### Reasonably foreseeable misuse

The following are prohibited:

- Using the microscope for any purpose not in accordance with the Declaration of Conformity (e.g. using it as an in-vitro diagnostic product according to Council Directive 98/79/EC or as a medical product according to Council Directive 93/42/EEC)
- Making use of auxiliary equipment to use the area between the microscope stage and objectives as a clamp or holder (working like a vice).
- Operating the microscope in an inclined position.
- Cleaning the microscope in a way other than specified in the manual.
- Modifying the built-in safety circuits in the instrument.
- Allowing unauthorized personnel to open the instrument.
- Using cables that Leica has not provided or permitted.
- Using combinations with non-Leica components that go beyond the scope of the manual

### Caution!

The manufacturer assumes no liability for damage caused by, or any risks arising from using the microscopes for other purposes than those for which they are intended or not using them within the specifications of Leica Microsystems CMS GmbH. In such cases the declaration of conformity shall cease to be valid.

### Caution!

These (IVD) instruments are not intended for use in the patient environment defined by DIN VDE 0100-710. Nor are they designed to be combined with medical instruments in accordance with EN 60601-1. If a microscope is electrically connected to a medical instrument in accordance with EN 60601-1, the requirements defined in EN 60601-1-1 shall apply.

Not suitable for examining potentially infectious specimens.

Only trained personnel may operate this type of device.

### Notes on handling laser devices

In the standard design without additional laser safety precautions, the microscopes are not suitable for coupling laser radiation (such as to the camera ports), since this radiation poses a hazard for the user (eye injury in particular).

For use of the microscope with lasers, Leica Microsystems offers special microscope variants with additional safety features. Laser couplings require corresponding safety devices that have to be inspected and installed by trained personnel.

For further information, please contact your authorized Leica Microsystems representative.

# 3. Safety Notes

## 3.1 General safety notes

These safety class 1 instruments were built and tested in accordance with the harmonized standards EN 61010-1 / IEC 61010-1 / UL 61010-1, Safety regulations for electrical measuring, control, and laboratory devices, and EN 61010-2-101 / IEC 61010-2-101, Safety regulations for electrical measuring, control, and laboratory devices, part 2 special requirements concerning in vitro diagnostics (IVD) medical instruments.

They also comply with EN 62471 / IEC 62471, Photobiological safety of lamps and lamp systems, and belong to the risk group 1 (low risk).

### Caution!

In order to maintain this condition and to ensure safe operation, the user must follow the instructions and warnings contained in this operating manual.

### Caution!

The instruments and accessories described in this manual have been safety-tested and checked for possible hazards.

The responsible Leica affiliate or the main plant in Wetzlar must be consulted whenever the instrument is altered, modified or used in conjunction with non-Leica components that are outside of the scope of this manual.

Unauthorized alterations to the device or noncompliant use shall void all rights to any warranty claims and product liability!

## 3.2 Electrical safety

### Microscope

For indoor use only.

Supply voltage:

Frequency:

Power input:

Fuses:

Ambient temperature:

Relative humidity:

Overvoltage category:

Pollution degree:

100-240 V AC

50/60 Hz

max. 14 VA

0.63 A, slow-blowing,  
breaking capacity H,  
250 VAC, size 5x20  
mm

15-35°C

max. 80% up to 30°C  
(non-condensing)

II

2

The power plug may only be plugged into an outlet equipped with a grounding contact.

Do not interfere with the grounding function by using an extension cord without a ground wire. Any interruption of the ground wire inside or outside of the instrument, or release of the protective conductor terminal, can cause the instrument to become hazardous. Intentional ground interruption is not permitted!

Use only original power cables or alternative cables with a VDE / HAR logo, which at least fulfill the requirement of  $3 \times 0.75 \text{ mm}^2$  and 10A/250V.

**Caution!**

Through connection to the grounding connection, ancillary equipment with its own and/or extra power supply may be brought to the same ground wire potential. For connections without a ground connector, Leica Service must be consulted.

**Caution!**

Never use any fuses as replacements other than those of the types and the current ratings listed here. Bypassing fuse holders is not permitted.  
The use of incorrect fuses may result in a fire hazard.

The microscope's electrical accessory components are not protected against water. Water can cause electric shock.

**Caution!**

Protect the microscope from excessive temperature fluctuations. Such fluctuations can lead to the accumulation of condensation, which can damage the electrical and optical components. Operating temperature: 15-35°C.

### 3. Safety Notes

Before exchanging the fuses or lamps, be absolutely certain to switch off the main power switch and remove the power cable.

Per definition, the connection between power supply cable and appliance inlet is the power separator of this device. Thus, the user is responsible for sufficient clearance and accessibility of this area at all times.

#### Caution!

Do not use the microscope in altitudes exceeding 2000 m ASL/NL.

#### Caution!

Do not use this instrument near sources of high electromagnetic radiation (for example, unshielded, intentionally operated ultra-high frequency sources), because these can disrupt proper operation.

We recommend assessing the electromagnetic environment before operation of this instrument and then giving corresponding instructions.

### 3.3 Transport and storage

#### Caution!

Transport and storage in a range of -20°C – +85°C and at a humidity not exceeding 90% (non-condensing).

### 3.4 Notes on handling light sources

#### CAUTION

Light sources pose a potential irradiation risk (glare, UV radiation, IR radiation). Therefore, lamps and LED light sources have to be operated in closed housings and only after being mounted.

Never look directly into the beam path (blinding hazard).

Always connect the light guide to the microscope first to prevent danger to the user from the high-energy light emitted by the Leica EL6000 compact light source.

Never look into the light escaping from the light guide!

Lamps and lamp housings may be hot!  
**They must be placed at least 10 cm away from the wall and away from flammable substances.**

Particularly, supply and data lines must not encounter lamp housings!

### 3.5 Notes on handling immersion oil

#### Caution!

When using immersion oil, take care to avoid skin contact! Ask the supplier for a safety data sheet!

### 3.6 Notes on handling acids and bases

For examinations using acids or other aggressive chemicals, particular caution must be taken.

#### Caution!

Be absolutely certain to avoid direct contact with these chemicals.

### 3.7 Disposal

Once the product has reached the end of its service life, please contact Leica Service or Sales about disposal.

Please observe and comply with the national and federal laws and regulations that are equivalent to EC directives such as WEEE.

#### Note!

Like all electronic devices, this instrument, its accessory components and consumables must never be disposed of with general household waste.

# 4. Overview of the Instrument

Specification	Leica DM IL LED
<b>Contrast methods</b>	<ul style="list-style-type: none"><li>• Transmitted light: Bright field, phase contrast, IMC</li><li>• Incident light: Fluorescence</li></ul>
<b>Transmitted light axis</b>	Integrated LED Illumination manual adjustment of <ul style="list-style-type: none"><li>• Brightness</li><li>• Aperture diaphragm</li><li>• Automatic shutoff (adjustable)</li></ul>
<b>Incident light axis (optional)</b>	Incident light fluorescence illuminator for up to eyepiece field number 20 with <ul style="list-style-type: none"><li>• Interchangeable slider for holding 3 filter systems</li><li>• Light trap for suppressing ambient light</li><li>• Shutter, selectable</li></ul>
<b>Tube</b>	optionally with <ul style="list-style-type: none"><li>• Fixed or variable viewing angle</li><li>• Up to 3 switching positions</li><li>• With or without camera output</li><li>• Ergotube with height-adjustable eye level and camera port</li></ul>
<b>Objective nosepiece</b>	<ul style="list-style-type: none"><li>• Manual</li><li>• 4-fold for objectives with M25 thread</li></ul>
<b>Stages</b>	<ul style="list-style-type: none"><li>• Fixed stage</li><li>• 3-plate mechanical stage</li><li>• Heating stage</li></ul>

Specification	Leica DM IL LED
<b>Condenser</b>	optionally with <ul style="list-style-type: none"><li>• S80/0.30 condenser</li><li>• S40/0.45 condenser</li></ul>
<b>Focus</b>	<ul style="list-style-type: none"><li>• Focus dial for coarse and fine focusing</li><li>• Height adjustment</li></ul>

**Note:**


Leica DM IL LED, Leica DM IL LED Cellfactory, Leica DM IL LED Fluo and Leica DM IL M LED microscopes are designed for an integrated LED lamp housing or for the LH115 LED lamp housing. Connection of lamp housings, other than LED type, are not supported by microscope outlets.

## 4. Overview of the Instrument

### Important assembly groups

The following overall views illustrate and list important assembly groups of the microscopes and its accessory components.


**Fig. 1** Right side of the Leica DM IL LED stand

**Fig. 1/2**

- 1 Binocular tube
- 2 Eyepieces
- 3 Eyepiece tube
- 4 Tube holder
- 5 Stabilizing plate
- 6 Brightness control
- 7 Fluorescence filter blocks
- 8 Main power switch
- 9 Fluorescence lamp housing
- 10 Empty slide or modulation or phase contrast slider
- 11 Integrated LED lamp housing
- 12 Holder for filter, Ø 32 mm
- 13 Aperture diaphragm
- 14 Condenser
- 15 Transmitted light illumination column
- 16 Transmitted light illumination carrier
- 17 Stop lever for condenser height adjustment
- 18 Dark stop
- 19 Coarse and fine focusing
- 20 Objective nosepiece and objectives
- 21 Empty slide or IMC modulator
- 22 C-mount video adapter
- 23 Specimen stage

**Fig. 2** Left side of the Leica DM IL LED stand

**Fig. 3** Front view of the Leica DM IL LED

### Stand

The Leica DM IL LED stand offers high stability due to the low center of gravity. When using the multiple discussion attachment\* or for long exposures in photomicrography, a stand stabilizing plate\* is available to improve the stability.

For incident light fluorescence applications, an incident light axis is integrated in a second stand version.

### Tube holder

The tube holder is the interface between the microscope stand and the tube. The tube holder allows the use of the DM ILB and DM ILT tubes and the IL/L tube adapter, which allows tubes from the upright microscope series to be used. (See tube adapter).

### Tube

The tube contains a 1x tube lens which creates the primary image in conjunction with the objective.

The binocular tube consists of a body, the binocular part, and the tube change ring.

The trinocular tube offers an additional documentation output to accommodate photo or video equipment. A switchable mirror directs 100% of the light to either the eyepieces or the photo output.

The tubes for the Leica DM IL LED are interchangeable and revolving.

### IL/L tube adapter

The tube adapter is for holding the tubes from the upright microscope series, the multiple discussion attachment\*, the magnification changer\* and an ErgoModule\*.

### Eyepieces

The eyepiece creates a magnified, virtual image of the real intermediate image that is projected by the objective. In this process, the eyepiece works as a magnifying lens.

### Brightness controller

The stand has built-in electronics for adjusting the brightness via the brightness controller.

### Main power switch

The illuminated main power switch is used to switch the power supply of the microscope on and off. The illumination allows you to tell immediately whether the microscope is switched on, even in darkened rooms.

### Transmitted light illumination unit

The transmitted light illumination unit consists of the transmitted light illumination carrier and the transmitted light illumination column. The transmitted light illumination carrier includes a pre-centered, high-intensity LED illuminator, a fixture for a diaphragm slide, a fixture for a light filter, a condenser and an aperture diaphragm.

Fig. 4 Leica DM IL LED illumination system

## 4. Overview of the Instrument

### Lamp housing

The Leica DM IL LED stand has an integrated lamp housing with an LED illuminator.

### Filters

The filters are generally used to improve the contrast of the specimen. They are firmly mounted in a spoon-shaped holder ( $\varnothing$  32 mm).

Different filters can be inserted into the filter holder of the transmitted light illumination unit.

### Aperture diaphragm

The aperture diaphragm defines the resolution, focus depth and contrast of the microscopic image. The best resolution is obtained when the apertures of the objective and the condenser are roughly the same.

#### Note:

The aperture diaphragm in the **illumination light path** is **not** for setting the image brightness. Only the rotary brightness adjustment knob or the neutral density filter should be used for this purpose.

### Condenser

The condenser is a lens system that collects the light and directs it to the specimen from the top. The condenser serves to utilize the numerical aperture in the objective.

### Stop lever for condenser height adjustment

The stop lever is used for condenser height adjustment by moving the transmitted light illumination carrier. The markings on the transmitted light column indicate the proper height to be set for the condenser used.

### Specimen stages and accessories

The specimen stage is used to support the specimens to be subjected to microscopic examination. For microscopy with the different specimens, multiple options are available, such as the attachable mechanical stage, retaining clips, 3-plate mechanical stage, heating stage, etc.

### Objective nosepiece and objectives

The objective nosepiece is designed to hold the objectives. Especially the L-type objectives with long working distance are able to compensate for the different thicknesses of flask bottoms.

All microscope objectives with magnifications 2.5 to 100 can be used. All objectives with a 25 mm thread are compatible. For additional information about objectives, refer to the "Technical Description" chapter or the respective objective lists.

(<http://www.leica-microsystems.com/products/light-microscopes/accessories/objectives/>)

### Coarse and fine focus adjustment

The coarse and fine focus adjustment enables you to adjust the microscopic image quickly and accurately. The focus is adjusted by moving the objective nosepiece vertically. The stroke length is 7 mm.

### **Incident light fluorescence device\***

The microscope version with incident light fluorescence device contains the integrated fluorescence axis and the lamp mount for mounting of a lamp housing.

### **Fluorescence filter block slide\***

The fluorescence filter block slide can hold up to three fluorescence filter blocks. The filter block slide can now be moved back and forth between the three switch positions.

One position of the slide can also be used as a bright-field position by leaving it empty (no filter block inserted).

### **Lamps and lamp housings\* for the incident light fluorescence device**

An additional illuminator is required for incident light fluorescence. Lamp housings of the 106z series can be used on the Leica DM IL LED series. Depending on the version used, the controls for the lamp housings are located on the right or left side.

The LH 106z lamp housing has a reflector which can be centered and focused. It also contains a multilens collector.

The following lamps can be used with their special sockets:

High-pressure mercury burner 100 W

For Leica DM IL M LED: The LH115 LED lamp housing is provided.

### **Supply unit\***

An external supply unit is used to regulate the lamp for incident light fluorescence and the corresponding lamp housings.

### **Modulation slider or phase contrast slider\*\***

The modulation contrast slider or phase contrast slider is part of a contrast method, either Integrated Modulation Contrast (IMC) or phase contrast.

For the S40/0.45 and S80/0.30 condensers, different sliders are used both for the phase contrast and for the IMC. In each case, three phase rings or three switch positions for IMC are integrated on the sliders.

If no phase slide or modulation slide is used, it is also possible to insert an empty slide into the corresponding holder on the condenser.

### **IMC modulator\*\***

For the Leica Integrated Modulation Contrast (IMC), the IMC modulator is offered in the stand.

All Leica DM IL LED stands come standard equipped with an empty slide.

\* Fluorescence stands

\*\* Integrated modulation contrast (optional)

## 4. Overview of the Instrument

**Fig. 5a** Rear view of the Leica DM IL LED microscope (Leica DM IL LED Cellfactory equivalent)

- 1 Connection for transmitted light illumination carrier with integrated LED illumination
- 2 Type label
- 3 Power supply connection
- 4 Earth (ground) terminal
- 5 Earth (ground) terminal symbol

**Fig. 5b** Rear view of the Leica DM IL LED Fluor microscope

- 1 Connection for transmitted light illumination carrier with integrated LED illumination
- 2 Type label
- 3 Power supply connection
- 4 Earth (ground) terminal
- 5 Earth (ground) terminal symbol
- 6 Lamp mount (for incident light fluorescence)

**Fig. 5c** Rear view of the Leica DM IL M LED microscope

- 1 Connection for transmitted light illumination carrier with integrated LED illumination
- 2 Type label
- 3 Power supply connection
- 4 Earth (ground) terminal
- 5 Earth (ground) terminal symbol
- 6 Lamp mount (for incident light fluorescence)

# 5. Unpacking

First, carefully remove all components from the transportation and packaging materials.

## Note:

If at all possible, avoid touching the lens surfaces of the objectives. If fingerprints do appear on the glass surfaces, remove them with a soft leather or linen cloth. Even small traces of finger perspiration can damage the surfaces in a short time. See the chapter on "Care of the Microscope" → p. 61 for additional instructions.

## Caution!

Do not connect the microscope and peripherals to the power supply at this point under any circumstances!

The following parts can be included in the delivery:

- Leica DM IL LED stand, including fixed stage
- Illumination and condenser holder
- Tube
- Eyepieces
- Objectives
- Condenser
- Dust cover
- Operating manual
- Optional components:
  - IL/L tube adapter
  - Phase slide
  - Focusing telescope
  - IMC slide
  - IMC module
  - Filters for transmitted light
  - Filter slide for fluorescence blocks
  - Fluorescence blocks
  - Lamp housing
  - High-pressure mercury burner
  - External supply unit
  - C-mount video adapter
  - Camera
  - Specimen stage accessories
  - Additional components from the upright microscope product line such as tubes, multiple discussion attachment, magnification changer, ErgoModule

**Fig. 6** Leica DM IL LED with discussion attachment

## 5. Unpacking

### Installation location

Work with the microscope should be carried in a dust-free room that is free of vapors (oil, chemicals etc.) and extreme humidity. At the workstation, large temperature fluctuations, direct sunlight and vibrations should be avoided. These conditions can distort measurements and micrographic images.

Permitted ambient conditions

Temperature	15–35°C
Relative humidity	maximum 80% up to 30°C non-condensing

Microscopes in warm and warm-damp climatic zones require special care in order to prevent fungus contamination.

See the chapter on "Care of the Microscope" → p. 61 for additional instructions.

#### Caution!

When installing the microscope, make sure the power inlet is freely accessible so that the instrument can be quickly disconnected from the mains if necessary.

#### Caution!

Electrical components must be assembled at least 10 cm from the wall and away from flammable substances.

### Transport

#### Caution:

For shipping or transporting the microscope and its accessory components, the original packaging must be used.

As a precaution to prevent damage from vibrations, the following components should be disassembled and packaged separately:

- Unscrew the objectives.
- Remove the condenser.
- Remove the lamp housings.
- Disassemble the burner of the 106z lamp housing.
- Remove all moving or loose parts.

### Storage

Put a dust cover over your microscope after use to protect it from dust.

### Weight

The weight of the microscope depends on the particular equipment.


Fully equipped, the microscope weighs about 13 kg. For transportation, the user has to take care of the corresponding actions.

#### Caution!


For transporting it is essential to remove all components listed under "Transport"!

Dimensions (specifications in mm)

+) for microscopes with 240 mm illumination stand, the value is correspondingly higher for microscopes with 480 mm illumination stand

Weight approx. 13 kg

# 6. Assembling the Microscope

## 6.1 Stand

- Place the base stand of the Leica DM IL LED on a sufficiently clear work table.
- Make certain that all four stand bases are pre-assembled and in place on the bottom on the stand.

### Caution!

Do not connect the stand to the power supply at this point under any circumstances!

- The microscope stage and transmitted light illumination axis are usually factory-assembled as shown in Fig. 7. If one has been delivered

Fig. 7 Stand with transmitted light illumination column

- 1 Screw for condenser collision protection
- 2 Stand bases

separately or another stage (3-plate mechanical stage or a heating stage) is retrofitted, fasten these stages using the three screw-on points. The stage is fastened using the guide pins and screwed on loosely at first. Then, all three screws are tightened once again.


### Caution!

The screws below the stage on the transmitted light illumination carrier must **not** be unscrewed. If these screws are unscrewed, the optical axis is moved.

- If your scope of delivery includes a stabilizing plate, screw it onto the bottom on the stand so that the two front stand feet fit into the notches (Fig. 8).

Tighten the screws and place the stand back on the table in the upright position.

Fig. 8 Stabilizing plate

### 6.2 Attaching the condensers

- Screw the S80/0.30 (9.2) or S40/0.45 (9.1) condenser from below into the condenser holder (10.2) of the transmitted light illumination carrier.

**Note:**

A condenser may not be required for the Leica DM IL LED Cellfactory version.

- Make sure that the transmitted light illumination carrier is locked in place.

**Note:**

For the stands with 3-plate mechanical stages, position the illumination carrier 25 mm below the marks, as the illumination axis is installed 25 mm higher on an adapter than on the fixed stage.

**Note:**

A screw (7.1) on the transmitted light illumination column prevents the condenser from colliding with the specimen stage.

If necessary, move the screw into the higher of the two openings, if using the S80/0.30 condenser.

### 6.3 Inserting the transmitted light illumination carrier

- Insert the transmitted light illumination carrier from above into the column by pressing and holding down the stop lever for condenser height adjustment (10.5).
- Position the transmitted light illumination carrier (10.3) on the transmitted light illumination column (10.4) depending on the condenser used (S40/0.45 or S80/0.30) and release the stop lever. (See note below.)


The marks (10.6) pertain to a liquid height of 15 mm. For stands with double marks, the upper and lower mark represent a range with different liquid heights.

**Fig. 10** Transmitted light illumination unit with condenser

- 1 Condenser
- 2 Condenser holder
- 3 Transmitted light illumination carrier
- 4 Transmitted light illumination column
- 5 Stop lever for condenser height adjustment
- 6 Marks

**Fig. 9** Condensers

- 1 S40/0.45 condenser
- 2 S80/0.30 condenser

## 6. Assembling the Microscope


### Changing the position of the transmitted light illumination carrier with the specimen stage

For the Leica DM IL LED Cellfactory version, the specimen stage must be rotated by 180° together with the transmitted light illumination carrier.

- Slacken the screws (11a.1) with a 3 mm Allen key and remove them.
- Rotate the stage with the transmitted light illumination unit through 180°.
- Insert the stage with the transmitted light illumination unit. Make sure the stage engages in the guide pins (11b.1) on the stand.
- Insert the screws (11c.1) and screw tight.
- Cable management clips (11c.2) have been supplied for laying the connecting cable. Stick these onto the underneath of the stage.


**Fig. 11a**

- 1** Screws for changing the position of the transmitted light illumination carrier

**Fig. 11c**

- 1** Screws for changing the position of the transmitted light illumination carrier  
**2** Cable management clips

**Fig. 11b**

- 1** Guide pins

### 6.4 Electrical connection of the transmitted light illumination carrier

- Use the connecting cable to connect the transmitted illuminator to the connection socket (12.1) on the rear panel of the instrument. The downward-bent plug faces down. Securely tighten the socket union.

**Note:**

The connecting cable for the Leica DM IL LED Cellfactory version is longer than usual. Fix the connecting cable as shown in Fig. 11c with the cable clips provided.

### 6.5 Inserting the tubes


The microscope comes standard with the DM ILB (13.2 binocular tube) or DM ILT (13.1 trinocular phototube). For additional information, refer to the "Technical Description" chapter.

#### DM ILB and DM ILT tube

- Unscrew the clamping screw (14.1) using a 3 mm hexagon socket wrench.
- Insert the tube (13.3) into the tube holder (14.2).
- Retighten the clamping screw.

**Fig. 13** Tubes

- 1 DM ILT trinocular phototube
- 2 DM ILB binocular tube


**Fig. 12** Rear view of the Leica DM IL LED  
Leica DM IL LED Fluo and DM IL M LED equivalent

- 1 Connection socket for connecting cable

**Fig. 14**

- 1 Clamping screw
- 2 Tube holder
- 3 Tube

## 6. Assembling the Microscope

- To set a new observation position, loosen the clamping screw (14.1), rotate the tube as desired and retighten the clamping screw.

### Tubes from the Leica DM1000-3000 product line

Instead of the standard DM IL tubes, you can also adapt tubes from the Leica DM1000-3000 product line.

Proceed as follows:

- Unscrew the clamping screw on the tube changer (14.1) using a 3 mm hexagon socket wrench.
- First, insert the IL/L tube adapter (Fig. 15) into the tube holder of the stand.
- Retighten the clamping screw.
- Unscrew the clamping screw on the IL/L tube adapter (15.1).
- Place a tube into the tube holder of the tube adapter.
- Retighten the clamping screw.
- To set a new observation position, loosen the clamping screw (15.1), rotate the tube as desired and retighten the clamping screw.

### Note:

To insert a multiple discussion tube\*, a magnification changer\* or an ErgoModule\*, refer to the chapter on "Assembling the Options".


**Fig. 15:** IL/L tube adapter  
1 Clamping screw

**Fig. 16** HC L1T tube

**Fig. 17** HC L1VT tube

## 6.6 Eyepieces and graticules\*

### Inserting the eyepieces

The eyepieces are inserted into the eyepiece tubes.

The following eyepieces are offered:

HC PLAN 10x/20 (M)

HC PLAN 12.5x/16 M

Widefield 16x/14 Br (M)<sup>+)</sup>

Widefield 25x/9.5 Br (M)<sup>+)</sup>

<sup>+)</sup>  also requires a spacer ring

Br) eyepiece for spectacle wearers

For information on the diameter, the viewable specimen area and the total magnification of the microscope, refer to the "Technical Description" chapter.

### Positioning the graticules\*

Graticules can be retrofitted by the user for the HC PLAN eyepieces listed above.

Basically, graticules can be used only with eyepieces with adjustable eyelens = Type M.

### ! Important:

Be extremely careful to maintain cleanliness. Otherwise, dust particles and fingerprints will appear in the field of view.

HC PLAN eyepieces have a uniform graticule diameter of 26 mm.

HC PLAN 10x/20 M and HC PLAN 12.5x/16 M eyepieces:

- Unscrew the retainer sleeve from the lower part of the eyepiece.
- Insert the graticule so that the coated side faces downward (toward the objective) and any marking appears laterally correct when viewed in the later observation direction.
- Screw the retainer sleeve back in.

## 6.7 Objectives

- Remove the screw caps on the objective threads.
- Screw the objectives into the nosepiece opening so that incremental change of the magnification levels is possible (e.g. in the sequence 4, 10, 20, 40).
- If objective threads are left unused, cover these with screw caps to protect the microscope optics against dust.

Fig. 18 Eyepiece pairs

- 1 HC PLAN 10x/20 ⚙
- 2 HC PLAN 10x/20 ⚙ M

Fig. 19 HI PLAN objective set

## 6. Assembling the Microscope

### 6.8 Inserting the filter

- Insert the filter (Fig. 20) into the filter holder (21.1) on the transmitted light illumination carrier.

**Fig. 20** Filters

**Fig. 21**

- 1 Filter receptacle

### 6.9 Specimen stage


#### Attaching the attachable mechanical stage

- Attach the attachable mechanical stage to the right side of the stage to accommodate brackets for different culture flasks (Fig. 23).
- Mount the attachable mechanical stage with a 3 mm hexagon socket wrench.

Some holders come with self-adhesive scale for reading the coordinate adjustment.

- Attach these in the recesses on the attachable mechanical stage.

**Fig. 23** Attachable mechanical stage with fixture for holders

6.10 Light source for the transmitted light axis

**Note:**

The Leica DM IL LED is equipped with integrated LED illumination. The service life of the LED is about 50000 hours. If, despite this, it should be necessary to change the LED, this task must be carried out by Technical Service only.

6.11 Electrical connection of the microscope

**Caution!**

Do not connect the microscope and supply unit\* to the power supply until all options are installed.

The power plug may only be plugged into an outlet equipped with a grounding contact. Do not interfere with the grounding function by using an extension cord without a ground wire.

- If you purchased additional options with the microscope, install these options first (see next chapter).

**Caution! Caution!**

For external lamp supply units, always carry out the power supply voltage setting as described in the manual which is supplied separately, or use a power transformer.

- Plug the power plug into the rear side of the microscope (24.1) and connect it to the wall socket.

**Caution!**

Make sure to follow the safety notes on pages 10-13!

**Fig. 24** Rear view of the Leica DM IL LED  
Leica DM IL LED Fluo and DM IL M LED equivalent  
**1** Power supply connection

# 7. Assembling the Options

**Note:**

These assembly tasks are omitted if no additional accessory components were purchased with the microscope.

## 7.1 Assembling the fluorescence filter blocks\*

**Note:**

For microscope with integrated incident light fluorescence device only.

The filter block slide (Fig. 25) can hold up to three fluorescence filter blocks.

- For assembly of the filter blocks, remove the slide cover (25.3).
- Insert the filter blocks (25.1) with the engraved side down and writing upside-down into the dovetail receptacle (25.2). Ensure that the filter blocks engage.
- Reinstall the slide cover.
- Check that the slide cover (25.3) is properly installed.
- Mark the filter positions using the enclosed adhesive labels.

### Inserting the filter block slide\*

- Hold the filter block slide so that the warning sticker is at the front left and insert it into the dovetail receptacle on the left side of the stand.
- The filter block slide can now be moved back and forth between the three switch positions.

**! Caution!**


The filter block slide is not protected against inadvertently sliding out.

When using the incident light beam splitter or polarization beam splitter in combination with fluorescence filter blocks, there is a blinding hazard if you switch over accidentally!

**Fig. 25**

- 1 Filter block
- 2 Dovetail receptacle
- 3 Slide cover
- 4 Lower part of the filter block slide

## 7.2 Inserting the phase contrast slider on the transmitted light illumination carrier\*

Depending on whether the S40/0.45 or S80/0.30 condenser is used, various light ring slides exist.

### ! Caution!

The light ring slide is not protected against inadvertently sliding out.

- Remove the dummy slide, where applicable (27.1).
- Hold the light ring slide (26.1, 26.2) so that the marking "TOP LEFT" (26.6) is at the top left and the marking of light rings 5, 10/20 and 40 is facing you. The click stops (26.3) are on the front longitudinal side of the slide.
- Insert the light ring slide laterally into the transmitted light illumination carrier (Fig. 27). The keyways should click in place when the slide is inserted.

Fig. 26

- 1 Slider for light rings (S80/0.30)
- 2 Slider for light rings (S40/0.45)
- 3 Click stop
- 4 Light rings
- 5 Bright-field position
- 6 "TOP LEFT" marking


Fig. 27 Transmitted light illumination carrier with holder for light ring slide and IMC slit diaphragm slider

- 1 Dummy slide

## 7. Assembling the Options

### 7.3 Inserting the IMC slit diaphragm slider\* on the transmitted light illumination carrier

Depending on whether the S40/0.45 or S80/0.30 condenser is used, various IMC slit diaphragm sliders exist.


**! Caution!**

The IMC slide is not protected against inadvertently sliding out.

- Remove the dummy slide, where applicable (27.1).
- Hold the IMC slide so that the marking "TOP LEFT" is at the top left (28.1). The click stops (28.2) are on the front longitudinal side of the slide.
- Insert the IMC slide laterally into the transmitted light illumination carrier (Fig. 27). You should feel the keyways click into place when the slide is inserted.

**Fig. 28**

- 1 "TOP LEFT" marking on the diaphragm slide
- 2 Click stops
- 3 IMC slit diaphragm slider for S80/0.30
- 4 IMC slit diaphragm slider for S40/0.45

### 7.4 Inserting the IMC modulator\*


- Remove the dummy slide, where applicable (Fig. 29, 31).
- Insert the IMC modulator so that the marking (30.3) faces forwards.
- Lock the slide in position BF (marking BF is visible) or in position IMC (marking IMC is visible).

**Fig. 29** Modulator dummy slide


**Fig. 30**

- 1 Click stops
- 2 IMC modulator
- 3 Marking of the IMC modulator

**Fig. 31**

- 1 Holder for IMC modulator or dummy slide

**7.5 Assembling the 106z lamp housing\* with mercury lamps**

**Note:**

For microscope with integrated incident light fluorescence device only.

**Note:**

If a Leica SFL100 is used for fluorescence illumination, please regard the relevant operating manual.

The 106z lamp housing is used with a variety of gas discharge lamps.

Light sources pose a potential irradiation risk (glare, UV radiation, IR radiation). Therefore, lamps have to be operated in closed housings and only after being mounted.

Ensure that the lamp housing has been disconnected from the power supply. Unplug the power plug and the power supply during assembly.

Never touch the glass parts of the burner with bare hands.

Never look directly into the beam path (blinding hazard).

**Caution!**

Make sure to follow the instructions and safety notes of the lamp supplier. Before changing lamps, allow at least 30 minutes for them to cool down!

The lamp may still be hot.

**Caution!**

First, read the user manual of the power supply (supply units)!

## 7. Assembling the Options

### Inserting the gas discharge lamps\* into the 106z lamp housing

Mercury lamps are powered by separate supply units.

Please also read the separate instruction manual provided with these supply units.


The following gas discharge lamps may be used and require different supply units and lamp mounts (Fig. 42):

Typ	Typical bulb life <sup>+) </sup>
100W high-pressure mercury burner (direct current)	200 hrs.
100W high-pressure mercury burner type 103 W/2 (direct current)	300 hrs.

<sup>+)</sup>  Please observe the data sheets of the lamp manufacturer.


**Fig. 37** Lamp housing 106z, open

- 1 Cover, in raised position
- 2 Collector
- 3 12 V 100 W halogen bulb or gas discharge lamp (see Fig. 42)
- 4, 9 Cover fastener
- 5 Reflector
- 6, 8 Adjusting screws for x and y centering of the lamp
- 7 Focus of the reflector
- 10 Fastening screws for lamp mount
- 11 Socket for disconnect plug

**Fig. 42** Lamp mount for gas discharge lamps

- 1 Upper clamping system
- 2 Lower clamping system
- 3 Cooling element

- To open the 106z lamp housing, unscrew the fastening screws (43.1) on the cover.
- Remove the transport anchorage (red plastic rod in place of the burner) in the lamp mount. To do so, remove the lower clamp (42.1). Pull up the cooling element (42.7) and turn it to the side. Detach the lower clamping system (42.3) and remove the transport anchorage.
- Install the burner in reverse order.

### Caution!

After installation, the marking must be upright.  
If a glass melt nipple is present, position it by turning the burner so that the nipple does not impede the beam path later, but instead is positioned sideways.

- Reinsert the lamp mount and retighten the fastening screw (37.10).
- Close the lamp housing and retighten the fastening screws.
- Place the lamp housing in the incident light lamp housing mount (5b.5 or 5c.5, p. 19) and fasten it with the clamping screw on the side.


**Fig. 43** 106 z lamp housing

- 1 Fastening screws for the cover

**Fig. 44** 106z L lamp housing with Hg 100 W lamp

- 1 Focus of the collector  
2 Vertical lamp adjustment  
3 Horizontal lamp adjustment  
4 Mercury lamp mount

## 7. Assembling the Options

### Note:

Hg 100 burner:

The burner holder of the Hg 100 lamp housing is higher than the free room below the lamp housing. Therefore, the microscope must be positioned 20 mm higher. For this purpose, a height compensation plate (Fig. 45) is offered.

The plate can also be used generally to increase the height for ergonomic reasons.

- Connect the lamp housing to the supply unit.

### Caution!

The burner must be adjusted immediately after being lit.

### Note:

An incorrect installation will reduce the lamp brightness to about 60% and will considerably limit the useful life of the lamp.

### Caution!

Be absolutely certain that the marks of the bulb socket and supply unit match. If, for example, the bulb socket is at L1 (or L2), the supply unit must also be at this setting in order to utilize the lamp fully and not shorten its service life.

### 7.6 Assembling the lamp housing LH115 LED

For the Leica DM IL M LED the lamp housing LH115 LED (Fig 48a) is provided.

- Place the lamp housing in the incident light lamp housing receptacle and fasten it with the clamping screw on the side.
- Connect the lamp housing to the socket (48b.1).


### Note:

A separate lamp power supply is necessary to connect other lamp housings.

**Fig. 46** Rear side of the ebq 100 supply unit

#### 1 Lamp connection

**Fig. 45** Height compensation plate

7.7 Leica EL6000\*

When using the Leica EL6000 compact light source, it is essential to observe the safety notes in the separate instructions.

Switch on light sources only if they are firmly connected to the microscope.  
Uncontrolled light output poses a blinding hazard!

For the Leica DM IL LED, you need an additional adapter to connect the Leica EL6000 compact light source (Fig. 49a).

Fig. 48b Rear side of DM IL M LED

- 1 Connection for LH115 LED
- 2 Lamp mount (for incident light)

Fig. 48a Lamp housing LH115 LED for DM IL M LED  
Lifetime: up to 50.000 hrs..


Fig. 49a Adapter for connecting the compact light source Leica EL6000

Fig. 49b Leica EL6000 compact light source

## 7. Assembling the Options

### 7.8 Adapting cameras to binocular photo tubes\*

You can adapt an image acquisition system, such as a video camera or digital camera, to each of the binocular photo tubes provided for the Leica DM IL LED.

#### Adapting

Various adapters are available for adapting cameras with a C-mount and B-mount objective mount. The adapters listed in the table below can be used on all trinocular photo tubes, although some tubes require an additional photo adapter tube. The image detail on the monitor depends on the adapter being used and the chip size of the camera.

### Calculating the magnification on the monitor

The magnification  $V_{TV}$  on the monitor can be calculated according to the following formula or measured using a stage micrometer and centimeter scale.

$$V_{TV} = \frac{\text{Objective magnification} \times \text{Factor magnification changer}^* \times \text{Magnification of adapter}^* \times \text{Screen diameter}}{\text{Chip diameter of the camera}}$$

	Recorded picture diagonal in mm with			
	1 inch camera	2/3 inch camera	1/2 inch camera	1/3 inch camera

#### Without zoom magnification, for 1-chip cameras only:

C-mount adapter 1 x HC	16	11	8	6
C-mount adapter 0.70 x HC	-	15.5	11.4	7.8
C-mount adapter 0.55 x HC	-	-	14.5	10.9
C-mount adapter 0.35 x HC	-	-	-	17.1

#### With zoom magnification (Vario TV adapter) for 1-3 chip cameras:

C-mount, 0.32-1.6 x HC	-	-	19 <sup>+)5</sup>	18-3.8
B-mount (ENG), 0.5-2.4 x HC (1/2 inch)	-	-	16-3.3	-

<sup>+) from zoom factor 0.42 x only!</sup>

#### Without zoom magnification, for 1-3 chip cameras:

C-mount adapter 1 x	-	-	16	12
B-mount adapter 1 x	-	-	16	12
B-mount adapter 1.25 x	-	17.5	-	-
F-mount adapter 1 x	-	-	16	12
F-mount adapter 1.25 x	-	17.5	-	-

**Additionally required for each:** TV optics 0.5 x HC

In principle, the same process is used for inserting the multiple discussion attachment, the magnification changer and the ErgoModule. The ErgoModule can either be assembled directly on the base stand (in combination with the DM ILB or DM ILT tubes) or on the DM IL/L tube adapter (in combination with the L tubes).

### 7.9 Inserting the multiple discussion attachment\*

#### Note:


When using the multiple discussion attachment, the stand stabilizing plate should be used (→ p. 23).

- Unscrew the clamping screw (50b.2) on the stand using a 3 mm hexagon socket wrench.
- Attach the IL/L tube adapter (50b.3).
- Retighten the clamping screw (50b.2).

- Unscrew the clamping screw (50b.4) on the tube adapter.
- Insert the multiple discussion attachment (50b.1) into the tube holder of the adapter.
- Tighten the clamping screw (50b.4).
- Detach the top part of the holding foot (50a.1). This part is not needed for the Leica DM IL LED.
- Unscrew the locking screw (50a.2).
- Set the height using the black knurled ring (50a.3).
- Fix the height using the screw.
- Unscrew the clamping screws (50b.5) on the multiple discussion attachment.
- Attach the tubes.
- Tighten the clamping screws (50b.5) on the multiple discussion attachment.

**Fig. 50a** Multiple discussion attachment holding foot

- 1 Upper part
- 2 Locking screw
- 3 Knurled ring for height adjustment

**Fig. 50b** Stand with multiple discussion attachment

- 1 Multiple discussion attachment
- 2 Stand clamping screw
- 3 IL/L tube adapter
- 4 Tube adapter clamping screw
- 5 Multiple discussion attachment clamping screws

## 7. Assembling the Options

### 7.10 Inserting the ErgoModule\*

- Unscrew the clamping screw (50b.2) on the stand using a 3 mm hexagon socket wrench.
- When using the DM IL/L tube adapter:
  - Attach the IL/L tube adapter (Fig 51).
  - Tighten the clamping screw (51.1).
  - Unscrew the clamping screw (50b.4) on the tube adapter.
- Insert the ErgoModule into the tube holder of the base stand or the IL/L tube adapter.
- Tighten the clamping screw (51.1).
- Unscrew the clamping screw on the ErgoModule.
- Attach the tube.
- Tighten the clamping screw on the ErgoModule.

**Fig. 51:** IL/L tube adapter

**1** Clamping screw

# 8. Operation

## Caution!

For examinations using acids or other aggressive chemicals, particular caution must be taken. Avoid direct contact of these substances with optical and mechanical parts.

## Caution!

The lamp housing must be mounted before operating the microscope!

## Caution!

After turning on the gas discharge lamp\*, the burner must be immediately adjusted. Therefore, **do not** turn on the supply unit\* yet. First, work in transmitted light in order to familiarize yourself with the microscope's controls.

Switch on light sources only if they are firmly connected to the microscope. Uncontrolled light output poses a blinding hazard!

## 8.1 Basic settings for transmitted light

### Switching on the LED transmitted light illumination

- Switch on the power switch (52.2).
- Adjust the brightness using the knob (52.1).

## Note:

By rotating the brightness controller once (up to three seconds after switching on the Leica DM IL LED) from the zero position to maximum brightness and back, you can activate or deactivate the automatic shutoff of the LED. This means that after two hours, the LED shuts off automatically and is re-enabled only by turning the brightness controller.

Fig. 52

- 1 Brightness control
- 2 Main power switch
- 3 Focus drive

## 8. Operation

The LED flashes briefly three times to indicate activation. The LED flashes twice for a longer time to indicate deactivation.

### Adjustment specimen

For initial setting of the microscope, we recommend using a specimen that has both high-contrast and low-contrast areas.

For incident light fluorescence of transparent specimens, we recommend carrying out an adjustment in transmitted light first.

### Focusing the specimen

- Position the desired objective. When doing this, the objective nosepiece should be lowered first. The objective is swung into place by turning the black knurled ring of the nosepiece. Be sure that each nosepiece audibly clicks into place.
- Bring the image into focus using the coarse and fine adjustment, which adjusts the height of the objective nosepiece. The level of the stage remains unchanged. The total adjustment is 7 mm. The focusing range extends (in air) from 1.0 mm below the stage surface to 6 mm above it.

### Caution!

Depending on the objective used, the objective nosepiece **must** be lowered before turning over the objective. Otherwise, there is a risk that the objective will collide with the stage.

### Adjusting the tubes and eyepieces

For microscopy work with glasses, the protection screen on the eyepieces must be removed or inverted, but it should always be used when observing without glasses.

- Adjust the interpupillary distance on the tube by pulling the eyepiece tube or pushing them together so that when viewing with both eyes, a congruent total image is seen and not a double image.
- Note your personal interpupillary distance.
- Additionally, for ergo tubes:  
Adjust the viewing angle ( $0^{\circ}$  –  $35^{\circ}$ ) by tilting the binocular eyepiece. To prevent symptoms of fatigue, vary the viewing angle from time to time if necessary.
- Close any unused tube openings, as otherwise stray light can interfere with observation.

Fig. 53

Always an unobstructed view of the specimen, even when using the trinocular phototube

**DM ILB binocular tube**

For eyepieces with inserted graticule only\*:

- Bring the specimen very much out of focus or remove it from the beam path.
- Focus the graticule by adjusting the eyelens with a relaxed eye. (The best way to relax your eye is to look for a moment at a distant object outside of the frame.)
- Focus the specimen only through the eyepiece using the graticule.
- Then, close that eye and focus on the object by adjusting the second eyepiece only.

Only if no graticule is inserted in both eyepieces:

- Bring the specimen very much out of focus or remove it from the beam path.
- Adjust the eyelens so that the limit of the field of view appears sharp. When you adjust the eyelens, a bright line is visible along the outer perimeter of the eyepiece body. It indicates the correct position of the eye lens for normal eyesight and for users wearing glasses during microscopy work with corrective glasses.

**Note:**

Eyeglasses with multifocal lenses (bifocals and smooth view glasses) must be removed while operating the microscope.

- Focus on the object through the eyepieces.

Only if one eyepiece does not have an adjustable eyelens:


- First focus the specimen through this eyepiece exactly (close the other eye).
- Then, likewise focus the image by adjusting the eyelens of the second eyepiece.

**Correction for vision problems**

- With your right eye, look through the right eyepiece and bring the specimen into sharp focus using the fine adjustment.
- Then, with your left eye, view the same position of the specimen and rotate the left eyepiece tube until this position is brought into sharp focus. While doing so, do not use the fine adjustment.
- If using eyepieces with adjustable eyelenses, do not compensate for vision problems by adjusting an eyepiece tube, but by adjusting the eyelens of the eyepiece.

**Fig. 54** Tubes

- 1 DM ILT trinocular phototube
- 2 DM ILB binocular tube

## 8. Operation

### DM ILT trinocular phototube

- Adjust the beam splitter to visual observation by moving the push rod. The meaning of the switch positions is illustrated by symbols on the side surface of the tube.  
Pulled-out push rod = photo  
Inserted push rod = visual
- Adjusting the eyepieces is identical to adjusting the binocular tube.
- Compensate for any vision problems by adjusting the eyelens of the eyepiece.

### 8.2 Objectives

#### Immersion objectives

**OIL:** Use optical immersion oil according to DIN/ISO standards only.

Cleaning → p. 62, marking → p. 63 ff.

#### Caution!

Follow safety notes for immersion oil!

**W:** Water immersion. The special water immersion objectives with ceramic front area can be used for all aqueous solutions.

**IMM:** Universal objective for water, glycerin, oil.

**Color-coding** of the objectives → p. 64.

#### Locking objectives

For some immersion objectives (with knurled wheel), the objective can be "shortened" (locked). Then, drops of immersion oil that have not been wiped off do not accidentally wet objectives and other specimens when the objective nosepiece is turned.

- Push the front part approx. 2 mm towards the nosepiece.
- Lock the objective in place with a small rotational movement.

#### Caution:

When you use the immersion objective again, be absolutely certain to release the lock, as otherwise the spring action to protect the specimen and objective is disabled and, furthermore, the other objectives are no longer parfocal to the immersion objective.

#### CORR objectives

These are special objectives that can be adjusted to the thickness of the cover slip.

- By turning the knurled ring, set the corrective mount roughly to a middle or estimated value.
- Optimize the focus on the specimen.
- Adjust the corrective mount until optimum contrast appears; refocus using fine adjustment if necessary.

### 8.3 Transmitted light

#### Bright-field illumination

Illumination methods in which the empty areas of the specimen make up the brightest areas of the image are called bright field. Bright-field observation requires absorbent specimen structures, i. e. most specimens will need specimen staining. Alternatives are optical contrast methods such as phase contrast or modulation contrast.

#### Adjusting the condenser

For correct height adjustment of the S80/0.30 and S40/0.45 condensers, the stand has marks (55.1). These marks pertain to a liquid height of 15 mm. For stands with double marks, the upper and lower mark represent a range with different liquid heights. Move the condenser holder while observing the specimen to attain an optimal image.

#### Note:

For the stands with 3-plate mechanical stages, position the illumination carrier 25 mm below the marks, as the illumination axis is installed 25 mm higher on an adapter than on the fixed stage.

- Press the stop lever (55.2) and adjust the transmitted light illumination carrier until the upper edge of the carrier is even with the corresponding condenser height mark.

#### Adjusting the aperture diaphragm

The aperture diaphragm (55.3) determines the resolution, the focus depth and the contrast of the microscope image. The best resolution is obtained when the apertures of the objective and the condenser are roughly the same.

When the aperture diaphragm is stopped down to be smaller than the objective aperture, resolution is reduced, but the contrast is enhanced. A noticeable reduction in the resolution is observed when the aperture diaphragm is stopped down to less than 0.6x of the objective aperture and should be avoided where possible.

- Adjust the aperture diaphragm subjectively, according to your opinion of the image.
- In principle, you can carry out a calibration yourself by comparing with the apertures of different objectives.
- You can visually compare the apertures of objective and condenser as follows:
  - Remove the eyepiece from the eyepiece tube or use a focusing telescope and focus.
  - Open or close the aperture diaphragm just enough for your image to be visible in the eyepoint (= illuminated circle) of the objective. This position is considered the normal position, i.e. condenser aperture = objective aperture.
  - Reinsert the eyepiece.

## 8. Operation

For specimens with lower contrast, you can close the aperture diaphragm further to make structure elements with lower contrast more clearly visible. In polarization microscopy, stopping down the aperture diaphragm generally results in more intense colors.

### Note:

The aperture diaphragm in the **illumination light path** is **not** for setting the image brightness. Only the rotary brightness adjustment knob or the neutral density filter should be used for this.

An aperture diaphragm in the **objective** is normally fully opened. The reduction in image brightness caused by stopping down results in:

- Greater focus depth
- Less cover slip sensitivity
- Dark field impression
- Change in contrast

### Possible faults

Incorrect thickness of the cover slip or incorrect objective. Specimen positioned with cover slip facing upwards rather than downwards.

Aperture diaphragm opened or closed too far.

Condenser in incorrect height position.

The fluorescence filters are in the beam path.


Light ring switched on by mistake.

IMC component switched on by mistake.

Optics dirty.

**Fig. 55** Transmitted light illumination unit with condenser

- 1 Marks
- 2 Stop lever for condenser adjustment
- 3 Aperture diaphragm

### 8.4 Phase contrast

Phase contrast is for contrasting unstained specimens.

- Adjust the condenser height.
- Insert the phase contrast light ring slide into the fixture (56.1).  
(Caution: There are different sliders for the condensers!)
- Swing the phase contrast objective (engraved PH) with the weakest magnification into place using the objective nosepiece.
- Open the aperture diaphragm (56.4), "PH" mark.
- Bring the specimen into focus using the coarse and fine adjustment. If you have difficulty finding the specimen plane: Stop down the aperture diaphragm temporarily or use a stained specimen. To do so, move the condenser disk into the BF position or pull out the light ring slide. Open the aperture diaphragm.
- Use the light ring of the light ring slide that corresponds to the magnification (5, 10/20 or 40).

The phase contrast light ring slide is encoded. As soon as the slider is taken out of the bright-field position (available opening) and engaged in a phase contrast position (any light ring), the brightness increases. Conversely, switching from phase contrast to bright field decreases the brightness.

#### Note:

The light rings do not need to be centered. The superimposition of light rings and phase rings is adjusted at the factory.


#### CAUTION

Do not look into the eyepieces when toggling the contrast method!

During the toggling procedure, the maximum radiant power of the light source may be present on the eyepieces for a short time and can temporarily blind the user!

Fig. 56

- 1 Holder for light ring and light segment slide
- 2 Transmitted light illumination column
- 3 Stop lever for condenser height adjustment
- 4 Aperture diaphragm
- 5 Filter holder,  $\varnothing$  32 mm

## 8. Operation

### Possible faults

Specimen: too thick, too thin or too brightly stained; refractive indexes of the mounting medium and specimen are identical, so that there is no phase jump.

#### Note:

Cover slip placed in a wedge-shaped position, so that the centering of the light and phase ring is no longer effective.

Incorrect light ring.

Aperture diaphragm not open.

Condenser in incorrect height position. Slide the condenser further up or down to obtain optimum phase contrast.

The fluorescence filters are in the beam path.

Incorrect light ring slide.

(There is one slide each for the S80/0.30 and S40/0.45 condenser.)

IMC modulator in IMC position.

The S40/0.45 condenser has been switched with the S80/0.30 condenser.

Optics dirty.

### 8.5 Integrated Modulation Contrast (IMC)

The Integrated Modulation Contrast is a special form of oblique illumination that is based on the principle of Hoffman modulation contrast.

In it, the phase gradients of an unstained specimen are converted into amplitude differences using a modulator.

The impression of a three-dimensional image results, similar to the image of a microscope with interference contrast. Unlike interference contrast, however, the specimen can also be viewed through birefringent plastic materials such as Petri dishes.

Additional advantages of this imaging method are:

- High contrast
- High resolution
- Halo-free, variable-contrast relief image
- Long working distance of the condenser
- Easy assembly and adjustment
- Use with stained and unstained specimens

#### Important!

The IMC can be used in combination with both the S40/0.45 and S80/0.30 condensers.

The standard bright field and phase contrast objectives can be used for the IMC. The magnification range from 10x to 40x can thus be covered.

The IMC is only possible with a standard stage assembly, not with a stage that has been rotated by 180°.

All objectives with the eyepoint **D** and **C** are suitable for the IMC. These include all C PLAN objectives (discontinued) and many HC PLAN objectives.

The following objectives are particularly well suited:

HI PLAN 20x  
 HI PLAN 40x  
 HI PLAN I 20x  
 HI PLAN I 40x  
 N PLAN 20x  
 N PLAN 40x  
 FL PLAN 20x  
 FL PLAN 40x  
 HC PL FLUOTAR 10x  
 HC PL FLUOTAR 20x  
 HC PL FLUOTAR 40x

And the corresponding phase contrast objectives.

The following objectives with eyepoint C are suitable:

HI PLAN I 10x  
 N PLAN L 20x  
 N PLAN L 40x  
 HCX PL FLUOTAR L 20x  
 HCX PL FL L 40x

And the corresponding phase contrast objectives.

(Refer also to "Possible faults.")

For the IMC, it is necessary to use the IMC modulator (Fig. 57) and the IMC slit diaphragm slider (Fig. 58).

**Fig. 57** IMC modulator

1 Adjusting screw on the IMC modulator slider

### Important!


The IMC modulator can be used in two different ways. If the marking faces forwards (the "IMC" marking is visible), all objectives with the eyepoint **D** can be used. If you use objectives with the eyepoint **C**, you have to turn the slide around its own axis and insert it so that the marking faces backward. This brings the modulator to another level and the IMC is optimized for objectives with eyepoint **C**.

### Adjusting the IMC modulator (slider)

- Remove the empty slide in the stand if present.
- Insert the IMC modulator from the right according to the eyepoint of the objective. (Refer to above note.)
- Engage the slide in the IMC position. The IMC modulator is flush with one side of the stand.

**Fig. 58**

- 1 "TOP LEFT" marking on the diaphragm slide
- 2 Click stops
- 3 IMC slit diaphragm slider for S80/0.30
- 4 IMC slit diaphragm slider for S40/0.45
- 5 Adjusting screw on the IMC slit diaphragm slider

## 8. Operation

### Inserting the IMC slit diaphragm slider on the transmitted light illumination carrier

- Remove the dummy slide (where applicable) or the phase contrast slider in the transmitted light illumination carrier.
- Hold the IMC slide so that the marking "TOP LEFT" is at the top left "TOP LEFT" (58.1) and the rest of the marking faces **forwards**. The click stops are on the front longitudinal side of the slide.
- Insert the diaphragm slide from the right into the transmitted light illumination carrier.
- Engage the slide in the IMC position. The round opening is the bright-field position.

The IMC slit diaphragm slider is encoded. As soon as the slider is taken out of the bright-field position (available opening) and engaged in the IMC position, the brightness increases. Conversely, switching from IMC to bright field decreases the brightness.

### Adjusting the aperture diaphragms

- Open the aperture diaphragm completely.
- Set the brightness to a medium setting, as otherwise the light gap will appear too bright.
- Any activated filters should be deactivated.
- Swing the objective with the lowest magnification into position, usually the 10x objective.
- Slide the modulator (according to whether the eyepoint is **C** or **D**) into the modulator holder, IMC position.
- To adjust the slit width, adjust the slider of the IMC slit diaphragm slider to the position allocated on the objective, e.g. to the 10x marking for the 10x objective.
- Remove one eyepiece and insert the focusing telescope.

- The light gap appears as a bright strip on the gray image of the modulator. Focus the light gap using the focusing telescope.
- Adjust the position of the light gap using the adjusting screw on the right of the IMC slit diaphragm slider (58.5). A suitable Allen key is provided for this purpose.

#### Caution:

Do not unscrew the screws on the top of the slide!

- The light gap must be located entirely on the gray field. For the 10x objective, the image of the modulator and the light gap are almost the same size. Adjust the aperture diaphragm until the edge of the bright gap is close to the darker edge.
- Now, swing the other objectives into position in ascending magnification and check the position of the light gap. For smaller deviations, split the difference in the position. While doing so, always ensure that the objective magnification corresponds to the position of the slider on the IMC slit diaphragm slider.

**Optimizing the IMC:**

When using the objective with the highest magnification, in some circumstances, optimum adjustment of the aperture diaphragm may not be possible (i.e. the light gap cannot be placed entirely on the gray field so that an offset is noticeable, either into the white or the dark zone). In this case, you can also carry out fine adjustment using the adjusting screw on the IMC modulator slider (57.1).

Using the Allen key, turn this screw in order to minimize this offset (in order to superimpose the gray field and the illumination gap). Afterwards, swing the 10x objective back into position (and the other objectives, one after the other) and adjust the modulator as outlined above. After several iterations, no offset should occur.

This setting usually has to be configured only once.

Once the IMC as adjusted optimally, remove the focusing telescope and reinsert the eyepiece.

**Possible faults**

Unsatisfactory image quality due to using objectives without eyepoint **D**.

Try to improve the image quality as follows:

Turn the IMC modulator around (marking towards the rear). Adjust the light gap so that the overlap margin is sufficient to prevent glare.

The position of the aperture diaphragm is not optimal.

The IMC modulator or IMC slit diaphragm is not engaged in the IMC position.

Condenser in incorrect height position. Slide the condenser further up or down to obtain optimum phase contrast.

The fluorescence filters are in the beam path.

Incorrect condenser slide. There is one slide for the S80/0.30 condenser and one for the S40/0.45 condenser.

Incorrect objective. Eyepoints A and B (see objective engraving) are not suitable.

The slides are not introduced correctly. Follow the instructions.

## 8. Operation

### 8.6 Incident light fluorescence

#### Note:

For microscope with integrated incident light fluorescence device only.

For incident light fluorescence of transparent specimens, we recommend carrying out an adjustment in transmitted light first.

- Open the light stop by actuating the lever (59.1).
  - = Light stop swung out
  - = Light stop swung in
- Slide the filter block into the beam path (59.2).
- Position the specimen and bring it into focus. The illuminated field diaphragm is installed precentered. It does not need to be adjusted.

**CAUTION**

The lamp housing and the lamp may still be hot!

Fig. 59

- 1 Light stop
- 2 Filter block slide

**CAUTION**

Never look directly into the beam path!

**CAUTION**

Light sources pose a potential irradiation risk (glare, UV radiation, IR radiation).

#### 8.6.1 Switching on and adjusting the halogen and mercury lamps in the 106z lamp housing\*

For the 106z lamp housing, the direct filament image (for halogen bulb) and direct image of the light arc (for gas discharge lamps) and their mirror image are focused and adjusted to each other separately.

- Switch on the lamp on the supply unit.
- Open the light stop.
- Close the aperture diaphragm, to avoid light reflection of the TL-illumination.
- Move the filter block into the beam path.
- Place a sheet of white paper on the specimen stage.
- Bring the surface into rough focus using a dry objective with weak to medium magnification.
- Move the filter system\* or the reflector\* into the beam path.
- Open the shutter (where applicable) and remove any diffusion filters\* from the beam path.
- Place a sheet of paper on the specimen stage and focus on the surface using a dry objective with weak to medium magnification.


- Set the illuminated field diaphragm and aperture diaphragm to a middle position.
- Using a pen, make a mark on the paper and move the mark into the center of the illuminated field.
- Remove the objective or select a position that is not occupied.

An image of the light source now appears on the paper. Adjust the lamp as follows while observing the light source.


### Centering the Hg 100 W mercury lamp\*

- On the paper, you see the direct image of the light arc and the mirror image, which are usually pushed together.
- Bring the direct image into focus using the collector (68.6).
- Using the adjusting knobs (68.2, 68.4) on the rear of the lamp housing, move the mirror image of the light arc to the side or out of the beam path entirely. The focused image of the light arc remains visible (Fig. 65).
- Using the adjusting knobs (68.1) and (68.5), position the direct image of the light arc in the middle of the centering surface; the bright tip of the light arc, the cathode spot, should be just off center (Fig. 66).
- Now, swing the mirror image of the light arc back into position using the adjusting knobs (68.2) and (68.4) and bring it into focus using the reflector (68.3).
- Align the mirror image so that is symmetrical to the direct image (Fig. 67). To do, again use the adjusting knobs (68.2) and (68.4).


**Fig. 65** Direct image of the light arc focused, but off-center (in reality, the image is blurrier).

**Fig. 66** Direct image of the light arc in the correct position (in reality, the image is blurrier)

**Fig. 67** Direct image of the light arc and mirror image in correct position (in reality, the image is blurrier).

## 8. Operation

The V-shaped images of the light arcs of the direct image and the mirror image can be superimposed.

### Caution!

The bright tips of the light arcs (the cathode spot) must never be projected superimposed, as in that case there is a danger of explosion from overheating.


### Caution!

For older lamps, the structure of the light arc is no longer clearly visible. The image is then more similar to that of an Hg 50 lamp. The image and mirror image can no longer be superimposed accurately. In this case, superimpose both images.

**Fig. 68** 106z lamp housing

- 1 Height adjustment of the lamp
- 2,4 Height and lateral adjustment of the mirror image
- 3 Focus of the reflector
- 5 Lateral adjustment of the lamp
- 6 Collector (focusing of the lamp image)

- Defocus the image via the collector using the knob (68.6) until the image of the light arc and the mirror image are no longer visible and the image is uniformly illuminated.

### Note:

For the Leica DM IL M LED the lamp housing LH115 LED is provided (→ p. 38).

### Note:

If a Leica SFL100 is used for fluorescence illumination, please regard the relevant operating manual.

### Possible faults

#### Weak fluorescence, too little brightness:

- Incorrectly stored, too old or bleached specimens.
- Quick bleaching of specimens (e. g. for FITC).
- Unspecific filter combination.
- Objective with too low of a numerical aperture.
- Eyepiece magnification too high.
- Spent lamp.
- Room where microscope is located too bright.

- Trinocular tube: incorrect beam splitter setting.
- Secondary light due to reflection at condenser.

### Low-contrast image due to:

- Excitation bandwidth too wide.
- Unspecific staining.
- Fluorescent mounting medium.
- Self-fluorescence of the objective or immersion oil.
- Dirty glass surfaces.

### Image background is not dark

- The fluorescence light is reflected by the TL-LED.  
Quick relief: a) close the aperture diaphragm completely b) put a dark stop into the TL filter slot c) if available move the phase contrast slider in position 40x d) cover your specimen with a dark lid.

# 9. Troubleshooting

Problem	Cause/remedy
---------	--------------

## Stand

The microscope does not respond.	<ul style="list-style-type: none"><li>▶ Ensure that the AC outlet has power.</li><li>▶ Make sure that the stand is connected to the power supply.</li><li>▶ Check the cable connections.</li></ul>
----------------------------------	--

## Illumination

The image is completely dark.	<p><u>Transmitted light:</u></p> <ul style="list-style-type: none"><li>▶ Make sure that the LED in the built-in transmitted light illuminator is not defective. Notify technical service if necessary.</li></ul> <p><u>Fluorescence:</u></p> <ul style="list-style-type: none"><li>▶ Make sure that the lamp housing is connected to the microscope and is not defective.</li><li>▶ Inform Leica Service and have the supply unit fuse checked.</li></ul>
-------------------------------	---

The image is not uniformly illuminated.	<p><u>Transmitted light and fluorescence:</u></p> <ul style="list-style-type: none"><li>▶ Remove all unneeded filters from the light path.</li><li>▶ Center the lamp (106z lamp housing)<ul style="list-style-type: none"><li>▶ Change out the old lamp.</li></ul></li></ul>
---	--

The illumination "flickers".	<p><u>Fluorescence:</u></p> <ul style="list-style-type: none"><li>▶ Be sure that there is no loose connection at the power supply.</li><li>▶ Change out the old lamp.</li></ul>
------------------------------	---

Problem	Cause/remedy
Fluorescence: The lamp does not light up immediately after being switched on.	<ul style="list-style-type: none"><li>▶ The external power supply must be switched on repeatedly.</li><li>▶ Allow hot mercury lamps to cool down before switching them off again.</li></ul>
<b>Focus</b>	
The specimen cannot be brought into focus.	<ul style="list-style-type: none"><li>▶ Use the correct immersion medium.</li><li>▶ Lay the specimen with the cover slip toward the bottom.</li><li>▶ Make sure that the cover slip thickness is correct and that it meets the specifications on the objective.</li></ul>
<b>Transmitted light</b>	
The image is not uniformly illuminated.	<ul style="list-style-type: none"><li>▶ Be sure that the right objective is being used.</li><li>▶ Aperture diaphragm opened or closed too far.</li><li>▶ Condenser in incorrect height position.</li><li>▶ Light ring or IMC component switched on by mistake.</li></ul>
Unwanted stray light.	<ul style="list-style-type: none"><li>▶ Clean the specimen and neighboring lens surfaces.</li></ul>

## 9. Troubleshooting

Problem	Cause/remedy
---------	--------------

### Phase contrast

Phase contrast cannot be adjusted.

- ▶ The specimen is too thick, too thin or too brightly stained.
- ▶ Refractive indexes of the mounting medium and specimen are identical, so that there is no phase jump.
- ▶ The cover slip is not placed uniformly.
- ▶ Check that the correct light ring is positioned.
- ▶ Check that the correct light slider is positioned.
- ▶ The S40/0.45 condenser has been switched with the S80/0.30 condenser.
- ▶ Open the aperture diaphragm completely.
- ▶ IMC modulator in IMC position.

### Integrated modulation contrast

IMC cannot be adjusted.

- ▶ Check that the correct objective is positioned (eyepoint C or D).
- ▶ Check the position of the aperture diaphragm.
- ▶ Check that the IMC modulator and the IMC slit diaphragm slider are inserted correctly and engaged in the IMC position.
- ▶ Check that the correct condenser (S40/0.45 or S80/0.30) is positioned and that the condenser height is correct.
- ▶ Switch off the fluorescence filter.
- ▶ Open the aperture diaphragm completely.
- ▶ IMC modulator in IMC position.

Problem	Cause/remedy
<b>Fluorescence</b>	
The image is completely dark (no fluorescence).	<ul style="list-style-type: none"> <li>▶ Open the shutter.</li> <li>▶ Check the antigen-antibody combination.</li> <li>▶ Insert a new lamp.</li> </ul>
The fluorescence is too weak.	<ul style="list-style-type: none"> <li>▶ Center the lamp.</li> <li>▶ Insert a new lamp.</li> <li>▶ Unsuitable specimen (stored incorrectly, too old, bleached).</li> <li>▶ Unspecific filter combination.</li> <li>▶ Objective with too low of a numerical aperture.</li> <li>▶ Eyepiece magnification too high.</li> <li>▶ Room where microscope is located too bright.</li> <li>▶ Incorrect beamsplit on the trinocular tube.</li> <li>▶ Secondary light due to reflection at condenser.</li> </ul>
The image has too little contrast.	<ul style="list-style-type: none"> <li>▶ Excitation bandwidth too wide.</li> <li>▶ Unspecific staining.</li> <li>▶ Fluorescent mounting medium.</li> <li>▶ Self-fluorescence of the objective or immersion oil.</li> <li>▶ Dirty glass surfaces.</li> </ul>
Image background is not dark.	<ul style="list-style-type: none"> <li>▶ Too much stray light from environment. : Room where microscope is located is too bright</li> <li>▶ The fluorescence light is reflected by the TL-LED. Quick relief: a) close the aperture diaphragm completely b) put a dark stop into the TL filter slot c) if available move the phase contrast slider in position 40x d) cover your specimen with a dark lid.</li> </ul>

# 10. Care of the Microscope

**Caution!**

Unplug the power supply before performing cleaning and maintenance work!  
Protect electrical components from moisture!

Microscopes in warm and warm-damp climatic zones require special care in order to prevent fungus contamination. The microscope should be cleaned after each use, and the microscope optics should be kept painstakingly clean.

## 10.1 Dust cover

**Note:**

To protect against dust, cover the microscope and accessories with the dust cover after each use.

**Caution!**

Let the microscope and lamp housings cool down before covering the stand with a dust cover. The dust cover is not heat-resistant. In addition, condensation may occur.

## 10.2 Cleaning

**Caution:**

Residual fiber and dust can create unwanted background fluorescence during fluorescence microscopy.

### Cleaning coated parts

Dust and loose dirt particles can be removed with a soft brush or lint-free cotton cloth.

Clinging dirt can be cleaned as necessary with a low-concentrated soap solution, petroleum ether or ethyl alcohol.

For cleaning coated parts, use a linen or leather cloth that is moistened with one of these substances.

**Caution!**

Thinners containing acetone, xylene or nitrogen can harm the microscope and thus must not be used.

Fig. 69 Microscope with dust cover

Test cleaning solutions of unknown composition on a less visible area of the unit first. Be sure that coated or plastic surfaces do not become matted or etched.

### Cleaning the specimen stage

Remove light-colored spots on the specimen stage by rubbing with paraffin oil or acid-free Vaseline.

### Cleaning Glass Surfaces and Objectives

Glass surfaces, and particularly objectives, are always to be cleaned as described in the brochure "Cleaning of Microscope Optics". You can download the information from

<http://www.leica-microsystems.com/products/light-microscopes/life-science-research/inverted-microscopes>

Select the DM IL LED microscope and switch to the "Download" page.

You can also contact our Technical Service with any questions.

### Removing immersion oil

#### Caution!

Follow safety notes for immersion oil!

First, wipe off the immersion oil with a clean cotton cloth, and then re-wipe the surface several times with ethyl alcohol.

### 10.3 Handling acids and bases

For examinations using acids or other aggressive chemicals, particular caution must be taken.

#### Caution:

Never allow the optics and mechanical parts to come into direct contact with these chemicals.

# 11. Technical Description

Due to basic principles of physics and the physiology of the eye, all methods – not just those in microscopy – have performance limits. Therefore, observe the following information for correct use of the microscope.

## Performance data of the objectives

The Leica DM IL LED microscope is based on the tube length  $\infty$  (infinity) and a tube lens focal length of  $f = 200$  mm.

### ! Caution:

Therefore, only objectives with the engraving  $\infty$  and the thread dimension M 25 may be used.

## Objective marking

Examples and meaning of the symbols:

$\infty / -$   
HI PLAN 10x/0.22

$\infty / 0.17$   
N PLAN 40x/0.65

$\infty / 0 / D$   
N PLAN 50x/0.75

- $\infty$  Objective for infinite tube length ( $\infty$ ).
- The objective can be used **with and without** cover slip.

**0.17** The objective may be used with a cover slip of the standard thickness 0.17 mm only. If the cover slip is missing or the thickness of the cover slip deviates greatly, significantly decreased performance will result, particularly for high objective apertures (see below).

**0** Application **without** cover slip, e.g. for cell smears, incident light. Not suited for inverted microscopes.

**D (or A, B, C)** Eyepoint of the objective (important for Integrated Modulation Contrast IMC, for example).

**L** Long working distance.

**10x/0.22** Magnification and aperture. The aperture determines the resolution, focus depth, contrast and brightness of the microscope image. Objectives with built-in iris diaphragm have the maximum and minimum aperture engraved, e.g. 0.85-0.55.

**Color-coding of the objectives**

In accordance with DIN/ISO standards, the magnification of each objective is indicated by a color ring along the outer edge:

100x 125x 150x 160x	63x	40x 50x	25x 32x	16x 20x	10x	6.3x	4x 5x	2.5x	1.6x
white	dark blue	light blue	dark green	bright green	yellow	orange	red	brown	gray

Immersion objectives are also marked by a second, lower color ring:

**Black** Oil or Imm (= universal objective for oil, water, glycerin)

**White** Water

**Orange** Glycerin

**Performance data of the eyepieces**

The following eyepieces are in the product line for the Leica DM IL LED:

Leica eyepiece type	Magnification/ Field number	Characteristics <sup>1)</sup>	
HC PLAN	10x/20	☞	M
HC PLAN	10x/20	☞	M
HC PLAN	12.5x/16	☞	M
HC PLAN	10x/20	☞	MF

Eyepiece tube diameter: 30 mm

- <sup>1)</sup> ☞ = With removable or invertible protection screen for eyeglass wearers and non-eyeglass wearers
- M = Adjustable eyelens (dioptric equalization) and fixture for graticules with diameter of 19 mm or 26 mm for HC eyepieces.
- MF = With illuminated graticule

**Fig. 70** HI PLAN objective set

**Fig. 71** Eyepiece pairs

- 1 HC PLAN 10x/20 ☞
- 2 HC PLAN 10x/20 ☞ M

## 11. Technical Description

### Eyepiece field number

For a certain microscope configuration, a certain eyepiece field number (see below), e.g. 20, must not be exceeded. If the maximum field number is exceeded, unwanted blurring may occur on the image field border and/or shading (vignetting) of the edge of the image, → following pages.

The eyepiece field number (FOV, field of view) designates the diameter of the intermediate image in the eyepiece in mm, i.e. the diameter of the circular diaphragm that limits the image and is located within the eyepiece.

This FOV is specified on the eyepiece after the magnification, e.g. 10x/20.

For the Leica DM IL LED microscope, the recommended maximum is **FOV 20**.

#### Note:

The maximum permitted eyepiece field number for a given equipment configuration is calculated from the following instrument data:

<b>Field performance of objectives</b>	<b>see below</b>
<b>Field performance of intermediate module(s)</b>	→ p. 66
<b>Field number of tube</b>	→ p. 67
<b>Condenser properties</b>	→ p. 69

The **lowest** occurring value is always decisive. If, for example, the intermediate modules allow a field number of only 20, but the objectives and tube allow 25, only eyepieces up to FOV 20 are permitted. In this case, eyepieces with FOV 25 can lead to vignetting. Specifically, the following applies:

The diameter of the **viewable specimen area** is calculated by dividing the diameter of the field of view by the magnification of the objective and the magnification factor of the stand optics.

Example:

10x/20 eyepiece

PLAN 4/0.10 objective

Magnification factor of the Leica DM IL LED

1x stand optics

Viewable specimen area

$$\frac{20 \text{ mm}}{4 \times 1} = \varnothing 5 \text{ mm}$$

The **total magnification** of the microscope is calculated by multiplying the eyepiece magnification by the magnification of the objective and the magnification factor of the stand optics.

Example:

10x/20 eyepiece

PLAN 4/0.10 objective

1x magnification factor

Total magnification  $10 \times 4 \times 1 = 40\times$

### Field performance of objectives

The field performance of objectives is not engraved on the objectives. It can fluctuate somewhat within a class, e.g. the low objective magnifications may well have slightly higher values than the guiding values listed below:

#### Objective series Max. recommended Eyepiece field number

	15	20	22	25
Achromats	████████			
HI PLAN achromats	████████			
APO L apochromats	████████			
N PLAN planachromats	████████		██	
PL FLUOTAR® semiapo.	████████	████████		
PL APO planapochromats	████████	████████	██	

You can always obtain an updated data sheet with all Leica objectives from your Leica representative.

**Field performance of intermediate modules**

The maximum permitted field performance intermediate modules is calculated from the model code listed both in the following table and on your invoice. This model code includes 2 values separated by a slash.

The first value is a relative measurement (vertical index) for the height of the module. If you multiply the vertical index by a factor of 15, you obtain the increase in height of the tube eyepiece or instrument height in **mm**. The second value is the maximum possible field number with this module.

- L 3/25 magnification changer
- L 3/20 discussion attachment (2 observers)

**Performance data of the filters**

Filters	Application
<b>Neutral density filters N</b>	Neutral density filters are used to attenuate the light without affecting the color temperature. The engraved value, e.g. N16, specifies the attenuation value. N16 thus means reduction to $1/16 = 100/16 = 6.25\%$ transmission.
<b>Green filter, GR panchromatic</b>	Contrast enhancement for black/white images.
<b>DLF</b>	Conversion filter (blue daylight filter, similar to the CB12) for color photography with daylight film, built into the filter magazine.
<b>ALF</b>	Artificial light filter for color photography with artificial light film to enhance the color contrast.

## 11. Technical Description

### Performance data of the tubes

The tube change is identical to the upright stands.

The tubes are rotatable and interchangeable.

### DM ILB binocular tube

The binocular tube consists of a body that carries the tube change ring on the bottom. The tube lens has a factor of 1x. The Siedentopf binocular head allows the interpupillary distance to be adjusted from 55 mm to 75 mm while maintaining a constant tube length. The viewing angle is 45°. The tube has an adjustable eyepiece tube that enables a field number of 20.

### DM ILT trinocular phototube

The trinocular tube consists of a body that carries the tube change ring on the bottom. The tube lens has a factor of 1x. The Siedentopf binocular head allows the interpupillary distance to be adjusted from 55 mm to 75 mm while maintaining a constant tube length. The viewing angle is 45°. The tube has an adjustable eyepiece tube that enables a field number of 20.

The side documentation output is operated with HC components only.


The tube contains a selectable mirror with two positions:

- a) 100% of the light to the binocular tube.
- a) 100% of the light to the photo tube.

The optical axis of the documentation output is shifted 88 mm to the left. This allows an unob-

**Fig. 72** Tubes

- 1 DM ILT trinocular phototube
- 2 DM ILB binocular tube

structured view of the specimen, even with this tube.

**Tubes from the upright microscope product line**

When using tubes from the Leica DM1000-3000 product line, you need a tube adapter. The DM IL/L tube adapter is a 60 mm long intermediate tube without optics for adjusting the eyepoint. A tube change ring is attached to the bottom, and the top has a tube change interface for the DM1000-3000 tubes.

**Field number of DM1000-3000 tubes**

The model code of the DM1000-3000 tubes also contains a number combination with concrete information about the maximum permitted eyepiece field number, e.g.

HC LB **0/3/4** binocular tube.

Here, the numbers **0/3/4** mean the maximum permitted height value of the intermediate modules (= vertical index) for the eyepiece field numbers **25, 22 and 20**.

In the above example, this means:

1. Number (0): Field of view 25 is possible only when the tube is adapted directly to the stand, without an intermediate system.
2. Number (3): Field number 22 is permitted only up to vertical index 3, e.g. L3/25 magnification changer.
3. Number (4): Field number 20 up to a max. vertical index of 4, e.g. 2 L 2/25 ErgoModules.

If the number is replaced by a dash -, e.g. LMP -/-/7 monocular tube, the tube cannot be used at all for the corresponding field number, i.e. in the

example, not for FOV 25 and 22, while FOV 20 is permitted up to index 7.

Exceeding the permitted values can cause vignetting (image shading around the edges).

**HC marking:** Only eyepieces of type **HC PLAN** and widefield 16x and 25x can be used. If the HC marking is missing, eyepieces of type **Leica L PLAN** may be used.

Additional examples:

0/4/4 Field number 25 is possible only when tube is adapted directly to stand (vertical index of the intermediate modules is 0), as long as the suitable objectives are used.

Field number 20 and 22 is possible up to a vertical index of 4, for example with the fluorescence equipment. Adapting another module would not be permitted; as a workaround, a tube with the following characteristic values could be used:

4/5/7 Field of view performance 25 is possible up to vertical index 4 (e. g. 2 L2/25 ErgoModules or L3/25 magnification changer). Field of view 22 is possible up to vertical index 5, field of view 20 is possible up to vertical index 7 (e.g. LRF 4/20 plus L3/245 magnification changer).

-/-/7 The tube enables a field of view up to 20 mm. If intermediate modules are used, the sum of its height values must not exceed 7.

## 11. Technical Description

### Performance data of the condensers

#### S80/0.30 condenser

For vials up to 80 mm in height and objectives with a numerical aperture up to 0.55. In the phase contrast method, the adjustment of the light ring to the phase ring is optimal up to a liquid height of 2 mm.

#### S40/0.45 condenser

For vials up to 55 mm in height and objectives with a numerical aperture up to 0.80. In the phase contrast method, the adjustment of the light ring to the phase ring is optimal up to a liquid height of 2 mm.

**Fig. 73** Condensers

- 1 S40/0.45 condenser
- 2 S80/0.30 condenser

Application options of the S80/0.30 and S40/0.45 condensers:

Illumination method	S80/0.30 objective	Light rings/ accessories	S40/0.45 objective	Light rings/ accessories
Bright field	2.5x 4x–100x	(with diffusion filter) –	2.5x 4x–100x	(with diffusion filter) –
Phase contrast	5x 10x–20x 40x	Phaco 0 Phaco 1 Phaco 2	5x 10x–20x 40x	Phaco 0 Phaco 1 Phaco 2
Integrated modulation contrast	All objectives with eyepoint C and D		All objectives with eyepoint C and D	IMC slide

**Incident light fluorescence illumination\***

The Leica DM IL LED microscope should preferably be equipped with mercury gas discharge lamps for incident light fluorescence due to the improved brightness of the image, but can also be operated with a 12 V 100 W halogen bulb.

**Performance data of the lamp housings**

**106z lamp housing\***

Lamp housing with a centerable and focusable reflector and 4 to 6-lens collector. A quartz collector is available by request. The following lamps can be used with their special sockets:

- 100 W high-pressure mercury burner (direct current, stabilized)
- 100 W high-pressure mercury burner (direct current, stabilized Type 103 W/2)

**107/2 lamp housing**

The shield connection of the lamp housing is connected to the equipotential bonding point of the 12 V 100 W supply unit. This lamp housing for incident light and transmitted light has a single-lens fixed collector and a fixed 12 V 100 W lamp.

**Leica EL6000 external compact light source**

For indoor use only.

Supply voltage:	100-240 V AC
Frequency:	50/60 Hz
Power input:	max. 200 W
Fuses:	5x20, 2,5 A, slow, breaking capacity H → manual EL6000
Ambient temperature:	0°-40°C
Relative humidity:	10-90% non-condensing
Overvoltage category:	II
Pollution degree:	2 (see separate manual)

**General specifications**

For indoor use only.

Supply voltage:	100-240 V AC
Frequency:	50/60 Hz
Power input:	max. 16 VA
Fuses:	0.63 A, slow-blowing, breaking capacity H, 250 VAC, size 5x20 mm
Ambient temperature:	15-35°C
Relative humidity:	max. 80% up to 30°C non-condensing
Overvoltage category:	II
Pollution degree:	2

Ordernumber fuses: 11 362 150 010 063

Typ	Typical bulb life <sup>+) </sup>
100W high-pressure mercury burner (direct current)	200 hrs.
100W high-pressure mercury burner type 103 W/2 (direct current)	300 hrs.

+ ) Please observe the data sheets of the lamp manufacturer.

# 12. Index

## A

Adjusting the aperture diaphragms 52  
 Adjusting the halogen and mercury lamps 54  
 Ambient conditions 22  
 Aperture diaphragm 18, 47, 49  
 Attachable mechanical stage 30

## B

Bright-field illumination 47  
 Brightness control 43  
 Brightness controller 17

## C

Camera 40  
 Cellfactory 25, 26, 27  
 Cleaning 62  
 Condenser 15, 18, 25, 47, 70  
 Condenser height adjustment 18, 49  
 Contrast methods 14

## D

Dimensions 23  
 Disposal 13

## E

EL6000 39, 71  
 Electrical connection 31  
 Electrical safety 10  
 ErgoModule 42  
 Eyepiece field number 66  
 Eyepieces 17, 29, 44, 65

## F

Filter block slide 32  
 Filter holder 49  
 Filters 18, 30

Fluorescence filter blocks 32  
 Fluorescence filter block slide 19  
 Focus 18  
 Focusing 44

## G

Gas discharge lamps 36  
 Graticule 29, 45

## H

HC L1T tube 28  
 Height compensation plate 38  
 Hg 100 W mercury lamp 71

## I

IL/L tube adapter 17, 28, 42  
 IMC modulator 19, 34  
 IMC slit diaphragm slider 34, 51, 52  
 Incident light axis 14  
 Incident light fluorescence 54  
 Incident light fluorescence device 19  
 Incident light fluorescence illumination 71  
 Integrated Modulation Contrast 50  
 Intermediate modules 67

## L

Lamp housing 18  
 Lamp housing 106z 54, 71  
 Lamp housing 107/2 71  
 Lamps 19  
 LED transmitted light illumination 43  
 Leica EL6000 39, 71  
 Leica SFL100 35, 56  
 Light ring 49

Light source 31

## M

Modulation slider 19  
 Multiple discussion attachment 41

## O

Objective marking 64  
 Objective nosepiece 14, 18  
 Objectives 18, 29, 64

## P

Phase contrast 49  
 Phase contrast slider 19, 33

## S

Safety notes 10  
 SFL100 35, 56  
 Specimen stage 18, 30  
 Stages 14  
 Storage 22  
 Supply unit 19, 38  
 Symbols 6

## T

Technical data 10, 71  
 Total magnification 66  
 Transmitted light 43, 47  
 Transmitted light axis 14  
 Transmitted light illumination carrier 25  
 Transmitted light illumination unit 17  
 Transport 12, 22  
 Tube 14, 17, 27, 44, 68

## W

Weight 22, 23

# 13. EU Declaration of Conformity

To download the EU Declaration of Conformity for your product use the link

<http://www.leica-microsystems.com/products/light-microscopes/life-science-research/inverted-microscopes>

or

<http://www.leica-microsystems.com/products/light-microscopes/industrial-materials/inverted-microscopes>

Select the type of your microscope and switch to the Download page.

- Administrative Measures on the Control of Pollution Caused by Electronic Information Products -

部件名称 Name of the part	有毒有害物质或元素 Hazardous substances						
	铅 (Pb)	汞 (Hg)	镉 (Cd)	六价铬 (Cr <sup>6+</sup> )	多溴联苯 (PBB)	多溴二苯醚 (PBDE)	
印刷电路板 printed circuit boards	X	0	0	0	0	0	
电子元器件 electronic components	X	0	0	0	0	0	
机械部件 mechanical parts	X	0	0	X	0	0	
光学元器件 optical components	X	0	X	0	0	0	
电缆 cables	0	0	0	0	X	X	
光源 light sources	0	X	0	0	0	0	

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Indicates that the concentration of the hazardous substance in all materials in the parts is below the relevant threshold of the SJ/T 11363-2006 standard.

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Indicates that the concentration of the hazardous substance of at least one of all materials in the parts is above the relevant threshold of the SJ/T 11363-2006 standard.

Note: The actual product may or may not include in all the part types listed above

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Copyright © Leica Microsystems CMS GmbH • Ernst-Leitz-Straße • 35578 Wetzlar • Germany 2012 • Tel. (06441)29-0 • Fax (06441)29-2599 • LEICA and the Leica logos are registered trademarks of Leica IR GmbH.