

# Technical note

## J.T.Baker® Binding affinity of PROchievA™ resin to IgG from different species



It is generally known that protein A, a 42 kDa surface protein of *Staphylococcus aureus*, has different binding affinities towards antibodies depending on the host species and isotypes due to amino acid sequence variation of the IgG H-chains [1]. BAKERBOND® PROchievA™ affinity resin has an agarose backbone modified with an engineered protein A ligand that is designed to provide high binding capacity and stability for use in mAbs and Fc fusion protein purification [2]. The objective of this study was to evaluate the binding profiles of PROchievA™ resin against IgG from different species to show feasibility of using the resin for different applications.

## MATERIALS AND METHODS

A frontal chromatographic technique was used to determine dynamic binding capacity (DBC) at 10% breakthrough according to the methods in Table 1 and 2. Purified IgG from five different species, human, rabbit, goat, bovine, and mouse were loaded onto the PROchievA™ column to determine dynamic binding capacity. Additionally, for human IgG capacity, three samples were tested, which included IgG<sub>4</sub> and two samples of IgG<sub>1</sub>. For purification of samples, 5 ml pre-packed PROchievA™ columns were used and 1 ml pre-packed PROchievA™ columns were used for DBC measurement.

To test the binding capacity during a purification process, human IgG<sub>4</sub> sample in clarified cell culture supernatant (3.57 mg/ml) was loaded onto the pre-packed 5 ml column and flowthrough fractions are analyzed to determine DBC at 10% breakthrough. The operating conditions are shown in Table 1 and 2.

Step	Buffer	Volume	Flow rate
Pre-equilibration wash	100mM HAc, pH 3.4	5 CV	1 ml/min
Equilibration	1X PBS, pH 7.4	10 CV	1 ml/min
Sample load	Purified samples (2 mg/ml) in 1X PBS, pH7.4	Until -15% breakthrough	0.13ml/min (8 min Residence Time)
Wash	1X PBS, pH 7.4	5 CV	1 ml/min
Elution	100mM HAc, pH 3.4	10 CV	1 ml/min
CIP	0.5N NaOH	5 CV	1 ml/min

**TABLE 1:** Operating condition for DBC measurement with 1 ml pre-packed column. (Samples: human IgG, rabbit, goat and bovine IgG)

Step	Buffer	Volume	Flow rate
Pre-equilibration wash	100mM HAc, pH 3.4	5 CV	1 ml/min
Equilibration	50mM Tris, 0.9M NaCl, pH 8.6	10 CV	1ml/min
Sample load	Purified mouse IgG, (1 mg/ml) in 50mM Tris, 0.9M NaCl, pH 8.6	Until -15% breakthrough	0.5 to 0.13 ml/min (2 - 8 min Residence Time)
Wash	50mM Tris, 0.9M NaCl, pH 8.6	5 CV	1 ml/min
Elution	100mM HAc, pH 3.4	10 CV	1 ml/min
CIP	0.5 N NaOH	5 CV	1 ml/min

**TABLE 2:** Operating conditions for DBC measurement of mouse IgG, with 1 ml pre-packed column (Sample: mouse IgG)




FIGURE 1: Chromatogram of DBC analysis of human IgG and rabbit IgG

### AFFINITY TOWARDS VARIOUS IgG

PROchievA™ showed varying degrees of binding affinity for IgG from different species. Chromatogram profile exhibiting such behavior is shown in Figure 1. The DBC values of human IgG<sub>1</sub> and rabbit IgG are significantly different at a same residence time. Table 3 shows the DBC data of PROchievA™ for human IgG, rabbit, goat, bovine and mouse IgG.

IgG	DBC (mg/mL)	Residence time (min)
Human IgG <sub>1</sub> (Sample 1)	68.8	8
Human IgG <sub>1</sub> (Sample 2)	68.9	8
Human IgG <sub>2</sub>	55.4	6
Rabbit IgG	36.3	8
Goat IgG	66.0	8
Bovine IgG	61.7	8
Mouse IgG <sub>1</sub>	31.6	8

TABLE 3: DBC of PROchievA™ towards various IgG from different species

### COMPARISON OF COMMERCIALY AVAILABLE PROTEIN A RESINS' DBC FOR MOUSE IgG<sub>1</sub>

Mouse IgG<sub>1</sub> has various applications, especially for diagnostics and assay development, while finding protein A resin with high capacity for mouse IgG<sub>1</sub> is challenging. Table 4 shows DBC of PROchievA™ at different residence times. To compare the performance, experiments were performed on competitor's new generation protein A resins at 2 min residence time.

Protein A resin	Residence time (min)	Mouse IgG <sub>1</sub> DBC (mg/mL)
PROchievA™	8	31.6
	6	30.3
	4	28.3
	2	20.3
Competitor resin 1	2	7.5
Competitor resin 2	2	3.7

TABLE 4: DBC of different agarose protein A resin for mouse IgG<sub>1</sub> sample

## CONCLUSION

The Fc-specific BAKERBOND® PROchievA™ showed high dynamic binding capacity for IgG antibodies from different species and two subclasses, IgG<sub>1</sub> and IgG<sub>4</sub>. For example, compared to the competitor's resins, the DBC for mouse IgG<sub>1</sub> was 3 to 5 times higher by PROchievA™ at the same process conditions.

## Materials used with ordering information

Materials	Avantor part number
PROchievA™ 1 ml column	C789-11
PROchievA™ 5 ml column	C789-18
Sodium Phosphate Dibasic anhydrous	3826
Sodium Phosphate Monobasic monohydrate	3802
Sodium Chloride	3625
Acetic acid	9526
Tris (Base)	4102
TrisHCl	4106
0.5N NaOH	0329



## Reference

1. Atkins KL, Burman JD, Chamberlain ES, Cooper JE, Poutrel B, Bagby S, Jenkins AT, Feil EJ, van den Elsen JM: S. aureus IgG-binding proteins SpA and Sbi: host specificity and mechanisms of immune complex formation. Mol Immunol. 2008 Mar;45(6):1600-11. doi: 10.1016/j.molimm.2007.10.021. Epub 2007 Dec 3. PMID: 18061675
2. Avantor: J.T.Baker® BAKERBOND® PROchievA™ Recombinant Protein A, Affinity Chromatography Resin. <https://us.vwr.com/store/product/30493609/j-t-baker-sup-sup-bakerbond-sup-sup-prochievatm-recombinant-protein-a-affinity-chromatography-resin>