




Chromatography consumables

Protecting from PFAS – together, we can do it

Chromatography and consumables for routine testing

Minimize risks of contamination with consumables

Due to the increased use of PFAS, there are many different risks of contamination to the analysis. Care must be taken during sample preparation, sample handling and analysis, particularly as detection limits are in the low to sub-ppt range. It is not only recommended to ensure risks of contamination is minimized in the laboratory but also with your chromatography solution. You can minimize the risk of contamination to your chromatographic equipment by using our PFAS kit (includes PFAS-free tubing, fittings), solvent filter inlets and our isolator column. Changing the tubing and using an isolator column as part of the solution, to separate out the analytical PFAS from the system PFAS, will result in more accurate quantitation.



PFAS columns overview

Isolator columns

The use of an isolator column is necessary when analysing PFAS, as it helps to separate background contaminants from sample analytes, allowing for a more accurate sample result. For this we recommend Thermo Scientific™ Accucore™ C18 column.

Analytical columns

For high retentivity and large volume direct injections a fully porous column such as Thermo Scientific™ Acclaim™ C18 column has excellent peak shape for PFAS analyses. A fully porous column with larger surface volume has more stationary phase to appropriately separate the analytes. For methods where smaller sample injection volumes are employed, Accucore RP-MS column is recommended. The Accucore column is a solid-core that provides excellent chromatographic separation and maintains robustness in challenging water matrices.

Guard columns

Use a guard column to prolong the lifetime of your analytical column. A guard column is a protective column installed between the injector and the analytical column. It serves to remove the impurities, particulate contaminants, suspended solids and highly absorptive compounds from samples from reaching the analytical column. Read more about prolonging the lifetime of analytical column [here](#). The Acclaim and Accucore guard columns are used for PFAS within the industry.

GC columns

For neutral (or volatile) PFAS, such as fluorotelomer alcohols (FTOHs), N-substituted fluoroalkylsulfonamides (FOSAs), the trifluoropropylmethyl polysiloxane mid-polarity phase of Thermo Scientific™ TraceGOLD™ TG-200MS provides a unique selectivity for fluorinated compounds for an adequate chromatographic separation of the targets combined with Gaussian peak shapes. Additionally, it also showed excellent peak shapes for ionic perfluoroalkylcarboxylic acids (PFCAs) and fluorotelomer carboxylic acids (FTCAs). The low polarity, arylene content of TraceGOLD TG-5SiIMS phase improves thermal stability and provides the lowest bleed, making it ideal for screening and unknown identification.

Recommended columns

Isolation columns

Accucore C18 column

Analytical columns

Acclaim C18 column

Accucore RP-MS column

Guard columns

Acclaim guard column

Accucore guard column

GC columns

TraceGOLD TG-5SiIMS column

TraceGOLD TG-200MS column

Columns solutions

Choose from one of the below products for your PFAS applications

Analytical and guard columns for U.S. EPA method 8327

Group	Description	Quantity	Cat. no
Column	Acclaim 120 C18 Column, 2.2 µm, 2.1 mm × 100 mm	Each	10039-116

Analytical and guard columns for U.S. EPA method 1633, 533 and 537

Group	Description	Quantity	Cat. no
Column	Acclaim 120 C18 Column, 2.2 µm, 2.1 mm × 100 mm	Each	10044-610
Column	Acclaim Polar Advantage HPLC Column, 2.2 µm, 2.1 mm × 100 mm	Each	10040-572
Guard cartridges	Acclaim 120 C18 Guard Cartridge Column, 5 µm, 2.1 mm × 10 mm	Pack of 2	10044-622
Guard holder	Acclaim Guard Holder and Coupler	Each	10046-750
Column	Accucore RP-MS Column, 2.6 µm, 2.1 mm × 100 mm	Each	10039-116
Column	Accucore RP-MS Column, 2.6 µm, 2.1 mm × 150 mm	Each	10039-108
Guard cartridge	Accucore RP-MS, 2.1 mm × 10 mm, Guard Cartridges	Pack of 4	10039-134
Guard holder	Thermo Scientific Universal Uniguard Holder for 2.1/3.0 mm ID	Each	10046-748

Analytical GC columns for volatile PFAS

Group	Description	Quantity	Cat. no
GC Column	TraceGOLD TG-200MS Column, 30 m, 0.25 mm, 1.0 µm Recommended for neutral or volatile PFAS analysis	Each	10056-332
GC Column	TraceGOLD TG-5SiIMS Column, 30 m, 0.25 mm, 0.25 µm Recommended for suspect screening and non-targeted analysis	Each	10056-602
GC liner	Thermo Scientific™ LinerGOLD™ SSL Liner, single taper, wool, 4 × 6.5 × 78.5 mm Recommended for splitless injection on S/SL inlet	Pack of 5	10171-652
GC liner	Thermo Scientific™ LinerGOLD™ Liner, concentric baffles, 2 × 2.75 × 120 mm Recommended for splitless injection on PTV inlet for Thermo Scientific branded GC	Pack of 5	10171-700
Gas filter	Thermo Scientific™ Super Clean™ Gas Cartridge Filter Recommended for removing moisture, oxygen, and contaminants from gas sources	Each	10055-390

SSL liner can be used for Thermo Scientific, Agilent, and PE branded GC. Gas cartridge is universal. More options and accessories are available.

Other column options are available depending on your type PFAS research, talk with your Avantor representative should you have further questions.



Sample handling solutions

PFAS sample handling solutions overview

Vials, caps, well plates and mats

To assure a robust method with accurate results use an absolute halogen free polypropylene vial or well plate, especially for the analysis of low concentrations of PFAS.

Caps and septa

To assure no contamination or background corruption for PFAS analysis, we have more than one recommended closure. You are able to use one piece polypropylene snap ring caps or polypropylene 9 mm screw thread caps with a thinned penetration area. When resealing for multiple injections is important, we offer snap ring and screw thread caps with a silicone/polyimide septum for trouble free analysis.

Note: Outer sample handling packaging qualified 100% PFAS free.



SureSTART vials and caps shown are in 100 pack, but are also available in 5,000 and 10,000 pack options.

Vials and caps

Group	Description	Quantity	Cat. no
Vials	Thermo Scientific™ SureSTART™ 0.4 mL clear PP 9 mm screw microvial, conical insert, Level 1	100/pack	TS6ESV9-04PP
Vials	SureSTART 1.5 mL clear PP 9 mm screw microvial, conical insert, Level 1	100/pack	TS6ESV9-1PP
Caps	SureSTART clear PP 9 mm screw cap, integral PP membrane, Level 1	100/pack	TS6ESC9-50
Caps	SureSTART 9 mm screw blue PP cap, white silicone/red polyimide septa, Level 1	100/pack	TS6ASC9SPI1
PFAS kit: vials/closures	SureSTART kit: 1.5. mL clear PP 9 mm screw vial, blue PP cap, white silicone/red polyimide septa 1.0 mm, specification certified, Level 2	100/pack	TS6AKSV9PP-15PI
PFAS kit: vials/closures	SureSTART kit: 1.5. mL clear PP 9 mm screw vial, clear short thread PP AVCS™ cap, clear polypropylene integral PP septum thinned penetration area, specification certified, Level 1	100/pack	TS6EKSV9PP-15

Sample handling solutions (continued)

Vials and caps for U.S. EPA method 1633

Group	Description	Quantity	Cat. no
Vials	SureSTART 0.4 mL clear PP, 9 mm screw microvial, conical insert, Level 1	100/pack	TS6ESV9-04PP
Vials	SureSTART 1.5 mL clear PP, 9 mm screw microvial, conical insert, Level 1	100/pack	TS6ESV9-1PP
Caps	SureSTART clear PP, 9 mm screw caps with integral PP membrane, Level 1	100/pack	TS6ESC9-50
Caps	SureSTART 9 mm screw blue PP cap, white silicone/red polyimide septa, Level 1	100/pack	TS6ASC9SPI1
PFAS kit: vials/closures	SureSTART kit: 1.5. mL clear PP 9 mm screw vial, blue PP cap, white silicone/red polyimide septa 1.0 mm, specification certified, Level 2	100/pack	TS6AKSV9PP-15PI
PFAS kit: vials/closures	SureSTART kit: 1.5. mL clear PP 9 mm screw vial, clear short thread PP AVCS cap, clear polypropylene integral PP septum thinned penetration area, specification certified, Level 1	100/pack	TS6EKS9PP-15

Vials and caps for U.S. EPA method 533 and 537

Group	Description	Quantity	Cat. no
Vials	SureSTART 0.4 mL clear PP 9 mm screw microvial, conical insert, Level 1	100/pack	TS6ESV9-04PP
Vials	SureSTART 1.5 mL clear PP 9 mm screw microvial, conical insert, Level 1	100/pack	TS6ESV9-1PP
Caps	SureSTART clear PP 9 mm screw caps with integral PP membrane, Level 1	100/pack	TS6ESC9-50
Caps	SureSTART 9 mm screw blue PP cap, white silicone/red polyimide septa, Level 1	100/pack	TS6ESC9-50
PFAS kit: vials/closures	SureSTART kit: 1.5. mL clear PP 9 mm screw vial, blue PP cap, white silicone/red polyimide septa 1.0 mm, specification certified, Level 2	100/pack	TS6ASC9SPI1
PFAS kit: vials/closures	SureSTART kit: 1.5. mL clear PP 9 mm screw vial, clear short thread PP AVCS cap, clear polypropylene integral PP septum thinned penetration area, specification certified, Level 1	100/pack	TS6EKS9PP-15
Vials	SureSTART 0.8 mL PP snap top microvials, flat bottom, Level 1 Everyday Analysis	100/pack	TS6ERV-118PPFB
Vials	SureSTART 0.8 mL PP snap top microvials, conical insert, Level 1	100/pack	TS6ERV11-08PPC
Vials	SureSTART 0.3 mL PP snap top microvials, conical insert, Level 1	100/pack	TS6ERV11-03PPC
Caps	SureSTART blue PE 11 mm snap caps with blue PP septa, Level 1	100/pack	TS6ERC11PE
Caps	SureSTART kit: 0.8 mL blue polyethylene 11 mm snap cap, white silicone/red polyimide 1.0 mm septa, Level 2	100/pack	TS6ARC11SPI1
PFAS kit: vials/closures	SureSTART kit: 0.8 mL clear PP snap vial 11 mm flat bottom insert, blue 11 mm blue polyethylene snap cap, white silicone/red polyimide 1.0 mm septa, specification certified, Level 2	100/pack	TS6AKRV11PP-08PI
PFAS kit: vials/closures	SureSTART kit: 0.8 mL clear PP snap vial 11 mm flat bottom insert, blue 11 mm blue polyethylene integral PE septum thinned penetration area, specification certified, Level 1	100/pack	TS6EKR11PP-08





Sample handling solutions (continued)

Vials and caps for U.S. EPA method 8327

Group	Description	Quantity	Cat. no
Vials	SureSTART 0.8 mL PP snap top microvials, flat bottom, Level 1	100/pack	TS6ERV-118PPFB
Vials	SureSTART 0.8 mL PP snap top microvials, conical insert, Level 1	100/pack	TS6ERV11-08PPC
Vials	SureSTART 0.3 mL PP snap top microvials, conical insert, Level 1	100/pack	TS6ERV11-03PPC
Caps	SureSTART blue PE 11 mm snap caps with blue PP septa, Level 1	100/pack	TS6ERC11PE
Caps	SureSTART kit: 0.8 mL blue polyethylene 11 mm snap cap, white silicone/red polyimide 1.0 mm septa, Level 2	100/pack	TS6ARC11SPI1
PFAS kit: vials/closures	SureSTART kit: 0.8 mL clear PP snap vial 11 mm flat bottom insert, blue 11 mm blue polyethylene snap cap, white silicone/red polyimide 1.0 mm septa, specification certified, Level 2	100/pack	TS6AKRV11PP-08PI
PFAS kit: vials/closures	SureSTART kit: 0.8 mL clear PP snap vial 11 mm flat bottom insert, blue 11 mm blue polyethylene integral PE septum thinned penetration area, specification certified, Level 1	100/pack	TS6EKR11PP-08

Well plates and mats

Group	Description	Quantity	Cat. no
Well plate	Thermo Scientific™ WebSeal™ 96-Well Deep Well Plates, non-coated plastic, round U-bottom, 1000 µL, Level 2	50/pack	10036-216
Well plate mat	WebSeal 96-Well Plate Sealing Mat, 7 mm diameter, round, Level 2	100/pack	76047-536
Well plate	WebSeal 96-Well Deep Well Plates, non-coated plastic, square V-bottom, 2000 µL, Level 2	50/pack	10036-220
Well plate mat	WebSeal 96-Well Plate Sealing Mat, 8 mm diameter, square, Level 2	100/pack	10036-182



Complete environmental workflow solutions

For environmental scientists pursuing lower detection limits and higher robustness in their analysis, selecting the correct workflow is imperative for success in delivering high-quality results.

We strive to create a better understanding of how to compose an optimal workflow allowing scientists to improve their throughput

and provide imperative environmental data faster and more consistently. The workflows in this brochure offer a sampling from Thermo Fisher Scientific.

For more information or other workflows, please [contact us](#).

Pesticides and Per- and polyfluoroalkyl substances (PFAS)


- Thermo Scientific™ HyperSep™ SPE products
- Thermo Scientific™ Accucore™ (U)HPLC columns
- Thermo Scientific™ Hypersil GOLD™ (U)HPLC columns
- Thermo Scientific™ Acclaim™ (U)HPLC columns
- Thermo Scientific™ TraceGOLD™ GC columns
- Gas chromatography (GC) parts and accessories
- Thermo Scientific™ SureSTART™ vials and caps

Volatile organic compounds (VOC) and semi-volatile organic compounds (SVOC)

- Solid-phase extraction products
- TraceGOLD GC columns
- GC parts and accessories
- SureSTART vials and caps

Persistent organic pollutants (POP)

- Solid-phase extraction products
- Accucore (U)HPLC columns
- TraceGOLD GC columns
- GC parts and accessories
- SureSTART vials and caps



Find more information about our products at [vwr.com/thermo_chrom](https://www.vwr.com/thermo_chrom)

 **avantor**[™]
delivered by **vwr**[™]

[VWR.COM](https://www.vwr.com)

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada unless otherwise noted, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit [vwr.com](https://www.vwr.com) to view our privacy policy, trademark owners, and additional disclaimers. © 2024 Avantor, Inc. All rights reserved.

For Research Use Only. Not for use in diagnostic procedures. © 2023, 2024 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. This information is presented as an example of the capabilities of Thermo Fisher Scientific products. It is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representatives for details. **Lit. No. 200326W 0424**