

GE Healthcare

# Multiphor II Electrophoresis System

User Manual





## Important user information

---

All users must read this entire manual to fully understand the safe use of Multiphor II Electrophoresis System.

### WARNING!

---



The WARNING! sign highlights instructions that must be followed to avoid personal injury. It is important not to proceed until all stated conditions are met and clearly understood.

### CAUTION!

---

The Caution! sign highlights instructions that must be followed to avoid damage to the product or other equipment. It is important not to proceed until all stated conditions are met and clearly understood.

### Note

---

The Note sign is used to indicate information important for trouble-free and optimal use of the product.

## CE Certifying

---

This product meets the requirements of applicable CE-directives. A copy of the corresponding Declaration of Conformity is available on request.

The **CE** symbol and corresponding declaration of conformity, is valid for the instrument when it is:

- used as a stand-alone unit, or
- connected to other CE-marked GE Healthcare instruments, or
- connected to other products recommended or described in this manual, and
- used in the same state as it was delivered from GE Healthcare except for alterations described in this manual.

### WARNING!

---

This is a Class A product. In a domestic environment this product may cause radio interference in which case the user may be required to take adequate measures.

## Recycling

---



This symbol indicates that the waste of electrical and electronic equipment must not be disposed as unsorted municipal waste and must be collected separately. Please contact an authorized representative of the manufacturer for information concerning the decommissioning of equipment.



# Contents

<b>1</b>	<b>Introduction</b>	
<b>2</b>	<b>Safety information</b>	
<b>3</b>	<b>Description of parts</b>	
<b>4</b>	<b>Installation</b>	
<b>5</b>	<b>Operation</b>	
5.1	Electrophoresis using ExcelGel SDS and buffer strips .....	20
5.2	Electrophoresis using buffer chambers.....	25
5.3	Isoelectric focusing using Ampholine PAGplate.....	28
5.4	Isoelectric focusing using CleanGel IEF .....	32
5.5	2-D electrophoresis using Immobiline DryStrip and ExcelGel SDS.....	34
5.6	Electrophoretic transfer Introduction .....	35
5.7	Stock solutions.....	42
5.8	Running conditions for precast gels.....	43
<b>6</b>	<b>Maintenance</b>	
6.1	Recycling .....	47
<b>7</b>	<b>Technical specifications</b>	
<b>8</b>	<b>Trouble shooting</b>	
<b>9</b>	<b>Multiphor II application kits and accessories</b>	
9.2	Buffer Strip Positioner.....	59
9.3	Immobiline DryStrip Reswelling Try.....	60
9.4	Immobiline DryStrip Kit.....	60
9.5	NovaBlot Kit.....	61
9.6	FilmRemover .....	62
9.7	Roller .....	63
9.8	Sample application accessories.....	63
<b>10</b>	<b>Ordering information</b>	
10.1	Multiphor II.....	67
10.2	MultiTemp III.....	69
10.3	EPS Power Supplies .....	69
10.4	Precast gels and buffer strips.....	70
10.5	Molecular weight and pI markers.....	71
10.6	Carrier ampholytes.....	71
10.7	PlusOne electrophoresis chemicals.....	73



# 1 Introduction

---



Multiphor™ II electrophoresis system is a versatile modular system for horizontal electrophoresis, isoelectric focusing, 2-D electrophoresis and electrophoretic transfer.

For ease of use and reproducible results, an innovative range of precast gels for all major electrophoretic techniques is available with Multiphor II:

Technique	Precast Gel
SDS and Native PAGE	ExcelGel™ SDS Gradient ExcelGel SDS Homogeneous
IEF	Ampholine™ PAGplate CleanGel™ IEF
2-D electrophoresis	Immobiline™ DryStrip ExcelGel SDS gradient Homogeneous

If laboratory cast gels are preferred an application kit and accessories can be added to the basic electrophoresis unit.

The following guide summarizes how you can expand and use Multiphor II Electrophoresis Unit with application kits and accessories.

<b>Application</b>	<b>Recommended Kit/Accessory</b>	<b>Code No.</b>
SDS and Native PAGE 0.5 x 125 x 260 mm homogeneous and gradient gel	SDS and Native PAGE, IEF Kit Gradient Maker	18-1102-45 18-1013-72
ExcelGel SDS	Buffer Strip Positioner	80-6442-90
IEF in polyacrylamide 0.5 x 125 x 260 mm	SDS and Native PAGE, IEF Kit	18-1102-45
2-D, first dimension: Immobiline DryStrip	Immobiline DryStrip Kit Reswelling Tray, 7-18 cm Reswelling Tray, 7-24 cm	18-1004-30 80-6371-84 80-6465-32
Electrophoretic transfer	NovaBlot Kit FilmRemover	18-1016-86 18-1013-75

This User Manual is comprised of the following sections:

1. "Introduction" includes a general description of Multiphor II system and dedicated precast gels, a guide to the application kits and the manual structure.
2. "Description of parts" describes in detail the components of Multiphor II Electrophoresis Unit.
3. "Installation" contains a detailed description of how to install Multiphor II Electrophoresis Unit and Multiphor II NovaBlot™ Unit.
4. "Operation" contains information on the operating procedure for SDS and native polyacrylamide gel electrophoresis, isoelectric focusing, 2-D electrophoresis and electrophoretic transfer.
5. "Maintenance" gives cleaning recommendations to help you maintain your Multiphor II unit.
6. "Trouble shooting" offers suggestions for correcting problems that may occur.
7. "Multiphor II application kits and accessories" describes in detail the contents, assembly and use of each Multiphor II application kit.
8. "Ordering information."
  - Multiphor II Electrophoresis Unit, application kits, accessories and replacement parts
  - MultiTemp™ III, EPS and GPS Power Supplies
  - Precast gels and buffer strips
  - Molecular weight and pI markers
  - Ampholine and Pharmalyte™ carrier ampholytes
  - Gel casting and electrophoresis chemicals



## 2 Safety information

---

To avoid any risk of injury, the instrument should be operated only by properly trained personnel and always in accordance with the instruction provided. Read this entire manual before using the instrument.



**WARNING!** The instrument is designed for indoor use only.



**WARNING!** Do not operate the system in extreme humidity (above 95% RH). Avoid condensing by equilibrating to ambient temperature, when taking the unit from a cooler to a warmer environment.



**WARNING!** Always check the wires for damage before using the unit.



**WARNING!** Always check that the electrodes are properly connected before using the lid.



**WARNING!** Always connect the lid according to the mounting instruction.



**WARNING!** Always connect the cables to the Power supply BEFORE turning the Power Supply ON.



**WARNING!** Always TURN OFF the Power Supply before removing the lid.



**WARNING!** Do NOT use concentrated acids, bases or halogenated and aromatic hydrocarbons.



**WARNING!** Only use water or coolant with high electrical resistance in the cooling plate.



**WARNING!** NEVER EXCEED the maximum pressure 0.5 bar in the cooling plate.



**WARNING!** NEVER EXCEED the maximum allowed voltage, current or power.



**WARNING!** The cooling plate is rated for operation at up to 3.5 kV (p-p).



**WARNING!** When using hazardous chemicals, take all suitable protective measures, such as wearing protective glasses and gloves resistant to the chemicals used. Follow local regulations and instructions for safe operation and maintenance of the system.



**WARNING!** When using hazardous chemicals, make sure that the entire system has been flushed thoroughly with bacteriostatic solution, e.g. NaOH, and distilled water before service and maintenance.

### 3 Description of parts

---

The Multiphor II Electrophoresis Unit includes; buffer tank with 4 levelling feet, cooling plate with accessories, safety lid, electrode holder with movable EPH/IEF electrodes (for buffer strips and electrode strips).

#### Unit contents – Code No. 18-1018-06

Designation	Code No
Buffer Tank	18-1122-25
Levelling Foot (4/pkg)	18-1026-40
Cooling Plate ceramic, 210 x 270 mm	18-1103-46
Grommet (2/pkg)	80-1106-58
Cooling Tubing, 8/12 mm, 4 m	80-1106-56
Tubing Connector Set 2 pcs female, 2 pcs male	18-1104-26
Hose Clamp (10/pkg)	18-1104-27
Safety Lid	18-1122-26
Electrode Holder	80-1106-55
EPH/IEF Electrode cathode	80-1121-52
EPH/IEF Electrode anode	80-1121-53
User Manual and Application Package	18-1103-44



The buffer tank is made of polypropylene, which is resistant to nearly all chemicals at room temperature.


The buffer tank contains four pin contacts. Viewing the buffer tank from the front, the cathode pins are located to the left and the anode pins to the right.

The larger pins are for connection to the safety lid and complete the electrical circuit when the lid is in position.

### 3 Description of parts

The small left hand pin is used to connect the EPH/IEF or card-mounted cathode electrodes. The small pin on the right connects to the card-mounted anode or the EPH/IEF cathode via the red lead mounted on the unit.


The buffer tank holds the four adjustable levelling feet, supports the cooling plate and is covered during electrophoresis with the safety lid. The buffer tank includes four buffer chambers, each with a 1 liter capacity, allowing the user to choose one of two orientations for electrophoretic runs.



The safety lid contains electrode leads, apertures for voltage measurement, and a safety interlock.

The well-recessed cathode connector (black) and anode connector (red) for connection to power supplies ensure safe operation at high voltages. The polycarbonate lid snugly fits the contours of the buffer tank. This makes it possible to reduce the atmospheric CO<sub>2</sub> content around the gel (important for IEF at basic pH intervals) and provides increased protection against condensation.

**Note:** *Polycarbonate is not resistant to concentrated acids and bases, or to halogenated and aromatic hydrocarbons.*



The ceramic (aluminium oxide) cooling plate measures 210 x 270 mm, supports the gel, and provides uniform temperature control. Aluminium oxide is an excellent heat conductor and electrical insulating material.



**WARNING!** The cooling plate is rated for operation at up to 3.5 kV (p-p).

To facilitate the correct positioning of electrophoresis gels, the surface of the cooling plate is screened with a template measuring 190 x 250 mm.


The two grommets are connected to the inlet and outlet tubes of the cooling plate which can then be connected to a thermostatic circulator such as MultiTemp III.



**WARNING!** Only use water with high electrical resistance as coolant and NEVER EXCEED maximum pressure of 0.5 bar.



The cooling tubing, tubing connector set and hose clamps provide a flexible and safe way to connect Multiphor II to a thermostatic circulator such as MultiTemp III.



The electrode holder holds the movable EPH/IEF electrodes. The holder keeps the electrodes away from the gel surface during alignment and then provides a uniform pressure over the buffer strips or electrode strips during the separation.

### 3 Description of parts

The electrode holder consists of a double strength glass plate with ground edges and four corner feet made of Rynite™ FR530. The electrode holder holds one anode and one cathode electrode.

The EPH/IEF electrodes consist of moulded polysulfone bars which support the platinum wire, held taut by stainless steel springs. The cables are spring reinforced for safety. The anode cable (red) carries the pin contact to be connected to the socket connector on the buffer tank.

The cathode cable (black) carries a female socket connector which fits to the buffer tank pin connector.


Clamping nuts located at each end of the electrode allow easy adjustment of the electrodes on the holder. The distance between the electrodes can be varied from 10 mm to 240 mm.

**Note:** *Polysulfone is not resistant to ketones, esters, halogenated and aromatic hydrocarbons.*

## 4 Installation

---

For all products, check the unassembled parts against the Packing List for the respective product to ensure that all items have been included.




For easy connection of Multiphor II Electrophoresis Unit to the cooling device, install tubing connectors on both inlet and outlet tubings. Use one male and one female connector on the Multiphor II unit side and thermostat unit side respectively. The cooling plate tubings and thermostatic circulator tubings can then be locked separately.




Two Multiphor II units can also be connected in series to one thermostatic circulator.

Place the rubber grommets on the cooling plate inlet and outlet tubes. Slide a short piece of tubing onto each tube and secure with hose clamps. Attach the male and female tubing connectors as described above. Repeat this process with the thermostatic circulator using longer pieces of tubing.

To lock the connectors, insert the male connector into the female and turn clockwise one-quarter turn until it clicks.






Screw one levelling foot into each corner of the buffer tank. Place the buffer tank on the lab bench where it will be used. Place the cooling plate on the unit, using the moulded guides to position it correctly. Fit the grommets into the cutouts in the back of the unit. Place a spirit level on the cooling plate and adjust the levelling feet until the unit is levelled.



To mount the EPH/IEF electrodes on the electrode holder, unscrew the clamping nut from each electrode.

When running electrophoresis across the width of the cooling plate, mount the electrodes as illustrated. Place the electrode under the electrode holder. Replace the nut and lightly tighten until the electrode is held firmly in place. To align the electrodes with the electrode strips, place the electrode holder with the electrodes onto the false holes on the buffer tank using the four corner feet.



When running electrophoresis along the length of the cooling plate, mount the electrodes as illustrated. Place the electrodes under the electrode holder.



Replace the nut and lightly tighten until the electrode is held firmly in place. To align the electrodes with the electrode strips, place the electrode holder with the electrodes onto the false holes on the buffer tank using the four corner feet.



Install the safety lid.



## 5 Operation



**WARNING!** When using hazardous chemicals, take all suitable protective measures, such as wearing protective glasses and gloves resistant to the chemicals used. Follow local regulations and instructions for safe operation and maintenance of the system.


This section, together with the information supplied with the precast gels, gives all the necessary information to run most analytical electrophoresis techniques using our precast gels. Running procedures for electrophoretic transfer are also included.

For laboratory cast gels, use the running conditions recommended in *Electrophoresis in Practice, A Guide to Theory and Practice* by Reiner Westermeier.


Multiphor II contains two alternative electrode configurations.

- EPH/IEF electrodes for use with buffer strips or electrode strips
- EPH electrodes for use with electrode wicks and buffer chambers

IEF, SDS-PAGE and native PAGE are most conveniently performed using EPH/IEF electrodes and buffer strips. The strips are applied on the gel edges with the electrodes on top.



ExcelGel



Electrode wicks

The buffer chambers are located below the cooling plate with the electrodes immersed in buffer solution. Paper wicks connect the buffer solution with the gel.

This method is used for SDS-PAGE and native PAGE.

The optional card-mounted EPH electrodes (18-1122-19 and 18-1122-20) – for electrophoresis using buffer chambers across the width of the cooling plate – are moulded from polypropylene and support the platinum wire. The anode cable (red) and cathode cable (black) carry female pin connectors for attachment to the male pins at the front of the buffer tank.

### 5.1 Electrophoresis using ExcelGel SDS and buffer strips

This section describes the running procedure for SDS PAGE using buffer strips. The running of ExcelGel SDS, gradient 8-18 using ExcelGel SDS buffer strips is chosen as an example, but the basic method is applicable to all SDS PAGE and native PAGE gels.

ExcelGel SDS, gradient 8-18, is a 0.5 mm-thin, precast polyacrylamide gel for horizontal electrophoresis of SDS denatured proteins. To facilitate handling, the gel is cast on a plastic support. During the run, the precast SDS buffer strips supply the gel with buffer ions. For further information, see the information supplied with ExcelGel SDS gels.




#### Sample preparation

Dissolve the samples in sample buffer B (for recipes, see Section 4.9 Stock solutions). Then heat the sample solution at 95 °C for 3 minutes. The sensitivity of your development technique and the volume of sample applied to the gel will determine the lower limit of your sample concentration. Generally, the sample must contain 200 to 500 ng of each component for Coomassie staining, and at least 10-25 ng of each component for silver staining. For molecular weight determination, we recommend the use of molecular weight calibration kits LMW and HMW/SDS.

#### Sample application

In horizontal electrophoresis there are three methods of applying the sample: application strips, paper pieces and sample wells. Sample application strips are put on the gel surface, forming sample slots. Silicone rubber sample application strips are specially designed for easy sample application.



The following application strips are available:

SDS application strips for up to 40 µl of sample in 26 slots. IEF/SDS application strips for up to 20 µl of sample in 52 slots.

Immobiline applicator strip for up to 5  $\mu\text{l}$  of sample in 52 slots. Immobiline applicator strip is designed to counteract lateral band spreading.

Sample application pieces hold approximately 20  $\mu\text{l}$  of sample. For smaller volumes, cut the paper pieces to an appropriate size. At least 24 application pieces size 5 x 10 mm and 50 application pieces size 2.5 x 5 mm can be placed on one gel. Apply the sample about 1 cm away from the cathodic buffer strip and 1 cm away from each short side of the gel. For the best results, remove the application pieces 15 min after electrophoresis has started.

ExcelGel SDS, ExcelGel Homogeneous 7.5, 12.5, 15 and CleanGel are available with various numbers of sample wells for various volumes. Samples are applied directly into the wells immediately prior to the run.

### Electrophoresis



Connect Multiphor II to MultiTemp II thermostatic circulator. Switch on MultiTemp II 15 minutes before starting the experiment and set the temperature to 15 °C.

Always wear clean gloves when working with polyacrylamide gels and buffer strips, particularly when using sensitive staining methods.





Remove one ExcelGel SDS from the package. Pipette 1 ml of insulating fluid (kerosene or light paraffin oil) onto the cooling plate. Place the gel with the stiff plastic film facing down in the middle of the cooling plate, making sure no air bubbles are trapped under the gel. Position the gel so that the polarity of the gel corresponds to that of the plate. Use the screened template on the cooling plate to centre the gel. Remove any excess solution with paper tissue. Remove the protective cover film from the gel.

Open the ExcelGel SDS strip packages and apply the cathodic and the anodic SDS buffer strips on the respective sides of the gel.



**Note:** *The narrowest side of the buffer strip should be placed on the gel surface.*

Choose an appropriate sample application method and apply the sample.




Place the electrode holder on the electrophoresis unit in the shallow depressions.



Align the EPH/IEF electrodes with the centre of the buffer strips by loosening the clamping nuts and sliding the electrodes to the appropriate position. Retighten the clamping nuts. Lift the electrode holder slightly and reposition the supporting feet over the deep holes. Lower carefully, so that the electrodes rest on the buffer strips. Connect the electrodes to the buffer tank.

During electrophoresis, the socket on the bridging cable **MUST** be attached to the pin connector at the front of the buffer tank as shown in the picture.

Connect the two electrodes to the buffer tank using the spring-loaded cables on the electrodes.



Connect the socket of the cathode electrode to the pin at the front of the unit and the anode pin to one of the sockets at the back.



Place the safety lid in position by matching the extensions on the back of the lid with the openings on the base unit. Using the extensions at the back as a hinge, connect the male and female banana plugs by pressing down firmly on the front of the lid.

Connect Multiphor II to the power supply. Follow the recommended electrical settings and running times given in the instructions supplied with the precast gel.

#### Running conditions

	Voltage (V)	Current (mA)	Power (W)	Time (min)
Run	600	50	30	75*

\* Approximate time, or until the Bromophenol Blue front reaches the anode buffer strip.

When the Bromophenol Blue front has reached the anodic buffer strip, electrophoresis is complete and should be stopped.

#### Ending the run

Turn off the power supply. Disconnect the Multiphor II unit from the power supply. Remove the safety lid from the unit. Carefully remove the electrode holder. Gently pull the strips from the gel and continue with detection techniques as required.



**WARNING!** Always TURN OFF the power supply before opening the safety lid.

### Detection

For automated silver and Coomassie™ staining of polyacrylamida gels see Protocol guide, Hoefer™ Automated Gel Stainer (80-6343-34).

#### *Coomassie staining*

On completion of the electrophoresis, immediately immerse the gel in Staining Tray 1 using the solutions and times indicated in table below.



The staining and destaining steps should be carried out on a shaking table. (See section 4.9 for stock solutions). Each step requires 250 ml of solution.

Step No	Solution	Time (min)	Temp °C
1	Fixing solution C	20	23
2	Destaining solution I	2	23
3	Staining solution K	10	60
4	Destaining solution I	20	23
5	Destaining solution I	30	23
6	Preserving solution L	10	23

The staining solution should be heated to 60 °C and poured over the gel. No further heating is necessary. Destain the gel using several changes of destaining solution (I) until the background is clear. Change the solution frequently (particularly at the beginning) in order to speed up the destaining. To preserve the gel, soak a cellophane sheet in preserving solution (L). Place it on the gel surface. Remove any air bubbles and wrap the excess cellophane around the glass plate. An additional glass plate may be placed underneath during drying to stop the cellophane from shifting or wrinkling. Leave the gel at room temperature until it is completely dry.

#### **Silver staining**

Silver staining is performed essentially as described by J. Heukeshoven and R. Dernick, *Electrophoresis Forum* 1986, 22-27. On completion of the electrophoresis, immediately immerse the gel in Staining Tray 1 using the solutions and times indicated in the table below. All steps should be carried out at room temperature in daylight, while gently shaking the solution. Use 250 ml of solution for each step. (See section 4.9 for stock solutions).



Time schedule for silver staining

Step No.	Solution	Time (min)
1	Fixing solution C	30
2	Incubation solution D	30
3	Distilled water	3 x 5
4	Silver solution E	40
5	Developing solution F	5-15 *
6	Stop solution G	10
7	Distilled water	3 x 5
8	Preserving solution L	20

\* Short development times will give a lightly stained gel. Long development times will give a dark gel.


To preserve the gel, soak a cellophane sheet in preserving solution (L) and lay it on the gel surface. Remove any air bubbles and wrap the excess around the glass plate. An additional glass plate may be placed underneath during drying to stop the cellophane from shifting or wrinkling. Leave the gel at room temperature until it is completely dry.

## 5.2 Electrophoresis using buffer chambers

This section describes the running procedure when using buffer chambers. To place the electrodes in the buffer chamber, remove the cooling plate from the buffer tank.

To reduce the effect of electrolysis products during the electrophoresis, the electrodes should be positioned as far as possible from the wicks.

Therefore, when performing electrophoresis across the large cooling plate, the electrodes should be placed in the grooves in the wall closest to the centre of the unit. The wicks lie at the outer edge of the buffer chamber. Place the cathode electrophoresis electrode in the left buffer chamber and the anode in the right buffer chamber.



Fill each chamber to the moulded line (indicating 1 liter volume) with buffer solution. (When running 120 x 250 mm gels, pour 1.2 liters of buffer into each chamber to ensure adequate buffer contact with the wicks.) Replace the cooling plate, making sure that the electrode socket connectors lie to the front and that the connecting cable is clear of the feet on the plate. Disconnect the anode bridging connector for isoelectric focusing and connect the electrodes to their respective pins.


When performing immunoelectrophoresis or agarose electrophoresis, center a small (84 x 94 mm) glass plate on the dry cooling plate and attach a strip of tape along the width of the cooling plate in alignment with the edges of the glass plate. These stop the small gels from shifting during application of the wicks and ensure that they are centred between the two buffer chambers during electrophoresis.

Switch on MultiTemp III thermostatic circulator and set the desired temperature (normally 10 °C for PAGE or agarose electrophoresis and 15 °C for SDS-PAGE) 15 minutes before starting the experiment. To ensure efficient heat transfer from the gel during electrophoresis, a uniform layer of a non-charged insulating fluid is applied under the gel. To do this, pour a few milliliters of kerosene or light paraffin oil towards one end of the cooling plate.

Starting with one end of the gel support, gradually lower the gel to the horizontal position, constantly checking for trapped air bubbles. If air becomes trapped, raise the gel just enough to release the air and then continue to lower it onto the cooling plate. Use the template markings to centre the gel on the cooling plate. Remove any excess solution with a tissue. Repeat this procedure for each gel. If voltage probe measurements are required, the gel(s) must be positioned with the direction of the current path across the width of the cooling plate.




Prepare the electrode wicks by aligning 8-10 pieces of filter paper for each buffer chamber. Starting at one end, slowly immerse the electrode wicks in the buffer, using capillary action to reduce the amount of air trapped in the paper.



When running large (195 x 250 mm) or medium sized gels (120 x 250 mm), place the wicks so that the long edge overlaps the gel by 15 mm over the entire length. The template markings are useful for checking that the alignment is correct.

For running small gels (84x94 mm) on the large cooling plate, place the wicks with the short edge overlapping the gel by 15 mm. In this case, one set of wicks is required for each gel.

Apply the samples as required. If voltage probe measurements are not required, remove the isoelectric focusing electrodes from the electrode holder.



Position the empty electrode holder so it lies directly on the electrode wicks. This will ensure even contact between the wicks and the gel and stop any moisture from condensing on the gel surface.



Replace the safety lid on the unit and connect the Multiphor II unit to the power supply. Set the power requirements, and start the experiment.

Typical power settings for agarose immunoelectrophoresis using buffer chambers are: constant voltage 20 V/cm and current and power set to maximum. Run time is 40-60 minutes. Bromophenol Blue is used as a tracking dye.

Typical power settings for SDS PAGE electrophoresis using buffer chambers.

Separation distance (cm)	Voltage (V)	Current (mA)	Power (W)	Time (min)	Temp (°C)
8	600	50	30	100	15
16	1200	50	30	165	15

### Ending the run

Turn off the power supply. Disconnect the Multiphor II unit from the power supply. Take off the safety lid from the unit. Carefully remove the electrode holder. Gently pull the wicks from the surface of the gel and continue with detection techniques as required.



**WARNING!** Always TURN OFF the power supply before opening the safety lid. Although user safety is not endangered, arcing may damage the contacts.

### 5.3 Isoelectric focusing using Ampholine PAGplate

This section describes the running procedure for IEF using Ampholine PAGplate. For further information, see the instructions supplied with the product.



Switch on MultiTemp II thermostatic circulator and set the temperature (normally 10 °C for polyacrylamide IEF) 15 minutes before starting the experiment.

Cover the holes in the safety lid with tape to limit the amount of CO<sub>2</sub> in contact with the gel, thereby improving the basic region of the pH gradient. If desired, 100 ml of 1 M sodium hydroxide solution may be poured into the buffer chambers to absorb CO<sub>2</sub> and further improve gradient stability.

Make up about 100 mls of the required electrode solutions. (See Table).

Electrode solutions for IEF using Ampholine PAGplate

pH range	Anode solution	Cathode solution
3.5–9.5	1 M Phosphoric acid	1.0 M NaOH
4.0–6.5	0.1 M Glutamic acid	0.1 M β-alanine
5.5–8.5	0.4 M HEPES	0.1 M NaOH
4.0–5.0	1 M Phosphoric acid	1.0 M Glycine
5.0–6.5	0.01 M Acetic acid	0.01 M NaOH

**Note:** Wear clean gloves when working with polyacrylamide gels.

Pipette 1 ml of insulating fluid (kerosene or light paraffin oil) onto the cooling plate. Open one package of Ampholine PAGplate and position the gel with the stiff plastic film facing down on the cooling plate. Make sure no air bubbles are trapped under the gel. Use the screened template on the cooling plate to centre the gel. Remove any excess solution with a paper tissue.

Soak the electrode strips evenly in the appropriate electrode solution (see table above), approximately 3 ml/strip. The surface of the strips should look wet. Remove excess solution with a tissue paper.




Apply the electrode strips to the long edges of the gel. Make sure the electrode strips are applied with the correct polarity. Use sharp scissors to cut off the strips which protrude beyond the ends of the gel.

### Sample application

Some proteins may be pH sensitive (they may precipitate at specific pH values) and better separation may be obtained when the sample is applied at a specific position in the pH gradient. This can be tested by applying the sample at various positions across the pH gradient.

There are three different methods for sample application. The method of choice depends primarily by the sample and the volume to be applied.

1. IEF/SDS applicator strip for 5–20  $\mu\text{l}$  sample volumes. This applicator strip makes sample loading quick and simple, especially when a large number of samples are to be applied. Check that the contact between the gel and the applicator strip is uniform. Leave the applicator strip on the gel during focusing.





2. Sample application pieces. Place the dry pieces on the Ampholine PAGplate surface at the desired position(s) in the gradient. Using a micropipette, apply 15–20  $\mu\text{l}$  volumes of sample solution on each piece. To apply larger volumes, use 2 or 3 pieces, stacked or placed end to end, for each sample applied. If smaller volumes are used, trim the paper proportionally before applying it to the gel. Remove the pieces after completing approximately half the total focusing time.
3. Very small sample volumes (e.g. 2  $\mu\text{l}$ ) can be applied as droplets directly onto the dry gel surface.

To determine the pH gradient, prepare and apply pI markers according to the instruction sheet included in the pI Marker Kit.


### Starting the IEF run

Place the electrode holder in the shallow depressions on the Multiphor II unit.



Align the electrodes with the centre of the electrode strips by loosening the clamping nuts and sliding the electrodes to the appropriate position. Retighten the clamping nuts.

Lift the electrode holder slightly and reposition the supporting feet over the deep holes. Lower carefully so that the electrodes rest on the electrode strips. As the anodic pin and socket connectors are different for electrophoresis using buffer chambers and isoelectric focusing, a red bridging cable has been provided.



During isoelectric focusing, the socket on the bridging cable MUST be attached to the pin connector.

Connect the two electrodes to the buffer tank using the spring-loaded cables on the electrodes.



Connect the socket of the cathode electrode to the pin at the front of the unit and the anode electrode pin to the socket at the back.



Place the safety lid in position by matching the extensions on the back of the lid with the openings on the base unit. Using the extensions at the back as a hinge, connect the male and female banana plugs by pressing down firmly on the front of the lid.

### Running conditions

Connect Multiphor II to the power supply. Set the running conditions given in the table below. Start the isoelectric focusing by turning on the power supply. Observe the time limits closely. If the isoelectric focusing is run too long, the pH gradient will begin to drift towards the cathode.

Suggested running conditions for Ampholine PAGplate

pH range	Voltage (V)	Current (mA)	Power (W)	Time (h)
3.5–9–5	1 500	50	30	1.5
4.0–6.5	2 000	25	25	2.5
5.5–8.5	1 600	50	25	2.5
4.0–5.0	1 400	50	30	3.0
5.0–6.5	2 000	15	20	3.0

**Note:** *If only half a gel is used, halve the current and power settings.*

Remove the sample application pieces with forceps after approximately half the focusing time has expired. When a sample application strip is used, let it remain on the gel during the focusing, but remove it before placing the gel in the fixing solution.

### Ending the run

Turn off the power supply and disconnect Multiphor II from the power supply. Remove the safety lid and electrode holder from the unit. Carefully remove the gel and proceed with staining or preparation for the second dimension as required.



**WARNING!** Always turn OFF the power supply before opening the safety lid



### Detection methods

For automated silver and coomassie staining of polyacrylamida gels see Protocol guide, Hoefer Automated Gel Stainer (80-6343-34).

#### *PhastGel Blue staining*

On completion of isoelectric focusing, remove the electrode strips from the gel. Immerse the gel in 250 ml of solution according to the schedule below (See section 4.9 for stock solutions). The staining and destaining steps should be carried out on a shaking table.

Step No.	Solution	Time (min)	Temp (°C)
1	Fixing solution N	20	23
2	Destaining solution I	2	23
3	Staining solution K	10	60
4	Destaining solution I	20	23
5	Destaining solution I	30	23
6	Preserving solution O	10	23

The staining solution should be heated to 60 °C and poured over the gel. No further heating is necessary. Destain the gel using several changes of destaining solution (I) until the background is clear. Change the solution frequently (particularly at the beginning) in order to speed up the destaining. To preserve the gel, soak a cellophane sheet in preserving solution (I) and lay it on the gel surface. Remove any air bubbles and wrap the excess cellophane around the glass plate. An additional glass plate may be placed underneath during drying to stop the cellophane from shifting or wrinkling. Leave the gel at room temperature until it is completely dry.


### 5.4 Isoelectric focusing using CleanGel IEF

CleanGel IEF is a washed and dried polyacrylamide gel optimized for analytical isoelectric focusing.



Prior to use, the dried gel is rehydrated in a solution containing Ampholine or/and Pharmalyte. Additives such as urea and/or detergents can also be added at this stage. CleanGel IEF is rehydrated in a flat tray (GelPool) to a thickness of 0.43 mm. Instructions are supplied with the product.





### Sample application

Use the same methods as for Ampholine PAGplate (see Section 4.4).

### Running conditions

In principle, CleanGel is run in the same way as Ampholine PAGplate except no electrode strips are used. The electrodes rest directly on the gel surface.

Recommended running conditions for one CleanGel IEF 3–10

Phase	Voltage (V)	Current (mA)	Power (W)	Time (min)
Prefocusing	700	12	8	20
Sample entrance	500	8	8	20
Isoelectric focusing	2000	14	14	90
Band sharpening	2500	14	18	10

**Note:** *If only half a gel is used, halve the current and power settings.*


### Detection Methods

Use the detection methods described for Ampholine PAGplate (see Section 4.4).

### **5.5 2-D electrophoresis using Immobiline DryStrip and ExcelGel SDS**

This chapter gives a brief description of the 2-D (two dimensional) electrophoresis method using Immobiline DryStrip and ExcelGel. For more detailed information on running conditions, please refer to the instructions supplied with Immobiline DryStrip Kit.

Isoelectric focusing using Immobiline DryStrip makes true isoelectric focusing possible and significantly improves the reproducibility of the spot distribution along the pH gradient axis of 2-D maps. Immobiline DryStrip also makes it possible to obtain distinct protein spots, even of basic proteins.



The Immobiline DryStrip Kit facilitates sample application, running and equilibration of Immobiline DryStrip for the first dimension of 2-D electrophoresis. The kit includes the accessories necessary to run up to 12 Immobiline DryStrip strips simultaneously on Multiphor II. Sample cup loading allows the application of up to 100  $\mu$ l on each Immobiline DryStrip.

See application Note 80-1443-47 for unmultiple miniformat 2-D lectrophoresis using ExcelGel 2-D, 12.5.

## 5.6 Electrophoretic transfer Introduction

The method of horizontal semi-dry electrophoretic transfer gives fast, even and efficient transfer of proteins from a gel to an immobilizing membrane. The resulting membranes may be used for a wide range of applications including general protein staining, identification of specific antigens or antibodies (immunoblotting) and glycoprotein detection. By using different electrode solutions and running conditions, it is possible to transfer proteins from SDS PAGE, native PAGE, agarose gels and isoelectric focusing gels.

The speed and efficiency of the electrophoretic transfer using NovaBlot system is dependant on:

- Characteristics of the immobilizing membrane
- Characteristics of the transfer buffer
- Molecular weight and charge of the protein
- Gel thickness and concentration of acrylamide and bisacrylamide
- Voltage, current and transfer time

The semi-dry transfer technique uses filter papers soaked in buffer as the only buffer reservoir. Both discontinuous and continuous buffer systems can be used in the filter paper layers. Methanol in the buffer solution prevents swelling of polyacrylamide gels. However, it may have the disadvantage of denaturing or fixing the proteins in the gel, resulting in a lower transfer efficiency. Conversely, methanol may increase the protein binding capacity of the nitrocellulose membrane by strengthening the hydrophobic interactions. The transfer speed and efficiency can also be increased by increasing the protein charge, i.e. adding 0.05% SDS in the transfer buffer.

The transfer is normally finished in about one hour. If it is necessary to transfer for more than 1 hour due to unusual sample characteristics, rewetting of the cathode filter paper is recommended. No cooling is necessary since negligible heat is produced during the transfer.

### Immobilizing membranes

The immobilizing membrane is an important factor affecting the efficiency of the electrophoretic transfer. The most important requirements for an immobilizing membrane are:

- High binding capacity for the molecules of interest
- Preservation of the biologic activity of the molecules of interest
- No interference with subsequent detection methods
- Chemical and mechanical stability to assay conditions
- Provision of an accurate reflection of the original separation

Nitrocellulose membranes are the standard medium for electrophoretic transfer of proteins and nucleic acids. This is due to their high binding capacity, versatility and easy use. Nitrocellulose membranes are available in various pore sizes, 0.45  $\mu\text{m}$  is most commonly used, however low molecular proteins may be lost. By using pore sizes of 0.2 or 0.1  $\mu\text{m}$ , most proteins are retained.




Nitrocellulose membranes can be probed several times. The membranes require no activation and the functional groups have an extended lifetime. Protein patterns on nitrocellulose membranes can be easily detected with most conventional stains, as well as by autoradiographic, immunoenzymatic and fluorescent methods.

Other membranes are: nylon-based membranes, diazobenzoyloxymethyl (DBM) and diazophenylthioether (DPT) papers.

**Transferring proteins from SDS polyacrylamide gels to nitrocellulose membrane.**

The support (or backing) film must be removed from all polyacrylamide and agarose gels before electrophoretic transfer. Using FilmRemover, the film is removed quickly and cleanly. Instructions for use are supplied with FilmRemover.



1. Saturate the graphite anode plate with distilled water and remove the excess water with absorbent paper. With the electrode lead to the front of the instrument, fit the lower anode plate onto the buffer tank.




Connect the anode socket (red lead) to the anode pin on the right side of the buffer tank. The transfer sandwich can now be assembled on the anode electrode.

**Note:** When assembling the transfer sandwich in NovaBlot unit, always wear gloves.



2. To ensure that the current passes only through the gel, cut the filter papers and the immobilizing and dialysis membranes to the same size as the gel to be transferred.  
When using a discontinuous buffer system, carefully soak the first layer of six filter papers in anode solution R (see Section 4.9 Stock solutions) by slowly immersing the papers under the surface of the electrode solution. Allow them to become wet by capillary action and avoid trapping air bubbles that may interfere with the transfer.



3. Carefully place the filter papers onto the anode electrode. When forming the first transfer sandwich, soak a further layer of three filter papers in anode solution S (see Section 4.9, Stock solutions), using the same method as above. Place them on top of the six filter papers on the anode electrode plate, again taking care to avoid trapping air bubbles. When using a continuous buffer system, all filter papers, cathodic and anodic are wetted in the same solution.

**Note:** To obtain optimal transfer of molecules from the gel, care should be taken to avoid trapping air bubbles at all stages of the assembly of the transfer sandwiches.




4. Cut the gel loose from the support film using FilmRemover. Do not move the gel, leave it on FilmRemover. Wet the immobilizing membrane in electrode solution S and carefully place it on top of the gel on the FilmRemover.


**Note:** Wear gloves to avoid contamination of the membrane.



- Loosen the support film from FilmRemover by pressing the handle and carefully lift the whole sandwich with the support film, immobilizing membrane and gel. Turn it over (support film up, immobilizing membrane down) and place it on the layer of three filter papers on the anode. Carefully remove the support film. If air bubbles become trapped under the gel, wet the surface of the gel with a few drops of electrode solution, and gently push out the bubbles.




- Immerse nine filter papers in cathode solution T. Place these filter papers on top of the gel to complete the transfer sandwich.




- Several gels of the same type and size can be transferred simultaneously. Two transfer sandwiches can be put on top of each other. The cellophane dialysis membrane placed between each transfer sandwich prevents crosscontamination between transfer sandwiches.

The maximum gel size is 200 x 250 mm. If small gels (125 x 250 mm) are to be transferred, NovaBlot will accept up to four gels for simultaneous transfer by assembling two transfer sandwich stacks side by side.

To ensure that the current passes through the gel, all components of the transfer sandwich are cut to the same size as the gels to be transferred.




8. Saturate the cathode electrode plate with distilled water and remove any excess with absorbent paper. Place the cathode on top of the transfer sandwich and connect the socket on the black cathode lead to the cathode pin in the Multiphor II base.



9. Close the Multiphor II safety lid and connect the unit to the power supply. It is recommended to run the transfer at a constant current of  $0.8 \text{ mA/cm}^2$ . A transfer time of approximately 1 hour is normal.

**Note:** The current is calculated using the surface area (total length  $\times$  width) of the transfer sandwiches, and this calculation applies irrespective of the number of transfer sandwiches in the stack.

**Note:** For transfer times longer than one hour turn off the power supply, remove the safety lid and carefully lift the cathode (top) electrode without disturbing the filter papers or gel. Carefully pour on additional transfer buffer to re-wetting the filter paper.





10. When the transfer is complete, turn off the power supply and disconnect NovaBlot from the power supply. Remove the safety lid and the upper cathode electrode. Carefully disassemble the transfer sandwiches and remove the immobilizing membranes for analysis. If necessary, save and stain the gel to monitor the transfer efficiency.

Clean the electrodes with distilled water.

**Note:** *Always turn off the power supply before opening the safety lid. Although user safety is not endangered, arcing may damage the contacts.*

#### *Detection methods*

Following electrophoretic transfer, the membrane can be stored, stained or probed immediately.

#### **Further reading**

Electrophoresis in Practice: A guide to theory and practice. Westermeier, R., Ed., (1993) VCH Verlagsgesellschaft mbH. Weinheim. Westermeier, R.

Beisiegel, U., *Electrophoresis*, 7, 1-18 (1986)

Kyhse-Andersen, J., *J. Biochem. Biophys. Meth.*, 10, 203 (1984)

Naaby-Hansen, S., Lihme, A. O. F., Bog-Hansen, T.C., Bjerrum, O.J., in *Lectins-Biology, Biochemistry, Clinical Biochemistry*, Walter de Gruyter & Co., Berlin & New York, 241 (1985)

Handbook of immunoblotting of proteins. ed. Bjerrum, O. J. & Heegaard, N. H. H., CRC Press, Florida, USA, Volume 1 Technical descriptions, Volume II experimental and clinical applications.

Hancock, K. and Tsang, V.S.W., *Anal. Biochem.* 133, 157-162 (1983)

Towbin, H., Staehlin, T., Gordon, J., *Proc. Natl. Acad. Sci. USA*, 76, 4350- 4354 (1979)

## 5.7 Stock solutions

### B. Sample buffer

0.050 mol/l Tris-HAc pH 7.5

Dissolve 0.3 g Tris in 40 ml distilled water. Carefully adjust to pH 7.5 with HAc (approximately 0.14 ml). Make up to 50 ml with distilled water. Add 0.4 g SDS and a few grains of Bromophenol Blue. Immediately before use add 40 mg of DTT.

### C. Fixing solution

Ethanol 400 ml

Acetic acid, HAc 100 ml

Make up to 1000 ml with distilled water.

### D. Incubation solution

Ethanol 75 ml

Sodium acetate 17.00 g

Glutaraldehyde (25% w/v) 1.25 ml

Sodium thiosulphate,  $\text{Na}_2\text{S}_2\text{O}_3 \times 5\text{H}_2\text{O}$  0.50 g

Make up to 250 ml with distilled water.

### E. Silver solution

Silver nitrate 0.25 g

Formaldehyde 50  $\mu\text{l}$

Make up to 250 ml with distilled water.

### F. Developing solution

Sodium carbonate 6.25 g

Formaldehyde 25  $\mu\text{l}$

Make up to 250 ml with distilled water.

### G. Stop solution

EDTA- $\text{Na}_2 \times 2\text{H}_2\text{O}$  3.65 g

Make up to 250 ml with distilled water.

### H. Preserving solution

Glycerol (87% w/w) 25 ml

Make up to 250 ml with distilled water.

### I. Destaining solution

Ethanol 250 ml

Acetic acid 80 ml

Make up to 1000 ml with distilled water.

### K. Coomassie solution

PhastGel Blue R 1 tablet

Make up to 400 ml with destaining solution.

Heat to 60 °C, stirring constantly, and filter before use.

### L. Preserving solution

Glycerol (87% w/w) 25 ml

Make up to 250 ml with destaining solution.

### N. Fixing solution

Trichloroacetic acid 100 g

Make up to 500 ml with distilled water.

Transfer buffers using a discontinuous buffer system

R. Anode solution 1, pH 10.4

Tris	36.3 g
Methanol	200 ml

Make up to 1000 ml with distilled water.

S. Anode solution 2, pH 10.4

Tris	3.03 g
Methanol	200 ml

Make up to 1000 ml with distilled water.

T. Cathode solution, pH 7.6

6-Amino-n-hexanoic acid	5.20 g
Methanol	200 ml

Make up to 1000 ml with distilled water

U. Transfer buffer using a continuous buffer system

Glycine	2.93 g
Tris	5.81 g
SDS	0.375 g
Methanol	200 ml

Make up to 1000 ml with distilled water.

**Note:** *In a continuous buffer system, this solution is used for both anode and cathode electrode solutions.*

## 5.8 Running conditions for precast gels

### ExcelGel SDS Gradient 8–18

Running conditions

Phase	Voltage (V)	Current (mA)	Power (W)	Time (min)	Temp °C
Run	600	50	30	75*	15

\* Or until the Bromophenol Blue front reaches the anode buffer strip.

### ExcelGel XL SDS Gradient 12-14

Running conditions

Phase	Voltage (V)	Current (mA)	Power (W)	Time (min)	Temp °C
Run	1000	40	40	165*	15

\* Or until the Bromophenol Blue front reaches the anode buffer strip.

### ExcelGel SDS Homogeneous 7.5, 12.5 and 15

Running conditions

ExcelGel SDS Homogeneous	Voltage (V)	Current (mA)	Power (W)	Time (min)	Temp °C
7.5 and 12.5	600	50	30	80*	15
15	600	30	30	140*	15

\* Or until the Bromophenol Blue front reaches the anode buffer strip.

**Ampholine PAGplate for IEF**

Running conditions for Ampholine PAGplate.

pH range	Voltage (V)	Current (mA)	Power (W)	Time (h)	Temp °C
3.5–9–5	1 500	50	30	1.5	10
4.0–6.5	2 000	25	25	2.5	10
5.5–8.5	1 600	50	25	2.5	10
4.0–5.0	1 400	50	30	3.0	10
5.0–6.5	2 000	15	20	3.0	10

If half a gel is used, halve the current and power settings.

**CleanGel IEF 3–10**

Running conditions for one CleanGel IEF 3–10.

Phase	Voltage (V)	Current (mA)	Power (W)	Time (min)	Temp °C
Prefocusing	700	12	8	20	10
Sample entry	500	8	8	20	10
Isoelectric focusing	2000	14	14	90	10
Band sharpening	2500	14	18	10	10

If half a gel is used, halve the current and power settings.

**2-D electrophoresis using Immobiline DryStrip and ExcelGel SDS****First dimension**

Option 1: EPS 3500 XL Power Supply, using a voltage gradient.

The parameters below may be used for up to 12 strips.

Programme for Immobiline DryStrip, pH 3–10, 110 mm. All steps are run at 10 °C.

Phase	Voltage (V)	Current (mA)	Power (W)	Time (h)	Vh
1	300	1	5	0.1	1
2	300	1	5	4.5	1350
3	2000	1	5	5	5750
4	2000	1	5	6.5	13000
Total				16	20100*

Programme for Immobiline DryStrip, pH 3–10 L, 180 mm. All steps are run at 10 °C.

Phase	Voltage (V)	Current (mA)	Power (W)	Time (h)	Vh
1	500	1	5	0.1	1
2	500	1	5	3	1500
3	3500	1	5	5	10000
4	3500	1	5	12.5	43750
Total				20.5	55250*

Programme for Immobiline DryStrip, pH 3–10 NL, 180 mm. All steps are run at 10 °C.

Phase	Voltage (V)	Current (mA)	Power (W)	Time (h)	Vh
1	500	1	5	0.1	1
2	500	1	5	5	2500
3	3500	1	5	5	10000
4	3500	1	5	9.5	32400
Total				19.5	44900*

Programme for Immobiline DryStrip, pH 4–7, 110 mm. All steps are run at 10 °C.

Phase	Voltage (V)	Current (mA)	Power (W)	Time (h)	Vh
1	300	1	5	0.1	1
2	300	1	5	6	1800
3	3500	1	5	5	9500
4	3500	1	5	5.5	19250
Total				16.5	30550*

Programme for Immobiline DryStrip pH 4–7, 180 mm. All steps are run at 10 °C.

Phase	Voltage (V)	Current (mA)	Power (W)	Time (h)	Vh
1	500	1	5	0.1	1
2	500	1	5	1	500
3	3500	1	5	5	10000
4	3500	1	5	10	35000
Total				16	45500*

\* The optimal total number of Volt-hours for these pH gradients depends on the type of sample, sample load ( $\mu\text{g}$ ) and sample volume.

#### Option 2: Using a Manual Power Supply

The power supply should run at constant voltage with the parameters set as below. All steps are run at 10 °C.

Running conditions for Immobiline DryStrip, pH 3–10, 110 mm.

Phase	Voltage (V)	Current (mA)	Power (W)	Time (h)	Vh
1	300	1	5	1	300
2	1400	1	5	14–15	20000

Running conditions for Immobiline DryStrip, pH 3–10 L, 180 mm.

Phase	Voltage (V)	Current (mA)	Power (W)	Time (h)	Vh
1	500	1	5	1	500
2	3500	1	5	15–16	55000

Running conditions for Immobiline DryStrip pH 3–10 NL 180 mm.

Phase	Voltage (V)	Current (mA)	Power (W)	Time (h)	Vh
1	500	1	5	1	500
2	3500	1	5	13	45000

Running conditions for Immobiline DryStrip, pH 4–7, 110 mm.

Phase	Voltage (V)	Current (mA)	Power (W)	Time (h)	Vh
1	300	1	5	1	300
2	2200	1	5	13.5	29700

Running conditions for Immobiline DryStrip, pH 4–7, 180 mm.

Phase	Voltage (V)	Current (mA)	Power (W)	Time (h)	Vh
1	500	1	5	1	500
2	3000	1	5	14.5	43500

### Second dimension

Running conditions for ExcelGel XL SDS gradient 12–14.

Step	Voltage (V)	Current (mA)	Power (W)	Time (min)	Temp. (° C)
1	1000	20	40	45*	15
2	1000	40	40	5**	15
3	1000	40	40	160***	15

Running conditions for ExcelGel SDS gradient 8–18.

Step	Voltage (V)	Current (mA)	Power (W)	Time (min)	Temp. (° C)
1	600	20	30	25–30*	15
2	600	50	30	3–5**	15
3	600	50	30	70***	15

\* When the Bromophenol Blue dye front has moved 4–6 mm for ExcelGel XL SDS gradient 12–14 and 1–2 mm for ExcelGel SDS, gradient 8–18 from Immobiline DryStrip, remove the strip and the application pieces.

\*\* When the front has moved a further 2 mm, move the cathodic buffer strip forward to cover the area of removed Immobiline DryStrip by 1–2 mm. Adjust the position of the cathodic electrode.

\*\*\* When the Bromophenol Blue front has just reached the anodic buffer strip, electrophoresis is continued for 5 min and should then be stopped. Remove the buffer strips.

Further information about the gels and running conditions are supplied with the products.

Running conditions for ExcelGel 2-D Homogeneous 12.5

Phase	Voltage (V)	Current (mA)	Power (W)	Duration (h:min)
1	600	20	30	~0:35 <sup>1</sup>
2	600	50	30	~1:15 <sup>2</sup>

## 6 Maintenance

---



**WARNING!** Remove liquid or dirt from the system surface using a cloth and, if necessary, a mild cleaning agent.



**WARNING!** When using hazardous chemicals, take all suitable protective measures, such as wearing protective glasses and gloves resistant to the chemicals used. Follow local regulations and instructions for safe operation and maintenance of the system.



**WARNING!** When using hazardous chemicals, make sure that the entire system has been flushed thoroughly with bacteriostatic solution, e.g. NaOH, and distilled water before service and maintenance.

A few standard measures are necessary to keep Multiphor II in full functioning order.

After isoelectric focusing, remove the electrodes from the electrode holder and rinse with distilled water to remove the strong acidic and basic solutions. Do not submerge the cable containing the pin or socket. Air dry or carefully dry with paper tissue. Check that the platinum wire is not damaged.

After electrophoresis using the buffer chambers, remove the electrodes. Rinse them in distilled water and air dry. Take care not to damage the platinum electrodes.

Rinse the buffer chambers with distilled water between buffer changes and after use. Do not immerse the socket connector. Air dry or carefully dry with a paper towel.

Following electrophoretic transfer, remove all remaining filter papers from the NovaBlot unit. Remove the anode and cathode plates and rinse them in distilled water. Do not immerse the electrode leads in water. Leave to air dry. For longer life of NovaBlot electrodes store them either by:

1. Placing 3 cm thick plastic foam between the electrodes as if for transfer
- or
2. Store the electrodes on "backs" without foam sandwiched between.

### 6.1 Recycling



This symbol indicates that the waste of electrical and electronic equipment must not be disposed as unsorted municipal waste and must be collected separately. Please contact an authorized representative of the manufacturer for information concerning the decommissioning of equipment.





## 7 Technical specifications

---

Maximum Voltage	3500 Vp-p (+/- 1750 V with reference to Ground)
Maximum Power	100 W
Max pressure cooling plate	0.5 bar
Dimensions	16 x 31 x 40 cm
Environment	+4 – +40°C, 20–95% relative humidity
Material of wetted parts	
Chemical resistance	The wetted parts are resistant to solvents commonly used in electrophoresis and solutions containing inorganic and organic acids, alkalis and alcohols.
Compliance with standards	The declaration of conformity is valid for the instrument only if it is: <ul style="list-style-type: none"> <li>• used in laboratory locations</li> <li>• used in the same state as it was delivered from GE Healthcare except for alterations described in the User Manual</li> <li>• connected to other CE labelled GE Healthcare modules or other products as recommended.</li> </ul>
Safety standards	This product meets the requirement of the Low Voltage Directive (LVD) 73/23/EEC through the following harmonized standards: <ul style="list-style-type: none"> <li>• EN 61010-1</li> <li>• IEC 61010-1</li> <li>• CAN/CSA-C22.2 No. 61010-1</li> <li>• UL61010-1</li> </ul>
EMC standards	This device meets the requirements of the EMC Directive 89/336/EEC through the following harmonized standards: <ul style="list-style-type: none"> <li>• EN 61326 (emission and immunity)</li> <li>• EN 55011, GR 2, Class A (emission)</li> <li>• This device complies with part 15 of the FCC rules (emission). Operation is subject to the following two conditions: <ol style="list-style-type: none"> <li>1 This device may not cause harmful interference.</li> <li>2 This device must accept any interference received, including interference that may cause undesired operation.</li> </ol> </li> </ul>



## 8 Trouble shooting



This symbol indicates that the waste of electrical and electronic equipment must not be disposed as unsorted municipal waste and must be collected separately. Please contact an authorized representative of the manufacturer for information concerning the decommissioning of equipment.



**WARNING!** Turn OFF the power supply before opening the lid.

### Trouble shooting guide to PAGE

Symptom	Cause	Remedy
No current reading	Safety plug improperly inserted in power supply outlet	Check the safety plug insertion
	Pin and socket connection from electrode to base incomplete	Check the pin and socket connections
	Anode bridging contact disconnected	Connect the bridging contact
	Banana plug connection in safety lid not completed	Press firmly on the safety lid
	Electrode holder not seated properly	Lower the electrode holder so that the electrodes are in contact with the electrode strips
	Poor contact between the electrodes and electrode strips	Check that the electrodes are clean and intact, and sit in the centre of the isoelectric focusing strips over the entire length
Uneven migration of the dye front	Bad electrical contact between the gel and the wicks and electrodes	Check the contact
	Poor cooling	Check the cooling
Burning at slots or accumulation of water in the slots	Polypeptide complexes are too big to enter the gel and cause electroendosmosis	If SDS PAGE, add DTT once again and boil the sample

### Trouble shooting guide to IEF

Symptom	Cause	Remedy
Current increases with time	Electrode strips applied incorrectly in relation to electrode polarity	Check the electrode polarity and the pH of the electrode strips
	Cathode and anode polarities reversed	Check the pin and socket connections, the gel orientation and the pH of the applied electrode strips.

Symptom	Cause	Remedy
Sparkling on the gel	Gel dried out, insufficient cooling	Check the temperature and flow of the cooling fluid. Lower the power
Water droplets on gel	Excessive condensation	Decrease condensation by adjusting the temperature of the cooling fluid. Wipe the electrode holder periodically to remove condensation
Drying out of the gel near the electrodes	Incorrect electrode solutions	Use the recommended electrode solution at the specified concentration.
Sparkling along edge of gel onto cooling plate	Excessive power setting	Check the power setting
	Excess moisture on gel or under cooling plate	Remove the excess moisture
	Electrode strips overhanging the ends of the gel	Cut the electrode strips short of the ends of the gel
Condensation over the entire surface of the glass electrode holder	Liquid expelled at sides of electrode strips due to electroendosmotic of water towards the cathode	Occasionally remove the excess flow fluid by blotting
	Excessive power setting	Check the power setting. When only a portion of the gel is used, reduce the power setting proportionally
Local condensation on the glass electrode holder	Insufficient cooling	Check the temperature and flow of the cooling fluid
	Local overheating due to a high salt concentration in the sample	Reduce the salt content of the sample by gel filtration using PD-10 columns pre-packed with Sephadex G-25
Excessive amount of condensation along electrode strips	Incorrect electrode solutions in relation to electrode polarity	Check the electrode polarity. Check the pH of the applied electrode strips
	Localized hot spots due to air bubbles under the gel	Use insulating fluid under the gel and check for air bubbles
	Cathodic drift may cause an electroendosmotic flow of water towards the cathode. Thus, cathode strip may become over saturated Reversed polarity of electrode strips (lower pH at cathode, higher pH at anode)	Cut the electrode strips shorter than the edge of the gel. If necessary, blot the pooled liquid.  Check pH of the strips and polarity of the plugs in the power supply. Reverse polarity if the strips have been incorrectly applied (should be acid at anode)

Symptom	Cause	Remedy
Skewed or wavy bands	Localized gradient disturbances due to excessive salt	Reduce the salt content of the sample by gel filtration using PD-10 columns pre-packed with Sephadex G-25. Salt content should be <50 mmol/l. Too much ammonium persulphate may also cause wavy bands
	Unevenly wetted electrode strips	Electrode strips must be evenly wetted and be neither too wet nor too dry
	Electrode strips too short	The strips should be cut just short of the edges of the gel

### Trouble shooting guide to electrophoretic transfer

Symptom	Cause	Remedy
Incomplete transfer	Gel concentration too high	Use reversible cross-linkers (e.g. DATD, BAC, DHEBA) and depolymerize gel before transfer. Lower monomer concentration. Convert molecule to smaller form by limited digestion with proteases (for proteins) or with nucleases or acid hydrolysis (for nucleic acids)
	Methanol present in transfer buffer	Remove methanol from transfer buffer
	Transfer time too short	Increase transfer time
	Field strength too low	Increase field strength
Poor transfer	Too low charge/mass ratio	Change transfer buffer pH further away from molecules pI. Add 0.1% SDS
	Air trapped between gel and membrane	Carefully push out all air bubbles from the layers of the transfer sandwich
Inefficient transfer	Too low binding efficiency (molecules migrate from the gel, but pattern is faint)	Use different immobilizing membrane (DEAE, NC, DBM, DPT). Immobilizing membrane needs to be activated. Remove interfering substances (denaturants, detergents). Raise/lower salt concentration. Raise/lower pH

**Symptom**

**Cause**

Field strength too high  
(e.g. low Mol. Wt. DNA)

Transfer time too long

Pore size too large

**Remedy**

Lower field strength

Shorten transfer time

With nitrocellulose, use smaller  
pore filters (0.1 or 0.7  $\mu\text{m}$ )


## 9 Multiphor II application kits and accessories

---

This section describes the contents of the Multiphor II application kits and accessories and provides instructions for assembly and use. For experimental details including preparation of samples and stock solutions, running conditions, staining and preserving procedures see Chapter 4. Operation. Further information can be found in "Electrophoresis in Practice" – Code No. 18-1104-12.

### 9.1 SDS and Native PAGE, IEF Kit

This kit is used for casting 0.5 mm homogeneous or gradient polyacrylamide gels. The gels are cast on a 1 mm thick glass plate (125 x 260 mm). Alternatively, casting can be done on GelBond™ PAGfilm (124 x 258 mm). Optional glass plates with U-frames allow casting of 1.0 and 2.0 mm thick gels. The kit includes sample application pieces and strips for applying the sample onto the gel surface. An optional template and tape allow preparation of a slot-former.




#### Kit contents – Code No. 18-1102-45

Designation	Code No.
Glass Plate, 125 x 260, 0.5 mm U-frame (2/pkg)	80-1106-89
Glass Plate, 125 x 260x1 mm (15/pkg)	80-1106-29
Glass Plate, 125 x 260x3 mm (2/pkg)	80-1106-99
FlexiClamp (6/pkg)	18-1013-73
IEF Electrode Strip (100/pkg)	18-1004-40
Electrophoresis Wick 104 x 253 mm (500/pkg)	80-1129-52
IEF Sample Application. Pieces (200/pkg)	80-1129-46
IEF/SDS Sample Application. Strip, 52 samples 5–20 $\mu$ l (5/pkg)	18-1002-26
Cellophane Sheets (50/pkg)	80-1129-38

**Optional accessories**


<b>Designation</b>	<b>Code No.</b>
Roller	80-1106-79
GelBond PAGfilm, 124 x 258 mm (50/pkg)	80-1129-36
Bind-Silane, 100 ml	17-1330-01
Repel-Silane, 500 ml	17-1331-01
SDS Sample Application Strip, 26 samples, 40 $\mu$ l	18-1002-74
Gradient Maker SG 100, 100 ml	80-6196-09
Glass Plate, 125 x 260 mm, 1,0 mm U-frame (2/pkg)	80-1106-91
Glass Plate, 125 x 260 mm, 2,0 mm U-frame (2/pkg)	80-1106-92
Tape, Dymo 0,25 x 9 mm, 3 m	80-1129-50
Template, 125 x 260 mm (10/pkg)	80-1129-55

The 3 mm thick glass plate is used as a support, either for the 1 mm glass plate or GelBond PAGfilm. The mould comprising the 3 mm glass plate, 1 mm glass plate or GelBond PAGfilm and glass plate with U-frame is clamped together using four FlexiClamps.



To prepare a slot-former for individual sample slots, a glass plate with Uframe, tape, 0,25 x 9 mm, and a template should be used. One or several layers of tape can be applied to the glass plate. For instance, 3 layers of 5 x 3 mm will make a sample slot for 10-20  $\mu$ l of sample. Wash the glass plate with detergent, rinse with distilled water and dry with a paper tissue. Using the template as a guide, apply the tape 30 mm from the open edge of the U-framed glass plate avoiding air bubbles. Check that all edges of the tape are cut perfectly even. Leave the slot former over night to ensure that the tape adheres completely.





To prevent the gel from sticking to the U-framed glass plate, coat the plate with Repel-Silane.


**Note:** For this operation use gloves and a fume hood.

Pour about 2 ml of Repel-Silane onto the glass plate and distribute it evenly with a tissue. Leave it to dry for a few minutes. Rinse the glass plate with distilled water and remove water drops by shaking or wiping lightly with a tissue. Leave the glass plate to dry.


When using the 1 mm thick glass plate as the gel support, simply lay it directly on top of the 3 mm thick glass plate. If the gel is to be permanently bound to the 1 mm thick glass plate, coat the plate with Bind-Silane, before preparing the mould.

**Note:** For this operation use gloves and a fume hood.


Pour about 2 ml of diluted Bind-Silane onto the glass plate and distribute it evenly with a tissue. Leave the glass plate to dry for a few minutes, rinse with distilled water and leave to dry.



When using GelBond PAGfilm, pour a few ml of water on to the 3 mm thick glass plate and lay the film over it with the hydrophilic side up (see Instructions supplied with the film). Centre the film on the glass plate. Beginning at one end, use the roller to apply even pressure over the film surface in order to eliminate air bubbles and seal the film to the plate with a minimum of water. Remove any excess water with a tissue.



Form the mould by placing the U-framed glass plate in position and clamp together using four FlexiClamps.



**Note:** *Gloves must be worn to protect the user from contact with the toxic acrylamide solution.*

Draw the gel solution into a syringe or a graduated pipette. Fill the mould, checking that air bubbles are not trapped along the rubber U-frame or around the slots.


When casting a gradient gel, position the mould horizontally using the Levelling Set and place the Gradient Maker as illustrated. Lay the end of the tubing from the Gradient Maker against the 1 mm glass plate or GelBond PAGfilm. The slot-former will otherwise disturb the flow of the solution.

To open the 1 mm glass plate mould, remove the four FlexiClamps. Carefully insert one or two thin-bladed spatulas between the gel surface and slot former on one of the short sides. Twist gently in order to introduce air across the whole of the short side. Twist more firmly to slowly separate the U-frame from the gel surface. Remove the U-frame. Carefully remove any unpolymerized acrylamide from the edge of the gel with a paper tissue. Separate the gel support from the thick glass plate. The gel is now ready to use.

To open the mould including GelBond PAGfilm, remove the four FlexiClamps and insert the spatula between the 3 mm thick glass plate and film. Remove the glass plate and dry the back of the film. Turn the mould upside down (with the glass plate with U-frame on top) and gently peel the film with gel away from the glass.

## 9.2 Buffer Strip Positioner

The Multiphor II Buffer Strip Positioner is a frame with slots that sits on top of an ExcelGel SDS gel on the Multiphor II cooling plate. The slots in the positioner facilitate placement of the buffer strips for electrophoresis and hold them securely in place. A locking cam secures the positioner on the cooling plate.



**Fig. 1.** Features of the Multiphor II Buffer Strip Positioner

**Slot# Use for placing**

- 1 Cathodic buffer strip, Phase 1 or entire run
- 2 Sample wells  
Immobiline DryStrip gels (IPG strips), Phase 1  
Cathodic buffer strip, Phase 2
- 3 Anodic buffer strip (with 11 x 25 cm ExcelGel SDS gels)
- 4 Anodic buffer strip (with 18 x 25 cm ExcelGel SDS gels)

**Designation**

Multiphor II Buffer Strip Positioner

**Code No.**

80-6442-90

### 9.3 Immobiline DryStrip Reswelling Try

The Immobiline DryStrip Reswelling Trays have twelve independent reservoir slots that can each hold a single IPG strip. Separate slots allow the rehydration of individual IPG strips in a minimal volume of solution.

Designation	Code No.
Immobiline DryStrip Reswelling Tray, for 7–18 cm IPG Strips	80-6371-84
Immobiline DryStrip Reswelling Tray, for 7–24 cm IPG Strips	80-6465-32

### 9.4 Immobiline DryStrip Kit

This kit is used for running IEF with Immobiline DryStrip for the first dimension in 2-D electrophoresis. Twelve strips can be focused simultaneously under a protective layer of silicone oil. The high sample capacity allows the application of up to 100 µl on each Immobiline DryStrip. Detailed instructions for use are available in the instruction manual provided with this kit.



#### Kit contents - Code No. 18-1004-30

Designation	Code No.
Tray and Electrode Holder	18-1004-31
DryStrip Aligner (4/pkg)	18-1004-34
DryStrip Kit Electrode, cathode	18-1018-67
DryStrip Kit Electrode, anode	18-1018-66
Sample Cup Bar	18-1004-33
Sample Cup (6 x 10/pkg)	18-1004-35
IEF Electrode Strip (100/pkg)	18-1004-40
IEF Sample Application. Piece (200/pkg)	80-1129-46
Instruction Manual	18-1038-63

#### Optional accessories

Designation	Code No.
Immobiline DryStrip Reswelling Tray, for 7–18 cm IPG Strips	80-6371-84
Immobiline DryStrip Reswelling Tray, for 7–24 cm IPG Strips	80-6465-32

## 9.5 NovaBlot Kit

This kit is used for electrophoretic transfer of proteins from polyacrylamide or agarose gels to an immobilizing membrane. The maximum gel size is 200 x 250 mm.

By building transfer sandwiches, simultaneous transfer from several gels of the same type can be achieved. Up to six transfer sandwiches can be stacked one on top of the other.

If 125 x 250 mm gels are to be transferred, NovaBlot accepts up to six gels for simultaneous transfer by assembling two transfer sandwiches side by side.

The operating procedures for NovaBlot Kit and FilmRemover are described and illustrated in Sections 4.8 and 7.13 respectively.



### Kit contents – Code No. 18-1016-86




Designation	Code No.
NovaBlot Electrode, cathode	18-1019-86
NovaBlot Electrode, anode	80-1257-87
Electrode Paper NovaBlot, 200 x 250 mm (500/pkg)	80-1106-19
Cellophane Sheets, 210 x 320 mm (50/pkg)	80-1129-38

### Optional accessories

Designation	Code No.
FilmRemover	18-1013-75
Nitrocellulose 0.20 µm, 150 x 200 mm (15/pkg)	80-1098-91

## 9.6 FilmRemover

FilmRemover is used for removing backing from a gel before electrophoretic transfer. Polyacrylamide or agarose gels with a thickness between 0.1 mm and 5.0 mm and a maximum gel size of 200 x 245 mm can be used.



Detailed instructions for the use of FilmRemover are available in the instruction manual provided with the product.

### Unit contents – 18-1013-75

Designation	Code No.
FilmRemover basic unit	80-1316-21
Lever and Wire Assembly (3/pkg)	18-1013-79
Instruction Manual	80-1316-37
Nitrocellulose 0.20 $\mu\text{m}$ , 150 x 200 mm (15/pkg)	80-1098-91
Nitrocellulose 0.45 $\mu\text{m}$ , 150 x 200 mm (15/pkg)	80-1098-90
ProBind 45 NC 0.45 $\mu\text{m}$ , roll 0.2 x 3.0 m	80-1247-86
GeneBind 45 nylon 0.45 $\mu\text{m}$ , roll 0.2 x 3.0 m	80-1247-87

### 9.7 Roller

For use when applying plastic support films onto glass plates with an interfacing fluid. The roller is used to provide even pressure over a large area, ensuring adhesion with a minimum amount of fluid and elimination of bubble formation.



Designation	Code No.
Roller	80-1106-79

Several accessories are available for simple and convenient sample application with Multiphor II.

### 9.8 Sample application accessories

Designation	Code No.
IEF Sample Application Pieces (200/pkg)	80-1129-46



The 5 x 10 mm sample application piece made of Paratex can be used for sample volumes in the range 15–20  $\mu$ l.

Up to 52 samples can be applied with this strip. Each well holds up to 20  $\mu$ l of sample. The applicator strip is made of flexible silicone and is applied directly onto the gel surface.

<b>Designation</b>	<b>Code No.</b>
SDS Sample Application Strip 26 samples, 40 $\mu$ l	18-1002-74



This strip is recommended for sample application on SDS gradient gels without preformed slots, e.g. ExcelGel SDS, gradient 8–18. Up to 26 samples can be applied and each well holds up to 40  $\mu$ l of sample. The strip is made of transparent flexible silicone and is applied directly onto the gel surface.

<b>Designation</b>	<b>Code No.</b>
IEF/SDS Sample Application Strip 52 samples, 5–20 $\mu$ l	18-1002-26



This applicator strip is recommended for use with PAGIEF gels, and SDS gradient gels without preformed slots, e.g. ExcelGel SDS, gradient 8–18. Up to 52 samples with a sample volume of 5–20  $\mu$ l can be applied in each well.

The applicator strip is made of flexible silicone and is applied directly on the gel surface.




<b>Designation</b>	<b>Code No.</b>
EPH/IEF Sample Application Foil 24 samples, 2-4 $\mu$ l	80-1129-47



This application foil with narrow slits is recommended for electrophoresis and IEF in agarose gels. Up to 24 samples can be applied, with a sample volume of 2-4  $\mu$ l in each slit. The foil is applied directly on the gel surface.

<b>Designation</b>	<b>Code No.</b>
EPH Electrode anode, long	80-1122-20
EPH Electrode cathode, long	80-1122-19



The electrophoresis electrodes are designed for use with the buffer vessels at the side of the buffer tank, allowing electrophoresis along the width of the cooling plate.



## 10 Ordering information



**WARNING!** Only spare parts approved or supplied by GE Healthcare may be used for maintaining and servicing of MULTIPHOR II

### 10.1 Multiphor II

Product	Quantity	Code No.
<b>Basic configuration</b>		
Multiphor II Electrophoresis Unit	1	18-1018-06
<b>Application kits and accessories</b>		
SDS and Native PAGE, IEF Kit	1	18-1102-45
Immobiline DryStrip Reswelling Tray, for 7–18 cm IPG strips	1	80-6371-64
Immobiline DryStrip Reswelling Tray, for 7–24 cm IPG strips	1	80-6465-32
Immobiline DryStrip Kit for running 1 to 12 Immobiline DryStrip gels (for use with Multiphor II only)	1	18-1004-30
NovaBlot Kit for electrophoretic transfer	1	18-1016-86
Multiphor II Buffer Strip Positioner, complete	1	80-6442-90
FlexiClamps	6	18-1013-73
FilmRemover for removing plastic gel backing before electrophoretic transfer	1	18-1013-75
Gradient Maker, 100 ml	1	80-6196-09
Roller	1	80-1106-79
Template, 125x260 mm	10	80-1129-55
Lever and Wire Assemblies (for Film-Remover)	3	18-1013-79
Levelling Feet	4	18-1026-40
<b>Cooling plates</b>		
Cooling Plate ceramic, 210 x 270 mm	1	18-1103-46
Grommets	2	80-1106-58
Cooling Tubing, 8/12 mm	4 m	80-1106-56
Tubing Connector Set, female and male	4	18-1104-26
Insulation for Cooling Tubing 14/27 mm	8 m	80-1116-11
Anode and Cathode Electrode	1	18-1037-44
Leads (for 18-1004-31 Immobiline DryStrip Kit Tray)		

Product	Quantity	Code No.
<b>Electrodes and electrode holders</b>		
EPH/IEF Electrode, anode	1	18-1121-53
EPH/IEF Electrode, cathode	1	18-1121-52
Electrode Holder (for 18-1106-60/-61)	1	80-1106-55
EPH Electrode (long), anode	1	18-1122-20
EPH Electrode (long), cathode	1	18-1122-19
Electrode (Immobiline DryStrip Kit), anode	1	18-1018-66
Electrode (Immobiline DryStrip Kit), cathode	1	18-1018-67
Tray and Electrode Holder for 18-1018-66/-67	1	18-1004-31
NovaBlot Electrode, anode	1	80-1257-87
NovaBlot Electrode, cathode	1	18-1019-86
<b>Glass plates and trays</b>		
125 x 260 x 3 mm	2	80-1106-99
125 x 260 x 1 mm	15	80-1106-29
125 x 260 mm, 0.5 mm U-frame	2	80-1106-89
125 x 260 mm, 1.0 mm U-frame	2	80-1106-91
125 x 260 mm, 2.0 mm U-frame	2	80-1106-92
<b>Glass plate treatment</b>		
Bind-Silane	100 ml	17-1330-01
Repel-Silane	500 ml	17-1331-01
<b>Gel support</b>		
GelBond PAGfilm, 124 x 258 mm	50	80-1129-36
<b>Paper electrode strip and wicks</b>		
IEF Electrode Strip	100	18-1004-40
EPH Electrode Wick, 82 x 130 mm	500	80-1129-53
EPH Electrode Wick 104 x 253 mm (also used as PEGG print paper and with agarose IEF)	500	80-1129-52
Electrode Paper NovaBlot, 200 x 250 mm	500	80-1106-19

Product	Quantity	Code No.
<b>Sample application</b>		
SDS Sample Application Strip, 26 samples, 40 µl	5	18-1002-74
IEF/SDS Sample Application Strip, 52 samples, 5–20 µl	5	18-1002-26
IEF Sample Application Pieces	200	80-1129-46
Sample Cups, Immobiline DryStrip Kit	60	18-1004-35
<b>Preserving</b>		
Cellophane Sheets, 210 x 320 mm	50	80-1129-38
Mylar™ Sheets, 125 x 260 mm		5080-1129-39
<b>Membranes, electrophoretic transfer</b>		
Nitrocellulose, 0.20 µm, 150x200 mm	15	80-1098-91
Nitrocellulose, 0.45 µm, 150x200 mm	15	80-1098-90
ProBind 45 NC 0.45 µm, roll 0.2x3.0 m	1	80-1247-86
GeneBind 45 nylon 0.45 µm, roll 0.2x3.0 m		80-1247-87
<b>Tapes</b>		
Dymo, 0.25x9 mm, 3 m	1	80-1129-50

## 10.2 MultiTemp III

Product	Quantity	Code No.
MultiTemp III thermostatic circulator, 100–120 V	1	18-1102-77
MultiTemp III thermostatic circulator, 220–220 V	1	18-1102-78
Cooling Tubing, 8/12 mm	4 m	80-1106-56
Tubing Connector Set, female and male	4	18-1104-26
Insulation for Cooling Tubing 14/27 mm	2 m	80-1116-11
3-way Valve Set	1	18-1106-39

## 10.3 EPS Power Supplies

Product	Quantity	Code No.
EPS 3501 XL 35-3500 V, 1-400 mA	1	18-1130-05
EPS 3501 35-3500 V, 1-150 mA	1	18-1130-04
EPS 1001 5-1000 V, 1-400 mA	1	18-1130-03
EPS 601 6-600 V, 1-400 mA	1	18-1130-02
EPS 301 5§-300 V, 1-400 mA	1	18-1130-01

### 10.4 Precast gels and buffer strips

Product	Quantity	Code No.
<b>SDS-PAGE and Native PAGE</b>		
ExcelGel SDS Homogeneous 7.5	6	80-1260-01
ExcelGel SDS Homogeneous 12.5	6	80-1261-01
ExcelGel SDS Homogeneous 15	6	80-1262-01
ExcelGel SDS, gradient 8–18	6	80-1255-53
ExcelGel XL SDS, gradient 12–14	3	17-1236-01
ExcelGel SDS Buffer Strips anode and cathode	6 each	17-1342-01
<b>IEF</b>		
CleanGel IEF	5	18-1035-32
GelPool for gel rehydration	1	18-1031-58
PaperPool for electrode strips	1	18-1031-59
Ampholine PAGplate pH 3.5–9.5	5	80-1124-80
Ampholine PAGplate pH 4.0–6.5	5	80-1124-81
Ampholine PAGplate pH 5.5–8.5	5	80-1124-82
Ampholine PAGplate pH 4.0–5.0	5	80-1124-83
<b>Immobiline DryStrip Gels</b>		
Immobiline DryStrip pH 3.5–4.5, 24 cm	12	17-6002-38
Immobiline DryStrip pH 4.0–5.0, 24 cm	12	17-6002-39
Immobiline DryStrip pH 4.5–5.5, 24 cm	12	17-6002-40
Immobiline DryStrip pH 5.0–6.0, 24 cm	12	17-6002-41
Immobiline DryStrip pH 5.5–6.7, 24 cm	12	17-6002-42
Immobiline DryStrip pH 6–9, 24 cm	12	17-6002-47
Immobiline DryStrip pH 3–7 NL, 24 cm	12	17-6002-43
Immobiline DryStrip pH 3–10, 24 cm	12	17-6002-44
Immobiline DryStrip pH 3–10 NL, 24 cm*	12	17-6002-45
Immobiline DryStrip pH 4–7, 24 cm	12	17-6002-46
Immobiline DryStrip pH 3.5–4.5, 18 cm	12	17-6001-83
Immobiline DryStrip pH 4.0–5.0, 18 cm	12	17-6001-84
Immobiline DryStrip pH 4.5–5.5, 18 cm	12	17-6001-85
Immobiline DryStrip pH 5.0–6.0, 18 cm	12	17-6001-86
Immobiline DryStrip pH 5.5–6.7, 18 cm	12	17-6001-87
Immobiline DryStrip pH 4–7, 18 cm	12	17-1233-01
Immobiline DryStrip pH 6–9, 18 cm	12	17-6001-88
Immobiline DryStrip pH 6–11, 18 cm	12	17-6001-97
Immobiline DryStrip pH 3–10 NL, 18 cm*	12	17-1235-01
Immobiline DryStrip pH 3–10, 18 cm	12	17-1234-01
Immobiline DryStrip pH 4–7, 13 cm	12	17-6001-13
Immobiline DryStrip pH 6–11, 13 cm	12	17-6001-96
Immobiline DryStrip pH 3–10 NL, 13 cm*	12	17-6001-15
Immobiline DryStrip pH 3–10, 13 cm	12	17-6001-14

Product	Quantity	Code No.
Immobiline DryStrip pH 4-7, 11 cm	12	18-1016-60
Immobiline DryStrip pH 6-11, 11 cm	12	17-6001-10
Immobiline DryStrip pH 3-10, 11 cm	12	18-1016-61
Immobiline DryStrip pH 4-7, 7 cm	12	17-6001-10
Immobiline DryStrip pH 6-11, 7 cm	12	17-6001-94
Immobiline DryStrip pH 3-10 NL, 7 cm*	12	17-6001-12
Immobiline DryStrip pH 4-7, 7 cm	12	17-6001-10
Immobiline DryStrip pH 3-10, 7 cm	12	17-6001-11
<b>IPG Buffer</b>		
IPG Buffer pH 3.5-5.0†	1 ml	17-6002-02
IPG Buffer pH 4.5-5.5	1 ml	17-6002-04
IPG Buffer pH 5.0-6.0	1 ml	17-6002-05
IPG Buffer pH 5.5-6.7	1 ml	17-6002-06
IPG Buffer pH 4-7§	1 ml	17-6000-86
IPG Buffer pH 6-11‡	1 ml	17-6001-78
IPG Buffer pH 3-10 NL*	1	17-6000-88
IPG Buffer pH 3-10	1 ml	17-6000-87

\* NL= increased resolution between pH 5-7

† Use IPG Buffer pH 3.5-5.0 for pH 3.5-4.5 and 4.0-5.0 IPG strips

‡ Use IPG Buffer pH 6-9 and pH 6-11 IPG strips.

§ Use IPG Buffer pH 4-7 for pH 3-7 IPG strips.

Product	Quantity	Code No.
<b>Second dimension</b>		
ExcelGel 2-D Homogeneous	6	17-6002-21
ExcelGel SDS, gradient 8-18 110 x 245 x 0.5 mm	6	80-1255-53
ExcelGel XL SDS, gradient 12-14 180 x 245 x 0.5 mm	3	17-1236-01
ExcelGel SDS Buffer Strips anode and cathode	6 each	17-1342-01

### 10.5 Molecular weight and pI markers

Product	Quantity	Code No.
MW range 2.512-16.949, 2 mg	1	80-1129-83
MW range 14.000-94.000, 200 µg/vial	10	17-0446-01
MW range 53.000-212.000, 200 µg/vial	10	17-0615-01
MW range 67.000-670.000, 200 µg/vial	10	17-0445-01
Broad pI kit pH 3.5-9.3		17-0471-01
Low pI kit pH 2.8-6.5		17-0472-01
High pI kit pH 5.2-10.3		17-0473-01
Carbamylate calibration kit		17-0582-01

## 10.6 Carrier ampholytes

Product	Quantity	Code No.
<b>Ampholine</b>		
Ampholine, preblended pH 3.5–9.5	25 ml	80-1127-15
Ampholine, preblended pH 4.0–6.5	25 ml	80-1127-17
Ampholine, preblended pH 5.0–8.0	25 ml	80-1127-19
Ampholine pH 3.5–10.0	25 ml	80-1125-87
Ampholine pH 3.5–5.0	25 ml	80-1125-89
Ampholine pH 4.0–6.0	25 ml	80-1125-90
Ampholine pH 5.0–7.0	25 ml	80-1125-91
Ampholine pH 5.0–8.0	25 ml	80-1125-92
Ampholine pH 6.0–8.0	25 ml	80-1125-93
Ampholine pH 7.0–9.0	25 ml	80-1125-94
<b>Pharmalyte</b>		
Pharmalyte pH 3–10	25 ml	17-0456-01
Pharmalyte pH 2.5–5	25 ml	17-0451-01
Pharmalyte pH 4–6.5	25 ml	17-0452-01
Pharmalyte pH 5–8	25 ml	17-0453-01
Pharmalyte pH 8–10.5	25 ml	17-0455-01
Pharmalyte pH 4.2–4.9	25 ml	17-0562-01
Pharmalyte pH 4.5–5.4	25 ml	17-0563-01
Pharmalyte pH 5–6	25 ml	17-0564-01
Pharmalyte pH 6.7–7.7	25 ml	17-0566-01
<b>Immobiline</b>		
Immobiline II pK 3.6	10 ml	80-1255-70
Immobiline II pK 4.6	10 ml	80-1255-71
Immobiline II pK 6.2	10 ml	80-1255-72
Immobiline II pK 7.0	10 ml	80-1255-73
Immobiline II pK 8.5	10 ml	80-1255-74
Immobiline II pK 9.3	10 ml	80-1255-75

Each bottle contains a ready-to-use 0.200±0.004 M solution.



## 10.7 PlusOne electrophoresis chemicals

### PlusOne electrophoresis chemicals

Product	Use	Quantity	Storage	Code No.
<b>Gel casting chemicals</b>				
Acrylamide IEF	IEF, PAGE, Sequencing	250 g	A	17-1300-01
Acrylamide IEF 40% solution	IEF, PAGE	1000 ml	D	17-1301-01
ReadyMix IEF	IEF	41.5 g1	C	17-1309-01
ReadySol IEF T40 C3	IEF	1000 ml	D	17-1310-01
Acrylamide PAGE	PAGE	250 g	A	17-1302-01
Acrylamide PAGE	PAGE	1000 g	A	17-1302-02
Acrylamide PAGE 40% Solution	PAGE	1000 ml	D	17-1303-01
ReadySol DNA PAGE T40 C5	PAGE, Sequencing	1000 ml	D	17-1308-01
N,N-Methylene bis-acrylamide	IEF, PAGE, Sequencing	25 g	C	17-1304-01
N,N-Methylene-bis-acrylamide	IEF, PAGE, Sequencing	100 g	C	17-1304-02
N,N-Methylene-bis-acrylamide 2% solution	IEF, PAGE, Sequencing	1000 ml	D	17-1306-01
Ammonium persulphate	IEF, PAGE, Sequencing	25 g	C	17-1311-01
TEMED	IEF, PAGE,	25 ml	C*	17-1312-01
Sequencing				
<b>Buffers</b>				
Tris	PAGE, Sequencing	500 g	A	17-1321-01
Boric acid	PAGE, Sequencing	500 g	A	17-1322-01
EDTA, di-sodium salt	PAGE, Sequencing	100 g	A	17-1324-01
Glycine	PAGE, Sequencing	500 g	A	17-1323-01
<b>Additives and sample treatment</b>				
Urea	IEF, PAGE, Sequencing	500 g	B	17-1319-01
Formamide	IEF, PAGE, Sequencing	250 ml	B	17-1320-01
Dithiothreitol	IEF, SDS-PAGE	1.0 g	F	17-1318-01
Dithiothreitol	IEF, SDS-PAGE	5 g	F	17-1318-02
Mercaptoethanol	IEF, SDS-PAGE	25 ml	B	17-1317-01
Glycerol 87%	IEF, PAGE, Sequencing	1000 ml	A	17-1325-01
<b>Detergents</b>				
Sodium dodecylsulphate	PAGE	100 g	A	17-1313-01
Triton™ X-100	IEF, PAGE, Sequencing	500 ml	G	17-1315-01
CHAPS IEF,	PAGE, Sequencing	1 g	F	17-1314-01
Tween™ 20	IEF, PAGE, Sequencing	500 ml	G	17-1316-01
<b>Stains</b>				
Silver Staining Kit, Protein	Protein detection	For 10–20 gels	D	17-1150-01
Silver Staining Kit, DNA	Nucleic and detection	For 10–20 gels	D	17-6000-30
Ethidium bromide solution 10 mg/ ml	DNA/RNA detection	10 ml	A	17-1328-01
Bromophenol Blue IEF,	PAGE, Sequencing	10 g	A	17-1329-01
<b>Glass plate treatment</b>				
Repel-Silane ES	IEF, PAGE, Sequencing	500 ml	C	17-1331-01
Bind-Silane	IEF, PAGE, Sequencing	25 ml	C	17-1330-01
Others				
DryStrip Cover Fluid 2-D Immobiline	DryStrip	1000 ml	G	17-1335-01
Amberlite IRN-150L	Purifying solutions	500 g	A	17-1326-01

**Storage:** A, room temp. B, dry at room temp. C, dry & dark at room temp. D, dark at 4 °C to 8 °C. E, dry & dark at 4 °C to 8 °C. F, dry at 4 to 8 °C. G, dark at room temp.  
\*Store well sealed.

1. Add 100 ml. 2. Add 500 ml.





www.gehealthcare.com

GE Healthcare Bio-Sciences AB  
Björkgatan 30  
751 84 Uppsala  
Sweden

Multiphor, Immobiline, CleanGel, Ampholine, MultiTemp, ExcelGel, Pharmalyte and NovaBlot are trademarks of GE Healthcare companies. GE, imagination at work and GE monogram are trademarks of General Electric Company.

Hoefer is a trademark of Hoefer Inc. Triton is a trademark of Union Carbide Chemicals and Plastics Co. Tween is a trademark of ICI Americas Inc. Rynite is a trademark of DuPont Company.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. GE Healthcare reserves the right, subject to any regulatory and contractual approval, if required, to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your local GE Healthcare representative for the most current information.

© 2006 General Electric Company – All rights reserved.

GE Healthcare Bio-Sciences AB a General Electric Company.

GE Healthcare Bio-Sciences AB  
Björkgatan 30, 751 84 Uppsala, Sweden

GE Healthcare Europe GmbH  
Munzinger Strasse 5, D-79111 Freiburg, Germany

GE Healthcare UK Ltd  
Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA, UK

GE Healthcare Bio-Sciences Corp  
800 Centennial Avenue, P.O. Box 1327, Piscataway, NJ 08855-1327, USA

GE Healthcare Bio-Sciences KK  
Sanken Bldg. 3-25-1, Hyakunincho, Shinjuku-ku, Tokyo 169-0073, Japan

Asia Pacific Tel: +852 2811 8693 Fax: +852 2811 5251 • Australasia Tel: +61 2 9899 0999 Fax: +61 2 9899 7511 • Austria Tel: 01/57606-1619 Fax: 01/57606-1627 • Belgium Tel: 0800 73 888 Fax: 03 272 1637 • Canada Tel: 800 463 5800 Fax: 800 567 1008 • Central, East, & South East Europe Tel: +43 1 982 3826 Fax: +43 1 985 8327 • Denmark Tel: 45 16 2400 Fax: 45 16 2424 • Finland & Baltics Tel: +358-09-512 39 40 Fax: +358 (0)9 512 39 439 • France Tel: 01 69 35 67 00 Fax: 01 69 41 96 77 • Germany Tel: 0761/4903-490 Fax: 0761/4903-405 • Italy Tel: 02 27322 1 Fax: 02 27302 212 Japan Tel: +81 3 5331 9336 Fax: +81 3 5331 9370 • Latin America Tel: +55 11 3933 7300 Fax: +55 11 3933 7304 • Middle East & Africa Tel: +30 210 9600 687 Fax: +30 210 9600 693 • Netherlands Tel: 0165 580 410 Fax: 0165 580 401 • Norway Tel: 815 65 555 Fax: 815 65 666 • Portugal Tel: 21 417 7035 Fax: 21 417 3184 • Russia & other C.I.S. & N.I.S Tel: +7 (095) 232 0250, 956 1137 Fax: +7 (095) 230 6377 • South East Asia Tel: 60 3 8024 2080 Fax: 60 3 8024 2090 • Spain Tel: 93 594 49 50 Fax: 93 594 49 55 • Sweden Tel: 018 612 1900 Fax: 018 612 1910 • Switzerland Tel: 0848 8028 12 Fax: 0848 8028 13 UK Tel: 0800 616928 Fax: 0800 616927 • USA Tel: 800 526 3593 Fax: 877 295 8102



imagination at work

User Manual 18-1103-43 AK 02/2006