



201PR-02

G-Biosciences ♦ 1-800-628-7730 ♦ 1-314-991-6034 ♦ [technical@GBiosciences.com](mailto:technical@GBiosciences.com)

A Geno Technology, Inc. (USA) brand name

# HOOK™ Agarose Coupling Kit (Amine Reactive)

For Agarose Affinity Columns  
Using Your Proteins or Peptides

(Cat. # 786-063)



think proteins! think G-Biosciences [www.GBiosciences.com](http://www.GBiosciences.com)

INTRODUCTION ..... 3

ITEM(S) SUPPLIED ..... 3

STORAGE CONDITION ..... 3

ADDITIONAL ITEMS REQUIRED ITEMS ..... 3

PREPARATION BEFORE USE ..... 4

PROTOCOL ..... 4

APPENDIX: GENERAL PROTOCOL FOR AFFINITY PURIFICATION OF PROTEIN ..... 6

    SAMPLE PREPARATION ..... 6

    ADDITIONAL MATERIALS REQUIRED ..... 6

    PROTOCOL ..... 6

RELATED PRODUCTS ..... 7

## INTRODUCTION

The amine reactive HOOK™ Agarose Coupling Kit is designed for the simple and efficient coupling of proteins and peptides to a solid agarose support.

This kit utilizes 6% agarose that has been activated to generate reactive aldehyde groups. The aldehyde groups of the agarose react spontaneously with primary amines, located at the N-terminus of proteins or in lysine residues, to form intermediate Schiff Base complexes. These, in turn, are selectively reduced by reductive amination to form stable amine linkages between the agarose and the ligand.

This kit allows for the generation of five stable and reusable columns (4ml slurry/column) for the affinity purification of antibodies, proteins and other protein interacting biomolecules.

## ITEM(S) SUPPLIED (Cat. # 786-063)

| Description   | Size     |
|---|----------|
| 5X Optimizer Buffer™ VI                             | 2 x 25ml |
| Sodium Cyanoborohydride (NaCNBH <sub>3</sub> ) [5M] | 0.5ml    |
| HOOK™ Activated Agarose slurry (Amine reactive)     | 20ml     |
| Empty Disposable Columns (Part# 069E-B)             | 5        |
| Quenching Buffer (3M)                               | 25ml     |

## STORAGE CONDITION

The kit is shipped at ambient temperature. On arrival store 5X Optimizer Buffer™ VI and HOOK™ Activated Agarose slurry (Amine reactive) at 4°C, Sodium Cyanoborohydride (NaCNBH<sub>3</sub>) desiccated at room temperature (RT) and all other components at RT. When stored and used properly the kit is stable for one year.

## ADDITIONAL ITEMS REQUIRED ITEMS

- NaOH (25mM)
- Phosphate Buffer Saline (PBS)
- Sodium Azide.

## PREPARATION BEFORE USE

Prepare 1X Optimizer Buffer™ VI (1ml 5X Optimizer Buffer™ VI in 4ml de-ionized water). The standard protocol uses 35ml 1X Optimizer Buffer™ VI.

**CAUTION:** No amine-containing buffers, Tris or glycine, should be used. If proteins are in amine-containing buffers, then dialyze the samples against 1X Optimizer Buffer™ VI to completely remove the amines. For easy and convenient dialysis, use G-Biosciences' Tube-O-Dialyzer™ (supplied separately), visit [www.GBiosciences.com](http://www.GBiosciences.com).

**CAUTION:** The pH of all buffers used should be greater than pH 7.0 to prevent the release of toxic gas. **DO NOT ACIDIFY SOLUTIONS.**

**WARNING:** Sodium cyanoborohydride is toxic, open tubes and prepare solutions in a fume hood.

## PROTOCOL

**NOTE:** To ensure complete coupling, we advise using the recommended quantities of reagents, as these provide a molar excess of protein or peptide to the reactive groups on the agarose. Ensure complete and thorough mixing during shaking of the agarose to prevent beads from drying out.

1. Allow the reagents to reach room temperature.
2. Gently resuspend the HOOK™ Agarose by inverting the bottle 2-3 times. Once a homogeneous slurry is formed, remove 4ml of the slurry and place in a 15ml tube. Centrifuge at 2000xg for 30-60 seconds. Remove the storage buffer with a pipette.  
**NOTE:** The binding capacity of the agarose is ~20mg protein/ml resin.
3. Equilibration Step: Add 10ml 1X Optimizer Buffer™ VI to the HOOK™ Agarose and place on a rocker for 3-5 minutes. Centrifuge at 2000xg for 30-60 seconds. Remove the buffer with a pipette. Repeat Step 3 once.
4. Protein/Peptide Solution: Make 2ml solutions of either 1-20mg/ml protein or 1-2mg/ml peptide in 1X Optimizer Buffer™ VI, H<sub>2</sub>O or phosphate buffer saline. Add the protein/peptide solution and an additional 3ml 1X Optimizer Buffer™ VI to the HOOK™ agarose from step 3.

**NOTE:** For peptides insoluble in aqueous solutions, we recommend the use of DMSO (<30%).

**OPTIONAL:** Remove an aliquot of the protein solution to be used as 'starting material' to determine coupling efficiency, if required.

**CAUTION: THE FOLLOWING STEP-5 AND 6 MUST BE PERFORMED IN A FUME HOOD**

5. Add 50µl Sodium Cyanoborohydride (results in ~ 50mM NaCNBH<sub>3</sub>) to the HOOK™ Agarose/protein/peptide solution from step 4. Mix with gentle shaking for 6 hours.
6. After this step, the reaction may be removed from the hood. Further incubation may be carried out at 4°C overnight.
7. Centrifuge at 2000xg for 30-60 seconds. Remove the supernatant with a pipette and save for testing coupling efficiency or discard.

**NOTE:** *The coupling efficiency is determined by measuring and comparing the concentration of the protein/peptide in 'starting material' (Step 4) and in the supernatant (Step 7). The concentration may be determined either by UV absorbance at 280nm or by protein assay (CB-X assay, Cat # 786-12X).*

8. Wash the coupled agarose twice with 5ml of 1X Optimizer Buffer™ VI as in Step 3.
9. Blocking Step: Add 5ml Quenching Buffer to the agarose and gently shake for one hour at room temperature.
10. Centrifuge at 2000xg for 30-60 seconds and remove the buffer with a pipette. Wash twice with 5ml PBS or buffer of choice as described in Step 3.
11. Transfer the coupled agarose to the column and store at 4°C in PBS or buffer of choice with 0.05% sodium azide. The column is now ready for use.

## APPENDIX: GENERAL PROTOCOL FOR AFFINITY PURIFICATION OF PROTEIN

### **Sample Preparation**

Dissolve 1-20mg protein or 1-2mg peptide to be immobilized in 2-3mL of 1X PBS or other suitable buffer.

For proteins already in solution either:

- Dilute sample 4-fold in 1X PBS or other suitable buffer
- Desalt or dialyze to buffer-exchange into 1X PBS or other suitable buffer

**NOTE:** *If the protein solution contains primary amines (e.g., Tris or glycine), these compounds must be thoroughly removed or they will compete with the intended protein-coupling reaction.*

### **Additional Materials Required**



- Binding/Wash Buffer: Phosphate Buffered Saline (PBS), Tris Buffered Saline (TBS) or other suitable buffer.
- Sample: See sample preparation above
- Elution Buffer: 0.1-0.2M glycine•HCl, pH 2.5-3.0
- Neutralization Buffer: 1M phosphate or 1M Tris; pH 9

### **Protocol**

1. Allow the affinity column to equilibrate to room temperature and allow the storage buffer to drain out.
2. Equilibrate column by adding 3 column volumes (CV) Binding/Wash Buffer and allowing it to drain from column.
3. Add Sample to column and allow it to flow into the resin bed. For samples <1CV, we recommend sealing the bottom of the column after the sample as entered the resin and incubating for up to 1 hour. For samples >1CV, seal the column after the addition of each CV and incubate to maximize binding.
4. Unseal the column and wash the column with 6CV of Binding/Wash Buffer.
5. Elute the bound protein by applying 4CV of Elution Buffer. Collect 0.5-1ml fractions. Monitor elution by absorbance at 280nm or a suitable protein assay. Pool fractions of interest and exchange into an appropriate storage buffer by desalting or dialysis.
6. The pH of each fraction can be adjusted to neutral by adding 50µl of Neutralization Buffer per 1ml of eluate.
7. Monitor elution by absorbance at 280nm. Pool fractions of interest and exchange into an appropriate storage buffer by desalting or dialysis.

## RELATED PRODUCTS

Download our Sample Preparation and Protein Purification Handbooks.



<http://info.gbiosciences.com/complete-protein-sample-preparation-handbook/>

<http://info.gbiosciences.com/complete-protein-purification-handbook>

For other related products, visit our website at [www.GBiosciences.com](http://www.GBiosciences.com) or contact us.

Last saved: 7/2/2015 CMH



[www.GBiosciences.com](http://www.GBiosciences.com)