



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

INSTRUCTION MANUAL

ZymoBIOMICS™ DNA Miniprep Kit Catalog Nos. **D4300T, D4300 & D4304**

Highlights

- Rapid, robust, and simple purification of high quality, inhibitor-free DNA from any sample including feces, soil, water, biofilms, swabs, saliva, body fluids, etc.
- ZymoBIOMICS™ innovative lysis system enables efficient and unbiased lysis of microbes including Gram-positive and negative bacteria, fungi, protozoans, algae, and viruses.
- Unbiased extraction of ultra-pure DNA makes the ZymoBIOMICS™ DNA Miniprep Kit ideal for 16S rRNA gene sequencing, shotgun metagenomic sequencing, arrays, PCR and other sensitive applications.

Contents

Product Contents & Specifications.....	1
Product Description	2-3
Protocol	4-5
Appendices	
A. Sample Collection	6
B. Application Notes	7-9
C. Standardize Sample Preparation with ZymoBIOMICS™ Microbial Standards.....	10-11
D. Troubleshooting	11-13
Ordering Information.....	14

Satisfaction of all Zymo Research products is guaranteed. If you should be dissatisfied with this product please call 1-888-882-9682.

¹ For optimal performance, add beta-mercaptoethanol to 0.5% (v/v) i.e., 500 µl per 100 ml.

² This equates to approximately 2×10^9 bacterial cells, 2×10^8 yeast cells and 2×10^7 mammalian cells.

³ For water samples, filter using desired filter (not provided). Cut the filter into small pieces and place into ZR BashingBead™ Lysis Tube (0.1 & 0.5 mm). Alternatively up to 250 µl water can be processed directly.

⁴ DNA/RNA Shield™ provides an accurate molecular signature of the sample at the time of collection by preserving nucleic acids at ambient temperature and inactivating organisms including infectious agents. See Appendix A for more information.

Product Contents

ZymoBIOMICS™ DNA Miniprep Kit (Kit Size)	D4300T (5 Preps.)	D4300 (50 Preps.)	D4304 (50 Preps.)	Storage Temperature
ZR BashingBead™ Lysis Tubes (0.1 & 0.5 mm)	5	50	-	Room Temp.
ZymoBIOMICS™ Lysis Solution	4 ml	40 ml	-	Room Temp.
ZymoBIOMICS™ DNA Binding Buffer ¹	6 ml	100 ml	100 ml	Room Temp.
ZymoBIOMICS™ DNA Wash Buffer 1	2 ml	50 ml	50 ml	Room Temp.
ZymoBIOMICS™ DNA Wash Buffer 2	5 ml	60 ml	60 ml	Room Temp.
ZymoBIOMICS™ DNase/RNase Free Water	3 ml	10 ml x 3	10 ml x 3	Room Temp.
Zymo-Spin™ IV Spin Filters (Orange Tops)	5	50	50	Room Temp.
Zymo-Spin™ IV- HRC Spin Filters (Green Tops)	5	50	50	Room Temp.
Zymo-Spin™ IIIC- Z Columns	5	50	50	Room Temp.
Collection Tubes	20	200	200	Room Temp.
Instruction Manual	1	1	1	-

Note - Integrity of kit components is guaranteed for up to one year from date of purchase. Reagents are routinely tested on a lot-to-lot basis to ensure they provide maximal performance and reliability.

Specifications


- **Sample Sources** – Bacterial, fungal, protozoan, algae, viral, mitochondrial, and host DNA is efficiently isolated from ≤ 200 mg of mammalian feces, ≤ 250 mg soil, and 50 – 100 mg (wet weight) of fungal bacterial cells², biofilms and water³.
- **Bead Beating System** – ZymoBIOMICS™ innovative lysis system enables complete homogenization/disruption of the microbial cells walls and accurate microbial DNA analysis, free of bias. To ensure unbiased lysis, calibration of each bead-beating device is recommended by using the ZymoBIOMICS™ Microbial Community Standard (see Appendix C for details).
- **DNA Purity** – High quality, inhibitor-free DNA is eluted with ZymoBIOMICS™ DNase/RNase Free Water and is suitable for all downstream applications including PCR and Next-Generation Sequencing.
- **DNA Integrity** – On average, post bead beating, genomic DNA is between 15-20 kb depending on the initial quality of the sample making it amenable to Next-Generation Sequencing platforms requiring high molecular weight DNA. For optimal DNA integrity, collect samples in DNA/RNA Shield™⁴.
- **DNA Recovery** – Up to 25 µg total DNA can be eluted into 100 µl (50 µl minimum) ZymoBIOMICS™ DNase/RNase Free Water.
- **Bioburden** – A single preparation is guaranteed to contain less than 3 bacterial genomic copies per 1 µl of eluate as determined by quantitative amplification of the 16S rRNA gene when eluted using 100 µl water.
- **Equipment** – Microcentrifuge, vortex/Disruptor Genie®, high speed cell disrupter (recommended).

ZYMO RESEARCH CORP.


Product Description

The ZymoBIOMICS™ DNA Miniprep Kit is designed for purifying DNA from a wide array of sample inputs (e.g. feces, soil, water, and biofilms), that is immediately ready for microbiome or metagenome analyses. The ZymoBIOMICS™ innovative lysis system eliminates bias associated with unequal lysis efficiencies of different organisms (e.g. Gram-negative/positive bacteria, fungi, protozoans, and algae)¹ making it ideal for microbial community profiling. Unbiased mechanical lysis of tough microbes is achieved by bead beating with the innovative ultra-high density BashingBeads™ and validated using the ZymoBIOMICS™ Microbial Community Standard²; as shown in Figure 4. In addition, the ZymoBIOMICS™ DNA Miniprep Kit is equipped with Zymo Research's Proprietary OneStep™ PCR Inhibitor Removal technology enabling PCR from the most PCR prohibitive environmental samples rich in humic and fulvic acids, tannins, melanin, and other polyphenolic compounds. Coupling state-of-the-art lysis technology with Zymo-Spin™ technology results in superior yields of ultra-pure DNA ideal for all downstream applications including PCR, arrays, 16S rRNA gene sequencing, and shotgun sequencing.

Innovation. Pure & Simple.™



Linear recovery with unparalleled sensitivity



Ultra-pure DNA from Inhibitor Rich Samples




Figure 1. The ZymoBIOMICS™ -96 Magbead DNA Kit produces linear recovery of DNA for sensitive applications, detecting pathogenic organisms such as *E. coli* and *H. pylori* in assays with up to a 1000x dilution factor. A dilution series was created using stool infected with 1 x 10⁶ shiga producing *E. coli* cells and stool infected with *H. pylori* cells. These stool samples were then extracted using the ZymoBIOMICS™ -96 Magbead DNA Kit, showing effective purification and qPCR amplification, even at 1000:1 dilution.

Figure 2. The ZymoBIOMICS™ DNA Miniprep Kit provides inhibitor-free DNA even when challenged with extremely inhibitor rich samples. Real-time PCR was used to evaluate eluates recovered using the ZymoBIOMICS™ DNA Miniprep Kit, and kits from Suppliers M, P, and Q. Reaction volumes consisted of either 10% or 35% of the eluate from each kit to detect the presence of PCR inhibitors. Each reaction contained 25 ng of *Brettanomyces* DNA. Delayed and/or no amplification indicates PCR inhibition from inefficient inhibitor removal.




Figure 3. The ZymoBIOMICS™ DNA Miniprep Kit provides superior yields when compared to Suppliers M, P, and Q. 80 mg of feces was processed using each kit according to the manufactures' recommended protocol. DNA was eluted using 100 µl ZymoBIOMICS™ DNase/RNase Free Water. 6 µl of each sample was analyzed in a 1.0% (w/v) agarose/ethidium bromide gel. Samples were processed in triplicate. L is a 1Kb ladder.

¹ Chemical, enzymatic, and inferior lysis matrices (beads) lead to unrealistic representation of organisms in downstream metagenomic analyses that is not reflective of actual abundance. To learn more about this topic see Figure 3.

² For more information on the ZymoBIOMICS™ Microbial Community Standard (D6300) & ZymoBIOMICS™ Microbial Community DNA Standard (D6305) see Appendix C.

³ DNA is predominately 15-20 kb and amenable to Next-Generation Sequencing techniques requiring high molecular weight DNA.

Zymo Research offers a full suite of ZymoBIOMICS™ Services for reliable, accurate microbial and metagenomic analyses.


Services include: Microbial Composition Profiling, Novel Microbe Identification, and Customizable Bioinformatics.

For details visit us at: <http://www.zymoresearch.com/services/metagenomics>

Or contact us at: services@zymoresearch.com

ZYMO RESEARCH CORP.

A) Bias Free Microbial DNA Extraction Using ZymoBIOMICS™ DNA Mini Kit Validated with the ZymoBIOMICS™ Microbial Community Standard



B) Bias Free Microbial DNA Extraction Using ZymoBIOMICS™ DNA Mini Kit From Human Stool




Figure 4. A) The ZymoBIOMICS™ DNA Miniprep Kit provides unbiased representation of the organisms extracted from the ZymoBIOMICS™ Microbial Community Standard. DNA was extracted from ZymoBIOMICS™ Microbial Community Standard using four different DNA extraction methods (ZymoBIOMICS™ DNA Miniprep Kit, Human Microbiome Project Protocol, Supplier M, and Supplier Q) and analyzed using 16S rRNA gene sequencing. 16S rRNA genes were amplified with primers targeting v3-4 region and the amplicons were sequenced on Illumina® MiSeq™ (2x250bp). Overlapping paired-end reads were assembled into complete amplicon sequences. The composition profile was determined based on sequence counts after mapping amplicon sequences to the known 16S rRNA genes of the eight different bacterial species.

B) The ZymoBIOMICS™ DNA Miniprep Kit reliably isolates DNA from even the toughest to lyse gram positive organisms, enabling unbiased analyses of microbial community compositions. There is a significant increase in yield and Gram-positive bacterial abundance when DNA was isolated using the ZymoBIOMICS™ DNA Miniprep Kit. Correlated with the results in Figure 3A it can be concluded that unbiased DNA isolation was achieved. DNA was extracted from 200 µl of human feces suspended in PBS (10 % m/v) using four different DNA extraction methods (ZymoBIOMICS™ DNA Miniprep Kit, Human Microbiome Project Protocol, Supplier M, and Supplier Q) and analyzed using 16S rRNA gene sequencing. 16S rRNA genes were amplified with primers targeting v3-4 region and the amplicons were sequenced on Illumina® MiSeq™ (2x250bp). Overlapping paired-end reads were assembled into complete amplicon sequences. Amplicon sequences were profiled with Qiime using Greengenes 16S rRNA gene database (gg_13_8).

Protocol

Before starting: For optimal performance, add beta-mercaptoethanol (user supplied) to the ZymoBIOMICS™ DNA Binding Buffer to a final dilution of 0.5% (v/v) *i.e.*, 500 µl per 100 ml.

For customers using **DNA/RNA Shield Lysis Tube (Microbe) [R1103]** proceed to Step 2 without adding Lysis Solution.

1. Add sample to a **ZR BashingBead™ Lysis Tubes (0.1 & 0.5 mm)**. Add 750 µl **ZymoBIOMICS™ Lysis Solution** to the tube and cap tightly.

Sample Type	Maximum Input
Feces	200 mg
Soil	250 mg
Liquid Samples ¹ and Swab Collections ²	250 µl
Cells (Suspended in DNA/RNA Shield™ or isotonic buffer, <i>e.g.</i> PBS)	50-100 mg (wet weight) (2 x 10 ⁹ bacterial, 2 x 10 ⁸ yeast cells, 2 x 10 ⁷ mammalian cells)
Samples in DNA/RNA Shield™ (10% v/v Sample) ³	250 µl

2. Secure in a bead beater fitted with a 2 ml tube holder assembly and process at maximum speed for ≥ 5 minutes.

Processing time will vary based on sample input and bead beater. Times may be as little as 5 minutes when using high-speed cell disrupters (FastPrep® -24) or as long as 20 minutes when using lower speeds (e.g., Disruptor Genie®).⁴

3. Centrifuge the **ZR BashingBead™ Lysis Tubes (0.1 & 0.5 mm)** in a microcentrifuge at ≥ 10,000 x *g* for 1 minute.
4. Snap off the base of the **Zymo-Spin™ IV Spin Filter (Orange Top)** and place in a **Collection Tube**. Transfer up to 400 µl supernatant to the Zymo-Spin™ IV Spin Filter and centrifuge at 8,000 x *g* for 1 minute.
5. Add 1,200 µl of **ZymoBIOMICS™ DNA Binding Buffer** to the filtrate in the Collection Tube from Step 4.
6. Transfer 800 µl of the mixture from Step 5 to a **Zymo-Spin™ IIIC-Z Column** in a Collection Tube and centrifuge at 10,000 x *g* for 1 minute.
7. Discard the flow through from the Collection Tube and repeat Step 6.
8. Add 400 µl **ZymoBIOMICS™ DNA Wash Buffer 1** to the Zymo-Spin™ IIIC-Z Column in a new Collection Tube and centrifuge at 10,000 x *g* for 1 minute. Discard the flow-through.
9. Add 700 µl **ZymoBIOMICS™ DNA Wash Buffer 2** to the Zymo-Spin™ IIIC-Z Column in a Collection Tube and centrifuge at 10,000 x *g* for 1 minute. Discard the flow-through.

For **Technical Assistance:**
1-888-882-9682 or E-mail
tech@zymoresearch.com

¹For water samples, filter using desired filter (not provided). Cut the filter into small pieces and place into ZR BashingBead™ Lysis Tubes (0.1 & 0.5 mm).

²Swabs can also be cut or broken and placed directly in bead beating tube. For more information on processing swab samples, see Appendix B.

³ See Appendix A for additional information on sample collection in DNA/RNA Shield™.

⁴For optimal lysis efficiency and unbiased profiling all bead beater devices beyond those validated by Zymo Research should be calibrated using the ZymoBIOMICS™ Microbial Community Standard. See Appendix C.

For **Technical Assistance:**
1-888-882-9682 or E-mail
tech@zymoresearch.com

⁵ See Appendix D for
additional elution
information.

⁶ In some cases a brown-
colored pellet may form at
the bottom of the tube after
centrifugation. Avoid this
pellet when collecting the
eluted DNA.

⁷ If fungi or bacterial cultures
were processed; the DNA is
now suitable for all
downstream applications.

⁸ For time savings, skip Step
12b and proceed to step 13.
This alternative method
does not affect downstream
performance, but can affect
accurate quantification using
spectrophotometry.

10. Add 200 μl **ZymoBIOMICS™ DNA Wash Buffer 2** to the Zymo-Spin™ IIC-Z Column and centrifuge at 10,000 x g for 1 minute.
11. Transfer the Zymo-Spin™ IIC-Z Column to a clean 1.5 ml microcentrifuge tube and add 100 μl ⁵ **ZymoBIOMICS™ DNase/RNase Free Water** directly to the column matrix and incubate for 1 minute. Centrifuge at 10,000 x g for 1 minute to elute the DNA^{6, 7}.
12. Preparing **Zymo-Spin™ IV-HRC Spin Filter (Green Top)**.⁸
 - a. Snap off the base of the Zymo-Spin™ IV-HRC Spin Filter (Green Top) and place into a clean Collection Tube. Centrifuge at 8,000 x g for 3 mins. Discard the flow-through.
 - b. Remove the cap and add 400 μl ZymoBIOMICS™ DNase/RNase Free Water to the Zymo-Spin™ IV-HRC Spin Filter. Loosely cap Zymo-Spin™ IV-HRC Spin Filter and centrifuge at 8,000 x g for 2 minute.
13. Transfer the eluted DNA (Step 11) to a prepared Zymo-Spin™ IV-HRC Spin Filter in a clean 1.5 ml microcentrifuge tube. Loosely cap the Zymo-Spin™ IV-HRC Spin Filter and centrifuge at exactly 8,000 x g for 1 minute.

Appendix A

Sample Collection

For high quality reproducible microbiomics data, **DNA/RNA Shield™** is recommended for sample collection to avoid bias or erroneous results due to compositional changes from nucleic acid degradation or microbial growth. DNA/RNA Shield™ provides an unbiased molecular snapshot of the sample at the time of collection by preserving nucleic acids at ambient temperature and inactivating organisms including infectious agents. Samples can be stored and transported easily and safely with DNA/RNA Shield™ and is ideal for applications such as PCR, 16S rRNA gene sequencing, and shotgun metagenomic sequencing. DNA/RNA Shield™ can preserve nucleic acids in nearly any sample including feces, soil, saliva, blood, and tissues.

DNA/RNA Shield™ - Lysis Tube (Microbe) – Simply add sample, seal and store at ambient temperature. The Lysis Tube is immediately ready for bead beating, thereby streamlining the collection to extraction transition. (Cat. No. **R1103**)

DNA/RNA Shield™ – Fecal Collection Tube – The collection device is specifically designed for easy collection and stabilization of feces. Includes a scoop built for collecting 1 gram of feces (or any other sample such as saliva or soil). (Cat. No. **R1101**)

DNA/RNA Shield™ – Swab Collection Tube – Easy collection of biological samples; swab has breakable tip to allow for easy sample collection and removes the need to dispose of a potentially biohazardous swab material. (Cat. No. **R1106 & R1107**)




Figure 5. A) Nucleic acids in stool are effectively stabilized in DNA/RNA Shield™ at room temperature. Graph shows spike-in DNA and RNA controls from stool purified at the indicated time points and analyzed by (RT)qPCR. Controls: HSV-1 and HIV (AcroMetrix™, Life Technologies).

B) Microbial composition of stool is unchanged after one month at ambient temperature with DNA/RNA Shield™. Stool samples suspended in DNA/RNA Shield™ and stored at room temperature were compared to stool without preservative for one month. They were sampled at the indicated time points and processed with ZymoBIOMICS™ DNA Miniprep Kit. The extracted DNA was then subjected to microbial composition profiling via 16S rRNA gene targeted sequencing. Graphs show both phylum composition (left) and genus composition (right). Samples stored with DNA/RNA Shield™ had a constant microbial composition while the samples stored without shifted dramatically.

C) Viruses, bacteria and yeast are effectively inactivated by DNA/RNA Shield™. Samples containing the infectious agent (viruses, bacteria, yeast) were treated with DNA/RNA Shield™ or mock (PBS) treated for 5 minutes. Titer (PFU) was subsequently determined by plaque assay. Validated by: Influenza A - D. Poole and Prof. A. Mehle, Department of Medical Microbiology and Immunology, University of Wisconsin, Madison; Ebola (Kikwit) - L. Avena and Dr. A. Griffiths, Department of Virology and Immunology, Texas Biomedical Research Institute; HSV-1/2.

Appendix B

Application Notes

DNA/RNA Shield™ Lysis Tubes (Microbe) (Cat. No. R1103)

Addition of ZymoBIOMICS™ Lysis Solution is not necessary for samples stored in DNA/RNA Shield™ Lysis Tubes (Microbe) and samples can be immediately processed via bead beating. Simply move to Step 2 of the protocol (page 4) and bead beat according to instructions provided. Proceed with the remaining protocol as written (page 4).

DNA Viruses

For unbiased metagenomic analysis of viruses, incorporating a Proteinase K digestion prior to bead beating is recommended. Add 5% v/v of Proteinase K (Cat. No. D3001-2-5) to the lysate after Step 2 (page 4) and incubate for 30 minutes at 55°C. Continue to Step 3 (page 4).

Cheese and Protein Rich Biofluids (e.g. Milk, Sputum, Saliva, Spinal Fluid, Blood, and Serum)

Substitute the following for Step 1 (page 4) in ZymoBIOMICS™ DNA Miniprep Kit protocol section:

1. Add 0.3-0.4 g of cheese or 200 µl of biofluid to the ZR BashingBead™ Lysis Tubes (0.1 & 0.5 mm). Add 650 µl of ZymoBIOMICS™ Lysis Solution.
2. Add 2% v/v of Proteinase K (cat. no. D3001-2-5) to the ZymoBIOMICS™ Lysis Tubes (0.1 & 0.5 mm) and incubate for 30 minutes at 55°C.
3. Continue on to Step 2 (page 4) in ZymoBIOMICS™ protocol for further lysis.

Plant Tissue (Leaves and other plant material)

Plant tissue such as leaves contain DNA sources within the host tissue that can overwhelm 16S rRNA gene targeted sequencing (from both mitochondria & chloroplast). Microbes must be removed from the plant material to exclude host tissue from the bead beating process.

Prior to Step 1 (Page 4), suspend plant tissue in PBS and gently sonicate with sonication bath for effective removal of microbes. Alternatively, place plant tissue in a submerging volume of PBS inside of a conical tube and vortex briefly. The plant tissue can then be removed and the microbes can be centrifuged at high speeds to concentrate. Alternatively, a filter can also be used to concentrate the microbes and water removal. The filter can subsequently be cut and placed directly into the ZR BashingBead™ Lysis Tubes (0.1 & 0.5 mm) for bead beating (Step 2).

Plant Root

Cut root into small pieces and place directly into ZR BashingBead™ Lysis Tubes (0.1 & 0.5 mm) with 750 µl of ZymoBIOMICS™ Lysis Buffer as stated in Step 2 (page 4). Lysis should be performed with a lower speed bead beating device (e.g. vortex adapter) to avoid the host tissue contamination. Proceed with the remaining protocol as written (page 4).

Water/Air Samples

Filter samples using desired filter (not provided) prior to Step 1 (page 4). Cut the filter into small pieces, place into ZR BashingBead™ Lysis Tubes (0.1 & 0.5 mm), and continue to Step 2.

Tissue and Insect Samples

Tissue and Insect samples can be processed three different ways, depending on the sample type and the equipment available. The recommendations are listed next to the options below:

(A) Proteinase K - *Tissue*

1. Add up to 15 mg of tissue to a 1.5 ml microcentrifuge tube, then add a solution of 95 µl water, 95 µl **Solid Tissue Buffer** (Cat. No. D4068-2-6) and 10 µl **Proteinase K** (Cat. No. D3001-2-5). Incubate for at least 1 hour at 55° C or until tissue clarifies (samples can be incubated overnight without affecting DNA quality).
2. Process the lysate by proceeding to Step 1, Page 4 (Liquid Sample).

(B) Bead beating - *Tissue and Insect*

1. Place up to 15 mg of tissue/insect sample in a **ZR BashingBead™ Lysis Tube (2.0 mm)** (Cat. No. S6003-50) with 750 µl of **ZymoBIOMICS™ Lysis Solution** (included).
2. Secure in a bead beater fitted with a 2 ml tube holder assembly and process at maximum speed for ≥ 5 minutes.

Processing time will vary based on sample input and bead beater. Times may be as little as 5 minutes when using high-speed cell disrupters (FastPrep® -24) or as long as 20 minutes when using lower speeds (e.g., Disruptor Genie®).

3. Transfer the entire lysate to the ZR BashingBead Lysis Tube (0.1 & 0.5 mm) and proceed to Step 2, Page 4.

(C) Mortar & Pestle - *Tissue and Insect*

1. Pre-homogenize up to 15 mg tissue/insect sample with a pestle and mortar while submersed in liquid nitrogen.
2. Proceed to Step 1, Page 4 and process the entire sample.

Samples Collected with Swabs

Place swab directly into the ZR BashingBead™ Lysis Tubes (0.1 & 0.5 mm) with 750 µl of ZymoBIOMICS™ Lysis Buffer. The swab can be cut at the height of the ZR BashingBead™ Lysis Tube (0.1 & 0.5 mm) and left inside for bead beating in Step 2 (page 4). Alternatively, vortex the swab in the ZR BashingBead™ Lysis Tube (0.1 & 0.5 mm) with the ZymoBIOMICS™ Lysis Solution for 30 seconds to transfer the microbes into solution. Remove the swab and proceed to bead beating in Step 2 (page 4).




Figure 6. Phylum composition of a simulated microbial community when bead beating was performed with and without the presence of a Puritan HydraFlock® sterile flocked collection device placed in a BeadBashing tube and processed at maximum speed (6.5 m/s) for 5 minutes. The extracted DNA was then subjected to microbial composition profiling via 16S rRNA gene targeted sequencing. Experiment was performed in technical duplicates

Urine

Pelleting of Cells by Centrifugation

1. Pellet the bacterial cells by centrifuging the urine at 15,000 x g for 10 minutes and then proceed to Step 4 below.

Alternatively

Separation of Cells by Centrifugation

2. Add 70 µl Urine Conditioning Buffer (Cat. No. D3061-1-140) for every 1 ml of urine and mix well by vortexing. Urine stabilized by the Urine Conditioning Buffer can be stored for up to 1 month at ambient temperature. When samples are ready to be processed, mix well by vortexing and proceed to Step 3.
3. Centrifuge at 3,000 x g for 15 minutes.
4. Without disturbing the pellet, slowly decant or pipette out the supernatant, leaving behind 100 – 400 µl of pellet.
5. Add ZymoBIOMICS™ Lysis Solution to a final volume of 800 µl and then transfer the mixture to a ZR BashingBead™ Lysis Tube (0.1 & 0.5 mm). Proceed with the remaining protocol as written (page 4), starting at Step 2.




Figure 7. Phylum composition of urine preserved in Urine Conditioning buffer™ (UCB™), which preserves the microbial composition of urine with simulated stool contamination for a month at room temperature. Urine with UCB™ added (Zymo Research, D3061-1-160) was stored at room temperature and analyzed over a month period. At the indicated time points (0 Days, 2 weeks, and 1 month), DNA was extracted using the ZymoBIOMICS™ DNA Miniprep Kit. The extracted DNA was then subjected to microbial composition profiling via 16S rRNA gene targeted sequencing. Experiment was performed in technical duplicates.

Appendix C

Standardize Sample Preparation with ZymoBIOMICS™ Microbial Community Standards

The **ZymoBIOMICS™ Microbial Community Standard (Cat. No. D6300)** is a mock microbial community of defined and well characterized composition making it the perfect control for all microbiome profiling and metagenomics analyses.

It is ideal for assessing bias of DNA extraction methods since it contains three easy-to-lyse Gram-negative bacteria (*e.g. Escherichia coli*), five tough-to-lyse Gram-positive bacteria (*e.g. Listeria monocytogenes*), and two tough-to-lyse yeasts (*e.g. Saccharomyces cerevisiae*).

Bead Beating Device Calibration Protocol:

Zymo Research suggests calibrating bead beating devices with the ZymoBIOMICS™ Microbial Community Standard in order to ensure bias free microbial extraction. For Disruptor Genie®, vortex adapters, and vortex lysis we suggest a time course ranging from 10-45 minutes with the vortex at maximum speed. For high speed cell disruptors such as the MP FastPrep -24® we suggest a time course at maximum speed with a range of 3-10 minutes. The resulting DNA should be evaluated by quantifying DNA yield and changes in microbial profile at each time point. The bead beating time that yields a profile that closely matches the theoretical composition should become standard operating procedure for the bead beating device.

ZymoBIOMICS™ Microbial Community DNA Standard (Cat. No. D6305) is a mixture of genomic DNA extracted from pure cultures of eight bacterial and two fungal strains. Genomic DNA from each culture was quantified before mixing. The ZymoBIOMICS™ Microbial Community Standard allows for assessment of bias from library preparation, sequencing, and bioinformatics analysis.

It serves perfectly as a microbial standard for benchmarking the performance of microbiomics or metagenomics analyses, including those provided by a 3rd party.

Figure 8. Accurate composition for reliable use to evaluate shotgun seq. and 16S rRNA seq.




Figure 8. Characterization of the microbial composition of the two ZymoBIOMICS™ standards with shotgun metagenomic sequencing (left panel) and 16S rRNA gene targeted sequencing (right panel). The measured composition of the two standards agrees with the theoretical/designed composition. "DNA Standard" represents ZymoBIOMICS™ Microbial Community DNA Standard (DNA version) and "Microbial Standard" represents ZymoBIOMICS™ Microbial Community Standard (cellular version). Genomic DNA composition by shotgun sequencing was calculated based on counting the amounts of raw reads mapped to each genome. 16S composition by 16S rRNA gene targeted sequencing was calculated based on counting the amount of 16S raw reads mapped to each genomes.

Figure 2. A) Use ZymoBIOMICS™ Microbial Standards for assessing GC-Bias in Shotgun Metagenomics




Figure 2. B) Perfect for tracking PCR Chimera in 16S rRNA Gene Sequencing




Figure 9.

A) Library preparation for shotgun metagenomic sequencing was performed in two different ways: one by supplier I and one by an in-house method. Shotgun sequencing was performed on Illumina® MiSeq™ with paired-end sequencing (2 x 150 bp). Raw reads were mapped to the 10 microbial genomes to evaluate the potential effect of GC content on sequencing coverage. Normalized coverage was calculated by normalization by the average sequencing coverage of each genome


B) PCR chimera increases with PCR cycle number in the library preparation process of 16S rRNA gene targeted sequencing. 20 ng ZymoBIOMICS™ Microbial Community Standard was used as a template. The PCR reaction was performed with ZymoBIOMICS™ PCR Premix and with primers that target v3-4 region of 16S rRNA gene. Chimera rate in percentage was determined with Uchime and using the 16S rRNA gene of the 8 bacterial strains in the standard as reference PCR.

Appendix D

Troubleshooting:

For **Technical Assistance**, please contact 1-888-882-9682 or E-mail tech@zymoresearch.com.

DNA Elution Guide



The Relationship between Elution Volume, DNA Yield, and DNA Concentration Using *Bacillus Subtilis* Culture. Using a smaller elution volume results in higher concentrations of DNA samples, but with reduced yields. Using a larger elution volume results in higher DNA yields, but at a reduced concentration. Choose an elution volume that best fits your individual application. Suggested elution volume is 100 µl.

Troubleshooting (Continued):

Problem	Possible Causes and Suggested Solutions
Background Contamination	<ul style="list-style-type: none"> • Clean workspace, centrifuge, and pipettes with 10% bleach to routinely to avoid contamination. • Use of kit in exposed environment without proper filtration. Check pipettes, pipette tips, microcentrifuge tubes, workspace, etc. for contamination. • Make sure bags of columns and buffer bottles are properly sealed for storage. Use of these outside a clean room or hood can result in contamination. <p><u>Lysis Methods</u></p> <ul style="list-style-type: none"> • When using a Disruptor Genie[®], vortex adapter, vortex, or similar processing times will vary. Suggested time is anywhere from 5-20 minutes. Calibrate bead beating times to your particular device and application by testing several different time points before using precious samples. (Suggested times to test: 10, 20, and 30 minutes.) See Appendix C for details. • When using FastPrep[™]-24 or similar devices run max speed for 5 minutes to ensure unbiased lysis. (6.5m/s on FastPrep[®]-24).
Low DNA Yield	<p><u>Incomplete Debris Removal</u></p> <ul style="list-style-type: none"> • For high density samples, ensure lysate is centrifuged properly to pellet insoluble debris following bead beating. Ensure that none of the debris is transferred to the Zymo-Spin[™] IV Spin Filter in the next step. <p><u>Input</u></p> <ul style="list-style-type: none"> • If the lysate does not pass through the column or is extremely viscous, use less input material. Too much sample input can cause cellular debris to overload the column and insufficient flow • Consult the Sample Type table on Page 4 for information on your particular input limit based on sample.

ZYMO RESEARCH CORP.

Binding Step

- Ensure that the ZymoBIOMICS™ DNA Binding Buffer is completely mixed with lysate before loading onto the column. Improperly mixed samples can lead to poor DNA recovery.

Low DNA Yield (Continued)

Elution Procedure

- Ensure the **ZymoBIOMICS™ DNase/RNase Free Water** hydrates the matrix for at least 1 minute before centrifugation.
- Make sure **Zymo-Spin™ IV–HRC** column cap is loosely tightened during the “wash” step. This can create a vacuum and cause lower elution volume.
- To increase yields, heat the **ZymoBIOMICS™ DNase/RNase Free Water** to 60°C before use. Additionally, users can reload the eluate onto the column matrix, incubate at room temperature for 3 minutes, and centrifuge again to increase yield without further dilution.

Ordering Information

Product Description	Catalog No.	Kit Size
ZymoBIOMICS™ DNA Microprep Kit	D4301	50 preps.
ZymoBIOMICS™ DNA Microprep Kit (Lysis Matrix Not Included)	D4305	50 preps.
ZymoBIOMICS™ DNA Miniprep Kit	D4300	50 preps.
ZymoBIOMICS™ DNA Miniprep Kit (Lysis Matrix Not Included)	D4304	50 preps.
ZymoBIOMICS™ -96 DNA Kit (Includes BashingBead™ Lysis Rack)	D4303	2x96 preps.
ZymoBIOMICS™ -96 DNA Kit (Includes BashingBead™ Lysis Tubes)	D4309	2x96 preps.
ZymoBIOMICS™ -96 DNA Kit (Lysis Matrix Not Included)	D4307	2x96 preps.
ZymoBIOMICS™ -96 Magbead DNA Kit (Includes BashingBead™ Lysis Rack)	D4302	2x96 preps.
ZymoBIOMICS™ -96 Magbead DNA Kit (Includes BashingBead™ Lysis Tubes)	D4308	2x96 preps.
ZymoBIOMICS™ -96 Magbead DNA Kit (Lysis Matrix Not Included)	D4306	2x96 preps.

For Individual Sale	Catalog No.	Amount
ZR BashingBead™ Lysis Tubes (0.1 & 0.5 mm)	S6012-50	50
ZymoBIOMICS™ Lysis Solution	D4300-1-40	40 ml
ZymoBIOMICS™ DNA Binding Buffer	D4300-2-100	100 ml
ZymoBIOMICS™ DNA Wash Buffer 1	D4300-3-50	50 ml
ZymoBIOMICS™ DNA Wash Buffer 2	D4300-4-60	60 ml
ZymoBIOMICS™ DNase/RNase Free Water	D4302-5-50	50 ml
Zymo-Spin™ IV Spin Filters (Orange Tops)	C1007-50	50
Zymo-Spin™ IV- HRC Spin Filters (Green Tops)	C1010-50	50
Zymo-Spin™ IIIC- Z Columns	C1006-50-G-ZB	50
Collection Tubes	C1001-50	50
	C1001-500	500
	C1001-1000	1,000

Sample Collection	Catalog No.	Amount
DNA/RNA Shield™ - Lysis Tube (Microbe)	R1103	50
DNA/RNA Shield™ – Fecal Collection Tube	R1101	10
DNA/RNA Shield™ – Swab and Collection Tube	R1106	10
	R1107	50
DNA/RNA Shield™	R1100-50	50 ml
	R1100-250	250 ml
DNA/RNA Shield™ (2X concentrate)	R1200-25	25 ml
	R1200-125	125 ml

ZymoBIOMICS™, Zymo-Spin™, and DNA/RNA Shield™ are trademarks of Zymo Research Corp. Disruptor Genie® is a registered trademark of Scientific Industries, Inc. FastPrep® is a registered trademark of Qbiogene, Inc. Illumina® MiSeq™, Illumina® Nextera® XT are trademarks or registered trademarks of Illumina Inc. AcroMetrix™ is a trademark of Thermo Fisher Scientific Inc. ZymoBIOMICS™ DNA Micro Kit is for research use only. ZymoBIOMICS™ DNA Micro Kit is not sold for use in diagnostic procedures. Reagents included with this kit are irritants. Follow the safety guidelines and rules enacted by your research institution or facility including the wearing of protective gloves and eye protection when using this kit.

ZYMO RESEARCH CORP.

Phone: (949) 679-1190 ▪ Toll Free: (888) 882-9682 ▪ Fax: (949) 266-9452 ▪ info@zymoresearch.com ▪ www.zymoresearch.com