

SpillEx

Super-absorbent. Super safe!



SpillEx

– the problem solver



Spillages on hard floors can cause considerable problems.

In cleanrooms and controlled environments, spills create slip hazards and are not easily removed. Spills usually require a bulky mop system or a large number of expensive wipes that may not absorb well. In Lab Animal Research facilities, animal urine on floors and benches require fast and sanitary removal.

Not to worry. There is a solution to spillage problems – SpillEx from Vileda Professional.

Why traditional floor cleaning doesn't work

When laundering microfiber mops that have been used to remove spills, you bring the contaminated mop in contact with up to 50 other mops. Not an ideal situation, even if you do a disinfection laundering process.

Similar problems occur when using a bucket and wringer. When dipping the used mop into the bucket, it will contaminate the clean water after the first rinse.

SpillEx takes care of the problem





SpillEx is a disposable, super-absorbing floor cloth that takes care of liquid spillages in minutes. One SpillEx absorbs and holds up to 1,200 ml of water and up to 500 ml of NaCl 0.9% (similar to urine) and it will not drip when lifting and moving, as SpillEx transforms the liquid spill into a gel.

SpillEx has a very low particle count. When tested against leading competitive products, SpillEx shows the lowest particle burden of items tested.


SpillEx is the ideal solution for when a spill occurs before the actual floor cleaning takes place. There is no risk of cross-contamination and no risk of soiling the bucket water with the spillage.

Easy as 1-2-3

1. Place SpillEx on the spillage with the absorbent unprinted side down.
2. Wait until the spillage is completely soaked up by SpillEx.
If the spillage area is bigger than the SpillEx, simply move and cover the remaining spillage area.
3. Pick up the used SpillEx and dispose of it.



Weight per cloth	1.6 oz/45 g
Size	20 x 15 in/51 x 37 cm
Typical spill uptake in application with tap water	41 oz/1,200 ml
Typical spill uptake in application with NaCl 0.9% (urine)	17 oz/500 ml



1.800.932.5000 | vwr.com

Prices and product details are current when published; subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | VWR, the VWR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks, of VWR International, LLC and its related companies. All other marks referenced are registered by their respective owner(s). | Visit vwr.com to view our privacy policy, trademark owners and additional disclaimers. ©2015 VWR International, LLC. All rights reserved.

0815 Lit. No. 061349W