



# HYDROGEN VS. HELIUM

VS.

# FOR GAS CHROMATOGRAPHY



## HYDROGEN

## HELIUM

The first element of the periodic table, hydrogen is a colorless, and insipid gas. It is the most flammable of the known substances and a powerful reductive agent, reacting with the oxides and chlorides of many metals.

Helium has many unique properties: low boiling point, low density, low solubility, high thermal conductivity and inertness, so it is used for any application which can exploit these properties.



### REACTIVE

Hydrogen reacts with a number of metals and non-metals and is a powerful reducer. It is highly combustible.



### INERT

Doesn't react with other elements and is non-combustible.



### NON-TOXIC

Although non-toxic, hydrogen can be absorbed through inhalation and cause oxygen deficiency in high concentrations.



### NON-TOXIC

Safe to use in a variety of applications.



### ATOMIC MASS

Low atomic mass: 1.007825 g·mol<sup>-1</sup>  
Low boiling point (1/c) - 252.8°C.



### LOWEST BOILING POINT - 268.9°C

Liquid at ultra-cool temp.



### SMALL MOLECULAR SIZE

Vanderwaals radius: 0.12 nm.




### SMALL MOLECULAR SIZE

Vanderwaals radius: 0.118 nm.

## SELECTING A GC CARRIER GAS

Historically, most gas chromatography labs have used helium as a carrier gas. Helium offers the benefits of being non-combustible, inert, and provides moderate speed of analysis. In contrast, Hydrogen, is flammable, and may be reactive under specific conditions. However, it produces high-speed analysis and generates sharper peak shapes. For those researchers who value time and increased sensitivity, hydrogen is a smart alternative. Further, given that helium is a limited resource and much more expensive than hydrogen, many labs are switching to hydrogen carrier gas.



## HYDROGEN VS. HELIUM

### ADVANTAGES

- Affordable
- Provides most time efficient separation
- Sharper peak shapes
- Efficient at high gas velocities (i.e., 60 cm/sec)
- Easily produced

### ADVANTAGES

- Very inert, will not react with analytes
- Non-flammable
- Provides time efficient separations

### DISADVANTAGES

- Can form explosive mixture with air
- Reactive under some circumstances

### DISADVANTAGES

- Expensive
- Non-replenishable resource

## 3 MYTHS OF HYDROGEN GAS USAGE

### HYDROGEN IS TOO DANGEROUS

- In spite of the fact that hydrogen has a history of unpleasant events in many labs being both flammable and explosive in air, a hazardous mixture can be easily avoided under standard lab conditions.

#### FACTS

- Being the lightest of the gasses, hydrogen rises quickly (45 miles/hr) and quickly dissipated in the lab environment.
- Hydrogen generators produce gas in an on demand and store less than 100 ml of hydrogen at low pressure—far less than a typical cylinder which holds 50 L at 200 atm.
- Flow rates from a hydrogen generator are controlled and far below the lower explosive limit for hydrogen in air when released in the oven of an average GC.
- Hydrogen generators are equipped with built-in leak detection sensors that shut down the system if a leak is detected.
- Most labs are already using hydrogen for a number of purposes, it is the fuel gas for FID detectors and is therefore already used in most GC labs.

### HYDROGEN IS TOO REACTIVE FOR MANY GC APPLICATIONS

- Some reducing chromatographers are deterred from using hydrogen as it is a reducing agent that can promote deterioration.

#### FACTS

- Hydrogenation is favored only at high temperatures and high pressures, and in the presence of a metallic catalyst such as nickel, platinum, or palladium
- Conditions for hydrogenation do not exist within regularly used open tubular fused silica columns
- Some precautions should be taken when using nickel or alumina (Al<sub>2</sub>O<sub>3</sub>) columns

### METHODS WON'T BE EASILY TRANSLATED FROM HELIUM TO HYDROGEN

- Changing carrier gas may involve some method redevelopment—however, this doesn't necessarily have to be difficult.

#### FACTS

- Not all methods have to be redeveloped. The likelihood that a method will need redevelopment depends on the complexity of the chromatogram.
- Method development software is available that can assist chromatographers when changing carrier gasses.
- For labs using regulated methods, consult the regulating agency—many new methods are available that utilize hydrogen as a carrier gas.

## BENEFITS OF HYDROGEN GENERATORS



### 1 SAFETY

Eliminates the risks associated with high pressure cylinders, which can cause significant injury if dropped, or physical injuries while moving or lifting heavy canisters.

### 2 CONSISTENCY

Gas purity may vary slightly from canister to canister impacting testing results. Further contamination can occur from piping joints, sealants, and leaks.

### 3 CONVENIENCE

Hydrogen generators eliminate inconvenient canister switches that interrupt work and cost time. Additional time spent ordering, transferring, inventorying, and monitoring canister stock can also be recovered.

### 4 COST-SAVINGS

As helium prices increase due to limited supply, switching to affordable hydrogen offers immediate savings. Further by producing hydrogen with a gas generator, indirect costs associated with purchasing, shipping and storing gas canisters can also be achieved.

### 5 IMPROVED RESULTS

Hydrogen offers faster separations with better peak resolution. Significantly decreasing analysis time greatly improved both throughput and productivity—in many cases results can be obtained in half the time.