

## PureCube Rho1D4 MagBeads

| Product  | Catalog No. | Package size  |
|--|-------------|---|
| PureCube Rho1D4 MagBeads (1 mL)  | 33201       | 1 x 1 mL  |
| PureCube Rho1D4 MagBeads (5 mL)  | 33205       | 1 x 5 mL  |
| PureCube Rho1D4 MagBeads (25 mL)   | 33225       | 1 x 25 mL   |
| Rho Starter Set 2: PureCube Rho1D4 MagBeads (1 mL) + Rho1D4 peptide (5 mg) | 33299       | 1 mL Rho1D4 magnetic beads +<br>1 x 5 mg Rho1D4 peptide |

### Product Description

PureCube Rho1D4 MagBeads were developed for the affinity purification of proteins with the rho1D4-tag (protein sequence TETSQVAPA), and are compatible with all prokaryotic and eukaryotic expression systems. The affinity matrix is based on spherical magnetic agarose beads, consisting of 6% cross-linked agarose. The material is highly porous to allow optimal protein interaction. Cross-linked agarose is also physically very stable, making it suitable for purification processes without deformation or destruction. Our magnetic beads are very homogeneous in size with a medium particle diameter of 30 µm, yielding a high degree of reproducibility between individual purification runs. The Rho1D4 antibody is coupled to the agarose resin in such a way as to obtain a matrix with highest binding capacity towards rho-tagged proteins and enhanced storage stability.

PureCube Rho1D4 MagBeads are delivered on ice as a 25% suspension. Therefore, 1 mL suspension will yield a 50 µL bed volume. The suspension contains 20% ethanol to prevent microbial growth.

### Protein Binding Capacity

The protein binding capacity is up to 3 mg protein per mL of settled beads, as determined by purification of a 35 kDa rho1D4-tagged membrane protein, and quantified via spectrophotometry.

### Compatibility

PureCube Rho1D4 MagBeads are compatible with low concentrations of most commonly used detergents.

### Shipping & Storage

|                      |  |
|----------------------|--|
| Shipment Temperature | Shipped on ice                         |
| Short-term Storage   | In equilibration buffer (see protocol) |
| Long-term Storage    | In 20% ethanol at 4 °C                 |

## **Additional Information**

For protein purification and cleaning protocols, as well as guidelines for fusing proteins to the rho1D4 tag, please visit our webpage at: [www.cube-biotech.com/protocols](http://www.cube-biotech.com/protocols). For purification of rho-tagged proteins with gravity flow columns and low pressure chromatography, we recommend using PureCube Rho1D4 Agarose. For affinity purification of GST-tagged, his-tagged or strep<sup>®</sup>-tagged proteins, Cube Biotech offers dedicated agarose resins, magnetic beads and prepacked cartridges. Also available are a range of ultrapure detergents and buffers for extraction and purification of proteins. See [www.cube-biotech.com/products](http://www.cube-biotech.com/products) for details.

Disclaimer: Our products are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.

Trademarks: Strep-tag<sup>®</sup> (IBA GmbH).