

Instruction Manual

Product	Size	Catalog Number
Keratinocyte Growth Medium 2 (Ready-to-use)	500 ml	C-20011
Keratinocyte Growth Medium 2 Kit	500 ml	C-20111
Keratinocyte Basal Medium 2	500 ml	C-20211
Keratinocyte Basal Medium 2, phenol red-free	500 ml	C-20216
Keratinocyte Growth Medium 2 SupplementMix	for 500 ml	C-39016
Keratinocyte Growth Medium 2 SupplementPack	for 500 ml	C-39011

Recommended for

- Normal Human Epidermal Keratinocytes (NHEK), juvenile foreskin
- Normal Human Epidermal Keratinocytes (NHEK), adult donor

Product Description

PromoCell Keratinocyte Growth Medium 2 is a serum-free medium that was developed for the *in vitro* cultivation of epidermal keratinocytes without feeder cells. The medium is optimized for primary human cells, but can also be used for horse, porcine, as well as mouse and rat keratinocytes when the calcium chloride (CaCl₂) concentration is adapted. In addition, it is suitable for rat corneal keratinocytes.

PromoCell Keratinocyte Growth Medium 2 is available as Medium (Ready-to-use) or as Medium Kit.

The Medium (Ready-to-use) consists of a 500 ml bottle of Basal Medium, one vial of CaCl₂ stock solution, and one

vial of SupplementMix. The Medium Kit consists of a 500 ml bottle of Basal Medium and the SupplementPack (a set of individual vials with pre-measured supplements) allowing the user full control over the media formulation.

Adding the Supplement solutions to the Basal Medium results in the complete Growth Medium.

Basal Medium (with or without phenol red), SupplementMix as well as SupplementPack can also be purchased separately.

Supplementation Details

PromoCell Keratinocyte Growth Medium 2 contains all the growth factors and supplements necessary for the optimal growth of human keratinocytes, e.g. Bovine Pituitary Extract, Epidermal Growth Factor, Insulin, and Transferrin (for details see table on page 2). PromoCell Keratinocyte Growth Medium 2 does not contain antibiotics or antimycotics and is formulated for use in an incubator with an atmosphere of 5% CO₂.

Note: CaCl₂ is provided as separate stock solution. This allows users to adjust the CaCl₂ concentration individually. To ensure the maximal growth performance of normal human Keratinocytes, PromoCell re-recommends adding the entire 60 µl CaCl₂ stock solution, which results in final 0.06 mM CaCl₂ concentration.

In case you need a higher CaCl₂ concentration (e.g. differentiation experiments), a 0.5 M CaCl₂-Solution can be ordered separately (see Related Products).

Preparation of the supplemented Medium for Use

Thaw the SupplementMix or SupplementPack at 15 to 25°C. Aseptically mix the supplement solutions by carefully pipetting up and down. Then, transfer the entire content of each supplement to the Basal Medium. Close the bottle and swirl gently until a homogenous mixture is formed.

Storage and Stability

Store the Basal Medium at 4 to 8°C in the dark and the SupplementMix, the CaCl₂ solution, or the SupplementPack at -20°C immediately after arrival. Do not freeze the Basal Medium. If stored properly, the products are stable until the expiry date stated on the label. After adding the supplements to the Basal Medium, the shelf life of the complete medium is 6 weeks at 4 to 8°C. Do not freeze the complete medium.

For use, pre-warm only an aliquot of the complete medium and keep the remaining medium refrigerated at 4 to 8°C.

Note: The SupplementMix, the CaCl₂ solution, and the SupplementPack are delivered thawed and can be frozen after arrival without losing any activity.

Quality Control

All lots of PromoCell Keratinocyte Growth Medium 2 are subjected to comprehensive quality control tests using primary human keratinocytes. Each lot is checked for growth promoting activity, adherence rate, and typical morphology of the tested keratinocytes. Approved in-house lots of media are used as a reference.

In addition, all lots of media have been tested for the absence of microbial contaminants (fungi, bacteria, mycoplasma).

Intended Use

The products are for *in vitro* use only and not for diagnostic or therapeutic proce-

dures. For safety precautions please see appropriate MSDS.

Note: Due to their low serum content or the absence of serum, PromoCell media are not suitable for trypsin neutralization (e.g. when splitting the cells). Instead we recommend using our DetachKit (C-41200, C-41210, C-41220), which contains HEPES BSS, Trypsin/EDTA and Trypsin Neutralizing Solution.

Related Products

- CaCl₂-Solution (0.5 M) (C-34006)

Final supplement concentrations (after addition to the medium)

■ Bovine Pituitary Extract	0.004 ml / ml
■ Epidermal Growth Factor (recombinant human)	0.125 ng / ml
■ Insulin (recombinant human)	5 µg / ml
■ Hydrocortisone	0.33 µg / ml
■ Epinephrine	0.39 µg / ml
■ Transferrin, holo (human)	10 µg / ml
■ CaCl ₂	0.06 mM

If you require special media modifications, we offer a Custom Media Service starting at 10 bottles per order. Please ask for details.

PromoCell GmbH

Sickingenstr. 63/65
69126 Heidelberg
Germany

Email: info@promocell.com
www.promocell.com

USA/Canada

Phone: 1 – 866 – 251 – 2860 (toll free)
Fax: 1 – 866 – 827 – 9219 (toll free)

Deutschland

Telefon: 0800 – 776 66 23 (gebührenfrei)
Fax: 0800 – 100 83 06 (gebührenfrei)

France

Téléphone: 0800 90 93 32 (ligne verte)
Téléfax: 0800 90 27 36 (ligne verte)

United Kingdom

Phone: 0800 – 96 03 33 (toll free)
Fax: 0800 – 169 85 54 (toll free)

Other Countries

Phone: +49 6221 – 649 34 0
Fax: +49 6221 – 649 34 40