



G-Biosciences ♦ 1-800-628-7730 ♦ 1-314-991-6034 ♦ [technical@GBiosciences.com](mailto:technical@GBiosciences.com)

A Geno Technology, Inc. (USA) brand name

# **XIT™ Genomic DNA from FFPE Tissue**

**For the isolation of genomic DNA from formalin fixed,  
paraffin embedded tissue**

**(Cat. #786-290)**



think proteins! think G-Biosciences [www.GBiosciences.com](http://www.GBiosciences.com)

Introduction..... 3

Items Supplied ..... 3

Storage Conditions ..... 3

ITEMS NEEDED BUT NOT SUPPLIED..... 3

Preparation Before Use ..... 3

PROTOCOL FOR FFPE FIXED TISSUE ..... 4

Related Products ..... 6

## INTRODUCTION

The *XIT™* Genomic DNA kit is designed for the isolation of genomic DNA from formalin fixed, paraffin embedded tissue. The *XIT™* kit uses solvent extraction, cell lysis, protein digestion and precipitation and finally DNA precipitation to isolate high quality genomic DNA.

*XIT™* Genomic DNA from FFPE Tissue kit is offered for the processing of a maximum of 0.25g of tissue. The purified DNA has a  $A_{260}/A_{280}$  ratio between 1.7 and 1.9, and is up to 200kb in size. The yield is 0.5-10µg per mg solid tissue.

## ITEM(S) SUPPLIED (Cat. # 786-290)

Description	Size
XIT™ Lysis Buffer	10ml
LongLife™ Proteinase K	0.5ml
XIT™ Protein Precipitation Buffer	2.5ml
Mussel Glycogen Solution	50µl
TE Buffer	1.5ml
LongLife™ RNase	0.5ml

## STORAGE CONDITIONS

The kit is shipped at ambient temperature. Upon arrival, store the *LongLife™* Proteinase K and *LongLife™* RNase at -20°C, all other kit components can be stored at room temperature. The kit components are stable for 1 year, if stored properly.

## ITEMS NEEDED BUT NOT SUPPLIED

Isopropanol, 70% ethanol, xylene.

## PREPARATION BEFORE USE

1. Read appropriate protocol and preheat waterbaths or heating blocks to appropriate temperatures.
2. Equilibrate TE Buffer to 50-60°C.


## PROTOCOL FOR FFPE FIXED TISSUE

1. Final chop <10mg formaldehyde fixed paraffin embedded (FFPE) tissue and transfer to a 1.5ml centrifuge tube.
2. Transfer 400µl xylene to the tube and incubate at room temperature with gentle shaking for 5 minutes.  
*NOTE: Wear gloves, safety goggles and lab coat when using xylene.*
3. Centrifuge the tube at 14,000g for 3 minutes to pellet the tissue. Carefully discard the supernatant.
4. Repeat steps 2 and 3 two more times.
5. Resuspend the tissue in 400µl 90% ethanol and incubate at room temperature with gentle shaking for 5 minutes.
6. Centrifuge the tube at 14,000g for 3 minutes to pellet the tissue. Carefully discard the supernatant.
7. Repeat steps 5 and 6.
8. Transfer 400µl *XIT*<sup>™</sup> Lysis Buffer to the tissue. Homogenize the sample until a homogeneous solution is obtained.  
*NOTE: For efficient grinding, we recommend G-Biosciences' EZ-Grind<sup>™</sup> (Cat. # 786-139), a high efficient grinding resin with matching pestle and tubes.*
9. Add 10µl *LongLife*<sup>™</sup> Proteinase K to the tube and mix by inverting the tube 20 times. Incubate at 55°C overnight for maximal yield. Invert the tube periodically during the incubation.
10. If tissue is not completely digested, add a further 10µl *LongLife*<sup>™</sup> Proteinase K and incubate at 55°C for 3 hours. Invert the tube periodically during the incubation.
11. Add 90µl *XIT*<sup>™</sup> Protein Precipitation Buffer to the sample and mix by inverting the tube 10-20 times.
12. Centrifuge at 14,000g for 5 minutes. Carefully, transfer the supernatant to a fresh tube.  
*NOTE: The precipitated protein should form a tight white pellet. If not, incubate the sample on ice for 5 minutes and repeat the centrifugation.*
13. Add 400µl isopropanol to the supernatant and mix by gently inverting the sample 30-50 times.  
*NOTE: If DNA concentrations is expected to be low (<10µg), add 1µl Mussel Glycogen Solution.*
14. Centrifuge at 14,000g for 5 minutes.
15. Discard the supernatant and use a pipette to carefully remove excess liquid.
16. Add 200µl 70% ethanol and invert the tube twice to wash the pellet.
17. Centrifuge at 14,000g for 2 minutes.
18. Discard the supernatant and drain the tube on a piece of clean absorbent paper. Allow to air dry for 15 minutes.
19. Add 50µl prewarmed TE buffer and 1µl *LongLife*<sup>™</sup> RNase to remove the RNA (if required).

20. Rehydrate the genomic DNA by incubating at 55-65°C for one hour, followed by an overnight incubation at room temperature to ensure complete genomic DNA hydration.
21. Store DNA at 4°C, for long term storage store at -20 or -80°C

## RELATED PRODUCTS

Download our Molecular Biology Handbook.



<http://info2.gbiosciences.com/complete-molecular-biology-handbook>

For other related products, visit our website at [www.GBiosciences.com](http://www.GBiosciences.com) or contact us.

Last saved: 3/24/2014 TNN





**[www.GBiosciences.com](http://www.GBiosciences.com)**