

Analysis of Acrylamide in Potato Chips by SPE and GC-MS

Anila I. Khan, Thermo Scientific, Runcorn, UK

Key Words

Hypercarb SPE, food, acrylamide, 2-propenamide, capillary GC, porous graphitic carbon (PGC), polyethylene glycol (PEG) GC column

Abstract


Acrylamide is an endogenous compound, formed when heating starchy or sugary foods. The production of potato chips can result in its formation. The method reported here detects acrylamide at the low ng/g levels at which it is produced. Potato chips were extracted using porous graphitic carbon for solid phase extraction (SPE). Analysis of acrylamide was performed using GC-MS on a polyethylene glycol phase GC column. A standard addition calibration curve was used to estimate the level of acrylamide in potato chips at 450 ng/g.

Introduction

Acrylamide (2-propenamide) is a potential human carcinogen. This toxic compound is usually formed as a by-product of Maillard reactions during the heating of carbohydrate-rich food. The World Health Organization (WHO) has set a safe limit of 500 ng/mL acrylamide in drinking water. Higher levels of 100–1000 ng/g are determined in some foods such as potato chips or french fries.

The extraction of acrylamide from potato chips is carried out using a Thermo Scientific™ HyperSep™ Hypercarb™ SPE cartridge. Hypercarb SPE material is 100% porous graphitic carbon (PGC) and offers retention of highly polar compounds that are not usually retained by traditional reversed phase C18 columns. HyperSep Hypercarb SPE can produce clean samples by removing potential matrix interferences.

The analysis of acrylamide was carried out using a GC-MS in electron ionization (EI) mode. Quantitative measurement in food can be difficult as matrix-derived ions can interfere with acrylamide fragment ions of m/z 71, 55, and 41 when using this mode. Acrylamide often requires derivatization to improve sensitivity on



a mass spectrometer. In this case, acrylamide is injected without derivatization onto a Thermo Scientific™ DSQ™ II mass spectrometer and an ultra low bleed Thermo Scientific™ TraceGOLD™ TG-WaxMS™ 30 m × 0.25 mm × 0.25 μm GC column.

Acrylamide is a highly polar water soluble compound having a logP value of -0.65 [1]. Such highly polar compounds are not readily amenable to GC, therefore a polar GC column is required. The TraceGOLD TG-WaxMS column is a polyethylene glycol-phase GC column that allows the analysis of polar compounds.

Experimental Details

Consumables		Cat. No.
Cartridge type:	HyperSep Hypercarb SPE cartridge, 500 mg/6 mL	10047-152
Column:	TraceGOLD TG-WaxMS, 30 m × 0.25 mm × 0.25 µm	10055-728
Septum:	Thermo Scientific BTO, 17 mm	10056-066
Liner:	Thermo Scientific™ Splitless FocusLiner™, 3 × 8 × 105 mm	101710-690
Column ferrules:	100% graphite ferrules for Thermo Scientific™ TRACE™ injector, 0.1–0.25 mm i.d.	10055-758
Column ferrules:	Graphite/Vespel® for transfer line 0.1–0.25 mm i.d.	10055-746
Vials and closures:	Assembled Advanced Vial Closure System vial and cap With ID Patch. Red PTFE/White Silicone	89239-442
Syringe filter:	Thermo Scientific™ Target2™ 30 mm GMF syringe filter membrane, 3.1 µm pore size	66064-840
Plastic syringe:	Thermo Scientific 3 mL plastic disposable syringes	66064-758

Sample Handling Equipment

HyperSep glass block manifold

Instrumentation

Thermo Scientific™ TRACE GC Ultra™ gas chromatograph

Thermo Scientific™ DSQ™ II single quadrupole mass spectrometer

Thermo Scientific™ TriPlus™ Autosampler

Chemicals and Reagents

HPLC grade water

HPLC grade methanol

Analytical grade formic acid

Sample Pretreatment

The potato chips were finely crushed with mortar and pestle and 1 g was weighed into a vial. A 1 g portion of the sample was spiked with 25, 50, 100, 250, 500, and 1000 ng/g of acrylamide standard in 2% formic acid / water. The sample was then filtered through a filter membrane.

Sample Preparation

Compounds:	Acrylamide and acrylamide-d ₃ (internal standard)
Matrix:	Potato chips
Conditioning stage:	Add 4 mL methanol, 4 mL water, and 4 mL 2% formic acid / water to the SPE cartridge.
Application stage:	Apply 1 mL of extract in 2% formic acid / water under vacuum at 1 mL/min to the SPE cartridge.
Washing stage:	Add 1 mL water to the SPE cartridge and dry for 20 min under vacuum.
Elution stage:	Apply 4 mL methanol to the SPE cartridge.
Additional stage:	Evaporate methanolic extract and reconstitute with 1 mL of 1 µg/mL of internal standard in methanol to the SPE cartridge.

Separation Conditions

Carrier gas:	Helium
Split flow:	50 mL/min
Column flow:	1.2 mL/min, constant flow
Oven temperature:	80 °C, 10 °C/min, 250 °C
Injector type:	Split/Splitless
Injector mode:	Splitless (1 min), constant septum purge
Injector temperature:	230 °C

Instrumentation

Transfer line temperature:	150 °C
Source temperature:	200 °C
Ionization conditions:	EI
Electron energy:	70 eV
SIM scan parameters:	m/z 71 for acrylamide and m/z 74 for acrylamide- d_3
Start time:	4.0 min
Dwell time:	0.1 s

Injection Conditions

Injection volume:	2 μ L
Pre- and post-needle injection dwell time:	0.5 s

Data Processing

Software:	Thermo Scientific™ Xcalibur™ software
-----------	---------------------------------------

Results

A standard addition calibration curve was constructed for acrylamide in matrix over the range 25–1000 ng/g. Standard addition calibration was chosen because acrylamide is endogenous in cooked foods and a suitable blank matrix was unavailable.

The amount of acrylamide present in the potato chips was calculated to be 450 ng/g. The chromatogram in Figure 1 shows the acrylamide peak in potato chips and acrylamide- d_3 internal standard spiked in potato chips.

The acrylamide concentration was calculated using the integrated response ratio of acrylamide/ acrylamide- d_3 (m/z 71/74). The acrylamide in the potato chips was calculated from the intercept of the x axis. An excellent linearity was demonstrated for this method with a coefficient of determination (R^2) of 0.999.

The accuracy of the back calculated concentrations for the amount of acrylamide spiked in potato chips was less than 10% (see Table 1).




Figure 1: TIC of chromatogram of 1 $\mu\text{g/mL}$ spiked acrylamide (m/z 71) and acrylamide- d_3 (m/z 74) extracted from potato chips

Specified Concentration ($\mu\text{g/mL}$)	Calculated Concentration	% Difference
0.25	0.225	-9.83
0.50	0.469	-6.22
1.00	0.983	-1.70
2.50	2.520	0.79
5.00	4.979	-0.41
10.0	10.037	0.37

Table 1: Accuracy data for the standard addition calibration curve for spiked acrylamide in potato chips



1.800.932.5000
vwr.com



Prices and product details are current when published; subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | VWR, the VWR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks, of VWR International, LLC and its related companies. All other marks referenced are registered by their respective owner(s). | Visit vwr.com to view our privacy policy, trademark owners and additional disclaimers. ©2015 VWR International, LLC. All rights reserved.

Thermo
SCIENTIFIC