EN

ISE F⁻ 9655
Refer to the meter manual for more information.

Caution
 ⚠ *Personal injury hazard. Glass components can break. Handle with care to prevent cuts.*
 ⚠ *Chemical exposure hazard. Refer to the current material safety data sheets (MSDS) for safety protocols.*

Introduction
This probe measures fluoride concentration in aqueous solutions. Use with a reference electrode.

Technical information

Range	0.02 mg/L to 2 g/L F ⁻
pH range	4 to 7
Linear range	0.1 mg/L to 2 g/L
Slope	59 mV/pF, 90 to 110% at 25 °C (77 °F)
Ionic strength adjuster (ISA)	Hach ISA 258999 pillow Hach TISAB solution 2829017
Operating temp	5 to 50 °C (41 to 122 °F)
Interferences	Fe ³⁺ , Al ³⁺ ; remove with TISAB

1 Connect the probe

2 Calibration
Add 1 ISA pillow or 5.0 mL of TISAB per 25 mL of standard.

3 Measurement
Add 1 ISA pillow or 5.0 mL of TISAB per 25 mL of sample.

Maintenance and storage
 Cleaning: rinse with distilled water between each standard and each sample.
 Storage: clean and dry the probe for storage.
 Damage or deposits on the membrane change the response of the electrode.

Troubleshooting
The ionic strength and temperature of samples and standards must be the same. Make sure the stir speed is the same for all standards and samples.

Warranty
6 months for manufacturing faults only. Damage from use is not covered.

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