

GE Healthcare

Amersham

Interleukin-4 [(h)IL-4] human,
Biotrak ELISA system

96 wells

Product booklet

Code: RPN2753



Page finder

1. Legal	3
2. Handling	4
2.1. Safety warnings and precautions	4
2.2. Storage	4
2.3. Expiry	4
3. Description	5
4. Introduction	6
5. Summary of the assay	7
6. Contents of the assay system	8
7. Assay methodology	9
7.1. Materials and equipment required	9
7.2. Sample preparation	9
7.3. Critical parameters	10
7.4. Assay procedure	11
7.5. Calculation of results	17
8. Additional Information	19
8.1. Specificity	19
8.2. Calibration	19
8.3. Reproducibility	19
8.4. Sensitivity	19
8.5. Recovery	20
8.6. Expected values	20
9. References	21
10. Related products	22

1. Legal

GE and GE monogram are trademarks of General Electric Company. Amersham and Biotrak are trademarks of GE Healthcare companies. GE Healthcare reserves the right, subject to any regulatory and contractual approval if required, to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation.

Contact your GE Representative for the most current information and a copy of the terms and conditions

© 2006 General Electric Company – All rights reserved.

<http://www.gehealthcare.com/lifesciences>

GE Healthcare UK Limited.

Amersham Place, Little Chalfont,
Buckinghamshire, HP7 9NA
UK

2. Handling

2.1. Safety warnings and precautions

Warning: *For research use only.*

Not recommended or intended for diagnosis of disease in humans or animals. Do not use internally or externally in humans or animals.

Warning: Contains methanol.

All chemicals should be considered as potentially hazardous. We therefore recommend that this product is handled only by those persons who have been trained in laboratory techniques and that it is used in accordance with the principles of good laboratory practice. Wear suitable protective clothing such as laboratory overalls, safety glasses and gloves. Care should be taken to avoid contact with skin or eyes. In the case of contact with skin or eyes wash immediately with water. See material safety

data sheet(s) and/or safety statement(s) for specific advice.

2.2. Storage

Store at 2-8°C.

2.3. Expiry

The expiry date is stated on the package and will normally be at least 4 weeks from the date of despatch.

3. Description

The Biotrak™ human interleukin-4 ELISA system from GE Healthcare provides a simple, specific, reliable and precise quantitative determination of (h)IL-4 in cell culture supernatants, serum, plasma and urine.

The assay system is based on a solid phase ELISA, which utilizes an antibody for (h)IL-4 bound to the wells of a microplate (12 x 8 well strip format) together with an antibody to (h)IL-4 and streptavidin conjugated to horseradish peroxidase. Although the Biotrak (h)IL-4 immunoassay contains recombinant (h)IL-4 and antibodies raised against recombinant (h)IL-4 it has been shown to quantitate accurately both natural (h)IL-4 and recombinant (h)IL-4.

(h)IL-4 can be measured in the approximate range 10–400 pg/ml (0.5–20 pg/well) in less than 4 hours using the protocol provided with the kit. Each pack contains sufficient material for 96 wells. If one standard curve is constructed, 42 unknowns can be measured in duplicate.

- High sensitivity - <2 pg/ml (0.05 pg/well)
- Same day protocol
- Pre-coated plate
- Specific for (h)IL-4

4. Introduction

Interleukin-4 (IL-4) is a multifunctional cytokine released by activated T cells that displays both stimulatory and inhibitory properties. IL-4 was first identified as a B cell growth factor capable of supporting the proliferation of B cells primed with antibodies to surface Ig. IL-4 has since been shown to play a pivotal role in regulating the immune response through its effects on B cells, T cells, monocytes and endothelial cells.

IL-4 production is limited to activated T cell subsets, T cell lines, a population of null spleen cells and, possibly, mast cell precursors. Human IL-4 is a 129 amino acid protein of 15–19 kDa (determined by SDS-PAGE) which displays variable N-linked glycosylation. IL-4 synthesis can be induced by agents that stimulate T cell activation, including antigen, anti-CD2, anti-CD3 or anti-TCR antibodies, phorbol esters, calcium ionophores and mitogenic lectins. The production of IL-4 is inhibited by TGF β and cyclosporin A.

In both humans and mice, the IL-4 receptor (IL-4R) consists of a single chain 140 kDa glycoprotein that binds IL-4 with high affinity ($2.6 \times 10^{-11}\text{M}$) in a species specific manner. IL-4R are expressed by both haematopoietic and non-haematopoietic cell types.

The treatment of mice with IL-4 antagonists has implicated IL-4 in the immune response to parasitic infections. IL-4 appears to promote the immune response to nematode infections. Finally, IL-4 produced in tumour cell lines transfected with IL-4 cDNA has been reported to inhibit solid tumour formation in mouse models. IL-4 has also been implicated in bone resorption.

The ability of IL-4 to inhibit solid tumour formation in mouse models suggest that it may have applications as an anti-neoplastic agent. IL-4 may also have therapeutic value in the management of inflammatory disease states by inhibiting the synthesis of inflammatory cytokines.

5. Summary of the assay

This assay employs the quantitative 'sandwich' enzyme immunoassay technique. An antibody specific for (h)IL-4 has been coated on the microplate provided in the kit. Samples are pipetted into the wells along with biotinylated antibody reagent. If present, the (h)IL-4 is bound by the immobilized antibody and the biotinylated antibody. After washing away any unbound sample proteins and biotinylated antibody, a streptavidin-HRP conjugate is added to the wells. Any (h)IL-4 which was bound by both the immobilized and the biotinylated antibody during the first incubation will be bound by the streptavidin conjugate. Following a wash to remove unbound conjugate, a substrate solution is added to the wells and colour develops in proportion to the amount of (h)IL-4 bound in the initial step.

In addition to the samples to be tested, a series of wells is prepared using known concentrations of the human IL-4 standard. A curve, plotting the optical density versus the concentration of the standard well, is prepared. By comparing the optical density of the samples to this standard curve, the concentration of the (h)IL-4 in the unknown samples is then determined.

6. Contents of the assay system

This pack contains the following assay components, sufficient material for 96 wells.

All reagents are stored refrigerated at 2–8°C. Refer to the expiry date on the kit box.

(h)IL-4 microplate - 96 well polystyrene microplate (12 strips of 8 wells) coated with an antibody against human IL-4.

Biotinylated antibody reagent - antibody against human IL-4 conjugated to biotin, with preservative, 8 ml.

(h)IL-4 standard - 2 vials of pre-diluted recombinant human IL-4, lyophilized.

Streptavidin-HRP concentrate - streptavidin conjugated to HRP with preservative, 50 µl.

Streptavidin-HRP dilution buffer - with preservative, 14 ml.

Standard diluent - with preservative, 12 ml.

Wash buffer concentrate - 30-fold concentrated solution, with preservative, 50 ml.

Pre-mixed TMB substrate reagent - with preservative and methanol, 12 ml. See safety data sheet supplied.

Stop solution - 0.18 M sulphuric acid, 15 ml.

Plate covers - 4 adhesive strips.

7. Assay methodology

Users are recommended to read this entire section before starting work.

7.1. Materials and equipment required

The following materials and equipment are required:

Pipettes or pipetting equipment with disposable tips (50 µl, 100 µl and 1.00 ml)*

Disposable polypropylene test tubes - do not use polystyrene, polycarbonate or glass

Measuring cylinder 2 L

Distilled or deionized water

Plate reader capable of reading at 450 nm

A centrifuge for preparing streptavidin-HRP solution

Optional equipment

Assays may be performed with commercially available microplate washers to aid convenience and assay throughput.

7.2. Sample preparation

Serum, plasma, cell culture supernatants and urine

Serum, plasma, urine and cell culture supernatant samples that are to be assayed within 24 hours should be stored at 2–8°C. Specimens to be stored for longer periods of time should be frozen at -70°C to avoid loss of biologically active cytokine. Avoid freezing and thawing samples more than once.

* GE Healthcare supplies a range of pipettes and disposable tips (see related products).

Dilution of test samples

If you suspect that the (h)IL-4 concentration of a sample exceeds the highest point of the standard curve, prepare one or more 5-fold dilutions of the test sample. Mix thoroughly between dilutions and before assaying.

7.3. Critical parameters

- 1) First select the number of strips to be run and allow samples and all reagents to reach room temperature prior to performing the assays. Do not use water baths to thaw samples or reagents.
- 2) Mix samples and all reagents thoroughly before use.
- 3) Avoid excessive foaming of reagents. Also avoid exposure of reagents to excessive heat or light during storage and incubation.
- 4) Avoid handling the tops of the wells both before and after filling.
- 5) Standards and samples should be assayed in duplicate.
- 6) Run a separate standard curve for each assay.
- 7) The total dispensing time for each plate should not exceed 20 minutes.
- 8) Use only coated wells from the same reagent batch for each assay. Also do not mix reagents from different kit lots.
- 9) It is important that the wells are washed thoroughly and uniformly. If using an automatic washer check operation of heads before starting. If washing by hand ensure that all wells are completely filled at each wash.
- 10) A small amount of precipitate may be present in some vials. It will not affect assay performance and should be ignored.

7.4. Assay procedure

Reagent preparation

Wash buffer concentrate

Any precipitate formed during storage will redissolve upon dilution. Dilute 30-fold to prepare 1500 ml of wash buffer. Store at 2–8°C until the expiry date of the kit. Do not use wash buffer if it becomes visibly contaminated on storage.

Streptavidin-HRP solution

Prepare the exact amount of streptavidin-HRP solution no more than 15 minutes prior to use.

The streptavidin-HRP concentrate may require spinning down to force the contents to the bottom of the vial. Add 30 µl of streptavidin-HRP concentrate per 12 ml of streptavidin-HRP dilution buffer in a plastic 15 ml tube and mix gently.

If running partial plates use 2.5 µl of streptavidin-HRP concentrate and 1 ml of streptavidin-HRP dilution buffer per strip being run.

Preparation of working standards

- a) Two vials of lyophilized standards are provided with this kit. Reconstitute and use one vial per partial plate.
- b) Prepare standards shortly before use. Use within one hour of reconstitution.

Do not store reconstituted standards.

- c) When **running culture supernatant samples**, reconstitute standard in distilled or deionized water. Reconstitution volume is stated on the standard vial label. The standard will take approximately 1 minute to dissolve. Mix by gently inverting the vial. Use your culture medium to prepare the dilutions of the standard curve, go to step e) opposite for further instructions. If running a partial plate, refer to step a) above.

- d) If running serum, plasma or urine samples**, reconstitute standard with distilled or deionized water. Reconstitution volume is stated on the standard vial label. The standard will take approximately 1 minute to dissolve. Mix by gently inverting the vial. Use the standard diluent provided to prepare the dilutions of the standard curve. If running a partial plate, refer to step a).
- e)** Label 6 tubes, one for each standard curve point: 400 pg/ml, 160 pg/ml, 64 pg/ml, 25.6 pg/ml, 10.24 pg/ml and 0 pg/ml. Then prepare 1:2.5 serial dilutions for the standard curve as follows:
- f)** Pipette 240 µl of appropriate diluent (see steps c) and d) above) into each tube.
- g)** Pipette 160 µl of the reconstituted standard into the first tube, 400 pg/ml and mix.
- h)** Pipette 160 µl of this dilution into the second tube labelled 160 pg/ml and mix.
- i)** Repeat serial dilutions three more times. These concentrations, 400 pg/ml, 160 pg/ml, 64 pg/ml, 25.6 pg/ml, 10.24 pg/ml and 0 pg/ml are your standard curve.

Running partial plates

This ELISA provides the flexibility to run two partial plates on separate occasions. Decide the number of strips you wish to run, leaving the strips to be used in the frame. Remove the unnecessary strips and store them in the foil pouch with the desiccant provided at 2–8°C, making sure the foil pouch is sealed tightly.

When adding the TMB substrate reagent, pour out from the bottle **ONLY** the amount needed to run the first half plate. Do not combine left over substrate with that reserved for the second half of the plate. Care must be taken to ensure that the remaining TMB substrate reagent is not contaminated. If the substrate reagent is bright blue prior to use, it has been contaminated. **DO NOT USE.**

Assay protocol

- 1) Prepare assay reagents and working standards as described in the previous sections.
- 2) Set up the microplate with sufficient wells to enable the running of all standards and samples as required (see figure 1).
- 3) Remove excess microplate strips from the frame and store in the resealable foil bag.
- 4) Add 50 μ l of biotinylated antibody reagent to each well that is to be used.
- 5) Add 50 μ l of standard or sample per well. Cover with adhesive strip provided and incubate for 2 hours at room temperature (20–25°C).
- 6) Aspirate or decant each well and wash, repeating the process twice for a total of three washes. Wash vigorously by filling each well with wash buffer using a washbottle, pipette or manifold dispenser. Complete removal of liquid at each step is essential for good performance. After the last wash, remove any remaining wash buffer by inverting the plate and blotting it against clean paper towelling.
- 7) Add 100 μ l of pre-diluted streptavidin-HRP conjugate. Cover with a new adhesive strip and incubate for 30 minutes at room temperature (20–25°C).
- 8) Repeat the aspiration/wash step as in step 6.
- 9) Add 100 μ l of TMB substrate solution into each well, incubate for 30 minutes at room temperature (20–25°C). If the substrate reagent is bright blue prior to use, do not use. **THE PLATE SHOULD BE DEVELOPED IN THE DARK.** Do not cover the plate with aluminium foil or an adhesive strip.

10) Add 100 μ l of stop solution to each well.

11) Determine the optical density of each well within 30 minutes, using a spectrophotometer set to 450 nm.

	1	2	3	4	5	6	7	8	9	10	11	12
A	0	0	S	S	S	S	S	S	S	S	S	S
B	10.24	10.24	S	S	S	S	S	S	S	S	S	S
C	25.6	25.6	S	S	S	S	S	S	S	S	S	S
D	64	64	S	S	S	S	S	S	S	S	S	S
E	160	160	S	S	S	S	S	S	S	S	S	S
F	400	400	S	S	S	S	S	S	S	S	S	S
G	S	S	S	S	S	S	S	S	S	S	S	S
H	S	S	S	S	S	S	S	S	S	S	S	S

Figure 1. Recommended positioning of standard (0–400 pg/ml) and sample wells (S).

Summary of assay protocol




Table 1. Assay protocol (all volumes are in microlitres)

	Zero standard (B_0)	Standards	Samples
Biotinylated antibody reagent	50	50	50
Standard	-	50	-
Standard diluent or cell culture media	50	-	-
Sample	-	-	50
Cover plate, incubate at room temperature for 2 hours.			
Aspirate/decant and vigorously wash all wells three times with wash buffer.			
Streptavidin-HRP conjugate	100	100	100
Cover plate, incubate at room temperature for 30 minutes			
Aspirate/decant and wash vigorously all wells			
three times with wash buffer			
Substrate	100	100	100
Incubate at room temperature for 30 minutes in the dark			
Stop solution	100	100	100
Determine optical density at 450 nm within 30 minutes.			

7.5. Calculation of results

Average the duplicate readings for each standard, control and sample and subtract the zero standard optical density.

Plot the optical density for the standards versus the concentration of the standards and draw the best curve. The data can be linearized by using a log/log plot and regression analysis may be applied to the log transformation.

Figure 2 shows such a plot of the data from table 2. The standard curve is provided for illustration only. A standard curve should be generated for each set of samples to be assayed. This allows for the measurements of 42 unknowns in duplicate.




Figure 2. Standard curve

Typical assay data

The following data (table 2) were obtained for a standard curve using the protocol provided.

Table 2. Typical assay data

Tube	Optical density subtracted	Zero standard
Zero standard	0.026	-
10.24 pg/ml standard	0.099	0.073
25.6 pg/ml standard	0.211	0.185
64 pg/ml standard	0.464	0.438
160 pg/ml standard	1.054	1.028
400 pg/ml standard	2.233	2.207

8. Additional Information

8.1. Specificity

This assay recognizes both natural and recombinant (h)IL-4. It does not cross react with (m)IL-4, (h)IL-2, (h)IL-6, (h)IL-10, (h)GM-CSF, (h)IFN gamma (h)TNF β or (h)TNF α .

8.2. Calibration

The standards in this ELISA are calibrated to the NIBSC reference lot 88/656.

One (1) pg of Biotrak standard = 6.5 NIBSC pg.

8.3. Reproducibility

Within-assay precision

The within-assay coefficient of variation of the ELISA has been determined to be <10%.

Between-assay precision

The between-assay coefficient of variation of the ELISA has been determined to be <10%.

8.4. Sensitivity

The minimum detectable dose of (h)IL-4 was determined to be <2 pg/ml (0.10 pg/well), by adding two standard deviations to the optical density value of zero and calculating the corresponding concentration from the standard curve.

8.5. Recovery

Recovery in the ELISA has been determined by spiking recombinant cytokine into cell culture media, normal human serum, plasma and urine samples. The values below are typical recoveries:

Spiked level (pg/ml)	30	100	200
Mean serum recovery	116	104	114
Mean plasma recovery	107	105	111
Mean urine recovery	100	103	98

8.6. Expected Values

Human IL-4 was not detected in 8 normal human serum samples, 8 normal human plasma samples and 5 normal human urine samples.

9. References

1. KUTZHALS, J.A. *et al.*, *J. Immunol. Methods*, **156**, p.239, 1992.
2. MAHANTY, S. *et al.*, *J. Immunol.*, **148**, p.3567, 1992.
3. VAN der POUW-KRAAN, T. *et al.*, *Eur. J. Immunol.*, **22**, p.1237, 1992.
4. BANCHEREAU, J. *et al.*, *Bull. Cancer*, **78**, p.299, 1991.
5. CAVALLO, M.G. *et al.*, *Clin. Exp. Immunol.*, **86**, p.256, 1991.
6. FUNAUCHI, M. *et al.*, *Ryumachi*, **31**, p.493, 1991.
7. BERGSTEDT-LINDQVIST, S. *et al.*, *Eur. J. Immunol.*, **18**, p.1073, 1988.
8. van KIMMENADE, A. *et al.*, *Eur. J. Biochem.*, **173**, p.109, 1988.
9. HU-LI, J. *et al.*, *J. Exp. Med.*, **165**, p.157, 1987.
10. YOKOTA, T. *et al.*, *Proc. Nat'l. Acad. Sci.*, **83**, p.5894, 1986.

10. Related products

Biotrak range of human cytokine ELISA systems

Interleukin-1 α [(h)IL-1 α]	RPN 2750
Interleukin-1 β [(h)IL-1 β]	RPN 2751
Interleukin-2 [(h)IL-2]	RPN 2752
Interleukin-6 [(h)IL-6]	RPN 2754
Interleukin-10 [(h)IL-10]	RPN 2755
Granulocyte-macrophage colony stimulating factor [(h)GM-CSF]	RPN 2756
Interferon-gamma [(h)IFN γ]	RPN 2757
Tumour necrosis factor, alpha [(h)TNF α]	RPN 2758
Interferon-alpha [(h)IFN α]	RPN 2759
Interleukin-8 [(h)IL-8]	RPN 2760

Biotrak range of mouse cytokine ELISA systems

Interleukin-1 α [(m)IL-1 α]	RPN 2719
Interleukin-1 β [(m)IL-1 β]	RPN 2720
Interleukin-2 [(m)IL-2]	RPN 2710
Interleukin-3 [(m)IL-3]	RPN 2711
Interleukin-4 [(m)IL-4]	RPN 2712
Interleukin-5 [(m)IL-5]	RPN 2713
Interleukin-6 [(m)IL-6]	RPN 2714
Interleukin-10 [(m)IL-10]	RPN 2715
Granulocyte-macrophage colony stimulating factor [(m)GM-CSF]	RPN 2716
Interferon-gamma [(m)IFN γ]	RPN 2717

Soluble intercellular adhesion molecule-1 [(m)ICAM-1]	RPN 2721
Tumour necrosis factor- α [(m)TNF- α]	RPN 2718

Range of unlabelled and radiolabelled growth factors and cytokines

Cell proliferation assay system and reagents

Cell proliferation assay system	RPN 210
Cell proliferation kit (for immunocytochemical/ immunohistochemical measurement)	RPN 20
Monoclonal anti-bromodeoxyuridine	RPN 202
Cell proliferation labelling reagent	RPN 201

Pipettes and pipette tips

Single channel, variable volume pipettes

Volume range

0.5–10 μ l	RPN 2340
5–50 μ l	RPN 2341
50–200 μ l	RPN 2342
200–1000 μ l	RPN 2343
1–5 ml	RPN 2344

Multi-channel, variable volume pipettes

8 channel, 5–50 μ l	RPN 2372
8 channel, 50–250 μ l	RPN 2373

GE Healthcare's range also includes fixed volume, 4 and 12 multi-channel variable volume pipettes, a range of pipette tips and related accessories.

GE Healthcare offices:

GE Healthcare Bio-Sciences AB
Björkgatan 30 751 84

Uppsala
Sweden

GE Healthcare Europe GmbH
Munzinger Strasse 5 D-79111
Freiburg
Germany

GE Healthcare UK Limited
Amersham Place
Little Chalfont
Buckinghamshire
HP7 9NA
UK

GE Healthcare Bio-Sciences
Corp
800 Centennial Avenue
P.O. Box 1327
Piscataway
NJ 08855-1327
USA

GE Healthcare Bio-Sciences KK
Sanken Bldg. 3-25-1
Hyakunincho Shinjuku-ku
Tokyo 169-0073
Japan

**GE Healthcare
regional office
contact numbers:**

Asia Pacific
Tel: +85 65 62751830
Fax: +85 65 62751829

Australasia
Tel: +61 2 8820 8299
Fax: +61 2 8820 8200

Austria
Tel: 01/57606-1613
Fax: 01/57606-1614

Belgium
Tel: 0800 73 890
Fax: 02 416 8206

Canada
Tel: 1 800 463 5800
Fax: 1 800 567 1008

**Central, East, & South
East Europe**
Tel: +43 1 972 720
Fax: +43 1 972 722 750

Denmark
Tel: 45 70 25 24 50
Fax: 45 45 16 2424

Eire
Tel: 1 800 709992
Fax: +44 1494 542010

Finland & Baltics
Tel: +358 9 512 3940
Fax: +358 9 512 39439

France
Tel: 01 69 35 67 00
Fax: 01 69 41 98 77

Germany
Tel: 0800 9080 711
Fax: 0800 9080 712

Greater China
Tel: +852 2100 6300
Fax: +852 2100 6338

Italy
Tel: 02 26001 320
Fax: 02 26001 399

Japan
Tel: +81 3 5331 9336
Fax: +81 3 5331 9370

Korea
Tel: 82 2 6201 3700
Fax: 82 2 6201 3803

Latin America
Tel: +55 11 3933 7300
Fax: +55 11 3933 7304

Middle East & Africa
Tel: +30 210 96 00 687
Fax: +30 210 96 00 693

Netherlands
Tel: 0800-82 82 82 1
Fax: 0800-82 82 82 4

Norway
Tel: +47 815 65 777
Fax: +47 815 65 666

Portugal
Tel: 21 417 7035
Fax: 21 417 3184

Russia, C.I.S. & N.I.S
Tel: +7 495 956 5177
Fax: +7 495 956 5176

Spain
Tel: 902 11 72 65
Fax: 935 94 49 65

Sweden
Tel: 018 612 1900
Fax: 018 612 1910

Switzerland
Tel: 0848 8028 10
Fax: 0848 8028 11

UK
Tel: 0800 515 313
Fax: 0800 616 927

USA
Tel: +1 800 526 3593
Fax: +1 877 295 8102

<http://www.gehealthcare.com/lifesciences>

GE Healthcare UK Limited

Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA
UK



imagination at work