

ExcelGel™ SDS Buffer Strips

ExcelGel SDS Buffer Strips for horizontal SDS electrophoresis

ExcelGel SDS precast anode and cathode buffer strips are to be used with:

- ExcelGel SDS Gradient 8-18 (Code No. 80-1255-53)
- ExcelGel XL SDS 12-14 (17-1236-01), and
- ExcelGel SDS Homogeneous 7.5, 12.5 and 15 (80- 1260-01, 80-1261-01 and 80-1262-01).

For complete information see Instructions provided with each package of ExcelGel SDS gel media.

During the run the precast SDS buffer strips supply the gel with buffer ions. The buffer system in the strips forms, together with the gel buffer, a dis-continuous buffer system. The buffer system and the polyacrylamide gradient gel jointly ensure good separation and sharp bands.

Method

Always wear clean gloves when working with polyacrylamide gels and buffer strips, particularly when using sensitive staining methods. Vinyl gloves are recommended because they have little tendency to stick to the buffer strip, contrary to gloves made of latex.

The anodic buffer strips are marked with the symbols +++ and the cathodic buffer strips are marked with - - -. The identification of the anodic buffer strips is facilitated by the addition of the dye Orange G

If you do not use vinyl gloves, you can prevent the gloves from sticking to the buffer strips by moistening the fingers with distilled water.

- 1 Open the SDS buffer strip by gently pulling away the aluminium foil. Apply pulling force perpendicular to the strip.
- 2 Loosen the strips at one of their ends by carefully putting a spatula under the gel.
- 3 Remove the strips by gripping them with both hands at least 2 cm away from the ends. (Should the strip break, you can still use it by carefully putting the ends together on the gel surface).
- 4 Apply the cathodic and the anodic SDS buffer strip on the respective side of the gel – **narrow side down**.




Fig. 1 Pull up carefully 2–3 cm at a time, keeping the peeling edge vertical to the plastic foil.

Technical data

ExcelGel SDS buffer strips	17-1342-01
Contents:	6 anode SDS buffer strips 6 cathode SDS buffer strips
Strip matrix:	polyacrylamide (T=12%, C=3%)
Dimension:	length 245 mm, height 4.5 mm
Buffer system Anode SDS buffer strip:	0.45 mol/L Tris/Acetate, pH 6.6, 4 g/L SDS and 0.05 g/L Orange G
Buffer system Cathode SDS buffer strip:	0.08 mol/L Tris, 0.80 mol/L Tricine and 6 g/L SDS, pH 7.1
Storage:	4°C to 8°C
Shelf life:	18 months



For local office contact information, visit
www.gelifesciences.com/contact

GE Healthcare Bio-Sciences AB
Björkgatan 30
751 84 Uppsala
Sweden

www.gelifesciences.com/protein-purification

GE, imagination at work and GE monogram are trademarks of General Electric Company.

ExcelGel is a trademark of GE Healthcare companies.

© 1993-2010 General Electric Company—All rights reserved.

Previously published Aug. 1993.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.

GE Healthcare UK Ltd
Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA, UK

GE Healthcare Bio-Sciences Corp
800 Centennial Avenue, P.O. Box 1327, Piscataway, NJ 08855-1327, USA

GE Healthcare Europe GmbH
Munzinger Strasse 5, D-79111 Freiburg, Germany

GE Healthcare Japan Corporation
Sanken Bldg. 3-25-1, Hyakunincho, Shinjuku-ku, Tokyo 169-0073, Japan



imagination at work