

# NxSeq UltraLow DNA Library Kit\* - 12 reactions

## Build high quality, Illumina-compatible DNA fragment libraries from 50 pg to 75 ng of DNA in only 3 hours

The NxSeq™ UltraLow DNA Library Kit and NxSeq Single Indexing Kits allow you to build complex, high quality DNA fragment libraries from extremely low DNA input amounts – as low as 50 pg (Table 1). If you have more DNA, no problem; you can use as much as 75 ng of input DNA with this system.

In order to build high quality DNA fragment libraries from limiting DNA inputs, the entire library prep workflow must be optimised starting with adaptor ligation efficiency (Figure 1) followed by un-biased, robust PCR amplification of the library to produce enough sequenceable material. This system is designed and built to produce the highest quality libraries possible (Table 1).




Figure 1. Efficiency of adaptor ligation to DNA library fragments. Two independent sets of libraries were prepped per kit/organism using manufacturer's recommended protocols and 1 ng of the same sheared genomic DNA input per library. Briefly, DNA fragments with adaptors ligated to both ends were measured using triplicate qPCR assays and an universal qPCR primer set, designed by LGC, Biosearch Technologies, that binds to and amplifies all adaptor-ligated DNA fragments independent of the kit used. Efficiency was determined by comparing the qPCR quantitation to fluorescence DNA quantitation. Efficiency data was averaged and then normalised to the corresponding NxSeq UltraLow Library Kit data (1.0) and plotted.

\* and NxSeq single indexing kits



High quality libraries and sequencing data start with high efficiency adaptor ligation

- **High quality data:** High efficiency adaptor ligation reactions produce complex libraries that yield improved sequencing depth, uniformity, and fewer “zero coverage” regions than other kits
- **Sensitive:** Construct DNA fragment libraries from as little as 50 pg to as much as 75 ng of sheared/fragmented DNA
- **Minimal bias:** Robust, uniform PCR amplification improves coverage and sequencing depth uniformity
- **Flexible:** Extensively tested in whole genome sequencing and re-sequencing, and compatible with other applications, such as exome-seq, ChIP-seq and sequencing of difficult FFPE DNA samples
- **Fast and easy-to-use:** Simplified, 3 hour protocol gets your samples on the sequencer quickly while decreasing the risk of handling errors
- **High value:** Cost-effective, high performance library and indexing kits

# NxSeq UltraLow DNA Library Kit, 12 reactions and NxSeq Single Indexing Kits specification



<b>Special instructions/requirements</b>	The NxSeq UltraLow DNA Library Kit and NxSeq Single Indexing Kits are only compatible with Illumina sequencers. <i>Please note that at least one NxSeq Single Indexing Kit is required to complete library prep and must be purchased along with a library kit.</i>
<b>Components</b>	The NxSeq UltraLow DNA Library Kit contains Enzyme Mix (EM), 2X Buffer (2XB), Ligase (LIG), 2X PCR Master Mix (MM) and Elution Buffer (EB). Each NxSeq Single Indexing Kit contains a Universal Adaptor and (12) different Primer Indexing Mixes with enough Universal Adaptor for 48 libraries and enough of each Primer Indexing Mix for 4 library amplification reactions (48 total reaction for all primer sets). Set A contains Primer Indexing Mixes 1-12 and Set B contains Primer Indexing Mixes 13-16, 18-23, 25, and 27. Each NxSeq Single Index equals the TruSeq® LT Index with the same number.

## Performance information

Sheared gDNA input	50 pg	250 pg	500 pg	75 ng
PCR cycles used	16	15	13	4
Reads sampled per library	312,500	312,500	312,500	312,500
No. mapped reads	303,442	304,135	304,894	304,935
Percent mapped reads	97.10%	97.32%	97.57%	97.58%
CLC mapped duplicates	0.84%	0.15%	0.08%	0.02%
Ave. coverage depth	9.62	9.72	9.74	9.84
Coverage depth standard dev.	3.75	3.71	3.71	3.6
No. of zero locations	109	91	91	56
Ave. zero region length (bp)	24	24	25	27
Total zero length (bp)	2651	2214	2175	1491

Table 1. Generation and sequencing of DNA fragment libraries from 50 pg to 75 ng of E. coli gDNA. Mechanically sheared (300 bp peak) E. coli K12 genomic DNA was diluted and used to make duplicate libraries starting with the indicated input amounts and number of PCR amplification cycles using the NxSeq UltraLow DNA Library Kit and a NxSeq Single Indexing Kit. The final libraries were quantitated, diluted to 2 nM, pooled and sequenced on a MiSeq using 2x150 bp chemistry. The average data from the duplicate libraries are presented.

High quality whole genome sequencing data from ≥50 pg of input DNA



	NxSeq UltraLow Kit	Competitor A	Competitor B
Total PF sampled reads	200,000,000	200,000,000	200,000,000
R1 mapped reads	96.52%	96.02%	96.10%
R2 mapped reads	95.51%	95.36%	94.97%
Fragment length (bp)	289 +/- 110	191 +/- 61	225 +/- 109
Duplicates (proper read pairs)	2.15%	3.73%	3.54%
Coverage depth	8.32	6.15	6.96
Autosomal coverage	98.64%	98.40%	98.73%


Table 2. HiSeq X Ten sequencing results with fragment libraries made with 10 ng of human genomic DNA. DNA fragment libraries were constructed using the indicated kits according to each manufacturer's recommended protocol. Briefly, human genomic DNA (NA15510, Coriell) was sheared (median size 300 bp), and 10 ng aliquots were used to make duplicate libraries with each kit. The libraries were PCR amplified for the following number of cycles: NxSeq UltraLow Kit, 8 cycles; Competitor A, 8 cycles; and Competitor B, 10 cycles. The amplified libraries were cleaned/size selected as recommended. The final libraries were sent to Hudson Alpha for sequencing, and each set of duplicates was pooled and sequenced on its own HiSeq X Ten lane using 2x150 bp chemistry. During analysis, the optical duplicates generated by the HiSeq X Ten were removed using the Clumpify program, and then the indicated number of reads were analysed. The averaged data from the replicate libraries are shown.

Improved sequencing performance on the HiSeq X Ten Instrument



## Ordering information

VWR cat no.	Size	Description
76180-800	12 rxn	NxSeq UltraLow DNA Library Kit, 12 reactions
76180-802	48 rxn (12x4 rxn)	NxSeq Single Indexing Kit, Set A
76180-804	48 rxn (12x4 rxn)	NxSeq Single Indexing Kit, Set B



**BIOSEARCH TECHNOLOGIES**  
GENOMIC ANALYSIS BY LGC

**Integrated tools.  
Accelerated science.**



For Research Use Only. Not for use in diagnostic procedures.

All trademarks and registered trademarks mentioned herein are the property of their respective owners. All other trademarks and registered trademarks are the property of LGC and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details. No part of this publication may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopying, recording or any retrieval system, without the written permission of the copyright holder. © LGC Limited, 2019. All rights reserved. GEN/618/SW/1219



800 932 5000 | VWR.COM

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2019 Avantor, Inc. All rights reserved.