

GE Healthcare

Amersham Rediprime II DNA Labeling System

Product Booklet

Codes: RPN1633
 RPN1634



Page finder

1. Legal	3
2. Handling	4
2.1. Safety warnings and precautions	4
2.2. Storage	4
2.3. Expiry	4
3. Components of the system	5
4. Description	6
5. Critical parameters	7
6. Additional equipment and reagents required	8
7. Rediprime II DNA Labeling System Protocol	9
8. Additional information	12
8.1. Quality control	12
8.2. Using the control DNA	12
8.3. Use of alternative reaction conditions	12
8.4. Monitoring the reaction and calculating the specific activity of labeled DNA	13
8.5. Removal of unincorporated nucleotides	14
9. Related products	16
10. References	17

1. Legal

GE, imagination at work and GE monogram are trademarks of General Electric Company.

Amersham, Rediprime, Rapid-hyb, Hybond, Hyperfilm, Hypercassette, Hyperscreen, Ready-To-Go, Sephadex and Megaprime are trademarks of GE Healthcare companies.

© 1995–2008 General Electric Company - All rights reserved.
Previously published 1995.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them.
A copy of these terms and conditions is available on request.

Contact your local GE Healthcare representative for the most current information.

<http://www.gehealthcare.com/lifesciences>

GE Healthcare UK Limited.

Amersham Place, Little Chalfont,
Buckinghamshire, HP7 9NA, UK

2. Handling

2.1. Safety warnings and precautions

Warning: For research use only.

Not recommended or intended for diagnosis of disease in humans or animals. Do not use internally or externally in humans or animals.

Caution: For use with radioactive material.

This product is to be used with radioactive material. Please follow the manufacturer's instructions relating to the handling, use, storage and disposal of such material.

All chemicals should be considered as potentially hazardous. We therefore recommend that this product is handled only by those persons who have been trained in laboratory techniques and that it is used in accordance with the principles of good laboratory practice. Wear suitable protective clothing, such as laboratory overalls,

Safety glasses and gloves. Care should be taken to avoid contact with skin or eyes. In the case of contact with skin or eyes, wash immediately with water. See material safety data sheet(s) and/or safety statement(s) for specific advice.

2.2. Storage

Store at 2–25°C.

2.3. Expiry

The kit components are stable for up to 6 months when stored under the recommended conditions.

3. Components of the system

Rediprime II DNA Labeling System

Component	RPN1633	RPN1634
Labeling reaction Buffered solution of dATP, dGTP, dTTP, exonuclease free Klenow enzyme and random primers in a dried, stabilized form	30 reactions	60 reactions
Control DNA 300 ng of lambda <i>Hind</i> III DNA in a dried, stabilized form	1 tube	2 tubes
Protocol summary card		

4. Description

Feinberg and Vogelstein (1, 2) introduced the use of random sequence hexanucleotides to prime DNA synthesis on denatured template DNA at numerous sites along its length. The primer-template complex is a substrate for the Klenow fragment of DNA polymerase I. By replacing a non-radioactive nucleotide with the radiolabeled equivalent in the reaction mixture, newly synthesized DNA is made radioactive.

Very small amounts of input DNA can be labeled, enabling very high specific activity probes to be produced with relatively small quantities of added nucleotides. These radioactive labeled fragments can then be used as sensitive hybridization probes for a wide range of filter based applications (3-6).

Traditional protocols for the random primer labeling of DNA required reaction times of at least 30 minutes. More recent procedures, such as that used in GE Healthcare popular Megaprime™ DNA Labeling system, allow the labeling of template DNA to the same specific activity but at a greatly accelerated rate. This rapid labeling at 37°C is achieved by the use of primers at high concentration giving more efficient priming from the template.

Rediprime™ II DNA Labeling System from GE Healthcare have been developed for extra convenience and performance. The systems provide individually dispensed reaction mixes which are dried in the presence of a stabilizer and a dye to make labeling probes easier for the user. There is no requirement to store the systems in a freezer, they can be stored in the fridge or on the laboratory bench ready for use. Rediprime reaction mixes have been formulated using an improved exonuclease-free Klenow to give probes with specific activity of 1.9×10^9 dpm/μg or greater after 10 minutes incubation at 37°C with the majority of DNA substrates.

5. Critical parameters

To get high labeling efficiencies, there are two critical factors:

- TE buffer must be used as a diluent for the DNA template.
- The DNA template must be boiled in a total volume of 45 μ l.
- Solutions which are too dilute to be used directly should be concentrated by ethanol precipitation. Redissolution in a volume ≤ 45 μ l of 10 mM Tris/HCl pH 8.0, 1 mM EDTA is recommended.

6. Additional equipment and reagents required

- Pipettes or pipetting equipment for 5, 50 and 100 μ l
- TE buffer
- Boiling water bath
- 37°C water bath or heating block
- 0.2 M EDTA
- Refrigerator
- Freezer

7. Rediprime II DNA Labeling System protocol

It is recommended that the protocol is read thoroughly before using the system.

Rediprime allows DNA from a variety of sources to be labeled to high specific activity using [α - 32 P]dCTP. The system has been designed for use with [α - 32 P]dCTP with a specific activity of 3000 Ci/mmol.

Each reaction tube can label up to 25 ng of DNA and, after incubation for 10 minutes at 37°C, probes with a specific activity of 1.9×10^9 dpm/ μ g or greater can be produced.

DNA prepared by standard minilysate methods may be used. DNA in restriction enzyme buffers may be added directly to the reaction, and the reaction can also be performed with DNA in low melting point agarose gel slices (see page 12).

Protocol	Notes
1. Dilute the DNA to be labeled to a concentration of 2.5–25 ng in 45 μ l of 10 mM Tris HCl pH8.0, 1 mM EDTA. (TE buffer)	1. To get high labeling efficiencies, it is important that TE buffer is used and that the final volume used to reconstitute the mix equals 45 μ l.
2. Denature the DNA sample by heating to 95–100°C for 5 minutes in a boiling water bath.	
3. Snap cool the DNA by placing on ice for 5 minutes after denaturation.	

Protocol	Notes
4. Centrifuge briefly to bring the contents to the bottom of the tube.	
5. Add the denatured DNA to the reaction tube.	5. Do not mix at this stage.
6. Add 5 μ l of [α - 32 P]dCTP and mix by pipetting up and down about 12 times, moving the pipette tip around in the solution.	
7. Incubate at 37°C for 10 minutes.	7. If desired, the labeling reaction can be allowed to proceed at room temperature, or may also be left to proceed overnight. Probes of high specific activity can be generated in 20–60 minutes at this temperature (see figure 1).
8. Stop the reaction by adding 5 μ l of 0.2 M EDTA. For use in hybridization, denature the labeled DNA by heating to 95–100°C for 5 minutes, then snap cool on ice for 5 minutes.	
9. Centrifuge the tube briefly and mix the contents of the tube well.	

Protocol	Notes
10. We recommend that if 25 ng of template was used in the labeling reaction, then 14 μl of this labeled probe is used per 5 ml of hybridization buffer.	10. Recommended probe concentration based on template quantity is 1.4 ng/ml. This is likely to give an actual probe concentration of 2.8 ng/ml due to the amount of probe synthesized during the reaction.

8. Additional Information

8.1. Quality control

Rediprime II DNA Labeling System is tested by our quality control group to ensure an incorporation rate greater than 55% after 10 minutes at 37°C using 50 μCi [α - ^{32}P]dCTP.

8.2. Using the control DNA

The performance of Rediprime systems can be checked by using the control DNA supplied. The vial contains 300 ng of lambda DNA which should be dissolved in TE buffer before use. Once reconstituted 5 μl of control DNA contains 25 ng. It may be stored at 2–8°C for up to 1 month, or for longer periods at -15°C to -30°C.

Protocol

1. Reconstitute the DNA by adding 60 μl of TE buffer and flick the tube until the DNA has dissolved.
2. Spin briefly and proceed with steps 1–8 of the Rediprime protocol.

8.3. Use of alternative reaction conditions

Labeling at room temperature

If required, labeling reactions can be carried out at room temperature. Typical results for labeling at this temperature are shown in figure 2. Room temperature should be used if reactions are to be left to incubate overnight.

Labeling of DNA fragments in low melting point agarose (7)

1. Cut out the DNA band from an ethidium bromide stained agarose gel and trim off excess agarose.
2. Weigh the agarose plug and add 3 ml of water for each gram of gel.

3. Place in a boiling water bath for 5–10 minutes to melt the gel and to denature the DNA. Store the DNA solution at 37°C until required.
4. Remove an aliquot and add 10 mM Tris HCl pH 8.0, 1 mM EDTA. (TE buffer) to give a final volume of 45 µl containing 2.5–25 ng DNA.
5. Proceed from step 4 of the Rediprime protocol, extending the labeling reaction time to 15–30 minutes at 37°C.

Labeling using less [α - ^{32}P]dCTP

For maximum probe specific activity (1.9×10^9 dpm/µg or greater) and sensitivity in hybridizations the recommended 50 µCi of [α - ^{32}P]dCTP should be used in all reactions. In certain cases, where a lower specific activity probe may be acceptable, the amount of [α - ^{32}P]dCTP may be reduced. For example, using 20 µCi [α - ^{32}P]dCTP will typically yield probes of specific activity up to 1×10^9 dpm/µg.

8.4. Monitoring the reaction and calculating the specific activity of labeled DNA

- It is possible to monitor the reaction using DE81 paper (8).
- It is also possible to monitor the reaction using precipitation with TCA (trichloroacetic acid) (9).
- The specific activity of labeled DNA can be calculated using the following formulae. First calculate the amount of DNA (template + probe) existing at the end of the reaction;

$$\text{Mass of DNA (ng)} = \frac{[\mu\text{Ci added}][13.2][\% \text{incorporation}] + \text{starting template (ng)}}{\text{Specific activity of } [\alpha\text{-}^{32}\text{P}]\text{dCTP}}$$

For example, if you have labeled 25 ng of DNA using 50 µCi of [α - ^{32}P]dCTP at a specific activity of 3000 Ci/mmol, and your incorporation is 70%;

$$\text{Mass of DNA} = \frac{[50 \mu\text{Ci}][13.2][70\%] + 25}{3000} = \mathbf{40.4 \text{ ng}}$$

Next, calculate the amount of radioactivity incorporated during the reaction in dpm;

$$=[50 \mu\text{Ci}][2.2 \times 10^4][70\%] = \mathbf{7.7 \times 10^7 \text{ dpm}}$$

The specific activity of the labeled DNA can now be calculated as:

$$\frac{[\text{dpm incorporated}][10^3]}{\text{DNA mass}} \quad \text{ie} \quad \frac{[7.7 \times 10^7 \text{ dpm}][10^3]}{40.4 \text{ ng}}$$

$$= \mathbf{1.9 \times 10^9 \text{ dpm}/\mu\text{g}}$$

8.5. Removal of unincorporated nucleotides

Removal of unincorporated nucleotides is sometimes desirable to reduce background during hybridization, particularly if incorporation is less than 50%.

It is also considered important to remove these free nucleotides if the probe is being kept for several days before being used. However, if GE Healthcare's new Rapid-hyb™ buffer is used, purification is not required unless the probe is to be stored for more than 24 hours before use. Probes can be purified by Sephadex™ chromatography or selective precipitation (9, 10).




Fig 1. Time course of incorporation of [α - 32 P]dCTP (17 pmols) in a Rediprime reaction at 37°C using the control DNA supplied with the system.




Fig 2. Time course of incorporation of [α - 32 P]dCTP (17 pmols) in a Rediprime reaction at room temperature using the control DNA supplied with the system.

9. Related products

Labeling systems

Code

Megaprime DNA Labeling System for use with any radiolabeled nucleotide	RPN1604, RPN1605
Megaprime DNA Labeling System for use with radiolabeled dCTP	RPN1606, RPN1607
Nick Translation Kit	N5000, N5500
Ready-To-Go™ DNA Labeling Beads (-dCTP)	27-9240-01

Hybridization buffers

Code

Rapid-hyb buffer	RPN1635, RPN1636
------------------	------------------

Hybridization products

Code

Hybond™ range of nylon and nitrocellulose blotting membranes	
Hybridization oven/shaker (220/240 V 50 Hz)	RPN2510
Hybridization oven/shaker (110/115 V 60 Hz)	RPN2511E

Autoradiography products

Hyperfilm™ high performance autoradiography films.
Hypercassette™ and Hyperscreen™ available from stock.

Safety products

Radiation safety products for safe handling and storage of ³²P.

See our current catalogues or contact your local GE Healthcare office for further details.

10. References

1. FEINBERG, A.P. and VOGELSTEIN, B., *Anal. Biochem.*, **132**, pp.6-13, 1983.
2. FEINBERG, A.P. and VOGELSTEIN, B., *Addendum Anal. Biochem.*, **137**, pp.266-267, 1984.
3. SOUTHERN, E.M., *J. Mol. Biol.*, **98**, pp.503-517, 1975.
4. THOMAS, P.S., *Proc. Nat. Acad. Sci., USA.*, **77**, pp.5201-5205, 1980.
5. MEINKOTH, J. and WAHL, G., *Anal. Biochem.*, **138**, pp.267-284, 1984.
6. GRUNSTEIN, M. and HOGNESS, D.S., *Proc. Natl. Acad. Sci. USA.*, **72**, pp.3961-3965, 1975.
7. SAMBROOK, J. FRITSCH, E.F. and MANIATIS, T., *Molecular Cloning, a laboratory manual* [second edition], p.10. 16, Cold Spring Harbor Laboratory, 1989.
8. SAMBROOK, J. FRITSCH, E.F. and MANIATIS, T., *Molecular Cloning, a laboratory manual* [second edition], p.E19, Cold Spring Harbor Laboratory, 1989.
9. SAMBROOK, J. FRITSCH, E.F. and MANIATIS, T., *Molecular Cloning, a laboratory manual* [second edition], p.E18, Cold Spring Harbor Laboratory, 1989.
10. MUNDY, C.R., CUNNINGHAM, M.W. and READ, C.A., *Essential Molecular Biology: A Practical Approach Vol 2* [T.A. Brown, ed]. pp.57-109, Oxford University Press, Oxford, 1991.

GE Healthcare offices:

GE Healthcare Bio-Sciences AB
Björkgatan 30, 751 84 Uppsala,
Sweden

GE Healthcare Europe GmbH
Munzinger Strasse 5, D-79111 Freiburg,
Germany

GE Healthcare Bio-Sciences Corp.
800 Centennial Avenue, P.O. Box 1327,
Piscataway, NJ 08855-1327,
USA

GE Healthcare Bio-Sciences KK
Sanken Bldg. 3-25-1, Hyakunincho,
Shinjuku-ku, Tokyo 169-0073,
Japan

For contact information for your local office,
please visit: www.gelifesciences.com/contact


GE Healthcare UK Limited
Amersham Place
Little Chalfont, Buckinghamshire,
HP7 9NA, UK

<http://www.gehealthcare.com/lifesciences>



imagination at work

The next two pages are the
protocol card, please add to the
back page as a tear off addition



Amersham

Rediprime II DNA Labeling System

Product protocol card

RPN1633PC

Add denatured DNA template in a final volume of 45 μ l.



Add 5 μ l of [α - 32 P] dCTP.



Pipette up and down to mix.



Incubate for 10 minutes at 37°C.



Warning: For research use only.

Not recommended or intended for diagnosis of disease in humans or animals. Do not use internally or externally in humans or animals.

Caution: For use with radioactive material.

This product is to be used with radioactive material. Please follow the manufacturer's instructions relating to the handling, use, storage and disposal of such material.

All chemicals should be considered as potentially hazardous. We therefore recommend that this product is handled only by those persons who have been trained in laboratory techniques and that it is used in accordance with the principles of good laboratory practice. Wear suitable protective clothing, such as laboratory overalls, safety glasses and gloves. Care should be taken to avoid contact with skin or eyes. In the case of contact with skin or eyes, wash immediately with water. See material safety data sheet(s) and/or safety statement(s) for specific advice.

Storage and stability

Store Rediprime™II DNA Labeling System between 2–25°C. The system is stable for up to 6 months when stored under these recommended conditions.

Using the control DNA

The performance of Rediprime II can be checked by using the control DNA supplied.

The vial contains 300 ng of lambda HindIII DNA which should be dissolved in 60 µl TE buffer before use by pipetting up and down. Once reconstituted 5 µl of control DNA contains 25 ng template.

Ordering information

Rediprime II DNA Labeling System

Code

RPN1633

GE, imagination at work and GE monogram are trademarks of General Electric Company.

Amersham and Rediprime are trademarks of GE Healthcare companies

© 1995–2008 General Electric Company - All rights reserved.

Previously published 1995.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them.

A copy of these terms and conditions is available on request.

Contact your local GE Healthcare representative for the most current information.

<http://www.gehealthcare.com/lifesciences>



imagination at work

GE Healthcare UK Limited.
Amersham Place, Little Chalfont,
Buckinghamshire, HP7 9NA, UK