



# Product specification

Before using this product, please read the instructions below for safe handling, storage and disposal

## ECL blocking agent

RPN 2125

### Safety warnings and precautions

**Warning:** For research use only. Not recommended or intended for diagnosis of disease in humans or animals. Do not use internally or externally in humans or animals.

We recommend that this product and components are handled only by those persons who have been trained in laboratory techniques and that it is used in accordance with the principles of good laboratory practice. As all chemicals should be considered as potentially hazardous, it is advisable when handling chemical reagents to wear suitable protective clothing, such as laboratory overalls, safety glasses and gloves. Care should be taken to avoid contact with skin or eyes. In case of contact with skin or eyes, wash immediately with water.

Quantity            40g

For blocking nitrocellulose and PVDF membranes to be used with ECL and ECL Plus<sup>TM</sup> Western blotting detection reagents.

### Storage and stability

On receipt store the container at 2-8°C. After first use store at room temperature, ensuring that the lid is tightly closed. It is stable for at least six months when stored under the recommended conditions.

**Caution:** Moisture will affect the solubility properties of the powder. This reagent has been selected to give effective blocking of membranes in virtually all blotting contexts, if used as directed.

### Quality control

This product is tested according to our quality control systems using an ECL or ECL Plus Western blotting application.

# Protocol

## 1) Block the membrane

Non-specific binding sites are blocked by immersing the membrane in a 5%(w/v) solution of the blocking agent in PBS or TBS Tween<sup>TM</sup> (0.1% is usually sufficient) for 1 hour at room temperature on an orbital shaker. These conditions have been found in our laboratories to be sufficient for most purposes, but may need some optimization.

## 2) Washing

PBS or TBS Tween (0.1%) (PBS-T/TBS-T) is recommended as the wash buffer. Briefly rinse the membrane twice, then wash 3 times, once for 15 minutes and twice for 5 minutes at room temperature with fresh changes of wash buffer on an orbital shaker.

## 3) Dilution of the primary antibody

During the washing step, dilute the primary antibody. Incubate the washed membrane in the diluted antibody. 1 hour at room temperature is often sufficient, but this should be optimized for each antibody.

## 4) Washing

Wash the membrane as detailed above

## 5) Detection

Detect the binding of the primary antibody according to usual protocols.

## Notes

1) It is essential to optimize both primary and secondary antibodies for results with high signal and low background when using ECL and ECL Plus, due to the extreme sensitivity of the detection reagents.

2) As a general rule, as large a volume of washing buffer as possible should be used each time. It may be necessary to adjust the blocking conditions for certain applications.

3) Do not use azide as a preservative for buffers to be used in immunodetection if HRP linked antibodies are to be used, as it is an inhibitor of horseradish peroxidase.

# Related products

## ECL Plus Western blotting detection reagents

- sufficient for 3000cm<sup>2</sup> membrane RPN 2133
- sufficient for 1000cm<sup>2</sup> membrane RPN 2132

## ECL Western blotting detection reagents

- sufficient for 6000cm<sup>2</sup> membrane RPN 2134
- sufficient for 4000cm<sup>2</sup> membrane RPN 2106
- sufficient for 2000cm<sup>2</sup> membrane RPN 2209
- sufficient for 1000cm<sup>2</sup> membrane RPN 2109

## ECL Western blotting analysis system

RPN 2108

For the detection of either mouse or rabbit membrane bound primary antibodies with ECL detection reagents. Sufficient for 1000cm<sup>2</sup> membrane

ECL streptavidin-HRP and blocking reagent RPN 2195

Mouse IgG, horseradish peroxidase linked whole antibody (from sheep) NA 931

Rabbit IgG, horseradish peroxidase linked whole antibody (from donkey) NA 934

Rat IgG, horseradish peroxidase linked whole antibody (from sheep) NA 932

Human IgG, horseradish peroxidase linked whole antibody (from sheep) NA 933

Mouse IgG, horseradish peroxidase lined F(Ab')<sub>2</sub> fragment (from sheep) NA 9310

Rabbit IgG, horseradish peroxidase lined F(Ab')<sub>2</sub> fragment (from donkey) NA 9340

Rat IgG, horseradish peroxidase lined F(Ab')<sub>2</sub> fragment (from sheep) NA 9320

Human IgG, horseradish peroxidase lined F(Ab')<sub>2</sub> fragment (from sheep) NA 9330

Streptavidin horseradish peroxidase conjugate RPN 1231

Streptavidin biotinylated horseradish peroxidase complex RPN 1051

## ECL protein molecular weight markers

RPN 2107

Rainbow<sup>TM</sup> coloured protein molecular weight markers RPN 755

(molecular weight range 2350 - 46000

Rainbow coloured protein molecular weight markers RPN 756

(molecular weight range 14300 - 200000)

Full Range Rainbow recombinant protein molecular weight markers RPN 800

(molecular weight range 10000 - 250000)

Hybond<sup>TM</sup>-ECL nitrocellulose membrane RPN 2020D

Hybond PVDf membrane RPN 2020P

Hyperfilm<sup>TM</sup>-ECL RPN 2103

For details of sizes, availability and ordering information, please contact your local sales office.

ECL, ECL device, ECL Plus, Rainbow, Hybond and Hyperfilm are trademarks of Amersham Biosciences Limited or its subsidiaries  
Amersham and Amersham Biosciences is a trademark of Amersham plc

Tween is a trademark of ICI Americas Inc  
©Amersham Biosciences UK Limited 1999 - All rights reserved

All goods and services are sold subject to the terms and conditions of sale of the company within the Amersham Biosciences group which supplies them. A copy of these terms and conditions is available on request.

<http://www.amershambiosciences.com>

**Amersham Biosciences UK Limited** Amersham Place

Little Chalfont Buckinghamshire England HP7 9NA

**Amersham Biosciences AB** SE-751 84 Uppsala Sweden

**Amersham Biosciences Inc** 800 Centennial Avenue PO Box 1327 Piscataway NJ08855 USA

**Amersham Biosciences Europe GmbH** Munzinger Strasse 9 D-79111 Freiburg Germany

RPN2125PS/AB