

**Index:**

<b>1</b>	<b>SCOPE</b>	<b>2</b>
<b>2</b>	<b>INTENDED USE</b>	<b>3</b>
<b>3</b>	<b>PRODUCT DESCRIPTION</b>	<b>3</b>
<b>4</b>	<b>PRODUCT CLASSIFICATION</b>	<b>4</b>
<b>5</b>	<b>GENERAL REQUIREMENTS</b>	<b>5</b>
<b>6</b>	<b>STERILISATION METHOD</b>	<b>6</b>
<b>7</b>	<b>PROPERTIES</b>	<b>6</b>
<b>8</b>	<b>APPENDICES</b>	<b>9</b>
<b>9</b>	<b>DOCUMENT ADMINISTRATION</b>	<b>10</b>
<b>9.1</b>	<b>Amendment information</b>	<b>10</b>

## 1 Scope

The Product specification applies to the following products:

Local

Article No.	Product Name	Designation	Drawing
NJ-4606027	NORM-JECT® Luer Solo	2 ML NORM-JECT BBRAUN	2-1470
NJ-4606051	NORM-JECT® Luer Solo	5 ML NORM-JECT BBRAUN	5-1470
NJ-4606067	NORM-JECT® Luer Solo	5 ML NORM-JECT ZENTRISCH BBRAUN	5-3920
NJ-4606108	NORM-JECT® Luer Solo	10 ML NORM-JECT BBRAUN	10-1470
NJ-4606110	NORM-JECT® Luer Solo	10 ML NORM-JECT ZENTRISCH BBRAUN	10-3920
NJ-4606205	NORM-JECT® Luer Solo	20 ML NORM-JECT BBRAUN	20-1470
NJ-4606701	NORM-JECT® Luer Lock Solo	2 ML NORM-JECT LL BBRAUN	2-3420
NJ-4606710	NORM-JECT® Luer Lock Solo	5 ML NORM-JECT LL BBRAUN	5-3420
NJ-4606728	NORM-JECT® Luer Lock Solo	10 ML NORM-JECT LL BBRAUN	10-3420
NJ-4606736	NORM-JECT® Luer Lock Solo	20 ML NORM-JECT LL BBRAUN	20-3420
NJ-9166017	NORM-JECT®-F Luer Solo	1 ML NORM-JECT-F BBRAUN	1-Z0160

OEM

Article No.	Product Name	Designation	Drawing
NJ-4606027-02	NORM-JECT® Luer Solo	2 ML NORM-JECT BBRAUN US	2-1470
NJ-4606051-02	NORM-JECT® Luer Solo	5 ML NORM-JECT BBRAUN US	5-1470
NJ-4606067-02	NORM-JECT® Luer Solo	5 ML NORM-JECT ZENTRISCH BBRAUN US	5-3920
NJ-4606108-02	NORM-JECT® Luer Solo	10 ML NORM-JECT BBRAUN US	10-1470
NJ-4606110-02	NORM-JECT® Luer Solo	10 ML NORM-JECT ZENTRISCH BBRAUN US	10-3920
NJ-4606205-02	NORM-JECT® Luer Solo	20 ML NORM-JECT BBRAUN US	20-1470
NJ-4606701-02	NORM-JECT® Luer Lock Solo	2 ML NORM-JECT LL BBRAUN US	2-3420
NJ-4606710-02	NORM-JECT® Luer Lock Solo	5 ML NORM-JECT LL BBRAUN US	5-3420
NJ-4606728-02	NORM-JECT® Luer Lock Solo	10 ML NORM-JECT LL BBRAUN US	10-3420
NJ-4606736-02	NORM-JECT® Luer Lock Solo	20 ML NORM-JECT LL BBRAUN US	20-3420
NJ-9166017-02	NORM-JECT®-F Luer Solo	1 ML NORM-JECT-F BBRAUN US	1-Z0160

## 2 Intended use


Single-use syringes, 2-piece

## 3 Product description

A syringe is a device that predominantly consists of a plunger which fits into a tube (“barrel”). The plunger can be pulled and pushed along inside the tube (“barrel”), allowing the syringe to withdraw or expel a fluid or gas. The open end of the syringe can be fitted with a hypodermic needle, a nozzle, or tubing to help direct the flow into and out of the barrel.

The mechanism of a syringe can be activated by applying a manually physical pressure on the plunger.

### 3.1 Exploded assembly drawing



### 3.2 Material

Individual components	Material		
<b>Syringe</b>	Barrel	Polypropylene	PP
	Plunger	Polyethylene	PE
<b>Packaging</b>	Paper	Medical Grade Paper	Paper
	Film	Polypropylene/Polyamide/Polyethylene	PP/PA/PE

**3.3 Basic Product Characteristics**

The single use syringes contain for their function and purpose the following elements:

- Sizes available: 1 mL - 20 mL
- Highly transparent barrel
- Permanent marking
- Good readability
- 1 mL with displacement spike: reduces the residual volume and minimizes unnecessary loss of medicament
- Safe plunger backstop
- Silicone oil free
- Latex-free
- According to ISO 7886-1
- Available with Luer Slip (1 mL - 20 mL) or Luer Lock (2 mL - 20 mL)

**3.4 Packaging**

Primary Packaging:	Individual single peel pack designated as a microbiological barrier to assure the sterility of the product.
Secondary Packaging:	Box containing a certain number of primary packaging. 100 pcs. per box.
Transport Packaging:	Case for dispatching and additional protection against mechanical damages during transportation.

**4 Product classification<sup>1</sup>**

Classification according to the Council Directive 93/42/EEC concerning medical devices Annex IX:

Non-invasive medical devices - Class I sterile, rule 2, additional control function  
(class Ism = class I with measurement function):

Single-use syringes, 2-piece (without needle)

- Standard syringes:  
Product groups: *Luer Solo*<sup>2</sup> (e.g. NORM-JECT® Luer Solo)  
*Luer Lock Solo* (e.g. NORM-JECT® Luer Lock Solo)
- Fine dosage syringes (for precise dosage of smallest volumes):  
Product groups: *F Luer Solo* (e.g. NORM-JECT®-F Luer Solo)

<sup>1</sup> The Classification has to occur from the party responsible for manufacture

<sup>2</sup> Solo: without needle

## 5 General requirements

The products fulfill the requirements of the standards in force at the time of manufacture.

### Biological requirements

EN ISO 10993-1	Biological evaluation of medical devices - Part 1: Evaluation and testing
EN ISO 10993-4	Biological evaluation of medical devices - Part 4: Selection of tests for interactions with blood
EN ISO 10993-5	Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity
ISO 10993-10	Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization
EN ISO 10993-11	Biological evaluation of medical devices - Part 11: Tests for systemic toxicity
EN ISO 11737-1	Sterilization of medical devices - Microbiological methods - Part 1: Determination of a population of microorganisms on products
DIN 58953-6	Sterilization - Sterile supply - Part 6: Microbial barrier testing of packaging materials for medical devices which are to be sterilized
Ph. Eur., chapter 2.6.8	European Pharmacopoeia Pyrogens
chapter 2.6.14	Bacterial Endotoxins

### Chemical requirements

ISO 7886-1	Sterile hypodermic syringes for single use - Part 1: Syringes for manual use
EN ISO 10993-18	Biological evaluation of medical devices - Part 18: Chemical characterization of materials
Ph. Eur., chapter 3.2.8	European Pharmacopoeia Sterile Single-Use Plastic Syringes

### Physical-technical requirements

ISO 7886-1	Sterile hypodermic syringes for single use - Part 1: Syringes for manual use
EN 1707	Conical fittings with a 6 % (Luer) taper for syringes, needles and certain other medical equipment - Lock fittings
EN 20594-1	Conical fittings with a 6 % (Luer) taper for syringes, needles and certain other medical equipment - Part 1: General requirements

**Packaging**

EN 868-5	Packaging for terminally sterilized medical devices Part 5: Sealable pouches and reels of porous materials and plastic film construction - Requirements and test methods
EN ISO 11607-1	Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems
EN ISO 11607-2	Packaging for terminally sterilized medical devices - Part 2: Validation requirements for forming, sealing and assembly processes
ISO 15223-1	Medical devices – Symbols to be used with medical device labels, labeling and information to be supplied - Part 1: General requirements
DIN 55529	Verpackung – Bestimmung der Siegelnahtfestigkeit von Siegelungen aus flexiblen Packstoffen
ASTM F 88	Standard Test Method for Seal Strength of Flexible Barrier Materials
ASTM F 1929	Standard Test Method for Detecting Seal Leaks in Porous Medical Packaging by Dye Penetration
ASTM D 4169	Standard Practice for Performance Testing of Shipping Containers and Systems
ASTM F 2096	Standard Test Method for Detecting Gross Leaks in Medical Packaging by Internal Pressurization (Bubble Test)

**Sterilization**

EN 556-1	Sterilization of medical devices - Requirements for medical device to be designated "STERILE - Part 1: Requirements for terminally sterilized medical devices
EN ISO 11135	Sterilization of health care products - Ethylene oxide - Requirements for development, validation and routine control of a sterilization process for medical devices

**Delivery**

The products are supplied in individual, sterile packs.

**Shelf life**

The product shelf life is 60 months.

**Storage conditions**

EN 1041	Medical devices; Information supplied by the manufacturer No specific storage conditions Storage as commonly used for medical products
---------	--

**Transport conditions**

EN 1041	Medical devices; Information supplied by the manufacturer No specific storage conditions Storage as commonly used for medical products
---------	--

**6 Sterilization method**

The product is EO sterilized.

The EO residuals are according to EN ISO 10993-7.

## 7 Properties

Properties	Standard
<u>Biology</u> - Haemolysis - Cytotoxicity - Sensitization - Intracutaneous Reactivity - Systemic toxicity (acute) - Microbial contamination - Bacterial endotoxins (LAL) - Pyrogens - Microbial barrier properties (in case of wetness) - Sterile - Sterilization residuals	EN ISO 10993-4 EN ISO 10993-5 ISO 10993-10 ISO 10993-10 EN ISO 10993-11 EN ISO 11737-1 Ph. Eur., chapter 2.6.14 Ph. Eur., chapter 2.6.8 EN ISO 11607-1 / DIN 58953-6 EN 556-1 / EN ISO 11135 EN ISO 10993-7
<u>Chemistry</u> - Extractable metals - Chemical characterization of materials - Appearance of solution - Acidity or alkalinity - Absorbance - Reducing substances - Transparency/Opaescence	ISO 7886-1 EN ISO 10993-18 Ph. Eur., chapter 3.2.8 Ph. Eur., chapter 3.2.8 Ph. Eur., chapter 3.2.8 Ph. Eur., chapter 3.2.8 Ph. Eur., chapter 3.2.8
<u>Technical</u> <u>Tightness</u> - Air tightness of syringe at vacuum - Fluid tightness of syringe with deflected plunger under compression - Water tightness of syringe and needle during injection - Air tightness of conical fitting assembly at vacuum - Fluid tightness of conical fitting assembly under compression <u>Tensile / Pressure Strength</u> - Breakaway force of syringe with water - Sliding force of syringe with water - Breakaway force of syringe with air - Sliding force of syringe with air - Secondary starting force of syringe with air - Separation force of plunger out of barrel - Separation force of fitting assembly - Flexural strength of barrel finger grips	ISO 7886-1, Annex B ISO 7886-1, Annex D -/- EN 20594-1 / EN 1707 EN 20594-1 / EN 1707 ISO 7886-1, Annex E ISO 7886-1, Annex E -/- -/- -/- -/- EN 20594-1 / EN 1707 -/-

Properties	Standard
<i>Visual</i>	
- Contamination of device	ISO 7886-1
- Damages of device	-/-
- Completeness of device	-/-
- Shaping of device	-/-
- Coloration of single parts and graduation	-/-
- Graduated scale of syringe	ISO 7886-1
- Wipe resistance of graduation	-/-
- Material projections of device	-/-
- Inclusions / flakes / bubbles / flow marks of single parts	-/-
- Transparency of syringe barrel	-/-
<i>Function</i>	
- Maximum dead space	ISO 7886-1, Annex C
- Tolerance on graduated capacity	ISO 7886-1
- Rolling behavior at 10°	ISO 7886-1
- Fit of plunger in barrel	ISO 7886-1
- Patency of syringe lumen	-/-
- Unscrewing torque of fitting assembly	EN 1707
- Ease of assembly of fitting assembly	EN 1707
- Resistance to overriding of fitting assembly	EN 1707
- Stress cracking of fitting assembly	EN 20594-1 / EN 1707
<i>Dimensional accuracy</i>	
- Luer conical fitting (male/female)	EN 20594-1 / EN 1707
- Position of scale - zero graduation line	ISO 7886-1
- Maximal useable capacity	-/-
- Minimum length of the plunger from the surface of the finger grips nearer to the push-button	ISO 7886-1
- Diameter of nozzle lumen	ISO 7886-1

<b>Properties</b>	<b>Standard</b>
<u>Packaging</u>	
<i>Tightness</i>	
- Air tightness of primary packaging (bubble test)	ASTM F 2096
- Fluid tightness seal seam of primary packaging (blue dye test)	ASTM F1929
<i>Tensile / Pressure Strength</i>	
- Strength of sealing seam, primary packaging	DIN 55529 / ASTM F 88
<i>Visual</i>	
- Labeling of all packages	ISO 7886-1
- Printing quality of all packages	ISO 7886-1
- Cleanliness of primary packaging	ISO 7886-1
- Contamination of all packages except primary packaging	ISO 7886-1
- Damages of all packages	-/-
- Closure of all packages	-/-
- Completeness of all packages	-/-
<i>Function</i>	
- Peelability of sealing seam, primary packaging	EN 868-5, Annex E
- Separation of primary packaging	-/-
<i>Dimensional accuracy</i>	
- Width of sealing seam, primary packaging	EN 868-5, Annex E

## 8 Appendices

none

– End of document –

## 9 Document administration

### 9.1 Amendment information

Version	Description of the changes
3.0	OEM devices (NORM-JECT® US) were added
2.0	ISO 868-5 was corrected into EN 868-5; 'Incompatibility with injection fluids' was deleted acc. to HC-CHC-ALMO-463; Feature 'Chemical characterization of materials' was moved from Biology into Chemistry; 'Cleanliness of device' was renamed into 'Contamination of device'; 'Cleanliness of all packages' was divided into 'Cleanliness of primary packaging' and 'Contamination of all packages except primary packaging'; Following features were deleted acc. to HC-CHC-ALMO-569 – harmonization with Product Risk Analysis: 'Flexural strength of Luer cone', 'Dimension according to drawing', 'Minimal film thickness of primary packaging'
1.0	New specification

Title: Two-piece Syringe (NORM-JECT®) - Product Specification Initiator: Martina ? Schreiber

This document is signed electronically in compliance with the B. Braun electronic signature policies and procedures by following persons:

UserName: Schreiber, Martina (schrmtde)  
Title: HC-PM-DE01 - Engineering  
Date: Wednesday, 06 November 2019, 10:36 W. Europe Daylight Time  
Meaning: Document signed as Author  
=====

UserName: Wende, Joerg (wendjede)  
Title: HC-RD-DE01 Head of developmentcenter syringes  
Date: Wednesday, 06 November 2019, 15:48 W. Europe Daylight Time  
Meaning: Approve Document  
=====

UserName: Henke, Gudrun (henkgude)  
Title: HC-RA Senior Manager Medical Devices CoE IVS  
Date: Friday, 08 November 2019, 18:34 W. Europe Daylight Time  
Meaning: Approve Document  
=====

UserName: Brand, Thomas (brantode)  
Title: HC-QM-DE08 Vice President QM for non-active Medical Devices  
Date: Tuesday, 12 November 2019, 09:29 W. Europe Daylight Time  
Meaning: Approve Document  
=====