



MessageBOOSTER™ cDNA Synthesis Kits

Do more with less - Produce cDNA from precious total RNA or cells

- **Fast:** Amplify RNA and produce cDNA in only 1 day
- **Unbiased Amplification:** Preserves the relative transcript abundance of a sample with a high fidelity, linear RNA amplification process
- **Increased Sample Usability:** Use samples you couldn't before, produce enough cDNA for multiple assays, and have enough left over for archiving
- **Sensitive:** Readily and reproducibly detect even low-abundance transcripts in RNA from a single cell or 10 pg of total RNA (CT values <35 cycles)

epicentre®
Exclusively available
thru Lucigen.

The MessageBOOSTER™ cDNA Synthesis for qPCR (Cat. No. MB060124) and the MessageBOOSTER™ cDNA Synthesis from Cells Lysates Kit (Cat. No. MBCL90310) enable users to perform sensitive qPCR assays using either limiting amounts of total RNA (as low as 10 pg) or from 1 to 500 cells, respectively.

The kits amplify the mRNA contained in small total RNA/cell samples and then converts the amplified RNA to cDNA that is ready for PCR or qPCR (Figure 1). The MessageBOOSTER cDNA Synthesis Kits both use oligo(dT) to prime cDNA synthesis, and as a result is best used with intact total RNA samples or cell samples.




Figure 1. Protocol overview of the MessageBOOSTER™ cDNA Synthesis Kits. Each MessageBOOSTER kit reaction amplifies the poly(A) RNA (mRNA) directly from either (A) purified total RNA or (B) crude cell lysates. The amplified RNA is then reverse-transcribed to cDNA that can be diluted up to 1,000-fold for qPCR.



Make Enough cDNA from Limiting Samples for Many qPCR Assays

Total RNA Amount or Cell Number Inputs	Number of qPCR Assays Possible	
	Low ^a to Medium ^b Abundance Transcripts	Medium ^b to High ^c Abundance Transcripts
10 pg (~1 cell)	≥10	≥100
100 pg (~10 cells)	≥100	≥1,000
500 pg (~50 cells)	≥500 – 1,000	≥5,000 – 10,000

Table 1. The number of qPCR assays that can be performed using cDNA produced by the MessageBOOSTER™ cDNA Synthesis Kits. The number of qPCR assays is dependent on the amount of total RNA used in each MessageBOOSTER reaction and the abundance of the target transcript(s). Similar numbers of assays are possible using the MessageBOOSTER cDNA Synthesis from Cell Lysates Kit and the indicated number of cells.

^a Low Abundance Transcripts = 1-1,000 copies per cell

^b Medium Abundance Transcripts = 1,000-10,000 copies per cell

^c High Abundance Transcripts = 10,000-100,000 copies per cell

Production of Sufficient cDNA for 1000's of qPCR Assays from a Single Cell




Figure 2. qPCR results from cDNA produced from a single NRK cell.

cDNA was made from a single NRK cell using the MessageBOOSTER™ cDNA Synthesis from Cell Lysates Kit and qPCR of the PBGD target was performed using undiluted (red), 1:10 diluted (green), 1:100 diluted (blue), and 1:1,000 diluted (purple) aliquots. The low-abundance PBGD transcript was readily detected.

Products	Size	Cat. No.
MessageBOOSTER™ cDNA Synthesis for qPCR	24 rxns	75927-956
MessageBOOSTER™ cDNA Synthesis from Cell Lysates Kit	10 rxns	75927-958

COMPONENTS

Both kits contain: MessageBOOSTER™ T7-Oligo(dT) Primer, RNase Inhibitor, MessageBOOSTER™ Reverse Transcription PreMix, MessageBOOSTER™ DNA Polymerase PreMix, MessageBOOSTER™ DNA Polymerase, MessageBOOSTER™ cDNA Finishing Solution, MessageBOOSTER™ Random Primers, MessageBOOSTER™ RNase H, MMLV Reverse Transcriptase, MessageBOOSTER™ In Vitro Transcription PreMix, MessageBOOSTER™ T7 RNA Polymerase, MessageBOOSTER™ T7 Transcription Buffer, RNase-Free DNase I, Forward and Reverse Control PCR Primers, DTT, RNase-Free Water, Poly(I), NRK Total RNA Control.

The MessageBOOSTER cDNA Synthesis from Cell Lysates Kit also contains: QuickExtract™ RNA Extraction Solution, Primer Mix, and Positive Control.

F209-VWR0917 For research or investigational use only. ISO 13485 Certified