

GE Healthcare

# Amersham Megaprime DNA Labeling Systems

Product Booklet

Codes:     RPN1604  
              RPN1605  
              RPN1606  
              RPN1607



# Page finder

1. Legal	3
2. Handling	4
2.1. Safety warnings and precautions	4
2.2. Storage and stability	4
2.3. Quality control	4
3. System components	6
3.1. Megaprime DNA Labeling Systems	8
4. Introduction	9
5. Megaprime DNA labeling protocols	11
5.1. Standard Megaprime protocol	11
5.2. New Megaprime protocol	15
5.3. Use of alternative reaction conditions	20
6. Appendices	27
6.1. Appendix I. Labeling of DNA fragments in low melting point agarose	27
6.2. Appendix II. Monitoring the reaction and calculating the specific activity of the labeled DNA	28
6.3. Appendix III. Removal of unincorporated nucleotides	32
6.4. Appendix IV. Additional equipment and reagents	34
7. Troubleshooting guide	35
8. References	38
9. Related products	39

# 1. Legal

GE, imagination at work and GE monogram are trademarks of General Electric Company.

Amersham, Megaprime, Hybond, Hyperfilm, Hypercassette, Hyperscreen, Rapid-hyb, Ready-To-Go, Rediprime, Sensitize, Sephadex and TrackerTape are trademarks of GE Healthcare companies.

All third party trademarks are the property of their respective owners.

© 2002–2010 General Electric Company – All rights reserved.  
First published 2002

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them.

A copy of these terms and conditions is available on request Contact your local GE Healthcare representative for the most current information.

<http://www.gelifesciences.com>

GE Healthcare UK Limited.  
Amersham Place, Little Chalfont,  
Buckinghamshire, HP7 9NA  
UK

## 2. Handling

### 2.1. Safety warnings and precautions

**Warning: For research use only.** Not recommended or intended for diagnosis of disease in humans or animals. Do not use internally or externally in humans or animals.

**Caution: For use with radioactive material.**

This product is to be used with radioactive material. Please follow the manufacturer's instructions relating to the handling, use, storage and disposal of such material.

All chemicals should be considered as potentially hazardous. We therefore recommend that this product is handled only by those persons who have been trained in laboratory techniques and that it is used in accordance with the principles of good laboratory practice. Wear suitable protective clothing such as laboratory overalls,

safety glasses and gloves. Care should be taken to avoid contact with skin or eyes. In the case of contact with skin or eyes, wash immediately with water. See material safety data sheet(s) and/or safety statement(s) for specific advice.

### 2.2. Storage and stability

Upon receipt of these systems components should be stored at -15°C to -30°C. The components are stable for at least 3 months when stored under recommended conditions.

### 2.3. Quality control

The Megaprime™ DNA Labeling Systems are tested by our quality control group to ensure an incorporation rate greater than 55% after 10 minutes at 37°C.

The performance of RPN1604/1605 is tested with the standard DNA provided

using 17 pmol/25 ng DNA of [ $\alpha$ - $^{32}$ P] labeled nucleotides, specific activity 3000 Ci/mmol and RPN 1606/1607 are tested using 17 pmol/25 ng DNA of [ $\alpha$ - $^{32}$ P]dCTP, 3000 Ci/mmol. Incorporations greater than 55% are achieved after 10 minutes incubation at 37°C, as assayed by thin-layer chromatography on PEI cellulose in 1.25M  $\text{KH}_2\text{PO}_4$ , pH3.4.

In addition components of the kits are checked for identity by HPLC and the DNA solutions for concentration by UV spectrophotometry.

### 3. System components

Megaprime DNA labeling	RPN1604	RPN1605	RPN1606	RPN1607
Primer solution: Random nonamer primers in an aqueous solution	150 µl	300 µl	150 µl	300 µl
Labeling buffer; dATP, dGTP and dTTP in Tris/HCl pH7.5, 2-mercaptoethanol and MgCl <sub>2</sub>	–	–	300 µl	600 µl
Nucleotide solutions				
(a) dATP	120 µl	240 µl	–	–
(b) cCTP	120 µl	240 µl	–	–
(c) dGTP	120 µl	240 µl	–	–
(d) dTTP	120 µl	240 µl	–	–
in Tris/HCl pH8.0, 0.5 mM EDTA				
Reaction buffer: A 10× concentrated buffer containing Tris/HCl pH7.5, 2-mercaptoethanol and MgCl <sub>2</sub>	150 µl	300 µl	–	–

<b>Megaprime DNA labeling</b>	<b>RPN1604</b>	<b>RPN1605</b>	<b>RPN1606</b>	<b>RPN1607</b>
Enzyme solution; 1 unit/ $\mu$ l DNA polymerase 1 Klenow fragment (cloned in 100 mM potassium phosphate pH6.5, 10 mM 2-mercapto- ethanol and 50% glycerol	60 $\mu$ l	120 $\mu$ l	60 $\mu$ l	120 $\mu$ l
Standard DNA solution; 5 ng/ $\mu$ l <i>Hind</i> III digested lambda DNA in 10 mM Tris/HCl pH 8.0, 1 mM EDTA	25 $\mu$ l	50 $\mu$ l	25 $\mu$ l	50 $\mu$ l
Carrier DNA solution; 500ng/ml sonicated fish sperm DNA in 10 mM Tris/HCl pH 8.0, 1 mM EDTA	1.25 ml	2.5 ml	1.25 ml	2.5 ml

### 3.1. Megaprime DNA Labeling Systems

30 standard labeling reactions –  
for use with any radioactive nucleotide RPN1604

60 standard labeling reactions –  
for use with any radioactive nucleotide RPN1605

30 standard labeling reactions –  
for use with radioactively labeled dCTP RPN1606


60 standard labeling reactions –  
for use with radioactively labeled dCTP RPN1607



## 4. Introduction

Feinberg and Vogelstein (1,2) introduced the use of random sequence hexanucleotides to prime DNA synthesis on denatured template DNA at numerous sites along its length. The primer-template complex is a substrate for the 'Klenow' fragment of DNA polymerase 1. By substituting a radiolabeled nucleotide for a non-radioactive equivalent in the reaction mixture newly synthesized DNA is made radioactive (see Figure 1). The absence of the 5'-3' exonuclease activity associated with DNA polymerase 1 ensures that labeled nucleotides incorporated by the polymerase are not subsequently removed as monophosphates. Very small amounts of input DNA can be labeled, enabling very high specific activity DNA probes to be produced with relatively small quantities of added nucleotides. These radioactive labeled fragments can then be used as sensitive hybridization probes for a wide range of filter based applications (3-6).

Previous protocols for the random primer labeling of DNA have required reaction times of at least 30 minutes. GE Healthcare's Megaprime DNA Labeling System allows the labeling of template DNA to the same high specific activity but at a greatly accelerated rate. Probes of specific activity  $1.9 \times 10^9$  dpm/ $\mu$ g can be produced with the majority of DNA substrates, using the standard protocol, after 10 minutes incubation at 37°C. This rapid labeling is achieved by the use of nonamer primers rather than the conventional hexamers (Figure 1). Nonamers allow for more efficient priming from the template DNA at 37°C, resulting in fast and efficient labeling of the DNA. A new alternative protocol has further reduced the variability in labeling which can occur with DNA template from a variety of sources. Both the standard Megaprime protocol and the new protocol are given as options in this booklet. The labeling of DNA in low melting point agarose takes only 15-30 minutes in contrast to conventional systems where overnight incubation are necessary.



**Figure 1.** Preparation of labeled probes using GE Healthcare's Megaprime DNA Labeling Systems.

## 5. Megaprime DNA labeling protocols

The Megaprime systems allow DNA from a variety of sources to be labeled *in vitro* to high specific activity with  $^{32}\text{P}$  and other radionuclides. The specific activity of the probes generated by these systems will vary according to the specific activity of the labeled dNTP used.

The standard Megaprime protocol is presented, together with a new protocol which reduces the variation in labeling efficiency that can occur with DNA template from a variety of sources.

The protocols given here are for use with 17 pmol[ $\alpha$ - $^{32}\text{P}$ ]dNTP, specific activity 3000 Ci/mmol. For alternative reaction conditions refer to page 20.

DNA prepared by standard minilysate methods may be used in either protocol. DNA solutions which are too dilute to be used directly should be concentrated by ethanol precipitation followed by redissolution in an appropriate volume of water or 10 mM Tris/HCl, pH 8.0, 1 mM EDTA. DNA in restriction enzyme buffers may be added directly to the reaction. The reaction can also be performed with DNA in agarose gel slices (see note 3 and Appendix 1).

### 5.1. Standard Megaprime protocol

Protocol	Notes
1. Dissolve the DNA to be labeled to a concentration of 2.5–25 ng/ $\mu\text{l}$ in either distilled water or 10 mM Tris/HCl, pH 8.0, 1 mM EDTA (TE buffer).	1. If desired, the labeling efficiency of a DNA sample can be compared with that of the standard DNA supplied with the kit. In this case 5 $\mu\text{l}$ of standard DNA should be used.

## Protocol

## Notes

---

2. Place the required tubes from the Megaprime system, with the exception of the enzyme, at room temperature to thaw. Leave the enzyme at  $-15^{\circ}\text{C}$  to  $-30^{\circ}\text{C}$  until required, and return immediately after use.
3. Place 25 ng of template DNA into a microcentrifuge tube and to it add 5  $\mu\text{l}$  of primers and the appropriate volume of water to give a total volume of 50  $\mu\text{l}$  in the final Megaprime reaction. Denature by heating to  $95\text{--}100^{\circ}\text{C}$  for 5 minutes in a boiling water bath.
4. Spin briefly in a microcentrifuge to bring the contents to the bottom of the tube.
5. Keeping the tube at room temperature, add the nucleotides and reaction buffer (RPN1604/5) or the labeling buffer (RPN1606/7) followed by the radiolabeled dNTP(s) and enzyme as follows:
3. When labeling DNA in low melting point agarose, first place the tube containing the stock DNA in a boiling water bath for 30 seconds to melt the agarose before removing the required volume. The volume of low melting point agarose DNA should not exceed 25  $\mu\text{l}$  in a 50  $\mu\text{l}$  reaction.
5. The reaction volume may be scaled up or down if more or less than 25 ng of DNA is to be labeled.

**Protocol****Notes**

Component	RPN1604/5	RPN1606/7
Labeling buffer		10 $\mu$ l
Unlabeled dNTPs	4 $\mu$ l of each omitting those to be used as label	-
Reaction buffer	5 $\mu$ l	-
Radiolabeled (dNTP)	5 $\mu$ l	5 $\mu$ l (dCTP)
Enzyme	2 $\mu$ l	2 $\mu$ l

6. Mix gently by pipetting up and down and cap the tube. Spin for a few seconds in a microcentrifuge to bring the contents to the bottom of the tube.

7. Incubate at 37°C for 10 minutes

6. Avoid vigorous mixing of the reaction mixture as this can cause severe loss of enzyme activity.

7. Purified DNA can be labeled to high specific activity in 10 minutes at 37°C but, if desired, can be labeled for up to 1 hour at this temperature. When labeling DNA in low melting point agarose, longer incubation of 15–30 minutes at 37°C are required for optimum labeling. Longer incubation

## Protocol

## Notes

7. Incubate at 37°C for 10 minutes *continued*.

8. Stop the reaction by the addition of 5 µl of 0.2 M EDTA. For use in a hybridization, denature the labeled DNA by heating to 95–100°C for 5 minutes, then chill on ice.

7. ***Continued.***

times (up to 60 minutes) are required when nucleotide analogues (e.g. [<sup>35</sup>S]dNTPαS) are used.

8. Labeled probe can be stored at -15°C to -30°C in a non frost-free freezer. Prolonged storage of <sup>32</sup>P-labeled probes can lead to substantial probe degradation (7). High specific activity probes should be stored for no longer than 3 days. Although probe purification is not usually necessary for most membrane applications, the removal of unincorporated nucleotide is sometimes useful to reduce background in filter hybridizations for probes >10<sup>9</sup> dpm/µg or when the reaction yields an incorporation of less than 50%. This procedure is described in Appendix III. Calculation of probe specific activity is described in Appendix II. Extensive experimentation with **Rapid-hyb buffer (RPN1635/6)** has shown that probe purification, even

Protocol	Notes
<b>8.</b> Stop the reaction by the addition of 5 $\mu$ l of 0.2 M EDTA. For use in a hybridization, denature the labeled DNA by heating to 95–100°C for 5 minutes, then chill on ice <i>continued</i> .	<b>8. <i>Continued</i></b> under the conditions given above is not required with the isotopes $^{32}\text{P}$ and $^{33}\text{P}$ . Purification of $^{35}\text{S}$ labeled probes is however required to reduce filter background.

## 5.2. New Megaprime protocol

Protocol	Notes
<b>1.</b> Dilute the DNA to a concentration of 5 ng/ $\mu$ l in either distilled water or 10 mM TE buffer.	<b>1.</b> DNA solutions at concentrations in the range 5–25 ng/ $\mu$ l can be used if desired. However the denaturing volume (step 3) should not be less than 10 $\mu$ l to maximize the efficiency of primer annealing. The labeling efficiency of a DNA sample can be compared with that of the standard DNA supplied with the kit. In this case 5 $\mu$ l of standard DNA should be used.
<b>2.</b> Place the required tubes from the Megaprime system with the exception of the enzyme at room temperature to thaw. Leave the enzyme at -15°C to -30°C until required, and return immediately after use.	

## Protocol

## Notes

3. Place 25 ng (5  $\mu$ l) of template DNA into a clean microcentrifuge tube and to it add 5  $\mu$ l of primers. Denature by heating to 95–100°C for 5 minutes in a boiling water bath.

4. Spin briefly in a microcentrifuge to bring the contents to the bottom of the tube.

5. Keeping the tube at room temperature add the nucleotides and 10 $\times$  reaction buffer (RPN1604/5) or the labeling buffer (RPN1606/7), water and enzyme:-

Component	RPN1604/5	RPN1606/7
Labeling buffer		10 $\mu$ l
Unlabeled dNTPs	4 $\mu$ l of each – omitting those to be used as label	
Reaction buffer	5 $\mu$ l	–
Enzyme	2 $\mu$ l	2 $\mu$ l
Water*	as appropriate for a final reaction volume of 50 $\mu$ l*	

3. If the volume of DNA and primers is less than 10  $\mu$ l make up to this volume with water. When labeling DNA in low melting point agarose first place the tube containing the stock DNA in a boiling water bath for 30 seconds to melt the agarose before removing the required volume. The volume of low melting point agarose DNA should not exceed 25  $\mu$ l in a 50  $\mu$ l reaction.

5. The enzyme can be added directly to the reaction mix or pipetted on to the side of the microcentrifuge tube and “washed” down with the water.

\*When calculating this volume remember to allow for the volume of radioactive nucleotide to be added.



## Protocol


## Notes

---

6. Cap the tube and spin for a few seconds in a microcentrifuge to bring the contents to the bottom of the tube.
  7. Add the radiolabeled dNTP, for example 5  $\mu$ l [ $\alpha$ - $^{32}$ P]dNTP, specific activity 3000 Ci/mmol. Mix by gently pipetting up and down. Spin for a few seconds in a microcentrifuge to bring the contents to the bottom of the tube.
  8. Incubate at 37°C for 10 minutes.
7. Avoid vigorous mixing of the reaction mixture as this can cause severe loss of enzyme activity.
  8. Purified DNA can be labeled to high specific activity in 10 minutes at 37°C but, if desired can be labeled for up to 1 hour at this temperature.  
When labeling DNA in low melting point agarose, longer incubation of 15–30 minutes at 37°C are required for optimum labeling. Longer incubation times (up to 60 minutes) are required when nucleotide analogues (e.g. [ $^{35}$ S]dNTP(S) are used.

Protocol	Notes
<p>9. Stop the reaction by the addition of 5 µl of 0.2 M EDTA. For use in a hybridization, denature the labeled DNA by heating to 95–100°C for 5 minutes, then chill on ice.</p>	<p>9. Labeled probe can be stored at -15°C to -30°C in a non frost-free freezer. Prolonged storage of <sup>32</sup>P-labeled probes can lead to substantial probe degradation(7). High specific activity probes should be stored for no longer than 3 days. Although probe purification is not usually necessary for most membrane applications the removal of unincorporated nucleotide is sometimes useful to reduce background in filter hybridizations for probes &gt;10<sup>9</sup> dpm/µg or when the reaction yields an incorporation of less than 50%. This procedure is described in Appendix III. Calculation of probe specific activity is described in Appendix II. Extensive experimentation with <b>Rapid-hyb buffer (RPN1635/6)</b> has shown that probe purification, even under the conditions given</p>

Protocol	Notes
9. Stop the reaction by the addition of 5 $\mu$ l of 0.2 M EDTA. For use in a hybridization, denature the labeled DNA by heating to 95–100°C for 5 minutes, then chill on ice <i>continued</i> .	9. <i>Continued</i> above is not required with the isotopes $^{32}\text{P}$ and $^{33}\text{P}$ . Purification of $^{35}\text{S}$ labeled probes is however required to reduce filter background.



**Figure 2.** Time course of incorporation of  $[\alpha\text{-}^{32}\text{P}]\text{dCTP}$  (17 pmoles) in a Megaprime reaction at 37°C. The DNA used was the standard DNA supplied with the system.

### 5.3. Use of alternative reaction conditions

#### a. Use of more than one labeled $[\alpha\text{-}^{32}\text{P}]\text{dNTP}$ .

Table 1 lists the results of a selection of standard reactions, using a variety of input labels under optimum conditions. Figure 3 gives more complete information on their use in Megaprime reactions. Reactions were carried out at 37°C for 5 minutes.

#### b. Use of alternative specific activity $[\alpha\text{-}^{32}\text{P}]\text{dNTPs}$ .

When using  $[\alpha\text{-}^{32}\text{P}]\text{dNTPs}$  of specific activity  $<3000$  Ci/mmol the incubation time should be extended to 1 hour at 37°C.

### **c. Use of [<sup>35</sup>S]dNTPαS.**

When using <sup>35</sup>S-labeled radionucleotides the incubation time should be extended to 1 hour at 37°C.

### **d. Labeling at room temperature.**

If desired, labeling reactions can be carried out at room temperature. Maximum incorporation occurs after an incubation time of 45–60 minutes. A decline in incorporation can be observed if reactions are left overnight.

### **e. Factors affecting the labeled DNA.**

#### **1. Specific activity**

Figure 3a should be used to ascertain the number and quantity of labeled dNTP's required in order to prepare a probe of the desired specific activity.

#### **2. Efficiency**

Figure 3b indicates the efficiency of the chosen reaction conditions, and thus permits a balance of specific activity and economy.

#### **3. Probe length**

Figure 3c gives a measure of mean probe lengths obtained under standard conditions. Probe lengths were measured by denaturing agarose gel electrophoresis followed by autoradiography with reference to molecular weight standards.

Probe length can be affected by the concentration of DNA, primer and nucleotide, the size of the template DNA and also radiolysis of the labeled probe. The data in the figure was obtained using linearized plasmid DNA, 4.5 Kb in length under the standard labeling conditions.

It is recommended that not less than 10 pmol and not more than 125 pmol of any labeled dNTP is used in the reaction and combinations shown offer optimum balance of stability, specific activity and economy.

**Table 1**

Compounds	Specific activity (see note a)	Formulation (see note b)	Quantity of each dNTP required				Specific activity of probe (see notes c,d and e)	
	TBq/mmol	Ci/mmol	MBq	μCi	μl	pmol	dpm/ug	
[α- <sup>32</sup> P]dCTP	~111	~3000	1	1.85	50	5	17	1.9 × 10 <sup>9</sup>
	~222	~6000	1	7.4	200	20	32	5.3 × 10 <sup>9</sup>
[α- <sup>32</sup> P]dATP	~111	~3000	1	1.85	50	5	17	3.4 × 10 <sup>9</sup>
[α- <sup>32</sup> P]dCTP	~111	~3000	1	1.85	50	5	17	
[α- <sup>32</sup> P]dATP	~111	~3000	1	1.85	50	5	17	3.7 × 10 <sup>9</sup>
[α- <sup>32</sup> P]dCTP	~111	~3000	1	1.85	50	5	17	
[α- <sup>32</sup> P]dGTP	~111	~3000	1	1.85	50	5	17	

**a.** At the specific activity reference date of the labeled nucleotide.


**b.** Formulation code 1 = 370 MBq/ml, 10 mCi/ml in stabilized aqueous solution.

**c.** The probe specific activities were calculated using observed incorporation levels which are similar to those found in figure 3b.

**d.** It is important to note that the specific activity of probes made from different amounts of labeled dNTP cannot be calculated on a proportional basis, because net DNA synthesis occurs.

**e.** Brackets enclose nucleotides used in combination.


### a) Specific activity



- i) One labeled dNTP
- ii) Two labeled dNTP
- iii) Three labeled dNTP

**Figure 3.** The use of  $[\alpha\text{-}^{32}\text{P}]\text{dNTPs}$  in the Megaprime DNA Labeling System (see notes on page 26).

## b) Incorporation efficiency




- i) One labeled dNTP
- ii) Two labeled dNTP
- iii) Three labeled dNTP

**Figure 3.** The use of [ $\alpha$ - $^{32}$ P]dNTPs in the Megaprime DNA Labeling System (see notes on page 26).



### c) Probe length



- i) One labeled dNTP
- ii) Two labeled dNTP
- iii) Three labeled dNTP

**Figure 3.** The use of  $[\alpha\text{-}^{32}\text{P}]\text{dNTPs}$  in the Megaprime DNA Labeling System (see below).

#### Notes to figure 3

- a. The results shown are the means of a number of experiments in which different nucleotides and combinations of nucleotides were used. Observed results may deviate  $\pm 10\%$  from those shown.
- b. As the number of different labeled nucleotides is increased, at a given level of total input label, the net synthesis of DNA is reduced.  
Although the overall incorporation efficiency is reduced the labeled product is of a higher specific activity.

- c. The data was generated using the standard labeling protocols.  
If dNTPs <3000 Ci/mmol are to be used, then the desired probe specific activity must be multiplied by a conversion factor, before determining the amount of input label.

For a single labeled dNTP:-

$$\text{Total input label (pmols)} = \frac{3000 \text{ Ci/mmol}}{\text{specific activity of dNTP to be used}} \times \frac{\text{required probe specific activity}}{\text{specific activity}}$$

For more than one labeled dNTP the mean specific activity of the labeled dNTP to be used should be inserted in the above calculation.

Having determined the required number of pmols of input label with reference to figure 3a, the required volume of each labeled dNTP can be calculated. Note that the figures give the **total** amount of input label required. If more than one labeled dNTP is to be used, this figure should be divided by the number of labeled dNTPs to be used to give the required number of pmols of each labeled dNTP.

Volume of each labeled dNTP required in  $\mu\text{l}$  =

$$\frac{\text{pmol of dNTP required} \times \text{specific activity of dNTP (Ci/mmol)} \times 10^{-3}}{\text{radioactive concentration of dNTP (mCi/ml)}}$$

## 6. Appendices

### 6.1. Appendix I. Labeling of DNA fragments in low melting point agarose

The DNA samples produced by the following protocol have been found to be labeled to approximately the same extent as purified DNA. 15–20 minutes at 37°C is optimum for labeling. The standard labeling protocol may be found to be more appropriate for labeling DNA in agarose as the volume of DNA to be added using the new protocol is limited to 5 µl, requiring a relatively high initial DNA concentration.

Protocol	Notes
<b>1.</b> Fractionate restriction endonuclease digested DNA in a suitable low melting point agarose gel containing 0.5 µg/ml ethidium bromide. Estimate the DNA content of the band by reference to a set of standards of known concentration on another track. 250 ng should allow 25 ng to be used in the standard labeling protocol without further concentration	<b>1.</b> A low melting point agarose of high purity for example SepRate-LMP is recommended for maximum labeling efficiency.
<b>2.</b> Excise the desired band cleanly, with the minimum of excess agarose and transfer to a pre-weighed 1.5 ml microcentrifuge tube.	<b>2.</b> It is recommended that the exposure to UV light is minimized, as prolonged exposure can damage the DNA.

Protocol	Notes
3. Add water to a ratio of 3 ml per gram of gel and place in a boiling water bath for 5 minutes to melt the gel and denature the DNA.	3. If the DNA is not to be used immediately divide the boiled samples into suitably sized aliquots and store at -15°C to -30°C in a non frost-free freezer.
4. If the DNA is to be used immediately remove the appropriate volume containing 25 ng, add to the primers as indicated in the labeling protocol (page 11, step 3). The volume of DNA should not exceed 25 µl for the standard labeling protocol.	4. When using DNA which has been previously boiled and then stored at -15°C to -30°C, first place the tube in a boiling water bath for 30 seconds to melt the agarose, before removing the required volume containing 25 ng. Do not reboil DNA aliquots more than twice.
5. Incubate the labeling reaction for 15–20 minutes at 37°C.	

## 6.2. Appendix II. Monitoring the reaction and calculating the specific activity of the labeled DNA

### A. Adsorption to DE81 paper

Monitoring of the progress of the labeling reaction and measurement of probe specific activity can be achieved by determining the proportion of the radionucleotide incorporated during the Megaprime reaction.

1. Remove a 1 or 2  $\mu\text{l}$  aliquot of the reaction mixture to a clean microcentrifuge tube containing 20  $\mu\text{l}$  of water or 10 mM Tris/HCl pH.8.0. 1 mM EDTA buffer. Mix well by pipetting up and down.
2. Spot, in quadruplicate, 5  $\mu\text{l}$  aliquots of this dilution on to Whatman DE81 chromatography paper squares (minimum size  $1 \times 1 \text{ cm}$ ), placed on a non-absorbent backing. These squares may be marked with a pencil for identification if required.
3. Take two of the filters and dry under a heat lamp. 10–15 minutes should be adequate.
4. Wash the remaining two filters twice for 5 minutes each, at room temperature in excess  $2 \times \text{SSC}$  (30 mM  $\text{Na}_3\text{citrate}$ , 300 mM  $\text{NaCl}$  pH7.0) using gentle agitation. Rinse briefly in distilled water and then once with ethanol for 5 minutes. Then dry the filters under a heat lamp.
4. In aqueous solution DE81 paper becomes fragile and care should be taken when handling. In order to stabilize the paper the squares are rinsed in ethanol.

## Protocol

## Notes

5. Place the squares in separate vials with at least 5 ml of scintillation fluid and count.

6. Efficiency of counting will vary, but the percentage incorporation can be used to calculate probe specific activity. Unlike the nick translation labeling reaction, Megaprime labeling leads to net DNA synthesis, and so the total amount of DNA at the end of the reaction must be calculated.

Total amount of DNA (A) ng =

$$\frac{\text{Total number of } \mu\text{Ci added} \times 13.2 \times \% \text{ incorporation}}{\text{Number of radioactive dNTPs added} \times \text{average specific activity of dNTPs added}} + 25$$

Number of radioactive dNTPs added  $\times$  average specific activity of dNTPs added

This assumes a 25% content of any one dNTP in the newly synthesized DNA, and 25 ng of template DNA.

\*13.2 equals four times the average molecular weight of the four dNTPs divided by 100.

5. Determination of the proportion of the  $^{32}\text{P}$  labeled nucleotide incorporated may be achieved using Cerenkov counting if desired in this case drying the filter is not necessary.

6. The mean value of the counts on the washed filter represents the proportion of the radionucleotide incorporated into the DNA probe, while the mean of the unwashed filters represents the total amount of radioactivity in the reaction mix, such that;

$$\% \text{ incorporation} = \frac{\text{mean counts on washed filters}}{\text{mean counts on unwashed filters}} \times 100$$

**6. Continued.**

The amount of radioactivity incorporated during the reaction (B) in dpm.

$B = \text{total number of } \mu\text{Ci added} \times 2.2 \times 10^4 \times \% \text{ incorporation}$

Thus the specific activity of the labeled DNA is

$\text{specific activity} = \frac{B}{A} \times 10^3 \text{ dpm per } \mu\text{g}$

**B. Precipitation with trichloroacetic acid**

Plastic or siliconized glass tubes must be used to avoid adsorption of DNA.

1. Dilute an appropriate aliquot of the reaction mixture as described in section A1.
2. Transfer 1–10  $\mu\text{l}$  of diluted reaction mixture to two duplicate tubes containing 200  $\mu\text{l}$  water or 0.2M EDTA and 50  $\mu\text{l}$  carrier DNA solution. Mix well. Use this mixture (less any set aside in step 3) for the TCA precipitation described in step 4 below.
3. Set aside an appropriate aliquot from each tube in step 2 for the determination of total input radioactivity.
4. To the diluted samples from step 2, add 2 ml ice-cold 10% trichloroacetic acid (TCA) solution, vortex, and allow to stand in an ice-bath for 10–15 minutes. The labeled and carrier DNA will co-precipitate. **Note that TCA is corrosive, and care should be taken in its handling.**
5. Collect the precipitated DNA by vacuum filtration on a glass fibre or nitrocellulose filter disc.

6. Wash the filter discs six times with 2 ml 10% TCA solution and dry the filter discs thoroughly, for example using an infra-red lamp. Avoid overheating and possible charring of the discs.
7. Count the dried filter discs by liquid scintillation or Cerenkov ( $^{32}\text{P}$ ) and count with the samples set aside in step 3.
8. Determine % incorporation and probe specific activity as in section A6.

### 6.3. Appendix III. Removal of unincorporated nucleotides

Removal of unincorporated nucleotides is sometimes desirable to reduce background produced by the probe during hybridization. It is considered important to remove these free nucleotides particularly if the radioactive probe is to be kept for several days before use or the incorporation is less than 50%. If  $^{32}\text{P}$  or  $^{33}\text{P}$ -labeled probes are to be used in combination with GE Healthcare's new Rapid-hyb buffer (RPN1635/6), purification is not required unless the probe is to be used more than 24 hours after preparation. Probes can be purified by Sephadex chromatography or selective precipitation (8,9).

#### A. Sephadex™G-50 spin columns

Probe reaction are passed through columns packed with Sephadex G-50, which retains the free nucleotides within the column matrix. A number of pre-packed columns are commercially available. However columns may also be prepared as indicated below:

1. Equilibrate Sephadex G-50 in TE buffer either overnight or at 65°C for 1–2 hours.
2. Plug a 1.0 ml syringe with a piece of siliconized glass wool.
3. Fill the syringe with the equilibrated Sephadex. Place in a 15 ml conical tube, in which a decapped 1.5 ml microcentrifuge tube has been inserted. Centrifuge at 1600 g for 5 minutes. Remove



any liquid from the microcentrifuge tube. Refill with Sephadex and centrifuge as before. Continue until the column is packed to a volume of 1 ml.

4. Add a volume of TE buffer equal to the reaction volume, to the top of the column and centrifuge, as in step 3. A minimum of 50  $\mu$ l should be applied to the column.
5. Repeat once more to ensure fractions of the correct size are collected from the column.
6. Place the column in a clean 15 ml conical tube containing a decapped 1.5 ml microcentrifuge tube.
7. Apply the DNA sample to the column. Centrifuge as before. The purified probe is collected in the microcentrifuge tube.

#### **B. Selective precipitation of labeled DNA**

The following protocol leads to precipitation of DNA greater than about 20 nucleotides in length with unincorporated nucleotides remaining in solution. Recovery of the labeled DNA by this method varies according to the DNA concentration and size, and may be as low as 50%.

1. Add one volume of 4 M ammonium acetate, pH 4.5 to the nick translation reaction, and mix gently by pipetting up and down.
2. Add four volumes of ethanol, mix by inversion. Chill the mixture for 15 minutes in a dry-ice ethanol bath or place at  $-70^{\circ}\text{C}$  for at least 30 minutes.
3. Thaw the mixture if necessary by placing at  $37^{\circ}\text{C}$  for 2 minutes.
4. Spin in a microcentrifuge for 15 minutes. Carefully aspirate and dispose of supernatant in a suitable manner.
5. Wash the pellet once in 0.5 ml of 0.67 M ammonium acetate, pH 4.5, 67% ethanol at room temperature by gentle inversion, centrifugation and aspiration.

6. Wash the pellet once in 90% ethanol, in the same manner. Dry the pellet.
7. Finally redissolve the DNA pellet in TE buffer for use as a probe and for storage.

## 6.4. Appendix IV. Additional equipment and reagents

TE buffer (10 mM Tris/HCl, pH 8.0, 1 mM EDTA)

0.2 M EDTA solution

Adjustable pipettes for example Pipetman™

Sterile pipette tips

Waterbaths at 37°C and 100°C

Polypropylene microcentrifuge tubes

Microcentrifuge

Gloves

Radiation safety equipment

DE81 ion-exchange chromatography paper (Whatman)

Trichloroacetic acid (TCA) solution: 10% (w/v) TCA in water

Filter discs; glass fibre or nitrocellulose

Plastic or siliconized glass tubes, capacity ~5 ml

Filtration apparatus

2× SSC (30 mM Na<sub>3</sub> citrate, 300 mM NaCl, pH 7.0)

## 7. Troubleshooting guide

If poor results are obtained, the following guide may help to determine the cause of the problem.

Problem	Possible cause	Remedy
<b>1. Low signal</b>	<b>1.</b> Incomplete denaturation of template DNA	<b>1.</b> Ensure denaturation protocol is followed.
	<b>2.</b> Low probe concentration	<b>2.</b> Accurately measure the concentration of template DNA used in the labeling reactions. Check recovery of probe if purification is performed to remove unincorporated nucleotide.
	<b>3.</b> Low probe specific activity	<b>3.</b> If the specific activity of the labeled DNA is lower than expected, a labeling reaction should be carried out using a sample of the control DNA supplied with the system. If this proceeds satisfactorily, check the concentration and purity of your DNA.

<b>Problem</b>	<b>Possible cause</b>	<b>Remedy</b>
	4. Loss of dNTP during evaporation	4. If the dNTP solution has been evaporated to dryness prior to use, handling losses may have occurred. Check this loss has not occurred during lyophilization of the solvent, during transfer of the reconstituted dNTP solution or by adsorption of the dNTP onto the walls of the tube. If necessary the reconstituted dNTP solution may be counted and an adjustment made before setting up the labeling reaction.
2. Non-specific background over whole of filter	1. Presence of unincorporated label	1. Unincorporated nucleotides can give high backgrounds. Remove by Sephadex G-50 spin columns or ethanol precipitation (see page 32 for protocol)

<b>Problem</b>	<b>Possible cause</b>	<b>Remedy</b>
	2. Concentrated probe has contacted membrane directly during probe addition	2. It is suggested that up to 1.0 ml of the buffer used for prehybridization is withdrawn for mixing with the probe. The mixture should then be added back to the hybridization container in an area away from the filter.
	3. Probe concentration is too high	3. Ensure measurement of template DNA concentration is accurate
	4. Probe not denatured	4. Non-denatured double-stranded probes often give high backgrounds.

## 8. References

1. FEINBERG, A.P. and VOGELSTEIN, B., *Anal. Biochem.*, **132**, pp.6-13, 1983.
2. FEINBERG, A.P. and VOGELSTEIN, B., *Addendum Anal. Biochem.*, **137**, pp.266-267, 1984.
3. SOUTHERN, E.M., *J.Mol.Biol.*, **98**, pp.503-517, 1975.
4. THOMAS, P.S., *Proc. Natl. Acad. Sci. USA.*, **77**, pp.5201-5205, 1980.
5. MEINKOTH, J. and WAHL, G., *Anal. Biochem.*, **138**, pp. 267-284, 1984.
6. GRUNSTEIN M. and HOGNESS, D.S., *Proc. Natl. Acad. Sci. USA.*, **72**, pp. 3961-3965, 1975.
7. HODGSON, C.P., FISK, R.Z. and WILLET, L.B., *Biotechniques*, **6**, pp.208-211.
8. SAMBROOK, J. FRITSCH, E.F. and MANIATIS, T., *Molecular Cloning*, a laboratory manual (second edition), Cold Spring Harbour Laboratory, 1989.
9. MUNDY, C.R., CUNNINGHAM, M.W. and READ, C.A., *Essential Molecular Biology; A Practical Approach Vol 2* (T.A. Brown, ed) Oxford University Press, Oxford, 1991. pp.57-109.

## 9. Related Products

### Labeling systems

Nick Translation Kits	N5000/5500
Rediprime II DNA Labeling System	RPN1633/4
Ready-To-Go DNA Labeling Beads (-dCTP)	27-9240-01

### Hybridization buffers

Rapid-hyb buffer	RPN1635/6
------------------	-----------

### Hybridization products

Hybond™ - Range of nylon and nitrocellulose blotting membranes	
Hybridisation oven/shaker (220/240V 50Hz)	RPN2510
Hybridisation oven/shaker (110/115V 60Hz)	RPN2511E

### Autoradiography products

Hyperfilm™ - high performance autoradiography films	
Hypercassettes™ and Hyperscreens™ - available from stock	
Sensitize™ preflash gun.	RPN2051
TrackerTape™	RPN2050

An adhesive waterproof tape that phosphoresces to give a permanent written image on autoradiography film

### Safety Products

Radiation safety products for safe handling and storage of  $^{32}\text{P}$ / $^{33}\text{P}$ / $^{35}\text{S}$  and  $^{125}\text{I}$ , liquid scintillation products

### Labeled dNTPs

See Table 2

**Table 2.** Labeled dNTPs and analogues available from GE Healthcare

Compound	Specific Activity		Product code
	TBq/mmol	Ci/mmol	
[8- <sup>3</sup> H]dATP	0.37–1.1	10–30	TRK 347
[1',2',5- <sup>3</sup> H]dCTP	1.85–3.14	50–85	TRK 625
[5- <sup>3</sup> H]dCTP	0.55–1.1	15–30	TRT 352
[8- <sup>3</sup> H]dGTP	0.185–0.740	5–20	TRK 350
[1'2,(- <sup>3</sup> H]dGTP	0.9–1.85	25–50	TRT 627
[methyl, 1',2'- <sup>3</sup> H]TTP	3.3–4.8	90–130	TRK 576
[methyl- <sup>3</sup> H]TTP		40–60	TRK 424
		30	TRK 354

See GE Healthcare Products catalogue for further details.









**GE Healthcare offices:**

GE Healthcare Bio-Sciences AB  
Björkgatan 30, 751 84 Uppsala,  
Sweden

GE Healthcare Europe GmbH  
Munzinger Strasse 5, D-79111  
Freiburg,  
Germany

GE Healthcare Bio-Sciences Corp.  
800 Centennial Avenue, P.O. Box 1327,  
Piscataway, NJ 08855-1327,  
USA

GE Healthcare Japan Corporation  
Sanken Bldg. 3-25-1, Hyakunincho,  
Shinjuku-ku, Tokyo 169-0073,  
Japan

For contact information for your local office,  
please visit: [www.gelifesciences.com/contact](http://www.gelifesciences.com/contact)

GE Healthcare UK Limited  
Amersham Place  
Little Chalfont, Buckinghamshire,  
HP7 9NA, UK

<http://www.gelifesciences.com/sampleprep>



imagination at work